

# The Guide to Nucleic Acid Labelling and Detection



# Contents

Introduction	page 2
Choosing a non-radioactive labelling system	page 3
AlkPhos Direct Labelling and Detection Systems	page 6
Gene Images Random-Prime Labelling and Detection Systems	page 10
Gene Images 3'-Oligolabelling and Detection System	page 12
ECL Direct Labelling and Detection System	page 14
Choosing a radioactive labelling system	page 16
Rediprime II DNA Labelling System	page 18
5'-End Labelling Kit	page 20
3'-End Labelling Kit	page 21
Ready-To-Go PCR Beads	page 22
RNA Labelling Kit	page 24
Rapid-hyb Buffer	page 25
Redivue Nucleotides	page 26
Hybond Membranes for Nucleic Acid Blotting	page 28
Autoradiography Films	page 30
<b>Amersham Biosciences Scanner</b>	<b>page 32</b>
Accessory products	page 34

# Introduction

For many years the use of radiolabels was the method of choice for nucleic acid labelling and detection. More recently, however, considerable research has lead to the development of a comprehensive range of quality systems using both radioactive and non-radioactive technologies. This extensive range now provides you with an excellent choice of systems for use with DNA, RNA and oligonucleotides in a variety of applications.

Significant developments in the area of non-radioactive labelling and detection include the use of light-based detection systems and direct labelling. The introduction of these techniques means that non-radioactive methods now present a viable alternative to the use of traditional radioactive techniques for many applications.

This breadth of choice can in itself cause problems. Choosing the most appropriate system for a particular application can be a challenge for even the most experienced researchers.

The choice is influenced by many factors, not all of which are directly related to the application. For example radioactive labelling has the advantage of high sensitivity and extremely robust protocols, however safety and long-term probe storage are serious considerations.

The aim of this guide is to give an overall introduction to the different radioactive and non-radioactive methodologies available, highlighting their advantages and drawbacks. Individual products are then presented in more detail allowing you to choose the most appropriate system for your application.

Factor	Considerations
Blotting system	Southern, Northern, dot/slot, library screen. Amount of target nucleic acid present
Label	Radioactive or non-radioactive Direct or indirect
Probe	DNA, RNA, oligonucleotide Size and purity of probe
Speed and convenience	Multiple exposures Time constraints
Probe storage	Stability
Results	Autoradiography film Image analysis with scanners
Membrane	Charged or uncharged nylon Nitrocellulose

Table 1 Factors affecting the choice of labelling and detection system

# Choosing a non-radioactive labelling system

**In the past the use of non-radioactive labelling and detection technologies has been limited by the relatively low levels of detection sensitivity as compared with radioactive detection. However, the further development of non-radioactive labelling methods, coupled with light-based detection has resulted in systems that now offer a real alternative to the use of radioisotopes in many applications where target levels are 50 fg and above.**

**Non-radioactive systems offer the following benefits:**

<b>Speed</b>	Exposure times with non-radioactive systems range from minutes to hours.
<b>Sensitivity</b>	Applications where target levels are 50 fg and above, e.g. single copy Southern and Northern blotting, colony/plaque blots.
<b>Stable probes</b>	Non-radioactively labelled probes are stable for at least 6 months, removing the need to label probes immediately before use.
<b>Ease of handling</b>	Eliminates handling and waste regulatory issues associated with the use of radioactivity.

Making the correct choice of non-radioactive system requires consideration of both the labelling and detection procedures.

## Labelling

Non-radioactive labelling methods are separated into two different approaches. Early non-radioactive systems utilized indirect labelling techniques, but these have now been largely superseded by the introduction of superior direct labelling methods, an example of which is AlkPhos Direct™.

**Direct labelling** Direct labelling\* using alkaline phosphatase is the most recent advance in the non-radioactive labelling of nucleic acids. It offers significant improvements in speed and convenience over indirect systems without compromising sensitivity.

With direct labelling systems, the enzyme molecule (alkaline phosphatase or horseradish peroxidase) is directly crosslinked to the nucleic acid probe in a simple 30-minute reaction.

Detection is rapid because direct labelling methods eliminate the need for an antibody conjugate incubation and the associated blocking and multiple washing steps.




Figure 1 Direct labelling

\*For licence information see page 37

# Choosing a non-radioactive labelling system

**Indirect labelling** Indirect labelling involves the introduction into the probe of nucleotides tagged with a hapten or 'reporter' molecule. Amersham Biosciences kits use fluorescein as the hapten molecule. The nucleotide analogues are readily incorporated into the probe by standard labelling methods and hybridization and stringency washing are carried out under standard conditions. Labelled probes are then detected with a highly specific anti-fluorescein antibody conjugated to either alkaline phosphatase (AP) or horseradish peroxidase (HRP) enzyme. This is followed by enzyme catalysed detection, using the appropriate substrate.




Figure 2 Indirect labelling

## Detection

Choice of detection method will depend on the enzyme used in the labelling stage as well as the detection medium required. Horseradish peroxidase can only be detected using a chemiluminescent substrate, but alkaline phosphatase can be detected using chemiluminescent and chemifluorescent substrates. Chemiluminescent detection systems are designed for use with autoradiography film whereas fluorescence based detection systems require the use of suitable scanning equipment.

## ECL

Based on the enhanced chemiluminescent reaction of luminol with horseradish peroxidase, ECL™ substrate can be used to detect probes which have been labelled either directly or indirectly with horseradish peroxidase. Rapid light output enables results to be achieved in 10 to 15 minutes. It is the substrate of choice for target amounts above 500 fg.

## CDP-Star

Based on the chemiluminescent breakdown of the dioxetane CDP-Star™ by alkaline phosphatase\*, this substrate can be used to detect probes which have been labelled either directly or indirectly with alkaline phosphatase. It is ideal for applications requiring high sensitivity detection such as genomic Southern and Northern blots.

Light output is rapid and continues for up to five days allowing exposure optimization and multiple exposures to be taken.

## ECF

Based on the breakdown of the ECF™ substrate by alkaline phosphatase\*, resulting in a highly fluorescent product which is localized at the site of hybridization. This chemifluorescent signal is detected using a suitable fluorescence scanner such as the Amersham Biosciences™ FluorImager™ or Storm™. ECF detection can be used to detect probes which have been either directly or indirectly labelled with alkaline phosphatase and is particularly suitable for applications where quantification is important.

Labelling and detection system	Sensitivity	Time from hybridization to detection	Duration of light output	Strip before re-probing	Quantification	Recommended application
AlkPhos Direct	0.06 pg	1 hour	5 days	yes	no	Single copy Southern and Northern
ECL Direct	0.5 pg	1 hour	1-2 hours	no	no	High target applications e.g. colony/plaques
Gene Images Random-Prime	0.1 pg	3 hours	5 days	yes	no	High sensitivity Northern
Gene Images 3'- end labelling with CDP-Star	0.1 pg	3 hours	5 days	yes	no	Oligo screening with stringency control
Gene Images 3'- end labelling with ECF	120 pg	3 hours	1-2 days	yes	yes	Quantification
ECL Random-Prime	0.5 pg	3 hours	1-2 hours	yes	no	Medium target Southern with DNA probes
ECL 3'-end labelling	0.2 pg	3 hours	1-2 hours	yes	no	Medium to high target Southern with oligo probes
ECF Random-Prime	0.25 pg	3 hours	1-2 days	yes	yes	Quantification

Table 2 Guide to the properties of non-radioactive labelling and detection systems

\*For licence information see page 37.

# AlkPhos Direct Labelling and Detection

**AlkPhos Direct™** combines the convenience of direct enzyme labelling (no blocking or antibody stages) with those of alkaline phosphatase detection (long light output and high sensitivity).

Labelling is complete within 30 minutes in a single tube protocol and the resulting probe is ready for use in hybridizations without further purification.

Due to the thermostable nature of the enzyme, hybridization stringency can be controlled by adjusting temperature as well as salt concentration.

The system can be used with either chemiluminescent or fluorescent detection reagents. The most sensitive result is achieved with **CDP-Star**, a chemiluminescent substrate with a fast light output. Alternatively, the use of **ECF** substrate generates a fluorescent signal suitable for use with fluorescence scanning devices such as Amersham Biosciences FluorImager or Storm instruments.

## Protocol Summary




Figure 3 Probe labelling

### Probe labelling

Denatured or single-stranded DNA or RNA probe is mixed with the labelling buffer and alkaline phosphatase. Formaldehyde is then added to covalently crosslink the enzyme to the probe. There is no need to purify the probe before hybridization.




Figure 4 Hybridization

### Hybridization and stringency

The probe is hybridized to the blot using the specially formulated AlkPhos Direct hybridization buffer, which stabilizes the activity of the enzyme. After a 15 minute pre-hybridization step, hybridization with probe is performed overnight at the required temperature. For higher target amounts, a 2 to 4 hour hybridization may be sufficient.

# Systems




Figure 5 Detection

Figure 5 Detection

## Detection

Chemiluminescence with CDP-Star:

Hybridized blots are detected with CDP-Star reagent. Following a two-minute incubation, blots are exposed to Hyperfilm™ MP for 1 to 2 hours or a light capture scanning device. Sensitivities as low as 60fg can be achieved using CDP-Star.

or

Chemifluorescence with ECF:

Hybridized blots are detected with ECF substrate. Non-fluorescent substrate is catalysed by alkaline phosphatase to produce a fluorescent signal which accumulates over time; low sensitivity applications yield results after 1 hour, high sensitivity applications usually require overnight incubation.

## Recommended Applications

AlkPhos Direct is a fast and easy-to-use system suitable for the majority of routine blotting applications using either DNA, RNA or oligonucleotide probes.

The versatility of the system, combined with the fact that both DNA and RNA probes are labelled equally efficiently, makes it the ideal choice for a busy laboratory carrying out a variety of different blotting and hybridization techniques.

For Southern and Northern blotting applications AlkPhos Direct, combined with CDP-Star detection, offers a highly sensitive system capable of detecting down to 60 fg of target DNA. Thus making it the ideal choice for the non-radioactive detection of single-copy genomic Southern blots.

Re-probing can be difficult with non-radioactive systems due to the need to remove both probe and antibody layers in an indirect system. The simplicity of the AlkPhos Direct system means that there are less components to be removed during the stripping procedure. Probe removal is therefore more effective with less damage to the membrane.

For labelling oligonucleotides, an optimized protocol has been developed which is capable of labelling probes down to 17 basepairs in length. Request TechTip number 174 for additional information.

Quantification of dot and slot blots is possible by using ECF detection in conjunction with AlkPhos Direct. Analysing results on a fluorescence scanner such as **Amersham Biosciences** Storm or FluorImager allows accurate quantification and analysis of results.



# AlkPhos Direct Labelling and Detection

**Speed** Time savings of 3 to 4 hours over conventional indirect methods are achievable as a result of the rapid and simple labelling reaction, and the absence of antibody incubations.

**Consistency** Quality control on each batch ensures consistent performance and provides confidence in probe labelling and concentration.

**Stringency control** Due to the thermostable nature of the enzyme, hybridization stringency can be controlled by elevating temperatures up to 75°C as well as by decreasing salt concentrations.

**Versatility** Labelling of DNA, RNA and oligonucleotide probes with detection by chemiluminescence or chemifluorescence offers a high degree of versatility.

**Stability** Probes are stable for up to six months if stored at -20°C in 50% glycerol.

**Efficiency** Nucleic acid probes are labelled efficiently and uniformly with enzyme.

**Accuracy** The crosslinking labelling reaction does not result in any net synthesis of probe and the amount of probe present before and after labelling does not change. Therefore probe concentrations during hybridization can be determined more accurately.




Figure 7  
Northern blot (Hybond-N+) of human skeletal muscle (total loadings of 0.5, 0.25, 0.125, 0.063 µg). Hybridized with GAPDH probe.  
(A) First detection  
(B) CDP-Star re-applied after stripping  
(C) CDP-Star re-applied after re-probing




Figure 6  
Southern blot  
Cosmid DNA digested with Not I and EcoR I, probed with a 1.1kb probe labelled with AlkPhos Direct and a competitor's hapten-based system  
(A) AlkPhos Direct  
(B) Competitor's hapten-based system  
Courtesy of Janet Bartels, Yale University

# Systems




Figure 8

Dot blots loaded with 5, 2, 1, 0.5, 0.1 and 0.01 ng of a 1.5kb N-ras insert on Hybond-N+ nylon membrane. The 7th dot was a control of 5ng carrier DNA. Results were obtained with CDP-Star detection reagent.




Figure 9

Control of stringency by wash temperature. Human genomic Southern blots (2, 1 and 0.5 µg loadings) hybridized with BCL2 at 5ng/ml. Washed at:

(A) 50°C (B) 55°C (C) 60°C  
0.2% blocking reagent included in the primary wash buffer.

## Ordering Information

AlkPhos Direct Labelling and Detection System with CDP-Star	25 reactions	RPN 3690
AlkPhos Direct Labelling and Detection System with CDP-Star	50 reactions	RPN 3691
AlkPhos Direct Labelling and Detection System with ECF	25 reactions	RPN 3692
AlkPhos Direct Labelling Module	25 reactions	RPN 3680
AlkPhos Direct Labelling Module	50 reactions	RPN 3681
AlkPhos Direct Hybridization Buffer	For 5000 cm <sup>2</sup> membrane	RPN 3688
CDP-Star Detection Reagent	For 2500cm <sup>2</sup> membrane	RPN 3682
CDP-Star Detection Reagent	For 5000cm <sup>2</sup> membrane	RPN 3683
ECF Substrate	For 2500cm <sup>2</sup> membrane	RPN 3685

## Related products

ECL Direct Nucleic Acid Labelling and Detection System	RPN3000
Hybond-N+ Membranes	see page 28
Hyperfilm MP	see page 30

# Gene Images Random-Prime Labelling and Detection

**Gene Images™ Random-Prime is designed to generate fluorescein labelled DNA probes using a random prime labelling reaction. The reaction requires only 50 ng of template DNA as the strand displacement properties of the enzyme, coupled with non-limiting concentrations of nucleotide, result in net synthesis of probe. The use of fluorescein allows the efficiency of the reaction to be checked before hybridization by a fluorescence assay.**

## Protocol Summary




Figure 10 Probe labelling

### Probe labelling

Fluorescein dUTP is incorporated into DNA in a random prime labelling reaction using Klenow DNA polymerase. Nonamer primers are used to ensure efficient hybridization to probe. Probes can be used without purification.

### Optional rapid labelling assay

The efficiency of probe labelling can be checked by spotting a small aliquot of labelled probe onto Hybond™-N+ membrane, washing to remove unincorporated fluorescein and visualizing on a UV transilluminator or fluorescent scanner.

### Hybridization

For optimum performance use the recommended hybridization buffer.

Sensitive detection is achieved by the use of a monoclonal antibody which is highly specific for alkaline phosphatase. 50 fg of target can be routinely detected in a human genomic Southern blot.

## Detection

Chemiluminescence with CDP-Star:

Hybridized blots are detected with CDP-Star reagent. Following a 2 minute incubation, blots are exposed to Hyperfilm MP for 1 to 2 hours. Sensitivities as low as 50 fg can be detected.

or

Chemifluorescence with ECF:

Hybridized blots are detected with ECF substrate. Non-fluorescent substrate is catalysed by alkaline phosphatase to produce a fluorescent signal which accumulates over time; low sensitivity applications yield results after 1 hour, high sensitivity applications usually require overnight incubation.




Figure 11 Indirect labelling

## Recommended Applications

Probes labelled with the Gene Images Random-Prime Labelling and Detection System can be used for most membrane applications especially in situations where there is insufficient probe for labelling by direct methods.


Gene Images Random-Prime Labelling and Detection Systems may be used with either the detection reagent CDP-Star for high sensitivity applications, or with ECF substrate for quantification on fluorescence scanning instrumentation.

Gene Images Random-Prime is especially recommended for high sensitivity Northern applications.


# Systems

**Confidence** The use of fluorescein-dUTP allows a simple and rapid (30 minute) probe labelling assay to be performed by monitoring the fluorescence of labelled probes on a UV transilluminator. This enables the efficiency of the labelling reaction to be checked before proceeding with the hybridization.

**Stability** Probes are stable for at least 9 months under recommended conditions minimizing the requirement for continual labelling of frequently used probes.



**Figure 12**  
5, 2, 1 and 0.5  $\mu$ g of human liver mRNA were blotted onto Hybond-N+, fixed using the Amersham crosslinker and hybridized with a fluorescein-labelled DNA probe, specific for p53, at 65°C. CDP-Star detection, 1 hour exposure on Hyperfilm-MP.



**Figure 13**  
Rapid labelling assay  
(a) Serial dilutions of fluorescein nucleotide 1/25, 1/50, 1/100, 1/250, 1/500, 1/1000. Membrane not washed  
(b) Aliquots of labelling reactions:  
(i) labelling reaction inhibited with EDTA at time 0  
(ii) labelled lambda DNA digested with Hind III  
(iii) labelled N-ras clone  
Membrane washed in 2xSSC at 60°C for 15 minutes. Visualized on a UV transilluminator

## Ordering Information

Gene Images Random-Prime Labelling and Detection System	30 reactions	RPN 3500
Gene Images Random-Prime Labelling and Detection System	60 reactions	RPN 3501
Gene Images Random-Prime Labelling Module	30 reactions	RPN 3540
Gene Images Random-Prime Labelling Module	60 reactions	RPN 3541
CDP-Star Detection Module	For 2500cm <sup>2</sup> membrane	RPN 3510
CDP-Star Detection Module	For 5000cm <sup>2</sup> membrane	RPN 3511
ECF Detection Module	For 2500cm <sup>2</sup> membrane	RPN 5750

## Related products

AlkPhos Direct Labelling and Detection System	RPN3690
Hybond-N+ membranes	see page 28

# Gene Images 3'-Oligolabelling and Detection

**Gene Images™ 3'-Oligolabelling system is based on the labelling of standard unmodified oligonucleotide sequences with fluorescein-dUTP in a reaction catalysed by terminal deoxynucleotidyl transferase. The optimized protocol results in a tail of 6-8 nucleotides, which provides good sensitivity with low background, and does not interfere with the specificity of hybridization. The use of fluorescein as the label allows the efficiency of the labelling reaction to be checked using a fluorescence assay before hybridization.**

## Protocol Summary




Figure 14 Probe labelling

### Probe labelling

Standard oligonucleotides are tailed with fluorescein-dUTP using the enzyme terminal deoxynucleotidyl transferase. There is no need to purify the probe before hybridization.

### Optional rapid labelling assay

Labelling efficiency can be checked by spotting a small aliquot of labelled probe onto Hybond-N+ membrane, washing to remove unincorporated fluorescein and visualizing on a transilluminator or fluorescence scanner.

## Hybridization

For optimum performance use the recommended hybridization buffer.

## Detection

Chemiluminescence with CDP-Star:

Hybridized blots are detected with CDP-Star reagent. Following a two minute incubation, blots are exposed to autoradiography film for 1 to 2 hours or detected using a light collection device.

or

Chemifluorescence with ECF:

Hybridized blots are detected with ECF substrate. Non-fluorescent substrate is catalysed by alkaline phosphatase to produce a fluorescent signal which accumulates over time; low sensitivity applications yield results after 1 hour; high sensitivity applications usually require overnight incubation.




Figure 15 Gene Images 3'-Oligolabelling with CDP-Star detection

# System

## Recommended Applications

The 3'-Oligolabelling system can be used to label probes for a wide variety of applications. The indirect labelling method and high sensitivity makes it especially suited to oligo screening applications requiring a high degree of stringency control.

**Ease of use** Simple protocols, and ready-to-use reagents ensures first-time success.

**Confidence** The use of fluorescein-dUTP allows a simple and rapid (30 minute) probe labelling assay to be performed by monitoring the fluorescence of labelled probes on a UV transilluminator. This enables the efficiency of the labelling reaction to be checked before proceeding with the hybridization. Use with any oligonucleotide, special synthesis is not required.

**Stability** Labelled probes are stable for at least 9 months, under recommended conditions allowing batches of probe to be prepared and stored until needed.



Figure 16

Dot blots of M13 mp8 (500 - 0.625 attomoles).  
3'-end labelled M13 17mer oligonucleotide probe.  
5 ng/ml, 2 hour hybridization and detected using  
CDP-Star. 1 hour exposure to Hyperfilm MP.

## Ordering Information

Gene Images 3'-Oligolabelling and CDP-Star Detection System	To label 1000 pmol oligo	RPN5776
Gene Images 3'-Oligolabelling Module	To label 1000 pmol oligo	RPN5770
CDP-Star Detection Module	For 2500 cm <sup>2</sup> membrane	RPN3510
CDP-Star Detection Module	For 5000 cm <sup>2</sup> membrane	RPN3511
ECF Detection Module	For 2500 cm <sup>2</sup> membrane	RPN5750

## Related products

AlkPhos Direct Labelling and Detection System with CDP-Star	RPN3690
---	---------

# ECL Direct Labelling and Detection System

The ECL Direct™ nucleic acid labelling and detection system is based on the direct labelling of DNA or RNA probes with the enzyme horseradish peroxidase (HRP) in a simple chemical reaction. Labelling takes only 10 minutes in a single tube reaction, and the resulting probe can be used directly in hybridization experiments without further purification. Detection of hybridized probe is achieved by generation of light via HRP catalysed oxidation of luminol with the use of an enhancer maximizing light output.




Figure 18

## Hybridization

The labelled probe is hybridized to the blot using the specially formulated ECL Gold Buffer, which stabilizes the activity of the enzyme. Following a 15 minute pre-hybridization, hybridization with probe is performed overnight. For colony or plaque screens a 2 to 4 hour hybridization is sufficient.

## Protocol Summary




Figure 17 Probe labelling

## Probe labelling

The labelling solution is added to the unlabelled, denatured single-stranded DNA or RNA probe. The positively charged HRP complex is bound electrostatically to the negatively charged DNA. Glutaraldehyde is added to covalently cross-link the HRP to the probe. There is no need to purify the probe before hybridization. Labelled probes can be stored in 50% glycerol for several months at -15°C to -30°C.




Figure 19 Detection

## Detection

The hybridized blot is soaked in the detection reagent. Where the probe is bound there is peroxidase-catalysed oxidation of luminol and subsequent enhanced chemiluminescence. The resulting light signal is detected on autoradiography film.

## Recommended Applications

ECL Direct™ is a fast and easy system for high to medium target amount applications such as colony/plaque screens, dot blots, genome mapping, and PCR® product analysis. It is particularly useful for applications which require a blot to be re-probed several times and screening applications requiring rapid throughput and results.




Figure 20  
Random YAC DNA  
digested with PvuII  
and hybridized with  
AluI human DNA  
repetitive probe  
(exposure time of  
30 minutes), using  
ECL Direct labelling.




Figure 21  
Comparison of 1st and 6th re-probing  
of Southern blotted human genomic  
DNA bound to Hybond-N+. Minimal  
reduction in signal intensity is observed.  
A) 6th re-probing  
B) 1st probing

**Ease of use** Simple protocol coupled with high initial signal output makes it the most rapid non-radioactive system available.

**Speed** Less than 4 hours from probe labelling to detection in high target amount applications with no lengthy antibody steps.

**Ease of re-probing** HRP is inactivated by the chemiluminescent reaction therefore no need to strip blots before reprobing. This saves time and avoids possible membrane damage.

**Consistency** Consistent results combining strong signals with very low backgrounds on nylon membranes.

**Convenience** Convenient, ready-to-use ECL Direct Gold hybridization buffer formulated with an exclusive rate enhancer\*.

**Efficiency** DNA or RNA probes are labelled equally efficiently.

**Stability** Probes can be stored for at least 3 months in 50% glycerol, avoiding the need for frequent labelling.

## Ordering Information

ECL Direct Labelling and Detection System	To label 5 µg nucleic acid	RPN 3000
ECL Direct Labelling and Detection System	To label 10 µg nucleic acid	RPN 3001
ECL Direct Detection Reagents	For 2000cm <sup>2</sup> membrane	RPN 3004
	For 4000cm <sup>2</sup> membrane	RPN 2105
ECL Direct Gold Buffer	For 4000cm <sup>2</sup> membrane	RPN 3006
ECL Direct Labelling System	To label 5 ng nucleic acid	RPN 3005

## Related products

AlkPhos Direct Labelling and Detection System	RPN3690
---	---------



# Choosing a radioactive labelling system

**Labelling with  $^{32}\text{P}$ -labelled nucleotides remains the most sensitive and robust technology available. Radiolabelled probes offer the following benefits:**

**Sensitivity**  $^{32}\text{P}$ -labelled probes offer the highest level of sensitivity, ideal for detection where there is a very low abundance of target.

**Prolonged emission** Radioactively labelled blots continue to emit signal for several weeks enabling multiple exposures and long exposure times to be carried out as required.

**Reliability** Radioactive labelling systems are very forgiving. Precise optimization of hybridization conditions is less critical, and hence less time is used determining ideal hybridization and stringency conditions.

**Flexibility** Incorporation of radiolabelled nucleotides into a probe can be controlled by the labelling conditions. Therefore probes can be produced with differing specific activities and incorporation levels can be measured accurately, allowing exact probe concentrations to be determined.

## Choice of radiolabel

The selection of radiolabel for a particular application depends on the level of sensitivity and resolution required. For most membrane hybridization techniques sensitivity is most important, so a high-emission energy label will provide sensitive results with shorter exposure times but resolution will be compromised.

In applications where resolution is most important, such as microsatellite analysis and sequencing, a lower energy-emitting isotope will give improved resolution, but longer exposure times will be required.

The three most widely used isotopes in the life science field are:

**$^{32}\text{P}$**  For high-sensitivity membrane applications  $^{32}\text{P}$  is the label of choice due to its high  $\beta$ -energy and the elevated specific activity of  $^{32}\text{P}$ -labelled probes. Labelled nucleotides with a specific activity of  $>110 \text{ TBq/mmol}$ ,  $>3000 \text{ Ci/mMol}$  are available for these applications.

**$^{33}\text{P}$**  With the advantage of a lower emission energy,  $^{33}\text{P}$  offers increased resolution as compared to  $^{32}\text{P}$ . It may be used for membrane hybridizations, but its lower energy emission means that longer exposure times are required. It is ideally suited to sequencing and *in situ* hybridization applications.

**$^{35}\text{S}$**  This label can be used for hybridizations but is not recommended because of its low emission energy. However, it provides high resolution and is suited to sequencing and *in situ* hybridization applications.

## Choice of labelling method

**Random prime** The most widely used method for the uniform labelling of DNA to provide high label density probes.

Ready-To-Go™ DNA Labelling Beads, Rediprime™ II, Megaprime™ and Multiprime™ are all based on the random prime technique which uses random sequence hexamers or nonamers to prime DNA synthesis on a denatured DNA template at numerous sites along its length.



Figure 22  
Ready-To-Go pre-dispensed  
reagent beads

The reaction uses the Klenow fragment of DNA polymerase I and leads to an efficient use of labelled nucleotides and net synthesis of probe. Therefore, very small amounts of input DNA are required, enabling high specific activity DNA probes to be produced with relatively small quantities of added label.

**Nick translation** This method provides uniformly labelled DNA probes by the introduction of random nicks into the template by the action of DNase I and subsequent filling in by DNA polymerase I. The existing nucleotide sequence of the DNA probe is renewed without net synthesis occurring. Nick translation can be used for the production of large amounts of high-specific activity probes.

**End labelling** This is the method of choice for labelling oligonucleotide probes. Either the 3' or 5' end can be labelled by the use of different techniques. End-labelled oligonucleotide probes may be used for membrane hybridizations. Probes labelled with  $^{35}\text{S}$  and  $^{33}\text{P}$  are suitable for use with *in situ* hybridization.

**RNA labelling** This technique produces radiolabelled RNA probes from inserts cloned into appropriate vectors. SP6 and T7 RNA polymerase may be used to produce asymmetric probes from the same or different vectors. Asymmetric probes are used in *in situ* hybridization applications to provide both positive results and negative controls.

Figure 23 Redivial™ showing removal of the splash guard which secures the radioactive solution at the bottom of the vial during ambient shipment. The splash guard is easily removed and disposed of as active waste.



Labelling system	Technology	Nucleotide	Amount of template	Labelling time	Probe specific activity (dpm/ $\mu\text{g}$ )	Recommended application
Rediprime II	Random-Prime	dCTP only	25 ng	10 minutes	$2 \times 10^9$	membrane hybridization
Ready-To-Go DNA labelling beads	Random-Prime	dCTP only	10 ng - 1 $\mu\text{g}$	5 minutes	$2 \times 10^9$	membrane hybridization
Megaprime	Random-Prime	any dNTP	25 ng	10 minutes	$2 \times 10^9$	membrane hybridization
Nick translation	Nick translation	any dNTP	1 $\mu\text{g}$	2-3 hours	$2 \times 10^9$	Production of large amounts of probe
RNA labelling	SP6/T7 RNA polymerase	UTP	1 $\mu\text{g}$	1-2 hours	$2 \times 10^9$	<i>in situ</i> hybridization
5'-end labelling	T4 Polynucleotide kinase	dATP	10 pmol ends	1 hour	$5 \times 10^6$	membrane hybridization <i>in situ</i> hybridization
3'-end labelling	Terminal deoxynucleotidyl transferase	any dNTP	10 pmol ends	30-60 minutes	$5 \times 10^6$	membrane hybridization <i>in situ</i> hybridization

Table 3 Guide to the properties of radioactive labelling systems

# Rediprime II

## DNA Labelling System

**Amersham Biosciences' premium radioactive labelling system consists of individually dispensed reaction mixes which are dried in the presence of a stabilizer and a dye to make labelling of probes easier. The system can be stored at 4°C or at room temperature ready for use.**

Rediprime II reaction mixes have been formulated using an improved exonuclease free Klenow to give probes with specific activities of  $1.9 \times 10^9$  dpm/mg or greater after 10 minutes incubation at 37°C with the majority of DNA substrates. When used with Redivue™ [ $^{32}\text{P}$ ]dCTP, Rediprime II reactions can be set up and completed to produce a DNA probe ready for hybridization in less than 15 minutes.

### Protocol Summary




Figure 24

Schematic diagram of the Rediprime II protocol

(A) Add denatured template to a final volume of 45 $\mu\text{l}$

(B) Flick tube and spin briefly

(C) Add 5ml Redivue [ $^{32}\text{P}$ ]dCTP, pipette up and down and spin briefly

(D) Incubate for 10 minutes at 37°C

## Recommended Applications

Labelling of DNA from a variety of sources to a high specific activity for use in Southern and Northern blot hybridizations. The system is designed for use with Redivue [ $^{32}\text{P}$ ]dCTP with a specific activity of 110TBq/mmol, 3000 Ci/mmol.

**Speed** Quick and convenient protocol requires the addition of template and labelled nucleotide only.

**Stability** Ambient temperature stable, therefore can be stored at room temperature.

**Flexibility** DNA can be labelled in the presence of Low Melting Point agarose or restriction enzyme buffers.

**Efficiency** Each labelling mix can label up to 25 ng of DNA to a specific activity of  $>10^9$  dpm/ $\mu\text{g}$ .




Figure 25

Northern blot probed with 3.8kb human EGFR cDNA fragments labelled with Rediprime. Result kindly supplied by J M Loughlin, Zeneca Pharmaceuticals, UK

## Ordering Information

Rediprime II DNA Labelling System	30 pre-mixed labelling reactions For use with radiolabelled dCTP	RPN1633
Rediprime II DNA Labelling System	60 pre-mixed labelling reactions For use with radiolabelled dCTP	RPN1634

### Related products

Redivue [ $\alpha$ - $^{32}\text{P}$ ]dCTP, 3000Ci/mmol	AA0005
Megaprime DNA Labelling System	RPN1604
Rapid-hyb Buffer	RPN1635
Hybond Nylon Membranes	see page 28

# 5'-End Labelling Kit

The 5'-End Labelling Kit exploits an optimized exchange buffer and T4 polynucleotide kinase to label oligonucleotides with or without a 5'-phosphate group. Labelling is complete within 30 minutes.

## Recommended Applications

Oligonucleotides labelled in this way are suitable for use in high target amount applications. For example as probes or primers to track PCR products, or in gel-shift and fragment analysis assays.

**Flexibility** Very small fragments of both DNA and RNA can be labelled.

**Reliability** Location of the labelled group is known.

## Protocol Summary




Figure 26  
The 5'-End Labelling reaction using T4 polynucleotide kinase


## Ordering Information

5'-End Labelling Kit	20 reactions	RPN1509
<b>Related products</b>		
Redivue $[\gamma\text{-}^{32}\text{P}]\text{dATP}$		AA0018
3'-End Labelling Kit		N4020

# 3'-End Labelling Kit

**Terminal deoxynucleotidyl transferase adds deoxyribonucleotides onto the 3'-ends of DNA fragments. It can be used in conjunction with  $^{32}\text{P}$ -,  $^{33}\text{P}$ -,  $^{35}\text{S}$ - or  $^3\text{H}$ -labelled nucleotides to label DNA for a variety of applications.**

## Protocol Summary



## Recommended Applications

One of the major applications for the 3'-End Labelling Kit is in the production of  $^{32}\text{P}$  end-labelled oligonucleotide probes for screening applications involving colony, plaque or PCR clones.  $^{35}\text{S}$  or  $^{33}\text{P}$  end-labelled probes are used for *in situ* hybridization applications.

**Flexibility** Template independent and all types of 3'-Ends can be labelled.

**Reliability** Location of the labelled group is known.

Figure 27  
The 3'-End Labelling reaction using terminal deoxynucleotidyl transferase

## Ordering Information

3'-End Labelling Kit	20 reactions	N4020
----------------------	--------------	-------

### Related products

Standard formulation radiolabelled nucleotides	RPN1509
5'-End Labelling Kit	RPN1509

# Ready-To-Go PCR Beads

**Ready-To-Go™ PCR Beads are pre-mixed, pre-dispensed, individual reactions designed for performing PCR amplification. Each dried room-temperature stable bead contains Taq® DNA polymerase, nucleotides and buffer, all optimized for standard PCR. Results are typically comparable with, or superior to, those using pre-formulated, aqueous 'master-mixes'.**



Figure 28  
Ready-To-Go pre-dispensed  
reagent beads

## Protocol Summary

Each pre-dispensed Ready-To-Go PCR Bead contains *Taq* DNA polymerase, nucleotides and buffer, all optimized for PCR.

Ready-To-Go PCR Beads are designed for 25 µl reactions. Simply add DNA template and primers and begin temperature cycling.

The beads are provided in either 0.5 ml or 0.2 ml tubes that are compatible with most thermocyclers. The 0.2 ml tubes come assembled in a convenient 96-well (8 x 12) plate format that allows individual strips of tubes to be easily removed. This flexibility allows the use of either the entire 96-well plate, strips of 8 or individual 0.2 ml tubes.

## Recommended Applications

Ready-To-Go PCR Beads can be used with a variety of templates including genomic DNA, viral DNA, plasmid DNA and cDNA. They provide an excellent source of probe template for both radioactive and non-radioactive labelling techniques. In addition low levels of radioactivity may be directly incorporated for use in applications such as single stranded conformational polymorphism (SSCP) assays.

\*For licence information see page 37

**Convenience** PCR components are provided in pre-formulated, single dose beads, reducing risk of pipetting errors and introducing contaminants. Pre-dispensed, individually packed PCR reactions, ensures a fresh reaction that dissolves quickly in a PCR compatible tube.

**Consistency** Application testing ensures each batch of beads delivers consistent results.

**Compatibility** Each reaction is provided in a thin-walled 0.5 ml or 0.2ml tube, which fit directly into most thermocyclers.




Figure 29  
Comparison of Ready-To-Go Beads with conventional PCR mixtures for the amplification of genes from interferon- $\gamma$  activated human aortic endothelial cells.  
Panel (A) amplification of TGF- $\alpha$   
Panel (B) amplification of  $\beta$ -actin.  
M = 100 Base Pair Ladder, beads = Ready-To-Go PCR Beads; mix = conventional PCR mixture; + = control reaction supplied with PCR beads.




Figure 30  
PCR results comparing performance of Ready-To-Go PCR Beads with pre-formulated, aqueous 'master mixes'. Both the PCR Beads and the 'master-mixes' contained all components for PCR except template and primer. 100 ng of human genomic template and primers specific for human aromatase, a single-copy gene, were added to each reaction to a final volume of 25  $\mu$ l. Reactions were subjected to 35 cycles of: 95°C for 1 minute, 55°C for 1 minute; 72°C for 2 minutes. An equal volume of each reaction was loaded onto an agarose gel. M = 100 base-pair ladder. Data courtesy of Denise Garvin.

## Ordering Information

Ready-To-Go PCR Beads (0.5 ml tubes)	100 reactions	27-9555-01
Ready-To-Go-PCR Beads (0.2 ml tubes/plate)	96 reactions	27-9553-01
Ready-To-Go-PCR Beads (0.2 ml tubes/plate)	5 x 96 reactions	27-9553-02

## Related products

Ready-To-Go DNA Labelling Beads	27-9240-01
---------------------------------	------------



# RNA Labelling Kit

The RNA Labelling Kit has been designed to give high-specific activity single-stranded RNA probes *in vitro*. This kit contains RNase-free DNase I for template removal and pre-prepared dithiothreitol (DTT) for added convenience. The kit comes with SP6 and T7 RNA polymerases and can also be used with T3 RNA polymerase, which must be purchased separately.



Figure 31  
Beta work tank

## Recommended Applications

For the generation of high specific activity RNA probes for membrane hybridization and *in situ* applications. Also full length RNA probes for RNase protection studies.

**Sensitivity** RNA probes show higher sensitivity than equivalent nick-translated probes in both Northern and *in situ* hybridizations.

**Stability** RNA:RNA hybrids are more stable than DNA:DNA or DNA:RNA hybrids, due to avoidance of reannealing and stronger hydrogen bonding.

**Specificity** Removal of non-specific bound probe by treatment with RNase A which is very specific for single-stranded RNA.

## Ordering Information

RNA Labelling Kit	20 labelling reactions up to 1 µg	RPN3100
-------------------	-----------------------------------	---------

### Related products

Redivue [ $\alpha$ - <sup>32</sup> P]UTP	AA0003
T3 RNA Polymerase	E70051Y

# Rapid-hyb Buffer

**Rapid-hyb™ Buffer is a rate enhancing hybridization buffer for rapid hybridization of radiolabelled nucleic acid probes to membrane-bound targets. In some Northern blotting experiments Rapid-hyb Buffer contributed to a five fold improvement in sensitivity.**

## Recommended Applications

Rapid-hyb Buffer is optimized for use in a wide range of applications, including Southern, Northern, dot/slot blots and colony/plaque lifts.

**Speed** Single-copy gene detection is possible after only a 2 hour hybridization with  $^{32}\text{P}$ -labelled probes.

**High signal to noise ratio** Inclusion of chemical blocking agents ensures low backgrounds.

**Stability** Stores at room temperature – ready to use without addition of carrier DNA.

**Compatibility** Compatible with DNA, RNA and oligonucleotide probes.

**Versatility** A wide range of hybridization temperatures (42°C – 70°C) can be used for optimal results.




Figure 32  
Northern blot analysis using:  
(A) Standard hybridization buffer (B) Rapid-hyb buffer  
Northern blots of HeLa cell total RNA (0.5 µg loadings). Linearized pHSP70 probe labelled with ( $\alpha$ - $^{32}\text{P}$ )dCTP using megaprime labelling system. Hybridizations were at 65°C for 1 hour using a probe concentration of 2 ng/ml. Exposure to Hyperfilm MP overnight.

## Ordering Information

Rapid-hyb Buffer	125 ml	RPN1635
Rapid-hyb Buffer	500 ml	RPN1636

## Related products

Rediprime II DNA Labelling System	RPN1633
5'-End Labelling Kit	RPN1509
3'-End Labelling Kit	N4020

# Redivue Nucleotides

Redivue™  $^{32}\text{P}$ ,  $^{33}\text{P}$  and  $^{35}\text{S}$  nucleotides contain a novel red dye and stabilizer solution, providing advantages in both handling and storage. Redivue is much easier to use being clearly visible in your pipette tip and in reaction mixes down to  $1\ \mu\text{l}$  in  $50\ \mu\text{l}$ . Redivue nucleotides can be used straight from the refrigerator in liquid form which saves time and avoids freeze thaw cycles.




Figure 34  
Redivue nucleotides can make experiments quicker to set up

## Recommended Applications

Redivue nucleotides can be substituted directly for standard formulation products in the majority of molecular biology applications.

**Ease of use** Addition of an intense red dye improves visibility and handling characteristics.

**Stability** Stabilized format allows storage at  $4^\circ\text{C}$  in a convenient liquid form.



Figure 33  
Use straight from the fridge. Dispensing is simpler with easy to use Redivue nucleotides

**Compatibility** Exhaustive testing of each Redivue product ensures equivalent performance to standard format radiolabelled nucleotides.

**Consistency** Elimination of repeated freeze thaw cycles ensuring consistent performance.




Figure 35  
Detection of B-myb mRNA by Northern hybridization using Redivue or standard  $[\alpha\text{-}^{32}\text{P}]\text{dCTP}$   
Results supplied by Dr E Lam, Ludwig Institute for Cancer Research, London

## Ordering Information

### Redivue Nucleotides

[ $\alpha$ - <sup>32</sup> P]dATP	~220 TBq/mmol, ~6000 Ci/mmol	AA0074
[ $\alpha$ - <sup>32</sup> P]dATP	~110 TBq/mmol, ~3000 Ci/mmol	AA0004
[ $\alpha$ - <sup>32</sup> P]dATP	~30 TBq/mmol, ~800 Ci/mmol	AA0084
[ $\alpha$ - <sup>32</sup> P]dATP	~15 TBq/mmol, ~400 Ci/mmol	AA0064
[ $\alpha$ - <sup>32</sup> P]dCTP	~220 TBq/mmol, ~6000 Ci/mmol	AA0075
[ $\alpha$ - <sup>32</sup> P]dCTP	~110 TBq/mmol, ~3000 Ci/mmol	AA0005
[ $\alpha$ - <sup>32</sup> P]dCTP	~30 TBq/mmol, ~800 Ci/mmol	AA0085
[ $\alpha$ - <sup>32</sup> P]dCTP	~15 TBq/mmol, ~400 Ci/mmol	AA0065
[ $\alpha$ - <sup>32</sup> P]dGTP	~110 TBq/mmol, ~3000 Ci/mmol	AA0006
[ $\alpha$ - <sup>32</sup> P]dGTP	~30 TBq/mmol, ~800 Ci/mmol	AA0086
[ $\alpha$ - <sup>32</sup> P]dGTP	~15 TBq/mmol, ~400 Ci/mmol	AA0066
[ $\alpha$ - <sup>32</sup> P]dTTP	~110 TBq/mmol, ~3000 Ci/mmol	AA0007
[ $\alpha$ - <sup>32</sup> P]dTTP	~30 TBq/mmol, ~800 Ci/mmol	AA0087
[ $\alpha$ - <sup>32</sup> P]dTTP	~15 TBq/mmol, ~400 Ci/mmol	AA0067
[ $\alpha$ - <sup>32</sup> P]dUTP	~110 TBq/mmol, ~3000 Ci/mmol	AA0003
[ $\gamma$ - <sup>32</sup> P]dATP	> 185 TBq/mmol, >5000 Ci/mmol	AA0018
[ $\gamma$ - <sup>32</sup> P]dATP	~110 TBq/mmol, ~3000 Ci/mmol	AA0068
[ $\gamma$ - <sup>33</sup> P]dATP	≥ 92.5 TBq/mmol, ≥ 2500 Ci/mmol	AH9968
[ $\alpha$ - <sup>33</sup> P]dATP	≥ 92.5 TBq/mmol, ≥ 2500 Ci/mmol	AH9904
[ $\alpha$ - <sup>33</sup> P]dCTP	≥ 92.5 TBq/mmol, ≥ 2500 Ci/mmol	AH9905
[ $\alpha$ - <sup>33</sup> P]dUTP	≥ 92.5 TBq/mmol, ≥ 2500 Ci/mmol	AH9903
[ <sup>35</sup> S]dATP $\alpha$ S	> 37 TBq/mmol, >1000 Ci/mmol	AG1000
[ <sup>35</sup> S]dATP $\alpha$ S	~22 TBq/mmol, ~600 Ci/mmol	AG1001
[ <sup>35</sup> S]dATP $\alpha$ S	~15 TBq/mmol, ~400 Ci/mmol	AG1002
[ <sup>35</sup> S]dCTP $\alpha$ S	> 37 TBq/mmol, >1000 Ci/mmol	AG0135
[ <sup>35</sup> S]dCTP $\alpha$ S	~22 TBq/mmol, ~600 Ci/mmol	AG0035
[ <sup>35</sup> S]dCTP $\alpha$ S	~15 TBq/mmol, ~400 Ci/mmol	AG0025

### Related products

Rediprime II DNA Labelling System	RPN1633
Ready-To-Go DNA Labelling Beads (dCTP)	27-9240-01
5'-End Labelling Kit	RPN1509
3'-End Labelling Kit	RPN4020
RNA Labelling Kit	RPN3100
Ready-To-Go PCR Beads	27-9555-01

# Hybond Membranes for Nucleic Acid Blotting

**Hybond™ membranes are available to support all major nucleic acid blotting applications.**

**Rely on Amersham Biosciences' technical experience in all types of blotting for clear guidance on the choice of membrane best suited for your application.**

**Individual technical advice is available on all aspects of blotting and hybridization via e:mail at**

**Hybond\_Help@eu.amershambiosciences.com**

## **Choosing the correct membrane**

The following table is designed to help in the choice of the best possible membrane for each application

<b>Application</b>	<b>Method</b>	<b>Recommended membrane</b>
Southern blotting	Radioactive detection	Hybond-XL
	Non-radioactive detection	Hybond-N+
	Alkali blotting	Hybond-XL / Hybond-N+
	Low volume hybridizations	Hybond-NX
	Rapid-hyb Buffer	Hybond-XL
	Fingerprinting	Hybond-Nfp
Northern blotting	Radioactive detection	Hybond-XL
	Non-radioactive detection	Hybond-N+
	Low volume hybridizations	Hybond-NX
	Rapid-hyb Buffer	Hybond-XL
Dot and slot blots	Radioactive detection	Hybond-XL
	Non-radioactive detection	Hybond-N+
Colony and plaque blots	Radioactive detection	Hybond-XL
	Non-radioactive detection	Hybond-N+

*Table 4. Choosing the correct membrane*

**Hybond-XL** Hybond-XL is a charged nylon membrane designed exclusively for radioactive nucleic acid transfer applications to achieve an improved signal-to-noise ratio. It has a greater concentration and a more even distribution of positive charge than other charged nylon membranes. Hybond-XL produces excellent results with a wide variety of applications and target/probe combinations. This membrane retains all the advantages of other nylon membranes such as high nucleic acid binding capacity and high tensile strength.

**Hybond-N+** Hybond-N+ is a positively charged nylon membrane which yields excellent sensitivity in both alkali blotting and conventional Southern blots. Nucleic acid samples may be fixed by simple alkali treatment or alkali blotting rather than UV exposure or baking, though UV fixation is recommended for maximum reproducibility.

**Hybond-N** Hybond-N is a neutral nylon membrane capable of high sensitivity in DNA and RNA blotting. It is a strong supported membrane that is inherently hydrophilic and requires no pre-wetting. It is easy to quickly crosslink nucleic acids to nylon using UV light. This makes Hybond-N ideal for any standard radioactive Southern or Northern blotting procedure, except for customers using low hybridization volumes where Hybond-NX is recommended.

**Hybond-NX** Hybond-NX is a neutral nylon membrane with all the properties of Hybond-N but has been specially developed for use with high throughput applications where low hybridization buffer volumes are used. Hybond-NX gives cleaner background than Hybond-N when low buffer volumes are used.




Figure 36  
Achieve cleaner results  
with Hybond-XL

**Hybond-Nfp** Hybond-Nfp is a neutral nylon membrane with all of the properties of Hybond-N, it was launched in response to the needs of researchers performing DNA fingerprinting. These customers preferred to validate a sample of a particular batch in their application before ordering the rest of the batch. Validation in this demanding technique can take up to 2 months, therefore to remove this inconvenience we offer Hybond-Nfp which is pre-validated by independent third party DNA fingerprinting experts.




Figure 37  
Human genomic Southern blot using:  
A) Hybond-XL B) leading competitor membrane  
Human genomic DNA, digested with Hind III 2, 1, 0.5, 0.1  $\mu$ g loadings.  
GAPDH random prime [ $^{32}$ P] labelled probe,  $2 \times 10^6$  cps/ml. Church  
hybridization buffer. Overnight exposure to Hyperfilm MP.

# Autoradiography Films

**Amersham Biosciences' expertise in the labelling and detection of biological molecules supports a range of autoradiography products designed specifically for life science applications.**

**This range includes the market leading Kodak™ X-Omat™ AR film and the high speed BioMax™ MR and BioMax MS films which utilize Kodak T-grain emulsion technology.**

Use the following table to help choose of the best possible film and exposure conditions for your application.

Application	Result required	Label	Method	Film
<b>Southern/Northern blots Colony/plaque blots</b>	High speed and sensitivity Accurate quantitation	<sup>32</sup> P	Pre-flash with screens at -70°C	Hyperfilm MP
	Maximum speed and sensitivity	<sup>32</sup> P	Screens at -70°C	BioMax MS
	Optimum resolution	<sup>33</sup> P/ <sup>35</sup> S	Direct autoradiography BioMax MR	Hyperfilm βmax
	Non-radioactive	ECL	Pre-flash and direct detection	Hyperfilm ECL
<b>Slot blots</b>	High speed and sensitivity accurate quantitation	<sup>32</sup> P	Pre-flash with screens at -70°C	Hyperfilm MP
	Maximum speed and sensitivity	<sup>32</sup> P	Screens at -70°C	BioMax MS
	non-radioactive	ECL	Pre-flash and direct detection	Hyperfilm ECL
<b>Dideoxy sequencing</b>	High speed	<sup>32</sup> P	Direct autoradiography	Hyperfilm MP BioMax MR
	Optimum resolution and maximum speed for difficult sequences	<sup>35</sup> S	Direct autoradiography	Biomax MR
	Good resolution and speed for routine work	<sup>35</sup> S	Direct autoradiography	Hyperfilm MP
<b>Cycle sequencing</b>	Optimum resolution and maximum speed for difficult sequences	<sup>33</sup> P	Direct autoradiography	BioMax MR
	Good resolution and speed for routine work	<sup>33</sup> P	Direct autoradiography	Hyperfilm MP
<b>In situ hybridization</b>	Macroscale <i>in situ</i> optimization of parameters	<sup>32</sup> P	Direct autoradiography	Hyperfilm MP




Figure 38 Rat kidney.

Total  $^{125}\text{I}$  ANP binding showing renal blood vessels, smooth muscle, renal papilla, glomeruli in cortex and the inner medulla. Reversed image from Hyperfilm  $^3\text{H}$ .

Photograph courtesy of John Wharton and Richard Rutherford. Department of Histochemistry, Royal Postgraduate School, London.

**Hyperfilm MP** An excellent multipurpose film. The clear base gives high contrast and good resolution while its double-sided construction and compatibility with intensifying screens means that it performs well in blotting applications.

**Hyperfilm  $\beta$ -max** A specialized film with high sensitivity, high resolution and a good linear range.

**Hyperfilm ECL** The recommended film for chemiluminescent applications. The clear base gives high contrast and good resolution while its sensitivity profile matches the emission spectrum of ECL and other chemiluminescent systems.

**BioMax MR** An excellent film for the direct autoradiography of  $^{35}\text{S}$ ,  $^{14}\text{C}$  and  $^{33}\text{P}$ . This clear, single-sided film gives high resolution and contrast and the special T-grain technology makes it twice as fast as X-Omat AR with weak emitters.

**BioMax MS** This film and screen combination offers exceptional sensitivity in  $^{32}\text{P}$  and  $^{125}\text{I}$  applications. The high performance film of choice for radioactive blotting.

**X-Omat AR** A general-purpose film for use with all commonly used isotopes. Coated with emulsion on both sides for both direct exposure and use with intensifying screens.

## Related autoradiography accessories

Hypercassette™ autoradiography cassettes  
Hyperprocessor™ automatic film processor  
Hyperscreen™ intensifying screens  
Hypertorch™  
Kodak BioMax TranScreen intensifying screens  
Safelights  
Sensitize™ preflash unit  
Trackertape™

## Ordering information

For further information please request 'The Complete Autoradiography Product Guide'.




Figure 39  
Cycle sequencing with  $^{33}\text{P}$  labelled primers. Detected on Hyperfilm MP.




# Amersham Biosciences Scanners

## Storm

The Amersham Biosciences™ Storm™ imaging system for gel, blot and macroarray analysis combines proven storage phosphor imaging technology for radioisotope detection with fluorescence and chemifluorescence imaging for non-radioactive detection. Researchers can make the move to non-radioactive gel and blot analysis without sacrificing conventional radioisotope techniques. Storm scans storage phosphor screens and fluorescent and chemifluorescent gels and blots in a large format 35 × 43cm scanning area or smaller. ImageQuant™ software for Windows NT™ or Macintosh™ is included. ECF detection reagents provide optimized chemifluorescent detection of Southern and Northern blots on the Storm system.

Place your storage phosphor screen, gel or blot on the sample area, and point and click to start your scan. The Storm system uses patented Variable Mode Imaging (VMI) architecture to activate the correct optical components for storage phosphor autoradiography, direct fluorescence or chemifluorescence.



*Figure 41*  
Large samples fit on the Storm systems 35 x 43cm scan area. Storm offers the high resolution needed for DNA base identification.



*Figure 40*  
Storm imaging system in use

## PhosphorImager SI

The **Amersham Biosciences PhosphorImager™ SI** offers advanced technology for autoradiography in a compact design. Storage phosphor screens capture images produced by x-rays and radioisotopes from radiolabelled blots, gels, arrays, TLC plates or tissues. Exposure times are reduced by at least 90% relative to an equivalent exposure to autoradiography film and the PhosphorImager SI can accurately quantitate both strong and weak signals from the same exposure. The PhosphorImager SI offers significant advantages in data throughput and analysis.




Figure 43  
Southern blotting with chemifluorescence. Bands were visualized by chemifluorescence in this Southern blot of human genomic DNA.

## FluorImager

FluorImager™ systems are multi-purpose fluorescent laser scanners for multiple colour analysis of gels, blots, arrays, TLC, and microplates. Samples are read directly, with no drying or exposure step required. Software for instrument control, quantitative image analysis and printing is included.




Figure 44  
Northern blotting with chemifluorescence. A fluoresceinated cDNA probe was used to identify p53 mRNA in A172 glioblastoma cells. Bands were visualized by chemifluorescence.



Figure 42  
FluorImager 595 provides fast, fluorescent gel and blot analysis

# Accessory products

## UV Crosslinker

The Amersham Biosciences™ UV Crosslinker reduces the time required to fix nucleic acids on membranes from hours to less than a minute.

**UV-fixed samples give consistently better hybridization signals than baked samples, especially after repeated re-probings. Nucleic acids can be UV-crosslinked to nylon, nitrocellulose and supported nitro-cellulose membranes.**

**Reliability** Short-wave UV light crosslinks DNA or RNA securely to membranes for repeated re-probings.

**Accuracy** Microprocessor controlled for accurate and reproducible UV dosage.

**Convenience** Keypad selection of pre-programmed or user-set energy or time exposures.

**Flexibility** Spacious chamber accepts membranes up to 25 x 30cm, but compact footprint requires only 35 x 40cm of bench space.



Figure 45  
The Amersham Biosciences UV  
Crosslinker

## Ordering Information

Ultraviolet Crosslinker 115 V 60 Hz Includes five 8 W short-wave UV lamps	80-6222-31
Ultraviolet Crosslinker 230 V 50 Hz Includes five 8 W short-wave UV lamps	80-6222-50
8 W, 254 nm replacement UV lamp	80-6223-64

## Vacuum Blotting Units

**Vacuum blotting is a quick and reproducible method for transferring nucleic acid to membranes. The specially designed vacuum pump draws the DNA/RNA through the agarose gel and onto the membrane, therefore reducing transfer times to less than an hour.**

**Speed** Rapidly transfer DNA or RNA from agarose gels up to 20 x 30 cm in size.

**Efficiency** Efficient and reproducible transfers improve blotting results.

**Convenience** Minimize gel handling by performing all pre-transfer steps in the unit.

**Accuracy** Use with VacuGene™ XL blotting pump ensures precise vacuum control.

## Ordering Information

VacuGene XL Vacuum Blotting System	13-2223-01
VacuGene XL Vacuum Blotting Unit	80-1266-24
VacuGene XL Vacuum Blotting Pump 110 V	80-1265-14
VacuGene XL Vacuum Blotting Pump 220 V	80-1265-15

## Submarine Electrophoresis Units

**Electrophoresis is an integral part of many molecular biology techniques, and as such requires equipment which is convenient, robust and adaptable.**

**The comprehensive range of submarine electrophoresis equipment from Amersham Biosciences can accommodate virtually any need for gel size, buffer volume and cooling requirement.**

### Amersham Biosciences HE 33 Mini Submarine Electrophoresis Unit

Coolant filled buffer chamber base provides passive cooling, allowing gels to be run at high voltages. 7 x 10cm agarose gels can be run in less than 20 minutes.

### Amersham Biosciences HE 99X Max Submarine Electrophoresis Unit

Versatile and durable, 3 different gel lengths and 12 different comb sizes allows customization of gel configuration. Gels are cast without the need for tape.

### Amersham Biosciences HE 100 SuperSub™ Submarine Electrophoresis Unit

High performance, large format electrophoresis unit. Contains a built-in magnetic stir bar for buffer circulation and an integral heat exchanger allowing high voltage or extended runs.

### EPS 301 Power Supply

Flexible, economical power supply for submarine electrophoresis units at up to 300V or 400mA with 80W power. This versatile unit also handles vertical gels and some electrophoretic blotting.



Figure 46  
EPS301 Power Supply

Figure 47

HE 33 Mini Submarine Electrophoresis Unit gives enhanced passive cooling for faster runs.





Figure 48

Separation of DNA size markers in HE 33 Submarine Unit. Lanes 1, 4, 7:  $\phi$ X-174 RF DNA-Hind II digest; lanes 2, 5, 8: lambda DNA-Hind III digest; lanes 3, 6:  $\phi$ X-174 RF DNA-Hae III digest. Markers were separated for 30 minutes at 200V in a 1% agarose gel.

## Ordering Information

Amersham Biosciences HE 33 Mini Submarine Electrophoresis Unit, complete	80-6052-45
Amersham Biosciences HE 99 Max Submarine Electrophoresis Unit, complete	80-6061-57
Amersham Biosciences HE 100 SuperSub Electrophoresis Unit, complete	80-6043-71
EPS 301 Power Supply	18-1130-01

# Accessory products

## Hybridization Oven/Shaker

The Hybridization Oven/Shaker rigidly controls the hybridization temperature and frequency of agitation, thus ensuring that stringency and washing conditions are maintained. Developed for use with both radioactive and non-radioactive labelling systems. A choice of hybridization modes is available by use of a variable speed rotisserie or platform shaker.

**Economy** Reduced hybridization buffer volumes minimizes probe usage.

**Safety** Acrylic/polycarbonate door minimizes exposure to beta radiation.

**Versatility** Interchangeable rotisserie and platform shaker.



Figure 49  
Hybridization Oven/Shaker

## Ordering Information

Hybridization Oven/Shaker 220/240 V 50 Hz	RPN2510
Hybridization Oven/Shaker 110/115 V 60 Hz	RPN2511
Hybridization Oven/Shaker 100 V 50/60 Hz	RPN2512

## Electrophoresis Chemicals

In gel electrophoresis the quality of the chemicals used for gel casting and gel staining is critical to obtain good results.

**Ethidium Bromide** Widely used for detecting DNA and RNA, ethidium bromide binds to double-stranded regions to form UV fluorescent complexes. As little as 50 ng of DNA can be detected with concentrations of 0.5 – 1.0 µg/ml of stain.

**Agarose** Agarose with low electroendosmosis is recommended for most DNA and RNA applications. To ensure maximum enzyme efficiency, a highly purified Agarose NA can be used. For preparative work, Agarose Prep with low gelling temperature should be used.

## Ordering Information

Ethidium Bromide Solution, 10 mg/ml	17-1328-01
Agarose NA, 100 g	17-0554-02
Agarose Prep, 50 g	80-1103-07

## Trademarks

AlkPhos Direct, ECF, ECL, ECL Direct, FluorImager, Gene Images, Amersham Biosciences, Hybond, Hypercassette, Hyperfilm, Hyperprocessor, Hyperscreen, Hypertorch, ImageQuant, Megaprime, **Amersham Biosciences**, Multiprime, PhosphorImager, Rapid-hyb, Ready-To-Go Rediprime, Redivial, Redivue, Sensitize, Storm, SuperSub, TrackerTape, TransVac, Vacublott, VacuGene are trademarks of Amersham **Biosciences** Limited or its subsidiaries.

Amersham and Amersham Biosciences are trademarks of Amersham plc

CDP-Star is a trademark of Tropix Inc

Kodak, BioMax, TranScreen, and X-Omat are trademarks of Eastman Kodak Co.

Windows NT is a trademark of Microsoft Corp.

Macintosh is a trademark of Apple Computers.

## Licence information

CDP-Star is protected by Tropix Inc. under one or more of US patent numbers 5326882 and 4931569

ECF substrate is manufactured for Amersham **Biosciences** Limited by JBL Scientific Inc, and is sold under licence from JBL Scientific Inc.

ECL Direct/AlkPhos Direct labelling reagents are covered by European patent number 120376B1 and sold under licence from European Molecular Biology Laboratory. ECL Direct (RPN3000/3001/3005) is covered by European patent numbers 373017 and 513040 and foreign equivalents and is sold under licence from Institut Pasteur.

Gene Images hybridization buffer is covered by US patent number 5512436 and foreign equivalents.

PCR: *Taq* DNA polymerase is sold under licensing arrangements with Roche Molecular Systems, F Hoffmann-La Roche Ltd and the Perkin Elmer Corporation. Purchase of these products is accompanied by a limited licence to use it in the polymerase Chain Reaction (PCR) process for research in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front licence fee, either by payment to Perkin-Elmer or as purchased, i.e. an authorized thermal cycler.

The Polymerase Chain Reaction (PCR) is covered by patents owned by Roche Molecular Systems and F Hoffmann-La Roche Ltd. A licence to use the PCR process for certain research and development activities accompanies the purchase of certain reagents from licensed suppliers such as Amersham **Biosciences** Limited and affiliates when used in conjunction with an authorized thermal cycler.

© Amersham **Biosciences** UK Limited 1998 - All rights reserved

All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham **Biosciences** group which supplies them. A copy of these terms and conditions is available on request.

**Amersham Biosciences UK Limited** Amersham Place Little Chalfont Buckinghamshire England HP7 9NA

**Amersham Biosciences AB** S-751 82 Uppsala Sweden

**Amersham Biosciences Inc** 800 Centennial Avenue PO Box 1327 Piscataway NJ 08855 USA

For more information about  
Nucleic Acid Labelling and Detection  
visit our website at:

**[www.amershambiosciences.com](http://www.amershambiosciences.com)**

**Asia Pacific**

Amersham Biosciences Asia Pacific Ltd

Tel: +852 2811 8693

Fax: +852 2811 5051

**Australasia**

Amersham Biosciences Australia Pty Ltd

Tel: +61 2 9894 5188

Fax: +61 2 9899 7511

**Austria**

MEDPRO

Tel: +43 1479 15 690

Fax: +43 1479 83 5217

**Belgium**

Amersham Biosciences Benelux

Tel: 0800 73888

Fax: +32 3272 1637

**Canada**

Amersham Biosciences, Inc.

Tel: 1 800 463 5800

Fax: 1 800 567 1008

**Central Europe**

Amersham Biosciences Export GmbH

Tel: +43 1 982 3826

Fax: +43 1 985 8327

**Denmark**

Amersham Biosciences

Tel: +33 4516 24000

Fax: +33 4516 2424

**Finland**

Amersham Biosciences Europe GmbH

Tel: 358 (9) 512 3940

Fax: 358 (9) 512 1710

**Russian Federation**

ZAO Amersham Biosciences

Tel: +7 (950) 232 0250/956 1137

Fax: +7 (950) 232 6377

**France**

Amersham Biosciences

Tel: +33 169 18 28 00

Fax: +33 169 18 28 28

**Germany**

Amersham Biosciences GmbH

Tel: +49 0761 4903 0

Fax: +49 0761 4903 247

**Italy**

Amersham Italia Srl

Tel: +39 02 27322 1

Fax: +39 02 27302 212

**Japan**

Amersham Biosciences KK

Tel: 81 3 5992 2828

Fax: 81 3 5992 5186

**Latin America**

Amersham Biosciences do Brasil Ltda

Tel: +55 11 3667 5700

Fax: +55 11 3667 58 99

**Middle East and Africa**

Amersham Biosciences Export GmbH

Tel: +30 (1) 96 00 687

Fax: +30 (1) 96 00 693

**Netherlands**

Amersham Biosciences Benelux

Tel: 0165 580 410

Fax: 0165 580 401

**Norway**

Amersham Biosciences

Tel: 63 89 23 10

Fax: 63 89 23 15

**South East Asia**

Amersham Biosciences AB

Tel: 60 3 724 2080

Fax: 60 3 724 2090

**South East Europe**

Amersham Biosciences GmbH

Tel: +43 (1) 982 3826

Fax: +43 (1) 985 8327

**Spain**

Amersham Biosciences Europe GmbH

Tel: 93 594 49 82

Fax: 93 594 49 55

**Sweden**

Amersham Biosciences AB

Tel: 018 164000

Fax: 018 712777

**Switzerland**

Amersham Biosciences Europe GmbH

Tel: 01 802 81 50

Fax: 01 802 81 51

**UK**

Amersham Biosciences UK Limited

Tel: 0800 515313

Fax: 0800 616927

**USA**

Amersham Biosciences Inc

Tel: +1 800 526 3593

Fax: +1 800 329 3593

