

MultiWash+ Microplate Washer

A 96- and 384-well automated washer

BENEFITS

- **Walk-up usability**
- **Start simply with wash and rinse bottles (included)**
- **Flexible system washes 96- or 384-well plates**
- **Automated rinse reduces clogging**
- **Quiet and efficient**

Introduction

MultiWash+ washer is an automated microplate washer that provides maximum flexibility with a minimum footprint. The MultiWash+ washer can be configured for both 96- and 384-well plates and comes with four wash/rinse bottles and one waste bottle for out-of-the-box use.

Wash protocols & modes


The on-board touch panel comes with the option of 20 different wash protocols, including up to 8 cycles within each protocol for walk-up usability and flexibility. Washing variations include adjustable speed and volume, adjustable aspiration speed and time, adjustable soak times and three modes of shaking are available to mix solutions.

Efficiency


Efficient washing is achieved with cross-wise aspiration, reducing residual volume within each well. The washer has a small footprint and is vacuum and pressure free, with on-board pumps that provide a quiet wash experience.

Mouse/Rat IL-22 Quantikine ELISA using the MultiWash+ washer

A Mouse/Rat IL-22 Quantikine ELISA from R&D systems was used to evaluate the MultiWash+ washer. An IL-22 standard curve was prepared and a sandwich ELISA was performed using the MultiWash+ plate washer in strip mode to wash the wells. After reading the ELISA plate on the EMax Plus Microplate Reader, the R2 value was 0.984. A control well was set up to verify the accuracy of the standard curve. The standard was within the range described with the kit (data not shown).



MultiWash+ washer can dispense in 50 µl increments and tubing can be inserted into a variety of containers since pressurized bottles are not required.



ELISA Assay

Technical specifications

Plate wash time (300 µl)	<ul style="list-style-type: none"> • 96-well plate, 1 wash: 1 min. 9 sec. • 96-well plate, 3 washes: 3 min. 35 sec.
Dispense volume	50 – 2000 µL in 50 µL increments
Residual volume	< 1 µL per well
Dead volume	15.62 mL
Precision	< 5% at 300 µL across the plate
User interface	2 X 16 characters backlit display keyboard with 5 function keys
Shaking	3 speeds
Microplates	<ul style="list-style-type: none"> • 96-well plates flat and round bottom (384-well plates with 16-manifold) • Fits plates with standard ANSI SBS footprint
Power requirement	90 – 250B, 60V auto-sensing, 50/60 Hz
Instrument dimensions (WxHxD)	<ul style="list-style-type: none"> • 21 x 21 x 44 cm, • 8.3" x 8.3" x 17.3"
Weight	8Kg (17.8lbs)

For more information, please contact your VWR representative.

Go to www.vwr.com for the latest news, special offers and details of your local VWR organisation.

Austria VWR International GmbH Graumannsgasse 7 1150 Vienna Tel.: +43 1 97 002 0 Email: info.at@vwr.com	Belgium VWR International bvba Researchpark Haasrode 2020 Geldenaaksebaan 464 3001 Leuven Tel.: +32 (0) 16 385 011 Email: vwr.be@vwr.com	Denmark VWR International A/S Tobaksvejen 21 2860 Søborg Tel.: +45 43 86 87 88 Email: info.dk@vwr.com	France VWR International S.A.S. Le Périgares – Bâtiment B 201, rue Carnot 94126 Fontenay-sous-Bois cedex Tel.: 0 825 02 30 30* (national) Tel.: +33 (0) 1 45 14 85 00 (international) Email: info.fr@vwr.com * 0,18 € TTC/min	Germany VWR International GmbH Hilpertstraße 20a D - 64295 Darmstadt Tel.: 0800 702 00 07* (national) Tel.: +49 (0) 6151 3972 0 (international) Email: info.de@vwr.com *Freecall
Ireland / Northern Ireland VWR International Ltd / VWR International (Northern Ireland) Ltd Orion Business Campus Northwest Business Park Ballycoolin Dublin 15 Tel.: +353 (0) 1 88 22 222 Email: sales.ie@vwr.com	The Netherlands VWR International B.V. Postbus 8198 1005 AD Amsterdam Tel.: +31 (0) 20 4808 400 Email: info.nl@vwr.com	Sweden VWR International AB Fagerstagatan 18a 163 94 Stockholm Tel.: +46 (0) 8 621 34 00 Email: kundservice.se@vwr.com	Switzerland VWR International GmbH Lerzenstrasse 16/18 8953 Dietikon Tel.: +41 (0) 44 745 13 13 Email: info.ch@vwr.com	UK VWR International Ltd Customer Service Centre Hunter Boulevard - Magna Park Lutterworth Leicestershire LE17 4XN Tel.: +44 (0) 800 22 33 44 Email: uksales@vwr.com