

GE Healthcare

Amersham Hybond-XL

a nylon membrane designed for nucleic acid transfer

Product Booklet

Codes:	RPN82S	82 mm
	RPN87S	87 mm
	RPN132S	132 mm
	RPN137S	137 mm
	RPN119S	11.9 x 7.8 cm
	RPN1210S	12 x 10 cm
	RPN1510S	15 x 10 cm
	RPN1520S	15 x 20 cm
	RPN2020S	20 x 20 cm
	RPN2222S	22 x 22 cm
	RPN3050S	30 x 50 cm
	RPN203S	20 cm x 3 m
	RPN303S	30 cm x 3 m



Page finder

1. Legal	3
2. Handling	4
2.1. Safety warnings and precautions	4
2.2. Storage	4
2.3. Stability	5
2.4. Expiry	5
3. Components	6
4. Other materials required	7
5. Critical parameters	8
6. Description	9
7. Solutions required for nucleic acid electrophoresis and blotting	10
8. Quality control	13
9. Blotting protocols	14
9.1. Protocol for capillary blotting	14
9.2. Southern blotting – gel treatments (Neutral transfer protocol)	16
9.3. Southern blotting - gel treatments (Alkali transfer protocol)	19
9.4. Northern blotting - gel preparation and treatment	20
9.5. Protocol for colony and plaque lifts (10)	23
9.6. Protocol for dot blotting (manual)	26
9.7. Hybridization protocols – hybridizations in bags and boxes	29
9.8. Hybridization in tubes	33
9.9. Stripping protocols	39
10. Additional information	41
10.1. Determination of the optimum UV crosslinking conditions using a UV transilluminator	41
10.2. Determination of the optimum alkali fixation	42
11. Related products	44
12. References	46

1. Legal

GE and GE monogram are trademarks of General Electric Company. Ficoll, Hybond, Hyperfilm, Megaprime, ProbeQuant, Rediprime, Redivue, SeedDNA and SepRate are trademarks of GE Healthcare companies.

AnalaR is a trademark of Merck

SybrGreen is a trademark of Molecular Probes Inc

Gilson is trademark of Gilson Medical Equipment

SaranWrap is a trademark of Dow Chemical Company

Kodak and BioMax are trademarks of Eastman Kodak Company.

© 2006 General Electric Company – All rights reserved.

GE Healthcare reserves the right, subject to any regulatory and contractual approval, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

Contact your GE Healthcare representative for the most current information and a copy of the terms and conditions.

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited.

Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA UK

2. Handling

2.1. Safety warnings and precautions

Warning: For research use

only. Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

We recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. As all chemicals should be considered as potentially hazardous, it is advisable when handling chemical reagents to wear suitable protective clothing, such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In case of contact with skin or eyes wash immediately with water.

Note that the procedures requires the use of:

Sodium Dodecyl Sulphate: irritant

Formaldehyde: toxic substance

Formamide: toxic substance

Ethidium Bromide: mutagenic substance

Sodium Hydroxide: corrosive

Hydrochloric Acid: corrosive

Diethylpyrocarbonate:

explosive, toxic substance

This product may also be used with radioactive materials.

Please follow the manufacturers' safety data sheets relating to the safe handling and use of these reagents.

2.2. Storage

Membranes should be stored in a clean, dry atmosphere away from excessive heat, light and noxious fumes. The membranes should be handled wearing gloves or using blunt ended forceps to prevent contamination.

2.3. Stability

Membranes are stable for up to three years. Membranes should be kept in the bags in which they are received. Performance is consistent when stored under the recommended conditions.

2.4. Expiry

Expiry information can be found on the external packaging.

3. Components

Hybond-XL membranes

RPN82S	82 mm diam, 50 discs
RPN87S	87 mm diam, 50 discs
RPN132S	132 mm diam, 50 discs
RPN137S	137 mm diam, 50 discs
RPN119S	11.9 x 7.8 cm, 50 sheets
RPN1210S	12 x 10 cm, 20 sheets
RPN1510S	15 x 10 cm, 20 sheets
RPN1520S	15 x 20 cm, 10 sheets
RPN2020S	20 x 20 cm, 10 sheets
RPN2222S	22 x 22 cm, 10 sheets
RPN3050S	30 x 50 cm, 5 sheets
RPN203S	20 cm x 3 m, 1 roll
RPN303S	30 cm x 3 m, 1 roll

4. Other materials required

Equipment

- Agarose gel electrophoresis apparatus, for example GE Healthcare. HE33 Mini or HE99X Max submarine gel electrophoresis systems
- Microwave
- Hybond blotting paper
- Absorbent paper towels
- Trays/dishes
- Glass plates
- 750 g weight
- Pipettes, for example, Gilson™ Pipetman P20, P200, P1000 and P5000.
- Assorted laboratory glassware
- Oven, at 80°C or UV transilluminator
- Orbital shaker
- SaranWrap™ or similar cling film

5. Critical parameters

- Storage

Membranes should be stored in a clean, dry atmosphere away from excessive heat, light and noxious fumes.

- Handling

The membranes should be handled wearing gloves or using blunted forceps to prevent contamination. All membranes should be cut using a clean sharp scissors to avoid damage to the membrane edges.

- Wettability

Nylon membranes are hydrophilic and do not require pre-wetting before use in blotting procedures. Wetting is however advised for large blots ($>100 \text{ cm}^2$) or when multiple blots are hybridized together. Wet the membrane first in water then equilibrate in an appropriate buffer.

- Fixation

The fixation procedure can significantly affect the eventual sensitivity of a system. Sub optimal fixation reduces the amount of available target sequences, particularly following stripping. Nucleic acid may be fixed using heat or UV crosslinking, or alkali. It is essential that fixation times, energy settings (UV) and concentration (alkali) are optimized.

6. Description

Hybond™-XL is a novel charged nylon membrane from GE Healthcare. It has been designed exclusively for nucleic acid transfer applications using radioactive detection to achieve an improved signal to noise ratio. The design process used ensures a robust product capable of producing excellent results with the wide variety of procedures, applications and target/probe combinations used in life science laboratories world-wide. This new membrane retains all of the advantages of other nylon membranes:

- high binding capacity
- superior binding of small fragments
- excellent handling properties
- effective probe removal
- efficiency of reprobing

The membrane, which is identical on both sides, is manufactured in long rolls known as “master rolls”. Production runs are carefully controlled and constantly monitored to ensure the most consistent product reaches the user. Samples are taken from the beginning, middle and end of each master roll and used in single copy gene detection (see Quality Control section page 13).

7. Solutions required for nucleic acid electrophoresis and blotting

All reagents should be of AnalaR™ grade where possible.

10 x nucleic acid loading dye mix

40 mg Bromophenol blue
40 mg Xylene cyanol
2.5 g Ficoll™ 400

Add approximately 8 ml of distilled water. Mix to dissolve. Make up to a final volume of 10 ml. Store at room temperature for up to 3 months.

50 x TAE (DNA electro-phoresis buffer)

242 g Trizma base
18.6 g Ethylenediaminetetra-acetic acid (EDTA), sodium salt

Add approximately 800 ml of distilled water. Mix to dissolve. Adjust to pH 8 with glacial acetic acid (~57 ml/l). Make up to a final volume of 1000 ml. Store at room temperature for up to 3 months.

Depurination solution

11 ml Hydrochloric acid, concentrated (HCl)
989 ml Distilled waterMix. Store at room temperature for up to 1 month.

Denaturation buffer

87.66 g Sodium Chloride (NaCl)
20 g Sodium hydroxide (NaOH)
Add approximately 800 ml of distilled water. Mix to dissolve. Make up to a final volume of 1000 ml. Store at room temperature for up to 3 months.

Neutralization buffer

87.66 g Sodium chloride NaCl
60.5 g Trizma base
Add approximately 800 ml of distilled water. Mix to dissolve. Adjust to pH 7.5 with concentrated hydrochloric acid. Make up to a final volume of 1000 ml. Store at room temperature for up to 3 months.

Nucleic acid transfer buffer (20 x SSC)

88.23 g Tri-sodium Citrate
175.32 g Sodium chloride (NaCl)
Add approximately 800 ml of distilled water. Mix to dissolve. Check the pH is 7–8. Make up to a final volume of 1000 ml. Store at room temperature for up to 3 months.

TE buffer

1.21 g Trizma base
0.372 g Ethylenediamine tetra-acetic acid (EDTA), sodium salt
Add approximately 800 ml of distilled water. Mix to dissolve. Adjust to pH 8 with concentrated Hydrochloric Acid. Make up to a final volume of 1000 ml. Store at room temperature for up to 3 months.

10 x MOPS buffer

41.2 g 3-(N-morpholino)propanesulphonic Acid (MOPS)
10.9 g Sodium Acetate, 3-hydrate
3.7 g Ethylenediamine tetra-acetic acid (EDTA), sodium salt
Add approximately 800 ml of nuclease free distilled water. Mix to dissolve. Adjust to pH 7 with Sodium hydroxide (NaOH) (prepared in nuclease free distilled water).

10 x MOPS buffer

Continued

Make up to a final volume of 1000 ml.
Filter sterilize. Store at room temperature protected from light. **Do not** use if the solution appears yellow in color.

100 x Denhardt's solution

2.0 g Bovine Serum Albumin
2.0 g Ficoll 400
2.0 g Polyvinylpyrrolidone
Add approximately 50 ml of distilled water.
Mix to dissolve. Make up to a final volume of 100 ml. Store at -15°C to -30°C for up to 3 months.

8. Quality control

Hybond-XL is QC tested using related GE Healthcare products and protocols to ensure maximum compatibility and optimum performance.

Date for individual batches and detailed QC protocols are always available on request.

Description:

Nylon hybridization transfer membrane.

Quality control assay:

Hind III restricted human genomic DNA, separated using SepRate™ -DNA agarose, is Southern blotted onto Hybond-XL and hybridized with N-*ras* proto-oncogene probe.

Specifications:

Detection of 0.1 pg of target DNA. Background below 0.15 OD. Hybridization volume 35 µl/cm²

Labelling and detection:

Performed using : Megaprime™ random prime labelling kit with Redivue™ ³²P-dCTP label, radioactive signal detected using Hyperfilm™ MP.

Storage instructions:

Store in a clean dry environment.

9. Blotting protocols

This instruction for use limits itself to the classical capillary blotting technique (1) used for the transfer of separated nucleic acid fragments from an agarose gel to a solid support (2, 3, 4, 5, 6), and is representative of the procedures used in GE Healthcare laboratories.

Details of the gel treatments required before the transfer of DNA or RNA may be found on pages 16 and 20. Figure 1 (page 17) is a diagrammatic representation of the transfer apparatus used in this technique.

9.1. Protocol for capillary blotting

Protocol

1. Prepare the gel for transfer.

Notes

1. **Details of gel treatments may be found on page 16 (Southern blotting) or page 20 (Northern blotting)**

2. **Cut a sheet of membrane to an appropriate size.**

2. The membrane should be cut with clean sharp scissors.

3. Half fill a tray or glass dish of a suitable size with the transfer buffer. Make a platform and cover with a wick made from three sheets of Hybond blotting paper saturated in transfer buffer.

3. At least 800 ml of buffer is required for a 20 x 20 cm gel and a dish 24 x 24 cm. Ensure the wick ends are immersed in the transfer buffer.

4. Place the treated gel on the wick platform. Avoid trapping any air bubbles between the gel and the wick. Surround the gel with cling film to

4. Air bubbles block the transfer of nucleic acid to the membrane. They may be removed at any stage by

Protocol

4. *Continued.*

prevent the transfer buffer being absorbed directly into the paper towels.

5. Place the membrane on top of the gel. Avoid trapping any air bubbles.

6. Place three sheets of Hybond blotting paper cut to size and saturated in transfer buffer, on top of the membrane. Avoid trapping any air bubbles.

7. Place a stack of absorbent towels on top of the Hybond blotting paper at least 5 cm high.

8. Finally, place a glass plate and a weight on top of the paper stack. Allow the transfer to proceed overnight. The weight should not exceed 750 g for a 20 x 20 cm gel.

9. After blotting, carefully dismantle the transfer apparatus. Before separating the gel and membrane, mark

Notes

4. *Continued.*

rolling a clean pipette or glass rod over the surface.

5. Do not attempt to move the membrane once it has touched the gel surface.

8. Small fragments (0.5–1.5 Kb) are rapidly transferred upwards in a few hours, larger fragments (>10 Kb) require at least overnight transfer. The efficiency of transfer of these larger fragments can be improved by depurination (fragmentation).

9. Rinsing the membrane following transfer is not advised. Extensive experimentation has shown

Protocol

9. *Continued.*
the membrane to allow identification of the tracks with a pencil or chinagraph pen.

Notes

9. *Continued.*
that rinsing the membrane before fixation produces blots of variable quality because nucleic acid is removed from the membrane during this step.

10. Fix the nucleic acid to the membrane by baking at 80°C for 2 hours or by using an optimized UV crosslinking procedure.

10. Details of an optimization procedure are given on page 39. The UVC 500 crosslinker from GE Healthcare has a pre-set UV exposure (70 000 micro-joules/cm²) which is suitable for Hybond-XL.

11. Blots may be used immediately. Blots must be thoroughly dried if storage is required.

11. Blots may be rinsed in 2 x SSC before storage or hybridization. Blots should be stored wrapped in SaranWrap desiccated at room temperature under vacuum.

9.2. Southern blotting – gel treatments (Neutral transfer protocol)

Protocol

1. Separate the DNA samples on a suitable neutral agarose gel.

Notes

1. Efficient separation of a range of DNA fragments may be achieved by varying the type and concentration of the

Protocol

Notes

1. *Continued.*

agarose in the gel. Ensure the optimum DNA concentration for detection is loaded into each track. 0.1 µg/ml Ethidium Bromide should be included in the gel for visualization.

2. Following electrophoresis visualize the DNA samples in the gel with UV light and photograph.

2. Minimize the exposure of the gel to UV light as this may cause excessive nicking of the nucleic acid.




Figure 1. Diagrammatic representation of a capillary blotting apparatus

Protocol

3. Process the gel for blotting, rinsing the gel in distilled water between each step.

3a. Depurination

Place in 0.125 M HCl so that the gel is completely covered in the solution. Agitate gently for 10 minutes. During this time the bromophenol blue dye present in the samples will change colour.

3b. Denaturation

Submerge the gel in sufficient denaturation buffer. Incubate for 30 minutes with gentle agitation. During this time the bromophenol blue dye will return to its original colour.

3c. Neutralization

Place the gel in sufficient neutralization buffer to submerge the gel. Incubate for 30 minutes with gentle agitation.

4. Set up the capillary blot as described on page 14.

Notes

- 3a. Depurination is not required for DNA fragments <10 Kb in size. Do not over depurinate, 10 minutes (or until the bromophenol blue turns yellow) is sufficient for most samples.

9.3. Southern blotting – gel treatments (Alkali transfer protocol)

Protocol

1. Separate the DNA samples on a suitable neutral agarose gel.

Notes

1. Efficient separation of a range of DNA fragments may be achieved by varying the type and concentration of the agarose in the gel. Ensure the optimum DNA concentration for detection is loaded into each track. 0.1 µg/ml Ethidium bromide should be included in the gel for visualization.

2. Following electrophoresis visualize the DNA samples in the gel with UV light and photograph

2. Minimize the exposure of the gel to UV light as this may cause excessive nicking of the nucleic acid.

3. Process the gel for blotting, rinsing the gel in distilled water between each step.

3a. Depurination

Place in 0.125 M HCl so that the gel is completely covered in the solution. Agitate gently for 10 minutes. During this time the Bromophenol blue dye present in the samples will change colour.

- 3a. Depurination is not required for DNA fragments <10 Kb in size. Do not over depurinate, 10 minutes (or until the Bromophenol blue turns yellow) is sufficient for most samples.

Protocol Notes

Notes

3b. Denaturation

Submerge the gel in sufficient denaturation buffer. Incubate for 30 minutes with gentle agitation. During this time the Bromophenol blue dye will return to its original colour.

4. Set up the capillary blot as described on page 14, using Denaturation solution as the transfer buffer.

9.4. Northern blotting – gel preparation and treatment

RNA is separated under denaturing conditions, the principle systems currently in use are the Glyoxal/Dimethylsulphoxide and the Formaldehyde/Formamide procedures. This booklet restricts itself to the latter. Successful Northern analysis (7, 8) depends on the quality of the reagents used as well as having pure undegraded RNA samples.

Avoid any contamination with RNases, use sterile disposable plastics wherever possible. Glassware may be treated by baking at 180°C overnight or incubating in 0.2% (v/v) Diethylpyrocarbonate (DEPC) followed by autoclaving or baking. Some plastics are also suitable for DEPC treatment.

Formaldehyde/Formamide protocol

Protocol

Notes

1. Prepare the MOPS/
1. The Agarose gel is ≈ 0.7 M with

Protocol

1. Continued.

Formaldehyde gel as follows:
Preheat 17.5 ml of
Formaldehyde and 30 ml
10 x MOPS buffer at 55°C.
Dissolve 3–4.5 g of Agarose
in 250 ml of nuclease free
water. Cool to 55°C. Add
the 10 x MOPS buffer and
Formaldehyde. Cast the gel
in an appropriate enclosure
and allow the gel to set.

Notes

1. Continued.

respect to Formaldehyde
and 1 x with respect to the
MOPS buffer. This formulation
can be scaled up or down as
appropriate for the size of gel
required.

SybrGreen™ or Ethidium
Bromide (0.01 µg/ml) may
be included in the gel for
visualization. RNA does not
stain as well with Ethidium
Bromide as the same amount
of DNA with Ethidium
bromide. Excessive amounts
of Ethidium Bromide will
also inhibit RNA transfer (9).
Other staining procedures
post electrophoresis for
example Ethidium Bromide
(10) or Acridine Orange (11), or
methods which stain the blot,
for example Methylene Blue
(12) may also be used.

2. Sample must be deproteinized.

Samples may be stored at
-15°C to -30°C for short
periods.

2. Prepare the RNA sample(s), using the table below:

RNA	Volume (µl)	final conc	V
Formaldehyde	5.5	2.2 M	
Formamide	15	50%	

Protocol

2. Continued.

10 x MOPS buffer	1.5	0.5 x
Water	8-V	
TOTAL	30	

Place the sample(s) at 55°C for 15 minutes to denature. After denaturation add 3 µl of 10 x nucleic acid dye loading buffer. Mix and load onto the Agarose gel.

3. Separate the RNA samples using 1 x MOPS buffer as the electrophoresis buffer.

Notes

2. Continued.

3. SybrGreen is recommended for visualization. When Ethidium Bromide is used for visualization the addition of Ethidium Bromide to the electrophoresis buffer (0.01 µg/ml) improves results. Nucleic acid loading buffer must be prepared using RNase free reagents/solutions .

4. Following electrophoresis, if appropriate, visualize the RNA within the gel with UV light and photograph.

4. The integrity of the RNA may be assessed by the absence of smearing and the fluorescent signal, the ratio of 28 S to 18 S RNA should be 2:1.

5. Place the gel in a suitable tray or dish and cover with distilled water. Incubate the gel with gentle agitation for 15 minutes.

Protocol

6. Discard the water and replace with sterile $10 \times \text{SSC}$.
Agitate for 15 minutes.
Repeat this step once more.
7. Set up the capillary blot as described on page 9 using the neutral transfer buffer.

Notes

7. $10 \times \text{SSC}$ or $20 \times \text{SSC}$ can be used as the transfer buffer.

9.5. Protocol for colony and plaque lifts (10)

Protocol

1. Plate out the cells or bacteriophage in the usual way. Incubate overnight at the required temperature.
2. Pre-cool the petri-dishes for at least 30 minutes at 4°C before taking the lift.

Notes

1. Do not allow the colonies to grow too large. A colony/plaque density of up to 200 per 83 mm plate is optimal for accurate selection of positive clones.
2. Pre-cooling prevents smearing of the colonies and separation of the top agar layer. Plates must be free of excess moisture.
3. The hydrophilic nature of nylon ensures accurate colony/plaque lifts. If desired the membrane may be pre-wet before use, for example on an unused agar plate or on a TE buffer saturated Hybond blotting pad. Excess liquid must be removed from the membrane before proceeding, this is achieved by placing the

Protocol

Notes

3. *Continued.*

disc on a dry sheet of Hybond blotting paper; once it has touched the agar surface.

4. Bend the membrane and place the resulting trough across the centre of the petri-dish. Release the trough and allow the membrane to sit on the surface. Mark the disc position on the plate at several positions using a pin to ensure correct orientation of the colonies/plaques in subsequent manipulations.

4. This procedure will prevent air being trapped under the membrane. Do not force the membrane down. As it unrolls, the membrane disc will flatten. Do not attempt to move the membrane disc

5. After 30–60 seconds remove the membrane from the petri-dish in one continuous movement using blunt ended forceps. Place colony/plaque side uppermost on a sheet of Hybond blotting paper.

5. Extending the time the membrane remains on the surface of the agar will cause diffusion of the colonies/plaques. Replicate filters can be prepared by placing a fresh membrane disc on top of this template membrane. Press the membrane firmly together using a replica planting tool, avoid any lateral movement. Mark the replica membrane. Replica filters should then be incubated on fresh agar plates under appropriate conditions until colonies of 0.5–1 mm

Protocol

6. The DNA must be liberated from the bacteria or bacteriophage, denatured and then fixed to the membrane following a neutralization step. This is achieved by placing the membrane discs colony/plaque uppermost on a series of solution saturated Hybond blotting paper pads:

- Denaturation step, denaturation buffer for 2–5 minutes.
- Neutralization step, neutralization buffer for 3 minutes. Repeat this step once more.

7. Finally, vigorously wash the membrane disc in 2 x SSC to remove the proteinous debris.

Notes

5. *Continued.*
diameter are obtained.
6. An initial (optional) lysis step, 10% (w/v) SDS for 1–3 minutes may be included in the colony lift procedure. The Hybond blotting paper pads should be moist, though not too wet as this will cause diffusion of the colonies/plaques. Timings should be optimized, prolonged incubations will cause diffusion of the target DNA making accurate selection of positive clones difficult. Avoid fluid reaching the upper surface of the membrane. When transferring membrane, remove as much fluid as possible from the underside of the membrane. This may be achieved by transferring briefly to dry Hybond blotting paper between treatments.
7. Adequate removal of cell debris from colony lifts is essential

Protocol	Notes
8. Transfer the disc, DNA side uppermost, to a pad of dry Hybond blotting paper, air dry.	
9. Fix the DNA to the membrane by baking for 2 hours at 80°C or by using an optimized UV crosslinking procedure.	9. Details of an optimization procedure using a UV transilluminator are given on page 39. The UV crosslinker has a pre-set UV exposure (70 000 micro-joules/cm ²) which is suitable for Hybond-XL.
10. Membranes may be used immediately or stored, once dry.	10. Membranes should be stored wrapped in Saran Wrap desiccated at room temperature under vacuum.

9.6. Protocol for dot blotting (manual)

The following is a general protocol for dot blotting target nucleic acids. A number of commercially available devices are also available, for example the PR648 Slot blot manifold available from GE Healthcare. These provide for a more consistent and even application of the sample than the manual procedure described below. This parameter is particularly important in those experiments requiring quantification.

Protocol	Notes
1. Cut the membrane to an appropriate size.	1. The membrane should be cut with clean sharp scissors.
2. Using a pencil, mark the membrane lightly with a grid	

Protocol

2. *Continued.*

or dots to guide subsequent sample application. There should be a minimum distance of 1 cm between samples applied in a volume 5 μ l or less.

3. Pre-wetting the membrane is not required.

4. Dilute the samples in an appropriate buffer to the required concentration. TE buffer or 2 x SSC may be used for DNA samples. RNA samples should be prepared using the information in the table on page 21. A sample size of 1–2 μ l is ideal for manual dot blotting.

5. Nucleic acid samples must be denatured to provide a suitable single-stranded target molecule for subsequent hybridizations. Denature the samples by

Notes

3. Membranes may be pre-wet if desired, see critical parameters page 8.

4. Carrier substance may be included in the diluent buffer to improve retention of very small amounts of target on the membrane. These include:

- sonicated herring sperm DNA for DNA samples
- tRNA for use with RNA samples.

Larger sample volumes of 50–200 μ l are common for commercial apparatus. This ensures an even application of the sample over the whole dot or slot.

5. Samples may also be prepared in 0.2 M NaOH incubated at 37°C for 15 minutes and dotted directly onto the membrane, in a denatured condition. Alternatively RNA samples may

Protocol

5. *Continued.*
heating in a boiling water bath for 5 minutes. Chill rapidly on ice, then centrifuge briefly to collect sample at the bottom of the tube.

6. Carefully apply the sample to the membrane, avoiding touching the membrane with the pipette tip. Leave the membrane to air dry.

7. Fix the nucleic acid sample to the membrane by UV crosslinking or baking at 80°C for 2 hours.

8. Blots may be used immediately or stored wrapped in SaranWrap desiccated at room temperature under vacuum.

Notes

5. *Continued.*
be preheated to 55°C for 15 minutes, see page 21
6. If the sample volume is greater than 2 μ l, then apply in successive 2 μ l aliquots to the same position on the membrane, allow the aliquot to dry between each application. This will prevent the sample from spreading.
7. Details of an optimization procedure are given on page 28. The GE Healthcare UV crosslinker has a pre-set UV exposure (70 000 micro-joules/cm²) which is optimum for Hybond-XL.

9.7. Hybridization protocols - hybridizations in bags and boxes

Protocol

1. Prepare the hybridization buffer, for example

Denhardt's Buffer

5 x SSC

5 x Denhardt's solution

0.5% (w/v) SDS

Modified Church and Gilbert Buffer (13)

0.5 M Phosphate buffer, pH 7.2

7% (w/v) SDS

10 mM EDTA(13)

Ensure the SDS is in solution before use. Gentle heating may be necessary.

Combine all the components, make up to the required volume.

Notes

1. This Denhardt's based buffer is that used in the quality control of all Hybond nylon membranes. A reduced concentration of SDS may lead to elevated backgrounds following hybridization. The Denhardt's hybridization buffer may be stored at -15°C to -30°C if required.

This modification of the Church and Gilbert buffer is routinely used in GE Healthcare's Laboratories. It has been shown to be suitable for Southern, Northern, dot blots and library screening applications. The hybridization buffer may be stored at room temperature. Ensure the SDS is fully dissolved before use. This may be achieved with gentle heating.

2. Prepare the radiolabelled probe using the appropriate procedure.
2. For radioactive applications use a probe concentration of $0.5\text{--}2 \times 10^6$ incorporated counts per ml of hybridization buffer for single copy gene

Protocol

Notes

2. *Continued.*

detection (i.e. high sensitivity application) or 0.125–0.5 x 10⁶ incorporated counts per ml of hybridization buffer for high target work (for example colonies/ plaques, PCR products etc). Probe purification, to remove unincorporated radioactive nucleotides, is strongly recommended

3. Preheat the required volume of hybridization buffer to the appropriate temperature.

4. Pre-wet the blot in a suitable buffer for example 5 x SSC or 0.5 M Phosphate buffer. Place the blot(s) in the hybridization buffer. 125 µl of hybridization buffer per cm² is a suitable volume. Prehybridize for at least 30 minutes with constant agitation, at the desired hybridization temperature (see Note 7).

5. When using labelled double stranded probes, pipette the required amount into a clean

4. Details of the pre-wetting procedures are given on page 6, Critical Parameters. Hybridization may be carried out in bags, or boxes, provided there is sufficient buffer for the container. Adequate circulation of the buffer is essential. When hybridizing several blots together, the blots should move freely within the buffer.

Protocol

5. *Continued.*

microcentrifuge tube. If the volume is less than 20 μ l, make up to this volume with water or TE buffer. Denature the probe by boiling for 5 minutes and snap cool on ice. Briefly centrifuge to draw the contents to the bottom of the tube.

6. Add the probe to the pre-hybridization buffer.

6. Avoid placing the probe directly on the blots, as this will cause excessive background.

7. Hybridize overnight with gentle agitation at the required temperature.

7. Hybridization temperature may vary with the probe. Lower temperatures achieve lower stringency. The temperature of hybridization used will depend on the degree of homology between the probe and the target. 65–68°C is suitable for most long probes (>100 bases). With short/oligo probes (<50 bases) hybridization temperature is usually defined as $T_m - 5^\circ\text{C}$:
 T_m (melting temperature) = $(4 \times \text{number of G+C bases}) + (2 \times \text{number of A+T bases}) / (14)$
Hybridization time can also

Notes

Protocol

Notes

7. *Continued.*

vary. Short hybridization times may be suitable for high target applications.

8. Prepare the stringency wash solutions. The wash solution should be used in excess, at least 1–5 ml/cm² of membrane.

For example:

Low stringency wash:

2 x SSC, 0.1% (w/v) SDS

Medium stringency wash:

1 x SSC, 0.1% (w/v) SDS

High stringency wash:

0.1 x SSC, 0.1% (w/v) SDS

8. Stringency washes will depend on the nature of the probe and target to be hybridized. Salt concentration and temperature should be taken into consideration. The lower the salt concentration, the greater the stringency. The higher the washing temperature, the greater the stringency.

Most commonly, stringency washes proceed from “high salt”/“low temperature”, for example 5 x SSC, 0.1% (w/v) SDS at room temperature, to “low salt/high temperature”, for example 0.1 x SSC, 0.1% (w/v) SDS at 65°C (nominal hybridization temperature).

9. After the hybridization, wash the blots by incubating twice, 5 minutes each, in 2 x SSC, 0.1% (w/v) SDS, followed by 1 x SSC, 0.1% (w/v) SDS for 15 minutes, and finally

9. Some procedures include room temperature washes under low stringency conditions. Do **not** allow the SDS to come out of solution during these washes, significant levels

Protocol	Notes
9. Continued. 0.1 x SSC, 0.1% (w/v) SDS for 2 x 10 minutes, at the hybridization temperature.	9. Continued. of background may result. Adequate circulation of the stringency buffer is essential when washing. Washing in boxes is advised.
10. Remove the blot from the last stringency wash, drain, wrap in SaranWrap and expose to X-ray film, for example Hyperfilm MP. Keep the blot moist if it is to be reprobbed. If reprobing is desired, it may be more suitable to seal the blot in a plastic bag.	10. The use of SaranWrap with ³⁵ S labelled probes will significantly increase exposure times. In this case the blot should be air dried before autoradiography, if reprobing is not required.

9.8. Hybridization in tubes

There are numerous commercially available rotisserie devices suitable for use as hybridization ovens for example GE Healthcare hybridization oven/shaker RPN2510E/2511E. These can accommodate 2 to 7 hybridization tubes. The major advantage of this approach to hybridization is the use of low volumes of hybridization buffer, and therefore minimal probe volumes. This is achieved because fluid is able to move continually over the membrane.

Protocol	Notes
1. Prepare the hybridization buffer, for example Denhardt's Buffer	1. There are a wide variety of hybridization buffers used by researchers. This Denhardt's

Protocol

1. *Continued.*

- 5 x SSC
- 5 x Denhardt's solution
- 0.5% (w/v) SDS

Modified Church and Gilbert

Buffer (13)

- 0.5 M Phosphate buffer,
- pH 7.2
- 7% (w/v) SDS
- 10 mM EDTA(13)

Ensure the SDS is in solution before use. Gentle heating may be necessary.

Combine all the components, make up to the required volume.

2. Prepare the radiolabelled probe using the appropriate procedure.

3. Preheat the required volume of hybridization buffer to the appropriate temperature.

Notes

1. *Continued.*

based buffer is that used in the quality control of all Hybond nylon membranes. A reduced concentration of SDS may lead to elevated backgrounds following hybridization. The Denhardt's hybridization buffer may be stored at -15°C to -30°C if required. This modification of the Church and Gilbert buffer is routinely used in GE Healthcare's Laboratories. It has been shown to be suitable for Southern, Northern, dot blots and library screening applications. The hybridization buffer may be stored at room temperature. Ensure the SDS is fully dissolved before use. This may be achieved with gentle heating.

3. Hybond-XL has been designed for use with very low volumes of hybridization buffer (30–70 $\mu\text{l}/\text{cm}^2$). High backgrounds will result if

Protocol

Notes

sub optimum volumes are used for the membrane and hybridization conditions.

4. Pre-wet the blot in a suitable dish, first in water then in an appropriate buffer. Ensure that the nucleic acid side is uppermost. Roll the blot along its length in such a way as to minimize overlap in the tube. Place inside the hybridization tube.

4. If there is significant overlap of the blot use of a nylon mesh should be considered. The mesh achieves separation of the blot layers allowing better probe access to these areas. It is strongly advised that hybridization volume should be increased to 70–125 $\mu\text{l}/\text{cm}^2$.

The nylon mesh should be at least 0.5 cm larger than the blot. Place the mesh in the pre-wetting solution before the blot, in subsequent manipulations treat as “one”.

The nylon mesh may be reused after washing in 10% (w/v) SDS and extensive rinsing in distilled water.

5. Add a small volume of appropriate buffer to the hybridization tube, cap the tube. Unroll the blot by rotating the tube in the opposite direction to the “rolled” blot.

5. It is important not to allow air to become trapped between the inner surface of the tube and the membrane. This can cause areas of no signal or background following hybridization.

Protocol

6. Drain the tube of excess liquid and replace with the appropriate volume of hybridization buffer.

7. Prehybridize for 30 minutes at the appropriate temperature. Ensure that the tube is placed in the correct orientation within the oven to avoid rolling up of the blot.

8. When using labelled double stranded probes, pipette the required amount into a clean microcentrifuge tube. If the volume is less than 20 μ l, make up to this volume with water or TE buffer. Denature the probe by boiling for 5 minutes and snap cool on ice. Briefly centrifuge to draw the contents to the bottom of the tube.

Notes

6. Hybond-XL has been designed for use with low volumes of hybridization buffer (30–70 μ l/cm²). High backgrounds will result if sub optimum volumes are used for the membrane and hybridization conditions.

8. For radioactive applications use a probe concentration of $0.5\text{--}2 \times 10^6$ incorporated counts per ml of hybridization buffer for single copy gene detection (i.e. high sensitivity application) or $0.125\text{--}0.5 \times 10^6$ incorporated counts per ml of hybridization buffer for high target work (for example colonies/plaques, PCR products etc). Probe purification, to remove unincorporated radioactive nucleotides, is strongly recommended.

Protocol

9. Add the probe to the pre-hybridization buffer.

10. Hybridize overnight at the required hybridization temperature.

Notes

9. Avoid placing the probe directly on the blot. Probe may be added to the hybridization while the tube is in a vertical position. If necessary probe may be mixed with a portion of the hybridization buffer and added to the tube in a larger volume.

10. Hybridization temperatures may vary with the probe. Lower temperatures achieve lower stringency. The temperature of hybridization used will depend on the degree of homology between the probe and the target. 65–68°C is suitable for most long probes (>100 bases). With short/oligo probes (<50 bases) hybridization temperature is usually defined as $T_m - 5^\circ\text{C}$:
 T_m (melting temperature) = $(4 \times \text{number of G+C bases}) + (2 \times \text{number of A+T bases})$
(14) Hybridization time can also vary. Short hybridization times may be suitable for high target applications .

Protocol

- 11.** Prepare the stringency wash solutions. The wash solution should be used in excess. Use a volume that occupies 30–50% of the tube. For example:
Low stringency wash:
2 x SSC, 0.1% (w/v) SDS
Medium stringency wash:
1 x SSC, 0.1% (w/v) SDS

High stringency wash:
0.1 x SSC, 0.1% (w/v) SDS

- 12.** After the hybridization wash the blot as follows:
a) rinse briefly in 2 x SSC, 0.1% (w/v) SDS
b) twice, 5 minutes each in 2 x SSC, 0.1% (w/v) SDS
c) twice, 10 minutes each in 1 x SSC, 0.1% (w/v) SDS
d) four times, 5 minutes each in 0.1 x SSC, 0.1% (w/v) SDS

- 13.** Remove the blot from the last stringency wash, drain and wrap in SaranWrap and expose to X-ray film, for example Hyperfilm MP. Keep the blot moist if it is to be reprobed.

Notes

- 11.** Washing in boxes is much more effective and is recommended if feasible. The inefficiencies of washing in tubes may be overcome by increasing the number of stringency washes while maintaining the same total wash time.

- 13.** The use of SaranWrap with ³⁵S labelled probes will significantly increase exposure times. In this case the blot should be air dried before autoradiography, if reprobing is not required.

9.9. Stripping protocols

For the successful removal of probes, membranes must never be allowed to dry out after hybridization and washing, since this process has the effect of “fixing” the hybrid.

The two most commonly used stripping procedures are described below. Some users have experienced variability in these methods. In these cases a combination of both of these procedures has proved suitable.

Hot SDS procedure

Protocol

Notes

1. Place the moist membrane in an appropriate sized tray.
2. Prepare a boiling solution of 0.1% (w/v) SDS, pour the solution onto the blot and allow to cool (15–30 minutes)
3. Rinse the blot briefly in 2 x SSC

2. This step may be repeated if the probe is particularly difficult to remove.

4. Check the removal of the probe using the appropriate procedure for the labelling and detection system used.
5. Hybridize overnight using the appropriate conditions.

Alkali procedure

Protocol

Notes

1. Place the moist membrane in an appropriate sized box.
2. Pre-heat 0.2 M NaOH to 42°C.
2. Higher concentrations of NaOH alone and/or in

Protocol

Notes

2. *Continued.*
combination with NaCl may also be used.
3. Incubation times may also vary. Optimization of the procedure for the probe and target in use is strongly advised.
3. Add the alkali solution to the blot and incubate at 42°C for 10 minutes with constant agitation.
4. Repeat step 3 with fresh alkali solution.

10. Additional information

10.1. Determination of the optimum UV crosslinking conditions using a UV transilluminator

Protocol

1. Produce five or six identical control blots, for example Lambda *Hind* III on the membrane of choice.
2. Protect the surface of the membrane by covering the transilluminator with SaranWrap. Expose each blot DNA side down on the transilluminator for a different length of time, for example 30 seconds to 5 minutes.

Notes

1. The type of blot should reflect the technique for which the calibration is being used.
2. The length of exposure required for optimum fixation will vary depending on the wavelength of the UV bulb and its age. The energy emitted from a UV bulb is reduced with use. Regular recalibration is advised if the apparatus is extensively used. This inconvenience may be overcome with the use of UV crosslinkers which are able to compensate for this effect, when used on the constant energy setting. A UV crosslinker with pre-set or manual energy and time settings is available from GE Healthcare. UVC 500 UV Crosslinker 115 VAC/230 VAC 80-6222-31/80-6222-50.

Protocol

3. Hybridize all the blots together with a suitably labelled probe.
4. Following autoradiography, the optimum UV exposure time will be indicated by selecting the blot showing the strongest signal.

10.2. Determination of the optimum alkali fixation

Protocol

1. Produce five or six identical control blots, for example Lambda *Hind* III on the membrane of choice.
2. Prepare a fresh solution of 0.4 M NaOH. Prepare and alkali soaked pad of Hybond blotting paper in an appropriate container.
3. Expose each blot DNA side uppermost on the alkali soaked Hybond blotting pad for a different length of time, for example 60 seconds to 10 minutes.
4. Hybridize all the blots together with a suitably labelled probe.

Notes

1. The type of blot should reflect the technique for which the calibration is being used.
2. Fresh alkali solution must be used. The soaked Hybond blotting paper must be drained of excess liquid.
3. For storage, blots must be neutralized. Following air drying rinse briefly in 2 x SSC, air dry again and store desiccated under vacuum.

Protocol

5. Following autoradiography, the optimum alkali incubation time will be indicated by selecting the blot showing the strongest signal.

11. Related products

Nucleotides

Application areas	Compound	Specific Activity Ci/mmol	Redivue code	Standard code
DNA labelling	[α - ³² P]dATP	~3000	AA0004	PB10204
	[α - ³² P]dCTP	~3000	AA0005	PB10205
	[α - ³² P]dGTP	~3000	AA0006	PB10206
	[α - ³² P]dTTP	~3000	AA0007	PB10207
RNA labelling	[α - ³² P]ATP	~3000		PB10200
	[α - ³² P]CTP	~3000		PB10202
	[α - ³² P]GTP	~3000		PB10201
	[α - ³² P]UTP	~3000		PB10203
End labelling of DNA	[γ - ³² P]ATP	~3000	AA0068	PB10168
		>5000	AA0018	PB10218
	[γ - ³² P]GTP	>5000		PB10244
	[α - ³² P]ddATP	~3000		PB10233

Only some of the wide variety of nucleotides codes available are shown, please refer to the catalogue for a comprehensive listing.

Labelling Kits

Megaprime DNA labeling system

Suitable for use with labeled dCTP

Suitable for use with any labeled dNTP

Ready to go labeling beads

Ambient stable single dose reaction

spheres for use with dCTP

RPN1606/1607

RPN1604/1605

27-9240-01

Rediprime™ II DNA labelling system	
premixed individually dispensed freeze dried reactions for use with dCTP only	RPN1633/1634
End labelling kits	
3' end labelling	RPN1509

Please contact your local customer service office for ordering details associated with GE Healthcare products.

Companion products

ProbeQuant™ G-50 Micro Column	27-5335-01
SeeDNA™, co-precipitant	RPN5200

A range of branded a autoradiography products

- Hyperfilm
- Hypercassettes
- Hyperscreens
- Kodak™ BioMax™
- Kodak X-Omat

Please refer to the catalogue for a comprehensive listing.

Please contact your local customer service office for ordering details associated with GE Healthcare products.

12. References

1. Southern, E.M., *J. Mol. Biol.* **98**, 503–517 (1975).
2. Lichenstein, A.V., Moiseer, V.L. and Zaboikin, M.M., *Anal. Biochem.* **191**, 187 (1990).
3. Smith, G.E., and Summer, M.D., *Anal. Biochem.* **109**, 123 (1980).
4. Chomczynski, P., *Anal. Biochem.* **201**, 134 (1992).
5. Olszewska, E and Jones, K., *Trends Genet.* **4**, 92 (1988).
6. Bittner, M., Kupferer, P and Morris, C.F. *Anal. Biochem.* **102**, 459 (1980).
7. Alwine, J.C., Kemp, D.J. and Stark, G.R., *Proc. Natl. Acad. Sci. USA* **74**, 5350–5354 (1977).
8. Alwine, J.C. *et al.*, *Methods Enzymol.* **68**, 220 (1979).
9. Thomas, P.S., *Methods Enzymol.* **100**, 255 (1983).
10. Sambrook, J. *et al.*, in *Molecular Cloning*, a laboratory manual, second edition, Cold Spring Harbor lab press, (1989).
11. Wilkinson, M., Doskow, J. and Lindey, S. *Nuc. Acid Res.* **19**, 679 (1991).
12. McMaster, G.M. and Carmichael, G.G., *Proc. Natl. Acad. Sci. USA* **74**, 4835 (1977).
13. Church, G.M. and Gilbert, W., *Proc. Natl. Acad. Sci. USA* **81**, 1991 (1984).
14. Wallace, R.B. and Miyada, C.G., *Meth. Enzymol.*, 152, p.438. 1990

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30 751 84

Uppsala

Sweden

GE Healthcare Europe GmbH

Münzinger Strasse 5 D-79111

Freiburg

Germany

GE Healthcare UK Limited

Amersham Place

Little Chalfont

Buckinghamshire

HP7 9NA

UK

GE Healthcare Bio-Sciences Corp

800 Centennial Avenue

P.O. Box 1327

Piscataway

NJ 08855-1327

USA

GE Healthcare Bio-Sciences KK

Sanken Bldg. 3-25-1

Hyakunincho Shinjuku-ku

Tokyo 169-0073

Japan

GE Healthcare

regional office contact numbers:

Asia Pacific

Tel: +85 65 6 275 1830

Fax: +85 65 6 275 1829

Australasia

Tel: +61 2 8820 8299

Fax: +61 2 8820 8200

Austria

Tel: 01 757606 1613

Fax: 01 757606 1614

Belgium

Tel: 0800 73 890

Fax: 02 416 82 06

Canada

Tel: 1 800 463 5800

Fax: 1 800 567 1008

Central & East Europe

Tel: +43 1 972720

Fax: +43 1 97272 2750

Denmark

Tel: +45 70 25 24 50

Fax: +45 16 24 24

Eire

Tel: 1 800 709992

Fax: +44 1494 542010

Finland & Baltics

Tel: +358 9 512 39 40

Fax: +358 9 512 39 439

France

Tel: 01 6935 6700

Fax: 01 6941 9877

Germany

Tel: 0800 9080 711

Fax: 0800 9080 712

Greater China

Tel: +852 2100 6300

Fax: +852 2100 6338

Italy

Tel: 02 26001 320

Fax: 02 26001 399

Japan

Tel: 81 3 5331 9336

Fax: 81 3 5331 9370

Korea

Tel: 82 2 6201 3700

Fax: 82 2 6201 3803

Latin America

Tel: +55 11 3933 7300

Fax: +55 11 3933 7304

Middle East & Africa

Tel: +30 210 9600 687

Fax: +30 210 9600 693

Netherlands

Tel: 0800 82 82 82 1

Fax: 0800 82 82 82 4

Norway

Tel: +47 815 65 777

Fax: +47 815 65 666

Portugal

Tel: 21 417 7035

Fax: 21 417 3184

Russia, C.I.S. & N.I.S

Tel: +7 495 956 5177

Fax: +7 495 956 5176

Spain

Tel: 902 11 72 65

Fax: 935 94 49 65

Sweden

Tel: 018 612 1900

Fax: 018 612 1910

Switzerland

Tel: 0848 8028 30

Fax: 0848 8028 11

UK

Tel: 0800 515 313

Fax: 0800 616 927

USA

Tel: +1 800 526 3593

Fax: +1 877 295 8102

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited

Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA

UK



imagination at work