
DOC022.97.93074

TitraLab® AT series workstations

09/2014, Edition 1

Basic User Manual
Manuel d'utilisation de base
Manual básico del usuario
Manual Básico do Usuário
基本用户手册
基本取扱説明書
기본 사용 설명서

English 3

Français 23

Español 45

Português 67

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Table of contents

Specifications on page 3	Startup on page 16
General information on page 3	Standard operations on page 17
Installation on page 8	Maintenance on page 19
Keypad on page 15	Troubleshooting on page 22

Additional information

Additional information is available on the manufacturer's website.

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	22 x 40 x 36 cm (8.7 x 15.7 x 14.2 in.)
Weight	4 kg (8.8 lb)
Power requirements	100–240 VAC, 50/60 Hz
Altitude	2,000 m (6,562 ft) maximum
Operating temperature	15 to 35 °C (59 to 95 °F)
Relative humidity	20 to 80%, non-condensing
Storage temperature	–5 to 40 °C (23 to 104 °F)
Installation category	II
Pollution degree	2
Certifications	Safety IEC/EN 61010-1; EMC IEC/EN 61326-1
Warranty	1 year (EU: 2 years)

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol indicates the presence of devices sensitive to Electro-static Discharge (ESD) and indicates that care must be taken to prevent damage with the equipment.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user.

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested


and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Product overview

The instrument operates with digital and analog sensors. Measurement applications are installed on the instrument to automate the measurement process. Instructions show on the display when user intervention is required. Refer to [Figure 1](#) for product features.

Figure 1 Product overview

1 Keypad	6 Beaker	11 Sensor holder
2 Display	7 Syringe protection cover	12 Pump 2 input/output
3 Sensor storage tubes	8 Syringe input/output	13 Pump 1 input/output
4 USB port	9 Tube clips	14 Pump access cover
5 Tube holder	10 Syringe	

Note: Depending on the model, there will be 1 or 2 syringes and syringe input/output ports, and 0, 1 or 2 pumps. Refer to [Table 1](#).


Table 1 Instrument configurations

Model	Syringes	Pumps
AT1102.98	1	0
AT1112.98	1	1
AT1122.98	1	2
AT1222.98	2	2

Instrument connections

Figure 2 shows the connections on the rear panel of the instrument. Use the USB port on the side of the instrument for the USB applications key supplied with the instrument. Use the USB port on the rear of the instrument to connect to a printer, mouse, keyboard or a USB hub.

Figure 2 Instrument connections


1 24 V external power supply port	4 External pump port	7 USB port
2 Sensor 1 port	5 External propeller port	8 Ethernet port
3 Sensor 2 port	6 Serial port	

Product components

Make sure that all components have been received. Refer to the packing list in the box. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.


Figure 3 Contents of the instrument box

1 Instrument	3 Sensor storage tubes (3x)	5 Power cord
2 Tube holder ¹	4 Power supply	

¹ 1 for each syringe position available on the instrument

Figure 4 Contents of the application box

1 Beakers (10 x 50 mL and 10 x 150 mL)	7 Syringe holding ring ⁴
2 Tube with anti-diffusion tip ¹	8 Syringe (refer to Table 1 on page 6 for quantity)
3 Conical adapters (2x)	9 USB applications key
4 Magnetic stir bars (10x)	10 Glass bottles ⁵
5 Legacy sensor adapter ²	11 Bottle caps (2 x GL45 and 1 x GL25)
6 Sensor ³	12 Empty desiccant cartridges (3x)

- ¹ If necessary for the application
- ² Not in all application kits
- ³ Type and quantity depends on application
- ⁴ 1 for each syringe
- ⁵ Not in all application kits

Installation

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

The instrument is available in different configurations (refer to [Table 1](#) on page 6). This manual supplies instructions for the installation of an instrument with one syringe and one pump. Adjust the installation procedure as applicable to accommodate the number of syringes and pumps in the instrument.

Installation guidelines

- This instrument is for indoor use only.
- The power supply connector on the rear panel must be easily accessible so the power can be disconnected quickly in case of emergency.
- Keep the instrument away from temperature extremes, including heaters, direct sunlight and other heat sources.
- Put the instrument on a stable and level surface in a well ventilated place.

- Make sure that there is at least 15 cm (6 in.) of space on all sides of the instrument to prevent electrical parts from overheating.
- Do not operate or keep the instrument in dusty, damp or wet locations.
- Always keep the surface of the instrument and all accessories dry and clean.

Connect to AC power

⚠ DANGER	
	Electrocution hazard. If this equipment is used outdoors or in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used for connecting the equipment to its main power source.
⚠ CAUTION	
	Electrical shock and fire hazards. Make sure that the supplied cord and non-locking plug meet the applicable country code requirements.
⚠ WARNING	
	Fire hazard. Use only the power supply that is specified for this instrument.

1. Connect the power cord to the power supply.
2. Connect the power supply to the instrument. Refer to [Figure 2](#) on page 6.
3. Connect the power cord to an electrical outlet.

Install the syringe


Before syringe installation, set the instrument power to on. Push the power button on the front of the instrument. Make sure that the startup sequence shows on the display. The syringe holder lowers to its operating position.







Note: Ignore any warning messages related to missing applications that show on the display.

The sensor holder has two positions: one over the magnetic stirrer and the second at 180° to the right. Move the sensor holder away from the instrument to the second position.

Refer to the illustrated steps that follow.

To install a second syringe, do steps 5 through 7 again.


Install the sensor storage tubes

Put the three sensor storage tubes into the holder that is on the side of the instrument. Refer to [Figure 1](#) on page 5. Keep the sensor in a storage tube when not in use.


Install the stir bar and the beaker

Add the stir bar to the beaker, and then attach the beaker to the sensor holder. Refer to the illustrated steps that follow.

Prepare the tubes

Remove any bends in the end of the tubes. Refer to the illustrated steps that follow.


Connect the tubes

Arrow symbols identify the inlet and outlet ports for the syringe and the pump connections. The “up” arrow is the outlet port. The “down” arrow is the inlet port. Turn the tube connectors on the inlet and outlet ports of the syringe and pump until they click.

The syringe outlet tube has a blue ring on it. If anti-diffusion tips are necessary, remove the pre-installed outlet tube from the syringe and install the tube from the application kit with the pre-installed anti-diffusion tip.

Push the outlet tubes into the tube holder slots so that they are correctly attached.

Refer to the illustrated steps that follow.


Install the sensor

Install the legacy adapter

Go to [Connect the sensor](#) on page 13 if no legacy adapter is included in the application kit.

1. Connect the measuring, reference and temperature sensors to the legacy adapter.
2. Connect the legacy adapter cable to a sensor socket on the rear panel of the instrument. Refer to [Figure 5](#).

Figure 5 Install the legacy adapter


1 Temperature sensor

2 Reference sensor

3 Measuring sensor

Connect the sensor

Use a conical adapter to hold the sensor tightly in the sensor holder. Connect the sensor to an available sensor port on the rear of the instrument. After the sensor is connected, make sure that the sensor icon shows in the banner at the top of the display. Refer to the illustrated steps that follow.

Install the titrant and the reagent

⚠ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.


Loosen the tube connector on the bottle cap.

Fill a desiccant cartridge with an applicable desiccant. Put the desiccant cartridge into the adapter on the bottle cap.


Push the inlet tube through the connector. Make sure that the end of the tube is at the bottom of the bottle. Tighten the connector on the bottle cap.

Refer to the "Application Note" on the USB applications key to identify the correct pump to connect to the reagent bottle. Refer to the illustrated steps that follow.


1

2


3

Tidy the work area

Attach the tubes to the instrument with the clips on the electrovalve and the sensor holder. Refer to the illustrated steps that follow.


User interface and navigation

Keypad

Figure 6 shows the keypad and gives the key functions.

Figure 6 Keypad

1 Power	4 Navigation keys
2 Printer	5 Home
3 Selection keys	

Key	Description
Power	Sets the instrument power to on or off. Push the key for 2 seconds to set the power to off.
Printer	Sends data to an attached printer. The printer key only operates if a printer is connected to the instrument.
Selection keys (contextual)	Shows the measurement options, selection and confirmation options. Use this key to exit the current menu display or open sub-menus. Available options show on the display above each key.
Navigation keys	Scrolls through menus and data, enter numbers and letters, enter checkbox settings and set options for the syringe and the pump.
Home	Goes to the main menu.

Startup

⚠ CAUTION	
	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

⚠ CAUTION	
Personal injury hazard. Never use the instrument without the syringe cover installed.	

Configure the instrument

1. From the main menu, select  **Settings**.
2. Select an option, then push **Select**.

Option	Description
Applications	Exports, changes, removes and makes copies of application data. Make sure that the duplication function does not make more than five applications for each syringe installed.
Operators	Adds, changes and removes operators.
Date + Time	Sets the instrument date and time.
Brightness	Sets the brightness of the display.
Sounds	Sets the sound options.
Language	Sets the language.
Info	Shows information about the instrument and the attached hardware.
Restore Defaults	Sets the instrument to the default configuration.
Options	Sets the application parameters view to expert mode, which lets the user customize the standard measurement applications. When the instrument is set to off, sets the syringe to empty into the titrant bottle. Changes the temperature display from Celsius to Fahrenheit. If a printer is connected, prints the measurement and derivative curves.

3. Push **Back**.

Install the applications

Use the supplied USB key to install the applications. The instrument can install a maximum of five applications for each syringe installed. For two syringes, the installed applications shown on the top line of the display refer to syringe one and the installed applications shown on the second line refer to syringe two. If any errors occur during installation, refer to [Troubleshooting](#) on page 22.

1. Push **Home** to go to the main menu.
2. Connect the USB key to the USB port on the side of the instrument. The applications on the USB key show on the display.
3. Push the arrow keys to highlight and select an application to install. Push the left or right arrow key to select it. Do this step again to select additional applications to install.
4. Push **Import** to install the selected applications.
5. Push **OK** to complete the installation. The installed applications show on the main menu.
Note: To install more applications, push **Home** to go to the main menu, then remove the USB key and reconnect it.

Prepare the instrument for measurement

1. From the main menu, select  **Purge**, then push **Start**. All attached devices are listed.
2. Select **All elements** to purge all the attached devices, or select one device to purge. Push **Select**. Air is removed from the device and filled with liquid from the bottle.
3. Push **OK** when the operation has completed.
4. Make sure that there are no air bubbles in the device. Do step 2 again if there are any air bubbles.
5. Select the next device to purge if individual devices are being selected.
6. Push **Exit** when all the tubes are filled with reagent and the device has no air bubbles.
Note: If a few small air bubbles can be seen on the inner wall and/or piston of the syringe, they can be left without effecting system performance.

Standard operations

▲ CAUTION	
	Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.
▲ CAUTION	
	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.
▲ CAUTION	
Personal injury hazard. Never use the instrument without the syringe protection cover in place.	
▲ CAUTION	
Chemical exposure hazard. Never remove the stir bar from the beaker before the end of a titration.	

Get a sample measurement

Use this option to get sample measurements with one of the installed applications.

1. From the main menu, select the measurement application, then push **Start**. Application information shows on the display.
2. Read the related "Application Note" from the USB applications key for more instructions.
3. If necessary, select an icon for more information or to change some data.
4. Fill a beaker with the recommended sample amount that shows on the display. If necessary, add sufficient deionized water to make sure the sensor is correctly installed in the sample.
5. Carefully put a magnetic stir bar into the beaker. Make sure there is no liquid spill.
6. Attach the beaker to the sensor holder.
7. Make sure that the icon at the bottom of the display  is highlighted. Do the instructions that show on the display adjacent to this icon. Refer to [Connect the sensor](#) on page 13 to make sure that the tubes and sensor are correctly aligned.
8. Push **Start** to start the measurement. Measurement data shows on the display.
9. If the default stirring speed needs to be adjusted, push the up and down arrow keys to increase or decrease the speed.
10. If the options are available during the measurement procedure, push **Skip** to ignore the current step or push **Stop** to abort the measurement.

11. When the measurement is complete, push the arrow keys to see the different measurement views.

12. Push **Next** for the options that follow:

Option	Description
Replicate sample	Use this option to start the same titration on the same sample. This is used to study the repeatability by successively analyzing several parts of the same sample. At the end of each measurement, a window shows the average value, the standard deviation and the relative standard deviation.
New sample	Use this option to start the same titration on a new sample. No standard deviation or relative standard deviation measurements will be done.

13. Push **Exit** to go back to the main menu.

Manage the data log

To select data to view, delete or export, specify data filters

1. From the main menu, select  **Data log**.
2. Select an option, then push **Select**.

Option	Description
View data log	Views measurement data. Select individual lines of data to view more content.
Delete data log	Removes measurement data from the system. Previews data selection before it is removed.
Export data log	Exports measurement data from the system to an external device. Preview data selection before it is exported. Make sure that an external device is connected to the instrument (e.g., a USB key, external hard drive, etc.).

3. Specify the data extraction parameters. Push the left and right arrow keys to make a selection. Push the up and down arrow keys to select an option.

Option	Description
Result type	Sets the type of result available.
Application	Sets the available applications.
Date	Sets the date range.
Operator	Sets the available operators.

4. Push **Preview** to preview the results.

Calibration

Calibrate the sensor

1. From the main menu, select  **Calibration**, then push **Electrode calibration**.
2. If more than one sensor is installed, push the arrow keys to highlight the sensor to calibrate, then push **Select**.
3. If more than one application includes calibration parameters for the sensor, push the arrow keys to choose the application, then push **Select**. Application information shows on the display.
4. If necessary, select an icon for more information or to change some data.
5. Do the instructions that show on the display, then push **Start** to start the calibration. Calibration data shows on the display.
6. If the default stirring speed needs to be adjusted, push the up and down arrow keys to increase or decrease the speed.
7. When the calibration is complete, push the arrow keys to see the different measurement views.

8. Push **Yes** to continue with the next calibration buffer solution.
9. When the calibration is complete, push **Yes** to accept the calibration or **No** to reject.
10. Push **Exit** to go back to the main menu.

Calibrate the titrant

1. From the main menu, select  **Calibration**, then push **Titrant calibration**. Related information shows on the display.
2. If necessary, select an icon for more information or to change some data.
3. Do the instructions that show on the display, then push **Start** to start the calibration. Calibration data shows on the display.
4. If the default stirring speed needs to be adjusted, push the up and down arrow keys to increase or decrease the speed.
5. When the calibration is complete, push the arrow keys to see the different measurement views.
6. Push **Continue** to continue with the calibration.
7. When the calibration is complete, push **Yes** to accept the calibration or **No** to reject.
8. Push **Exit** to go back to the main menu.

Purge

Use this procedure to remove air bubbles from the system. Refer to [Prepare the instrument for measurement](#) on page 17 for instructions.

Maintenance

⚠ CAUTION	
	Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

NOTICE
Do not disassemble the instrument for maintenance. If the internal components must be cleaned or repaired, contact the manufacturer.

Clean the instrument

NOTICE
Never use flammable or corrosive solvents to clean any part of the instrument. Use of these solvents can degrade the environmental protection of the instrument and may void the warranty.

Clean the exterior surface with a moist cloth or with a mixture of water and mild detergent. Dry with a soft cloth.

Maintenance menu

Select  **Maintenance** from the main menu.

Syringe activation

Do a check of the syringe. Make sure the syringe fills and empties correctly.

1. From the maintenance menu, push **Syringe activation**.

Option	Description
Fill	Fills the syringe with titrant solution. The process stops automatically when the syringe is full.

Option	Description
Empty to bottle	Discards the contents of the syringe into the titrant bottle.
Empty to beaker	Discards the contents of the syringe into the beaker. Make sure that the tube from the outlet port of the syringe is inside the beaker.
Stop	Stops the operation.

2. If a second syringe is installed, push **Toggle** to change to the second syringe.

Syringe replacement

To replace the syringe, select **Syringe replacement** from the maintenance menu. Obey the instructions on the display. Refer to [Install the syringe](#) on page 9.

Pump activation

This option only applies to instruments with pumps installed. Do a check of the pump. Make sure the pump fills and empties correctly.

1. From the maintenance menu, push **Pump activation**.



Option	Description
Start	Starts the pump. The reagent is pumped through the tubes into the beaker. Make sure that the tube from the outlet port of the pump is in the beaker.
Stop	Stops the operation.


2. If more than one pump is installed, push **Toggle** to change to another pump.

Pump cassette replacement

This option only applies to instruments with pumps installed.

To replace the pump cassette, refer to the illustrated steps that follow and obey the instructions on the display .

Other maintenance options

1. Select an option, then push **Select**.

Option	Description
Stirring activation	Do a check of the magnetic stirrer. Push the up and down arrow keys to increase or decrease the stirring speed.
Live measure	This option is not available for all sensors. The option shows continuous measurement data with connected sensors to quickly check measurements. The installed applications and the automatic additions to the sample are set to off. Continuous measurements are not temperature compensated, so measurement differences may occur in the same sample between continuous measurements and measurements that use installed applications with temperature compensation.
Maintenance summary	See the number of days remaining for maintenance tasks. After doing a task, push Reset to set the number of days remaining to the default value.
Maintenance schedule	See the list of maintenance tasks. Push Edit to change the default value.

Option	Description
Reagent replacement	Use this option to replace the reagents. Obey the instructions on the display.
Clean PtPt electrode	This option is only available if a PtPt sensor is connected to the instrument. Use a beaker for all the cleaning liquids. Obey the instructions on the display.

Troubleshooting

Refer to the following table for common problem messages or symptoms and possible causes.

Error/Warning	Description
Stop requested	The operator has pushed the stop button
Measurement unstable	The measurement is unstable
Measurement out of range	The measurement is out of range
Temperature out of range	The temperature measurement is out of range
Slope out of range	The calibration slope is outside accepted limits
Offset out of range	The calibration offset is outside accepted limits
Cell constant out of range	The calibration value of the cell constant is outside accepted limits
Equiv. point not found	The titration fails to determine the equivalence point
Out of range	The computed result is outside accepted limits
Electrode calibration has expired	The Calibration date has expired
Titrant calibration has expired	The Calibration date has expired
Maximum number of applications is reached	The maximum number of applications per line has been reached
Requires at least one titrant that cannot be installed	The application is not compatible with applications already installed as it uses a different titration
Requires too many titrants for this instrument model configuration	The instrument hardware configuration is not compatible with this application
Buffer unstable	The value of the buffer is unstable
The connected electrode(s) cannot be calibrated or there is no application that contains calibration parameters!	The type of sensor cannot be calibrated (e.g., PtPt, ORP, etc.) or there is no installed application that contains calibration parameters
Standard unstable	The standard is unstable
Same buffer	The buffer has already been used
Calibration solution already used	The calibration solution has already been used
No titrant information	System failure
No measure received	System failure
Burette fails to deliver	The syringe cannot be emptied
Burette fails to fill	The syringe cannot be filled
Burette failed to read delivered volume	System failure

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)