Algae & Cyanobacteria

Conditions for Customer Ownership

We hold permits allowing us to transport these organisms. To access permit conditions, <u>click here</u>. *Never purchase living specimens without having a disposition strategy in place*.

• There are currently no USDA permits required for any algae or cyanobacteria. Do not release into the environment.

Primary Hazard Considerations

Wash your hands before and after handling algae and cyanobacteria. It is not safe to release your algae into the environment because it could disrupt the normal ecosystem.


Availability

- Most algae and cyanobacteria are cultured in our labs and are available year-round.
- The following algae are collected, so shortages may occur:
 Nostoc colonial: Shortages possible December

 February and July

 –August.
 Porphyra: Shortages possible October

 –March.



Your algae culture will arrive in a plastic jar or sterile glass test tube filled with media. Immediately upon receipt, loosen the jar cover or test tube cap to allow gas exchange and store in a cool area (15–20°C) with dim light. Do not store your culture in direct sunlight or at temperatures above 30 °C. In its packing container, the culture can retain its high quality for 4–7 days at room temperature. If the culture is not going to be used in this time frame, it should be subcultured to prevent death and overcrowding.

Care

Habitat:

- Your cultures arrive in a habitat that is suitable for short-term use in classrooms.
- If you wish to subculture, a more specialized environment is required. Algae require specific lighting and media in order to flourish. Generally, cultures should be grown in a 16-hour light period alternating with an eight-hour dark period. Ideally, the cultures should be illuminated by 40-watt cool-white fluorescent tubes on a timer. A 40-watt fluorescent tube at a distance of about 15 cm will provide roughly 500 foot candles of illumination. Freshwater algal cultures should be grown under a light intensity of 400 to 500 foot candles. At this light intensity, cultures will reach optimum growth or density for microscopic viewing in 7–14 days, depending on the species and condition of the initial algal inoculum. After this period, reduce the light to 50–100 foot candles. Marine algae grow best in slightly lower intensities than those required by freshwater algae: 200–300 foot candles. The media that algal cultures are grown in varies based upon the requirements of each individual culture (see tables below, links at end of information section).

Subculturing:

Flasks, tubes, bottles, or Petri dishes can be used as culture containers. If using a 250 mL Erlenmeyer flask, fill the flask to approximately 150 mL with freshly prepared media and sterilize the media (autoclave for 20 minutes at 121°C and 15 psi). Allow the media to cool overnight before adding any algae. Add a small amount (2–5 mL) of inoculum from the culture provided, handling both containers in a sterile manner. Always prepare more than one subculture in case one of the new cultures becomes contaminated. Make new subcultures from the freshest cultures.



Basic Types of Algae:

Blue-Green Algae (Cyanobacteria)

The cyanophytes are the only prokaryotic algae. They are found in virtually every type of environment including terrestrial, freshwater, and marine habitats. Since cyanobacteria are prokaryotes, they lack membrane bound organelles. However, the external structure can range from unicellular or colonial, to branched or unbranched, and filamentous. The life-cycle is binary fission with the majority undergoing fragmentation to split. Spores/akinetes are produced as protected cells when unfavorable conditions exist. Like the rhodophytes, the cyanophytes possess no flagellated or ciliated cells at any stage of their lifecycle, although simple movements such as bending and swaying are made possible by internal pressure changes exerted on the cell wall. They are heavily pigmented with chlorophyll a, beta carotene, and several xanthophylls. The presence of several phycobiliproteins gives the cyanophyta their unique blue-green coloration. Food is stored in the form of glycogen.

Organism		Media (captive care)	Reproduction	Wild habitat	Shape/Characteristic (Special Notes)
Anabaena	86-1800	Basic	Trichome fragmentation	Freshwater	filamentous; unbranched; heterocysts
Cyanophora	86-1840	Basic	Binary fission	Freshwater	symbiont in protozoa
Cylindrospermum	86-1930	Basic	Trichome fragmentation	Freshwater	filamentous; heterocysts & akinetes terminal
Fischerella	86-1960	Basic	Trichome fragmentation	Freshwater	true branching, dimorphic structures
Glaucocystis	86-1980	Basic	Binary fission	Freshwater	symbiont, lives within colorless host
Gloeocapsa	86-2000	Basic	Mitosis and colony fragmentation	Freshwater	colonial; concentric layers of gelatinous sheath
Gloeotrichia	86-2010	Basic	Trichome fragmentation	Freshwater	filamentous; trichomes radiating from basil heterocysts
Lyngbya	86-2040	Basic	Trichome	Freshwater	filamentous; unbranched trichome, thinly sheathed
Merismopedia	86-2100	Soil-water	Binary fission	Freshwater	colonial; sheet of many cells one cell thick
Nostoc, colonial	86-2150	Basic	Trichome fragmentation	Freshwater	
				or Soil	macroscopic hollow balls; collected
Nostoc, filamentous	86-2155	Basic	Trichome fragmentation	Freshwater	filamentous; unbranched,
					contorted filaments
Oscillatoria	86-2300	Basic	Trichome fragmentation	Freshwater	filamentous; unbranched trichome,
					thinly sheathed
Scytonema	86-2350	Basic	Trichome fragmentation	Freshwater	filamentous; displays false branching
Spirulina	86-0700	Marine	Trichome fragmentation	Marine	brackish; helically coiled filament

Green Algae

Chlorophytes are a diverse group and are common in fresh water, salt water, and soil. They are very similar to plants, and most botanists agree the ancestor of higher plants can be found somewhere within this group. Chlorophyte reproduction varies greatly, from asexual division to isogamy and heterogamy to oogamy. Cell walls are constructed of cellulose and pectin. The food storage product is true starch, the same as plants. This can be demonstrated by staining with IKI, which turns the starch in the algae blue-black. Green algae possess true chloroplasts, which contain the same pigments found in higher plants: chlorophyll a and b, alpha and beta carotene, and many xanthophylls.

Organism		Media (captive care)	Wild habitat	Shape/Characteristic (Special Notes)
Ankistrodesmus	86-0010	Proteose agar	Freshwater	needle-like cells
Carteria	86-0080	Basic	Freshwater	large cells, four flagella
Chlamydomonas reinhardtii (+)	86-0102	Basic/Proteose agar	Freshwater	small, round, biflagellated
Chlamydomonas reinhardtii (-)	86-0103	Basic/Proteose agar	Freshwater	small, round, biflagellated
Chlamydomonas moewusii (+)	86-0104	Basic/Proteose agar	Freshwater	small, round, biflagellated
Chlamydomonas moewusii (-)	86-0105	Basic/Proteose agar	Freshwater	small, round, biflagellated
Chlorella	86-0126	Basic/Proteose agar	Freshwater	small, round, non-motile
Cladophora	86-0145	Erdschreiber's	Marine	filamentous; freely branched
Closterium ehrenbergii	86-0160	Basic	Freshwater	elongated desmid
Closterium littorale	86-0164	Basic	Freshwater	homothallic strain
Cosmarium	86-0165	Basic	Freshwater	circular desmid
Dunaliella	86-0171	Erdschreiber's	Marine	halophilic flagellate

Ward's science

– continued –

		Media		Shape/Characteristic
Organism		(captive care)	Wild habitat	(Special Notes)
Eremosphaera	86-0175	Basic	Freshwater	circular, large; macroscopic
Eudorina	86-0180	Basic	Freshwater	colonial; with spherical, biflagellated cells
Fritschiella	86-0181	Basic	Freshwater	filamentous; displays differentiation
Gonium	86-0185	Basic	Freshwater	colonial; 4-32 cells, flattened colony
Hydrodictyon	86-0210	Basic	Freshwater	colonial; net-like, multinucleated cells
Micrasterias	86-0270	Basic	Freshwater	circular, large desmid
Microspora	86-0275	Basic	Freshwater	filamentous; H-shaped cell walls
Mougeotia	86-0280	Basic	Freshwater	filamentous; relative to Spirogyra
Netrium	86-0290	Basic	Freshwater	oblong desmid
Nitella	86-0300	Spring water	Freshwater	macroalga; used for cytoplasmic streaming; collected
Oedogonium cardiacum (+)	86-0350	Basic	Freshwater	filamentous; heterothallic strain
0edogonium cardiacum (-)	86-0351	Basic	Freshwater	filamentous; heterothallic strain
Oedogonium foveolatum	86-0355	Basic	Freshwater	filamentous; homothallic strain
Pandorina	86-0400	Basic	Freshwater	colonial; elliptical; made up of 16-32 cells
Pediastrum	86-0425	Basic	Freshwater	colonial; polygonal cells
Protococcus	86-0495	None	Tree bark	colonial; on wood bark; collected
Scenedesmus	86-0600	Basic	Freshwater	colonial; four cells with spines
Selenastrum	86-0620	Proteose agar	Freshwater	lunate
Spirogyra	86-0650	Soil-water	Freshwater	filamentous; spiral chloroplasts
Staurastrum	86-0660	Basic	Freshwater	triangular desmid
Stigeoclonium	86-0690	Basic	Freshwater	filamentous; produces zoospores
Ulothrix	86-0750	Basic	Freshwater	filamentous; unbranched filaments
Ulva	86-0770	Sea water	Marine	macroalga, membranous thallus; collected
Volvox aureus	86-0800	Basic/soil	Freshwater	colonial; cytoplasmic connections about the size of flagella
Volvox globator	86-0805	Basic/soil	Freshwater	colonial; thickest cytoplasmic collections
Zygnema	86-0900	Basic	Freshwater	filamentous; 2 stellate chloroplasts per cell

Diatoms

Bacillariophytes, which occur in fresh water, salt water, and terrestrially, date back to the Cretaceous Period. They are single-celled algae with shells constructed of two overlapping valves composed of pectin and impregnated with silica; these shells can be quite ornate. Although the diatoms are single-celled organisms, they can form colonies and filaments. The group comprises two main types: centric and pennate. Centric diatoms are radially symmetrical and contain numerous plasmids, while pennate diatoms are bilaterally symmetrical and contain fewer plasmids. Many diatoms have conspicuous oil droplets within the cell, which is the photosynthetic food reserve, chrysolaminarin. The plastids of diatoms contain the pigments chlorophyll a and b, alpha and beta carotene, and several xanthophylls.

Organism		Media (captive care)	Wild habitat	Shape/Characteristic (Special Notes)
Cyclotella	86-0168	Marine	Marine	small, centric
Navicula	86-1210	Basic	Freshwater	pennate, with bilateral symmetry
Synedra	86-0692	Basic	Freshwater	pennate, with bilateral symmetry
Thalassiosira	86-0695	Marine	Marine	brackish, centric



Dinoflagellates

Dinoflagellates are mostly marine organisms and they compose nearly all marine plankton. They occur as free-living flagellates, sessile unicells, colonies, and filamentous forms. The fossil record of the dinoflagellates can be dated back to the Cambrian period, with some evidence suggesting they existed even earlier. The term dinoflagellate actually refers to the twirling motion exhibited by the pair of whip-like undulapodia (flagella). These flagella originate in the sulcus, or groove, of the organism. Some dinoflagellates have thecal plates embedded in their cytoplasmic membrane and are called armored, others lack these plates and are called naked. Food is stored in the form of true starch and oils. Dinoflagellates contain the pigments chlorophyll a and c, beta carotene, and several xanthophylls that often give these organisms a brownish color. Some dinoflagellates produce powerful toxins with potentially dangerous results. When "blooms" occur, the water can take on a pinkish or red hue known as a red tide. This often causes massive fish kills and can be dangerous to humans as well. Some, such as Noctiluca, are bioluminescent, and can cause ocean waves to glow at night. This is the only example of bioluminescence in the algae kingdom.

Organism		Media (captive care)	Wild habitat	Shape/Characteristic (Special Notes)
Amphidinium	86-1560	Erdschreiber's	Marine	Lower cell half is longer and broader than upper
Peridinium	86-2900	Soil-water	Freshwater	Thick, armored plates
Prorocentrum	86-2905	Marine	Marine	Prominent anterior spine

Euglenoids

Typically green and unicellular, euglenoid flagellates live in fresh water. They have characteristics of both plants and animals yet are distinct in many ways. Most are photosynthetic, but many, lacking chloroplasts, are heterotrophs. Most do not reproduce sexually. Euglenoids lack a cellulose cell wall; instead, they have a proteinaceous pellicle just inside the plasmalemma. The plastids contain chlorophyll a and b, beta carotene, and xanthophylls. If placed in the dark over the course of several divisions, the chloroplasts of *Euglena gracilis* will become colorless. When returned to the light, the plastid structure is reformed and the green color returns.

		Media		Shape/Characteristic
Organism		(captive care)	Wild habitat	(Special Notes)
Euglena acus	86-2600	Euglena or soil– water	Freshwater	large; slow swimmer
Euglena gracilis	86-2650	Euglena agar	Freshwater	Z strain; used in vitamin B ₁₂ assay
Euglena sp.	87-0100	Euglena	Freshwater	good for general flagellates study
Phacus	86-2750	Soil-water	Freshwater	heart-shaped
Trachelomonas	86-2800	Soil-water	Freshwater	loricated flagellate

Brown Algae

Multicellular and structurally complex, with no colonies or simple, unbranched filaments, the Phaeophytes, or brown seaweed, are primarily marine algae; less than one percent occur in fresh water. They are most abundant and reach their maximum development in the colder water of the oceans. While some species of Sargassum are found floating in enormous numbers in the Atlantic, the algae are usually firmly attached to a substrate by means of elaborate holdfast structures. Food is stored as soluble carbohydrates such as lamarin, fats and the alcohol mannitol. The plastids of the brown algae contain pigments chlorophyll a and c, c-carotene, and xanthins; an accessory pigment, fucoxanthin, gives the algae their characteristic dark brown or olive green color. The Phaeophytes are an economically important resource, used for alginic acid, fertilizer, and food. The brown algae are all collected from marine sources.

Organism		Media (captive care)	Wild habitat	Shape/Characteristic (Special Notes)
Dictyota	86-1940	Seawater	Marine	thin, flat forking branches; macroalgae; collected
Ectocarpus	86-1950	Seawater	Marine	filamentous; branched thallus, zoospores
Fucus	86-3921	Seawater	Marine	bladderlike floats, disk-shaped holdfasts for clinging to rocks; macroalgae; collected
Laminaria	86-3920	Seawater	Marine	long, leathery laminae and large in size; kelp; collected



Golden Algae

Chrysophytes are a large and complex group characterized by plastids containing distinctive golden yellow pigments. The group is diverse in form, yet all feature this yellow color, permitting easy identification. Chrysophytes are usually found in cold freshwater lakes and ponds, although some marine forms are common. *Synura*, existing in colonies in fresh water, can cause a fishy odor in reservoirs even in low concentrations, but is not harmful.

Organism		Media (captive care)	Wild habitat	Shape/Characteristic (Special Notes)
Coccolithophora	86-2915	Marine	Marine	coccoid, from coccoliths
Ochromonas	86-1250	Ochromonas	Freshwater	two unequal flagella
Synura	86-1350	Soil-water	Freshwater	colonial; radically arranged colonies

Red Algae

Although Rhodophytes, the most abundant type of seaweed, are widely distributed in the oceans, most occur in tropical and subtropical littoral zones. Of the 4,000 species, the vast majority are marine. Rhodophytes are not mobile –they possess no flagellated or ciliated cells at any stage of their life cycle — yet all reproduce sexually. Many red algae, such as *Corallina*, are calcified and encrusted appearing much like coral. This calcification has made it possible to trace the Rhodophytes to the Paleozoic Period. Single-celled forms such as *Porphyridium* are a rarity. Rhodophytes are characterized by reddish plastids, called rhodoplasts, which contain the pigments chlorophyll a and d, alpha and beta carotene, some xanthophylls and phycobiliproteins.

Organism		Media (captive care)	Wild habitat	Shape/Characteristic (Special Notes)
Bangia	86-2510	Erdschreiber's	Marine	filamentous; unbranched thallus
Batrachospermum	86-2825	Basic	Freshwater	filamentous
Callithamnion	86-1568	Erdschreiber's	Marine	filamentous
Corallina	86-3933	Seawater	Marine	macroalga; pink, coral like appearance; collected
Gigartina	86-3926	Seawater	Marine	macroalga, common high-intertidal alga; collected
Polysiphonia	86-2910	Seawater	Marine	filamentous; polysiphonious construction; collected
Porphyra	86-3929	Seawater	Marine	macroalga, long, irregular lobed; collected
Porphyridium	86-2850	Erdschreiber's	Marine	thin, gelatinous, blood red layer

Yellow-Green Algae

Xanthophytes are highly successful in fresh water and terrestrial environments, although some marine forms also exist. The yellow-green algae have pectin-rich cellulose walls. Starch is absent and food is stored in the form of oils. Xanthophytes are characterized by yellow-green plastids (xanthoplasts) which contain pigments chlorophyll a and c, several xanthins, and beta carotene. *Vaucheria* is a large, macroscopic, filamentous form that was classified a chlorophyte until pigment analysis showed the absence of chlorophyll b and true starch. *Tribonema* is a typical freshwater, unbranched, filamentous form which clearly demonstrates overlapping walls.

Organism		Media (captive care)	Wild habitat	Shape/Characteristic (Special Notes)
Tribonema	86-1400	Basic	Freshwater	filamentous; H-shaped overlapping cell wall
Vaucheria	86-1500	Soil-water	Freshwater	filamentous; oogamous



Media for Algae and Cyanobacteria

See the label or tables above for the appropriate media for each specimen.

Basic Culture Solution Concentrate

For culture of freshwater algae. Sterile. Makes one liter.

88 W 3250-100 mL bottle

Basic Culture Solution, Working Solution

For culture of freshwater algae. Sterile.

88 W 3251-1 L bottle

Euglena Medium

Sterile.

88 W 5200-1 L bottle

Erdschreiber's Medium

For culture of marine algae. Sterile.

88 W 4351-1 L bottle

Marine Algae Culture Medium

For culture of macroscopic and microscopic marine algae. Sterile. Makes 1 liter.

88 W 4155—125 mL bottle

Soil-Water Medium

For culture of freshwater forms. Sterile.

88 W 4150-1 L bottle

Spring Water

To maintain *Hydra*, *Planaria*, protists, and algae.

88 W 7000—1 gal. jug

Water, distilled

For defined media.

88 W 7005—1 gal.

Pond Water

Non-sterile; may contain organisms.

88 W 7010—1 gal. jug

Seawater

From our marine tanks. Used in culturing saltwater invertebrates or in media recipes. Non-sterile.

88 W 7011—1 gal. jug

Disposition

- Please dispose of excess living material in a manner to prevent spread into the environment. Consult with your school to identify their preferred method of disposal.
- You can safely use one of the following methods:
 - Treat culture with a 10% bleach solution for 24 hours (1 part bleach to 9 parts culture medium or water culture medium removed). Then rinse bleach solution down the drain with water until you can no longer smell bleach. Rinse remaining materials and containers with water and dispose of them in a general garbage container.
 - Carefully wrap specimens and their containers in a biohazard bag (without containing anything sharp that might puncture the bag) and tie closed (a twist tie works well). Autoclave the bag for 30 minutes at 121 degrees C and at a pressure of 15 lbs. per square inch. Dispose of autoclaved bag as your school recommends.

