

Molecular Cloning

TECHNICAL GUIDE

UPDATE
2016

VWR  | vwr.com
We Enable Science

 NEW ENGLAND
BioLabs[®] Inc.

be INSPIRED
drive DISCOVERY
stay GENUINE



Molecular Cloning Overview

Molecular cloning refers to the process by which recombinant DNA molecules are produced and transformed into a host organism, where they are replicated. A molecular cloning reaction is usually comprised of two components:

1. The DNA fragment of interest to be replicated.
2. A vector/plasmid backbone that contains all the components for replication in the host.

DNA of interest, such as a gene, regulatory element(s), operon, etc., is prepared for cloning by either excising it out of the source DNA using restriction enzymes, copying it using PCR, or assembling it from individual oligonucleotides. At the same time, a plasmid vector is prepared in a linear form using restriction enzymes (REs) or Polymerase Chain Reaction (PCR). The plasmid is a small, circular piece of DNA that is replicated within the host and exists separately from the host's chromosomal or genomic DNA. By physically joining the DNA of interest to the plasmid vector through phosphodiester bonds, the DNA of interest becomes part of the new recombinant plasmid and is replicated by the host. Plasmid vectors allow the DNA of interest to be copied easily in large amounts, and often provide the necessary control elements to be used to direct transcription and translation of the cloned DNA. As such, they have become the workhorse for many molecular methods such as protein expression, gene expression studies, and functional analysis of biomolecules.

During the cloning process, the ends of the DNA of interest and the vector have to be modified to make them compatible for joining through the action of a DNA ligase, recombinase, or an *in vivo* DNA repair mechanism. These steps typically utilize enzymes such as nucleases, phosphatases, kinases and/or ligases. Many cloning methodologies and, more recently kits have been developed to simplify and standardize these processes.

This technical guide will clarify the differences between the various cloning methods, identify NEB® products available for each method, and provide expert-tested protocols and FAQs to help you troubleshoot your experiments.

TABLE OF CONTENTS

| | |
|-------|--|
| 3–5 | Cloning Workflow Descriptions |
| 3 | Seamless Cloning |
| 4 | Traditional Cloning |
| 4 | PCR Cloning |
| 5 | Ligation Independent Cloning (LIC) |
| 5 | Recombinational Cloning |
| 6–7 | Cloning Workflow Comparison |
| 8 | Migrate to Monarch |
| 9–18 | DNA Preparation |
| 9 | cDNA Synthesis |
| 10 | Restriction Enzyme Digestion |
| 10 | Protocol |
| 10 | Optimization Tips |
| 11–16 | Performance Chart |
| 17–18 | PCR |
| 17 | Protocols |
| 17 | Optimization Tips |
| 18 | Product Selection |
| 19–22 | Common DNA End Modifications |
| 19 | Phosphorylation |
| 19 | Protocol |
| 19 | Optimization Tips |
| 20 | Dephosphorylation |
| 20 | Protocol |
| 20 | Optimization Tips |
| 20 | Product Selection |
| 21 | Blunting/End-repair |
| 21 | Protocol |
| 21 | Optimization Tips |
| 21 | Product Selection |
| 22 | A-tailing |
| 22 | Protocol |
| 22 | Product Selection |
| 22 | Activity in CutSmart® Buffer |
| 23–24 | Vector and Insert Joining |
| 23–24 | DNA Ligation |
| 23 | Protocol |
| 23 | Optimization Tips |
| 24 | Product Selection |
| 25 | Transformation |
| 25 | Protocol |
| 25 | Optimization Tips |
| 25 | Product Selection |
| 26 | DNA Analysis |
| 26 | Product Selection |
| 27–31 | Cloning & Mutagenesis |
| 27 | NEBuilder® HiFi DNA Assembly |
| 27 | Overview |
| 28 | Protocol/Optimization Tips |
| 28 | Gibson Assembly® |
| 29 | NEB Golden Gate Assembly |
| 29 | Overview |
| 30 | NEB PCR Cloning Kit |
| 30 | Overview/Protocols |
| 31 | Q5® Site-Directed Mutagenesis Kit |
| 31 | Protocols/Optimization Tips |
| 32 | DNA Assembly Selection Chart |
| 33–34 | Traditional Cloning Quick Guide |
| 35–37 | Troubleshooting Guide |
| 38–39 | Ordering Information |



Cloning Workflow Descriptions

There are several methods that can be used to generate DNA constructs, each of which is described below. A comparison of the various workflows discussed can be found on page 6.

Seamless Cloning/Gene Assembly

The group of cloning methods we refer to as “seamless cloning” typically combine attributes from more established cloning methods to create a unique solution to allow sequence-independent and scarless insertion of one or more DNA fragments into a plasmid vector. Various commercial systems, such as NEBuilder HiFi DNA Assembly, NEB Gibson Assembly and In-Fusion® employ PCR to amplify the gene of interest, an exonuclease to chew back one strand of the insert and vector ends, and either a ligase, recombination event, or *in vivo* repair to covalently join the insert to the vector through a true phosphodiester bond. The ability to quickly join a single insert to a plasmid at any sequence in the vector, without a scar, makes these technologies very appealing cloning methods. Additionally, the ability to join 5–10 fragments in a predetermined order, with no sequence restrictions or scars, provides a powerful technique for synthetic biology endeavors, such as moving whole operons for metabolic engineering or whole genome reconstructions.


ADVANTAGES

- No sequence constraints
- Efficient assembly of multiple fragments
- High cloning efficiency
- Exquisite control of higher-order gene assembly


DISADVANTAGES

- Cost, relative to traditional methods
- PCR primers for vector and insert must be designed and ordered

Overview of the NEBuilder HiFi DNA Assembly cloning method



NEB Golden Gate Assembly workflow



In its simplest form, Golden Gate Assembly requires a Type IIS recognition site, in this case, Bsal (GGTCTC), added to both ends of a dsDNA fragment. After digestion, these sites are left behind, with each fragment bearing the designed 4-base overhangs that direct the assembly.

Golden Gate Assembly is another method of seamless cloning that exploits the ability of Type IIS restriction enzymes (such as Bsal) to cleave DNA outside of the recognition sequence. The inserts and cloning vectors are designed to place the Type IIS recognition site distal to the cleavage site, such that the Type IIS restriction enzyme can remove the recognition sequence from the assembly. The advantages of such an arrangement are three-fold: 1. the overhang sequence created is not dictated by the restriction enzyme, and therefore no scar sequence is introduced; 2. the fragment-specific sequence of the overhangs allows orderly assembly of multiple fragments simultaneously; and 3. the restriction site is eliminated from the ligated product, so digestion and ligation can be carried out simultaneously. The net result is the ordered and seamless assembly of DNA fragments in one reaction.



Traditional Cloning

Traditional Cloning usually refers to the use of restriction endonucleases to generate DNA fragments with specific complementary end sequences that can be joined together with a DNA ligase, prior to transformation. This typically involves preparing both a DNA fragment to be cloned (insert) and a self-replicating DNA plasmid (vector) by cutting with two unique restriction enzymes that flank the DNA sequence, and whose cut sites are present at the preferred site of insertion of the vector, often called the multiple cloning site (MCS). By using two different REs, two non-compatible ends are generated, thus forcing the insert to be cloned directionally, and lowering the transformation background of re-ligated vector alone. Directional cloning is useful to maintain open reading frames or another positional requirement with *cis*-acting regulatory elements. Non-directional cloning can also be performed with compatible ends generated by a single restriction enzyme; in this case the clones will need to be screened to determine that the gene orientation is correct. Typically the vector needs to be de-phosphorylated to prevent self-ligation, which directly competes with the insert and lowers the efficiency of the cloning reaction.

In the early years of cloning, genomic DNA was often cloned into plasmid vectors using DNA adaptors to add the required restriction sites to a sequence of interest, prior to ligation. Additionally, genes or other DNA elements were swapped between vectors using compatible ends contained by both vectors. More recently, PCR is used as an upstream step in a cloning protocol to introduce the necessary restriction sites for directional cloning prior to preparation of the vector and insert by restriction digests, followed by fragment purification, fragment ligation, and transformation into an *E. coli* cloning strain for plasmid amplification. Transformed colonies, now resistant to an antibiotic due to a resistance gene harbored by the plasmid, are screened by colony PCR or restriction digest of plasmid DNA for the correct insert. Direct sequencing of the recombinant plasmid is often performed to verify the sequence integrity of the cloned fragment.

PCR Cloning

PCR cloning differs from traditional cloning in that the DNA fragment of interest, and even the vector, can be amplified by PCR and ligated together without the use of restriction enzymes. PCR cloning is a rapid method for cloning genes, and is often used for projects that require higher throughput than traditional cloning methods can accommodate. It also allows for the cloning of DNA fragments that are not available in large amounts. Typically, a PCR reaction is performed to amplify the sequence of interest and then it is joined to the vector via a blunt or single-base overhang ligation prior to transformation. Early PCR cloning often used *Taq* DNA Polymerase to amplify the gene. This results in a PCR product with a single template-independent base addition of an adenine (A) residue to the 3' end of the PCR product, through the normal action of the polymerase. These "A-tailed" products are then ligated to a complementary T-tailed vector using T4 DNA Ligase, followed by transformation. High-fidelity polymerases are now routinely used to amplify DNA sequences with the PCR product containing no 3' extensions. The blunt-end fragments are joined to a plasmid vector through a typical ligation reaction or by the action of an "activated" vector that contains a covalently attached enzyme, typically Topoisomerase I, that facilitates the vector:insert joining. PCR cloning with blunt-end fragments is non-directional. Some PCR cloning systems contain engineered "suicide" vectors that include a toxic gene into which the PCR product must be successfully ligated to allow propagation of the strain that takes up the recombinant molecule during transformation. A typical drawback common to many PCR cloning methods is that a dedicated vector must be used. These vectors are typically sold by suppliers, like NEB, in a ready-to-use, linearized format and can add significant expense to the total cost of cloning. Also, the use of specific vectors restricts the researcher's choice of antibiotic resistance, promoter identity, fusion partners, and other regulatory elements.

ADVANTAGES

- Low cost
- Versatile
- Many different vector choices
- Directional cloning can be easily done

DISADVANTAGES

- Possible sequence constraints due to presence and/or translation of restriction site

ADVANTAGES

- High efficiency, with dedicated vectors
- Amenable to high throughput

DISADVANTAGES

- Higher cost
- Multi-fragment cloning is not straight forward
- Directional cloning is difficult



Ligation Independent Cloning (LIC)

Ligation Independent Cloning (LIC) is a technique developed in the early 1990s as an alternative to restriction enzyme/ligase cloning. Inserts are usually PCR amplified, and vectors are made linear either by restriction enzyme digestion or by PCR. This technique uses the 3'→5'-exo activity of T4 DNA Polymerase to create overhangs with complementarity between the vector and insert. Incorporation of only dGTP in the reaction limits the exonuclease processing to the first complementary C residue, which is not present in the designed overlap, where the polymerization and exonuclease activities of T4 DNA Polymerase become “balanced”. Joined fragments have 4 nicks that are repaired by *E. coli* during transformation. This technique allows efficient creation of scarless recombinant plasmids at many, but not all, positions in a vector.

More recently, the technique has evolved to include many useful variations. One in particular, Sequence and Ligation Independent Cloning (SLIC), has been adopted by many researchers. In this variation, all dNTPs are initially excluded from the reaction with T4 DNA Polymerase. This allows the exo activity of T4 DNA Polymerase to proceed and generate the complementary overlaps between insert and vector. After the overlap is generated, dCTP is added back to the reaction, which shifts the enzyme back into a polymerase. It then stalls due to the lack of a complete set of dNTPs in the buffer, and the complementary overlap is retained. The product contains 4 nicks, just like the original LIC product, and is repaired by *E. coli* during transformation. This modification of the protocol allows a scarless and sequence-independent insertion into nearly any vector.

ADVANTAGES

- Low cost
- Many different vector choices

DISADVANTAGES

- Some types of sequence modifications not possible

Recombinational Cloning

Recombinational cloning became popular with the introduction of three cloning systems: Gateway[®], Creator[™], and Echo Cloning[™] systems. These systems use a site-specific recombinase (Integrase in Gateway and Cre Recombinase in Creator and Echo) to allow the reliable transfer of a fragment from one vector to another without using restriction enzymes and ligases. Typically, a researcher would clone a sequence of interest into a holding vector (“Entry” for Gateway and “Donor” for Creator) using traditional cloning methods. Once the new clone is made, it is easily shuttled to many different “destination” or “acceptor” vectors that contain the appropriate sequence recognized by the recombinase (attachment sites *attB* and *attP* with Gateway and *loxP* with Creator/Echo). Higher throughput is possible with these systems and they have become a useful tool for screening many different expression hosts for protein expression projects or for multiple reporter vectors for functional analysis studies. At this time, only the Gateway system is still commercially supported, although NEB does sell Cre Recombinase, an essential reagent for the *in vitro* recombination step used by the Creator and Echo Cloning systems.


ADVANTAGES

- Allows high-throughput vector creation
- Widely available ORF collections

DISADVANTAGES

- Cost relative to traditional methods
- Vector sets typically defined by supplier
- Proprietary enzyme mixes often required

Cre/*loxP* Site-specific Recombination







Cloning Workflow Comparison



The figure below compares the steps for the various cloning methodologies, along with the time needed for each step in the workflows.

INSERT PREPARATION



* Note that times are based on estimates for moving a gene from one plasmid to another. If the source for gene transfer is gDNA, add 2 hours to calculation for the traditional cloning method. Total time does not include transformation, isolation, or analysis.

** 70 minutes for recombination occurs on second day.







SELECTION CHARTS & PROTOCOLS

Need help with locating product selection charts & protocols?


- 10 cDNA Synthesis
- 11 Restriction Enzymes
- 17 PCR
- 19 Phosphorylation
- 20 Dephosphorylation
- 21 Blunting
- 22 A-tailing
- 23 Ligation
- 25 Transformation
- 26 DNA Analysis
- 27 Cloning & Mutagenesis

VECTOR PREPARATION



+

ional (Univector)
Up
est (Standard) n (Time Saver)
Up
pecific combination
ase K Treatment
min **
Endpoint vector
5 hr 20 min



Migrate to Monarch[®].

Environmentally-friendly Nucleic Acid Purification Kits from New England Biolabs

Want to feel good about your choice in DNA purification? With our fast and reliable Monarch Nucleic Acid Purification Kits, you can achieve optimal purification while creating less waste. Available for plasmid minipreps, DNA gel extraction and enzymatic cleanup (including PCR), our products use up to 44% less plastic and are packaged using responsibly-sourced, recyclable materials. Make the change and migrate to Monarch today.



“ *These kits might be the best I have used for the price. The best part is that it uses less plastic for production!! Thank you for caring about our environmental impacts, NEB!!!* ”

– NEB customer



DNA Preparation

The first step in any of the cloning workflows mentioned previously is the preparation of DNA. In most cases, this involves the preparation of the vector backbone and insert. When starting with DNA, restriction enzyme digestion or PCR is performed. When starting with RNA, a cDNA synthesis step is performed using a reverse transcriptase.

cDNA Synthesis

When RNA is used as starting material, a reverse transcriptase can be used to generate cDNA, which can then be used as template for any of the cloning methods listed previously. Depending on which workflow is being followed, the resulting DNA may require a clean-up step. This can be performed using a spin column or by gel extraction.

Protocol: cDNA Synthesis

| | DENATURATION PROTOCOL |
|--|---|
| Total RNA | 1–6 μ l (up to 1 μ g) |
| d(T)₂₃ VN (50 μM) | 2 μ l |
| Nuclease-free Water | to a total volume of 8 μ l |
| Incubation | 65°C for 5 minutes spin briefly and put on ice |

| | SYNTHESIS PROTOCOL |
|----------------------|--------------------------------------|
| Denatured RNA | 8 μ l |
| Reaction Mix | 10 μ l |
| Enzyme Mix | 2 μ l |
| Incubation | 80°C for 5 minutes store at –20°C |

cDNA Synthesis Selection Chart

| cDNA SYNTHESIS | FEATURES |
|--|---|
| KITS | |
| ProtoScript® II First Strand cDNA Synthesis Kit | Generates cDNA at least 10 kb in length |
| | Contains ProtoScript II Reverse Transcriptase, an enzyme with increased thermostability and reduced RNase H activity |
| | Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix |
| ProtoScript First Strand cDNA Synthesis Kit | Generates cDNA at least 5 kb in length |
| | Contains M-MuLV Reverse Transcriptase |
| | Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix |
| STANDALONE REAGENTS | |
| ProtoScript II Reverse Transcriptase An alternative to SuperScript® II | RNase H ⁻ mutant of M-MuLV Reverse Transcriptase with increased thermostability and reduced RNase H activity |
| | Increased reaction temperatures (37–50°C) |
| M-MuLV Reverse Transcriptase | Robust reverse transcriptase for a variety of templates |
| | Standard reaction temperatures (37–45°C) |
| AMV Reverse Transcriptase | Robust reverse transcriptase for a broad temperature range (37–52°C) |
| | Can be used for templates requiring higher reaction temperatures |

TIPS FOR OPTIMIZATION

STARTING MATERIAL

- Intact RNA of high purity is essential for generating cDNA for cloning applications.
- Total RNA or mRNA can be used in the reverse transcription reaction. Total RNA is generally sufficient for cDNA synthesis reactions. However, if desired, mRNA can be easily obtained using a PolyA Spin mRNA Isolation Kit or Magnetic mRNA Isolation Kit.
- The amount of RNA required for cDNA cloning depends on the abundance of the transcript-of-interest. In general, 1 ng to 1 μ g total RNA or 0.1–100 ng mRNA are recommended.

PRODUCT SELECTION

- Streamline your reaction setup by using the ProtoScript II First Strand cDNA Synthesis Kit. This kit combines ProtoScript II Reverse Transcriptase, a thermostable M-MuLV (RNase H⁻) Reverse Transcriptase, and recombinant RNase Inhibitor in an enzyme Master Mix, along with a separate Reaction Mix containing dNTPs. Additionally, the kit contains two optimized reverse transcription primer mixes.

YIELD

- ProtoScript II Reverse Transcriptase is capable of generating cDNA of more than 10 kb up to 48°C. We recommend 42°C for routine reverse transcription.
- You can increase the yield of a long cDNA product by doubling the amount of enzyme and dNTPs.

ADDITIVES

- For most RT-PCR reactions, RNase H treatment is not required. But for some difficult amplicons or sensitive assays, add 2 units of *E.coli* RNase H to the reaction and incubate at 37°C for 20 minutes.



Restriction Enzyme Digestion


Restriction enzyme sites that are unique to both the insert and vector should be chosen. Unidirectional cloning is achieved using two different restriction enzymes, each with unique recognition sites at an end of the insert. Depending on the RE chosen, ends can be blunt or sticky (cohesive). Restriction enzyme digestion is generally used in traditional cloning.

Protocol: Restriction Enzyme Reactions

| | STANDARD PROTOCOL | TIME-SAVER® PROTOCOL |
|-------------------------------|-------------------|----------------------|
| DNA | up to 1 µg | up to 1 µg |
| 10X NEBuffer | 5 µl (1X) | 5 µl (1X) |
| Restriction Enzyme | 10 units* | 1 µl |
| Total Volume | 50 µl | 50 µl |
| Incubation Temperature | enzyme dependent | enzyme dependent |
| Incubation Time | 60 minutes | 5–15 minutes** |

*Sufficient to digest all types of DNAs.

**Time-Saver qualified enzymes can also be incubated overnight with no star activity.



EcoRI-HF shows no star activity in overnight digests, even when used at higher concentrations. 50 µl reactions were set up using 1 µg of Lambda DNA, the indicated amount of enzyme and the recommended reaction buffer. Reactions were incubated overnight at 37°C. Marker M is the 1 kb DNA Ladder.

TIPS FOR OPTIMIZATION

ENZYME

- Keep on ice when not in the freezer
- Should be the last component added to reaction
- Mix components by pipetting the reaction mixture up and down, or by “flicking” the reaction tube. Follow with a quick (“touch”) spin-down in a microcentrifuge. Do not vortex the reaction.
- In general, we recommend 5 – 10 units of enzyme per µg DNA, and 10 – 20 units per µg of genomic DNA in a 1 hour digest

STAR ACTIVITY

- Unwanted cleavage that can occur when an enzyme is used under sub-optimal conditions, such as:
 - Too much enzyme present
 - Too long of an incubation time
 - Using a non-recommended buffer
 - Glycerol concentrations above 5%
- Star activity can be reduced by using a High-Fidelity (HF) enzyme, reducing the number of units, reducing incubation time, using a Time-Saver enzyme or increasing reaction volume

DNA

- Should be free of contaminants such as phenol, chloroform, alcohol, EDTA, detergents or excessive salts. Extra wash steps during purification are recommended.
- Methylation of DNA can effect digestion with certain enzymes.

BUFFER

- Use at a 1X concentration
- BSA is included in NEBuffer 1.1, 2.1, 3.1 and CutSmart Buffer. No additional BSA is needed.
- Restriction enzymes that do not require BSA for optimal activity are not adversely affected if BSA is present in the reaction

REACTION VOLUME

- A 50 µl reaction volume is recommended for digestion of up to 1 µg of substrate. This helps maintain salt levels introduced by miniprep DNA low enough that they don't affect enzyme activity.
- Enzyme volume should not exceed 10% of the total reaction volume to prevent star activity due to excess glycerol
- Additives in the restriction enzyme storage buffer (e.g., glycerol, salt), as well as contaminants found in the substrate solution (e.g., salt, EDTA, or alcohol), can be problematic in smaller reaction volumes

| | RESTRICTION ENZYME* | DNA | 10X NEBUFFER |
|-------------|---------------------|--------|--------------|
| 10 µl rxn** | 1 unit | 0.1 µg | 1 µl |
| 25 µl rxn | 5 units | 0.5 µg | 2.5 µl |
| 50 µl rxn | 10 units | 1 µg | 5 µl |

* Restriction Enzymes can be diluted using the recommended diluent buffer when smaller amounts are needed

** 10 µl rxns should not be incubated for longer than 1 hour to avoid evaporation

INCUBATION TIME

- Incubation time for the Standard Protocol is 1 hour. Incubation for the Time-Saver Protocol is 5–15 minutes.
- It is possible, with many enzymes, to use fewer units and digest for up to 16 hours.

STORAGE AND STABILITY

- Storage at –20°C is recommended for most restriction enzymes. For a few enzymes, storage at –80°C is recommended.
- 10X NEBuffers should be stored at –20°C
- The expiration date is found on the label
- Long term exposure to temperatures above –20°C should be minimized whenever possible



Performance Chart for Restriction Enzymes

New England Biolabs supplies > 210 restriction enzymes that are 100% active in a single buffer, CutSmart. This results in increased efficiency, flexibility and ease-of-use, especially when performing double digests.

This performance chart summarizes the activity information of NEB restriction enzymes. To help select the best conditions for double digests, this chart shows the optimal (supplied) NEBuffer and approximate activity in the four standard NEBuffers for each enzyme. Note that BSA is included in all NEBuffers. In addition, this performance chart shows recommended reaction temperature, heat-inactivation temperature, recommended diluent buffer, methylation sensitivity, whether the enzyme is Time-Saver qualified (cleaves substrate in 5–15 minutes under recommended conditions, and can be used overnight without degradation of DNA), and whether the enzyme works better in a substrate with multiple sites.

Chart Legend

| | | | |
|----------|--|------------|---|
| U | Supplied with a unique reaction buffer that is different from the four standard NEBuffers. The compatibility with the four standard NEBuffers is indicated in the chart. | SAM | Supplied with a separate vial of S-adenosylmethionine (SAM). To obtain 100% activity, SAM should be added to the 1X reaction mix as specified on the product data card. |
| | Recombinant | | dcm methylation sensitivity |
| | Time-Saver qualified | | CpG methylation sensitivity |
| | Engineered enzyme for maximum performance | | Indicates that the restriction enzyme requires two or more sites for cleavage |
| | dam methylation sensitivity | | |

Activity Notes (see last column)

FOR STAR ACTIVITY

- Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.
 - Star activity may result from extended digestion.
 - Star activity may result from a glycerol concentration of > 5%.
- * May exhibit star activity in this buffer.

FOR LIGATION AND RECUTTING

- Ligation is less than 10%
- Ligation is 25% – 75%
- Recutting after ligation is < 5%
- Recutting after ligation is 50% – 75%
- Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

NEBuffer Compositions (1X)

| | |
|---------------------|--|
| NEBuffer 1.1 | 10 mM Bis Tris Propane-HCl, 10 mM MgCl ₂ , 100 µg/ml BSA (pH 7.0 @ 25°C). |
| NEBuffer 2.1 | 10 mM Tris-HCl, 10 mM MgCl ₂ , 50 mM NaCl, 100 µg/ml BSA (pH 7.9 @ 25°C). |
| NEBuffer 3.1 | 50 mM Tris-HCl, 10 mM MgCl ₂ , 100 mM NaCl, 100 µg/ml BSA (pH 7.9 @ 25°C). |
| CutSmart | 20 mM Tris-acetate, 10 mM magnesium acetate, 50 mM potassium acetate, 100 µg/ml BSA (pH 7.9 @ 25°C). |
| Diluent A | 50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA (pH 7.4 @ 25°C). |
| Diluent B | 300 mM NaCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 500 µg/ml BSA, 50% glycerol (pH 7.4 @ 25°C). |
| Diluent C | 50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 0.15% Triton X-100, 200 µg/ml BSA 50% glycerol (pH 7.4 @ 25°C). |

| | ENZYME | SUPPLIED NEBUFFER | % ACTIVITY IN NEBUFFERS | | | | INCUB. TEMP. (°C) | INACTIV. TEMP. (°C) | DILUENT | UNIT SUBSTRATE | METHYLATION SENSITIVITY | NOTE |
|--|---------|-------------------|-------------------------|------|------|----------|-------------------|---------------------|---------|----------------|-------------------------|---------|
| | | | 1.1 | 2.1 | 3.1 | CUTSMART | | | | | | |
| | AatII | CutSmart | < 10 | 50* | 50 | 100 | 37° | 80° | B | Lambda | | |
| | AbaSI | CutSmart | 25 | 50 | 50 | 100 | 25° | 65° | C | T4 wt Phage | | e |
| | AccI | CutSmart | 50 | 50 | 10 | 100 | 37° | 80° | A | Lambda | | |
| | Acc65I | 3.1 | 10 | 75* | 100 | 25 | 37° | 65° | A | pBC4 | | |
| | AcII | CutSmart | < 10 | 25 | 100 | 100 | 37° | 65° | A | Lambda | | |
| | AcII | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | No | B | Lambda | | |
| | AcuI | CutSmart + SAM | 50 | 100 | 50 | 100 | 37° | 65° | B | Lambda | | 1, b, d |
| | AfeI | CutSmart | 25 | 100 | 25 | 100 | 37° | 65° | B | pXba | | |
| | AflII | CutSmart | 50 | 100 | 10 | 100 | 37° | 65° | A | phiX174 | | |
| | AflIII | 3.1 | 10 | 50 | 100 | 50 | 37° | 80° | B | Lambda | | |
| | AgeI | 1.1 | 100 | 75 | 25 | 75 | 37° | 65° | C | Lambda | | |
| | AgeI-HF | CutSmart | 100 | 50 | 10 | 100 | 37° | 65° | A | Lambda | | |
| | AhdI | CutSmart | 25 | 25 | 10 | 100 | 37° | 65° | A | Lambda | | a |
| | AleI | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | 80° | B | Lambda | | |
| | AluI | CutSmart | 25 | 100 | 50 | 100 | 37° | 80° | B | Lambda | | b |
| | AlwI | CutSmart | 50 | 50 | 10 | 100 | 37° | No | A | Lambda dam- | | 1, b, d |
| | AlwNI | CutSmart | 10 | 100 | 50 | 100 | 37° | 80° | A | Lambda | | |
| | ApaI | CutSmart | 25 | 25 | < 10 | 100 | 25° | 65° | A | pXba | | |



DNA PREPARATION – RESTRICTION ENZYME DIGESTION

| ENZYME | SUPPLIED NEBUFFER | % ACTIVITY IN NEBUFFERS | | | | | INCUB. TEMP. (°C) | INACTIV. TEMP. (°C) | DILUENT | UNIT SUBSTRATE | METHYLATION SENSITIVITY | NOTE |
|----------|-------------------|-------------------------|------|------|----------|-----|-------------------|---------------------|---------------------|----------------|-------------------------|------|
| | | 1.1 | 2.1 | 3.1 | CUTSMART | | | | | | | |
| ApaLI | CutSmart | 100 | 100 | 10 | 100 | 37° | No | A | Lambda HindIII | CpG | | |
| ApeKI | 3.1 | 25 | 50 | 100 | 10 | 75° | No | B | Lambda | CpG | | |
| ApoI | 3.1 | 10 | 75 | 100 | 75 | 50° | 80° | A | Lambda | | | |
| ApoI-HF | CutSmart | 10 | 100 | 10 | 100 | 37° | 80° | B | Lambda | | | |
| AscI | CutSmart | < 10 | 10 | 10 | 100 | 37° | 80° | A | Lambda | CpG | | |
| Asel | 3.1 | < 10 | 50* | 100 | 10 | 37° | 65° | B | Lambda | | 3 | |
| AsiSI | CutSmart | 50 | 100 | 100 | 100 | 37° | 80° | B | pXba (Xho digested) | CpG | 2, b | |
| AvaI | CutSmart | < 10 | 100 | 25 | 100 | 37° | 80° | A | Lambda | CpG | | |
| AvaII | CutSmart | 50 | 75 | 10 | 100 | 37° | 80° | A | Lambda | dcm CpG | | |
| AvrII | CutSmart | 100 | 50 | 50 | 100 | 37° | No | B | Lambda HindIII | | | |
| BaeI | CutSmart + SAM | 50 | 100 | 50 | 100 | 25° | 65° | A | pXba | CpG | e | |
| BaeGI | 3.1 | 75 | 75 | 100 | 25 | 37° | 80° | A | Lambda | | | |
| BamHI | 3.1 | 75* | 100* | 100 | 100* | 37° | No | A | Lambda | | 3 | |
| BamHI-HF | CutSmart | 100 | 50 | 10 | 100 | 37° | No | A | Lambda | | | |
| BanI | CutSmart | 10 | 25 | < 10 | 100 | 37° | 65° | A | Lambda | dcm CpG | 1 | |
| BanII | CutSmart | 100 | 100 | 50 | 100 | 37° | 80° | A | Lambda | | 2 | |
| BbsI | 2.1 | 100 | 100 | 25 | 75 | 37° | 65° | B | Lambda | | | |
| BbvI | CutSmart | 100 | 100 | 25 | | 37° | 65° | B | pBR322 | | 3 | |
| BbvCI | CutSmart | 10 | 100 | 50 | 100 | 37° | No | B | Lambda | CpG | 1, a | |
| BccI | CutSmart | 100 | 50 | 10 | 100 | 37° | 65° | A | pXba | | 3, b | |
| BceAI | 3.1 | 100* | 100* | 100 | 100* | 37° | 65° | A | pBR322 | CpG | 1 | |
| BcgI | 3.1 + SAM | 10 | 75* | 100 | 50* | 37° | 65° | A | Lambda | dam CpG | e | |
| BciVI | CutSmart | 100 | 25 | < 10 | 100 | 37° | 80° | C | Lambda | | b | |
| BclI | 3.1 | 50 | 100 | 100 | 75 | 50° | No | A | Lambda dam- | dam | | |
| BcoDI | CutSmart | 50 | 75 | 75 | 100 | 37° | No | B | Lambda | CpG | | |
| BfaI | CutSmart | < 10 | 10 | < 10 | 100 | 37° | 80° | B | Lambda | | 2, b | |
| BfuAI | 3.1 | < 10 | 25 | 100 | 10 | 50° | 65° | B | Lambda | CpG | 3 | |
| BfuCI | CutSmart | 100 | 50 | 25 | 100 | 37° | 80° | B | Lambda | CpG | | |
| BglI | 3.1 | 10 | 25 | 100 | 10 | 37° | 65° | B | Lambda | CpG | | |
| BglII | 3.1 | 10 | 10 | 100 | < 10 | 37° | No | A | Lambda | | | |
| BlpI | CutSmart | 50 | 100 | 10 | 100 | 37° | No | A | Lambda | | d | |
| BmgBI | 3.1 | < 10 | 10 | 100 | 10 | 37° | 65° | B | Lambda | CpG | 3, b, d | |
| BmrI | 2.1 | 75 | 100 | 75 | 100* | 37° | 65° | B | Lambda HindIII | | b | |
| BmtI | 3.1 | 100 | 100 | 100 | 100 | 37° | 65° | B | pXba | | 2 | |
| BmtI-HF | CutSmart | 50 | 100 | 10 | 100 | 37° | 65° | B | pXba | | | |
| BpmI | 3.1 | 75 | 100 | 100 | 100 | 37° | 65° | B | Lambda | | 2 | |
| Bpu10I | 3.1 | 10 | 25 | 100 | 25 | 37° | 80° | B | Lambda | | 3, b, d | |
| BpuEI | CutSmart + SAM | 50* | 100 | 50* | 100 | 37° | 65° | B | Lambda | | d | |
| BsaI | CutSmart | 75* | 75 | 100 | 100 | 37° | 65° | B | pXba | dcm CpG | 3 | |
| BsaI-HF | CutSmart | 50 | 100 | 25 | 100 | 37° | 65° | B | pXba | dcm CpG | | |
| BsaAI | CutSmart | 100 | 100 | 100 | 100 | 37° | No | C | Lambda | CpG | | |
| BsaBI | CutSmart | 50 | 100 | 75 | 100 | 60° | 80° | B | Lambda dam- | dam CpG | 2 | |
| BsaHI | CutSmart | 50 | 100 | 100 | 100 | 37° | 80° | A | Lambda | dcm CpG | | |
| BsaJI | CutSmart | 50 | 100 | 100 | 100 | 60° | 80° | A | Lambda | | | |
| BsaWI | CutSmart | 10 | 100 | 50 | 100 | 60° | 80° | A | Lambda | | | |
| BsaXI | CutSmart | 50* | 100* | 10 | 100 | 37° | No | B | Lambda | | e | |
| BseRI | CutSmart | 100* | 100 | 75 | 100 | 37° | 80° | A | Lambda | | d | |
| BseYI | 3.1 | 10 | 50 | 100 | 50 | 37° | 80° | B | Lambda | CpG | d | |
| BsgI | CutSmart + SAM | 25 | 50 | 25 | 100 | 37° | 65° | B | Lambda | | d | |
| BsiEI | CutSmart | 25 | 50 | < 10 | 100 | 60° | No | A | Lambda | CpG | | |
| BsiHKAI | CutSmart | 25 | 100 | 100 | 100 | 65° | No | B | Lambda | | | |
| BsiWI | 3.1 | 25 | 50* | 100 | 25 | 55° | 65° | B | phiX174 | CpG | | |
| BsII | CutSmart | 50 | 75 | 100 | 100 | 55° | No | A | Lambda | dcm CpG | b | |



| RESTR. ENZ. (RR) | ENZYME | SUPPLIED NEBUFFER | % ACTIVITY IN NEBUFFERS | | | | INCUB. TEMP. (°C) | INACTIV. TEMP. (°C) | DILUENT | UNIT SUBSTRATE | METHYLATION SENSITIVITY | NOTE |
|------------------|-----------------------|-------------------|-------------------------|------|------|----------|-------------------|---------------------|---------|----------------|-------------------------|---------|
| | | | 1.1 | 2.1 | 3.1 | CUTSMART | | | | | | |
| RR | BsmI | CutSmart | 25 | 100 | < 10 | 100 | 65° | 80° | A | Lambda | | |
| RR | BsmAI | CutSmart | 50 | 100 | 100 | 100 | 55° | No | B | Lambda | CpG | |
| RR | BsmBI | 3.1 | 10 | 50* | 100 | 25 | 55° | 80° | B | Lambda | CpG | |
| RR | BsmFI | CutSmart | 25 | 50 | 50 | 100 | 65° | 80° | A | pBR322 | dcm CpG | 1 |
| RR | BsoBI | CutSmart | 25 | 100 | 100 | 100 | 37° | 80° | A | Lambda | | |
| RR | Bsp1286I | CutSmart | 25 | 25 | 25 | 100 | 37° | 65° | A | Lambda | | 3 |
| RR | BspCNI | CutSmart + SAM | 100 | 75 | 10 | 100 | 25° | 80° | A | Lambda | | b |
| RR | BspDI | CutSmart | 25 | 75 | 50 | 100 | 37° | 80° | A | Lambda | dam CpG | |
| RR | BspEI | 3.1 | < 10 | 10 | 100 | < 10 | 37° | 80° | B | Lambda dam- | dam CpG | |
| RR | BspHI | CutSmart | < 10 | 50 | 25 | 100 | 37° | 80° | A | Lambda | dam | |
| RR | 2*site BspMI | 3.1 | 10 | 50* | 100 | 10 | 37° | 65° | B | Lambda | | |
| RR | BspQI | 3.1 | 100 | 100 | 100 | 100 | 50° | 80° | B | Lambda | | 3 |
| | BsrI | 3.1 | < 10 | 50 | 100 | 10 | 65° | 80° | B | phiX174 | | b |
| RR | BsrBI | CutSmart | 50 | 100 | 100 | 100 | 37° | 80° | A | Lambda | CpG | d |
| RR | BsrDI | 2.1 | 10 | 100 | 75 | 25 | 65° | 80° | A | Lambda | | 3, d |
| RR | e BsrFcoI | CutSmart | 25 | 25 | 0 | 100 | 37° | No | C | pBR322 | CpG | |
| RR | BsrGI | 2.1 | 25 | 100 | 100 | 25 | 37° | 80° | A | Lambda | | |
| RR | e BsrGI-HF | CutSmart | 10 | 100 | 100 | 100 | 37° | 80° | A | Lambda | | |
| RR | BssHII | CutSmart | 100 | 100 | 100 | 100 | 50° | 65° | B | Lambda | CpG | |
| RR | e BssS ^o I | CutSmart | 10 | 25 | < 10 | 100 | 37° | No | B | Lambda | | |
| RR | BstAPI | CutSmart | 50 | 100 | 25 | 100 | 60° | 80° | A | Lambda | CpG | b |
| RR | BstBI | CutSmart | 75 | 100 | 10 | 100 | 65° | No | A | Lambda | CpG | |
| RR | BstEII | 3.1 | 10 | 75* | 100 | 75* | 60° | No | A | Lambda | | 3 |
| RR | e BstEII-HF | CutSmart | < 10 | 10 | < 10 | 100 | 37° | No | A | Lambda | | |
| RR | BstNI | 3.1 | 10 | 100 | 100 | 75 | 60° | No | A | Lambda | | a |
| | BstUI | CutSmart | 50 | 100 | 25 | 100 | 60° | No | A | Lambda | CpG | b |
| RR | BstXI | 3.1 | < 10 | 50 | 100 | 25 | 37° | 80° | B | Lambda | dcm | 3 |
| RR | BstYI | 2.1 | 25 | 100 | 75 | 100 | 60° | No | A | Lambda | | |
| RR | BstZ17I | CutSmart | 75 | 100 | 100 | 100 | 37° | No | B | Lambda | CpG | 3, b |
| RR | e BstZ17I-HF | CutSmart | 100 | 100 | 10 | 100 | 37° | No | A | Lambda | CpG | |
| RR | Bsu36I | CutSmart | 25 | 100 | 100 | 100 | 37° | 80° | C | Lambda HindIII | | b |
| RR | BtgI | CutSmart | 50 | 100 | 100 | 100 | 37° | 80° | B | pBR322 | dcm | |
| RR | BtgZI | CutSmart | 10 | 25 | < 10 | 100 | 60° | 80° | A | Lambda | CpG | 3, b, d |
| RR | e Bts ^o I | CutSmart | 100 | 100 | 25 | 100 | 55° | No | A | Lambda | | |
| RR | e BtsIMutI | CutSmart | 100 | 50 | 10 | 100 | 55° | 80° | A | pUC19 | | b |
| RR | BtsCI | CutSmart | 10 | 100 | 25 | 100 | 50° | 80° | B | Lambda | | |
| | Cac8I | CutSmart | 50 | 75 | 100 | 100 | 37° | 65° | B | Lambda | CpG | b |
| RR | Clal | CutSmart | 10 | 50 | 50 | 100 | 37° | 65° | A | Lambda dam- | dam CpG | |
| RR | 2*site CspCI | CutSmart + SAM | 10 | 100 | 10 | 100 | 37° | 65° | A | Lambda | | e |
| RR | CviAII | CutSmart | 50 | 50 | 10 | 100 | 25° | 65° | C | Lambda | | |
| RR | CviKI-1 | CutSmart | 25 | 100 | 100 | 100 | 37° | No | A | pBR322 | dcm | 1, b |
| RR | CviQI | 3.1 | 75 | 100* | 100 | 75* | 25° | No | C | Lambda | | b |
| RR | Ddel | CutSmart | 75 | 100 | 100 | 100 | 37° | 65° | B | Lambda | | |
| RR | DpnI | CutSmart | 100 | 100 | 75 | 100 | 37° | 80° | B | pBR322 | CpG | b |
| RR | DpnII | U | 25 | 25 | 100* | 25 | 37° | 65° | B | Lambda dam- | dam | |
| RR | Dral | CutSmart | 75 | 75 | 50 | 100 | 37° | 65° | A | Lambda | | |
| RR | e DralII-HF | CutSmart | < 10 | 50 | 10 | 100 | 37° | No | B | Lambda | CpG | b |
| | DrdI | CutSmart | 25 | 50 | 10 | 100 | 37° | 65° | A | pUC19 | CpG | 3 |
| RR | EaeI | CutSmart | 10 | 50 | < 10 | 100 | 37° | 65° | A | Lambda | dcm CpG | b |
| RR | EagI | 3.1 | 10 | 25 | 100 | 10 | 37° | 65° | C | pXba | CpG | |
| RR | e EagI-HF | CutSmart | 25 | 100 | 100 | 100 | 37° | 65° | B | pXba | CpG | |
| RR | EarI | CutSmart | 50 | 10 | < 10 | 100 | 37° | 65° | B | Lambda | CpG | b, d |

1. Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

2. Star activity may result from extended digestion.

3. Star activity may result from a glycerol concentration of > 5%.

* May exhibit star activity in this buffer.



DNA PREPARATION – RESTRICTION ENZYME DIGESTION

| RR | ENZYME | SUPPLIED NEBUFFER | % ACTIVITY IN NEBUFFERS | | | | CUTSMART | INCUB. TEMP. (°C) | INACTIV. TEMP. (°C) | DILUENT | UNIT SUBSTRATE | METHYLATION SENSITIVITY | NOTE |
|----|------------|-------------------|-------------------------|------|------|-----|----------|-------------------|---------------------|------------------------|----------------|-------------------------|---------|
| | | | 1.1 | 2.1 | 3.1 | | | | | | | | |
| RR | EcoI | CutSmart | 100 | 50 | 50 | 100 | 37° | 65° | A | Lambda | CpG | 2 | |
| RR | Eco53kl | CutSmart | 100 | 100 | < 10 | 100 | 37° | 65° | A | pXba | CpG | 3, b | |
| RR | EcoNI | CutSmart | 50 | 100 | 75 | 100 | 37° | 65° | A | Lambda | | b | |
| RR | EcoO109I | CutSmart | 50 | 100 | 50 | 100 | 37° | 65° | A | Lambda HindIII | dcm | 3 | |
| RR | EcoP15I | 3.1 + ATP | 75 | 100 | 100 | 100 | 37° | 65° | A | pUC19 | | e | |
| RR | EcoRI | U | 25 | 100* | 50 | 50* | 37° | 65° | C | Lambda | CpG | | |
| RR | EcoRI-HF | CutSmart | 10 | 100 | < 10 | 100 | 37° | 65° | C | Lambda | CpG | | |
| RR | EcoRV | 3.1 | 10 | 50 | 100 | 10 | 37° | 80° | A | Lambda | CpG | | |
| RR | EcoRV-HF | CutSmart | 25 | 100 | 100 | 100 | 37° | 65° | B | Lambda | CpG | | |
| RR | FatI | 2.1 | 10 | 100 | 50 | 50 | 55° | 80° | A | pUC19 | | | |
| RR | FauI | CutSmart | 100 | 50 | 10 | 100 | 55° | 65° | A | Lambda | CpG | 3, b, d | |
| RR | Fnu4HI | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | No | A | Lambda | CpG | a | |
| RR | FokI | CutSmart | 100 | 100 | 75 | 100 | 37° | 65° | A | Lambda | dcm | CpG | 3, b, d |
| RR | FseI | CutSmart | 100 | 75 | < 10 | 100 | 37° | 65° | B | Adenovirus-2 | dcm | CpG | |
| RR | FspI | CutSmart | 10 | 100 | 10 | 100 | 37° | No | C | Lambda | CpG | b | |
| RR | FspEI | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | 80° | B | pBR322 | dcm | | 2, e |
| RR | HaeII | CutSmart | 25 | 100 | 10 | 100 | 37° | 80° | A | Lambda | CpG | | |
| RR | HaeIII | CutSmart | 50 | 100 | 25 | 100 | 37° | 80° | A | Lambda | | | |
| RR | HgaI | 1.1 | 100 | 100 | 25 | 100 | 37° | 65° | A | phiX174 | CpG | 1 | |
| RR | HhaI | CutSmart | 25 | 100 | 100 | 100 | 37° | 65° | A | Lambda | CpG | | |
| RR | HincII | 3.1 | 25 | 100 | 100 | 100 | 37° | 65° | B | Lambda | CpG | | |
| RR | HindIII | 2.1 | 25 | 100 | 50 | 50 | 37° | 80° | B | Lambda | | 2 | |
| RR | HindIII-HF | CutSmart | 10 | 100 | 10 | 100 | 37° | 80° | B | Lambda | | | |
| RR | Hinfi | CutSmart | 50 | 100 | 100 | 100 | 37° | 80° | A | Lambda | CpG | | |
| RR | HinPII | CutSmart | 100 | 100 | 100 | 100 | 37° | 65° | A | Lambda | CpG | | |
| RR | HpaI | CutSmart | < 10 | 75* | 25 | 100 | 37° | No | A | Lambda | CpG | 1 | |
| RR | HpaII | CutSmart | 100 | 50 | < 10 | 100 | 37° | 80° | A | Lambda | CpG | | |
| RR | HphI | CutSmart | 50 | 50 | < 10 | 100 | 37° | 65° | B | Lambda | dam | dcm | b, d |
| RR | Hpy99I | CutSmart | 50 | 10 | < 10 | 100 | 37° | 65° | A | Lambda | CpG | | |
| RR | Hpy166II | CutSmart | 100 | 100 | 50 | 100 | 37° | 65° | C | pBR322 | CpG | | |
| RR | Hpy188I | CutSmart | 25 | 100 | 50 | 100 | 37° | 65° | A | pBR322 | dam | | 1, b |
| RR | Hpy188III | CutSmart | 100 | 100 | 10 | 100 | 37° | 65° | B | pUC19 | dam | CpG | 3, b |
| RR | HpyAV | CutSmart | 100 | 100 | 25 | 100 | 37° | 65° | | Lambda | CpG | | 3, b, d |
| RR | HpyCH4III | CutSmart | 100 | 25 | < 10 | 100 | 37° | 65° | A | Lambda | | b | |
| RR | HpyCH4IV | CutSmart | 100 | 50 | 25 | 100 | 37° | 65° | A | pUC19 | CpG | | |
| RR | HpyCH4V | CutSmart | 50 | 50 | 25 | 100 | 37° | 65° | A | Lambda | | | |
| RR | I-CeuI | CutSmart | 10 | 10 | 10 | 100 | 37° | 65° | B | pBHS Scal-linearized | | | |
| RR | I-SceI | CutSmart | 10 | 50 | 25 | 100 | 37° | 65° | B | pGPS2 NottI-linearized | | | |
| RR | KasI | CutSmart | 50 | 100 | 50 | 100 | 37° | 65° | B | pBR322 | CpG | 3 | |
| RR | KpnI | 1.1 | 100 | 75 | < 10 | 100 | 37° | No | A | pXba | | 1 | |
| RR | KpnI-HF | CutSmart | 100 | 25 | < 10 | 100 | 37° | No | A | pXba | | | |
| RR | LpnPI | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | 65° | B | pBR322 | | 2, e | |
| RR | MboI | CutSmart | 75 | 100 | 100 | 100 | 37° | 65° | A | Lambda dam- | dam | CpG | |
| RR | MboII | CutSmart | 100* | 100 | 50 | 100 | 37° | 65° | C | Lambda dam- | dam | | b |
| RR | MfeI | CutSmart | 75 | 50 | 10 | 100 | 37° | No | A | Lambda | | 2 | |
| RR | MfeI-HF | CutSmart | 75 | 25 | < 10 | 100 | 37° | No | A | Lambda | | | |
| RR | MluI | 3.1 | 10 | 50 | 100 | 25 | 37° | 80° | A | Lambda | CpG | | |
| RR | MluI-HF | CutSmart | 25 | 100 | 100 | 100 | 37° | No | A | Lambda | CpG | | |
| RR | MluCI | CutSmart | 100 | 10 | 10 | 100 | 37° | No | A | Lambda | | | |
| RR | MlyI | CutSmart | 50 | 50 | 10 | 100 | 37° | 65° | A | Lambda | | b, d | |
| RR | MmeI | CutSmart + SAM | 50 | 100 | 50 | 100 | 37° | 65° | B | phiX174 | CpG | b, c | |
| RR | MnII | CutSmart | 75 | 100 | 50 | 100 | 37° | 65° | B | Lambda | | b | |

a. Ligation is less than 10%
b. Ligation is 25% – 75%

c. Recutting after ligation is <5%
d. Recutting after ligation is 50% – 75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.



| ENZYME | SUPPLIED NEBUFFER | % ACTIVITY IN NEBUFFERS | | | | INCUB. TEMP. (°C) | INACTIV. TEMP. (°C) | DILUENT | UNIT SUBSTRATE | METHYLATION SENSITIVITY | NOTE | |
|--------|-------------------|-------------------------|------|------|----------|-------------------|---------------------|---------|----------------|-------------------------|---------|---------|
| | | 1.1 | 2.1 | 3.1 | CUTSMART | | | | | | | |
| RR | MscI | CutSmart | 25 | 100 | 100 | 100 | 37° | 80° | B | Lambda | dcm | |
| RR | MseI | CutSmart | 75 | 100 | 75 | 100 | 37° | 65° | A | Lambda | | |
| RR | MspI | CutSmart | 50 | 50 | < 10 | 100 | 37° | 80° | A | Lambda | | |
| RR | MspI | CutSmart | 75 | 100 | 50 | 100 | 37° | No | A | Lambda | | |
| RR | MspA1I | CutSmart | 10 | 50 | 10 | 100 | 37° | 65° | B | Lambda | CpG | |
| RR | MspJI | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | 65° | B | pBR322 | | 2, e |
| RR | MwoI | CutSmart | < 10 | 100 | 100 | 100 | 60° | No | B | Lambda | CpG | |
| RR | 2*site NaeI | CutSmart | 25 | 25 | < 10 | 100 | 37° | No | A | pXba | CpG | b |
| RR | 2*site NarI | CutSmart | 100 | 100 | 10 | 100 | 37° | 65° | A | pXba | CpG | |
| RR | Nb.BbvCI | CutSmart | 25 | 100 | 100 | 100 | 37° | 80° | A | pUB | | e |
| RR | Nb.BsmI | 3.1 | < 10 | 50 | 100 | 10 | 65° | 80° | A | pBR322 | | e |
| RR | Nb.BsrDI | CutSmart | 25 | 100 | 100 | 100 | 65° | 80° | A | pUC19 | | e |
| RR | e Nb.BssSI | 3.1 | 10 | 100 | 100 | 25 | 37° | No | B | pUC19 | | |
| RR | Nb.BtsI | CutSmart | 75 | 100 | 75 | 100 | 37° | 80° | A | phiX174 | | e |
| RR | NciI | CutSmart | 100 | 25 | 10 | 100 | 37° | No | A | Lambda | CpG | b |
| RR | NcoI | 3.1 | 100 | 100 | 100 | 100 | 37° | 80° | A | Lambda | | |
| RR | e NcoI-HF | CutSmart | 50 | 100 | 10 | 100 | 37° | 80° | B | Lambda | | |
| RR | NdeI | CutSmart | 75 | 100 | 100 | 100 | 37° | 65° | A | Lambda | | |
| RR | 2*site NgoMIV | CutSmart | 100 | 50 | 10 | 100 | 37° | No | A | Adenovirus-2 | CpG | 1 |
| RR | NheI | 2.1 | 100 | 100 | 10 | 100 | 37° | 65° | C | Lambda HindIII | CpG | |
| RR | e NheI-HF | CutSmart | 100 | 25 | < 10 | 100 | 37° | 80° | C | Lambda HindIII | CpG | |
| RR | NlaIII | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | 65° | B | phiX174 | | |
| RR | NlaIV | CutSmart | 10 | 10 | 10 | 100 | 37° | 65° | B | pBR322 | dcm CpG | |
| RR | 2*site NmeAIII | CutSmart + SAM | 10 | 10 | < 10 | 100 | 37° | 65° | B | phiX174 | | c |
| RR | NotI | 3.1 | < 10 | 50 | 100 | 25 | 37° | 65° | C | pBC4 | CpG | |
| RR | e NotI-HF | CutSmart | 25 | 100 | 25 | 100 | 37° | 65° | A | pBC4 | CpG | |
| RR | NruI | 3.1 | < 10 | 10 | 100 | 10 | 37° | No | A | Lambda | dam CpG | b |
| RR | e NruI-HF | CutSmart | 0 | 25 | 50 | 100 | 37° | No | A | Lambda | dam CpG | |
| RR | NsiI | 3.1 | 10 | 75 | 100 | 25 | 37° | 65° | B | Lambda | | |
| RR | e NsiI-HF | CutSmart | < 10 | 20 | < 10 | 100 | 37° | 80° | B | Lambda | | |
| RR | NspI | CutSmart | 100 | 100 | < 10 | 100 | 37° | 65° | A | Lambda | | |
| RR | Nt.AIwI | CutSmart | 10 | 100 | 100 | 100 | 37° | 80° | A | pUC101 dam-dcm- | dam | e |
| RR | Nt.BbvCI | CutSmart | 50 | 100 | 10 | 100 | 37° | 80° | A | pUB | CpG | e |
| RR | Nt.BsmAI | CutSmart | 100 | 50 | 10 | 100 | 37 | 65° | A | pBR322 | CpG | e |
| RR | Nt.BspQI | 3.1 | < 10 | 25 | 100 | 10 | 50° | 80° | B | pUC19 | | e |
| RR | Nt.BstNBI | 3.1 | 0 | 10 | 100 | 10 | 55° | 80° | A | T7 | | |
| RR | Pacl | CutSmart | 100 | 75 | 10 | 100 | 37° | 65° | A | pNEB193 | | |
| RR | PaeR7I | CutSmart | 25 | 100 | 10 | 100 | 37° | No | A | Lambda HindIII | CpG | |
| RR | PciI | 3.1 | 50 | 75 | 100 | 50* | 37° | 80° | B | pXba | | |
| RR | PfiFI | CutSmart | 25 | 100 | 25 | 100 | 37° | 65° | A | pBC4 | | b |
| RR | PfiMI | 3.1 | 0 | 100 | 100 | 50 | 37° | 65° | A | Lambda | dcm | 3, b, d |
| RR | PI-PspI | U | 10 | 10 | 10 | 10 | 65° | No | B | pAKR XmnI | | |
| RR | PI-Scel | U | 10 | 10 | 10 | 10 | 37° | 65° | B | pBSvdeX XmnI | | |
| RR | 2*site PfiI | CutSmart | 25 | 50 | 25 | 100 | 37° | 65° | A | Lambda | CpG | b, d |
| RR | 2*site PfuTI | CutSmart | 100 | 25 | < 10 | 100 | 37° | 65° | A | pXba | CpG | b |
| RR | PmeI | CutSmart | < 10 | 50 | 10 | 100 | 37° | 65° | A | Lambda | CpG | |
| RR | PmlI | CutSmart | 100 | 50 | < 10 | 100 | 37° | 65° | A | Lambda HindIII | CpG | |
| RR | PpuMI | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | No | B | Lambda HindIII | dcm | |
| RR | PshAI | CutSmart | 25 | 50 | 10 | 100 | 37° | 65° | A | Lambda | CpG | |
| RR | PsiI | CutSmart | 10 | 100 | 10 | 100 | 37° | 65° | B | Lambda | | 3 |
| RR | PspGI | CutSmart | 25 | 100 | 50 | 100 | 75° | No | A | T7 | dcm | 3 |
| RR | PspOMI | CutSmart | 10 | 10 | < 10 | 100 | 37° | 65° | B | pXba | dcm CpG | |

1. Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

2. Star activity may result from extended digestion.
3. Star activity may result from a glycerol concentration of > 5%.

* May exhibit star activity in this buffer.



DNA PREPARATION – RESTRICTION ENZYME DIGESTION

| RR | ENZYME | SUPPLIED NEBUFFER | % ACTIVITY IN NEBUFFERS | | | | CUTSMART | INCUB. TEMP. (°C) | INACTIV. TEMP. (°C) | DILUENT | UNIT SUBSTRATE | METHYLATION SENSITIVITY | NOTE |
|----|----------|-------------------|-------------------------|------|------|------|----------|-------------------|---------------------|---------------------|----------------|-------------------------|------|
| | | | 1.1 | 2.1 | 3.1 | | | | | | | | |
| RR | PspXI | CutSmart | < 10 | 100 | 25 | 100 | 37° | No | B | Lambda HindIII | CpG | | |
| RR | PstI | 3.1 | 75 | 75 | 100 | 50* | 37° | 80° | C | Lambda | | | |
| RR | PstI-HF | CutSmart | 10 | 75 | 50 | 100 | 37° | No | C | Lambda | | | |
| RR | PvuI | 3.1 | < 10 | 25 | 100 | < 10 | 37° | No | B | pXba | CpG | | |
| RR | PvuI-HF | CutSmart | 25 | 100 | 100 | 100 | 37° | No | B | pXba | CpG | | |
| RR | PvuII | 3.1 | 50 | 100 | 100 | 100* | 37° | No | B | Lambda | | | |
| RR | PvuII-HF | CutSmart | < 10 | < 10 | < 10 | 100 | 37° | No | B | Lambda | | | |
| RR | RsaI | CutSmart | 25 | 50 | < 10 | 100 | 37° | No | A | Lambda | CpG | | |
| RR | RsrII | CutSmart | 25 | 75 | 10 | 100 | 37° | 65° | C | Lambda | CpG | | |
| RR | SacI | 1.1 | 100 | 50 | 10 | 100 | 37° | 65° | A | Lambda HindIII | | | |
| RR | SacI-HF | CutSmart | 10 | 50 | < 10 | 100 | 37° | 65° | A | Lambda HindIII | | | |
| RR | SacII | CutSmart | 10 | 100 | 10 | 100 | 37° | 65° | A | pXba | CpG | | |
| RR | SalI | 3.1 | < 10 | < 10 | 100 | < 10 | 37° | 65° | A | Lambda HindIII | CpG | | |
| RR | SalI-HF | CutSmart | 10 | 100 | 100 | 100 | 37° | 65° | A | Lambda HindIII | CpG | | |
| RR | SapI | CutSmart | 75 | 50 | < 10 | 100 | 37° | 65° | B | Lambda | | | |
| RR | Sau3AI | 1.1 | 100 | 50 | 10 | 100 | 37° | 65° | A | Lambda | CpG | b | |
| RR | Sau96I | CutSmart | 50 | 100 | 100 | 100 | 37° | 65° | A | Lambda | dcm CpG | | |
| RR | SbfI | CutSmart | 50 | 25 | < 10 | 100 | 37° | 80° | A | Lambda | | 3 | |
| RR | SbfI-HF | CutSmart | 50 | 25 | < 10 | 100 | 37° | 80° | B | Lambda | | | |
| RR | Scal-HF | CutSmart | 100 | 100 | 10 | 100 | 37° | 80° | B | Lambda | | | |
| RR | ScrFI | CutSmart | 100 | 100 | 100 | 100 | 37° | 65° | C | Lambda | dcm CpG | 2, a | |
| RR | SexAI | CutSmart | 100 | 75 | 50 | 100 | 37° | 65° | A | pBC4 dcm- | dcm | 3, b, d | |
| RR | SfiNI | 3.1 | < 10 | 75 | 100 | 25 | 37° | 65° | B | phiX174 | CpG | 3, b | |
| RR | SfiI | CutSmart | 75 | 50 | 25 | 100 | 37° | 65° | B | Lambda | | 3 | |
| RR | SfiI | CutSmart | 25 | 100 | 50 | 100 | 50° | No | C | Adenovirus-2 | dcm CpG | | |
| RR | SfoI | CutSmart | 50 | 100 | 100 | 100 | 37° | No | B | Lambda HindIII | dcm CpG | | |
| RR | SgrAI | CutSmart | 100 | 100 | 10 | 100 | 37° | 65° | A | Lambda | CpG | 1 | |
| RR | SmaI | CutSmart | < 10 | < 10 | < 10 | 100 | 25° | 65° | B | Lambda HindIII | CpG | b | |
| RR | SmlI | CutSmart | 25 | 75 | 25 | 100 | 55° | No | A | Lambda | | b | |
| RR | SnaBI | CutSmart | 50 | 50 | 10 | 100 | 37° | 80° | A | T7 | CpG | 1 | |
| RR | SpeI | CutSmart | 75 | 100 | 25 | 100 | 37° | 80° | C | Adenovirus-2 | | | |
| RR | SpeI-HF | CutSmart | 25 | 50 | 10 | 100 | 37° | 80° | C | pXba | | | |
| RR | SphI | 2.1 | 100 | 100 | 50 | 100 | 37° | 65° | B | Lambda | | 2 | |
| RR | SphI-HF | CutSmart | 50 | 25 | 10 | 100 | 37° | 65° | B | Lambda | | | |
| RR | SrfI | CutSmart | 10 | 50 | 0 | 100 | 37° | 65° | B | pNEB193-SrFI | CpG | | |
| RR | SspI | U | 50 | 100 | 50 | 50 | 37° | 65° | C | Lambda | | | |
| RR | SspI-HF | CutSmart | 25 | 100 | < 10 | 100 | 37° | 65° | B | Lambda | | | |
| RR | StuI | CutSmart | 50 | 100 | 50 | 100 | 37° | No | A | Lambda | dcm | | |
| RR | StyD4I | CutSmart | 10 | 100 | 100 | 100 | 37° | 65° | B | Lambda | dcm CpG | | |
| RR | StyI | 3.1 | 10 | 25 | 100 | 10 | 37° | 65° | A | Lambda | | b | |
| RR | StyI-HF | CutSmart | 25 | 100 | 25 | 100 | 37° | 65° | A | Lambda | | | |
| RR | Swal | 3.1 | 10 | 10 | 100 | 10 | 25° | 65° | B | M13mp19 RFI | | b, d | |
| RR | TaqI | CutSmart | 50 | 75 | 100 | 100 | 65° | 80° | B | Lambda | dam | | |
| RR | TfiI | CutSmart | 50 | 100 | 100 | 100 | 65° | No | C | Lambda | CpG | | |
| RR | TseI | CutSmart | 75 | 100 | 100 | 100 | 65° | No | B | Lambda | CpG | 3 | |
| RR | Tsp45I | CutSmart | 100 | 50 | < 10 | 100 | 65° | No | A | Lambda | | | |
| RR | TspMI | CutSmart | 50* | 75* | 50* | 100 | 75° | No | B | pUCAdeno | CpG | d | |
| RR | TspRI | CutSmart | 25 | 50 | 25 | 100 | 65° | No | B | Lambda | | | |
| RR | Tth111I | CutSmart | 25 | 100 | 25 | 100 | 65° | No | B | pBC4 | | b | |
| RR | XbaI | CutSmart | < 10 | 100 | 75 | 100 | 37° | 65° | A | Lambda HindIII dam- | dam | | |
| RR | XcmI | 2.1 | 10 | 100 | 25 | 100 | 37° | 65° | C | Lambda | | 2 | |
| RR | XhoI | CutSmart | 75 | 100 | 100 | 100 | 37° | 65° | A | Lambda HindIII | CpG | b | |
| RR | XmaI | CutSmart | 25 | 50 | < 10 | 100 | 37° | 65° | A | Adenovirus-2 | CpG | 3 | |
| RR | XmnI | CutSmart | 50 | 75 | < 10 | 100 | 37° | 65° | A | Lambda | | b | |
| RR | ZraI | CutSmart | 100 | 25 | 10 | 100 | 37° | 80° | B | Lambda | CpG | | |

a. Ligation is less than 10%
b. Ligation is 25% – 75%

c. Recutting after ligation is <5%
d. Recutting after ligation is 50% – 75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.



PCR/Amplification

Amplification can be performed to generate a blunt insert, or to have a 1-base overhang, depending on the polymerase used. Additionally, primers can be used to incorporate RE recognition sites. After amplification, the insert can be used directly or cloned into a holding vector, or RE digestion can be performed to generate cohesive ends. Amplification is often the first step for PCR cloning, seamless cloning, ligation independent cloning and recombinational cloning.

Protocol: High-Fidelity PCR with Q5

| | 25 μ l REACTION | 50 μ l REACTION | FINAL CONCENTRATION |
|--|---------------------|---------------------|---------------------------|
| 5X Q5 Reaction Buffer* | 5 μ l | 10 μ l | 1X |
| 10 mM dNTPs | 0.5 μ l | 1 μ l | 200 μ M |
| 10 μ M primers (forward and reverse) | 1.25 μ l | 2.5 μ l | 0.5 μ M |
| Template DNA | variable | variable | < 1 μ g |
| Nuclease-free water | to 25 μ l | to 50 μ l | |
| Q5 High-Fidelity DNA Polymerase** | 0.25 μ l | 0.5 μ l | 0.02 units/50 μ l rxn |

* Q5 High GC Enhancer can be used for difficult amplicons.

** For amplicons > 6 kb, up to 2 units/50 μ l rxn can be added.

Protocol: Routine PCR with OneTaq[®]

| | 25 μ l REACTION | 50 μ l REACTION | FINAL CONCENTRATION |
|--|---------------------|---------------------|---------------------------|
| OneTaq Standard 5X Reaction Buffer* | 5 μ l | 10 μ l | 1X |
| 10 mM dNTPs | 0.5 μ l | 1 μ l | 200 μ M |
| 10 μ M primers (forward and reverse) | 0.5 μ l | 1 μ l | 0.2 μ M |
| Template DNA | variable | variable | < 1 μ g |
| Nuclease-free water | to 25 μ l | to 50 μ l | |
| OneTaq DNA Polymerase** | 0.125 μ l | 0.25 μ l | 1.25 units/50 μ l rxn |

* If reaction buffer is 5X, volume should be doubled.

** Amount of polymerase added will depend on polymerase used.

| | CYCLES | TEMP. | TIME |
|-----------------------|--------|----------|----------------------|
| Initial denaturation: | 1 | 98°C | 30 seconds |
| Denaturation | 30 | 98°C | 5–10 seconds |
| Annealing | | 50–72°C* | 10–30 seconds |
| Extension | | 72°C | 20–30 seconds per kb |
| Final extension: | 1 | 72°C | 2 minutes |
| Hold: | 1 | 4–10°C | |

* T_m values should be determined using the NEB T_m calculator.

Please note that Q5 and Phusion[®] annealing temperature recommendations are unique.

| | CYCLES | TEMP. | TIME |
|-----------------------|--------|----------|-----------------|
| Initial denaturation: | 1 | 94°C | 30 seconds |
| Denaturation | 30 | 94°C | 15–30 seconds |
| Annealing | | 45–68°C* | 15–60 seconds |
| Extension | | 68°C | 1 minute per kb |
| Final extension: | 1 | 68°C | 5 minutes |
| Hold: | 1 | 4–10°C | |

* T_m values should be determined using the NEB T_m calculator.

TIPS FOR OPTIMIZATION

When switching from a *Taq* product to a high-fidelity polymerase, remember to use:

- Higher annealing temps
- Higher denaturation temps – particularly beneficial for difficult templates
- Higher primer concentrations
- Shorter cycling protocols

DNA TEMPLATE

- Use high-quality, purified DNA templates whenever possible. Refer to specific product information for amplification from unpurified DNA (i.e., colony or direct PCR).
- For low-complexity templates (i.e., plasmid, lambda, BAC DNA), use 1 μ g – 10 ng of DNA per 50 μ l reaction
- For higher complexity templates (i.e., genomic DNA), use 1 ng – 1 μ g of DNA per 50 μ l reaction
- Higher DNA concentrations tend to decrease amplicon specificity, particularly for high numbers of cycles

PRIMERS

- Primers should typically be 20–40 nucleotides in length, with 40–60% GC content
- Primer T_m values should be determined with NEB's T_m Calculator
- Primer pairs should have T_m values that are within 5°C

- Avoid secondary structure (i.e., hairpins) within each primer and potential dimerization between the primers
- Higher than recommended primer concentrations may decrease specificity
- When engineering restriction sites onto the end of primers, 6 nucleotides should be added 5' to the site

ENZYME CONCENTRATION

- Optimal concentration is specific to each polymerase
- Master mix formulations already contain optimal enzyme concentrations for most applications

MAGNESIUM CONCENTRATION

- Most PCR buffers provided by NEB already contain sufficient levels of Mg²⁺ at 1X concentrations
- Excess Mg²⁺ may lead to spurious amplification; insufficient Mg²⁺ concentrations may cause reaction failure

DEOXYNUCLEOTIDES

- Ideal dNTP concentration is typically 200 μ M each
- The presence of uracil in the primer, template, or deoxynucleotide mix will cause reaction failure when using archaeal PCR polymerases. Use OneTaq or Taq DNA Polymerases for these applications.

STARTING REACTIONS

- Unless using a hot start enzyme, assemble all reaction components on ice
- Add the polymerase last, whenever possible

- Transfer reactions to a thermocycler that has been pre-heated to the denaturation temperature. Pre-heating the thermocycler is not necessary when using a hot start enzyme (e.g., Q5 Hot Start or OneTaq Hot Start).

DENATURATION

- Avoid longer or higher temperature incubations unless required due to high GC content of the template
- NEB's aptamer-based hot start enzymes do not require additional denaturation steps to activate the enzymes

ANNEALING

- Primer T_m values should be determined using the NEB T_m Calculator
- Non-specific product formation can often be avoided by optimizing the annealing temperature or by switching to a hot start enzyme (e.g., Q5 Hot Start High-Fidelity DNA Polymerase or OneTaq Hot Start DNA Polymerase)

EXTENSION

- Extension rates are specific to each PCR polymerase. In general, extension rates range from 15–60 s/kb.
- Longer than recommended extension times can result in higher error rates, spurious banding patterns and/or reduction of amplicon yields



PCR Polymerase Selection Chart for Cloning

For almost 40 years, New England Biolabs, Inc. has been a world leader in the discovery and production of reagents for the life science industry. NEB offers a wide range of DNA polymerases, and through our commitment to research, ensures the development of innovative and high quality tools for PCR and related applications. The following table simplifies the selection of a polymerase that best suits your cloning experiment.

| | STANDARD PCR | | HIGH-FIDELITY PCR | | SPECIALTY PCR |
|---------------------------|--------------------------------|------------------------|--------------------|---|--|
| | OneTaq/ OneTaq Hot Start | Taq / Hot Start Taq | Highest Fidelity | | Long Amplicons |
| | | | Q5/Q5 Hot Start | Phusion ^{®(1)} / Phusion ⁽¹⁾ Flex | LongAmp [®] / LongAmp Hot Start Taq |
| PROPERTIES | | | | | |
| Fidelity vs. Taq | 2X | 1X | > 100X | > 50X | 2X |
| Amplicon Size | < 6 kb | ≤ 5 kb | ≤ 20 kb | ≤ 20 kb | ≤ 30 kb |
| Extension Time | 1 kb/min | 1 kb/min | 6 kb/min | 4 kb/min | 1.2 kb/min |
| Resulting Ends | 3' A/Blunt | 3' A | Blunt | Blunt | 3' A/Blunt |
| 3' → 5' exo | Yes | No | Yes | Yes | Yes |
| 5' → 3' exo | Yes | Yes | No | No | Yes |
| Units/50 µl Reaction | 1.25 | 1.25 | 1.0 | 1.0 | 5.0 |
| Annealing Temperature | Tm ⁻⁵ | Tm ⁻⁵ | Tm ⁺³ | Tm ⁺³ | Tm ⁻⁵ |
| APPLICATIONS | | | | | |
| Routine PCR | ★ | ● | ● | ● | ● |
| Colony PCR | ★ | ● | ● | ● | ● |
| Enhanced Fidelity | ● | ● | ★ | ● | ● |
| High Fidelity | ● | ● | ★ | ● | ● |
| High Yield | ★ | ● | ★ | ● | ● |
| Fast | ● | ● | ★ | ● | ● |
| Long Amplicon | ● | ● | ★ | ● | ★ |
| GC-rich Targets | ★ | ● | ★ | ● | ● |
| AT-rich Targets | ★ | ● | ★ | ● | ● |
| High Throughput | ● | ● | ● | ● | ● |
| Multiplex PCR | ● | ★ ⁽²⁾ | ● | ● | ● |
| Site-directed Mutagenesis | ● | ● | ★ | ● | ● |
| FORMATS | | | | | |
| Hot Start Available | ● | ● | ● | ● | ● |
| Kit | ● | ● | ● | ● | ● |
| Master Mix Available | ● | ● | ● | ● | ● |
| Direct Gel Loading | ● | ● ⁽³⁾ | ● | ● | ● |

(1) Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific.

(2) Use Multiplex PCR 5X Master Mix.

(3) Use Quick-Load 2X Taq Master Mix

★ indicates recommended choice for application
ND indicates not determined

GETTING STARTED

- When choosing a polymerase for PCR, we recommend starting with OneTaq or Q5 DNA Polymerases (highlighted to the left in orange). Both offer robust amplification and can be used on a wide range of templates (routine, AT- and GC-rich). Q5 provides the benefit of maximum fidelity, and is also available in a formulation specifically optimized for next generation sequencing.



Common DNA End Modifications

Modification of the termini of double-stranded DNA is often necessary to prepare the molecule for cloning. DNA ligases require a 5' monophosphate on the donor end, and the acceptor end requires a 3' hydroxyl group. Additionally, the sequences to be joined need to be compatible, either a blunt end being joined to another blunt end, or a cohesive end with a complementary overhang to another cohesive end. End modifications are performed to improve the efficiency of the cloning process, and ensure the ends to be joined are compatible.

Phosphorylation

Vectors and inserts digested by restriction enzymes contain the necessary terminal modifications (5' phosphate and 3' hydroxyl), while ends created by PCR may not. Typical amplification by PCR does not use phosphorylated primers. In this case, the 5' ends of the amplicon are non-phosphorylated and need to be treated by a kinase, such as T4 Polynucleotide Kinase, to introduce the 5' phosphate. Alternatively, primers for PCR can be ordered with 5' phosphate to avoid the need to separately phosphorylate the PCR product with a kinase.

Protocol: Phosphorylation with T4 Polynucleotide Kinase

| | STANDARD PROTOCOL |
|---------------------------------------|---------------------------------|
| DNA | 1–2 µg |
| 10X Polynucleotide Kinase Buffer | 5 µl |
| 10 mM Adenosine 5'-Triphosphate (ATP) | 5 µl (1 mM final concentration) |
| T4 Polynucleotide Kinase (PNK) | 1 µl (10 units) |
| Nuclease-free water | to 50 µl |
| Incubation | 37°C, 30 minutes |

TIPS FOR OPTIMIZATION

ENZYME

- T4 Polynucleotide Kinase and T4 DNA Ligase can be used together in the T4 DNA Ligase Buffer.
- T4 Polynucleotide Kinase is inhibited by high levels of salt (50% inhibition by 150 mM NaCl), phosphate (50% inhibition by 7 mM phosphate) and ammonium ions (75% inhibited by 7 mM $(\text{NH}_4)_2\text{SO}_4$).
- If using T4 Polynucleotide Kinase and working with 5'-recessed ends, heat the reaction mixture for 10 min at 70°C, chill rapidly on ice before adding the ATP (or Ligase Buffer containing ATP) and enzyme, then incubate at 37°C.

ADDITIVES

- The addition of PEG 8000 (up to 5%) can improve results.



Dephosphorylation

Dephosphorylation is a common step in traditional cloning to ensure the vector does not re-circularize during ligation. If a vector is linearized by a single restriction enzyme or has been cut with two enzymes with compatible ends, use of a phosphatase to remove the 5' phosphate reduces the occurrence of vector re-closure by intramolecular ligation and thereby reduces the background during subsequent transformation. If the vector is dephosphorylated, it is essential to ensure the insert contains a 5' phosphate to allow ligation to proceed. Each double-strand break requires that one intact phosphodiester bond be created before transformation (and *in vivo* repair).

Protocol: Dephosphorylation using Quick Dephosphorylation Kit

| | STANDARD PROTOCOL |
|---------------------|---------------------|
| DNA | 1 pmol of ends |
| 10X CutSmart Buffer | 2 µl |
| Quick CIP | 1 µl |
| Nuclease-free water | to 20 µl |
| Incubation | 37°C for 10 minutes |
| Heat Inactivation | 80°C for 2 minutes |

Phosphatase Selection Chart

| | Recombinant Shrimp Alkaline Phosphatase (rSAP) | Antarctic Phosphatase | Alkaline Phosphatase Calf Intestinal (CIP) | Quick Dephosphorylation Kit |
|-------------------------------|--|-----------------------|--|-----------------------------|
| FEATURES | | | | |
| 100% heat inactivation | 5 minutes/65°C | 2 minutes/80°C | No | 2 minutes/80°C |
| High specific activity | • | | • | • |
| Improved stability | • | | | • |
| Works directly in NEB buffers | • | • | • | • |
| Requires additive | | • (Zn ²⁺) | | |
| Quick Protocol | | | | • |

TIPS FOR OPTIMIZATION

ENZYME

- When trying to dephosphorylate a fragment following restriction enzyme digest, if the restriction enzyme(s) used are heat-inactivable, then a clean-up step prior to the addition of the phosphatase is not needed. Alternatively, if the restriction enzyme(s) used are not heat inactivable, a DNA clean step is recommended prior to the dephosphorylation step. For this we recommend the Monarch PCR & DNA Cleanup Kit.
- When working with the Quick Dephosphorylation Kit, rSAP or AP, which are heat-inactivable enzymes, a DNA clean-up step after dephosphorylation is not necessary prior to the ligation step. However, when using CIP, a clean-up step (e.g., Monarch PCR & DNA Cleanup Kit) prior to ligation is necessary.


ADDITIVES

- AP requires the presence of Zn²⁺ in the reaction, so don't forget to supplement the reaction with 1X Antarctic Phosphatase Reaction Buffer when using other NEBuffers.

Migrate to Monarch[®]

Monarch PCR & DNA Cleanup Kit

- Purify DNA from a variety of enzymatic reactions
- Elute in as little as 6 µl
- Prevent buffer retention and salt carry-over with optimized column design
- Save time with fast, user-friendly protocol
- Uses less plastic & recyclable packaging



The improved design of the columns supplied with the Monarch Gel Extraction and PCR & DNA Cleanup Kits enables elution in as little as 6 µl, and eliminates buffer retention



Blunting/End-repair

Blunting is a process by which the single-stranded overhang created by a restriction digest is either “filled in”, by adding nucleotides on the complementary strand using the overhang as a template for polymerization, or by “chewing back” the overhang, using an exonuclease activity. Vectors and inserts are often “blunted” to allow non-compatible ends to be joined. Sequence information is lost or distorted by doing this and a detailed understanding of the modification should be considered before performing this procedure. Often, as long as the sequence being altered is not part of the translated region or a critical regulatory element, the consequence of creating blunt ends is negligible. Blunting a region of translated coding sequence, however, usually creates a shift in the reading frame. DNA polymerases, such as the Klenow Fragment of DNA Polymerase I and T4 DNA Polymerase, included in our Quick Blunting Kit, are often used to fill in (5′→3′) and chew back (3′→5′). Removal of a 5′ overhang can be accomplished with a nuclease, such as Mung Bean Nuclease.

Protocol: Blunting using the Quick Blunting Kit

| | STANDARD PROTOCOL |
|----------------------------|--|
| DNA | up to 5 µg |
| 10X Blunting Buffer | 2.5 µl |
| 1 mM dNTP Mix | 2.5 µl |
| Blunt Enzyme Mix | 1 µl |
| Nuclease-free water | to 25 µl |
| Incubation | room temperature; 15 min for RE-digested DNA; 30 min for sheared/nebulized DNA or PCR products* |
| Heat Inactivation | 70°C, 10 minutes |

* PCR generated DNA must be purified before blunting by using a purification kit, phenol extraction/ethanol precipitation, or gel extraction.

Blunting Selection Chart

| | T4 DNA Polymerase* | DNA Polymerase I, Large (Klenow) Fragment | Quick Blunting Kit | Mung Bean Nuclease |
|-------------------------|--------------------|---|--------------------|--------------------|
| APPLICATION | | | | |
| Fill in of 5′ overhangs | • | • | • | |
| Removal of 3′ overhangs | • | • | • | • |
| Removal of 5′ overhangs | | | | • |

* T4 DNA Polymerase has a strong 3′→5′ exo activity.

TIPS FOR OPTIMIZATION

ENZYME

- Make sure that you choose the correct enzyme to blunt your fragment. The Quick Blunting Kit, T4 DNA Polymerase and DNA Polymerase I, Large (Klenow) Fragment will fill 5′ overhangs and degrade 3′ overhangs. Mung Bean Nuclease degrades 5′ overhangs.
- T4 DNA Polymerase and DNA Polymerase I, Large (Klenow) Fragment are active in all NEBuffers. Please remember to add dNTPs.

CLEAN-UP

- When trying to blunt a fragment after a restriction enzyme digestion, if the restriction enzyme(s) used are heat inactivable, then a clean-up step prior to blunting is not needed. Alternatively, if the restriction enzyme(s) used are not heat inactivable, a DNA clean-up step is recommended prior to blunting.
- When trying to blunt a fragment amplified by PCR, a DNA clean-up step (e.g., Monarch PCR & DNA Cleanup Kit) is necessary prior to the blunting step to remove the nucleotides and polymerase.
- When trying to dephosphorylate a fragment after the blunting step, you will need to add a DNA clean-up step (e.g., Monarch PCR & DNA Cleanup Kit) after the blunting and before the addition of the phosphatase.

TEMPERATURE

- When trying to blunt a fragment with Mung Bean Nuclease, the recommended temperature of incubation is room temperature, since higher temperatures may cause sufficient breathing of the dsDNA ends that the enzyme may degrade some of the dsDNA sequence. The number of units to be used and time of incubation may be determined empirically to obtain best results.

HEAT INACTIVATION

- Mung Bean nuclease reactions should not be heat inactivated. Although Mung Bean Nuclease can be inactivated by heat, this is not recommended because the DNA begins to “breathe” before the Mung Bean Nuclease is inactivated and undesirable degradation occurs at breathing sections. Purify DNA by phenol/chloroform extraction and ethanol precipitation or spin column purification.



A-tailing

Tailing is an enzymatic method to add a non-templated nucleotide to the 3' end of a blunt, double-stranded DNA molecule. Tailing is typically done to prepare a T-vector for use in TA cloning or to A-tail a PCR product produced by a high-fidelity polymerase (not *Taq* DNA Polymerase) for use in TA cloning. TA cloning is a rapid method of cloning PCR products that utilizes stabilization of the single-base extension (adenosine) produced by *Taq* DNA Polymerase by the complementary T (thymidine) of the T-vector prior to ligation and transformation. This technique does not utilize restriction enzymes and PCR products can be used directly without modification. Additionally, PCR primers do not need to be designed with restriction sites, making the process less complicated. One drawback is that the method is non-directional; the insert can go into the vector in both orientations.

TIPS FOR OPTIMIZATION

- If the fragment to be tailed has been amplified with a high-fidelity polymerase, the DNA needs to be purified prior to the tailing reaction. For this we recommend the Monarch PCR & DNA Cleanup Kit. Otherwise, any high-fidelity polymerase present in the reaction will be able to remove any non-templated nucleotides added to the end of the fragments.

Protocol: A-tailing with Klenow Fragment (3' → 5' exo-)

| | STANDARD PROTOCOL |
|--------------------------------|-----------------------|
| Purified, blunt DNA | 1–5 µg* |
| NEBuffer 2 (10X) | 5 µl |
| dATP (1 mM) | 0.5 µl (0.1 mM final) |
| Klenow Fragment (3' → 5' exo-) | 3 µl |
| H ₂ O | to 50 µl |
| Incubation | 37°C, 30 minutes |

* If starting with blunt-ended DNA that has been prepared by PCR or end polishing, DNA must be purified to remove the blunting enzymes.

A-tailing Selection Chart

| | Klenow Fragment (3' → 5' exo-) | <i>Taq</i> DNA Polymerase |
|----------------------|-----------------------------------|---------------------------|
| FEATURES | | |
| Reaction temperature | 37°C | 75°C |
| Heat inactivated | 75°C, 20 minutes | No |
| Nucleotide cofactor | dATP | dATP |

Activity of DNA Modifying Enzymes in CutSmart Buffer

A selection of DNA modifying enzymes were assayed in CutSmart Buffer, in lieu of their supplied buffers. Functional activity was compared to the activity in its supplied buffer, plus required supplements. Reactions were set up according to the recommended reaction conditions, with CutSmart Buffer replacing the supplied buffer.

| ENZYME | ACTIVITY IN CUTSMART | REQUIRED SUPPLEMENTS |
|---|----------------------|---------------------------|
| Alkaline Phosphatase (CIP) | +++ | |
| Antarctic Phosphatase | +++ | Requires Zn ²⁺ |
| <i>Bst</i> DNA Polymerase | +++ | |
| CpG Methyltransferase (M. SssI) | +++ | Requires SAM |
| DNA Polymerase I | +++ | |
| DNA Polymerase I, Large (Klenow) Fragment | +++ | |
| DNA Polymerase Klenow Exo | +++ | |
| DNase I (RNase free) | +++ | Requires Ca ²⁺ |
| <i>E. coli</i> DNA Ligase | +++ | Requires NAD |
| Endonuclease III (Nth), recombinant | +++ | |
| Endonuclease VIII | +++ | |
| Exonuclease III | +++ | |
| GpC Methyltransferase (M. CviPI) | + | Requires DTT |
| McrBC | +++ | |

+++ full functional activity ++ 50–100% functional activity + 0–50% functional activity

| ENZYME | ACTIVITY IN CUTSMART | REQUIRED SUPPLEMENTS |
|---|----------------------|---------------------------|
| Micrococcal Nuclease | ++ | Requires Ca ²⁺ |
| phi29 DNA Polymerase | +++ | |
| RecJ | +++ | |
| Shrimp Alkaline Phosphatase (rSAP) | +++ | |
| T3 DNA Ligase | +++ | Requires ATP + PEG |
| T4 DNA Ligase | +++ | Requires ATP |
| T4 DNA Polymerase | +++ | |
| T4 Phage β-glucosyltransferase (T4-BGT) | +++ | |
| T4 Polynucleotide Kinase | +++ | Requires ATP + DTT |
| T4 PNK (3' phosphatase minus) | +++ | Requires ATP + DTT |
| T7 DNA Ligase | +++ | Requires ATP + PEG |
| T7 DNA Polymerase (unmodified) | +++ | |
| T7 Exonuclease | +++ | |
| USER™ Enzyme, recombinant | +++ | |



Vector and Insert Joining

DNA Ligation

Ligation of DNA is a critical step in many modern molecular biology workflows. The sealing of nicks between adjacent residues of a single-strand break on a double-strand substrate and the joining of double-strand breaks are enzymatically catalyzed by DNA ligases. The formation of a phosphodiester bond between the 3' hydroxyl and 5' phosphate of adjacent DNA residues proceeds in three steps: Initially, the ligase is self-adenylated by reaction with free ATP. Next, the adenylyl group is transferred to the 5' phosphorylated end of the “donor” strand. Lastly, the formation of the phosphodiester bond proceeds after reaction of the adenylylated donor end with the adjacent 3' hydroxyl acceptor and the release of AMP. In living organisms, DNA ligases are essential enzymes with critical roles in DNA replication and repair. In the lab, DNA ligation is performed for both cloning and non-cloning applications.

Molecular cloning is a method to prepare a recombinant DNA molecule, an extra-chromosomal circular DNA that can replicate autonomously within a microbial host. DNA ligation is commonly used in molecular cloning projects to physically join a DNA vector to a sequence of interest (“insert”). The ends of the DNA fragments can be blunt or cohesive and at least one must contain a monophosphate group on its 5' ends. Following the mechanism described above, the covalent bonds are formed and a closed circular molecule is created that is capable of transforming a host bacterial strain. The recombinant plasmid maintained in the host is then available for amplification prior to downstream applications such as DNA sequencing, protein expression, or gene expression/functional analysis.

Protocol: Ligation

| | Quick Ligation Kit | T4 DNA Ligase | Instant Sticky-End Master Mix | Blunt/TA Master Mix |
|----------------------------|--------------------------|-------------------------------|-------------------------------|---------------------|
| Format | Kit | Enzyme | Master Mix | Master Mix |
| Vector (3 kb) | 50 ng | 50 ng | 50 ng | 50 ng |
| Insert (1 kb) | 50 ng | 50 ng | 50 ng | 50 ng |
| Buffer | 2X Quick Ligation Buffer | T4 DNA Ligase Reaction Buffer | 5 µl (Master Mix) | 5 µl (Master Mix) |
| Ligase | 1 µl | 1 µl | N/A | N/A |
| Nuclease-free water | to 20 µl | to 20 µl | to 10 µl | to 10 µl |
| Incubation | 25°C, 5 minutes | 25°C, 2 hrs; 16°C, overnight* | N/A, instant ligation | 25°C, 15 minutes |

* For sticky-end ligation, the incubation time can be shortened to 25°C for 10 minutes.

TIPS FOR OPTIMIZATION

REACTION BUFFERS

- T4 DNA Ligase Buffer should be thawed on the bench or in the palm of your hand, and not at 37°C (to prevent breakdown of ATP).
- Once thawed, T4 DNA Ligase Buffer should be placed on ice.
- Ligations can also be performed in any of the four standard restriction endonuclease NEBuffers or in T4 Polynucleotide Kinase Buffer supplemented with 1 mM ATP.
- When supplementing with ATP, use ribo-ATP. Deoxyribo-ATP will inhibit ligation.
- Before ligation, completely inactivate the restriction enzyme by heat inactivation, spin column (e.g., Monarch PCR & DNA Cleanup Kit) or Phenol/EtOH purification.

DNA

- Use purified DNA preparations without high EDTA or salt concentrations.
- Either heat inactivate (AP, rSAP) or remove phosphatase (CIP, BAP or SAP) before ligation.
- Keep total DNA concentration between 5–10 µg/ml.
- Vector:Insert molar ratios between 1:1 and 1:10 are optimal for single insertions.
- For cloning more than one insert, we recommend the NEBuilder HiFi DNA Assembly Master Mix or Cloning Kit.
- If you are unsure of your DNA concentration, perform multiple ligations with varying ratios.

LIGASE

- For most ligations (blunt or cohesive), the Quick Ligation Kit or the master mixes are recommended.
- For large inserts, reduce insert concentration and use concentrated ligase at 16°C overnight.
- T4 DNA Ligase can be heat inactivated at 65°C for 20 minutes.
- Do not heat inactivate if there is PEG in the reaction buffer because transformation will be inhibited.
- Electroporation is recommended for large constructs (> 10,000 bp). If planning to electroporate, we recommend ElectroLigase for your ligation step.

TRANSFORMATION

- Add between 1–5 µl of ligation mixture to competent cells for transformation.
- Extended ligation with PEG causes a drop off in transformation efficiency.
- Do not heat inactivate when using the Quick Ligation Buffer or Ligase Master Mixes, as this will inhibit transformation.



DNA Ligase Selection Chart for Cloning

| | Instant Sticky-end Ligase Master Mix | Blunt/TA Ligase Master Mix | ElectroLigase® | T4 DNA Ligase | Quick Ligation Kit | T3 DNA Ligase | T7 DNA Ligase | Tag DNA Ligase |
|---------------------------------|--------------------------------------|----------------------------|----------------|---------------|--------------------|---------------|---------------|----------------|
| DNA APPLICATIONS | | | | | | | | |
| Ligation of sticky ends | ●●● | ●● | ●● | ●● | ●●● | ●● | ●● | ● |
| Ligation of blunt ends | ● | ●●● | ●● | ●● | ●●● | ●● | | |
| T/A cloning | ● | ●●● | ●● | ●● | ●● | ● | ● | |
| Electroporation | | | ●●● | ●● | | | | |
| Ligation of sticky ends only | | | | | | | ●●● | |
| Repair of nicks in dsDNA | ●●● | ●●● | ●●● | ●●● | ●●● | ●●● | ●●● | ●●● |
| High complexity library cloning | ●● | ●● | ●● | ●●● | ●● | | | |

| FEATURES | | | | | | | | |
|---|---|---|---|---|---|---|---|---|
| Salt tolerance (> 2X that of T4 DNA Ligase) | | | | | | ✓ | | |
| Ligation in 15 min. or less | ✓ | ✓ | | ✓ | ✓ | ✓ | ✓ | ✓ |
| Master Mix Formulation | ✓ | ✓ | | | | | | |
| Thermostable | | | | | | | | ✓ |
| Recombinant | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |

KEY

- Recommended product(s) for selected application
- Works well for selected application
- Will perform selected application, but is not recommended

GETTING STARTED

For traditional cloning, follow the ligation guidelines specified by the ligase supplier. If they suggest a 3:1 molar ratio of insert to vector, try this first for the best result. Using a 3:1 mass ratio is not the same thing (unless the insert and vector have the same mass). Ligation usually proceeds very quickly and, unless your cloning project requires the generation of a high-complexity library that benefits from the absolute capture of every possible ligation product, long incubation times are not necessary.



Transformation

Transformation is the process by which an organism acquires exogenous DNA. Transformation can occur in two ways: natural transformation and artificial transformation. Natural transformation describes the uptake and incorporation of naked DNA from the cell's natural environment. Artificial transformation encompasses a wide array of methods for inducing uptake of exogenous DNA. In cloning protocols, artificial transformation is used to introduce recombinant DNA into host bacteria. The most common method of artificial transformation of bacteria involves use of divalent cations (e.g., calcium chloride) to increase the permeability of the bacterium's membrane, making them chemically competent, and thereby increasing the likelihood of DNA acquisition. Another artificial method of transformation is electroporation, in which cells are shocked with an electric current, to create holes in the bacterial membrane. With a newly-compromised cell membrane, the transforming DNA is free to pass into the cytosol of the bacterium. Regardless of which method of transformation is used, outgrowth of bacteria following transformation allows repair of the bacterial surface and selection of recombinant cells if the newly acquired DNA conveys antibiotic resistance to the transformed cells.

Protocol: High Efficiency Transformation

| | STANDARD PROTOCOL |
|---------------------------------|--|
| DNA | 1–5 µl containing 1 pg – 100 ng of plasmid DNA |
| Competent <i>E. coli</i> | 50 µl |
| Incubation | On ice for 30 minutes |
| Heat Shock | Exactly 42°C for exactly 30 seconds* |
| Incubation | On ice for 5 minutes Add 950 µl room temperature SOC 37°C for 60 minutes, with shaking |

* Follow specific heat shock recommendations provided for the *E. coli* competent cell strain being used.

Competent Cell Selection Chart

| | NEB 5-alpha Competent <i>E. coli</i> | NEB Turbo Competent <i>E. coli</i> | NEB 5-alpha F' I ⁺ Competent <i>E. coli</i> | NEB 10-beta Competent <i>E. coli</i> | <i>dam</i> ⁻ / <i>dcm</i> ⁻ Competent <i>E. coli</i> | NEB Stable Competent <i>E. coli</i> |
|---|--|--|--|--|--|---|
| FEATURES | | | | | | |
| Versatile | • | | | • | | • |
| Fast growth (< 8 hours) | | • | | | | |
| Toxic gene cloning | | • | • | | | • |
| Large plasmid/BAC cloning | | | | • | | • |
| Dam/Dcm-free plasmid growth | | | | | • | |
| Retroviral/lentiviral vector cloning | | | | | | • |
| FORMATS | | | | | | |
| Chemically competent | • | • | • | • | • | • |
| Electrocompetent | • | • | | • | | |
| Subcloning | • | | | | | |
| 96-well format* | • | | | | | |
| 384-well format* | • | | | | | |
| 12x8-tube strips* | • | | | | | |

* Other strains are available upon request.

TIPS FOR OPTIMIZATION

THAWING

- Cells are best thawed on ice.
- DNA should be added as soon as the last trace of ice in the tube disappears.
- Cells can be thawed by hand, but warming above 0°C decreases efficiency.

DNA

- Up to 5 µl of DNA from a ligation mix can be used with only a 2-fold loss of efficiency.

INCUBATION & HEAT SHOCK

- Incubate on ice for 30 minutes. Expect a 2-fold loss in transformation efficiency (TE) for every 10 minutes this step is shortened.
- Both temperature and time are specific to the transformation volume and vessel. Typically, 30 seconds at 42°C is recommended.

OUTGROWTH

- Outgrowth at 37°C for 1 hour is best for cell recovery and for expression of antibiotic resistance. Expect a 2-fold loss in TE for every 15 minutes this step is shortened.
- SOC gives 2-fold higher TE than LB medium.
- Incubation with shaking or rotation results in 2-fold higher TE.

PLATING

- Selection plates can be used warm or cold, wet or dry with no significant effects on TE.
- Warm, dry plates are easier to spread and allow for the most rapid colony formation.

DNA CONTAMINANTS TO AVOID


| CONTAMINANT | REMOVAL METHOD |
|------------------------|--|
| Detergents | Ethanol precipitate |
| Phenol | Extract with chloroform and ethanol precipitate |
| Ethanol or Isopropanol | Dry pellet before resuspending |
| PEG | Column purify (e.g., Monarch PCR & DNA Cleanup Kit) or phenol/chloroform extract and ethanol precipitate |



DNA Analysis

Agarose-gel electrophoresis is the standard method used for separation, identification and purification of DNA fragments. DNA is visualized on a gel after soaking or pre-impregnating the gel with ethidium bromide, a DNA intercalating agent that fluoresces under UV illumination. Using the marker or ladder as a reference, it is possible to determine the size and relative quantity of the DNA of interest. The original DNA markers were made of genomic DNAs digested with a restriction enzyme to exhibit a banding pattern of known fragment sizes. Later, markers were made of fragments with evenly-spaced sizes and the resulting banding pattern resembles a ladder. The bands are visible under UV illumination; since the bands of the marker/ladder are not visible under normal lighting conditions. To track the progress of the gel as it runs, the marker contains a dye or combination of dyes that identify the leading edge of well contents, also called the dye front.


The following DNA Ladders are Now Available in Quick-Load Purple Format



* Available in Quick-Load® and TriDye™ formats

*** Free Loading Dye included

Additional DNA Ladders from New England Biolabs



*** Free Loading Dye included




Cloning & Mutagenesis Kits


NEBuilder HiFi DNA Assembly

NEBuilder HiFi DNA Assembly enables virtually error-free joining of DNA fragments, even those with 5'- and 3'-end mismatches. Available with and without competent *E. coli*, this flexible kit enables simple and fast seamless cloning utilizing a new proprietary high-fidelity polymerase. Make NEBuilder HiFi your first choice for DNA assembly and cloning.

Overview of the NEBuilder HiFi DNA Assembly cloning method



NEBuilder HiFi DNA Assembly offers improved efficiency and accuracy with lower amounts of DNA by increasing overlap length



Reactions were set up in a 4-fragment assembly reaction according to recommended reaction conditions. Amount of DNA and size of overlap is shown.

RECOMMENDED PRODUCTS

NEBuilder HiFi DNA Assembly Cloning Kit

NEBuilder HiFi DNA Assembly Master Mix

NEBuilder HiFi DNA Assembly Bundle for Large Fragments

- Simple and fast seamless cloning in as little as 15 minutes
- Use one system for both "standard-size" cloning and larger gene assembly products (up to 12 fragments and 20 kb)
- DNA can be used immediately for transformation or as template for PCR or RCA
- Adapts easily for multiple DNA manipulations, including site-directed mutagenesis
- Enjoy less screening/re-sequencing of constructs, with virtually error-free, high-fidelity assembly
- Use NEBuilder HiFi in successive rounds of assembly, as it removes 5'- and 3'-end mismatches
- Bridge two ds-fragments with a synthetic ssDNA oligo for simple and fast construction (e.g., linker insertion or gRNA library)
- No licensing fee requirements from NEB for NEBuilder products
- NEBuilder HiFi DNA Assembly Cloning Kit includes the NEBuilder HiFi DNA Assembly Master Mix and NEB 5-alpha Competent *E. coli*



Cloning & Mutagenesis Kits (Cont.)

Optimization Tips for NEBuilder HiFi DNA Assembly

Assembly Reaction

- When directly assembling fragments into a cloning vector, the molar concentration of assembly fragments should be 2–3 times higher than the concentration of vector.
- For multiple (4–12) fragment assembly, design 25–30 bp overlap regions between each fragment to enhance assembly efficiency. Use 0.05 pmol of each fragment in the assembly reaction.
- For assembly of 1–3 fragments, 15 minute incubation times are sufficient. For assembly of 4–6 fragments, 60 minute incubation times are recommended. Reaction times less than 15 minutes are generally not recommended.

Primer Design

- For help with primer design, we recommend using NEBuilder Assembly Tool.

Transformation

- The NEBuilder HiFi DNA Assembly Cloning Kit and the Gibson Assembly Cloning Kit include NEB 5-alpha Competent *E. coli*. NEB recommends using the competent cells provided with the kit because of their high efficiency. The components of the master mix may inhibit the functionality of competent cells from other companies if not diluted. The NEBuilder HiFi DNA Assembly Bundle for large fragments includes NEB 10-beta Competent *E. coli*, ideal for assembling larger fragments (>15 kb).

Protocol: Assembly

Before use, thaw and vortex the master mix thoroughly and keep on ice.

1. Set up the following reaction on ice.

| | REACTION SETUP | | |
|---|-------------------------------|-----------------------------|-----------------------|
| | 2–3 Fragment Assembly | 4–6 Fragment Assembly | Positive Control** |
| PCR Fragment(s) + linearized vector | X µl (0.02– 0.5 pmols)* | X µl (0.2– 1 pmols)* | 10 µl |
| Assembly Master Mix (2X) | 10 µl | 10 µl | 10 µl |
| Deionized H ₂ O | 10–X µl | 10–X µl | 0 |
| Total Volume | 20 µl*** | 20 µl*** | 20 µl |

* Optimized cloning efficiency requires about 50–100 ng of vector and at least 2-fold excess inserts. Use 5X more insert if the size is less than 200 bps.

** Control reagents are provided for five reactions.

*** If greater numbers of fragments are assembled, additional Gibson Assembly Master Mix may be required.

2. Incubate samples in a thermocycler at 50°C for 15–60 minutes depending on number of fragments being assembled. Following incubation, store samples on ice or at –20°C for subsequent transformation.

Protocol: Transformation with NEB 5-alpha cells

| | STANDARD PROTOCOL |
|-----------------------------|--|
| DNA | 2 µl |
| Competent <i>E. coli</i> | 50 µl |
| Incubation | On ice for 30 minutes |
| Heat Shock | Exactly 42°C for exactly 30 seconds |
| Incubation | On ice for 5 minutes Add 950 µl room temperature SOC 37°C for 60 minutes, with shaking |

Gibson Assembly Cloning Kit

Gibson Assembly enables multiple, overlapping DNA fragments to be joined in a single-tube isothermal reaction, with no additional sequence added (scar-less). The Gibson Assembly Master Mix includes three different enzymatic activities that perform in a single buffer (described below). The assembled, fully-sealed construct is then transformed into NEB 5-alpha competent *E. coli*. The entire protocol, from assembly to transformation, takes just under two hours.




NEB Golden Gate Assembly

The efficient and seamless assembly of DNA fragments, commonly referred to as Golden Gate assembly (1,2), has its origins in 1996 when, for the first time, it was shown that multiple inserts could be assembled into a vector backbone using only the sequential (3) or simultaneous (4) activities of a single Type IIS restriction enzyme and T4 DNA Ligase. This method can be accomplished using Type IIS restriction enzymes, such as BsaI, and can also be used for the cloning of single inserts. The method is efficient and can be completed in one tube in as little as 10 minutes for single inserts, or can utilize cycling steps for multiple inserts (see page 3 for workflow).

The NEB Golden Gate Assembly Mix incorporates digestion with BsaI and ligation with T4 DNA Ligase into a single reaction, and can be used to assemble up to 10 fragments in a single step.

NEB Golden Gate Assembly Mix offers improved assembly



Assembly reactions were set up using either precloned inserts or PCR amplicons directly. Reaction conditions were set up according to manufacturer, and are shown above. Two separate experiments are shown for each reaction type.

RECOMMENDED PRODUCTS

NEB Golden Gate Assembly Mix

- Seamless cloning – no scar remains following assembly
- Ordered assembly of up to 10–20 fragments in a single reaction
- Efficient with regions with high GC content and areas of repeats
- Compatible with a broad range of fragment sizes (< 100 bp to > 15 kb)

References:

1. Engler, C. et al. (2008) *PLoS ONE*, 3: e3647.
2. Engler, C. et al. (2009) *PLoS ONE*, 4: e5553.
3. Lee, J.H. et al. (1996) *Genetic Analysis: Biomolecular Engineering*, 13; 139-145.
4. Padgett, K.A. and Sorge, J.A. (1996) *Gene*, 168, 31-35.

Time for change.

Introducing Monarch Nucleic Acid Purification Kits


- Optimized for maximum performance and minimal environmental impact
- Unique column design eliminates buffer retention and offers elution in lower volumes
- Fast user-friendly protocols
- Available for plasmid miniprep, DNA gel extraction and reaction cleanup





NEB PCR Cloning Kit

The NEB PCR Cloning Kit (with or without competent cells) enables quick and simple cloning of all your PCR amplicons, regardless of the polymerase used. This kit utilizes a novel mechanism for background colony suppression – a toxic minigene is generated when the vector closes upon itself – and allows for direct cloning from your reaction, with no purification step. The NEB PCR Cloning Kit is supplied with the pMiniT 2.0 vector, which allows *in vitro* transcription from both SP6 and T7 promoters, features more unique restriction sites for subcloning (including four 8-base cut sites) and can be used for Golden Gate Assembly as the plasmid has no internal BsaI sites.



TIPS FOR OPTIMIZATION

- For first time use of the kit, prepare a positive control reaction containing 2 µl (30 ng) of the 1 kb amplicon cloning control included with the kit.
- 3:1 insert:vector ratio is best, but ratios from 1:1 to 10:1 can also be utilized.

Protocol: Ligation

| | STANDARD PROTOCOL |
|---|--------------------|
| Linearized pMiniT 2.0 Vector (25 ng/µl) | 1 µl |
| Insert + H ₂ O | 4 µl |
| Cloning Mix 1 | 4 µl |
| Cloning Mix 2 | 1 µl |
| Incubation | 5–15 minutes, 25°C |


Protocol: Transformation

| | STANDARD PROTOCOL |
|--------------------------|--|
| Ligation Reaction | 2 µl |
| Competent <i>E. coli</i> | 50 µl |
| Incubation | On ice for 20 minutes |
| Heat Shock | 42°C for exactly 30 seconds |
| Incubation | On ice for 5 minutes Add 950 µl room temperature SOC 37°C for 60 minutes, with shaking |

Protocol: Plating

1. Mix cells thoroughly by flicking or inversion and spread 50 µl of the 1 ml outgrowth onto 37°C pre-warmed agar plates containing 100 µg/ml ampicillin. If a 15 minute ligation time was used, also plate 50 µl of a 1:10 dilution prepared with SOC.
2. Invert plate and incubate overnight at 37°C or for 24 hours at 30°C. Do not use room temperature growth as the slow growth rate will interfere with selection of constructs with inserts.
3. After colonies appear, use the plate with well separated colonies for screening.

Features within Sequence Flanking the Toxic Minigene/Cloning Site




Top map shown above displays the construct formed if no insert is present. Unique restriction sites are shown in bold. Additional restriction sites that can be used for subcloning are also shown. Expanded box below shows location of cloning analysis primers for cloning PCR or sequencing, restriction sites for subcloning or linearization for *in vitro* transcription, RNA Polymerase promoter sequences and placement of insertion site within the toxic minigene.



Q5 Site-Directed Mutagenesis Kit

The Q5 Site-Directed Mutagenesis Kit (with or without competent cells) enables rapid, site-specific mutagenesis of double-stranded plasmid DNA in less than 2 hours. The kit utilizes Q5 Hot Start High-Fidelity DNA Polymerase, along with custom mutagenic primers to create substitutions, deletions and insertions in a wide variety of plasmids. Transformation into high-efficiency NEB 5-alpha Competent *E. coli* cells ensures robust results with plasmids up to, at least, 14 kb in length.

Overview of Q5 Site-Directed Mutagenesis Kit



TIPS FOR OPTIMIZATION

- No purification of your plasmid is necessary, either before or after the KLD reaction.
- You can expect a high frequency of your desired mutation (> 90%).
- While the Q5 SDM Kit is supplied with high-efficiency, NEB competent *E. coli*, you can use your own chemically competent cells for cloning; results will vary, according to the quality and efficiency of the cells.

Protocol: Assembly

Before use, thaw and vortex the master mix thoroughly and keep on ice.

1. Exponential Amplification

| | 25 μ L RXN | FINAL CONC. |
|--|----------------|-------------|
| Q5 Hot Start High-Fidelity 2X Master Mix | 12.5 μ l | 1X |
| 10 μ M Forward Primer | 1.25 μ l | 0.5 μ M |
| 10 μ M Reverse Primer | 1.25 μ l | 0.5 μ M |
| Template DNA (1–25 ng/ μ l) | 1 μ l | 1–25 ng |
| Nuclease-free water | 9.0 μ l | |

2. KLD Reaction

| | VOLUME | FINAL CONC. |
|------------------------|-----------|-------------|
| PCR Product | 1 μ l | |
| 2X KLD Reaction Buffer | 5 μ l | 1X |
| 10X KLD Enzyme Mix | 1 μ l | 1X |
| Nuclease-free Water | 3 μ l | |

Protocol: Transformation with NEB 5-alpha

| | STANDARD PROTOCOL |
|--------------------------|---|
| KLD Mix | 5 μ l |
| Competent <i>E. coli</i> | 50 μ l |
| Incubation | On ice for 30 minutes |
| Heat Shock | Exactly 42°C for exactly 30 seconds |
| Incubation | On ice for 5 minutes Add 950 μ l room temperature SOC 37°C for 60 minutes, with shaking |



DNA Assembly Selection Chart

New England Biolabs now offers several products that can be used for DNA assembly and cloning. Use this chart to determine which product would work best to assemble your DNA.

| | NEBuilder HiFi DNA Assembly | Gibson Assembly | NEB Golden Gate Assembly Mix | USER™ Enzyme |
|--|-----------------------------|-----------------|------------------------------|--------------|
| PROPERTIES | | | | |
| Removes 5' or 3' End Mismatches | ★★★ | ★ | N/A | N/A |
| Assembles with High Fidelity at Junctions | ★★★ | ★★ | ★★★ | ★★★ |
| Tolerates Repetitive Sequences at Ends | ★ | ★ | ★★★ | ★★★ |
| Generates Fully Ligated Product | ★★★ | ★★★ | ★★★ | NR |
| Joins dsDNA with Single-stranded Oligo | ★★★ | ★★ | NR | NR |
| Assembles with High Efficiency with Low Amounts of DNA | ★★★ | ★★ | ★★ | ★★ |
| Accommodates Flexible Overlap Lengths | ★★★ | ★★★ | ★ | ★★ |

| | | | | |
|---|-----|-----|-----|-----|
| APPLICATIONS | | | | |
| Simple Cloning (1-2 Fragments) | ★★★ | ★★★ | ★★★ | ★★★ |
| 4-6 Fragment Assembly | ★★★ | ★★★ | ★★★ | ★★★ |
| >6 Fragment Assembly | ★★★ | ★★ | ★★★ | ★★★ |
| Template Construction for <i>In vitro</i> Transcription | ★★★ | ★★★ | ★★★ | ★ |
| Synthetic Whole Genome Assembly | ★★★ | ★ | ★ | ★ |
| Multiple Site-directed Mutagenesis | ★★★ | ★★ | ★★ | ★★ |
| Library Generation | ★★ | ★★ | ★★ | ★★ |
| Pathway Engineering | ★★★ | ★★ | ★★ | ★★★ |
| TALENs | ★★ | ★★ | ★★★ | ★★ |
| Short Hairpin RNA Cloning (shRNA) | ★★★ | ★★ | ★ | ★ |
| gRNA Library Generation | ★★★ | ★★ | ★ | ★ |
| Large Fragment (>10 kb) Assembly | ★★★ | ★★★ | ★★★ | ★★ |
| Small Fragment (<100 bp) Assembly | ★★★ | ★ | ★★★ | ★★★ |
| Use in Successive Rounds of Restriction Enzyme Assembly | ★★★ | ★ | NR | ★ |

KEY

| | | | |
|-----|--|-----|------------------------------------|
| ★★★ | Works best for selected application | N/A | Not applicable to this application |
| ★★ | Suitable for selected application, but other product(s) perform better | NR | Not recommended |
| ★ | Will perform selected application, but is not recommended | | |



Traditional Cloning Quick Guide

Preparation of insert and vectors

Insert from a plasmid source

- Digest plasmid with the appropriate restriction enzymes to produce a DNA fragment that can be cloned directly into a vector. Unidirectional cloning is achieved with restriction enzymes that produce non-compatible ends.

Insert from a PCR product

- Design primers with appropriate restriction sites to clone unidirectionally into a vector
- Addition of 6 bases upstream of the restriction site is sufficient for digestion with most enzymes
- If fidelity is a concern, choose a proofreading polymerase such as Q5 High-Fidelity DNA Polymerase
- Purify PCR product by running the DNA on an agarose gel and excising the band or by using a spin column (e.g., Monarch DNA Gel Extraction Kit, Monarch PCR & DNA Cleanup Kit)
- Digest with the appropriate restriction enzyme

Standard Restriction Enzyme Protocol

| | |
|------------------------|--|
| DNA | 1 µg |
| 10X NEBuffer | 5 µl (1X) |
| Restriction Enzyme | 10 units is sufficient, generally 1 µl is used |
| Nuclease-free Water | To 50 µl |
| Incubation Time | 1 hour* |
| Incubation Temperature | Enzyme dependent |

* Can be decreased by using a Time-Saver qualified enzyme.

Time-Saver Restriction Enzyme Protocol

| | |
|------------------------|------------------|
| DNA | 1 µg |
| 10X NEBuffer | 5 µl (1X) |
| Restriction Enzyme | 1 µl |
| Nuclease-free Water | To 50 µl |
| Incubation Time | 5–15 minutes* |
| Incubation Temperature | Enzyme dependent |

* Time-Saver qualified enzymes can also be incubated overnight with no star activity

Insert from annealed oligos

- Annealed oligos can be used to introduce a fragment (e.g., promoter, polylinker, etc.)
- Anneal two complementary oligos that leave protruding 5' or 3' overhangs for ligation into a vector cut with appropriate enzymes
- Non-phosphorylated oligos can be phosphorylated using T4 Polynucleotide Kinase

Typical Annealing Reaction

| | |
|----------------------|---|
| Primer | 1 µg |
| 10X T4 Ligase Buffer | 5 µl |
| Nuclease-free Water | To 50 µl |
| Incubation | 85°C for 10 minutes, cool slowly (30–60 min.) |

Vector

- Digest vector with appropriate restriction enzymes. Enzymes that leave non-compatible ends are ideal as they prevent vector self-ligation.

Dephosphorylation

- Dephosphorylation is sometimes necessary to prevent self ligation. NEB offers four products for dephosphorylation of DNA:
- The Quick Dephosphorylation Kit, Shrimp Alkaline Phosphatase (rSAP) and Antarctic Phosphatase (AP) are heat-inactivatable phosphatases. They work in all NEBuffers, but AP requires supplementation with Zn²⁺. The Quick Dephosphorylation Kit is optimized for fast and robust dephosphorylation in 10 minutes, and is heat inactivated in 2 minutes.
- Calf Intestinal Phosphatase (CIP) will function under many different conditions and in most NEBuffers. However, CIP cannot be heat inactivated and requires a purification step (such as with Monarch PCR & DNA Cleanup Kit) before ligation.

Dephosphorylation of 5' ends of DNA using the Quick Dephosphorylation Kit

| | |
|---------------------|---------------------|
| DNA | 1 pmol of DNA ends |
| 10X CutSmart Buffer | 2 µl |
| Quick CIP | 1 µl |
| Nuclease-free Water | to 20 µl |
| Incubation | 37°C for 10 minutes |
| Heat Inactivation | 80°C for 2 minutes |

Note: Scale larger reaction volumes proportionally.

Blunting

- In some instances, the ends of the insert or vector require blunting
- PCR with a proofreading polymerase will leave a predominantly blunt end
- T4 DNA Polymerase or Klenow will fill in a 5' overhang and chew back a 3' overhang
- The Quick Blunting Kit is optimized to blunt and phosphorylate DNA ends for cloning in less than 30 minutes
- Analyze agarose gels with longwave UV (360 nm) to minimize UV exposure that may cause DNA damage

Blunting with the Quick Blunting Kit

| | |
|-----------------------|--|
| DNA | Up to 5 µg |
| Blunting Buffer (10X) | 2.5 µl |
| dNTP Mix (1 mM) | 2.5 µl |
| Blunt Enzyme Mix | 1 µl |
| Nuclease-free Water | To 25 µl |
| Incubation | 15 minutes for RE-digested DNA/sheared or 30 minutes for nebulized DNA or PCR products |
| Heat Inactivation | 70°C for 10 minutes |

* PCR-generated DNA must be purified before blunting using a commercial purification kit, phenol extraction/ethanol precipitation or gel electrophoresis (e.g., Monarch PCR & DNA Cleanup Kit)



Traditional Cloning Quick Guide (Cont.)

Phosphorylation

- For ligation to occur, at least one of the DNA ends (insert or vector) should contain a 5' phosphate
- Primers are usually supplied non-phosphorylated; therefore, the PCR product will not contain a 5' phosphate
- Digestion of DNA with a restriction enzyme will always produce a 5' phosphate
- A DNA fragment can be phosphorylated by incubation with T4 Polynucleotide Kinase

Phosphorylation With T4 PNK

| | |
|---------------------|-------------------------|
| DNA (20 mer) | 1–2 µg |
| 10X T4 PNK Buffer | 5 µl |
| 10 mM ATP | 5 µl (1 mM final conc.) |
| T4 PNK | 1 µl (10 units) |
| Nuclease-free Water | To 50 µl |
| Incubation | 37°C for 30 minutes |

Purification of Vector and Insert

- Purify the vector and insert by either running the DNA on an agarose gel and excising the appropriate bands or by using a spin column, such as Monarch DNA Gel Extraction Kit
- DNA can also be purified using β-Agarase I with low melt agarose or an appropriate spin column or resin
- Analyze agarose gels with longwave UV (360 nm) to minimize UV exposure that may cause DNA damage

Ligation of Vector and Insert

- Use a molar ratio of 1:3 vector to insert
- If using T4 DNA Ligase or the Quick Ligation Kit, thaw and resuspend the Ligase Buffer at room temp. If using Ligase Master Mixes, no thawing is necessary.
- The Quick Ligation Kit is optimized for ligation of both sticky and blunt ends
- Instant Sticky-end Ligase Master Mix is optimized for instant ligation of sticky/cohesive ends
- Blunt/TA Ligase Master Mix is optimized for ligation of blunt or single base overhangs, which are the more challenging type of ends for T4 DNA Ligase
- Following ligation, chill on ice and transform
- DO NOT heat inactivate when using the Quick Ligation Buffer or Ligase Master Mixes, as this will inhibit transformation
- Electroligase is optimized for ligation of both sticky and blunt ends and is compatible with electroporation (i.e., no cleanup step required)

Ligation with the Quick Ligation Kit

| | |
|--------------------------|--------------------------------|
| Vector DNA (3 kb) | 50 ng |
| Insert DNA (1 kb) | To 50 ng |
| 2X Quick Ligation Buffer | 10 µl |
| Quick T4 DNA Ligase | 1 µl |
| Nuclease-free Water | 20 µl (mix well) |
| Incubation | Room temperature for 5 minutes |

Ligation with Instant Sticky-end Ligase Master Mix

| | |
|---------------------|----------|
| Vector DNA (3 kb) | 50 ng |
| Insert DNA (1 kb) | 50 ng |
| Master Mix | 5 µl |
| Nuclease-free Water | To 10 µl |
| Incubation | None |

Ligation with Blunt/TA Ligase Master Mix

| | |
|---------------------|---------------------------------|
| Vector DNA (3 kb) | 50 ng |
| Insert DNA (1 kb) | 50 ng |
| Master Mix | 5 µl |
| Nuclease-free Water | To 10 µl |
| Incubation | Room temperature for 15 minutes |

Transformation

- To obtain transformants in 8 hrs., use NEB Turbo Competent *E. coli*
- If recombination is a concern, then use the RecA⁻ strains NEB 5-alpha Competent *E. coli* or NEB-10 beta Competent *E. coli* or NEB Stable Competent *E. coli*
- NEB-10 beta Competent *E. coli* works well for constructs larger than 5 kb
- NEB Stable Competent *E. coli* can be used for constructs with repetitive sequences such as lentiviral constructs
- If electroporation is required, use NEB 5-alpha Electrocompetent *E. coli* or NEB 10-beta Electrocompetent *E. coli*
- Use pre-warmed selection plates
- Perform several 10-fold serial dilutions in SOC for plating

Transformation with NEB 5-alpha Competent *E. coli*

| | |
|--------------------------|--|
| DNA | 1–5 µl containing 1 pg – 100 ng of plasmid DNA |
| Competent <i>E. coli</i> | 50 µl |
| Incubation | On ice for 30 minutes |
| Heat Shock | Exactly 42°C for exactly 30 seconds |
| Incubation | On ice for 5 minutes Add 950 µl room temperature SOC 37°C for 60 minutes, with shaking |



Troubleshooting Guide for Cloning

We strongly recommend running the following controls during transformations. These controls may help troubleshoot which step(s) in the cloning workflow has failed.

- 1 Transform 100 pg – 1ng of uncut vector to check cell viability, calculate transformation efficiency and verify the antibiotic resistance of the plasmid.
- 2 Transform the cut vector to determine the amount of background due to undigested plasmid. The number of colonies in this control should be < 1% of the number of colonies in the uncut plasmid control transformation (from control #1).
- 3 Transform a vector only ligation reaction. The ends of the vector should not be able to re-ligate because either they are incompatible (e.g., digested with two restriction enzymes that do not generate compatible ends) or the 5' phosphate group has been removed in a dephosphorylation reaction (e.g., blunt ends treated with rSAP). This control transformation should yield the same number of colonies as control #2.
- 4 Digest vector DNA with a single restriction enzyme, re-ligate and transform. The ends of the vector DNA should be compatible and easily joined during the ligation reaction, resulting in approximately the same number of colonies as control #1.

The cloning workflow often benefits from an accurate quantitation of the amount of DNAs that are being worked with. We recommend quantification of DNAs whenever possible.

| PROBLEM | CAUSE | SOLUTION |
|-----------------------------|---|---|
| Few or no transformants | Cells are not viable | <ul style="list-style-type: none"> Transform an uncut plasmid (e.g., pUC19) and calculate the transformation efficiency of the competent cells. If the transformation efficiency is low (<10⁶) re-make the competent cells or consider using commercially available high efficiency competent cells. |
| | Incorrect antibiotic or antibiotic concentration | <ul style="list-style-type: none"> Confirm antibiotic and antibiotic concentration |
| | DNA fragment of interest is toxic to the cells | <ul style="list-style-type: none"> Incubate plates at lower temperature (25–30°C). Transformation may need to be carried out using a strain that exerts tighter transcriptional control over the DNA fragment of interest (e.g., NEB-5-alpha F' /^h Competent <i>E. coli</i>) |
| | If using chemically competent cells, the wrong heat-shock protocol was used | <ul style="list-style-type: none"> Follow the manufacturer's specific transformation protocol (Note: going above the recommended temperature during the heat shock can result in competent cell death) |
| | If using electrocompetent cells, PEG is present in the ligation mix | <ul style="list-style-type: none"> Clean up DNA by drop dialysis prior to transformation Try NEB's ElectroLigase |
| | If using electrocompetent cells, arcing was observed or no voltage was registered | <ul style="list-style-type: none"> Clean up the DNA prior to the ligation step Tap the cuvette to get rid of any trapped air bubbles Be sure to follow the manufacturer's specified electroporation parameters |
| | Construct is too large | <ul style="list-style-type: none"> Select a competent cell strain that can be transformed efficiently with large DNA constructs (≥ 10 kb, we recommend trying NEB 10-beta Competent <i>E. coli</i>) For very large constructs (> 10 kb), consider using electroporation |
| | Construct may be susceptible to recombination | <ul style="list-style-type: none"> Select a Rec A- strain such as NEB 5-alpha or NEB 10-beta Competent <i>E. coli</i> or NEB Stable Competent <i>E. coli</i> |
| | The insert comes directly from mammalian or plant DNA and contains methylated cytosines, which are degraded by many <i>E. coli</i> strains | <ul style="list-style-type: none"> Use a strain that is deficient in McrA, McrBC and Mrr, such as NEB 10-beta Competent <i>E. coli</i> |
| | Too much ligation mixture was used | <ul style="list-style-type: none"> Use < 5 µl of the ligation reaction for the transformation |
| | Inefficient ligation | <ul style="list-style-type: none"> Make sure that at least one fragment being ligated contains a 5' phosphate moiety Vary the molar ratio of vector to insert from 1:1 to 1:10 Purify the DNA to remove contaminants such as salt and EDTA ATP will degrade after multiple freeze-thaws; repeat the ligation with fresh buffer Heat inactivate or remove the phosphatase prior to ligation Ligation of single base-pair overhangs (most difficult) may benefit from being carried out with Blunt/TA Master Mix, Quick Ligation Kit or concentrated T4 DNA Ligase Test the activity of the ligase by carrying out a ligation control with Lambda-HindIII digested DNA |
| Inefficient phosphorylation | <ul style="list-style-type: none"> Purify the DNA prior to phosphorylation. Excess salt, phosphate or ammonium ions may inhibit the kinase. If the ends are blunt or 5' recessed, heat the substrate/buffer mixture for 10 minutes at 70°C. Rapidly chill on ice before adding the ATP and enzyme, then incubate at 37°C. ATP was not added. Supplement the reaction with 1mM ATP, as it is required by T4 Polynucleotide Kinase Alternatively, use 1X T4 DNA Ligase Buffer (contains 1 mM ATP) instead of the 1X T4 PNK Buffer | |



Troubleshooting Guide for Cloning (cont.)

| PROBLEM | CAUSE | SOLUTION |
|---|---|--|
| Few or no transformants | Inefficient blunting | <ul style="list-style-type: none"> Heat inactivate or remove the restriction enzymes prior to blunting Clean up the PCR fragment prior to blunting Sonicated gDNA should be blunted for at least 30 minutes Do not use > 1 unit of enzyme/μg of DNA Do not incubate for > 15 minutes Do not incubate at temperatures > 12°C (for T4 DNA Polymerase) or > 24°C (for Klenow) Make sure to add a sufficient amount of dNTPs to the reaction (33 μM each dNTP for DNA Polymerase I, Large (Klenow) Fragment and 100 μM each dNTP for T4 DNA Polymerase). When using Mung Bean Nuclease, incubate the reaction at room temperature. Do not use > 1 unit of enzyme/μg DNA or incubate the reaction > 30 minutes. |
| | Inefficient A-Tailing | <ul style="list-style-type: none"> Clean up the PCR prior to A-tailing. NEB recommends the Monarch PCR & DNA Cleanup Kit. High-fidelity enzymes will remove any non-templated nucleotides. |
| | Restriction enzyme(s) didn't cleave completely | <ul style="list-style-type: none"> Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence Use the recommended buffer supplied with the restriction enzyme Clean up the DNA to remove any contaminants that may inhibit the enzyme. NEB recommends the Monarch PCR & DNA Cleanup Kit. When digesting a PCR fragment, make sure to have at least 6 nucleotides between the recognition site and the end of the DNA molecule |
| Colonies don't contain a plasmid | Antibiotic level used was too low | <ul style="list-style-type: none"> Increase the antibiotic level on plates to the recommended amount Use fresh plates with fresh antibiotics |
| | Satellite colonies were selected | <ul style="list-style-type: none"> Choose large, well-established colonies for analysis |
| Colonies contain the wrong construct | Recombination of the plasmid has occurred | <ul style="list-style-type: none"> Use a RecA⁻ strain such as NEB 5-alpha, or NEB 10-beta Competent <i>E. coli</i>, or NEB Stable Competent <i>E. coli</i> |
| | Incorrect PCR amplicon was used during cloning | <ul style="list-style-type: none"> Optimize the PCR conditions Gel purify the correct PCR fragment. NEB recommends the Monarch DNA Gel Extraction Kit. |
| | Internal recognition site was present | <ul style="list-style-type: none"> Use NEBcutter[®] to analyze insert sequence for presence of an internal recognition site |
| | DNA fragment of interest is toxic to the cells | <ul style="list-style-type: none"> Incubate plates at lower temperature (25–30°C) Transformation may need to be carried out using a strain that exerts tighter transcriptional control of the DNA fragment of interest (e.g., NEB 5-alpha F['] I['] Competent <i>E. coli</i>) |
| | Mutations are present in the sequence | <ul style="list-style-type: none"> Use a high-fidelity polymerase (e.g., Q5 High-Fidelity DNA Polymerase) Re-run sequencing reactions |
| Too much background | Inefficient dephosphorylation | <ul style="list-style-type: none"> Heat inactivate or remove the restriction enzymes prior to dephosphorylation |
| | Kinase is present/active | <ul style="list-style-type: none"> Heat inactivate the kinase after the phosphorylation step. Active kinase will re-phosphorylate the dephosphorylated vector. |
| | Restriction enzyme(s) didn't cleave completely | <ul style="list-style-type: none"> Check the methylation sensitivity of the restriction enzyme(s) to be sure it is not inhibited by methylation of the recognition sequence Use the recommended buffer supplied with the restriction enzyme Clean up the DNA to remove contaminants (e.g., too much salt). NEB recommends the Monarch PCR & DNA Cleanup Kit). |
| | Antibiotic level is too low | <ul style="list-style-type: none"> Confirm the correct antibiotic concentration |
| Ran the ligation on a gel and saw no ligated product | Inefficient ligation | <ul style="list-style-type: none"> Make sure at least one DNA fragment being ligated contains a 5' phosphate Vary the molar ratios of vector to insert from 1:1 to 1:10 Purify the DNA to remove contaminants such as salt and EDTA. NEB recommends the Monarch PCR & DNA Cleanup Kit. ATP will degrade after multiple freeze-thaws; repeat the ligation with fresh buffer Heat inactivate or remove the phosphatase prior to ligation Ligation of single base-pair overhangs (most difficult) may benefit from being carried out with Blunt/TA Master Mix, Quick Ligation Kit or concentrated T4 DNA Ligase Test the activity of the ligase by carrying out a ligation control with Lambda-HindIII digested DNA |
| The ligated DNA ran as a smear on an agarose gel | The ligase is bound to the substrate DNA | <ul style="list-style-type: none"> Treat the ligation reaction with Proteinase K prior to running on a gel |
| The digested DNA ran as a smear on an agarose gel | The restriction enzyme(s) is bound to the substrate DNA | <ul style="list-style-type: none"> Lower the number of units Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA |
| | Nuclease contamination | <ul style="list-style-type: none"> Use fresh, clean running buffer Use a fresh agarose gel Clean up the DNA. NEB recommends the Monarch PCR & DNA Cleanup Kit. |
| Incomplete restriction enzyme digestion | Cleavage is blocked by methylation | <ul style="list-style-type: none"> DNA isolated from a bacterial source may be blocked by Dam and Dcm methylation DNA isolated from eukaryotic source may be blocked by CpG methylation Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence If the enzyme is inhibited by Dam or Dcm methylation, grow the plasmid in a <i>dam-1 dcm-</i> strain |
| | Salt inhibition | <ul style="list-style-type: none"> Enzymes that have low activity in salt-containing buffers (NEBuffer 3.1) may be salt sensitive, so clean up the DNA prior to digestion. NEB recommends the Monarch PCR & DNA Cleanup Kit. DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. To prevent this, DNA solution should be no more than 25% of total reaction volume. |
| | Inhibition by PCR components | <ul style="list-style-type: none"> Clean up the PCR fragment prior to restriction digest. NEB recommends the Monarch PCR & DNA Cleanup Kit. |
| | Using the wrong buffer | <ul style="list-style-type: none"> Use the recommended buffer supplied with the restriction enzyme |
| | Too few units of enzyme used | <ul style="list-style-type: none"> Use at least 3–5 units of enzyme per μg of DNA |
| | Incubation time was too short | <ul style="list-style-type: none"> Increase the incubation time |
| | Digesting supercoiled DNA | <ul style="list-style-type: none"> Some enzymes have a lower activity on supercoiled DNA. Increase the number of enzyme units in the reaction. |



| PROBLEM | CAUSE | SOLUTION |
|--|--|--|
| Incomplete restriction enzyme digestion | Presence of slow sites | <ul style="list-style-type: none"> Some enzymes can exhibit slower cleavage towards specific sites. Increase the incubation time, 1–2 hours is typically sufficient. |
| | Two sites required | <ul style="list-style-type: none"> Some enzymes require the presence of two recognition sites to cut efficiently |
| | DNA is contaminated with an inhibitor | <ul style="list-style-type: none"> Assay substrate DNA in the presence of a control DNA. Control DNA will not cleave if there is an inhibitor present. Mini prep DNA is particularly susceptible to contaminants. NEB recommends the Monarch PCR & DNA Cleanup Kit. Clean DNA with a spin column, resin or drop dialysis, or increase volume to dilute contaminant |
| Extra bands in the gel | If larger bands than expected are seen in the gel, this may indicate binding of the enzyme(s) to the substrate | <ul style="list-style-type: none"> Lower the number of units in the reaction Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the substrate |
| | Star activity | <ul style="list-style-type: none"> Use the recommended buffer supplied with the restriction enzyme Decrease the number of enzyme units in the reaction Make sure the amount of enzyme added does not exceed 10% of the total reaction volume. This ensures that the total glycerol concentration does not exceed 5% v/v Decrease the incubation time. Using the minimum reaction time required for complete digestion will help prevent star activity. Try using a High-Fidelity (HF) restriction enzyme. HF enzymes have been engineered for reduced star activity. |
| | Partial restriction enzyme digest | <ul style="list-style-type: none"> Enzymes that have low activity in salt-containing buffers (e.g., NEBuffer 3.1) may be salt sensitive. Make sure to clean up the DNA prior to digestion. NEB recommends the Monarch PCR & DNA Cleanup Kit. DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. To prevent this, DNA solution should be no more than 25% of total reaction volume Clean-up the PCR fragment prior to restriction digest. NEB recommends the Monarch PCR & DNA Cleanup Kit. Use the recommended buffer supplied with the restriction enzyme Use at least 3–5 units of enzyme per μg of DNA Digest the DNA for 1–2 hours |
| No PCR fragment amplified | Used the wrong primer sequence | <ul style="list-style-type: none"> Double check the primer sequence |
| | Incorrect annealing temperature | <ul style="list-style-type: none"> Use the NEB Tm calculator to determine the correct annealing temperature |
| | Incorrect extension temperature | <ul style="list-style-type: none"> Each polymerase type has a different extension temperature requirement. Follow the manufacturer's recommendations. |
| | Too few units of polymerase | <ul style="list-style-type: none"> Use the recommended number of polymerase units based on the reaction volume |
| | Incorrect primer concentration | <ul style="list-style-type: none"> Each polymerase has a different primer concentration requirement. Make sure to follow the manufacturer's recommendations. |
| | Mg ²⁺ levels in the reaction are not optimal | <ul style="list-style-type: none"> Titrate the Mg²⁺ levels to optimize the amplification reaction. Follow the manufacturer's recommendations. |
| The PCR reaction is a smear on a gel | Difficult template | <ul style="list-style-type: none"> With difficult templates, try different polymerases and/or buffer combinations |
| | If bands are larger than expected it may indicate binding of the enzyme(s) to the DNA | <ul style="list-style-type: none"> Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA |
| Extra bands in PCR reaction | Annealing temperature is too low | <ul style="list-style-type: none"> Use the NEB Tm calculator to determine the annealing temperature of the primers |
| | Mg ²⁺ levels in the reaction are not optimal | <ul style="list-style-type: none"> Titrate the Mg²⁺ levels to optimize the amplification reaction. Make sure to follow the manufacturer's recommendations. |
| | Additional priming sites are present | <ul style="list-style-type: none"> Double check the primer sequence and confirm it does not bind elsewhere in the DNA template |
| | Formation of primer dimers | <ul style="list-style-type: none"> Primer sequence may not be optimal. Additional primers may need to be tested in the reaction. |
| | Incorrect polymerase choice | <ul style="list-style-type: none"> Try different polymerases and/or buffer combinations |



ORDERING INFORMATION

Selected Products for PCR

| PRODUCT | VWR CAT.NO. | SIZE |
|--|-------------|---------------|
| HIGH-FIDELITY DNA POLYMERASES | | |
| Q5 High-Fidelity DNA Polymerase | 102500-138 | 100 units |
| | 102500-136 | 500 units |
| Q5 Hot Start High-Fidelity DNA Polymerase | 102500-146 | 100 units |
| | 102500-144 | 500 units |
| Q5 High-Fidelity 2X Master Mix | 102500-134 | 100 reactions |
| | 102500-132 | 500 reactions |
| Q5 Hot Start High-Fidelity 2X Master Mix | 102500-142 | 100 reactions |
| | 102500-140 | 500 reactions |
| Q5 High-Fidelity PCR Kit | 102855-120 | 50 reactions |
| | 102855-122 | 200 reactions |
| Phusion High-Fidelity PCR Master Mix with HF Buffer | 101641-018 | 100 reactions |
| | 101641-016 | 500 reactions |
| Phusion High-Fidelity PCR Master Mix with GC Buffer | 101641-022 | 100 reactions |
| | 101641-020 | 500 reactions |
| Phusion Hot Start Flex 2X Master Mix | 102500-126 | 100 reactions |
| | 102500-124 | 500 reactions |
| Phusion High-Fidelity PCR Kit | 101640-804 | 50 reactions |
| | 101640-802 | 200 reactions |
| Phusion High-Fidelity DNA Polymerase | 101641-014 | 100 units |
| | 101641-012 | 500 units |
| Phusion Hot Start Flex High-Fidelity DNA Polymerase | 102500-130 | 100 units |
| | 102500-128 | 500 units |
| DNA POLYMERASES | | |
| One <i>Taq</i> DNA Polymerase | 101640-970 | 200 units |
| | 101640-968 | 1,000 units |
| | 101640-972 | 5,000 units |
| One <i>Taq</i> Hot Start DNA Polymerase | 101640-976 | 200 units |
| | 101640-974 | 1,000 units |
| | 101640-978 | 5,000 units |
| One <i>Taq</i> 2X Master Mix with Standard Buffer | 101640-982 | 100 reactions |
| | 101640-980 | 500 reactions |
| One <i>Taq</i> 2X Master Mix with GC Buffer | 101640-986 | 100 reactions |
| | 101640-984 | 500 reactions |
| One <i>Taq</i> Quick-Load 2X Master Mix with GC Buffer | 101641-002 | 100 reactions |
| | 101641-005 | 500 reactions |
| One <i>Taq</i> Quick-Load 2X Master Mix with Standard Buffer | 101640-998 | 100 reactions |
| | 101640-996 | 500 reactions |
| One <i>Taq</i> Hot Start 2X Master Mix with Standard Buffer | 101640-990 | 100 reactions |
| | 101640-988 | 500 reactions |
| One <i>Taq</i> Hot Start 2X Master Mix with GC Buffer | 101640-994 | 100 reactions |
| | 101640-992 | 500 reactions |
| One <i>Taq</i> Hot Start Quick-Load 2X Master Mix with Standard Buffer | 101641-006 | 100 reactions |
| | 101641-004 | 500 reactions |
| One <i>Taq</i> Hot Start Quick-Load 2X Master Mix with GC Buffer | 101641-010 | 100 reactions |
| | 101641-008 | 500 reactions |
| <i>Taq</i> DNA Polymerase with ThermoPol™ Buffer | 101227-638 | 400 units |
| | 101227-636 | 2,000 units |
| | 200067-186 | 4,000 units |
| | 102877-582 | 20,000 units |
| <i>Taq</i> DNA Polymerase with Standard <i>Taq</i> Buffer | 101227-642 | 400 units |
| | 101227-640 | 2,000 units |
| | 200067-188 | 4,000 units |
| <i>Taq</i> DNA Polymerase with Standard <i>Taq</i> (Mg-free) Buffer | 101447-600 | 400 units |
| | 101447-598 | 2,000 units |
| <i>Taq</i> PCR Kit | 101294-872 | 200 reactions |
| Quick-Load <i>Taq</i> 2X Master Mix | 200064-426 | 500 reactions |
| <i>Taq</i> 2X Master Mix | 101294-868 | 500 reactions |
| <i>Taq</i> 5X Master Mix | 101444-748 | 500 reactions |
| Multiplex PCR 5X Master Mix | 101446-500 | 100 reactions |

Selected Products for PCR (Cont.)

| PRODUCT | VWR CAT.NO. | SIZE |
|--|-------------|-----------------|
| DNA POLYMERASES (CONT'D) | | |
| Hot Start <i>Taq</i> DNA Polymerase | 102715-950 | 200 units |
| | 102715-948 | 1,000 units |
| Hot Start <i>Taq</i> 2X Master Mix | 102715-954 | 100 reactions |
| | 102715-952 | 500 reactions |
| Vent DNA Polymerase | 101228-292 | 200 units |
| | 101228-290 | 1,000 units |
| Vent (exo-) DNA Polymerase | 101228-300 | 200 units |
| | 101228-298 | 1,000 units |
| Deep Vent DNA Polymerase | 101227-258 | 200 units |
| | 101227-256 | 1,000 units |
| Deep Vent (exo-) DNA Polymerase | 101227-262 | 200 units |
| | 101227-260 | 1,000 units |
| LongAmp <i>Taq</i> DNA Polymerase | 101444-774 | 500 units |
| | 101444-772 | 2,500 units |
| LongAmp Hot Start <i>Taq</i> DNA Polymerase | 102500-072 | 500 units |
| | 102500-070 | 2,500 units |
| LongAmp <i>Taq</i> 2X Master Mix | 101444-754 | 100 reactions |
| | 101444-752 | 500 reactions |
| LongAmp Hot Start <i>Taq</i> 2X Master Mix | 102500-068 | 100 reactions |
| | 102500-066 | 500 reactions |
| LongAmp <i>Taq</i> PCR Kit | 101444-720 | 100 reactions |
| PCR CLONING & MUTAGENESIS | | |
| NEB PCR Cloning Kit | 102877-564 | 20 reactions |
| NEB PCR Cloning Kit (Without Competent Cells) | 102969-096 | 20 reactions |
| Q5 Site-Directed Mutagenesis Kit | 102855-188 | 10 reactions |
| Q5 Site-Directed Mutagenesis Kit (Without Competent Cells) | 102877-562 | 10 reactions |
| dNTPs | | |
| Deoxynucleotide (dNTP) Solution Set | 101227-334 | 25 µmol of each |
| Deoxynucleotide (dNTP) Solution Mix | 101228-418 | 8 µmol of each |
| | 101228-416 | 40 µmol of each |

Products for cDNA Synthesis

| PRODUCT | VWR CAT.NO. | SIZE |
|--|-------------|---------------|
| ProtoScript II First Strand cDNA Synthesis Kit | 102885-124 | 30 reactions |
| | 102885-126 | 150 reactions |
| ProtoScript First Strand cDNA Synthesis Kit | 101640-910 | 30 reactions |
| | 101640-908 | 150 reactions |
| ProtoScript II Reverse Transcriptase | 102500-080 | 4,000 units |
| | 102500-078 | 10,000 units |
| M-MuLV Reverse Transcriptase | 102500-082 | 40,000 units |
| | 101228-288 | 10,000 units |
| AMV Reverse Transcriptase | 101228-286 | 50,000 units |
| | 101417-976 | 200 units |
| | 101417-974 | 1,000 units |
| | 101417-978 | 500 units |



Products for Restriction Digestion

| PRODUCT | VWR CAT.NO. | SIZE |
|--|-------------|--------------|
| HIGH-FIDELITY (HF®) RESTRICTION ENZYMES | | |
| AgeI-HF | 101641-202 | 300 units |
| | 101641-200 | 1,500 units |
| ApoI-HF | 102969-126 | 1,000 units |
| | 102969-124 | 5,000 units |
| BamHI-HF | 101447-638 | 10,000 units |
| | 101641-128 | 50,000 units |
| | 101641-132 | 10,000 units |
| | 101641-130 | 50,000 units |
| BmtI-HF | 102500-026 | 300 units |
| | 102500-024 | 1,500 units |
| Bsal-HF | 101641-198 | 1,000 units |
| | 101641-196 | 5,000 units |
| BsrGI-HF | 102902-478 | 1,000 units |
| | 102902-476 | 5,000 units |
| BstEII-HF | 102500-042 | 2,000 units |
| | 102500-040 | 10,000 units |
| | 102500-044 | 10,000 units |
| BstZ171-HF | 103258-572 | 1,000 units |
| | 103258-570 | 5,000 units |
| DraIII-HF | 101641-194 | 1,000 units |
| | 101641-192 | 5,000 units |
| EagI-HF | 101447-648 | 500 units |
| | 101641-188 | 2,500 units |
| | 101641-190 | 2,500 units |
| EcoRI-HF | 101447-634 | 10,000 units |
| | 101641-102 | 50,000 units |
| | 101641-106 | 10,000 units |
| | 101641-104 | 50,000 units |
| EcoRV-HF | 101447-646 | 4,000 units |
| | 101641-178 | 20,000 units |
| | 101641-182 | 4,000 units |
| | 101641-180 | 20,000 units |
| HindIII-HF | 101641-112 | 10,000 units |
| | 101641-108 | 50,000 units |
| | 101641-114 | 10,000 units |
| | 101641-110 | 50,000 units |
| KpnI-HF | 101641-152 | 4,000 units |
| | 101641-148 | 20,000 units |
| | 101641-150 | 20,000 units |
| MfeI-HF | 101447-650 | 500 units |
| | 101447-204 | 2,500 units |
| MluI-HF | 102902-474 | 1,000 units |
| | 102902-472 | 5,000 units |
| NcoI-HF | 101447-644 | 1,000 units |
| | 101641-174 | 5,000 units |
| | 101641-176 | 5,000 units |
| NheI-HF | 101444-816 | 1,000 units |
| | 101641-120 | 5,000 units |
| | 101641-122 | 5,000 units |
| NotI-HF | 101447-642 | 500 units |
| | 101641-170 | 2,500 units |
| | 101641-172 | 2,500 units |
| NruI-HF | 102902-470 | 1,000 units |
| | 102902-468 | 5,000 units |
| NsiI-HF | 102902-466 | 1,000 units |
| | 102902-464 | 5,000 units |
| PstI-HF | 101641-144 | 10,000 units |
| | 101641-140 | 50,000 units |
| | 101641-146 | 10,000 units |
| | 101641-142 | 50,000 units |

Products for Restriction Digestion (Cont.)

| PRODUCT | VWR CAT.NO. | SIZE |
|--|-------------|--------------|
| PvuI-HF | 101641-156 | 500 units |
| | 101641-154 | 2,500 units |
| PvuII-HF | 101444-820 | 5,000 units |
| | 101641-158 | 25,000 units |
| | 101641-160 | 25,000 units |
| Sacl-HF | 101447-640 | 2,000 units |
| | 101641-162 | 10,000 units |
| | 101641-164 | 10,000 units |
| Sall-HF | 101444-818 | 2,000 units |
| | 101641-134 | 10,000 units |
| | 101641-138 | 2,000 units |
| SbfI-HF | 101641-136 | 10,000 units |
| | 101447-652 | 500 units |
| Scal-HF | 101641-206 | 2,500 units |
| | 101444-814 | 1,000 units |
| SpeI-HF | 101641-116 | 5,000 units |
| | 101641-118 | 5,000 units |
| | 102500-164 | 500 units |
| SphI-HF | 102500-162 | 2,500 units |
| | 102500-166 | 2,500 units |
| | 101444-822 | 500 units |
| SspI-HF | 101641-166 | 2,500 units |
| | 101641-168 | 2,500 units |
| | 101447-636 | 1,000 units |
| StyI-HF | 101641-124 | 5,000 units |
| | 101641-126 | 5,000 units |
| | 101641-186 | 3,000 units |
| | 101641-184 | 15,000 units |
| OTHER POPULAR RESTRICTION ENZYMES | | |
| AclI | 101229-226 | 500 units |
| | 101229-224 | 2,500 units |
| AvtI | 101228-918 | 100 units |
| | 101228-916 | 500 units |
| BglII | 101228-786 | 2,000 units |
| | 101228-782 | 10,000 units |
| | 101228-784 | 10,000 units |
| BsaI | 101229-150 | 1,000 units |
| | 101229-148 | 5,000 units |
| BsmBI | 101229-308 | 200 units |
| | 101229-306 | 1,000 units |
| DpnI | 101228-926 | 1,000 units |
| | 101228-924 | 5,000 units |
| MluI | 101229-010 | 1,000 units |
| | 101229-008 | 5,000 units |
| NcoI | 101228-984 | 1,000 units |
| | 101228-980 | 5,000 units |
| | 101228-986 | 1,000 units |
| NdeI | 101228-982 | 5,000 units |
| | 101228-646 | 4,000 units |
| NheI | 101228-644 | 20,000 units |
| | 101228-714 | 1,000 units |
| | 101228-710 | 5,000 units |
| PacI | 101228712 | 5,000 units |
| | 101229-194 | 250 units |
| PmeI | 101229-192 | 1,250 units |
| | 101229-234 | 500 units |
| SmaI | 101229-232 | 2,500 units |
| | 101228-770 | 2,000 units |
| SpeI | 101228-768 | 10,000 units |
| | 101228-726 | 500 units |
| | 101228-722 | 2,500 units |
| | 101228-724 | 2,500 units |



Products for Restriction Digestion (Cont.)

| PRODUCT | VWR CAT.NO. | SIZE |
|--------------------------------------|---------------------------------|--------------|
| XhoI | 101228-798 | 5,000 units |
| | 101228-796 | 25,000 units |
| | 101418-166 | 25,000 units |
| XbaI | 101228-792 | 3,000 units |
| | 101228-788 | 15,000 units |
| | 101228-794 | 3,000 units |
| | 101228-790 | 15,000 units |
| | 101228-940 | 500 units |
| XmaI | 101228-936 | 2,500 units |
| | 101228-938 | 2,500 units |
| | FEATURED GEL LOADING DYE | |
| Gel Loading Dye, Purple (6X) | 102877-610 | 4 ml |
| Gel Loading Dye, Purple (6X), no SDS | 102877-816 | 4 ml |

Products for End Modification

| PRODUCT | VWR CAT.NO. | SIZE |
|---|-------------|---------------|
| Quick Dephosphorylation Kit | 103258-516 | 100 reactions |
| | 103258-556 | 500 reactions |
| Shrimp Alkaline Phosphatase (Recombinant) | 102855-172 | 500 units |
| | 102855-174 | 2,500 units |
| Antarctic Phosphatase | 101228-350 | 1,000 units |
| | 101228-348 | 5,000 units |
| Alkaline Phosphatase, Calf Intestinal (CIP) | 101228-354 | 1,000 units |
| | 101228-352 | 5,000 units |
| T4 DNA Polymerase | 101228-186 | 150 units |
| | 101228-184 | 750 units |
| DNA Polymerase I, Large (Klenow) Fragment | 101228-212 | 200 units |
| | 101228-208 | 1,000 units |
| | 101228-210 | 1,000 units |
| Quick Blunting Kit | 101417-942 | 20 reactions |
| | 101417-940 | 100 reactions |
| Mung Bean Nuclease | 101228-276 | 1,000 units |
| | 101228-274 | 5,000 units |
| T4 Polynucleotide Kinase | 101228-174 | 500 units |
| | 101228-172 | 2,500 units |
| Klenow Fragment (3' → 5' exo) | 101228-222 | 200 units |
| | 101228-218 | 1,000 units |
| | 101228-220 | 1,000 units |
| β-Agarase I | 101228-396 | 100 units |
| | 101228-394 | 500 units |

Products for Ligation

| PRODUCT | VWR CAT.NO. | SIZE |
|--------------------------------------|-------------|---------------|
| Blunt/TA Ligase Master Mix | 102715-940 | 50 reactions |
| | 102715-938 | 250 reactions |
| | 102715-946 | 50 reactions |
| Instant Sticky-End Ligase Master Mix | 102715-944 | 250 reactions |
| | 102715-942 | 50 reactions |
| ElectroLigase | 101228-180 | 20,000 units |
| | 101228-176 | 100,000 units |
| | 101228-182 | 20,000 units |
| | 101228-178 | 100,000 units |
| | 101227-656 | 30 reactions |
| Quick Ligation Kit | 101227-654 | 150 units |
| | 102500-168 | 100,000 units |
| T3 DNA Ligase | 102500-170 | 750,000 units |
| | 102500-172 | 100,000 units |
| T7 DNA Ligase | 102500-174 | 750,000 units |
| | 101228-154 | 2,000 units |
| Taq DNA Ligase | 101228-152 | 10,000 units |

Products for Transformation

| PRODUCT | VWR CAT.NO. | SIZE |
|---|-------------|---------------------|
| <i>dam-/dcm-</i> Competent <i>E. coli</i> | 200064-394 | 20 x 0.05 ml/tube |
| | 101417-894 | 6 x 0.2 ml ml/tube |
| | 101417-898 | 20 x 0.05 ml/tube/ |
| NEB 5-alpha Competent <i>E. coli</i> (High Efficiency) | 101417-900 | 6 x 0.2 ml/tube/ |
| | 102855-190 | 1 x 96 well plate/ |
| | 103218-956 | 1 x 384 well plate/ |
| | 103218-958 | 12 x 8 tube strips |
| NEB 5-alpha Competent <i>E. coli</i> (Subcloning Efficiency) | 200067-176 | 6 x 0.4 ml/tube |
| NEB 5-alpha Electrocompetent <i>E. coli</i> | 101417-902 | 6 x 0.1 ml/tube |
| NEB 5-alpha F' <i>lacZ</i> Competent <i>E. coli</i> (High Efficiency) | 200067-178 | 20x 0.05/ |
| | 101417-904 | 6x 0.2 ml |
| NEB 10-beta Competent <i>E. coli</i> (High Efficiency) | 200067-180 | 20 x 0.05 ml/tube/ |
| | 101417-920 | 6 x 0.2 ml ml/tube |
| NEB 10-beta Electrocompetent <i>E. coli</i> | 101417-922 | 6 x 0.1 ml/tube |
| NEB Turbo Competent <i>E. coli</i> (High Efficiency) | 200064-396 | 20 x 0.05 ml/tube/ |
| | 101417-896 | 6 x 0.2 ml/tube |
| NEB Turbo Electrocompetent <i>E. coli</i> | 101447-584 | 6 x 0.1 ml/tube |
| NEB Stable Competent <i>E. coli</i> | 102877-612 | 20 x 0.5 ml/tube/ |
| | 102877-560 | 6 x 0.1 ml/tube |

Products for Nucleic Acid Purification

| PRODUCT | VWR CAT.NO. | SIZE |
|--------------------------------------|-------------|-----------|
| Monarch Plasmid Miniprep Kit | 102971-698 | 50 preps |
| | 102971-696 | 250 preps |
| Monarch DNA Gel Extraction Kit | 102971-670 | 50 preps |
| | 102971-668 | 250 preps |
| Monarch PCR & DNA Cleanup Kit (5 µg) | 102971-674 | 50 preps |
| | 102971-672 | 250 preps |

Products for DNA Analysis

| PRODUCT | VWR CAT.NO. | SIZE |
|-------------------------------------|-------------|-----------------|
| 1 kb DNA Ladder | 101228-494 | 200 gel lanes |
| | 101228-492 | 1,000 gel lanes |
| TriDye 1 kb DNA Ladder | 101228-508 | 125 gel lanes |
| Quick-Load 1 kb DNA Ladder | 101228-426 | 125 gel lanes |
| | 101449-372 | 375 gel lanes |
| 100 bp DNA Ladder | 101228-490 | 100 gel lanes |
| | 101228-488 | 500 gel lanes |
| TriDye 100 bp DNA Ladder | 101227-144 | 125 gel lanes |
| Quick-Load 100 bp DNA Ladder | 101227-136 | 125 gel lanes |
| | 101449-370 | 375 gel lanes |
| 2-Log DNA Ladder (0.1 - 10.0 kb) | 101228-486 | 200 gel lanes |
| | 101228-484 | 1,000 gel lanes |
| TriDye 2-Log DNA Ladder | 101227-142 | 250 gel lanes |
| Quick-Load 2-Log DNA Ladder | 101228-484 | 250 gel lanes |
| Quick-Load Purple 2-Log DNA Ladder | 102877-598 | 250 gel lanes |
| | 103258-568 | 750 gel lanes |
| 50 bp DNA Ladder | 101228-506 | 200 gel lanes |
| | 101228-504 | 1,000 gel lanes |
| Quick-Load Purple 50 bp DNA Ladder | 103218-930 | 250 gel lanes |
| Quick-Load 1 kb Extend DNA Ladder | 102500-150 | 125 gel lanes |
| Quick-Load Purple 1 kb DNA Ladder | 102877-812 | 1.25 ml |
| Quick-Load Purple 100 bp DNA Ladder | 102877-810 | 1.25 ml |
| Low Molecular Weight DNA Ladder | 101228-498 | 100 gel lanes |
| | 101228-496 | 500 gel lanes |



Products for DNA Analysis (Cont.)

| PRODUCT | VWR CAT.NO. | SIZE |
|---|-------------|---------------|
| Quick-Load Purple Low Molecular Weight DNA Ladder | 103218-932 | 125 gel lanes |
| Fast DNA Ladder | 101641-050 | 200 gel lanes |
| PCR Marker | 101228-502 | 100 gel lanes |
| | 101228-500 | 500 gel lanes |

Products for Seamless Cloning

| PRODUCT | VWR CAT.NO. | SIZE |
|--|-------------|---------------|
| NEBuilder HiFi DNA Assembly Cloning Kit | 102877-808 | 10 reactions |
| NEBuilder HiFi DNA Assembly Master Mix | 102877-804 | 10 reactions |
| | 102877-802 | 50 reactions |
| NEBuilder HiFi DNA Assembly Bundle for Large Fragments | 103218-906 | 20 reactions |
| Gibson Assembly Cloning Kit | 102715-912 | 10 reactions |
| | 102500-052 | 10 reactions |
| Gibson Assembly Master Mix | 102500-054 | 50 reactions |
| | 102902-438 | 15 reactions |
| NEB Golden Gate Assembly Mix | 102902-438 | 15 reactions |
| BioBrick® Assembly Kit | 101449-400 | 50 reactions |
| BbsI | 101229-162 | 300 units |
| | 101229-160 | 1,500 units |
| BsaI | 101229-150 | 1,000 units |
| | 101229-148 | 5,000 units |
| BsaI-HF | 101641-198 | 1,000 units |
| | 101641-196 | 5,000 units |
| BsmBI | 101229-308 | 200 units |
| | 101229-306 | 1,000 units |
| T4 DNA Polymerase | 101228-186 | 150 units |
| | 101228-184 | 750 units |
| Taq DNA Ligase | 101228-154 | 2,000 units |
| | 101228-152 | 10,000 units |
| T4 DNA Ligase | 101228-180 | 20,000 units |
| | 101228-176 | 100,000 units |
| | 101228-182 | 20,000 units |
| | 101228-178 | 100,000 units |
| T5 Exonuclease | 101710-280 | 1,000 units |
| | 101710-278 | 5,000 units |
| USER™ Enzyme | 101228-400 | 50 units |
| | 101228-398 | 250 units |

Products for Recombinational Cloning

| PRODUCT | VWR CAT.NO. | SIZE |
|-----------------|-------------|-----------|
| Cre Recombinase | 101228-370 | 50 units |
| | 101228-368 | 250 units |
| | 101641-516 | 250 units |



One or more of these products are covered by patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

Gibson Assembly[®] is a registered trademark of Synthetic Genomics, Inc. Phusion[®] is a registered trademarks of Thermo Fisher Scientific. Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. SuperScript[®] II, Gateway[®] and Geneart[®] are registered trademark of Life Technologies, Inc. In Fusion[®] is a registered trademark of Clontech Laboratories, Inc. DH5[™] and DH10B[™] are trademarks of Invitrogen. ECHO[™] is a trademark for Life Technologies, Inc. BIOBRICK[®] is a trademark of The BioBricks Foundation.

VWR  **1.800.932.5000**
We Enable Science **VWR.COM**



Prices and product details are current when published; subject to change without notice. [Certain products may be limited by federal, state, provincial, or local regulations.] VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. [VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s).] Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

