

Amersham ECL Anti-rabbit IgG, Horseradish Peroxidase- Linked Species-Specific Whole Antibody (from donkey)

Product Specification Sheet

Code: NA934

Warning

For research use only.

Not recommended or intended for diagnosis of disease in humans or animals.

Do not use internally or externally in humans or animals.

Storage

Store at 2–8°C. Do not freeze. Under these conditions, the product is stable for at least 6 months from the date of despatch.

Expiry

See outer packaging.

Safety warnings and precautions

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

Components

Horseradish Peroxidase conjugated antibody is supplied in Phosphate Buffered Saline (Sodium Phosphate 0.1 M, NaCl 0.1 M) pH 7.5, containing 1% (w/v) Bovine Serum Albumin and an anti-microbial agent.

Description

Purification to ensure species-specificity

The antibody is prepared by hyper-immunizing donkeys with purified immunoglobulin fractions from normal rabbit serum to produce high affinity antibodies. The pooled antiserum is used to produce an immunoglobulin preparation which is then affinity adsorbed to remove cross-reacting antibodies towards rat, human and mouse immunoglobulins. These activities are thoroughly depleted to ensure species-specificity.

Finally, to select for specific binding to rabbit IgG, the antibodies are purified using an affinity column of rabbit IgG. After washing to remove non-specific serum components and low affinity antibodies, the species-specific antibodies are eluted using carefully selected, mild conditions which minimize aggregation and preserve immunological activity, yet which will elute high affinity antibodies.

Preparation of labelled antibody

The enzyme Horseradish Peroxidase is attached to the immunoglobulin molecules using an adaptation of the periodate oxidation technique (1). This method has been found not to affect the effective binding of the antibody to the antigen or the activity of the enzyme.

Quality control

For every batch of enzyme-linked antibody that is produced the antibody titre is determined in an ELISA. The substrate used for the peroxidase is 2,2'-Azinobis[3-Ethylbenzothiazoline Sulphonate, diammonium salt], ABTST[™].

Every batch is also QC tested in a Western blotting system. This is performed using Hybond[™] ECL[™] membrane containing serially diluted Beta-galactosidase protein and immunodetected with primary antibody Anti-Beta-galactosidase and secondary antibody NA934, anti-rabbit HRP. Blots are detected using ECL and ECL Plus[™] detection systems.

Applications

1. Protein blotting

a) Detection with ECL (2) Western blotting reagents

This reagent has been shown to be suitable for use in ECL Western blotting applications. The control system used was the detection of Anti-Beta Galactosidase.

We have found in our laboratories that dilutions of 1:5000 for Anti-Beta Galactosidase and 1:50 000 for NA934 are suitable for the detection of 6 ng of beta-galactosidase on Hybond ECL membrane, exposed to Hyperfilm[™] ECL for 5 minutes.

To achieve the same sensitivity level on Hybond-P PVDF, concentrations would typically be Anti-Beta Galactosidase - 1:5000 and NA934 - 1:100 000.

b) Detection with ECL Plus (3, 4) Western blotting reagents

ECL Plus Western blotting reagent is highly sensitive, giving an increase, for this antibody, of 4 to 20 fold over ECL detection.

This property can be utilized in 2 ways:

- Preservation of antibodies that are rare or costly
- Increase in detectable sensitivity levels

The control system used was the same as for ECL.

The suitable antibody dilutions, to detect 6 ng of Beta-Galactosidase on Hybond ECL membrane are Anti-Beta Galactosidase - 1:10 000 and NA934 - 1:100 000. For Hybond-P



PVDF antibody dilutions are typically Anti-Beta Galactosidase-1:20 000 and NA934 - 1:200 000.

c) Colorimetric detection

A dilution of 1:300 is recommended.

2. ELISA

If this reagent is to be used to detect rabbit immunoglobulins, we have found in our laboratories that a dilution of 1:9000 is suitable for the detection of 1 µg of IgG. For greater sensitivity (for example down to 300 pg) the reagent should be diluted rather less (for example 1:5000). A suitable diluent is Phosphate-Buffered Saline containing 0.05% (v/v) Tween™ 20.

3. Immunocytochemistry

When using the reagent as a second antibody in immunocytochemistry on sections of formalin-fixed wax-embedded tissue the antibody can be typically diluted 1:100 in Phosphate-Buffered Saline. The user may wish to adjust this to obtain the required sensitivity for the tissue under investigation. If frozen sections are used, acceptable staining may be obtained using even higher dilutions of the reagent.

Protocol recommendations

Membranes

Nitrocellulose and PVDF membranes are suitable for use with both detection systems. PVDF membrane is highly recommended for use with ECL Plus detection reagents.

For high quality results the following guidelines should be followed:

Blocking: Use enough blocking agent to block all non-specific sites. A typical block 5% non-fat dried milk in PBS Tween or TBS Tween. See 'Tech-Tips' No. 136 available from GE Healthcare for further details.

Washing: The volume of wash buffer (eg PBS-T or TBS-T) must be sufficient to cover the membrane completely.

Optimization of primary and secondary antibodies

ECL detection

ECL Western blotting is a very sensitive technique. As such it is essential to optimize the system under study for high signal and low background for both primary and secondary antibodies.

Dot blots are a quick and effective method of determining the optimum dilutions required for primary and secondary antibodies. Optimization details are set out in the RPN2106/2108/2109/2209/2134 booklets and 'Tech-Tips' No. 129 available from GE Healthcare.

ECL Plus detection

Due to the improved sensitivity of ECL Plus compared to ECL, optimization details as set out in the RPN2132/2133 booklets and 'Tech-Tips' No. 169 available from GE Healthcare are recommended.

Typical anti-mouse secondary antibody dilution ranges:

ECL for nitrocellulose membrane 1:5000 to 1:50 000
ECL Plus for nitrocellulose membrane 1:10 000 to 1:100 000

For PVDF membrane the use of higher dilutions may be necessary. The exact concentration of the secondary antibody will always be dependent upon the primary antibody used and the sensitivity and exposure times required.

Detection: Ensure any excess ECL or ECL Plus detection reagents are sufficiently drained prior to exposure.

Exposure times:

ECL - exposure times of 1 to 15 minutes are suggested.
ECL Plus - initial exposure times of 1 to 5 minutes are suggested.

Signal can still be obtained up to 24 hours after the application of ECL Plus reagents, and for this exposure times of 1 to 2 hours may be required.

Related products

ECL Western blotting detection reagents	RPN2106/2108/2109/2209/2134
ECL Plus Western blotting detection system	RPN2132/2133
Hybond ECL membrane	RPN2020D
Hybond-P PVDF membrane	RPN2020F
Hyperfilm ECL membrane	28-9068-35/28-9068-36/ 28-9068-37/28-9068-38/ 28-9068-39/28-9068-40/ 28-9068-41
ECL protein molecular weight markers	RPN2107
ECL Plus Western Blotting Reagent Pack	RPN2124
ECL Blocking Agent	RPN2125
Rainbow™ Molecular Weight Markers	RPN800E/756E/755E

References

1. NAKANE, P.K. and KAWAOI, A., *Journal of Histochemistry and Cytochemistry*, **22**, pp.1084-1091, 1974.
2. WHITEHEAD, T.P. *et al.*, *Clin. Chem.*, **25**, pp.1531-1546, 1979.
3. AKHAVEN-TAFTI, H. *et al.*, *Clin. Chem.*, **41**, pp.1368-1369, 1995.
4. AKHAVEN-TAFTI, H. *et al.*, *Biolom. And Chemilum. Fundamentals and Applied Aspects*, pp.199-202, Chichester, 1994.

Legal

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

Amersham, Hybond, ECL, ECL Plus and Hyperfilm are trademarks of GE Healthcare companies.

All third party trademarks remain the property of their respective owners.

ECL Plus Western blotting detection reagents are manufactured for GE Healthcare by Lumigen Inc. This component is covered by US patent numbers 5491072 and 5593845 and foreign equivalents and is sold under licence from Lumigen Inc.

This component is covered by US Patent Nos: 5491072 and 5593845 and is sold under licence from Lumigen Inc.

© 2006-2008 General Electric Company – All rights reserved
Previously published 2006

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited
Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA UK

GE Healthcare Bio-Sciences AB
Björkgatan 30 751 84,
Uppsala Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5 D-79111
Freiburg Germany

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue PO Box 1327
Piscataway NJ 08855-1327 USA

GE Healthcare Bio-Sciences KK
Sanken Bldg 3-25-1, Hyakunincho Shinjuku-ku
Tokyo 169-0073 Japan