

GE Healthcare



illustra NICK Columns

Gravity flow columns for the removal of unincorporated radiolabeled nucleotides from DNA labeling reactions

Product booklet

Codes: 17-0855-01 (20 purifications)
17-0855-02 (50 purifications)



Page finder

1. Legal	3
2. Handling	4
2.1. Safety warnings and precautions	4
2.2. Storage conditions	4
2.3. Expiry	4
3. Components	5
3.1. Kit contents	5
3.2. Materials to be supplied by user	5
3.3. Equipment needed	5
4. Description	6
4.1. Introduction	6
4.2. The basic principle	7
4.3. Product specifications	9
5. Protocol	11
5.1. Preparation of working solutions	11
5.2. Protocol for the removal of unincorporated radiolabeled nucleotides from DNA labeling reactions	12
6. Appendices	15
6.1. Related products	15
 Quick Reference Protocol Card	Back Cover
Tear off sheet containing protocols for the experienced user for the removal of unincorporated radiolabeled nucleotides from DNA labeling reactions	

1. Legal

Product use restriction

The **illustra™ NICK™ Columns** and components have been designed, developed and sold **for research purposes only**. They are suitable **for *in vitro* use only**. No claim or representation is intended for their use to identify any specific organism or for clinical use (diagnostic, prognostic, therapeutic, or blood banking).

It is the responsibility of the user to verify the use of the **illustra NICK Columns** for a specific application, as the performance characteristics of this product have not been verified for any specific organism.

GE, imagination at work and GE monogram are trademarks of General Electric Company.

NICK, MicroSpin, Rediprime, Sephadex, Ready-To-Go, ProbeQuant, Hybond, Redivue, Megaprime, Hyperfilm, Hypercassette and **illustra** are trademarks of GE Healthcare companies.

Kathon is a trademark of Rohm and Haas Company.

© 2006–2007 General Electric Company-All rights reserved.
Previously published 2006. Date of revision Apr. 2007.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available upon request. Contact your GE Healthcare representative for the most current information (see back cover for contact details).

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited.
Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA UK

2. Handling

2.1. Safety warnings and precautions

Warning: For research use only.

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

All chemicals should be considered as potentially hazardous. Only persons trained in laboratory techniques and familiar with the principles of good laboratory practice should handle these products. Suitable protective clothing such as laboratory overalls, safety glasses and gloves should be worn. Care should be taken to avoid contact with skin or eyes; if contact should occur, wash immediately with water. (see Material Safety Data Sheet(s) and/or Safety Statements(s) for specific recommendations).

2.2. Storage conditions

Store at ambient temperature (4°C–30°C). **Do not freeze.**

2.3. Expiry

For expiry date please refer to outer packaging label.

3. Components

3.1. Kit contents

Identification	Pack size Cat. No.	20 purifications 17-0855-01	50 purifications 17-0855-02
	illustra™ NICK™ Columns	20	50

Refer to the Certificate of Analysis for a complete list of kit components.

3.2. Materials to be supplied by user

Disposables:

DNase-free collection tubes

A suitable receptacle to catch waste buffer flow through

Chemicals:

Buffer 1

Any nuclease-free buffer is suitable, including water or Tris/EDTA (TE); both are available from GE Healthcare. 5 ml **Buffer 1** per purification is ample.

Buffer 2

This buffer solution should be the same as the one in which your un-purified probe is solved. This may be the same as **Buffer 1**, but does not have to be. **Buffer 2** is used for the Column Equilibration, Sample Application and Elution steps and it is important to use the same buffer for each of these steps. 5 ml **Buffer 2** per purification is ample.

3.3. Equipment needed

A column support

4. Description


4.1. Introduction

illustra NICK Columns are single-use disposable columns pre-packed with Sephadex™ G-50 DNA Grade resin and require only gravity to run. They allow DNA purification by the process of gel filtration. Molecules larger than the largest pores in the matrix are excluded from the matrix and elute first. Intermediate size molecules penetrate the matrix to varying extents, depending on their size. Penetration of the matrix retards progress through the column; very small molecules elute last. The volume required to elute these small molecules is dependent on the volume available both inside and outside the pores (i.e. the bed volume).

illustra NICK Columns are designed for the rapid and convenient separation of nick-translated DNA from unincorporated ^{32}P -labeled nucleotides and similar separations. They can be used for any DNA greater than 20 bases in length. They will not remove or denature enzyme.

4.2. The basic principle

Use of **illustra NICK Columns** Columns involves the following steps:



Step	Comments	Component
1. Column Preparation	The storage buffer is poured away and the top of the column rinsed with buffer chosen by user (Buffer 1).	
2. Column Equilibration	Column is equilibrated with Buffer 2 (buffer in which un-purified probe is solved).	
3. Sample Application	Sample is applied to column Additional Buffer 2 is applied to the column.	
4. Elution	Purified sample is eluted from the column with Buffer 2 .	

4.3. Product specifications

illustra NICK Columns are recommended for the removal of unincorporated radiolabeled nucleotides from DNA labeling reactions. A summary of the product specifications are given in Table 1 below:

Table 1. illustra NICK Column specifications

Sample Type:	Oligonucleotides and small DNA fragments > 20 bases in length
Principle	Gel filtration
Column matrix	Sephadex G-50 DNA grade
Input sample volume	1–100 µl
Column buffer	Distilled water containing 0.15 % Kathon™ CG/ICP Biocide as preservative
Yield/recovery of DNA	> 90%
Purity of recovered DNA	Typically < 0.2 % salt contamination
Length of labeled DNA recovered	> 20 bp (N.B. there is no maximum length of oligonucleotide that can be purified)
Volume of eluted purified sample	400 µl
Major subsequent applications	Hybridization
Gel bed dimensions	0.9 x 2.0 cm
Column capacity (maximum amount of DNA that can be loaded onto column)	100 µg Do not load a DNA solution of a concentration greater than 1 mg/ml. Higher concentrations reduce column resolution and give lower yield due to increased viscosity. Samples at a concentration greater than 1 mg/ml should be diluted with Buffer 1 prior to loading.

For the quantitative removal of unincorporated labeled nucleotides from a DNA labeling reaction using spin-column chromatography, we recommend use of illustra ProbeQuant™ G-50 Micro Columns. These columns are suitable for purification of labeled DNA greater than 20 bases in length, and are also suitable for purification of biotinylated probes. Please note that these columns are supplied in 150 mM STE buffer containing 0.05 % Kathon.

For purification of labeled DNA less than 20 bases in length, we recommend use of illustra MicroSpin™ G-25 Columns. These columns are designed for the rapid purification of DNA, and will remove unincorporated nucleotides from end-labeled oligonucleotides. They can be used for any DNA greater than 10 bases in length, and are therefore ideal for the purification of oligonucleotides or very small DNA fragments. They will not remove or denature enzyme. Please note that these columns are supplied in double-distilled water containing 0.05% Kathon.

Although tests have shown low levels of RNase activity for illustra NICK Columns, they are not specifically treated to be RNase free. Customers do regularly use illustra NICK columns for purification of RNA and are very satisfied. However, as the columns have not been tested for the absence of RNases we cannot guarantee that the RNA will not be degraded-but in principle it works well.

Biotinylated probes can be purified using illustra NICK Columns (as it is the size of the probe and not the modification that is relevant).

5. Protocol

Use of icons

The Key below describes the purpose of the icons used throughout the protocol booklet.



This icon is used to highlight particularly critical steps within the protocol that must be adhered to. If this advice is not followed it will have a detrimental impact on results.



This icon is used to highlight technical tips that will enhance the description of the step. These tips may indicate areas of flexibility in the protocol or give a recommendation to obtain optimum performance of the kit.

5.1. Preparation of working solutions

See section 3.2 and 3.3 for Materials & Equipment to be supplied by user.

Buffer 1

Any nuclease-free buffer is suitable, including water or Tris/EDTA (TE); both are available from GE Healthcare. 5 ml **Buffer 1** per purification is ample.

Buffer 2

This buffer solution should be the same as the one in which your unpurified probe is solved. This may be the same as **Buffer 1**, but does not have to be. It is important to use the same buffer in the Column Equilibration, Sample Application and Elution steps. 5 ml **Buffer 2** per purification is ample.

5.2 Protocol for the removal of unincorporated radiolabeled nucleotides from DNA labeling reactions



Note: Prior to commencing, ensure that the NICK Columns have equilibrated to room temperature (20–25°C).

1. Column Preparation

a. Remove the **top** cap from the NICK Column and pour off the excess storage liquid.

b. Rinse the top of the column once with **Buffer 1** and pour off the excess.



Note: Any nuclease-free buffer is suitable, including water or Tris/EDTA (TE); both are available from GE Healthcare.

c. Remove the **bottom** cap. Support the column over a suitable waste receptacle.



Rinse top of column with Buffer 1. Pour off excess.



2. Column Equilibration

a. Equilibrate the column with 3 ml of **Buffer 2**.




Note: **Buffer 2** should be the same as the one in which your un-purified probe is solved. This may be the same as **Buffer 1**, but does not have to be.



Note: This volume corresponds to 1 complete refill of the column.

b. Allow **Buffer 2** to completely enter the gel bed by gravity flow. Do not apply positive pressure.




3 ml Buffer 2




3. Sample Application

- a. Add the sample to the top-centre of the resin column in any volume ranging from 1–100 μl . Allow the sample to enter the gel bed completely.



Note: Concentration of the DNA sample should be less than 1 mg/ml, as higher concentrations tend to reduce resolution and give lower yields due to increased viscosity. If necessary, dilute sample with **Buffer 1** prior to loading, but do not exceed maximum sample volume of 100 μl per column.


- b. Add 400 μl **Buffer 2** to the column and allow to enter the gel bed completely.



4. Elution

- a. Place an appropriately sized collection tube under the column.

- b. Elute the purified sample with 400 μl **Buffer 2**.



- c. Store the purified sample at -20°C .



6. Appendix

6.1. Related products

A full range of Molecular Biology reagents can be found on the GE Healthcare web site and in the catalog <http://www.gehealthcare.com/illustra>

A full range of Detection Products and available pack sizes can be found on the GE Healthcare web site and in the catalog <http://www.gehealthcare.com/newhyperfilm>

If you need further information, GE technical services are happy to assist (world-wide phone numbers can be found on the back cover).

Application	Product	Product code	Pack size
Buffer 1	TE Buffer, 50 ×	US75834	100 ml
	Water, nuclease-free	US70783	500 ml
Blotting	Hybond™-N+ (82 mm)	RPN82B	50 discs
	Hybond-N+ (15 × 20 cm)	RPN1520B	10 sheets
	Hybond-NX (82 mm)	RPN82T	50 discs
	Hybond-NX (15 × 20 cm)	RPN1520T	10 sheets
	Hybond-N (82 mm)	RPN82N	50 discs
	Hybond-N (15 × 20 cm)	RPN1520N	10 sheets

Application	Product	Product code	Pack size
	Hybond-XL (82 mm)	RPN82S	50 discs
	Hybond-XL (15 × 20 cm)	RPN1520S	10 sheets
	Hybond blotting paper (20 × 20 cm)	RPN6101M	100 sheets
Radioactive labeling	Redivue™ nucleotides	AA0085	250 UCI – 1 MCI
	Rediprime™ II DNA Labeling System	RPN1633	30 reactions
	Ready-To-Go™ DNA Labeling Beads (-dCTP)	27-9240-01	1 kit
	Megaprime™ DNA Labeling System, dNTP	RPN1604	30 reactions
	Megaprime DNA Labeling System, dCTP	RPN1606	30 reactions
	Nick Translation Kit, dNTP	N5000	20 reactions
	5'-End Labeling Kit	RPN1509	20 reactions
Detection	Hyperfilm™ MP (18 × 24 cm)	28-9068-43	50 sheets
	Hyperfilm MP Enveloped (18 × 24 cm)	28-9068-50	50 sheets

Application	Product	Product code	Pack size
	Hypercassette™	RPN11642	1
Purification	illustra MicroSpin G-50 Columns	27-5330-01	50 purifications
	illustra MicroSpin G-25 Columns	27-5325-01	50 purifications
	illustra ProbeQuant G-50 MicroColumns	28-9034-08	50 purifications

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30 751 84
Uppsala
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5 D-79111
Freiburg
Germany

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire
HP7 9NA
UK

GE Healthcare Bio-Sciences
Corp
800 Centennial Avenue
P.O. Box 1327
Piscataway
NJ 08855-1327
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1
Hyakunincho Shinjuku-ku
Tokyo 169-0073
Japan

**GE Healthcare
regional office
contact numbers:**

Asia Pacific:
Tel: +65 6275 1830
Fax: +65 6275 1829

Australasia
Tel: +61 2 9899 0999
Fax: +61 2 9899 7511

Austria
Tel: 01/57606-1619
Fax: 01/57606-1627

Belgium
Tel: 0800 73 888
Fax: 02 416 82 06

Canada
Tel: 1 800 463 5800
Fax: 1 800 567 1008

**Central, East, & South
East Europe**
Tel: +43 1 972720
Fax: +43 1 97272 2750

Denmark
Tel: 45 16 2400
Fax: 45 16 2424

Eire
Tel: 01494 544000
Fax: 0044 1494 542010

Finland & Baltics
Tel: +358-(0)9-512 39 40
Fax: +358 (0)9 512 39 439

France
Tel: 01 6935 6700
Fax: 01 6941 9677

Germany
Tel: (089) 96281 660
Fax: (089) 96281 620

Greater China
Tel: +852 2100 6300
Fax: +852 2100 6338

Italy
Tel: 02 27322 1
Fax: 02 27302 212

Japan
Tel: +81 3 5331 9336
Fax: +81 3 5331 9370

Latin America
Tel: +55 11 3933 7300
Fax: +55 11 3933 7304

Middle East & Africa
Tel: +30 210 9600 687
Fax: +30 210 9600 693

Netherlands
Tel: 0800 82 82 82 1
Fax: 0800 82 82 82 4

Norway
Tel: 815 65 555
Fax: 815 65 666

Portugal
Tel: 21 417 7035
Fax: 21 417 3184

Russia & other C.I.S. & N.I.S
Tel: +7 (495) 956 5177
Fax: +7 (495) 956 5176

Spain
Tel: 93 594 49 50
Fax: 93 594 49 55

Sweden
Tel: 018 612 1990
Fax: 018 612 1910

Switzerland
Tel: 0848 8028 12
Fax: 0848 8028 13

UK
Tel: 0800 616928
Fax: 0800 616927

USA
Tel: +1 800 526 3593
Fax: +1 877 295 8102


<http://www.gehealthcare.com/illustra>

GE Healthcare UK Limited
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA
UK



imagination at work

The next two pages are a
protocol card.
If required please add to the
back page as a tear off addition



If not then delete these three
pages.


Quick Reference Protocol Card

illustra™ NICK™ Columns

17-0855-01 (20 purifications)
17-0855-02 (50 purifications)

A. Protocol for the removal of unincorporated radiolabeled nucleotides from DNA labeling reactions

- Ensure suitable **Buffer 1** is available
- Ensure suitable **Buffer 2** is available



: Add

1. Column preparation

- Remove top cap
- Rinse top of column with **Buffer 1** and pour off excess
- Remove bottom cap




2. Column equilibration

- ➕ 3 ml **Buffer 2**
- Allow to completely enter gel bed by gravity flow




3. Sample application

- ➕ 1–100 µl sample
- Allow to completely enter gel bed by gravity flow
- ➕ 400 µl **Buffer 2**
- Allow to completely enter gel bed by gravity flow



4. Elution

- Place an appropriate collection tube under column
-  400 µl Buffer 2
- Collect eluate by gravity flow
- Store purified sample at -20°C



GE, imagination at work and GE monogram are trademarks of General Electric Company.
illustra and NICK are trademarks of GE Healthcare Companies.

© 2006-2007 General Electric Company-All rights reserved, previously published 2006.
All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available upon request.

Contact your GE Healthcare representative for the most current information.

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited,
Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA UK



imagination at work

17-0855-01PC Rev-B 06/2007