product codes:

RPN 8510, RPN 8511

# Nucleon PhytoPure

#### DNA from plant tissue

Nucleon<sup>TM</sup> PhytoPure<sup>TM</sup> system can produce excellent yields of high quality DNA in a fraction of the time taken by conventional methods. Not only is it more efficient but it is considerably simpler. Nucleon Phytopure is:

- Novel: utilizes a revolutionary proprietary resin incorporating borate chemistries to ensure polysaccharide-free DNA preparations
- Fast: enables extraction of DNA in less than 1 hour
- Safe: eliminates the need to use phenol
- Simple: easy-to-use protocol requires only one centrifugation step prior to DNA precipitation
- Pure: DNA is of high quality and suitable for RFLP, RAPD and AFLP analyses

Polysaccharides are common contaminants in plant DNA extracts and can inhibit further enzymatic analysis of DNA. Nucleon PhytoPure DNA extraction system has been developed specifically to solve this problem.

After breaking of the cell wall, the cells are lysed in reagents containing potassium SDS which is known to complex with proteins and polysaccharides. Chloroform is added along with the Nucleon PhytoPure proprietary resin (Fig 1).

This resin contains free boric acid (B(OH)<sub>2</sub>) groups which covalently bind polysaccharides (Fig 2), thus removing them from the sample. The resin forms a semi-solid stratum during partitioning with chloroform facilitating DNA recovery, ensuring good recovery of high quality DNA (Fig 3).

Nucleon PhytoPure has been used successfully on a wide range of plant species (Table 1). DNA can be extracted from fresh, frozen or freeze dried material. The kit is available in two sizes and the protocols can accommodate large and small sample preparations.




Fig 1. Principle steps in the Nucleon Phytopure protocol




Fig 2. Polysaccharide binding mechanism




Fig 3. Formation of semi-solid stratum after addition of PhytoPure resin and chloroform

Arabidopsis	Brassica oleracea	Brassica napus	Capsicum annuum
Capsicum frutescens	Cereals (barley, maize, rye, wheat)	Cocos nucifera	Helianthus annus
Helianthus tuberosus	Hevea braziliensis	Irvingia gabonensis	Lotus japonicus
Lycopersicon esculentum	Musa spp	Malus spp	Nicotiana
Phaseolus vulgaris	Pinus sylvestris	Pisum sativum	Rhododendron spp
Salix spp	Solanum tuberosum	Swietenia macrophylla	Sphagnum, bog moss

Table 1. Plant species from which DNA has been successfully extracted using PhytoPure



Research carried out at IARC - Long Ashton, UK compared Nucleon PhytoPure against a standard SDS/phenol extraction method using fresh leaf material from four plant species (willow, maize, rye and rhododendron). The DNA extracted was of sufficiently high purity to use in restriction digests, RAPD analysis (Fig 4) and AFLP analysis (Fig 5). High quality DNA was recovered without phenol and in less than half the time taken by the conventional method.




Fig 4. RAPD analysis of 4 species using operon primer OPA 10. DNA diluted to approx. 25 ng per track and run on a 1.2% agarose in TBE.

Tracks 1 & 10: 1kb ladder

Track 2: willow DNA extracted by Nucleon PhytoPure

Track 3: willow DNA extracted by SDS/phenol

Track 4: maize DNA extracted by Nucleon PhytoPure

Track 5: maize DNA extracted by SDS/phenol

Track 6: rye DNA extracted by


Nucleon PhytoPure Track 7: rye DNA extracted by

SDS/phenol

Track 8: rhododendron DNA extracted by Nucleon PhytoPure

Track 9: rhododendron DNA extracted by SDS/phenol

Note the failure of the amplification reactions with willow DNA extracted by the SDS/phenol method.



Pst | 1 & Mse | 2

Fig 5. AFLP analysis of DNA from four plant species using two primer sets.

First primer set: Pst I.1 & Mse I.2 Second primer set: Pst 1.2 & Mse 1.2

Track 1: willow DNA extracted by Nucleon PhytoPure

Track 2: willow DNA extracted by SDS/phenol

Track 3: maize DNA extracted by Nucleon PhytoPure

Track 4: maize DNA extracted by SDS/phenol

Track 5: rye DNA extracted by Nucleon PhytoPure

Track 6: rye DNA extracted by

SDS/phenol Track 7: rhododendron DNA

extracted by Nucleon PhytoPure Track 8: rhododendron DNA extracted by SDS/phenol

AFLP's run on a 6% acrylamide gel. Mse 1.2 primer was labelled with 33P.

Note the complete failure of the AFLP analysis in track 2 when SDS/phenol was used to extract the DNA

RPN 8502

## Ordering Information

Nucleon PhytoPure for plant DNA extraction 50 preparations of 0.1 g	RPN 8510	
<b>Nucleon PhytoPure</b> for plant DNA extraction 50 preparations of 1.0 g	RPN 8511	
Kits also available :		
Nucleon MiY for yeast mini preps DNA extraction 75 preparations	RPN 8518	
Nucleon HT for hard tissue & paraffin sections 50 preps of up to 25 mg per prep	RPN 8509	

### **Nucleon BACC1**

50 preps of 1ml whole blood or cultured cells x 10<sup>6</sup> to 3 x 10<sup>6</sup>) RPN 8501

#### cleon BACC2

50 preps of 10 ml whole blood or cultured cells x 10<sup>6</sup> to 1 x 10<sup>7</sup>)

100 preps of 1 ml whole blood and cultured cells x 10<sup>6</sup> to 3 x 10<sup>6</sup>)

#### cleon BACC3

50 preps of 10 ml whole blood RPN 8512 s kit contains all reagents necessary for 50 x 10 ml blood extractions)

## Product information

Asia Pacific Tet: +852 2811 8693 Fax: +852 2811 5251 Australasia Tet: +61 2 9894 5152 Fax: +61 2 9899 7511 Austria Tet: 01 576 0616 25 Fax: 01 576 0616 27 Belgium Tet: 0800 73 888 Fax: 03 272 1637 Canada Tel: 1 800 463 8800 Fax: 1 800 567 1008 Central, East, and South East Europe Tel: 43 1 982 3826 Fax: +43 1 985 8327 Denmark Tel: 45 16 2400 Fax: 45 16 2424 Finland Tel: 09 512 3940 Fax: 09 512 1710 France Tel: 01 6935 6700 Fax: 01 6941 9677 Germany Tel: 0761 4903 106 Fax: 0761 4903 405 Italy Tel: 02 27322 1 Fax: 02 27302 212 Japan Tel: 81 3 5331 9334 Fax: 81 3 5331 9369 Latin America Tel: +55 11 3667 5700 Fax: 0165 580 410 Fax: 09 Middle East and Africa Tel: +230 (1) 96 00 697 Fax: -30 (1) 96 00 693 Netherlands Tel: 01 6935 80 410 Norway Tel: 2318 5800 Fax: 60 3 724 2090 Spain Tel: 93 594 49 50 Fax: 93 594 49 55 Sweden Tel: 018 612 0000 Fax: 018 612 1200 Switzerland Tel: 01 802 81 50 Fax: 01 802 81 51 UK Tel: 0800 616928 Fax: 0800 616927 USA Tel: +1 800 526 3593 Fax: +1 877 295 8102

Manufactured for Amersham Biosciences by Teonel Life Sciences plc. Nucleon and PhytoPure are trademarks of Teonel Life Sciences plc. \*The Polymerase Chain Reaction (PCR) is covered by patents owned by Roche Molecular Systems Manuactured for Amersham Biosciences by Lepnel Life Sciences pic. Nucleon and PhytoPure are trademarks of Lepnel Life Sciences pic. \*1 Ine Polymerase Chain Reaction (PCR) is covered by patients owned by Roche Molecular Systems and F Hoffmann-La Roche Lifd. A license be to use the PCR process for certain research and development activities accompanies the purchase of certain in reagents from licensed suppliers such as Amersham Biosciences Limited and affiliates when used in conjunction with an authorized thermal cycler. Amersham is a trademark of Nycomed Amersham pic. Amersham Biosciences UK Limited,
Amersham Place, Little Chalfont, Buckinghamshire HP7 9NA, England. Amersham Biosciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden. Amersham Biosciences Inc., 800 Centennial Avenue, PO Box 1327, Piscataway, NJ 08855, USA Amersham Biosciences Biosciences September 1, Freiburg, German, All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Biosciences group which supplies them. A copy of these terms and conditions is available on request. © Amersham Biosciences UK Limited 2000 - All rights reserved.

