

VWRbioMarke™

DECODING THE GENE

VWR is now part of Avantor
Learn more on page 31

| Page 4 Rapid detection of gene-editing events and genotyping | Page 14 Capture your genetic snapshot | Page 21 Simplifying immune oncology research | Page 24 Efficient transient gene expression in neurosphere | Page 32 High-yield, scalable library preparation |



DECODING THE GENE

Sometimes our greatest advancements come from the most unlikely of places. Such is the case with Genetics, which finds its roots in the pea plants of the Austrian Scientist Gregor Mendel. Mendel's work on the inheritance patterns observed in *Pisum sativum* went largely unnoticed in his own time but would ultimately serve as the basis for the modern field of genetics.

In the century and a half since Mendel published his work, we have seen our knowledge grow by leaps and bounds; highlighted by Watson and Crick's landmark discovery of the double helix structure of DNA and the completion of the Human Genome. Today we find ourselves at the dawn of a new era. One in which we have the capacity to sequence a genome in hours and the potential to engineer that sequence to our liking. The potential applications of our Decoding the Genome is revolutionizing all aspects of biological research and medicine. To think, it all started with the humble pea.



BIOMARKERS



Biomarker Service Offerings



OnTarget ctDNA Analysis

Detect and monitor circulating tumor DNA in blood with OnTarget enrichment technology.

Biomarker Discovery

Measure the progression of diseases, the effectiveness of therapies or the diagnosis of conditions.

Liquid Biopsy

Test and detect either circulating tumor cell-associated biomarkers or circulating tumor DNA.

Clinical Chemistry


Measure analytes in serum, plasma, cells, tissue lysates, or other biological materials.

Multiple Biomarker Types

- Genomics
- Transcriptomics
- Proteomics
- Metabolomics
- Lipidomics
- And many more!

Sample Suppliers

- Covance
- Eurofins-Cerep-Panlabs
- CiToxLab
- Metabolon
- Human Metabolome Tech
- Personal Genome Diagnostics
- And many more!



For more info visit VWR.com/VWRSciencePortal or email VWR@scientist.com

This service is not available in Canada. Please contact your VWR Sales Rep to learn about similar options available in your region.



CONTENTS

Rapid Detection of Gene-Editing Events and Genotyping Using DNA Isolated with QuickExtract™ DNA Extraction Solution 4

Demonstration of FlashGel™ Recovery with Sheared DNA..... 10

Capture Your Genetic Snapshot 14

How Can a Chromatography Bead Free Up Your Time?..... 19

Simplyfing Immune Oncology Research 21

Efficient Transient Gene Expression in Neurosphere with jetMESSENGER™ mRNA Transfection Reagent..... 24

High-Yield, Scalable Library Preparation with the NEBNext® Ultra™ II FS DNA Library Prep Kit..... 32

SUPPLIER	PAGE	SUPPLIER	PAGE	SUPPLIER	PAGE
Agilent Technologies.....	7, 13	GE Healthcare	19	PeptoTech.....	18
Analytik Jena.....	9	Greiner Bio-One	28	Polyplus-transfection	24–26
Avantor.....	31	Hardy Diagnostics	29	Quantabio.....	8, 30
Azure Biosystems.....	20	Lonza	10–12, 13	Scientist.com	2, 23
Beckman Coulter.....	27	Lucigen	4–7	StemBioSys.....	16, 17
Biotium.....	8, 30	Mirus Bio	17	VWR Life Science Seradigm	18
Corning	21–22	New England BioLabs	12, 16, 32–35	Zymo Research.....	14–15

Rapid Detection of Gene-Editing Events and Genotyping Using DNA Isolated with QuickExtract™ DNA Extraction Solution

By Ken Doyle

Gene-editing techniques have attracted considerable attention in recent years due to their therapeutic potential and application in diverse areas of research, such as disease modeling, drug discovery, and agrigenomics. A fundamental requirement common to all gene-editing methods is validation of the results. The QuickExtract DNA Extraction Solution has emerged as a popular reagent to extract genomic DNA for PCR-based validation, due to its speed, convenience, and ease of automation. It has also been used in PCR-based genotyping with a broad range of sample types.

INTRODUCTION

When isolating nucleic acids for downstream applications, researchers generally choose between two types of methods: purification, which is a multi-step, longer procedure that provides high-quality DNA or RNA, and extraction, which is a rapid procedure designed specifically for subsequent PCR or RT-PCR. The QuickExtract family of products provides a simple, rapid method for the extraction of PCR-ready DNA from a variety of sources. For most sample types, the procedure can be completed in as little as 8 minutes. The QuickExtract products are easy to use, with a single-tube protocol that does not require the use of spin columns, which can cause considerable sample loss.

The QuickExtract workflow avoids the use of toxic chemicals and is easily adapted for multiwell plates with robotic automation systems. Because of the speed of the QuickExtract method, it is ideal for applications where the results from PCR amplification are required to proceed with further experiments, such as screening and genotyping assays.

QuickExtract DNA Extraction Solution is suitable for a broad range of sample types, from hair follicles to cultured cells. Due to the more challenging nature of some sample types, additional QuickExtract kits




Fig. 1. The QuickExtract DNA Extraction Solution isolates PCR-ready DNA in 8 minutes or less.

were developed to handle these difficult samples. The QuickExtract FFPE DNA Kit eliminates the need for tedious extraction procedures that use organic solvents for formalin-fixed, paraffin-embedded (FFPE) samples. The QuickExtract Plant DNA Extraction Solution enables rapid extraction of plant genomic DNA from leaf tissues for both endpoint and real-time PCR. This white paper provides a brief

overview of published applications that represent the diverse uses for the QuickExtract DNA Extraction Solution.

CRISPR-CAS GENE EDITING

The ability to make precise, targeted changes to genes within living organisms has undergone a revolution in recent years, largely thanks to the rapid development of CRISPR-Cas technology. This technology has its origins in a bacterial immune defense mechanism that uses clustered regularly interspaced short palindromic repeats (CRISPR) together with CRISPR-associated (Cas) nucleases to recognize and destroy foreign DNA from invading viruses. Several types of CRISPR-Cas systems have been identified, of which the CRISPR-Cas9 system is the most commonly used.

Typically, foreign DNA is fragmented within bacteria and incorporated into CRISPR sequences. When the CRISPR locus is transcribed, the transcripts are processed into smaller CRISPR RNAs (crRNAs). In conjunction with bacterial trans-activating CRISPR RNAs (tracrRNA), these crRNAs can target and cleave complementary sequences in bacteriophage DNA during future infections.

Genome-editing techniques take advantage of the relatively small number of components in the CRISPR-Cas9 system.

The sequences of the crRNA and tracrRNA are often combined into a single guide RNA (sgRNA). Sequences coding for sgRNA and Cas9 can be incorporated into a single vector that is used to transfect or transduce the host cell. Alternatively, purified Cas9 and transcribed guide RNA can be complexed and delivered to the host cell as a ribonucleoprotein (RNP). Once inside the cell, the guide will target Cas9 activity to a specific locus. Cas9 activity results in a double-stranded break (DSB), which will be repaired by the error-prone nonhomologous end-joining (NHEJ) pathway. This type of repair often produces insertion and/or deletion (indel) mutations at the target site. If the sgRNA and Cas9 sequences are delivered to the cell with a "repair template," then the cell uses an alternative repair mechanism called homology-directed repair (HDR). A repair template contains a DNA sequence to be inserted so that a locus can be modified in a specific manner. The modified target sequences are then validated, typically by PCR-based screening or next-generation sequencing (NGS).

Screening for Mutations

An important part of any gene-editing method is the screening of mutated and wild-type cells to confirm the success and efficiency of the procedure. Several assays have been developed for this purpose. Of these, enzyme mismatch cleavage (EMC) assays have been widely used in CRISPR-Cas9 based methods. The first EMC assays used the bacteriophage resolvase enzymes T4 endonuclease VII¹ and T7 endonuclease I². Subsequently, plant-based S1 nucleases were introduced, and the technology is currently commercialized in the

line of Surveyor mutation-detection assays from Integrated DNA Technologies (IDT).

In a typical EMC assay, genomic DNA is isolated from a population of cells, containing wild-type and mutant alleles. The targeted region is amplified by PCR and then denatured and reannealed. This process results in mismatched duplexes at the targeted sites, due to cross-annealing of the wild-type and mutated sequences. Digestion of the amplicons with T7 endonuclease I, which selectively cleaves distorted duplex DNA, is followed by analysis of the DNA fragments by polyacrylamide gel electrophoresis.

The use of QuickExtract DNA Extraction Solution to prepare genomic DNA for an EMC assay was originally reported to screen for mutations introduced by zinc-finger nucleases (ZFNs)³. The assay is sensitive down to ~1% gene modification in a population, and the authors cite its speed, low cost, and convenience as advantages over other mutation-detection methods for ZFN-based gene editing.

Subsequently, the laboratory of Dr. Feng Zhang at the Massachusetts Institute of Technology developed a standard EMC protocol for mutation detection following CRISPR-Cas9 gene editing, using genomic DNA isolated with QuickExtract DNA Extraction Solution⁴. Extract workflow contribute to this specific, scalable, and cost-efficient assay for genome editing in mammalian cells.

Genome-Scale Knockout

Joung et al.⁵ described a system for performing pooled genome-scale

knockout or transcriptional activation screening using CRISPR-Cas9 in a variety of human and mouse cell lines. They constructed plasmids containing variable sgRNAs that, when pooled, make up an sgRNA library for genome-scale targeting. The system format can either encode the sgRNA, Cas9, and other necessary components on the same plasmid or on separate plasmids. The sgRNA libraries were packaged into lentivirus and transduced into the cell lines being studied. After applying selective pressure to the cells, the resulting population was screened by NGS for changes in sgRNA distribution and, from the results, a number of candidate target genes were selected for validation. Following cloning and lentiviral packaging of validation sgRNAs targeting the candidate genes of interest, these select sgRNAs were transduced into cells. Once the cells reached confluency, genomic DNA was isolated using QuickExtract DNA Extraction Solution in 96-well plates, and the regions containing the sgRNA sequences were amplified by a two-step PCR for initial validation by targeted NGS on an Illumina[®] platform. The researchers report that although their method is designed for loss-of-function and gain-of-function screening, it can be used for other types of screening as well.

Improved Gene-Editing Efficiency

Research using human pluripotent stem cells (hPSCs) has enabled an array of novel technologies in the study of human development, as well as disease progression and therapeutic intervention. The ability of CRISPR-Cas9 technology to introduce specific, targeted mutations in hPSCs has added a new dimension to this rapidly expanding field of research. Xie et al.⁶

describe a novel, episomal vector-based CRISPR-Cas9 system (epiCRISPR) for the generation of knockout mutations in hPSCs. The vector includes components from Epstein-Barr virus that allow it to replicate once per cell division in eukaryotes, allowing persistent expression of Cas9 and sgRNAs. Further, the transfected cells can be enriched by puromycin selection. The epiCRISPR system was tested in separate experiments for gene knockout and genomic deletions, using custom-designed sgRNA panels. After transfection of hPSCs and puromycin selection, genomic DNA was isolated using the QuickExtract DNA Extraction Solution. The target region was amplified by PCR, purified, and analyzed by restriction fragment length polymorphism (RFLP). Off-target effects were also examined using NGS. The authors reported high efficiency of gene knockouts and up to 100% efficiency for indel generation with the epiCRISPR system, with no off-target effects.

Allele Exchange

While the majority of CRISPR-Cas9 methods are used for gene knock-out or knock-in experiments, the exchange of an entire allele is an attractive prospect for potential therapeutic applications. Kelton et al.⁷ provide a proof-of-concept method to reprogram the polymorphic major histocompatibility complex (MHC) locus in murine-derived antigen-presenting cell lines, using a CRISPR-Cas9 system. The authors used Cas9 with sgRNAs that introduced DSBs flanking the MHC-I allele (~3.4kb), enabling its replacement with a MHC donor cassette, which was introduced as a double-stranded DNA template. To confirm insertion of MHC donor DNA, candidate cells were treated with QuickExtract DNA Extraction Solution, and the genomic DNA was amplified using primers covering the putative cleavage sites at the MHC locus for analysis. The authors conclude that their approach demonstrates

the feasibility of replacing large MHC alleles at the native locus, suggesting its future utility for correcting MHC mismatches in allogenic cellular transplantations.

TALENS

Concern about the off-target effects of CRISPR-derived gene-editing methods has spurred renewed interest in TALENs, which typically offer high targeting specificity. Wang et al.⁸ optimized a TALEN design and tested their method in human cultured cell lines, as well as human embryonic stem cells (hESCs). The method combined an optimized TALEN with a single-strand oligodeoxynucleotide (ssODN) for high-efficiency gene editing. The cutting efficiency of the optimized TALEN was validated in HEK293T cells after transfection with TALEN plasmid by treating the cells with QuickExtract DNA Extraction Solution and analyzing the percentage of gene modification using a Surveyor nuclease assay. Subsequent experiments established human tumor cell lines and hESCs with homozygous deletions of two microRNA genes, miR-21 and miR-9-2. The authors suggest that the versatility of the TALEN-ssODN method could contribute to its utility in therapeutic applications.

GENOTYPING

The speed of the QuickExtract DNA workflow has made it popular for PCR-based genotyping in a variety of organisms, from bivalves to mice.

Mussels

Ferguson et al.⁹ characterized population structure, dispersal potential, and reproductive strategies in the freshwater mussel *Lampsilis cardium* from three watershed sites in Ohio, USA. Genomic DNA was extracted from adult mussels using the PureGene kit (Gentra) and from ethanol-preserved glochidia (mussel larvae) using QuickExtract DNA Extraction Solution. The QuickExtract procedure was modified to

extend the time for both heating steps: 30 minutes at 65°C followed by 7 minutes at 98°C. The extracted DNA was used in PCR with primers for 12 microsatellite loci to assess local population structure relative to within-population patterns of relatedness and parentage. The study concluded that long-distance fertilization observed in this mussel species may have implications for improved recovery of imperiled freshwater mussel populations.

Mice

The G protein-coupled trace amine-associated receptor 1 (TAAR1) is stimulated by neurotransmitters and metabolites, and the receptor has multiple physiological roles, including predisposition to drug abuse and related effects. Shi et al.¹⁰ studied the functional effects of single-nucleotide polymorphisms (SNPs) in the TAAR1 gene in mice. Genomic DNA from several mouse strains was extracted from ear or tail tissue using QuickExtract DNA Extraction Solution and amplified by PCR with primers flanking the SNP-containing regions. The amplicons were purified and sequenced. The results showed that TAAR1 mutants resulted in expression of functional, subfunctional, and nonfunctional receptors, with widely differing responses to methamphetamine. The authors suggest that SNPs in TAAR1 could provide a useful screening tool for determining predisposition to a variety of human diseases, as well as for individualizing treatments using TAAR1-specific therapies.

Transgenic Mice

Tumor-specific antigens (TSAs) constitute attractive therapeutic targets in a wide range of cancers. Unlike TSAs, tumor-associated antigens (TAAs) are expressed at varying levels in normal cells as well as tumors, making therapeutic strategies more challenging. Yong et al.¹¹ studied the TAA human epidermal growth



factor 2 (Her2), which is expressed in normal breast tissue but upregulated in 15%-30% of breast cancers. They used a transgenic mouse model to examine therapies that specifically targeted Her2 antigen in tumor cells. While heterozygous Her2^{+/-} mice appear to develop in a similar manner to wild-type mice (Her2^{-/-}), it has proven difficult to generate homozygous Her2^{+/+} mice, potentially due to embryonic lethality. Genomic DNA was extracted from tail snips of Her2^{+/+}, Her2^{+/-} and Her2^{-/-} mice using QuickExtract DNA Extraction Solution, followed by PCR to examine the integration site of the Her2 transgene.

the Her2 transgene was indeed responsible for embryonic lethality, due to its insertion into the gene encoding precocious dissociation of sisters (Pds5b), which is implicated in embryonic development.

CONCLUSION

The speed and convenience of the QuickExtract methods make them ideal for PCR validation during intermediate steps in gene-editing and genotyping workflows. The QuickExtract procedure allows for a large number of samples to be screened in

a rapid, automation-friendly process. The utility of the QuickExtract DNA Extraction Solution in applications requiring PCR-based screening has been recognized by its incorporation into standard CRISPR-Cas9 gene-editing protocols provided by IDT (Alt-R™ Genome Editing Detection Kit, Protocol CRS-10056-PR) and New England Biolabs (EnGen® Mutation Detection Kit, Technical Manual E3321).

REFERENCES

For a complete list of references, visit vwr.com/vwrbioMarke.

Subsequent whole-genome sequencing confirmed that the integration site of

Description	Cat. No.	Unit
QuickExtract™ DNA Extraction Solution	76081-768	5 mL
QuickExtract™ DNA Extraction Solution	76081-766	50 mL

Agilent QuikChange Lightning Kits

Discover the Fastest and Latest Generation of the Market-Leading Site-Directed Mutagenesis Kits

- Fast, reliable, and easy QuikChange mutagenesis protocol
- 75% reduction in time for both single site and multiple site mutations
- >80% mutation efficiency for both short and long templates (up to 14 kb)

QuikChange kits provide researchers with a fast, easy, and efficient non-PCR method to reliably perform site-directed mutagenesis. The QuikChange method uses a linear amplification strategy with only the parental strand serving as the DNA template. Combining this with our highest fidelity polymerases leads to a significant reduction in unwanted second-site errors.



Description	Size	Cat. No.	Unit
QuikChange HT Protein Engineering System 150nt	10 Sites	76193-646	Each
QuikChange HT Protein Engineering System 150nt	20 Sites	76193-648	Each
QuikChange HT Protein Engineering System 150nt-Academic	10 Sites	76193-654	Each
QuikChange HT Protein Engineering System 150nt-Academic	20 Sites	76193-656	Each
QuikChange HT Protein Engineering System 200nt	10 Sites	76193-650	Each
QuikChange HT Protein Engineering System 200nt	20 Sites	76193-652	Each
QuikChange HT Protein Engineering System 200nt-Academic	10 Sites	76193-658	Each
QuikChange HT Protein Engineering System 200nt-Academic	20 Sites	76193-660	Each

For Research Use Only. Not for use in diagnostic procedures.

These products will be available soon in Canada. Please contact your VWR Sales Rep to learn about similar options available in your region in the meantime.



High-Performance EvaGreen® qPCR Master Mixes



Forget-Me-Not™ EvaGreen qPCR Master Mixes Offer Superior Performance and Sensitivity

- EvaGreen dye for superior real-time qPCR and melt curve analysis — higher signal, better sensitivity, no PCR inhibition, and environmentally friendly
- Fast-activating, chemically-modified, hot-start Cheetah™ Taq for excellent specificity and speed
- Color-tracking to streamline sample set-up and minimize pipetting errors
- Unrivalled performance at unbeatable prices! As low as 20¢/rxn

Forget-Me-Not EvaGreen qPCR Master Mix is ideally suited for qPCR and HRM applications. The unique features of EvaGreen offer several advantages over SYBR® Green dye-based mixes. The master mixes support regular and fast-cycling PCR protocols, and reagent color-tracking helps minimize errors, saving precious samples, reagents and time. PCR products can be directly analyzed by gel electrophoresis, without staining.



Description	Size	Cat. No.	Unit
Forget-Me-Not EvaGreen qPCR Master Mix			
Low ROX	5 mL	76221-616	Each
Low ROX	20 mL	76204-718	Each
High ROX	1 mL	76221-722	Each
High ROX	20 mL	76204-720	Each
Forget-Me-Not EvaGreen qPCR Master Mix with 2-Color Tracking			
No ROX	1 mL	10001-796	Each
No ROX	20 mL	76204-724	Each
ROX Provided Separately	1 mL	10001-800	Each
ROX Provided Separately	20 mL	76204-722	Each



A Faster, Smaller, Better Way to qPCR

- **Ultra-Fast Data Acquisition** — 35 cycles in 25 minutes
- **Unrivalled Performance** — Detect 2-fold expression level differences
- **Portable and Compact** — 4.5 lbs; easy transportation without recalibration
- **Scalable and Wireless** — Connect up to 10 instruments
- **Magnetic Induction Technology** — Eliminate variability vs block-based cyclers

The Q PCR cycler uses a patented magnetic induction technology to rapidly heat and cool samples in order to acquire data in as little as 25 minutes. A unique spinning rotor provides superior temperature uniformity eliminating edge effects and variability versus traditional block-based cyclers, allowing you to generate data that is superior in reproducibility and repeatability, while enabling the detection of 2-fold gene expression level differences and difficult class IV SNP's.




Description	Size	Cat. No.	Unit
Q 2-Channel qPCR Instrument	—	76175-392	Each
Q 4-Channel qPCR Instrument	—	76175-394	Each
Accessories			
Q Tubes and Caps	20 Racks of 48 Tubes	76202-252	Each

With a Studio in Your Lab, Your Research is Unlimited



UVP ChemStudio™ touch

- Advanced application possibilities with flexibility in choice of camera and transilluminators, and optional overhead NIR module
- Prioritized data integrity with unmodified raw image data
- Designed for convenience and streamlined workflow, featuring a 13.3" integrated touch screen monitor, fold down door, and Thin-Line Transilluminator on a slide tray

analytikjena
An Endress+Hauser Company



Description	Camera & Lens	Software Control	Cat No.	Unit
UVP ChemStudio Touch	8.1MP, f/0.95	Integrated Touch Ccreen	75873-962	Each
UVP ChemStudio	8.1MP, f/0.95	PC-Controlled	75873-966	Each



Demonstration of FlashGel™ Recovery with Sheared DNA

By Christopher Veilleux and Hugh White, Lonza Rockland, Inc.

Introduction

More and more molecular biologists are working with sheared DNA. Sheared DNA is critical for constructing genomic libraries, and is utilized in DNA reassociation and hybridization analysis. Published methods for DNA fragmentation can be classified into four categories: sonication; enzymatic digestion; hydrodynamic shearing; and nebulization. Regardless of the fractionation method, gel electrophoresis is one of the best ways to estimate mean size and distribution of the DNA material.¹

The FlashGel™ system for recovery is a fast and simple tool for separating and recovering DNA, including fragmented DNA.² With the recent addition of a double-tier 2.2% concentration cassette format, it is possible to separate and recover a narrow window within a sheared DNA sample. DNA may be recovered in as little as 5 minutes, without the need for UV light or downstream purification.

Capabilities of the FlashGel™ System for Recovery of Sheared DNA

The sensitivity of the FlashGel™ system allows for a wide range of starting material concentrations. As with most gel recovery systems, the higher the concentration of starting sample, the higher the recovery yield. Intensity of recovered DNA from FlashGel™ cassettes scales well with load levels of starting DNA, (Figures 1A and 1B). The FlashGel™




Figure 1A: Sheared DNA on the FlashGel™ system
A dilution series of sheared DNA 600–37.5ng in amount and approximately 200–4,000bp in size were run on a double-tier 1.2% FlashGel™ recovery cassette at 275V for 4 minutes. Lanes QL contain FlashGel™ quantladder.




Figure 1B: Recovered samples from sheared DNA
4µL aliquots of Figure 1A samples post recovery were run on a single-tier 1.2% FlashGel™ cassette at 275V for 5 minutes. Lanes QL contain the FlashGel™ quantladder. Lane C contains 50ng of the original sheared DNA sample as a reference. Numeric lane labels indicate amount of starting material correlated from Figure 1A.




Figure 2: Recovered samples of various size selections
Recovered samples of various size selections from fragmented genomic lambda DNA. Lanes M contain the FlashGel™ 100–4,000bp marker, lanes QL contain the FlashGel™ quantladder, lane C contains 50ng of sheared DNA, lanes 1–5 contain DNA extractions from several experiments taken from 1.2% FlashGel™ recovery cassettes. Samples were run on a single-tier 1.2% FlashGel™ cassette at 275V for 7 minutes. Lane 1: 280–550bp; lane 2: 650–1,190bp; lane 3: 425–800bp; lane 4: 425–800bp (smaller load volume of Lane 3); lane 5: 800–3,000bp M QL C 1 2 3 4.




Figure 3: Sheared DNA samples following size selection recovery

Lower MW fragments from two sheared DNA samples were recovered from a double-tier 2.2% FlashGel™ recovery cassette. Several seconds of 275V current were applied after the selected areas were removed from the samples. The images show the area where the recovered samples were removed. Samples of higher molecular weight are ready to be recovered. Lanes QL contain FlashGel™ quantladder.

recovery system is capable of recovering various size ranges and various size range windows within the sheared DNA (Figure 2).

Recovery steps may be performed in parallel lanes to recover samples from multiple lanes, and multiple samples may be recovered from a single lane. Simply run the smallest of the desired size range in to the recovery well, stop the voltage, add FlashGel™ recovery buffer, and collect the sample. Then continue running the unwanted DNA through the recovery well, until the next desired size range reaches the well. Then stop the voltage, add more recovery buffer, and collect the next sample. As long as current is applied, DNA that is not recovered will migrate through the well and leave no trace contamination (Figure 3). Only when the current is stopped, and the

FlashGel™ recovery buffer is added, will DNA remain in the recovery well. Recovering DNA using the FlashGel™ system is highly efficient. Figure 3 shows clear voids of samples removed during recovery.

Capabilities of 1.2% and 2.2% FlashGel™ Cassettes

Extraction range windows are predefined by selecting and using the 1.2% or the 2.2% recovery cassettes. Optimal separation range for the 1.2% FlashGel™ recovery cassettes is 50–4,000bp, and optimal separation range for the 2.2% recovery cassettes is 10–1,000bp. The 2.2% cassettes provide a larger spread of bands in the 10–400bp range, while the 1.2% cassettes show a wider spread in the larger size range




Figure 4: Demonstration of size separation Sheared DNA consisting of fragments 50–3,000bp run on 1.2% and 2.2% FlashGel™ cassettes. Lane QL is FlashGel™ quantladder, lane 1 is the DNA.

(Figure 4). The 2.2% cassettes will provide a tighter and more defined range than the 1.2% cassettes within the same size range (Figure 5). The 1.2% cassettes may be better suited for a wider size selection range in a single extraction cycle.

Summary

The FlashGel™ system is an efficient tool for fragmented DNA size selection. The two recovery cassette concentrations (1.2% and 2.2%) are optimized to provide narrow separation over a wide size range. With the FlashGel™ system, DNA fragments may be separated, recovered and photographed in as little as 5 minutes, without the use of UV light. Recovered samples are immediately ready for downstream applications — eliminating agarose gel preparation, band excision and purification.

REFERENCES

For a complete list of references, visit vwr.com/vwrbiomarke.




Figure 5: Demonstration of size recovery Lanes labeled 1 and 3 are extractions of samples through recovery well window from 1.2% and 2.2% cassettes respectively. Lanes 2 and 4 are slightly higher load volumes of lanes 1 and 3 respectively. Lane QL contains the FlashGel™ quantladder.



Lonza FlashGel™ System for Rapid Genomic Analysis

Lonza

Say Goodbye to Gel Preparation, Band Excision, Purification, and UV Light

- Separations complete in as little as 2–3 minutes
- Monitor in real time and capture gel images at the bench with no UV light
- Recover samples directly from the FlashGel cassette

Complete separation, recovery, and documentation safely, at the bench, in minutes with the FlashGel system (running dock and cassettes required, camera optional). FlashGel is the ideal preparatory tool for isolation and/or screening of DNA or RNA samples for genome editing, sequencing, and other applications. FlashGel cassettes contain ready-to-use precast gels with buffer and stain; simply prepare your samples, load, and run.



Description	Cat. No.	Unit
FlashGel DNA System	95045-604	Each
FlashGel Device Kit	89400-686	Each
DNA Cassette, 1.2%, 12+1 Well, Single Tier	95015-618	Pack of 9
DNA Cassette, 2.2%, 16+1 Well, Double Tier	95015-624	Pack of 9
RNA Cassette, 1.2%, 12 + 1 Well, Single Tier	95015-626	Pack of 9
Recovery Cassette, 1.2% Agarose, 8+1 (18-Well), Double Tier	95053-310	Pack of 9
DNA Marker (100 bp–4 kb), 500 µL	95015-632	Each
Loading Dye, 5 x 1 mL, 5X Concentrate	95015-630	Each

Accelerate Your Cloning Experiments with Reagents from New England Biolabs



Race to NEB for the Highest Quality Products That Save Valuable Time!

- **Quick Dephosphorylation Kit** — Complete dephosphorylation in 10 minutes with rapid, irreversible heat inactivation
- **Instant Sticky-end Ligase Master Mix** — Specially formulated to instantly ligate 2–4 base pair overhangs

While cloning has typically been a time-consuming process, advances in both technologies and reagents have enabled the process to now take much less time. For a limited time, try our Fast Sample Pack, containing samples of our Instant Sticky-end Ligase Master Mix and Quick CIP from our Quick Dephosphorylation Kit. Contact your VWR representative to request a sample.



Description	Size	Cat. No.	Unit
Quick Dephosphorylation Kit	100 Reactions	103258-556	Each
Quick Dephosphorylation Kit	500 Reactions	103258-516	Each
Instant Sticky-end Ligase Master Mix	50 Reactions	102715-946	Each
Instant Sticky-end Ligase Master Mix	250 Reactions	102715-944	Each

Agilent CRISPRa/i Libraries



Agilent Technologies

Designed with Content from UCSF

- Synthesized with industry's lowest error rates
- Full control of targets and guide sequence
- Available as Ready-to-Amplify and Ready-to-Clone libraries in catalog or custom options

Agilent synthesizes the highest fidelity, longest oligos in the industry, and we are now applying this technology to make CRISPR gRNA libraries. CRISPRa/i libraries with gene targets from UCSF are available in genome-wide and prevalidated subsets relevant to disease-related research.



Description	Cat. No.	Unit
Genome Wide Libraries		
SureGuide Human CRISPRa - Genome Wide	76204-942	Each
SureGuide Human CRISPRi - Genome Wide	76204-940	Each
SureGuide Mouse CRISPRa - Genome Wide	76204-946	Each
SureGuide Mouse CRISPRi - Genome Wide	76204-944	Each
Genome Wide Libraries		
SureGuide Human CRISPRa - Cancer/Apoptosis	76204-950	Each
SureGuide Human CRISPRa - Drug/Kinase/Phosphatase	76204-954	Each
SureGuide Human CRISPRa - Stress/Proteostasis	76204-958	Each
SureGuide Human CRISPRi - Gene Expression	76204-966	Each
SureGuide Human CRISPRi - Membrane Proteins	76204-970	Each
SureGuide Human CRISPRi - Mito/Traffic/Motility	76204-960	Each

The above six sublibraries are all available for Human and Mouse, Activation and Interference. For Research Use Only. Not for use in diagnostic procedures.

These products will be available soon in Canada. Please contact your VWR Sales Rep to learn about similar options available in your region in the meantime.

Take the 20 Minute Route to Mycoplasma Detection

Lonza

Get Reliable Results within 20 Minutes Using the Lonza MycoAlert™ PLUS Assay

- Bioluminescence-based technology with no DNA extraction necessary
- Simply add two reagents to the culture supernatant and perform two luminescence readings
- Detects all common mollicute contaminations, except Ureaplasma (Mycoplasma, Acholeplasma, Entomoplasma, and Spiroplasma)

The MycoAlert PLUS Assay is a selective biochemical test that exploits the activity of mycoplasmal enzymes, which are not present in eukaryotic cells. If viable mycoplasmas are present, the enzymes react with the MycoAlert PLUS Substrate and generate ATP. The ATP is then transferred into a light signal via the luciferase enzyme in the MycoAlert PLUS Reagent to show that mycoplasma is present. The MycoAlert Assay Control Set (does NOT contain mycoplasma) is sold separately.

Visit vwr.com/lonzamycoalert for more information!



Description	Size	Cat. No.	Unit
MycoAlert™ Plus Mycoplasma Detection Kit	10 Tests	75860-358	Each
MycoAlert™ Plus Mycoplasma Detection Kit	30 Tests	75860-360	Each
MycoAlert™ Plus Mycoplasma Detection Kit	50 Tests	75860-362	Each
MycoAlert™ Plus Mycoplasma Detection Kit	100 Tests	75866-212	Each
MycoAlert™ Assay Control Set	10 Tests	75870-450	Each

Capture Your Genetic Snapshot

By Shaun Veran, Sr. Research Associate, Zymo Research

The Importance of Sample Collection & Preservation

Sample collection is the first step in the purification of all workflows that involve nucleic acids. In this critical step, it is imperative that steps are taken to ensure proper storage and preservation of the samples as these methodologies can dramatically affect downstream analyses. We have all heard the saying “garbage in, garbage out.” Compositional changes and bias can occur due to nucleic acid degradation, cellular growth/decay, and the logistics of collection.

Current preservation methods use cold-chain logistics to prevent or slow these changes to a sample. However, these systems are very costly and without proper storage conditions, can result in a misrepresentation of an analyte’s abundance, bias results, reduce assay sensitivity, and decrease reproducibility.

These steps are most critical for RNA which is especially vulnerable to degradation from the prevalence of RNases and an unstable phosphodiester bond of the backbone. However, DNA can also be prone to rapid degradation and complete signal loss. Because of this, it is necessary to reliably preserve these nucleic acids.

Stabilization of Nucleic Acids

Things change, make sure your results do not! Standardized sample collection in the research and clinical setting is becoming more important. But rather than using costly cold-chain logistics, Zymo Research has reimagined the sample collection process with a single preservation reagent, DNA/RNA Shield™.

DNA/RNA Shield™ stabilizes both DNA and RNA without the need for cold storage (-80°C), protecting samples even at ambient




Figure 1. Whole blood samples +/- DNA/RNA Shield™ were subjected to > 2 freeze thaw cycles. Total RNA was subsequently purified using the Quick-RNA™ Whole Blood kit.




Figure 2. RNA from tissue stored in DNA/RNA Shield™ (included with the Quick-RNA™ Miniprep Plus Kit) is preserved at ambient temperature. RNA from muscle tissue (mouse) was purified using the Quick-RNA™ Miniprep Plus Kit and analyzed by RT-PCR.

temperatures for up to 30 days. This eliminates the need for costly storage conditions and specialized sample collection procedures while maintaining a preserved snapshot of the DNA/RNA profile, even through multiple freeze-thaw cycles.

RNA samples are very susceptible to degradation. Time, temperature, and freeze-thaw cycling all contribute to issues with sample stabilization. DNA/RNA Shield™ provides worry-free protection against these changes. As seen in Figure 1, the DNA/RNA Shield™ effectively prevents degradation after whole blood samples underwent multiple freeze-thaw cycles before purification.

DNA/RNA Shield™ protects samples kept at ambient temperatures over time. Analysis of RNA purified from tissue stored in DNA/RNA Shield™ shows there is no change in the average C_t values of muscle tissue stored in DNA/RNA Shield™ for up to 30 days at ambient temperatures (Figure 2).

Maintain Nucleic Acid Profiles

Shifting nucleic acid profiles can also be observed for unprotected samples. This is most exemplified in the microbial composition of a stool sample (Figure 3). Without proper stabilization, high levels of variability and bias are introduced. The compositional proportion of DNA for specific species changes drastically over time in unprotected samples. Conversely, DNA/RNA Shield™-protected samples change little compared to the original collected sample.



SAMPLE COLLECTION PRODUCT GUIDE

1. Choose a Collection Device:

	<p>Swabs</p> <p>8mL tube with breakable sterile swab</p> <ul style="list-style-type: none"> • Mouth • Nose • Throat • Fluid 		<p>Stool</p> <p>15mL tube with screwcap scoop</p> <ul style="list-style-type: none"> • Virus • Microbe • Host
	<p>Blood</p> <p>10mL blood tube</p> <ul style="list-style-type: none"> • Fresh • EDTA • Citrate • Heparin 		<p>Swabs</p> <p>2mL tube with high density beads</p> <ul style="list-style-type: none"> • Animal • Plant • Insect • Microbe

2. Choose a Purification Kit:

SAMPLE TYPE	DNA	RNA	BOTH
Blood Collection Tube			
Entire Blood Tube	Quick-DNA/RNA™ Blood Tube Kit		
200µL Whole Blood Input	Quick-DNA™ Plus Kits	Quick-RNA™ Whole Blood Kit	Quick-DNA/RNA™ Kit
Fecal or Soil Samples	Quick-DNA™ Plus Kits	Quick-RNA™ Plus Kits	—
Cell or Soft Tissue	Quick-DNA™ Plus Kits	Quick-RNA™ Plus Kits	Quick-DNA/RNA™ Kit
Microbiomic Samples	ZymoBIOMICS® DNA Kits	ZymoBIOMICS® RNA Kits	ZymoBIOMICS® DNA/RNA Miniprep Kits
Viral Samples	Quick-DNA™ Viral Kits	Quick-RNA™ Plus Kits	Quick-DNA/RNA™ Viral Kit



Description	Size	Cat. No.	Unit
DNA/RNA Shield Blood Collection Tube	50 Tubes	76020-434	Each
DNA/RNA Shield Fecal Collection Tube	10 Tubes	76020-422	Each
DNA/RNA Shield Lysis Tube (Microbe)	50 Tubes	76020-426	Each
DNA/RNA Shield Lysis Tube (Tissue)	50 Tubes	76020-098	Each
DNA/RNA Shield Collection Tube with Swab (1 mL fill)	10 Tubes	76020-100	Each
DNA/RNA Shield Collection Tube with Swab (2 mL fill)	10 Tubes	76020-430	Each
DNA/RNA Shield Reagent	50 mL	76020-420	Each
DNA/RNA Shield Reagent (2X Concentrate)	25 mL	76020-444	Each

These specific products are not available in Canada. Please contact your VWR Sales Rep to learn about similar options available in your region.




Figure 3. Microbial composition (indicated by different colors for different species) from a stool sample is unchanged after one month at ambient temperature with DNA/RNA Shield™. Stool samples suspended in DNA/RNA Shield™ and stored at room temperature were compared to stool without preservative for one month. Samples were collected at the indicated time points and DNA was subsequently extracted with the ZymoBIOMICS® DNA Miniprep Kit. The DNA was then subjected to microbial composition profiling via 16S rRNA gene targeted sequencing. Samples stored with DNA/RNA Shield™ had a constant microbial composition while the samples stored without shifted dramatically.

It provides an easy and effective method to stabilize both DNA and RNA over time at various temperatures, despite freeze-thaw cycles.

Simplified Purification

DNA/RNA Shield™ drastically streamlines the process of sample collection and nucleic acid purification. It is compatible with all Zymo Research purification kits, and unlike other preservatives such as RNAlater®, it does not need to be removed prior to sample processing. Simply add binding reagent to a sample protected in DNA/RNA Shield™ and purify.



To help determine which kit would work best for your needs, please refer to the Sample Collection Product Guide listed at the end of this article. DNA/RNA Shield™ is offered in a variety of collection vessels to best suit the needs of your workflow: Blood Collection Tube, Fecal Collection Tube, Lysis Tube, or Swab & Collection Tube. Once a vessel is chosen, the appropriate extraction kit can be determined from the starting sample type.



Monarch® Kits for DNA and RNA Purification



Sustainably-Designed Purification Kits for Total RNA Extraction, Isolation of Plasmids, and Cleanup of DNA Fragments

- **Excellent performance** — Consistently high nucleic acid yields and purity
- **Ease of use** — Enjoy fast, simple protocols
- **Flexibility and convenience** — In addition to convenient kit formats, buffers and columns are available separately

Monarch kits enable quick and easy purification of high-quality DNA and RNA, suitable for direct use in a variety of downstream applications.

Recover pure, intact DNA and RNA in minutes with fast, user-friendly protocols and optimized buffer systems, and focus your time on the experiments that will drive your research forward. Request a free sample from your VWR representative today.



Description	Size	Cat. No.	Unit
Monarch Total RNA Miniprep Kit	50 Preps	103529-148	Each
Monarch Plasmid Miniprep Kit	50 Preps	102971-698	Each
Monarch Plasmid Miniprep Kit	250 Preps	102971-696	Each
Monarch DNA Gel Extraction Kit	50 Preps	102971-670	Each
Monarch DNA Gel Extraction Kit	250 Preps	102971-668	Each
Monarch PCR and DNA Cleanup Kit	50 Preps	102971-674	Each
Monarch PCR and DNA Cleanup Kit	250 Preps	102971-672	Each
Monarch DNA Cleanup Columns (5 µg)	100 Columns	102971-676	Each

CELLvo™ Cells



The Cell of Choice for Discerning Scientists. Never Exposed to Artificial Substrates.

- More proliferative
- More natural response to experimental variables
- Less phenotypic drift

CELLvo Cells are primary cells isolated and grown on CELLvo Matrix, a cell-derived substrate that mimics the natural extracellular microenvironment containing over 150 proteins. CELLvo Cells are never exposed to artificial substrates and when grown on CELLvo Matrix express a more natural gene profile compared to cells grown on standard tissue culture plastic.



Description	Size	Cat. No.	Unit
CELLvo Human Bone-Marrow Mesenchymal Stem Cells	>1,000,000	76203-012	Each
CELLvo Human Cord Blood Endothelial Progenitor Cells from a Single Donor Cryopreserved at P4	>500,000	76169-452	Each
CELLvo Human Mesenchymal Stem Cell – Wharton’s Jelly Cryopreserved at P2	>500,000	76200-216	Each
CELLvo Primary H-Chondrocytes (HC-a) P0	>1,000,000	76203-918	Each
CELLvo H-Chondrocytes (HC-a) P1	>1,000,000	76203-926	Each

Mirus Bio TransIT® Lentivirus System

An Optimized Transfection Reagent and Lentivirus Packaging Mix for High, Functional Titers with a Streamlined Workflow

- Produce 2–3X higher functional titers than Lipofectamine® 2000 and 3000
- Achieve higher titers without the need to concentrate
- No media change required, single harvest

Combined with an efficient host-genome integration mechanism capable of infecting both dividing and non-dividing cells, recombinant lentivirus has become a central gene delivery tool for robust and stable transgene expression in target cells. The lentivirus platform is often limited by insufficient viral titers requiring concentration before use. To address this limitation, Mirus Bio developed the TransIT® Lentivirus System for enhanced delivery and expression of the essential lentivirus genes required for higher-titer lentivirus production in adherent or suspension HEK 293 cell types. The TransIT® Lentivirus System is composed of the TransIT®-Lenti Transfection Reagent, Lentivirus Packaging Mix Powered by MISSION® Genomics, and TransduceIT™ Reagent.



Description	Cat. No.	Unit
TransIT® Lentivirus System RXN5	76210-950	Each
TransIT® Lentivirus System RXN34	76210-952	Each
TransIT®-Lenti Transfection Reagent 0.3mL	75814-982	Each
TransIT®-Lenti Transfection Reagent 1.5mL	75814-980	Each
Lentivirus Mix Powered by MISSION® RXN5	76218-672	Each
Lentivirus Mix Powered by MISSION® RXN34	76211-034	Each
TransIT®-Lenti Transfection Reagent 5 x 1.5 mL (7.5mL)	75814-986	Each
TransIT®-Lenti Transfection Reagent 10 x 1.5 mL (15 mL)	75814-994	Each

These specific products are not available in Canada. Please contact your VWR Sales Rep to learn about similar options available in your region.

CELLvo™ Matrix – The Next Evolution in Cell Culture™



Use in Mammalian, Adherent Cell Culture to Avoid Unnatural Cell-Matrix Interactions

- Cells grown on CELLvo Matrix express a more natural gene profile and reduced donor-donor variability
- Cells grown on CELLvo Matrix exhibit more natural responses to experimental variables
- CELLvo Matrix enhances isolation of primary cells

CELLvo Matrix is an extracellular matrix of proteins synthesized *in vitro* by bone marrow stromal cells and contains more than 150 proteins that were secreted and assembled by the cells during production. The final product is cell free with only the ECM attached to the surface of the dish providing a native three-dimensional microenvironment for rapid expansion of high quality adherent cells...There is no place like home!



Description	Growth Area	Cat. No.	Unit
CELLvo Matrix			
6-Well Flat Bottom Plates	9.5 cm ²	75853-938	Pack of 5
T75 Vented Flasks	75 cm ²	75853-856	Pack of 5
T150 Vented Flasks	150 cm ²	75853-858	Pack of 5
CELLvo XF Matrix			
XF-Matrix 6-Well Flat Bottom Plates	9.5 cm ²	76183-420	Pack of 5
XF-Matrix T75 Vented Flasks	75 cm ²	76183-422	Pack of 5
XF-Matrix T150 Vented Flasks	150 cm ²	76183-424	Pack of 5

VWR Life Science Seradigm FB Essence

Proven Performance and Consistency

- Nutritionally rich, cost-effective alternative to Fetal Bovine Serum (FBS)
- Proven effective across a broad range of cell types and origins
- Long term stability of supply and pricing
- 100% US Origin

FB Essence is a nutritionally rich, cost-effective alternative to Fetal Bovine Serum (FBS) containing FBS, Bovine Calf Serum, Equine Serum, and a proprietary blend of supplements, vitamins, minerals, and growth factors.

FB Essence offers proven performance and consistency, low endotoxin and hemoglobin levels, and can produce equivalent results to FBS in many applications. FB Essence is 100% US origin.



Description	Size	Cat. No.	Unit
FB Essence	50 mL	10805-184	Each
FB Essence	500 mL	10803-034	Each
FB Essence, Heat Inactivated	50 mL	10799-384	Each
FB Essence, Heat Inactivated	500 mL	10799-390	Each

Gamma irradiation available upon request.

Reagents that Advance Scientific Discovery



PeproTech Cytokines and Growth Factors

- RUO and animal-free grade offered in human, murine, and rat
- Flexible sizes
- Lot-to-lot consistency

PeproTech offers a comprehensive line of full length and fully biologically active recombinant proteins that are developed in-house by PeproTech's experienced protein scientists.



Description	Size	Cat. No.	Unit
Human FGF-basic (146aa)	10 µg	10771-944	Each
Human FGF-basic (154aa)	10 µg	10778-898	Each
Human GM-CSF	20 µg	10780-398	Each
Human IL-2	50 µg	10779-566	Each
Human IL-4	20 µg	10779-606	Each
Human M-CSF	10 µg	10773-706	Each
Human SCF	10 µg	10780-450	Each
Human TPO	10 µg	10773-600	Each



Distributor
GE Healthcare

How Can a Chromatography Bead Free Up Your Time?

By Emma Lind, Product Manager for Ion Exchange Chromatography at GE Healthcare

Have you ever considered how a resin bead in your chromatography column can affect your ion exchange (IEX chromatography) run time?

Today's protein researcher is pressed to deliver results faster than ever before. The less time spent on purification, the more time you will have to focus on studying your protein.

Ion exchange chromatography can be used in the first steps of protein purification for both native proteins and untagged recombinant proteins. While protein purification has become almost routine, new technology and materials always provide an opportunity to look for improvements. That is the case when choosing a Capto™ ImpRes resin for use with viscous samples.

Example with 50 mL *E. coli* lysate run in cold room



Speeding Up Run Time for Viscous Samples

Viscous samples, such as *E. coli* extract, can increase back pressure on the column and force you to reduce the flow rate, resulting in longer load times. Utilizing a robust, modern, and rigid resin — such as Capto Q ImpRes — reduces sample run time significantly by allowing high flow rate with low column backpressure, as seen in the example above.

HiTrap Capto ImpRes contains Capto ImpRes high-flow agarose resin, which has much better pressure/flow properties than Sepharose™ High Performance resin. (Capto ImpRes resin has maximum flow rate of 300cm/h vs 150cm/h for Sepharose High Performance.) This allows more flexibility — even in the HiTrap format — for viscous samples such as *E. coli* lysate, and allows them to be loaded at higher flow rates in a HiTrap Capto ImpRes than a HiTrap HP column. In this example, the IEX total run time is 30% shorter with HiTrap Capto Q ImpRes, than a HiTrap Q HP, with an otherwise similar performing column.

HiTrap Capto ImpRes: a Versatile Tool

HiTrap Capto ImpRes columns can be used for viscous samples as well as for other high-resolution IEX chromatography procedures with a broad range of samples. Since the resolution is similar to Sepharose High Performance, Capto ImpRes can be used in these applications as well, and there is no need to have both columns on hand.

Capto ImpRes resin has been used in diverse applications, including:

- Analyzing the influence of a mitochondrial protein in mammalian cytosolic Fe-S repair
- Researching the mode of action of antimalarial drug primaquine

HiTrap™ Capto™ ImpRes can serve as a universal column for routine work in a research lab. Use HiTrap Capto Q ImpRes for anion exchange chromatography, and HiTrap Capto SP ImpRes for cation exchange chromatography.

REFERENCES

For a complete list of references, visit vwr.com/vwrbioMarke.



Product Description	Pack Size	Cat. No.	Unit
HiTrap™ Capto™ SP ImpRes Prepack Columns	5 x 1 mL	97067-776	Each
HiTrap Capto SP ImpRes Prepack Columns	5 x 5 mL	97067-778	Each
HiTrap Capto SP ImpRes Resin	25 mL	97067-934	Each
HiTrap Capto SP ImpRes Resin	100 mL	97067-932	Each
HiTrap Capto Q ImpRes Prepack Columns	5 x 1 mL	97067-780	Each
HiTrap Capto Q ImpRes Prepack Columns	5 x 5 mL	97067-782	Each
HiTrap Capto Q ImpRes Resin	25 mL	97067-938	Each
HiTrap Capto Q ImpRes Resin	100 mL	97067-936	Each



Image Now! Who Has Time to Wait?

The c300 Imager is Ideal for an Individual Lab's Western Blotting Needs

- Chemiluminescent performance that's affordable
- Push button imaging and easy quantitation
- Lab workhorse

The c300 is a powerful yet budget-conscious system for fast-moving scientists who would like to image their chemiluminescent blot and protein/DNA gels now, and not when a colleague down the hall is done using the shared imaging system. Additionally, it's fully upgradable to a fluorescent imaging system if your experimental needs change.

**Description**

c300: The Darkroom Replacer Imaging System

Cat. No.

10147-220

Unit

Each

Simplifying Immune Oncology Research

CORNING

Novel 3D Tumor Model Streamlines High Throughput Immunotherapy Testing

Content provided by Corning

Immunotherapy has become a topic of increased interest for scientists studying cancer treatment. By using a patient's own immune cells, immunotherapy has shown promise in attacking cancer cells with best-case clinical results of total cancer remission. But immunotherapy is not a one-size-fits-all approach. Treatment strategy and efficacy can vary across patients and cancer types.

In order to more carefully study immunotherapy as a means to treat cancer, scientists need to be able to utilize efficient, *in vitro* screening models to examine targeted therapy options. This type of research has historically been conducted using cell culture models that allow scientists to study immune cell migration, and that can be easily adapted for high throughput screening. However, these models typically utilize two-dimensional (2D) cell culture monolayers which do not accurately reflect a tumor's complexity in the human body, leading to a growing trend to use more *in vivo*-like three dimensional (3D) models. While current 3D models are more *in vivo*-like, in order to incorporate the immune cell component, these 3D models require directly adding immune cells to the 3D spheroids. This may not be an accurate reflection of how immune cells migrate towards, and further invade a tumor within a patient suffering from cancer.

According to Hannah Gitschier, Corning Life Sciences Applications Lab Manager in Kennebunk, ME, "With immunotherapy, it is imperative that researchers have solid models to best recapitulate the *in vivo* microenvironment and interactions between the engineered immune cells and patient tumor samples. These models must be easy to use and transferrable to a high throughput screening environment, enabling rapid screening of conditions to best predict clinical outcomes. Traditionally, 2D assays were the method of choice given the high proliferation rate of cells in 2D culture, as well as the reproducibility of assays. However, with emerging technologies and the development of novel lab consumables, 3D assays are becoming easier to adapt and more relevant for researchers end goals."

Creating a Better Way

Hilary Sherman, a Corning Applications Scientist created a 3D culture model to observe immune cell and tumor cell interactions by combining the Corning® 96-well spheroid microplate with the



high throughput screening (HTS) Transwell®-96 tissue culture system. Coated with Corning's unique Ultra-Low Attachment surface, the spheroid microplate allows for highly reproducible, single multicellular tumor spheroids to form in each well. The HTS Transwell permeable support systems are commonly used in migration and invasion studies, with a flat bottom receiver plate for monolayer culture. This novel combination of both products allows researchers to observe changes to tumor behavior, and by adding immune cells to the Transwell insert, to observe their rates of migration through the membrane and towards the tumor spheroids in the receiver plate. By combining both technologies one can study how these two cell types interact in a more *in vivo*-like environment.

Adds Gitschier, "We frequently utilize both HTS Transwell systems and the spheroid microplates in our Applications lab. When




Figure 1. HTS Transwell 96-well permeable supports were placed into 96-well spheroid microplates for 3D immune cell migration and tumor invasion assays.

Hilary [Sherman] realized the two products fit together perfectly and could be used for many applications requiring a migration/invasion component, as well as a 3D cell culture model component, it was very clear this would be very useful. This product synergy will provide a more realistic model to study more *in vivo*-like cell behavior and interactions.”

Rather than using the flat-bottom Transwell® receiver plate, Sherman placed the spheroid microplate under the HTS Transwell insert plate. This enabled her to more accurately observe tumor cell activity in a 3D system while facilitating migration of immune cells towards the tumor cell model. Doing so allows for the investigation of immune cell homing, tumor toxicity, and tumor immune evasion in an easy-to-use, reproducible, 3D high throughput assay.

Study Methods

To form the tumor spheroid models, cells from a cancerous lung tissue cell line (A549 cells) were seeded into 96-well spheroid microplates. The following day, Sherman placed the HTS Transwell 96-well permeable supports in the spheroid microplates. To assess immune cell homing and migration, natural killer (NK) cells were

added into the inserts and were allowed to migrate for 24 hours. (See Figure 1.)

Implications on Future Research

It is well understood that in order to effectively research the body’s response to targeted immunotherapy, a model must be used that successfully allows for the immune cells to reach the tumor cells of interest. Sherman’s first-of-its-kind research model confirms that not only are the NK cells in this study able to reach the target tumor cells, but they are also able to infiltrate the 3D spheroid structure and lead to the desired cytotoxic outcome. This novel system is a more streamlined model to observe all necessary components of this immunotherapy research in a single, high throughput, reproducible *in vivo*-like format.

“This model allows researchers to study immune system–tumor interactions in a more comprehensive way that also has the benefit of being high throughput and easy to use. The system can also be used to study other immune and cancer cells as well as the ability to form more complex models by adding endothelial or glial cells to the Transwell membrane to create more advanced blood brain barrier models,” says Gitschier.

Learn more about this new research model, request the Corning application note, CLS-AN-425 from your VWR Life Science Specialist.

Offering Solutions for Cell Culture Success

Since the development of the first cell culture flask, Corning has been a leader in cell culture innovation. And our dedication to this innovation has not wavered. From the development of Corning Matrigel® matrix 30 years ago, to breakthrough technologies such as our automation-friendly HTS Transwell inserts or our spheroid microplate with Ultra-Low Attachment surface, enabling uniform and reproducible 3D multicellular spheroid formation, it is our goal to create innovative 3D cell culture solutions. These solutions allow you to produce optimal environments for growing cells that exhibit *in vivo*-like behaviors and functionality — across all aspects of research.

Description	Size	Cat. No.	Unit
Ultra-Low Attachment Surface Spheroid Microplate, Black/Clear	96-Well	10185-094	Case of 5
Ultra-Low Attachment Surface Spheroid Microplate, Black/Clear	384-Well	10185-096	Case of 5
HTS Transwell Permeable Support Systems and Plates, PC, 0.3 µm	96-Well	89089-932	Case of 8
HTS Transwell Permeable Support Systems and Plates, PC, 0.4 µm	96-Well	89089-924	Case of 5
Corning® Matrigel® Matrix, Growth Factor Reduced	10 mL	47743-720	Case of 1
Corning® DMEM [+/-] 4.5 g/L Glucose, L-Glutamine, Sodium Pyruvate	1 L	45000-306	Pack of 6



VWR SCIENCE PORTAL

YOUR GATEWAY TO RESEARCH RESOURCES



The Lab is Always Open

- Help available 24/7 from the Research Concierge™
- View latest tools, technologies, and services

Simplified Ordering

- Minimize your search time
- Consolidate thousands of suppliers
- Reduce approval process time

Over 3,500 Research Services Available

- Assay Development
 - Cell & Tissue Histology
 - Custom Antibody Production
 - Expression Profiling
 - Next Generation RNA/DNA Sequencing
 - Biochemical or Cell-Based Assays
 - Chemical Analysis
 - Custom Cell Line Creation
 - Gene Editing, CRISPR/Cas9
 - Peptide Synthesis
 - Biomarker Discovery
 - Compound Synthesis
 - Epigenetics
 - Genotyping
 - Protein Expression
- And many more!

How It Works

- Visit vwr.com/vwrscienceportal and click "Continue" to begin your search for the services or products needed for your research
- Define your scope of work and receive quotes from suppliers
- Choose the supplier that matches your criteria and add the services to your basket

Access these unique features of Scientist.com through the VWR Science Portal:

 <p>Research Concierge™ Get help from MS/Ph.D.scientists and sourcing experts</p>	 <p>Innovation Hub™ Browse the latest cutting-edge tools, technologies, and services</p>
 <p>Supplier Ratings View supplier ratings and reviews</p>	 <p>Simplification Obtain approvals quickly (IACUC, BioBank, IP, etc)</p>

Research Concierge™ and Innovation Hub™ are Trademarks of scientist.com



Get started today at vwr.com/vwrscienceportal



Efficient Transient Gene Expression in Neurosphere with jetMESSENGER™ mRNA Transfection Reagent

By Barbara Demeneix, Karine Le Blay, and Sylvie Remaud [Muséum National d'Histoire Naturelle (MNHN), Paris, France], Maxime Dumont, Patrick Erbacher, Géraldine Guérin-Peyrou, and Alengo Nyamay'Antu (Polyplus-transfection®, Illkirch, France)




Fig. 1: Comparison between mRNA transfection using jetMESSENGER™ and DNA transfection using jetPRIME®.

Introduction

In light of their remarkable properties of both self-renewal and differentiation, adult neuronal stem cells are a unique cellular system that has demonstrated promising potential in developmental biology, regenerative medicine, and disease therapy. Cultivating adult neuronal stem cells *in vitro* as non-adherent free-floating spheroids is a frequently used 3D cell culture model to mimic the *in vivo* niche. Growing neural stem cells as neurospheres is an established approach to studying molecular mechanisms that drive

self-renewal and differentiation. Gene expression strategies to examine the role of proteins in modulating these cellular mechanisms rely on the ability to induce ectopic or endogenous gene overexpression.

As neurospheres are composed of both stem cells and progenitor cells, it is important to find a gene delivery method that will target all cell types, irrespective of their self-renewal capacity and consequently of their division rate. As an alternative to viral-based transfection methods which could give rise to

undesirable random genomic integration and subsequent potential pathogenicity, Polyplus-transfection® has developed an mRNA transfection reagent, jetMESSENGER, for difficult-to-transfect cells. Messenger RNA transfection is a growing alternative transfection method to plasmid DNA delivery which allows reproducible and highly efficient transfection in a wide variety of cells by circumventing the need for nuclear import of DNA, hindered in cells that are dividing slowly or have stopped altogether¹. This process has many advantages over DNA transfection, including a high percentage

of transfected cells irrespective of their division rate and earlier and more controlled protein expression with no risk of genome integration. Concomitantly, toxicity and immune response usually induced by mRNA entry into the cytosol, has been reduced with the latest generation of modified mRNA.

Here, we propose a mRNA transfection based-method to reach high transfection efficiencies in neuronal stem cells grown in spheroids, based on our combined expertise on neuronal stem cells and the novel mRNA transfection reagent jetMESSENGER™ developed by Polyplus-transfection®. The aim of this study was two-fold. Firstly, to test mRNA as an alternative method to DNA transfection for neuronal stem cells. Secondly, to determine the transfection efficiency of the heterogeneous population of neurosphere-forming stem cells including neuroblasts, oligodendrocyte progenitor cells, as well as other cell types, such as astrocytes and glial cells (Fig 1.).

Material & Methods

Neurospheres were transfected with either EGFP mRNA using jetMESSENGER™ (Polyplus-transfection®) or Lipofectamine®-2000 (Thermo Fisher Scientific®), or with DNA using jetPRIME® (Polyplus-transfection®) according to the manufacturers' recommendations. Neurospheres were gently dissociated and plated together at equal densities of 40,000 cells per well of a 24-well plate in complete medium. Twenty four hours post-seeding of cells, transfection was performed by adding 0.5µg of mRNA. The GFP encoding mRNA was used with respectively a mRNA: reagent ratio of





Fig. 2: Comparative transfection efficiency of mRNA versus DNA transfection reagents

1:2.4 and 1:3 ratio for jetMESSENGER™ and Lipofectamine®. For DNA transfection, 0.5µg of EGFP plasmid DNA was complexed with 1.5µL of jetPRIME® per well. EGFP expression in transfected cells was examined by fluorescence microscopy from 24 hours and up to 7 days post-transfection.

Results


First, we analyzed the transfection efficiency of both jetMESSENGER™ and Lipofectamine®-2000 transfection reagents. To this end, transfected cells were fixed at different time points following transfection (from 24 hours to 7

days). For both transfection reagents, we detected EGFP expression 24 hours post-transfection. Both continued to provide high levels of EGFP expression up to 7 days post-transfection (Fig 2A).

To evaluate more precisely the transfection efficiency, we quantified the percentage of EGFP+ transfected cells 48 hours post-transfection. We observed a 4-fold increase in transfected cells using either jetMESSENGER™ (40.93 ± 1.57) or Lipofectamine®-2000 (44.56 ± 2.51) compared to jetPRIME (13.15 ± 1.59) (Fig 2B), showing that mRNA transfection is more efficient for transfection of

differentiated neuronal cells from neurospheres than DNA transfection. Since neurospheres contain a mix of cell types, we studied which cell type is preferentially targeted by each

transfection reagent. To this end, immunohistochemistry using antibodies directed against markers expressed in NSCs (GFAP and SOX2) or astrocytes (GFAP) was performed. For each reagent, the EGFP signal was mostly detected in the cytoplasm of (i) NSC co-expressing SOX2 and GFAP and (ii) SOX2+ progenitors (Fig 3A).



jetMESSENGER™ transfects more efficiently both NSCs and progenitors (53.37 ± 2.55 and 23.04 ± 2.81) compared to either Lipofectamine®-2000 (42.55 ± 2.32 and 13.44 ± 2.09) or jetPRIME® (39.42 ± 4.73 and 18.53 ± 3.67). Moreover, jetMESSENGER™ privileges transfection of NSCs and progenitors as opposed to astrocytes, GFAP+ astrocytes showing only 16.92 ± 1.9 transfection efficiency as opposed to either Lipofectamine® 2000 (29.46 ± 2.2) or jetPRIME® (28.42 ± 4.49) (Fig 3B). Thus, jetMESSENGER™ shows more specific targeting of NSC/progenitors than Lipofectamine® 2000.

Conclusion & Perspectives

Here we show that mRNA transfection is an efficient technique to express exogenous proteins in neurospheres. Indeed, mRNA transfection leads to significantly higher transfection efficiency than DNA transfection, by a 3-fold increase.

To target different cell types within the neurosphere, such as neuronal stem cells and progenitors, mRNA delivery using jetMESSENGER™ was the most efficient approach.

In the future, we will use an mRNA delivery approach for our gene expression experiments in neurospheres using jetMESSENGER™ as it is the most efficient way to target our cells of interest, neuronal stem cells, and progenitors.

Reagent Size	Cat No.	Each
1 mL	75806-306	10 mL
0.75 mL	75806-224	60 mL
1.5 mL	75806-226	2 x 60 mL

Fig. 3: mRNA transfection using jetMESSENGER™ targets preferentially both NSCs and progenitor-derived neurospheres.

VWRbioMarke Supports Life Science Research


The VWRbioMarke Program features premier Suppliers and Brands that support Life Science applications and workflows in Cellular Science, Genomics, Proteomics and Microbiology. These suppliers have a proven track record of delivering industry leading technologies and products. They are committed to working closely with VWR Representatives and Life Science Specialists to provide you, our customers, with technical support, early access to new technologies, and value offerings.

Find more information on the VWRbioMarke Suppliers, their product offerings, technical articles, and more at vwr.com/vwrbiomark.

GO VIOLET AND GO SMALL.



HARNESS THE POWER OF 405 NM LASER FOR NANOPARTICLE DETECTION.



The boundaries for flow cytometry are being pushed by the need to measure and evaluate characteristics of smaller particles. The CytoFLEX Flow Cytometer has the resolution to detect particles down to 200nm while also collecting information in the fluorescent channels to help you define cellular origin.

Description	Cat. No.	Unit
CytoFLEX System, B4-RO-VO	76183-344	Each
CytoFLEX System, B3-R1-VO	76183-346	Each
CytoFLEX System, B2-RO-V2	76183-348	Each
CytoFLEX System, B2-R2-VO	76183-350	Each




CELLSTAR® Cell-Repellent Surface

Cell Culture Vessels for Suspension and Spheroids Used in 3D Cell Culture Models

Product Applications:

- Spheroid cultures of tumor cells and non-cancerous cells
- Aggregation of stem cells
- Suspension culture of semi-adherent and adherent cell lines



Descriptions	Cat. No.	Unit
250 mL Cell Culture Flask, Cell-Repellent Surface, PS, Sterile, White Filter Screw Cap	10836-541	Case of 15
100x20 mm Dish, Cell-Repellent Surface, PS, Clear, Sterile, with Lid	30618-024	Case of 5
6-Well Plate, Cell-Repellent Surface, PS, Clear, Sterile, with Lid	30618-022	Case of 5
24-Well Multiwell Plate, Cell-Repellent Surface, PS, Clear, Sterile, with Lid	10014-320	Case of 5
96-Well Plate, Cell-Repellent Surface, PS, Round U Bottom, Clear, Sterile, with Lid	30618-026	Case of 6
96-Well Plate, Cell-Repellent Surface, PS, Flat Bottom, Chimney Style, Clear, Sterile, with Lid	30618-028	Case of 6
96-Well Plate, Cell-Repellent Surface, PS, Flat Bottom, Black µClear, Sterile, with Lid	10014-318	Case of 32
384-Well Plate, Cell-Repellent Surface, PS, U Bottom, Clear, Sterile, with Lid	76199-198	Case of 8

CELLCOAT® — Protein Coated Cell Culture Vessels

Now Available in Small Packaging Sizes!

- Poly-D-Lysine — Synthetic protein
- Collagen type I — Extracellular matrix protein

Applications:

- Improved adhesion
- Improved cell proliferation of primary cell lines
- Reduced-serum or serum-free cultivation




Product Descriptions	Cat. No.	Unit
CELLCOAT® TC-Flask, 250 mL, 75 CM ² , Poly-D-Lysine	75796-006	Pack of 5
CELLCOAT® TC-Flask, 250 mL, 75 CM ² , Collagen I	75795-996	Pack of 5
CELLCOAT® TC-Plate, 6-Well, Poly-D-Lysine	75796-008	Pack of 5
CELLCOAT® TC-Plate, 12-Well, Poly-D-Lysine	75795-982	Pack of 5
CELLCOAT® TC-Plate, 24-Well, Poly-D-Lysine	75795-992	Pack of 5
CELLCOAT® TC-Plate, 6-Well, Collagen I	75795-998	Pack of 5
CELLCOAT® TC-Plate, 12-Well, Collagen I	75796-000	Pack of 5
CELLCOAT® TC-Plate, 24-Well, Collagen I	75796-002	Pack of 5

Compact Dry™

There's Always a Better Way

- Rigid transparent frame with removeable lid keeps your colonies intact when subculturing
- Self-diffusing dry media transforms into gel with the 1mL sample — no need for an extra step
- 18-month shelf life, room temperature storage, and stackable

The Compact Dry plate line presents a variety of organism-specific and total count tests, packaged up nicely in a sophisticated and sturdy cassette. This product line was designed for the end user with the intent to save technicians time while performing high-throughput colony counts on food, beverage, cannabis, cosmetic, environmental, meat, and raw material samples.



Description	Application	Cat. No.	Unit
Compact Dry EC	EC — <i>E. coli</i> and Coliforms	10145-954	Pack of 100
CompactDry ETB	ETB — <i>Enterobacteriaceae</i>	10753-986	Pack of 100
CompactDry LS	LS — <i>Listeria spp.</i>	10789-408	Pack of 100
CompactDry X-SA	X-SA — <i>Staphylococcus aureus</i>	10145-970	Pack of 100
CompactDry SL	SL — <i>Salmonella spp.</i>	10789-458	Pack of 100
CompactDry YM	YM — Yeast and Mold	10145-972	Pack of 100
CompactDry YMR	YMR — Yeast and Mold Rapid Test	76076-714	Pack of 100
CompactDry TC	TC — Total Plate Count	10145-968	Pack of 100

Dilucup™ Elegance

Dilucup Saves Time and Money Spent on Preparing Your Own Dilution Vials!

- The Dilushaker™ uses orbital rotations to make sure each diluent is thoroughly vortexed the same way, every time
- Simplify your serial dilution process with attached, pre-filled vials that are easy to open
- Save bench space with this small unit and refrigerator space with Dilucups that are easily stackable and store at room temperature

This automated, vortexing system maximizes lab efficiency with pre-filled, sterile, and ready-to-use dilution cups. Bringing automation to your vortexing process virtually eliminates variable results from tech to tech. The Dilucup line also reduces the occurrence of repetitive motion injuries with its ergonomic design.



Description	Cat. No.	Unit
Dilushaker III Digital, 3 Rows, LED, 21 Cups, 7x3	75860-388	Each
Dilushaker III Digital, 6 Rows, LED, 42 Cups, 7x6	75860-392	Each
Dilucup Elegance Solutions		
Dilucup Elegance BPW, 32 Trays, 7 x 3 Cups, 9 mL	75860-416	Each
Dilucup Elegance BPW, 16 Trays, 7 x 6 Cups, 9 mL	75860-418	Each
Dilucup Elegance MRD, 32 Trays, 7 x 3 Cups, 9 mL	75860-406	Each
Dilucup Elegance MRD, 16 Trays, 7 x 6 Cups, 9 mL	75860-408	Each
Dilucup Elegance NaCl, 32 Trays, 7 x 3 Cups, 9 mL	75860-412	Each
Dilucup Elegance NaCl, 16 Trays, 7 x 6 Cups, 9 mL	75860-414	Each

These specific products are not available in Canada. Please contact your VWR Sales Rep to learn about similar options available in your region.

sparQ DNA Frag and Library Prep Kit

Cost-Effective Enzymatic Library Prep Solution with Industry-Leading Performance

- High quality libraries under 3 hours from 1ng– 1µg of input DNA
- Tunable and reproducible frag profiles to accommodate different applications
- Streamlined workflow with exceptional sensitivity and high yield
- Novel library chemistry and high-fidelity amplification minimizing bias

sparQ DNA Frag and Library Prep Kit provides rapid enzymatic fragmentation and construction of DNA libraries from double-stranded DNA for sequencing on Illumina NGS platforms. Using a simple two-step protocol, high-quality libraries can be completed in less than 3 hours with minimal hands-on time. The optimized chemistry coupled with high fidelity amplification results in better library yield, lower duplication bias, and consistent coverage across a broad range of GC content.



Description	Size	Cat. No.	Unit
sparQ DNA Frag and Library Prep Kit	24 Reactions	76183-242	Each
sparQ DNA Frag and Library Prep Kit	96 Reactions	76183-244	Each
sparQ Adapter Barcode, Set A	12 Single Index Barcodes for 96 Reactions	76169-154	Each
sparQ Adapter Barcode, Set B	12 Single Index Barcodes for 96 Reactions	76169-156	Each

dsDNA Quantitation Kits for Qubit® Fluorometers and Fluorescent Plate Readers



Ultrasensitive, Accurate, and Versatile Options for the Widest Range of DNA Concentrations

- **Versatile** — Choose from kits with different detection limits, for ultra-high sensitivity or broad dynamic range
- **Superior sensitivity** — AccuBlue® NextGen detects as little as 1–5 pg DNA in microplate assay, outperforming PicoGreen®
- **Economical** — Seamlessly switch from Qubit dsDNA HS and BR assay kits to the equivalent AccuGreen™ kits and SAVE (as much as 25%)
- **Robust** — Tolerant to common contaminants including free nucleotides, protein, detergents, salts, and contaminating RNA

AccuBlue, AccuGreen, and AccuClear® assays use highly sensitive and selective fluorescent dsDNA binding dyes for precise quantitation of purified DNA with greater accuracy than absorbance-based methods. Easy-to-use kits are available for use on any instrument, with options for a wide range of DNA concentrations, and are compatible with dsDNA from a variety of sources.



Description	Size	Cat. No.	Unit
AccuGreen High Sensitivity dsDNA Quantitation Kit (for Qubit)	100 Assays	75845-718	Each
AccuGreen High Sensitivity dsDNA Quantitation Kit (for Qubit)	500 Assays	75845-716	Each
AccuGreen Broad Range dsDNA Quantitation Kit (for Qubit)	500 Assays	76204-710	Each
AccuClear Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards	200 Assays	89493-596	Each
AccuClear Ultra High Sensitivity dsDNA Quantitation Kit with DNA Standards	1000 Assays	89427-078	Each
AccuBlue NextGen dsDNA Quantitation Kit with DNA Stds.	500 Assays	76204-714	Each
AccuBlue High Sensitivity dsDNA Quantitation Kit with DNA Standards	1000 Assays	89139-008	Each
AccuBlue Broad Range dsDNA Quantitation Kit with DNA Standards	1000 Assays	89139-010	Each



We Set Science in Motion to Create a Better World.

From breakthrough discovery to agile delivery of innovative products and services, Avantor® and VWR are one company.

Our shared vision: to be your trusted partner, providing discovery-to-delivery solutions, high-performance materials and expert support for your unique personal protection and facility requirements.

Learn more about our new company at
www.SettingScienceinMotion.com

High-Yield, Scalable Library Preparation with the NEBNext® Ultra™ II FS DNA Library Prep Kit

Improving Performance, Ease of Use and Reliability of Enzymatic DNA Fragmentation

Lynne Apone, Vaish Panchapakesa, and Karen McKay
New England Biolabs, Inc.

Introduction

The continued expansion of the use of next gen sequencing depends in large part on overcoming the limitations and bottlenecks associated with high-quality library preparation. The requirement for numerous steps and expensive equipment can result in sample loss, errors, and limited throughput. To address these issues, we have built upon our NEBNext Ultra II DNA library prep workflow to create a kit that includes a fragmentation system: the NEBNext Ultra II FS DNA Library Prep Kit. This kit integrates a new enzymatic DNA fragmentation reagent into the library prep workflow, such that fragmentation is combined with end repair and dA-tailing. The combination of these reactions eliminates the need for equipment to shear DNA mechanically and also reduces the number of sample transfers and losses. Subsequently, adaptor ligation is also carried out in the same vial, after which a single cleanup step is performed. For low input samples, PCR amplification is performed prior to sequencing.

Importantly, the enzymatic shearing of DNA with the FS kit does not introduce bias into the library, and the kit is suitable for input DNA from the full range of GC content. The reduced sample loss and increased efficiencies of the workflow enable use of lower input amounts, with a range of 100pg–0.5µg, and insert sizes of 100bp to 1kb can be generated.

Workflow

The NEBNext Ultra II FS DNA Library Prep Kit protocol, including

	Fragmentation/ End Repair/ dA-Tailing	Adaptor Ligation	Clean Up/ Size Selection	Amplification	Clean Up	Total Workflow
Hands-On	2 min.	1 min.	5 min.	0–1 min.	0–5 min.	8–14 min.
Total	37–62 min.	16–31 min.	27–37 min.	0–34 min.	0–27 min.	1.3–3.2 hr.

Figure 1. NEBNext Ultra II FS DNA workflow.

fragmentation, is fast (~2.5 hours) and simple, and can accommodate 100pg–0.5µg of input DNA (Figure 1). Fragmentation, End Repair, and dA-Tailing reagents are combined, and there is no clean-up before adaptor ligation. The same fragmentation protocol is followed for all input amounts and for all GC contents. The kit also includes options for PCR-free workflows. The protocol is compatible with adaptors and

primers from the NEBNext product line (“NEBNext Oligos”) or from other sources.

Increased Library Yields

The use of enzymatic fragmentation can result in higher library yields than workflows incorporating mechanical shearing of DNA, due to both reductions in sample loss and decreased DNA damage. The NEBNext Ultra II FS kit further




Figure 2. NEBNext Ultra II FS DNA Produces the Highest Yields From a Range of Input Amounts

Libraries were prepared from Human NA19240 genomic DNA using the input amounts and numbers of PCR cycles shown. For NEBNext Ultra II FS, a 20-minute fragmentation time was used. For Kapa™ HyperPlus libraries, input DNA was cleaned up with 3X beads prior to library construction, as recommended, and a 20-minute fragmentation time was used. Illumina® recommends 50ng input for Nextera®, and not an input range; therefore, only 50ng was used in this experiment. “Covaris” libraries were prepared by shearing each input amount in 1X TE Buffer to an insert size of ~200bp using a Covaris instrument, followed by library construction using the NEBNext Ultra II DNA Library Prep Kit. Error bars indicate standard deviation for an average of 3–6 replicates performed by 2 independent users.




Figure 3. Consistent and Reliable Library Preparation with NEBNext Ultra II FS DNA

Libraries were prepared from Human NA19240 genomic DNA using the input amounts shown. NEBNext Ultra II FS libraries were prepared using a 20-minute fragmentation time. For Kapa HyperPlus, input DNA was cleaned up with 3X beads prior to library construction, as recommended, and a 20-minute fragmentation time. Library size was assessed using the Agilent Bioanalyzer. Low input (1ng and below) libraries were loaded on the Bioanalyzer without a dilution. High input libraries were loaded with a 1:5 dilution in 0.1X TE.

increases library yields through the integration of the fragmentation reagent with end repair and dA-tailing reagents (thereby minimizing loss during transfer steps), by not requiring a clean-up step before adaptor ligation, and through the high reaction efficiencies of each step in the workflow.

Achieving sufficient library yields for high quality sequencing from very low input amounts can be especially challenging with mechanical shearing of DNA, a situation compounded by the preference to amplify libraries using as few PCR cycles as possible. The NEBNext Ultra II FS kit overcomes this low-input challenge, and users can now obtain higher library yields with input amounts as low as 100pg of human

genomic DNA, as shown in Figure 2 (below, left).

Robustness of DNA Fragmentation

Consistent and reliable fragmentation of DNA is critical for an enzymatic method to be adopted. Fragmentation using the NEBNext Ultra II FS kit is time-dependent, with fragmentation times amenable for manual or automated library preparation. Importantly, the final library size reflects the initial fragmented DNA, demonstrating that sufficient fragments of the desired size were produced during the fragmentation step.

Since in practice the exact quantity of DNA in a sample, and the exact GC content, may be unknown, it can be challenging when

different fragmentation protocols are required for different input amounts and GC contents. Additionally, the ability to use input DNA in a range of buffers greatly simplifies the very start of the workflow, especially in situations where the input DNA buffer composition may be unknown or uncertain. The Ultra II FS kit addresses all of these issues by requiring a single fragmentation protocol for the full range of input amounts (100pg–0.5µg) (Figure 3) and for the full range of GC content (see next section). Additionally, DNA can be in water, Tris, 0.1X TE or 1X TE.

Uniform GC Coverage

While sufficient yield of a library is required for successful sequencing, the quality of a library is also critical, regardless of the




Figure 4. NEBNext Ultra II FS DNA Provides Uniform GC Coverage for Microbial Genomic DNA Over a Broad Range of GC Composition. Libraries were prepared using 1ng of a mix of genomic DNA samples from *Haemophilus influenzae*, *Escherichia coli* (K-12 MG1655), *Rhodopseudomonas palustris* and the library prep kits shown, with 9 PCR cycles for consistency across samples, and sequenced on an Illumina MiSeq. NEBNext Ultra II FS libraries were prepared using a 20-minute fragmentation time. For Kapa HyperPlus libraries, input DNA was cleaned up with 3X beads prior to library construction, as recommended, followed by a 25-minute fragmentation time. "Covaris" libraries were prepared by shearing 1ng of DNA in 1X TE Buffer to an insert size of ~200bp using a Covaris instrument, followed by library construction using the NEBNext Ultra II DNA Library Prep Kit. Reads were mapped using Bowtie 2.2.4 and GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library.

input amount or GC content of the sample DNA. A high-quality library will have uniform representation of the original sample, including even coverage across the GC spectrum.

For enzymatic DNA fragmentation methods, randomness with regard to GC content can be of concern, and ensuring uniform GC coverage in libraries produced with the NEBNext Ultra II FS kit was of great importance. The Ultra II FS kit shows not only consistent uniformity of GC coverage at the full range of input amounts, but uniformity of GC coverage superior to mechanical shearing workflows

and alternative enzymatic shearing methods (Figure 4, next page).

Library Quality

As noted above, an ideal library will represent completely and proportionally the sequence of the sample DNA, regardless of the input amount and GC content of the sample DNA. Especially when input amounts for a library are low, there is a risk that the resulting library will lack this diversity, and that some sequences will be over- or under-represented. This could be due to preferential amplification of some sequences during PCR of the final library, or non-random DNA fragmentation

at the beginning of the library prep process. Comparison of the level of sequence coverage, in 10kb intervals, achieved for libraries with different input amounts is a useful measure of the diversity of a library. This is especially useful when the comparison is with a PCR-free library made with a higher input amount, as such a library is typically not affected by the factors that lead to bias in libraries. We see good correlation of a PCR-free library prepared with 100ng of Covaris-sheared human genomic DNA, and amplified libraries prepared with 1ng of Covaris-sheared DNA or with the NEBNext Ultra II FS kit.

PCR-Free Libraries

Construction of a library using a PCR-free workflow removes the risk of incorporation of bias during library amplification. However, the input amounts required to produce sufficient amounts of a high-diversity library without an amplification step are necessarily higher, and this can be limiting. With the NEBNext Ultra II FS protocol, by performing fragmentation, end repair and dA-tailing in the same reaction vial, sample loss due to transfer or clean-up steps after DNA fragmentation are eliminated, and yields are higher. Additionally, the use of Ultra II FS enzymatic DNA fragmentation rather than mechanical shearing reduces DNA damage, and consequently increases library yields, making the omission of the PCR step now feasible for lower nanogram level input amounts.

Variant Detection

Enzymatic DNA fragmentation enables increased library yields in comparison to mechanical shearing methods, in part due to the reduced DNA damage that occurs during shearing. Additionally, sequence markers indicative of oxidative damage to DNA that have been found in libraries constructed with mechanically-sheared



DNA, are absent in libraries constructed using the NEBNext Ultra II FS kit, and greater differences are seen with lower input amounts. This highlights the higher quality of libraries constructed with the Ultra II FS kit compared to Covaris-sheared DNA libraries, especially at low input amounts.

Conclusion

The NEBNext Ultra II FS DNA Library Prep Kit for Illumina provides a simple and reliable solution for DNA fragmentation integrated with library construction. The

kit enables production of high quality libraries from a broad range of input amounts and GC contents with a single fragmentation protocol, greatly simplifying the scalability of library construction:

- Perform fragmentation, end repair, and dA-tailing with a single enzyme mix
- Obtain robust fragmentation with a single protocol, regardless of DNA input amount or GC content
- Generate high quality libraries from 100pg–0.5µg input DNA
- Use with DNA in standard buffers (TE, Tris-HCl) and water

- Save time with a streamlined workflow: (~ 2.5 hours), with < 15 minutes hands-on time
- Experience reliable fragmentation, even with very low input amounts
- Generate high yields with increased reaction efficiencies and minimized sample loss
- Vary incubation time to generate a wide range of insert sizes

References

For a complete list of references, visit vwr.com/vwrbiomarke.



Description	Size	Cat. No.	Unit
NEBNext® Ultra™ II FS DNA Library Prep Kit for Illumina®	24 rxns	103529-140	Each
NEBNext® Ultra™ II FS DNA Library Prep Kit for Illumina®	96 rxns	103529-168	Each
NEBNext® Ultra™ II FS DNA Library Prep with Sample Purification Beads	24 rxns	103529-166	Each
NEBNext® Ultra™ II FS DNA Library Prep with Sample Purification Beads	96 rxns	103529-164	Each
NEBNext® Ultra™ II FS DNA Module	24 rxns	103529-144	Each
NEBNext® Ultra™ II FS DNA Module	96 rxns	103529-142	Each

These products will be available soon in Canada. Please contact your VWR Sales Rep to learn about similar options available in your region in the meantime.


VWR® TRANSFECTION SELECTOR

Search hundreds of transfection reagents for bioproduction, *in vivo* work, viral production, DNA and siRNA transfection, CRISPR/Cas 9, and more!

SEARCH for specialized transfection reagents

ADDRESS rare cell cultures, neurons, and primary cells

FIND additional products such as electroporation and reagents



vwr.com/transfectionselection

ATTN: POST OFFICE: Electronic Service Requested



Prices and product details are current when published, subject to change without notice. Certain products may be limited by federal, state, provincial, or local regulations. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2018 VWR International, LLC. All rights reserved.

ATTN:
MAILROOM:
If addressee is no longer with the organization, please deliver to laboratory products buyer. Thank you.

0818 35M Lit. No. 930424



IgG Antibody Selector...

Search | Select | Simple



vwr.com/antibody