PUREZERO* HG5 NITRILE GLOVES

12 in. / 30.5 cm

Ambidextrous

Accelerator Free*

Low Dermatitis Potential

HALYARD* PUREZERO* HG5 White Nitrile Gloves

- 12 in./30.5 cm Length
- Accelerator Free* • Low Dermatitis Potential
- Ambidextrous
- Textured
- For the Cleanroom Environment
- For Industrial Use Only

NOTICE: THIS INSERT SHOULD BE FURNISHED OR MADE AVAILABLE TO THE USERS OF THESE GLOVES AS A SAFETY PRECAUTION.

Product certified to the Regulation (EU) 2016/425 ("PPE Regulation") by Notified Body 2797 and to the Regulation 2016/425 on personal protective equipment, as amended to apply in GB by Approved Body 0086.

Warning: Gloves offer protection against chemicals (Splash) and micro-organisms.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in the physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly.

For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Before use, inspect the gloves for defects or imperfections such as holes and tears. If the gloves are ripped or puncture during use, dispose of them immediately.

Always refer to your organization's procedures for donning, doffing, and usage. First, wash or disinfect hands. Make sure they dry completely prior to donning. To don glove, insert hand into cuff and carefully pull up over hand and wrist. Repeat with the second glove. Adjust as needed for comfort and dexterity. To remove, grasp the exterior of the glove carefully pulling the glove off the hand, turning it inside-out. Put the balled up glove in your opposite hand, then repeat the previous step. Discard appropriately. When donned/doffed correctly, gloves are meant to protect the user, and no interface issues with other PPE should exist. Please refer to halyardhealth.com/information for additional recommendations for proper donning/doffing practices.

*Not formulated with the following vulcanizing chemicals and accelerators: Sulfur, Thiurams, Thiazoles, Guanidines, Carbamates

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ISO 21420:2020 Dexterity Classification = 5

HALYARD* PUREZERO* HG5	
	White Nitrile Gloves
XS (5.5)	CLN30512XS
S (6.0 - 6.5)	CLN30512SM
M (7.0 - 7.5)	CLN30512MD
L (8.0 - 8.5)	CLN30512LG
XL (9.0-10.0)	CLN30512XL

<https://halyardhealth.info/cleanroom>

CE 2797
(PPE Cat. III)

UK
CA
0086

O&M Halyard, Inc., 9120 Lockwood Blvd., Mechanicsville, VA 23116; In USA, please call 1-844-425-9273 • halyardhealth.com

CE REP ArcRoyal Unlimited Company, Virginia Road, Kells, Co Meath, Ireland

BE O&M Halyard Belgium BV, Berkenlaan 88, 1831 Machelen (Brab.), Belgium

UK RP: RLD Quality Limited Centenary House, Peninsula Park, Rydon Lane, Exeter, EX2 7XE, United Kingdom

Product certified to the Regulation (EU) 2016/425 by Notified Body 2797 BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, and to the Regulation (EU) 2016/425 on personal protective equipment, as amended to apply in Great Britain by Approved Body BSI Assurance UK Ltd (0086) Kitemark Court, Davy Avenue, Milton Keynes, MK5 8PP.

Sponsored in Australia by O&M Halyard Australia Pty Ltd., Tenancy 1, Level 1, Trinititi 3, 39 Delhi Road, North Ryde NSW 2113.

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*Registered Trademark or Trademark of O&M Halyard, Inc., or its affiliates. ©2018. All rights reserved.
Patent and/or Pending Application. See www.halyardhealth.com/patents. 2022-07-27

15-HM-040-0-01