VWR shiroGEL

Vertical 20

INSTRUCTION MANUAL

European Catalogue Number:

700-0353 VE 20 ShiroGel Vertical

700-0354 VE 20 Dual Casting System

Version: 1

Issued: 05.21.2014

Legal Address of Manufacturer

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011

Country of origin

UK

Contents:

	<u>Page</u>
Safety Information	1
Package Contents	2
Usage Guidance and restrictions:	3
Care and Maintenance	4
Introduction	5
Components	6
Key Features	7
Gel Casting	10
Set-up	11
Converting from 2 to 4 gel configuration	14
4 gel assembly for 1.5 and 2mm thick gels	16
Gel Preparation	17
Gel Selection	18
Gel Pouring	19
Sample Prep and Loading	20
Gel Running	23
Stock Solutions	24
References	26
Combs	27
Connecting to a Chiller Unit	28
Technical Service and Warranty	29
Equipment Disposal	30

Warning

When used correctly, these units pose no health risk. However, these units can deliver dangerous levels of electricity and are to be operated only by qualified personnel following the guidelines laid out in this instruction manual.

Anyone intending to use this equipment should read the complete manual thoroughly.

Safety Information

To avoid electrical shock:

The unit must never be used without the safety lid correctly in position.

To avoid electrical shock:

The unit should not be used if there is any sign of damage to the external tank or lid.

These units comply with the statutory CE safety directives:

Low Voltage Directive: 2006/95/EC	RoHS Directive: 2011/65/EU
EMC Directive: 2004/108/EC	WEEE 2012/19/EU
IEC 61010-1:2010 (Third Edition)	IEC 61326-1:2005

Avoiding Damage to the Instrument

The units should never come into contact with the following cleaning agents, these will cause irreversible and cumulative damage:-

Acetone, Phenol, Chloroform, Carbon tetrachloride, Methanol, Ethanol, Isopropyl alcohol, Alkalies.

Package Contents

VWR shiroGEL Vertical 20:

Units include tank, lid, internal module and electrodes and include the following accessories:

	Glass Plates	Combs	Casting base	Cooling Coil	Cables
shiroGEL Vertical 20	700-0361– Notched, Pk/2 700-0364– Plain with bonded 1mm spacers, Pk/2 700-0369 – Dummy Plate	2 of 700-0388 1mm thick, 24 sample	1 of 700-0356 (includes 1x 700-0359)	700-0358	1 set
shiroGEL Vertical 20 Complete	700-0361– Notched, Pk/2 700-0364– Plain with bonded 1mm spacers, Pk/2 700-0369 – Dummy Plate	2 of 700-0388 1mm thick, 24 sample	1 of 700-0356 (includes 1x 700-0359) + External CASTING STAND (no electrodes) 700-0355	700-0358	1 set

The packing lists should be referred to as soon as the units are received to ensure that all components have been included. The unit should be checked for damage when received.

Please contact your supplier if there are any problems or missing items.

Usage Guidance and restrictions:

- Maximum altitude 2,000m.
- Temperature range between 4°C and 65°C.
- Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.
- Not for outdoor Use.

Care and Maintenance:

Cleaning shiroGEL Vertical 20

Units are best cleaned using warm water and a mild detergent. **Water at temperatures above 60° C can cause damage to the unit and components.**

The unit should be thoroughly rinsed with warm water or distilled water to prevent build up of salts but care should be taken not to damage the plate electrodes and vigorous cleaning is not necessary or advised.

Air drying is preferably before use.

The units should only be cleaned with the following:

- *Warm water with a mild concentration of soap or other mild detergent.*

Compatible detergents include dishwashing liquid, Hexane and Aliphatic hydrocarbons

The units should not be left to in detergents for more than 30 minutes.

The units should never come into contact with the following cleaning agents, these will cause irreversible and accumulative damage:-

Acetone	Phenol	Chloroform	Carbon tetrachloride
Ethanol	Methanol	Isopropyl alcohol	Alkalis

RNase Decontamination

This can be performed using the following protocol:-

Clean the units with a mild detergent as described above.

Wash with 3% hydrogen peroxide (H₂O₂) for 10 minutes.

Rinsed with 0.1% DEPC- (diethyl pyrocarbonate) treated distilled water,

Caution: DEPC is a suspected carcinogen. Always take the necessary precautions when using. RNaseZAP™ (Ambion) can also be used. Please consult the instructions for use with acrylic gel tanks.

Instructions:

Introduction

The new shiroGEL VERTICAL 20 Vertical Electrophoresis System is the latest product innovation for large-format vertical gel electrophoresis. Designed to perform a variety of separations, including first- and second-dimension SDS-PAGE, native, preparative, gradient and high-resolution nucleic acid electrophoresis, plus capillary tube gel IEF and electroblotting, the shiroGEL VERTICAL 20 is one of the most versatile large-format vertical systems available.

By introducing innovative, new vertical screw-clamp technology within the PAGE insert only four screws are now necessary to secure as many 20x20cm gels. This gives the shiroGEL VERTICAL 20 the selective advantage of a much faster set up speed compared to competitor products whose traditional clamping configurations require as many as 24 screws to secure just two glass plates. In addition, the shiroGEL VERTICAL 20's vertical screw-clamp configuration distributes pressure evenly along the height of the gel rather than in the center to eliminate plate bowing and gel compression, but still maintains a leak-proof seal during casting; while the ergonomic shiroGEL VERTICAL 20 design of the PAGE insert, aids both handling and set up.


Whatever your requirements are the shiroGEL VERTICAL 20 can be made to meet them. Regardless of whether it is running 2 or 4 gels, electroblotting, and IEF using capillary tube gels or IPG strips, all of these techniques may be performed using the same versatile unit while retaining the benefits of large format electrophoresis, such as extended separation distances, greater sample throughput and superior resolution.

Before using the SHIROGEL VERTICAL 20 we recommend assembling and disassembling the unit before using it for electrophoresis, and familiarizing yourself with the component parts. Please see Figure 1 below.

Figure 1. shiroGEL VERTICAL 20

Vertical Component Parts

shiroGEL VERTICAL 20 Component Parts
1. Lid and power cables
2. PAGE insert
3. Vertical screw-pin
4. Clamping bars
5. Glass plates
6. Inner buffer chamber
7. Gasket
8. Detachable cooling coil
9. Outer tank
10. Cam-pin caster
11. Ultra-soft casting mat
12. Combs

Key Features

- **Vertical screw-clamp technology** - using only 4 screw-pins - dramatically reduces assembly time; also act in conjunction with color-coded clamping bars (see below) to distribute pressure evenly along the height of the gel to prevent plate bowing or compression (Fig.1, number 3)
- **A built-in inner buffer chamber within the PAGE insert** - allows set up to be completed without inclusion of a top tank or upper buffer chamber (Fig.1, number 6)
- **Color-coded clamping bars** – use green clamping bars to run up to two 2-mm-thick gels in the standard shiroGEL VERTICAL 20 configuration, or by using the thinner, yellow clamping bars to run up to four 2-mm-thick gels in the shiroGEL VERTICAL 20. (Fig.1, number 4; also see section on

Converting the shiroGEL VERTICAL 20 from a 2- to 4-gel Configuration)


- **Glass plates** - compress gently against a flat, level gasket to prevent current leakage from the inner buffer chamber during electrophoresis (Fig.1, number 5); notched glass plate is placed innermost and compresses directly against the gasket to allow buffer to fill the upper wells within the gel, thereby facilitating electrical connectivity between the negative electrode and the top of the gel, whereas the plain glass plate with bonded spacers is outermost to prevent buffer leakage from the inner chamber
- **Detachable inner cooling coil** - connects to the laboratory water supply or a recirculating chiller to provide uniform, smile-free electrophoresis, while allowing runs to be performed at higher voltage (Fig.1, number 8)
- **Deep gel tank** - with adequate clearance beneath the glass plates to allow a magnetic stirrer to maintain buffer recirculation and uniform pH (Fig.1, number 9)
- **Cam-caster** - with very forgiving ultra-soft silicone mat compensates for glass plate misalignment to ensure leak-free casting (Fig.1, number 10 & 11)

Other Features

- Notched glass plates with bonded spacers supplied with shiroGEL VERTICAL 20 systems allow gel capacity to be doubled from 2 to 4, while dummy plate allows single gels to be run (See Fig. 2)

Fig. 2. Glass plate set-up for 2-gel sandwich using yellow spacers


1. Plain glass plate with bonded spacer
2. Notched glass plate with bonded spacer
3. Notched glass plate

- Bonded spacers and combs color coded for thickness (not shown)
- Asymmetric lid design (Fig. 1, number 1) and color-coded screw pins in PAGE insert prevent polarity reversal (See Fig. 3)


Fig. 3. Color-coded screw pins prevent polarity reversal

1. Red vertical screw pin (1a.) corresponding to positive electrode plug (1b.) color-coded with **positive** power cable
2. Black vertical screw pin corresponding to negative electrode plug color-coded with **negative** power cable

Leak-free Casting Using Vertical Screw-Pin Technology

Fig. 4. The vertical screw-clamp technology of the shiroGEL VERTICAL 20 PAGE insert facilitates fast, leak-proof gel casting.

Gel Casting for 2 Gels Using the shiroGEL VERTICAL 20

General rules

1. Clean a set of glass plates for each gel first with distilled water and then with 70 % ethanol. One set of glass plates constitutes one notched glass plate and one plain glass plate with bonded spacers. When using a 2-gel glass plate sandwich in the shiroGEL VERTICAL 20 to cast 2 gels either side of the PAGE insert, two notched glass plates are required - one with and the other without bonded spacers - as well as a set of plain glass plates with bonded spacers (Please also see section on **Converting the shiroGEL VERTICAL 20 from a 2- to 4-gel Configuration**). The plain glass plate is positioned outermost, followed by a notched glass plate with bonded spacers and then a second notched glass plate (Please see **Fig. 2. Glass plate set-up for 2-gel sandwich using yellow spacers**).

N.B. - All glass plates, modules and casting base accessories must be completely dry during set – up. Wet components are more likely to misalign and cause leaks.

2. Assemble the glass plates so that the bottom of the glass plates and the spacers are perfectly aligned; the standard glass plates supplied with the shiroGEL VERTICAL 20 or shiroGEL VERTICAL 20 systems, are supplied with bonded spacers and do not require manual alignment. However, if using free spacers, which are not included as standard, these need to be perfectly aligned. This is best performed using a small spacer or comb to push the spacers apart.

NOTE: The glass plates with bonded spacers have an arrow in the top of the spacers which are slightly longer than the glass plate to indicate the top.

Please follow the illustrations below for detailed set up of shiroGEL VERTICAL 20 system


1. Assemble each gel cassette on a flat level surface, by placing the plain glass plate down with the spacers facing upwards followed by the notched glass plate.

2. Loosen the vertical screw-pins in the PAGE insert to release the locking mechanism, allowing the gel clamps to sit in the resting slots.


3. Insert a gel cassette into each side of the inner buffer chamber in the PAGE insert, and begin tightening the vertical screw-pins.


4. Continue to tighten the screw-pins until the gel clamps glide out of the resting slots and fix firmly against the glass plates pushing them upright.


5. Check the bottom of the glass plates to ensure that they are flush together, and place the PAGE insert on the casting base; make sure that the cams are facing downwards.

6. Insert cams and turn until tightened, drawing the PAGE insert onto the casting to form a leak-proof seal.

7. Pour in the gel solution, insert the combs and allow the wells to polymerize. Gel pouring is best performed using a 25ml or 50ml pipette.

8. Transfer the PAGE insert to the gel tank. Fill the inner and outer buffer chambers before loading samples.

9. Replace the lid, connect to the power supply and run.

Converting the shiroGEL VERTICAL 20 from a 2- to 4-gel Configuration

The thicker green sliding gel clamps are recommended to secure up to 2 gels (i.e. 1 gel either side of the PAGE insert) for gels up to a maximum thickness of 2mm. For 4 gels (i.e. 2 gels either side of the PAGE insert: made using 1 plain glass plate and 1 notched glass plate, both with bonded spacers, and 1 notched plate without spacers), the thinner YELLOW sliding gel clamps **must** be used.

To convert the shiroGEL VERTICAL 20 from a 2- to 4-gel configuration please adhere to the following instructions.

1. To replace the green sliding gel clamps, begin by unscrewing the color-coded vertical screw pins. There should be no glass plates within the PAGE insert, in the side(s) being unscrewed. Once the screw pins are unscrewed sufficiently the green clamp should sit in the resting slot as shown.

2. Gently push the green sliding clamp horizontally towards the core of PAGE insert until it can move no further.


3. Push the green sliding gel clamp out of the PAGE insert by gently pressing the protruding pin as shown.

4. Once the green sliding gel clamp is removed the PAGE insert is ready to accept the thinner yellow sliding gel clamps for 2-gel sandwiches, either side of the PAGE insert, to convert the shiroGEL VERTICAL 20 to a 4-gel configuration.

5. Insert the yellow sliding gel clamp into the hole closest to the core of the PAGE insert. Once inserted gently withdraw the sliding clamp outwards away from the core of the PAGE insert as shown.

6. Withdraw the sliding gel clamps so that they sit suspended in the resting slots as shown. Repeat steps 1-6 to replace the remaining green sliding gel clamps. Once complete the PAGE insert is ready for 4-gel assembly.

4-gel Assembly for 1.5- and 2-mm-thick Gels


Once converted into the 4-gel configuration using the yellow sliding gel clamps the shiroGEL VERTICAL 20 is ready to accept and run a mum of four gels ranging in thickness from 1mm to 2mm.

1-mm-thick 2-gel sandwiches may be inserted into the PAGE insert from the top as shown in Figures 1-9 in the ***Detailed set up of the shiroGEL VERTICAL 20 insert*** section on pages 13-14. However, for 1.5- and 2-mm-thick 2-gel sandwiches, it may be necessary to load the shiroGEL VERTICAL 20 by positioning the PAGE insert on its side in order to overcome any resistance that may be posed by the gasket. The instructions below illustrate how best to undertake this.

1. With the PAGE insert upright in its normal position, unscrew the screw pins sufficiently so that that the sliding gel clamps sit suspended within the resting slots. Turn the PAGE insert on its side and push the sliding clamps until they lie completely flat within the PAGE insert parallel with the bench surface.

2. Form a 2-gel sandwich on an even bench surface, as described in **General Rules** on page 10 in the Vertical Gel Casting Section. Push the 2-gel sandwich with the notched plate uppermost into the PAGE insert. This notched plate will be closest to the core of the PAGE insert when it is upright. The uppermost notched glass plate may encounter some resistance from the gasket. This may be overcome by pushing the plain and notched glass beneath it into the PAGE insert (**encircled**).

3. Once the resistance is overcome push the uppermost notched glass plate into the PAGE insert. Restore the PAGE to its normal upright position and begin securing the glass plates for vertical electrophoresis, as described in Steps 3-9 of the **Detailed set up of the PAGE insert** on pages 13-14.

Gel Preparation:

1. It is always advisable to work using stock solutions which allow added convenience and save time when it comes to gel pouring. Page 22 lists stock solutions for SDS PAGE gels which should be pre-made beforehand. For native gel formula and running conditions, please consult a laboratory manual. The protocol below is given for use of the standard stock solutions advised. This should be adjusted if you are using different stock solutions or gel formulas.
2. Table 1 below shows the total volume of gel solution required. In subsequent tables, amounts of gel and solutions are given for two 1mm thick gels so adjustments are needed for when running single or more than two gels and for 0.75, 1.5 or 2mm thick spacers.

Table 1.

shiroGEL VERTICAL 20	
	Total Gel volume for a 1mm thick gel.
For different thicknesses of gel, multiply the below amounts by the spacer thickness.	
Single – one gel, one dummy plate	35ml
Double – two gels	70ml
Using a Triple Plate 2-gel sandwich – four gels	140ml

Gel Selection:

Care should be taken when selecting the pore size of the gel to be used.

The pore size or % of gel determines the resolving ability given different sizes of protein.

See Table 2 below which details which percentage of gel to use to separate the sizes of proteins indicated.

Table 2.

Acrylamide Percentage	Separating Resolution
5 %	60 - 220 KD
7.5 %	30 - 120 KD
10 %	20 - 75 KD
12%	17 – 65 KD
15 %	15 -45 KD
17.5%	12 – 30 KD

3. Prepare gel solutions as per tables below. These give the volumes of solutions from the standard stock solutions. These should be gently mixed avoiding generation of bubbles which will inhibit polymerization by removing free radicals.

Table 3: Preparation of the separating gel solution for two 20 x 20cm shiroGEL VERTICAL 20 gels using 1 mm spacers.

Solution	5 %	7.5%	10 %	12%	15 %	17.5%
Distilled Water	41ml	35.25ml	29.6ml	24.7ml	17.6ml	11.7ml
30 % Stock Acrylamide Solution	11.7ml	17.6ml	23.5ml	28.2ml	35.25ml	41.1ml
4 X Resolving Tris Solution	17.6ml	17.6ml	17.6ml	17.6ml	17.6ml	17.6ml
10 % Ammonium Persulphate	700µl	700µl	700µl	700µl	700µl	700µl

Gel Pouring:

For gels with stacking layers:

1. Insert the comb into the glass plates and mark a point on the glass plates 1cm below where the comb teeth finish. This indicates where to add the resolving gel to.
2. Add 70 μ l of TEMED to the gel solution and mix well, but avoid generating air bubbles.
3. Fill the glass plates again avoiding generating any air bubbles. Filling must be performed quickly before the TEMED causes the gel to become too viscous.
4. Overlay the gel extremely carefully with 1 ml of Isobutanol, Isopropanol or distilled water. When using distilled water extra care must be taken to ensure there is no mixing with the gel solution.
- 5 Let the resolving gel polymerize. Usually this takes around 15 minutes but this can vary due to the freshness of the reagents used. If polymerization is taken a lot longer than this, use fresher stock solutions or add more APS and TEMED.
6. Prepare the stacking gel using Table 4 below as a guide. Again stock solutions are given on pages 22 and 23.

Table 4.

Solution	shiroGEL VERTICAL 20
Distilled Water	16.8ml
30 % Stock Acrylamide Solution	2.6ml
4 X Stacking Gel Tris Solution	6.4ml
10 % Ammonium Persulphate	178 μ l

7. Carefully mix the stacking gel solution, avoiding generating air bubbles.
8. Pour off the overlay liquid and rinse the gel with distilled water.
9. For shiroGEL VERTICAL 20 gels add 26.8 μ l of TEMED. Mix well. Use a Pasteur pipette to fill the glass plates up to the top with stacking gel solution.

10. Carefully insert the comb making sure that no air bubbles get trapped under the ends of the comb teeth as these will inhibit sample progression.
11. Allow the stacking gel polymerize for 30 minutes.

For gels without stacking layers:

1. Add 70µl of TEMED for shiroGEL VERTICAL 20 gels and mix well but avoid generating air bubbles.
2. Fill the glass plates again avoiding generating any air bubbles. Filling must be performed quickly before the TEMED causes the gel to become too viscous.
3. Carefully insert the comb making sure that no air bubbles get trapped under the ends of the comb teeth as these will inhibit sample progression.
4. Let the gel polymerize. Usually this takes around 15 minutes but this can vary due to the freshness of the reagents used. If polymerization is taken a lot longer than this, use fresher stock solutions or add more APS and TEMED.

Sample Preparation & Loading

Preparation of denatured protein samples for loading:

The instructions given below are for denatured samples. For Native samples, please consult a laboratory handbook.

1. Prepare the protein samples for loading. The volume of sample depends on the capacity of the wells (See Comb specifications page 25).
2. Using a 0.5 ml micro-centrifuge tube or other convenient receptacle, combine the protein sample and 4 X sample buffer. It is always advisable to use protein markers in one of the end lanes to indicate sizes of bands. These should be prepared according to the manufacturer's instructions.
3. Heat the samples in a water bath or heating block for 2 minutes to denature the samples.
4. Centrifuge the samples in a micro-centrifuge for 20 seconds at 12,000 rpm. The protein samples are now ready to load.

Loading the samples:

1. Transfer the Inner gel module containing cast gels into the main tank in the correct orientation as indicated (+) on the module aligned with (+) on the tank, (-) on the module aligned with (-) on the tank.
2. Fill the outer tank with 1 x reservoir buffer. See page 22 for recommended running buffer solution. Table 6 shows the volume of buffer required.
3. Load the samples into the wells using a pipette tip taking care not to damage the wells or induce any air bubbles.
4. Fill any unused wells with 1 X sample buffer.
5. It is a good idea to note the orientation and order the samples were loaded in. This can be done by noting which samples were loaded adjacent to each electrode.

Buffer Volume:

Table 5.

Buffer Volume	shiroGEL VERTICAL 20
Minimum – Inner tank is filled to above the wells. Outer Tank is filled to just flood the bottom of the glass plates. Cooling potential is at a minimum which may affect resolution.	Total Inner Buffer Chamber:: 640ml Outer Tank: 1L
mum – Inner tank is filled to above the wells. Outer Tank is filled to the mum fill line. Cooling is high offering good resolution of samples. This may be further enhanced by using the cooling coil	Total Inner Buffer Chamber:: 640ml Outer Tank for 2 gels: 5.4L Outer Tank for 4 gels: 4.8L

Gel Running:

1. Fit the lid and connect to a power supply.
2. Consult Table 6, below, for details on recommended power supply voltage settings.
3. Turn the power supply off when the loading dye reaches the bottom of the gel, sooner if your proteins are below 4Kd in size.
4. Remove the gel running module, first emptying the inner buffer into the main tank. Buffer can be re-used but this may affect run quality if continued.
5. Unscrew the glass plates and gently pry apart the glass plates. The gel will usually stick to one of the plates and can be removed by first soaking in buffer and then gently lifting with a spatula.
6. The gel is now ready to be stained with Coomassie or silver stain or the proteins in the gel can be transferred to a membrane by electroblotting for specific band identification and further analysis.

Table 6.

Recommended Voltages and Resultant Current for 1mm thick, 12% gels.	shiroGEL VERTICAL 20
2-4 gels	35mA (constant); up to 350V maximum Runtime to 5h – no cooling; 4h with cooling

Stock Solutions for SDS PAGE gels:

Stock 30% Acrylamide Gel Solution:-

30.0 g acrylamide

0.8 g methylene bisacrylamide

Distilled Water to 100ml

Stock 4 X Resolving Gel Tris (1.5 M Tris HCl pH8.8, 0.4 % SDS)

To 110ml Distilled Water add 36.4 g of Tris base

Add 8ml of 10 % SDS

Adjust pH to 8.8 with 1N HCl

Adjust the final volume to 200ml with Distilled Water.

Stock 4 X Stacking Tris (0.5 M Tris HCL pH6.8, 0.4 % SDS)

To 110ml Distilled Water add 12.12 g of Tris base

Add 8ml of 10 % SDS

Adjust pH to 6.8 with 1N HCl

Add Distilled Water to a final volume of 200ml

Stock 4 X Tris-glycine tank buffer - SDS

36 g Tris base

172.8 g glycine

Distilled Water to 3 L

1 x Tris-glycine tank buffer - SDS

750ml of 4 X Tris-glycine reservoir buffer - SDS

30ml of 10 % SDS

Distilled Water to 3L

10 % AP (ammonium persulphate solution)

0.1 g ammonium persulphate

0.2 1ml Distilled Water

TEMED**Stock 4 X Sample Buffer**

4ml glycerol

2ml 2-mercaptoethanol

1.2 g SDS

5ml 4 X Stacking Tris

0.03 g Bromophenol blue

Aliquot into 1.5ml microcentrifuge tubes. Store at -20°C.

References:

1. Sambrook, Fritsch, and Maniatis, **Molecular Cloning A Laboratory Manual**, Second Edition, Cold Spring Harbor Laboratory Press, 1989.
2. **Current Protocols in Molecular Biology**, Greene Publishing Associates and Wiley-Interscience, 1989.

Combs:

(MC denotes multi-channel pipette compatible)

shiroGEL VERTICAL 20

Code	Description	Sample Volume μ l Per Well
700-0376	Comb 1 Prep, 1 Marker, 0.75mm thick	1100
700-0377	Comb 5 sample, 0.75mm thick	160
700-0378	Comb 10 sample, 0.75mm thick	80
700-0379	Comb 18 sample MC, 0.75mm thick	40
700-0380	Comb 24 sample, 0.75mm thick	30
700-0381	Comb 30 sample, 0.75mm thick	25
700-0382	Comb 36 sample MC, 0.75mm thick	20
700-0383	Comb 48 sample, 0.75mm thick	15
700-0384	Comb 1 Prep, 1 Marker, 1mm thick	1500
700-0385	Comb 5 sample, 1mm thick	200
700-0386	Comb 10 sample, 1mm thick	100
700-0387	Comb 18 sample, 1mm thick	50
700-0388	Comb 24 sample, 1mm thick	40
700-0389	Comb 30 sample, 1mm thick	35
700-0390	Comb 36 sample MC, 1mm thick	25
700-0391	Comb 48 sample, 1mm thick	20
700-0392	Comb 1 Prep, 1 Marker, 1.5mm thick	2200
700-0393	Comb 5 sample, 1.5mm thick	320
700-0394	Comb 10 sample, 1.5mm thick	160
700-0395	Comb 18 sample, 1.5mm thick	80
700-0396	Comb 24 sample, 1.5mm thick	60
700-0397	Comb 30 sample, 1.5mm thick	50
700-0398	Comb 36 sample MC, 1.5mm thick	40
700-0399	Comb 48 sample, 1.5mm thick	30
700-0400	Comb 1 Prep, 1 Marker, 2mm thick	3000
700-0401	Comb 5 sample, 2mm thick	400
700-0402	Comb 10 sample, 2mm thick	200
700-0403	Comb 18 sample, 2mm thick	100
700-0404	Comb 24 sample, 2mm thick	80
700-0405	Comb 30 sample, 2mm thick	70
700-0406	Comb 36 sample MC, 2mm thick	50
700-0407	Comb 48 sample, 2mm thick	40

Connecting to a Chiller Unit

1. Vertical PAGE may sometimes result in high currents. High current may cause the build-up of heat which can affect sample migration and gel resolution as a result. To counteract the adverse effects of heat generation, the manufacturer recommends using an effective chiller unit, preset to 4°C.
2. For active temperature regulation attach two short lengths of hose from the inlet and outlet ports of the chiller unit to the respective outlet and inlet connectors of the cooling coil sitting within the shiroGEL VERTICAL 20 outer tank.
3. Ensure that the flow rate of the chiller unit is set at a maximum 1L/minute. Do not exceed this limit as damage to the cooling platform might occur.
4. If using a chiller unit which exceeds this limit, a T-connector may be used. The T-connector can divert some coolant volume back to the recirculating chiller reducing the pressure incident within the coil.

The unit is now ready to be used with active cooling.

Technical service

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

EQUIPMENT DISPOSAL

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle-end by handing it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland

VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
E-mail: info@pl.vwr.com

Portugal

VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43-3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey

Pro Lab Laboratuvar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İşmerkezi No.32/1
34896 Pendik - İstanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

Australia

VWR International, Pty Ltd.
Unit 1/31 Archimedes Place
Murarrie, Queensland 4172
Tel.: 1300 727 696
Fax: 1300 135 123

China

VWR International China Co., Ltd
Rm.219, 2100 Dongming Road
Pudong New District
Shanghai 200123
Tel.: +86-21-5898 6888
Fax: +86-21-5855 8801
E-mail: info_china@vwr.com

India

VWR Lab Products Private Limited
135/12, Brigade Towers, 2nd Floor
Front wing, Brigade Road,
Bengaluru, India – 560 025
Tel.: +91-80-41117125/26 (Bengaluru)
Tel.: +91-2522-647911/922 (Mumbai)
Fax: +91-80-41117120
E-mail: vwr_india@vwr.com

New Zealand

Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

**GO TO VWR.COM FOR THE LATEST
NEWS, SPECIAL OFFERS AND
DETAILS OF YOUR LOCAL VWR
DISTRIBUTOR**

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