Title: Technical Data Sheet Rovers® Cervex-Brush® (Sterile)
(Reference: 380100331)

Date: September 16th, 2021
Doc. ID: 2021259LEA
Version: 1.3

DUTCH DESIGN
MADE IN THE NETHERLANDS

THE INFORMATION CONTAINED HEREIN IS CONFIDENTIAL AND THE PROPERTY OF ROVERS MEDICAL DEVICES B.V., LEKSTRAAT 10, 5347 KV, OSS, THE NETHERLANDS. THE RECIPIENT IS RESPONSIBLE FOR SAFE-KEEPING, AND PREVENTION OF UNAUTHORIZED APPROPRIATION, USE, DISCLOSURE, AND COPYING.

1 Introduction

Rovers® Cell Sampling Devices (Sterile) for professional and home use are based on "thin-hair technology". Cell material collected with Rovers® Cell Sampling Devices can be processed with conventional as well as liquid based cytology(LBC) and also be used for bacteriological testing, virus testing (e.g. for Sexually Transmitted Diseases or HPV) and DNA analysis.

Cell collection with Rovers® Cell Sampling Devices is quick, simple, painless and reliable.

In accordance with the European Medical Device Directive 93/42/EEC (as amended by Directive 2007/47/EC) the Rovers® Cell Sampling Devices (Sterile) are Class **Is** Medical Devices with a low risk profile.

Author:

L. Evers,
QA/RA Manager

Signature:

A handwritten signature in blue ink, appearing to be "L. Evers", written over a light blue horizontal line.

Approved:

H. Vissers,
Senior QA/RA Manager

Signature:

A handwritten signature in blue ink, appearing to be "H. Vissers", written over a light blue horizontal line.

Table of Contents

1	INTRODUCTION	2
	TABLE OF CONTENTS	3
2	ABBREVIATIONS	4
3	BACKGROUND INFORMATION	4
3.1	MANUFACTURER PROFILE / IDENTIFICATION	4
3.2	COMPANY ADDRESS	5
3.1	WEBSITES	5
4	ROVERS® CERVEX-BRUSH®	6
4.1	MARKET PRESENTATION	6
4.2	EU DEVICE CLASSIFICATION	6
4.3	QUALITY CERTIFICATES	7
4.4	GMDN	7
4.5	UMDNS	7
4.6	THE ROVERS® CERVEX-BRUSH® IN RESEARCH STUDIES	7
4.7	DOWNLOAD LINK PRODUCT BROCHURE	7
4.8	INTENDED USE	7
4.9	INTENDED USER	7
4.10	INTENDED PLACE OF USE / USE ENVIRONMENT	8
4.11	PRECAUTIONS (CONTRA-INDICATIONS)	8
4.12	RESTRICTIONS	8
4.13	WARNINGS	8
4.14	SPECIAL CARE, STORAGE AND TRANSPORT CONDITIONS	8
4.15	MATERIAL INFORMATION ROVERS® CERVEX-BRUSH®	8
4.16	PRIMARY PACKAGING INFORMATION	9
4.17	SECONDARY PACKAGING PARTICULARS	10
4.18	TERTIARY PACKAGING PARTICULARS	11
5	INSTRUCTIONS FOR USE	12

2 Abbreviations


AAMI	Association for the Advancement of Medical Instrumentation
EEC	European Economic Community
EN	European Norm
ETO	Ethylene Oxide
EU	European Union
HPV	Human Papilloma Virus
hrHPV	high-risk Human Papilloma Virus
LBC	liquid based cytology
RMD	Rovers Medical Devices B.V.

3 Background Information

3.1 Manufacturer Profile / Identification

Rovers Medical Devices B.V. (RMD) designs, develops and manufactures and worldwide markets medical devices (cancer screening devices) for cytological analysis, bacteriological-, HPV-, viral- and DNA testing.

The company was founded in 1982.

Our main product line is Rovers® Cell Sampling Devices. The "thin-hair technology" used for this product has been developed by RMD.

Cell material collected with our sampling devices can be processed with conventional as well as liquid based cytology(LBC) and also be used for bacteriological testing, virus testing (e.g. for Sexually Transmitted Diseases or HPV) and DNA analysis.

Our Rovers® products are being used in more than 60 countries all over the world. For the development of new products we work together with medical specialists in gynaecology, cytology, pathology etc.

3.2 Company Address

Rovers Medical Devices B.V. (RMD) Lekstraat 10 5347 KV, Oss The Netherlands	 A photograph of the Rovers Medical Devices building, a modern three-story structure with a grey facade and blue accents around the windows and entrance. A large tree stands in the foreground, and a black metal fence runs along the front of the property. The sky is blue with some clouds.
--	--

3.1 Websites

<http://www.roversmedicaldevices.com/>

3.2 Contact person

H. Vissers - Quality Assurance / Regulatory Affairs Manager
T: +31 412 648870
F: +31 412 623835
Email: hvissters@roversmedicaldevices.com

4 Rovers® Cervex-Brush®

“BECAUSE THE HAIRS ARE FLEXIBLE, THE MUCOSAL SURFACE IS WELL SAMPLED BUT NOT DAMAGED, THUS REDUCING THE NUMBER OF UNSATISFACTORY SMEARS.”

Prof. G.P. Vooy's

The Cervex-Brush® enables simultaneous collection of ectocervical, endocervical and transformation-zone cells with a single device. The brush is made up of well-defined semi-circular soft, flexible hairs for the optimal sample collection. The Cervex-Brush® can be used for HPV testing, conventional cytology and liquid-based cytology. The hydrophobic material of the brush facilitates the release of the cell material into the fluid or on to a glass slide.

Benefits

- Sweeps the portio and reaches well beyond the transformation-zone into the cervical canal for a rich cell yield.
- Significantly reduces the number of inadequate smears.¹
- Avoids bleeding and pain which are normally associated with rigid spatulas and less considerate endocervical devices.
- The Rovers® Cervex-Brush® can be used for liquid-based as well as conventional cytology. The hydrophobic material of the brush facilitates the release of the cell material into the fluid.
- Optional detachable brush head can be supplied depending on the LBC method used.

4.1 Market Presentation

Product Name:	Article Number
Rovers® Cervex-Brush® (STERILE)	380100331

Table 2 : Article # Rovers® Cervex-Brush®

- 100 PCS / Box + 1 IFU

Product Name:	UDI
UDI-DI	08719189246353
Primary packaging	08719189246353
Secondary packaging	08719189246056
Tertiary packaging	18719189246053

4.2 EU Device Classification

In accordance with the European Medical Device Directive 93/42/EEC (as amended by Directive 2007/47/EC) the Rovers® Cervex-Brush® (STERILE) is a Class **Is** Medical Device.


EU market access is granted by:

Notified Body #: 0044, TÜV NORD CERT GmbH, Langemarckstr. 20, 45141 Essen, Germany.

4.3 Quality Certificates

Copies of quality certificates can be downloaded via:
<https://www.roversmedicaldevices.com/about-rovers/quality-certificates/>

4.4 GMDN

Term: Cytology scraper, single-use

GMDN Code: **42537**

Definition: A wooden, plastic, or metal instrument used to collect a scraping of cells from tissue for the purposes of a diagnostic examination. This is a single-use device.

4.5 UMDNS

UMDNS Code:

15-018, Cell Sampling Devices

4.6 The Rovers® Cervex-Brush® in research studies

<https://www.roversmedicaldevices.com/cell-sampling-devices/cervex-brush/>

4.7 Download Link Product Brochure

https://www.roversmedicaldevices.com/wp-content/uploads/2016/07/ROVER0140_FOLDER-CERVEX-BRUSH.pdf

4.8 Intended Use

The Rovers® Cervex-Brush® is a single-use cell sampler intended to collect simultaneous cell material from the cervical os, ectocervical and endo cervical in women.

4.9 Intended User

They are for use by healthcare professionals (qualified medical personnel).

4.10 Intended Place of Use / Use Environment

Clinical setting.

4.11 Precautions (Contra-indications)

The Rovers® Cervex-Brush® should not be used on patients after the first 10 weeks of pregnancy.

4.12 Restrictions

To be used only by qualified medical personnel.

4.13 Warnings

Do not re-use.

Do not use if package is damaged.

Caution, consult instructions for use.

Intended for single use only.

Reuse causes contamination and misdiagnosis.

To be used only by qualified medical personnel.

4.14 Special Care, storage and transport conditions

The primary packaged products are all protected for transport and storage by a secondary and tertiary carton box. The expiry period is three years after date of production, provided that the products are stored under ambient storage conditions and protected from sunlight.

The white cardboard boxes (Secondary Packaging) are packaged according a predefined stacking pattern in a brown cardboard box (Tertiary Packaging, also called the shipping container) on a wooden pallet. Hereafter the shipping container is secured with ratchet straps and foil wrapped and ready for transport. For smaller amounts the quantity of packaging and protection materials will be adjusted accordingly.

4.15 Material Information Rovers® Cervex-Brush®


The blue coloured handle is made of Polypropylene, the white coloured brush is made of Polyethylene. After assembling/manufacturing the product is packaged in plastic pouch and sealed with a breathable medical paper. Sterilization is performed with a validated EO process.

Part description	Material
Brush	Polyethylene / White Colorant
Handle	Polypropylene / Blue Colorant
Pouch	Polyamide / Polyethylene or Polypropylene / Polyamide / Polyethylene
Paper lidding	Medical Paper

Product Drawing:

The Rovers® Cervex-Brush® product consists of a blue handle (length 175 mm, diameter approx. 4 mm(B)) and a white brush (l x w x h = 32 mm x 19 mm(C) x 6 mm) with hairs. The total length is 198(A) mm.

Total Length Rovers® Cervex Brush®:

By design the total length of the Rovers® Cervex Brush® is determined by the individual dimensions of the handle and the brush. During manufacturing the appearance of the handle, brush and assembled product is monitored.

Product Weight:

– 2 g

(The provided product weight is indicative and has not been determined as an in-process and/or release control parameter).

4.16 Primary Packaging Information

Dimensions	26,5×4.2 cm
------------	-------------

Information as available on the primary packaging:

Rovers® Cervex-Brush®, Cervical Cell Sampler

 (1)	 (2)	 (3)	 Liquid-Based Pap Preparation (A)	 Traditional Pap Preparation (B)	 Do not reuse. Do not use if package is damaged.	 CE 0044	 LOT
					 30 cm	 Attention, see instructions for use.	 MD

Rovers Medical Devices B.V., Lekstraat 10, 5347 KV, Oss, The Netherlands REF 380100331

4.17 Secondary Packaging Particulars

	White coloured carton box
	100 PCS of product + 1 Instructions For Use
	29,5×16,5×10,0 cm


Information as available on the label:

Rovers® Cervex-Brush® 100 pcs

LOT 41xxx **REF** 380100331

 **2021-03-01** **C €0044**

 **2024-03-01** **STERILE EO**

www.roversmedicaldevices.com    (01)08719189246056
Rovers Medical Devices BV, Lekstraat 10, 5347 KV Oss, The Netherlands (17)240301 (10)41XXX **MD**

Rovers® Cervex-Brush® is a patented product. Registered design.

Rovers® and Cervex-Brush® are registered trademarks of Rovers Medical Devices B.V. V2103

4.18 Tertiary Packaging Particulars

Tertiary Packaging	Brown coloured carton box
	5 layers of 4 box's (Secondary Packaging)
	71.5×30.5×54.5 cm

For smaller amounts the quantity of packaging and protection materials will be adjusted accordingly.

Rovers® Cervex-Brush®

REF 380100331

LOT 41xxx

 2021-03-01


 2024-03-01

MD

Rovers®
medical devices

(01)08719189246216
(17)240301
(10)41XXX

20 x 100 pcs

STERILE EO

C €0044

www.roversmedicaldevices.com
Rovers Medical Devices BV, Lekstraat 10, 5347 KV Oss, The Netherlands
Rovers® Cervex-Brush® is a patented product. Registered design. V2103
Rovers® and Cervex-Brush® are registered trademarks of Rovers Medical Devices B.V.

11/13

5 Instructions For Use

Instructions For Use information as supplied with every Secondary Packaging:

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