

DIGE Gel and DIGE Buffer Kit

2-D Fluorescence Difference Gel Electrophoresis (2-D DIGE) is a key technology for analyzing differences in protein abundance in biological samples. DIGE Gel and DIGE Buffer Kit are tools to increase the quality, reproducibility, and convenience of 2-D DIGE using Ettan™ DIGE System. DIGE Gel is a precast polyacrylamide gel in a low fluorescent glass cassette for second dimension 2-D electrophoresis. The DIGE Buffer Kit is a specially designed two-buffer system for running up to twelve DIGE Gels simultaneously (Fig 1).

DIGE Gel and DIGE Buffer Kit offer:

- *High quality protein spot maps:* DIGE Gel provides high sensitivity and high resolution protein spot maps of equal quality compared to homecast polyacrylamide gels.
- *Reproducibility:* DIGE Gel run with DIGE Buffer Kit result in spot maps with higher reproducibility compared to homecast polyacrylamide gels, reducing technical variation between gels.
- *Convenience:* DIGE Gel and DIGE Buffer Kit are ready-to-use products with a long shelf-life that saves you time by eliminating the need to make buffers, cast gels and wash glass plates, leaving more time for more important work.

Description

DIGE Gel is a 12.5% precast, low fluorescent polyacrylamide gel cast in a low fluorescent glass cassette specially developed for 2-D DIGE analysis. DIGE Gel should be used with the DIGE Buffer Kit, which consists of concentrated running buffers and Sealing Solution for attaching Immobiline™ DryStrip Gels (IPG strips) to the top of the polyacrylamide gel. The capacity of the buffer system is similar to the commonly used Laemmli (Tris-Glycine) buffer system, and the separation performance of DIGE Gel is comparable to 12.5% Laemmli gels (1). Therefore, DIGE Gel and DIGE Buffer Kit offer the high quality and performance of the Laemmli system with the added advantages of higher reproducibility, convenience, and time savings. The main characteristics of DIGE Gel and DIGE Buffer Kit are listed in Table 1.



Fig 1. DIGE Gel and DIGE Buffer Kit for differential protein analysis using the Ettan DIGE System.

Table 1. Characteristics of DIGE Gel and DIGE Buffer Kit

DIGE Gel	
Gel composition	T = 12.5%, C = 3% (12.125% acrylamide, 0.375% bisacrylamide)
Separation range	12 000 to 120 000 M _r , relative molecular mass
Gel dimensions	255 × 196 × 1 mm
Buffer in gel	Buffer based on piperidinopropionamide (PPA)
Gel cassette	Low fluorescent glass
Shelf life	12 months
Storage	4°C to 8°C
DIGE Buffer Kit	
Anode Buffer (2 bottles)	Buffer based on PPA
Cathode Buffer (4 bottles)	0.25 M Tris, 1.92 M glycine, 1% (w/v) SDS
Sealing Solution	Gel Buffer with 0.5% agarose and 0.002% bromophenol blue
Shelf life	Estimated 12 months
Storage	4°C to 8°C



DIGE Gel is intended for use with the DIGE Buffer Kit and Ettan DALT Large Vertical System. Ettan DALT_{twelve} and Ettan DALT_{six} instruments are designed for large-format, second dimension 2-D electrophoresis using either 18 or 24 cm Immobiline DryStrip Gels. With Ettan DALT_{twelve}, up to 12 DIGE Gels can be run simultaneously with separations complete in 4–6 h. The DIGE Buffer Kit contains the amount of reagents necessary for a single run of up to 12 DIGE Gels using the Ettan DALT_{twelve} electrophoresis unit, or two runs of up to six gels using the Ettan DALT_{six} unit.

High quality protein spot maps

Homecast polyacrylamide gels give high quality separation and resolution, but the manual step of casting second dimension gels introduces gel-to-gel variation. Figure 2 shows that DIGE Gels provide high sensitivity and high resolution protein spot maps of equal quality compared to homecast polyacrylamide gels.

It is important to note that protein samples run on homecast polyacrylamide gels, DIGE Gel or DALT Gel 12.5 (12.5%) will give rise to differences in spot maps due to varying gel composition and buffer system. DIGE Gel and DALT Gel 12.5 result in very similar protein spot maps, and the DIGE Buffer Kit is compatible with both types of precast gels. In Figure 3, DIGE Gels were run with conventional Laemmli (Tris-Glycine) buffer which can result in protein spot maps distorted in the low molecular weight region (Fig 3A). Figure 3B shows that using the DIGE Buffer Kit, together with DIGE Gels, results in superior protein spot maps with focused and clearly resolved protein spots, even in the low molecular weight region.

Sample labeling: *E.coli* protein sample (50 µg) labeled with CyDye DIGE Fluor, Cy5 minimal dye

Sample buffer: 7 M urea, 2 M thiourea, 2% CHAPS, 130 mM DTT, and 2% IPG Buffer pH 4–7.

First dimension IEF: Immobiline DryStrip pH 4–7, 24 cm rehydrated overnight at room temperature in DeStreak Rehydration Solution + 2% IPG Buffer pH 4–7. Separation using Ettan IPGphor™ 3 IEF unit, anodic cup loading, 60–62 kVh.




Fig 2. Second dimension separation using Ettan DALT_{twelve} with (A) 12.5% homecast polyacrylamide gels or (B) 12.5% precast DIGE Gels. DIGE Gels provide high resolution protein spot maps of equal quality compared to homecast polyacrylamide gels.

Sample labeling: *E.coli* protein sample (50 µg) labeled with CyDye DIGE Fluor, Cy5 minimal dye

Sample buffer: 7 M urea, 2 M thiourea, 2% CHAPS, 130 mM DTT, and 2% IPG Buffer pH 4–7.

First dimension IEF: Immobiline DryStrip pH 4–7, 24 cm rehydrated overnight at room temperature in DeStreak Rehydration Solution + 2% IPG Buffer pH 4–7. Separation using Ettan IPGphor™ 3 IEF unit, anodic cup loading, 60–62 kVh.




Fig 3. Second dimension separation using 12.5% precast DIGE Gels run on Ettan DALT_{twelve} with (A) Tris-Glycine buffer or (B) DIGE Buffer Kit. With conventional Tris-Glycine buffer, protein spot maps can be distorted in the low molecular weight region (red circle). Using DIGE Buffer Kit, together with DIGE Gels, gives clearly resolved spot maps even in the low molecular weight region.

Technical replicates were run using either DIGE Gels or homecast polyacrylamide gels. Spot matching was performed, and the number of matched protein spots in the replicate analyses were plotted. Figure 4 shows that DIGE Gel gives spot maps with a higher reproducibility compared to homecast gels. The figure shows that more protein spots were consistently matched from the technical replicates using DIGE Gels with DIGE Buffer Kit.




Fig 4. Higher reproducibility with DIGE Gel: DIGE Gels run with DIGE Buffer Kit or homecast polyacrylamide gels run with Laemmli buffer. Eight technical replicates of DIGE Gels, and six technical replicates of homecast polyacrylamide gels. Log abundance for all matched spots exported from DeCyder Biological Variation Analysis (BVA) module, recalculated to covariance (%) and plotted.

References

1. Laemmli, U.K. Cleavage of structural proteins during the assembly of the head of bacteriophage T4. *Nature* **227**, 680–685 (1970).

Synchronized solutions for better 2-D electrophoresis - improved 2-D DIGE

2-D electrophoresis is a powerful technology for comparing complex protein mixtures from biological samples in proteomics research. With the aim of simplifying the process and improving quantitation, GE Healthcare has made significant improvements across the entire workflow, from sample prep and separation to data analysis. 2-D DIGE

(Difference Gel Electrophoresis) allows you to analyze differences in low abundance proteins with unparalleled statistical confidence. Many of our 2-D electrophoresis products improve data quality from traditional 2-D experiments as well, and the products can eventually be integrated into a complete solution to maximally improve your 2-D results.



Sample Preparation	Protein Labeling	Protein Separation	Image Acquisition	Image Analysis	Validation
2-D Protein Extraction Buffer	DIGE Trial Pack CyDye™ DIGE Fluor 2 nmol kit CyDye DIGE 400 nmol sat dye CyDye DIGE repackaging	IPGbox Immobiline™ DryStrip gels IPG Buffer DIGE Gel and DIGE Buffer Kit	Etan™ DIGE Imager Typhoon™ Variable Mode Imager	DeCyder™ 2-D software v7.0	ECL Plex™ ImageQuant™ TL 7.0 with ImageQuant TL SecurITy

For more information, visit www.gelifesciences.com/2DE

Ordering information

Product	Quantity	Code No	Related Products	Quantity	Code No
DIGE Gel	3	28-9374-51	<i>First dimension electrophoresis</i>		
DIGE Buffer Kit	2 × 125 ml anode buffer 4 × 125 ml cathode buffer 12 tubes agarose sealing solution	28-9374-52	Etan IPGphor3 Isoelectric Focusing Unit	1	11-0033-64
			IPGphor Manifold ceramic tray, complete	1	80-6498-38
			IPGbox	1 IPGbox + IPGbox Kit	28-9334-65
			IPGbox Kit	10 Reswell Trays + IPGbox Insert	28-9334-92
			IPG Buffer pH 3.5-5.0	1 ml	17-6002-02
			IPG Buffer pH 5.5-6.7	1 ml	17-6002-06
			IPG Buffer pH 4-7	1 ml	17-6000-86
			IPG Buffer pH 6-11	1 ml	17-6001-78
			IPG Buffer pH 7-11 NL	1 ml	17-6004-39
			IPG Buffer pH 3-10 NL	1 ml	17-6000-88
			IPG Buffer pH 3-10	1 ml	17-6000-87
			IPG Buffer pH 3-11 NL	1 ml	17-6004-40
			DeStreak Rehydration Solution	5 × 3 ml	17-6003-19
			DeStreak Reagent	1 ml	17-6003-18
Related Products	Quantity	Code No			
<i>Sample labeling</i>					
DIGE Trial Pack, incl. CyDye DIGE Fluor minimal dye labeling kit (2 nmol) + 1 month free DeCyder 2-D trial license	1	28-9373-73			
CyDye DIGE Fluor minimal dye labeling kit	2 nmol	28-9345-30			
CyDye DIGE Fluor minimal dye labeling kit	5 nmol	25-8010-65			
CyDye DIGE Fluor Cy2 minimal dye	5 nmol	25-8010-82			
CyDye DIGE Fluor Cy3 minimal dye	5 nmol	25-8010-83			
CyDye DIGE Fluor Cy5 minimal dye	5 nmol	25-8010-85			
CyDye DIGE Fluor Cy2 minimal dye	10 nmol (2 × 5)	25-8008-60			
CyDye DIGE Fluor Cy3 minimal dye	10 nmol (2 × 5)	25-8008-61			
CyDye DIGE Fluor Cy5 minimal dye	10 nmol (2 × 5)	25-8008-62			
CyDye DIGE Fluor Cy2 minimal dye	25 nmol (5 × 5)	RPK0272			
CyDye DIGE Fluor Cy2 minimal dye	25 nmol (5 × 5)	RPK0273			
CyDye DIGE Fluor Cy2 minimal dye	25 nmol (5 × 5)	RPK0275			
CyDye DIGE Fluor Labeling Kit for Scarce Samples	1	25-8009-83			
CyDye DIGE Fluor Labeling Kit for Scarce Samples and Preparative Gel Labeling	1	25-8009-84			
CyDye DIGE Fluor Preparative Gel Labeling for Scarce Samples	400 nmol	28-9366-83			

Immobiline DryStrip Gels for IEF

pH range	Quantity	Code No	
		18 cm	24 cm
3.5-4.5	12/pack	-	17-6002-38
3-7 NL	12/pack	-	17-6002-43
4-7	12/pack	17-1233-01	17-6002-46
6-9	12/pack	17-6001-88	17-6002-47
6-11	12/pack	17-6001-97	-
3-10	12/pack	17-1234-01	17-6002-44
3-10 NL	12/pack	17-1235-01	17-6002-45
3-5.6 NL	12/pack	17-6003-56	17-6003-57
5.3-6.5	12/pack	17-6003-61	17-6003-62
6.2-7.5	12/pack	17-6003-66	17-6003-67
7-11 NL	12/pack	17-6003-71	17-6003-72
3-11 NL	12/pack	17-6003-76	17-6003-77

Related Products	Quantity	Code No
Second dimension electrophoresis		
Ettan DALT ^{twelve} Separation Unit & Power Supply/Control Unit, 230 V	1	80-6466-27
Ettan DALT ^{twelve} Separation Unit & Power Supply/Control Unit, 115 V	1	80-6466-46
Ettan DALT ^{six} Electrophoresis Unit, 220 V	1	80-6485-27
Ettan DALT ^{six} Electrophoresis Unit, 115 V	1	80-6485-08

Related Products	Quantity	Code No
Image capture and analysis		
Typhoon Trio™ with ImageQuant™ TL and PC	1	63-0055-88
Typhoon Trio+ with ImageQuant TL and PC	1	63-0055-90
Typhoon 9400 with ImageQuant TL and PC	1	63-0055-79
Typhoon 9410 with ImageQuant TL and PC	1	63-0055-81
Ettan DIGE Imager	1	63-0056-42
DeCyder 2-D v7.0 Software Package	1	28-9435-83
DeCyder 2-D v7.0 Oracle 10gRZ, 5 user license	1	28-9435-88
DeCyder 2-D v7.0, 1 user license	1	28-9442-75

Related literature	Code No
2-D Electrophoresis: Principles and Methods, Handbook	80-6429-60
Ettan DIGE System User Manual	18-1173-17

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Cy, CyDye, ECL Plex, Ettan, DeCyder, ImageQuant, Immobiline, IPGphor, Typhoon, and Typhoon Trio are trademarks of GE Healthcare companies.

2-D Fluorescence Difference Gel Electrophoresis (2-D DIGE) technology is covered by US patent numbers 6,043,025, 6,127,134 and 6,426,190 and equivalent patents and patent applications in other countries and exclusively licensed from Carnegie Mellon University. CyDye: This product or portions thereof is manufactured under an exclusive license from Carnegie Mellon University under US patent numbers 5,569,587, 5,627,027 and equivalent patents in other countries.

The cyanine dyes in the product are manufactured under an exclusive license from Carnegie Mellon University under US patent numbers 5,268,486 and equivalent patents in the US and other countries.

The purchase of CyDye products includes a limited license to use the CyDye products for internal research and development but not for any commercial purposes. A license to use the CyDye products for commercial purposes is subject to a separate license agreement with GE Healthcare. Commercial use shall include:

1. Sale, lease, license or other transfer of the material or any material derived or produced from it.
2. Sale, lease, license or other grant of rights to use this material or any material derived or produced from it.
3. Use of this material to perform services for a fee for third parties, including contract research and drug screening.

If you require a commercial license to use this material and do not have one, return this material unopened to GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden and any money paid for the material will be refunded.

DIGE Gel and DIGE Buffer Kit: The buffer system in this gel and buffer kit is covered by patent application WO9616724 granted in US, EP and JP.

DeCyder: This release of DeCyder (software) is provided by GE Healthcare to the customer under a non-exclusive license and is subject to terms and conditions set out in the 2-D Differential Gel Electrophoresis Technology Access Agreement. Customer has no rights to copy or duplicate or amend the Software without the prior written approval of GE Healthcare.

All third party trademarks are the property of their respective owners.

© 2009 General Electric Company—All rights reserved
First published Mar. 2009

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Limited Amersham Place
Little Chalfont
Buckinghamshire, HP7 9NA
UK

GE Healthcare Europe, GmbH
Munzinger Strasse 5
D-79111 Freiburg
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327
Piscataway, NJ 08855-1327
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg., 3-25-1, Hyakunincho
Shinjuku-ku, Tokyo 169-0073
Japan

For local office contact information, visit
www.gelifesciences.com/contact

www.gelifesciences.com/2DE

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden



imagination at work