Laerdal ALS Baby Trainer

Directions for Use

Page 2 - 7

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Seite 8 - 13

Mode d'emploi

Page 14 - 19

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
Páginas 20 - 25

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Pagine 26 - 31

Handleiding

Blz. 32 - 37

Laerdal
helping people live

CONTENTS

Cautions and Warnings	2
Limited Warranty	2
Introduction	3
Instructions for Use - Manikin Features	4
Maintenance	6
Specifications	38

Cautions and warnings*Latex*

The manikin face, airway as well as the leg replacement pads contain latex. *Users who suffer from latex allergy should take precautions while using or handling the latex parts by wearing non-latex protective gloves.*

ECG connectors

The ECG connectors on the chest skin are intended for monitoring only! The Laerdal ALS Baby Trainer must not be defibrillated. Doing so may be potentially hazardous to the operator and may also damage the ECG simulator.

Alcohol

Do not use alcohol to clean any of the latex parts (face, airway or leg replacement pads). Also do not use Laerdal Airway Lubricant to lubricate the airways, as it contains ethanol. Instead, use the lubricant jelly supplied with the ALS Baby.

Staining

Do not allow the manikin's skin to come in direct contact with ink or photocopied paper, as this can cause permanent stains on the skin. Using protective gloves when handling the manikin will also reduce the risk of staining the face skin and the intraosseous replacement pads. However, avoid coloured vinyl or latex gloves as these may cause discolouration of the manikin's skin.

Upper airways

Always make sure that the airway and tube have been properly lubricated before practising intubation on the manikin. This is important both to provide the proper realism and to reduce the risk of damaging the manikin's upper airways.

Limited Warranty

Laerdal Medical warrants to the purchaser that its product(s) is/are free from defects in material and workmanship for a period of one (1) year from the date of purchase by the original user. During the designated one (1) year period, Laerdal Medical will, upon receipt of a product found to be defective due to materials or workmanship from the purchaser and notification in writing of the defect at its option repair or replace any parts found to be defective or the entire product.

Products found to be defective and notification of defects may also be sent to the authorized Laerdal Medical dealer from whom the product was purchased. All postage, shipping, or handling charges shall be the sole responsibility of the purchaser.

Laerdal Medical is responsible for the effects of safety, reliability and performance of its product(s) only if:

- service, repair, readjustment or modification is carried out by Laerdal Medical or persons authorized by Laerdal Medical.
- the product is used in the proper manner in strict compliance with its Directions for Use.

Laerdal Medical shall not be liable under this warranty for incidental or consequential damages if unauthorized repairs or modifications have been made or attempted or when the product, or any part thereof, has been damaged by accident, misuse or abuse.

This warranty does not cover batteries, bulbs, fuses, normal wear and tear, staining, discoloration or other cosmetic irregularity which does not impede or degrade product performance.

Some states in the USA do not allow the exclusion or limitation of incidental or consequential damages, so those limitations or exclusions may not apply to you.


There are no other express or implied warranties, whether of merchantability, fitness or purpose, or otherwise, on the product, its parts and accessories.

Service

All manikin parts are replaceable. Should any part break or wear out, it can be identified by referring to the exploded view drawing and parts list. Spare parts are available from Laerdal or the local Laerdal distributor.

Since the Laerdal ALS Baby Trainer is easy to disassemble and reassemble, many repairs can be accomplished by the manikin's owner. Should a problem arise, however, we will be glad to assist you.

Introduction

The Laerdal ALS Baby Trainer comes with ECG monitoring chest skin.

It is also possible to upgrade the Laerdal Resusci® Baby to a Laerdal ALS Baby Trainer by means of an Update Kit. The Laerdal ALS Baby Trainer simulates a 3 month old, 5 kilogram / 11 pound infant.

It is designed to provide the following common diagnosis and treatment modalities:

Airway management

- Ventilation via bag-valve-mask
- Endotracheal and nasotracheal intubation
- Auscultation of breath sounds
- Bilateral chest movement and stomach distention
- Oral/Nasal Airways
- Oral/Nasal Gastric Tube Insertion

Vascular access practice


- scalp vein simulation
- realistic intraosseous access and drug therapies.

ECG monitoring and arrhythmia recognition

Assembly - Getting started

Before using your ALS Baby, perform the following preparations:


Lubrication

During this procedure, wear protective gloves to avoid contamination. Take approximately 1 cm of lubricant jelly (supplied with manikin) on your fingertip, an applicator or a cotton swab and apply gently to the manikin's tongue, the back wall of the mouth, and around the tracheal opening. Also lubricate the tip of the endotracheal tube and laryngoscope blade with lubricant jelly before performing endotracheal or nasotracheal intubation.


Placement of IO pads

- Unpack one of the intraosseous access (IO) pads.
- Hold the pad with the blood reservoir pointing upwards.

- Insert the blood reservoir end of the pad into the round cavity in the knee.

- Carefully bend the foot downward, and while holding it down: Slip the tab at the lower end of the IO pad into the square cavity inside the foot.

- Let the foot bend back to normal position.


ECG connection

The three connectors on the lead from the manikin are snapped into place on the corresponding mount on the rhythm simulator. On the Laerdal Heartsim® simulators as well as on the Laerdal Cardiac Rhythm Simulator (model previous to the Heartsim 2000), these connectors are located on the back of the unit.

Instructions for Use - Manikin Features


Ventilation

Airway open/closed - Stomach distension

In the flexed position, with the head resting forward onto the chest, the airways are closed. If an attempt is made to ventilate the manikin via bag-mask between this position and the neutral position, air will not enter the lungs. If an air pressure of more than 10 mm Hg is applied, stomach distention can be observed on the manikin.

Head tilt

In the neutral position, the airways are open. If ventilated via bag-mask in this position, air will enter the lungs. If a too rapid ventilation is performed, causing an air pressure in excess of 10 mm Hg, stomach distention will occur.


Hyperextension of the head/neck will not cause the airways to close.

Jaw thrust

The manikin also allows jaw thrust to be performed. Performed in the flexed position, the airways will open and allow ventilation by bag-mask.

Intubation


The manikin can be intubated via the oral and nasal routes.

A 3,5 mm tube is recommended to avoid backward air leakage. A No. 1 laryngoscope blade is recommended.

If intubated too deeply, the tube will pass into the right mainstem bronchus, causing ventilation to the right lung only.

Important: Always make sure that the airway and tube have been properly lubricated before practising intubation on the manikin. This is important both to provide the proper realism and to reduce the risk of damaging the manikin's upper airways. Use the lubricant jelly supplied with the Laerdal ALS Baby Trainer. Do not use the Laerdal Airway Lubricant.

The cricoid pressure technique (Sellick manoeuvre) can be realistically performed on the ALS Baby Trainer.

Lungs, stomach

The manikin contains two separate lungs and a stomach. They are designed to provide realistic movement and sounds during ventilation of either lung or distention of the stomach. To expel air from the stomach, simply press downward on the manikin's stomach.

The manikin will expel all air back through the mouth and nose.

Airway hygiene

Please note that if mouth-to-mouth or mouth-to-nose ventilation has been performed, a thorough cleaning of the upper airways, and a replacement of the lower airways are necessary. See section "Cleaning upper airways". However, we recommend that this manikin not be used for mouth-to-mouth or mouth-to-nose ventilation.

Vascular access

Intravenous access via scalp vein


A graphical representation of a scalp vein is indicated on the manikin's face skin. This allows the student to simulate the procedure of locating and attempting to enter a vein. However the scalp vein does not simulate the feel of the vein or flashback of blood.

Intraosseous access via tibia bone

The left leg is designed with a recess that houses a replaceable pad containing a simulated tibia with artificial bone marrow.

For pad placement, refer to section "Assembly - Getting Started".

With the pad in place, intraosseous access can be practised according to clinical protocol.

To confirm that the needle is correctly inserted in the tibia bone, simulated bone marrow can be aspirated into a syringe.

When the needle is withdrawn, the sleeve surrounding the tibia bone will seal around the entrance point so that


further punctures can be performed without excessive leakage. A permanent hole, however, will be left in the tibia bone after each puncture. Up to 15 punctures can normally be performed before pad replacement will be necessary. However, this will vary depending on how densely the punctures are performed.

To replace an intraosseous pad, carefully bend the left foot downward and release the lower part of the pad from the recess in the foot.

Pull out the pad and discard.

Place a new pad as described in the section "Assembly - Getting Started".

Cardiac features

A pulse is simulated via a manually operated pulse bulb.


ECG can be monitored using the three ECG electrodes on the chest skin if the manikin is connected to the 3-lead outlet of the Laerdal Heartsim 2000 or other simulator with compatible outlets as described in the section "Assembly - Getting Started".

Warning!

The Laerdal ALS Baby Trainer must not be defibrillated. The monitoring electrodes are not designed to absorb high voltage / high energy shocks. Application of such shocks may cause a potential shock hazard to the operator and may also damage the connected cardiac rhythm simulator.

Maintenance


Lubrication

It is essential that the upper airway is properly lubricated in order to obtain the intended realism and durability.

Take care to lubricate the airway, endotracheal tube and laryngoscope blade as described in section "Lubrication of the airway" on page 3 before using the manikin.

Replacement of IO pads

- Carefully bend the foot downwards, and while holding it down.
- Lift out the tab at the lower end of the IO pad from the square cavity inside the foot.
- Lift the pad out of the recess and pull out from the round cavity in the knee.
- Insert a new pad as described in section "Placement of IO pads" on page 3.

Replacing or cleaning the face skin

To remove

- Lift off the back of the head.

- Lift the face skin off the tabs at the chin.
- Lift face skin off.


To clean

- To clean face skin, wash by hand in lukewarm, mild soapy water. Then rinse with water and let dry. Alternatively replace with a new face skin (Cat. No. 08 21 30).

To install

- Make sure inside of mouth is clean, see the section "Cleaning Upper Airways".
- Put face skin back on skull and fasten chin tabs.

- Pull face skin approximately 1 cm over the ridge at the back of the skull.


- Pull back of head down over skin until it locks in place.

Replacing the lungs/stomach

To remove


- Unbutton and lift off chest skin.

- Lift out chest plate.

- Remove stomach pressure valve band.

- Unsnap the two lung connectors.

- Unsnap the silicone stomach tube from the oesophagus end of the head/airway assembly.
- Discard lung/stomach.

To install

- Unpack new lung/stomach (Cat. No. 08 20 20).
- Connect stomach tube to oesophagus end of head/airway assembly.
- Place new airway/stomach with the lung connectors into the two holes in the chest plate.
- Connect the two lung tubes to the two lung connectors in chest plate.
- Mount the stomach pressure valve band over the two tabs. Put chest plate in place over the chest compression slide.
- Put chest skin back on. Make sure the cables are placed inside the chest skin so that they do not restrict chest movement.

Cleaning upper airways (Requires cleaning kit, optionally available as Cat. No. 08 15 00. See parts list with illustration.)

To remove

- Open up manikin and disconnect airway/stomach as described in "To remove" in section "Replacing the lungs/stomach".


- Press down tab on inner neck connector and pull head away from body.

- Remove back of head and face skin as described in "To remove" in section "Replacing or cleaning the face skin".

To clean

- Open and empty the cleaning kit basin.
- Place the manikin head face down in the cleaning basin.

- Connect the two lung tube connectors to the double Y-piece in the tube system.


- Connect the oesophagus to the single connection in the tube system.

- Insert the cleaning syringe into the valve opening and place the free end of the tube system into the basin.

- Fill the basin with liquid to just underneath the edge. Circulate liquid through the manikin airways by pumping the syringe plunger.

- Use this setup to go through the following four steps, changing liquid in the container for each new step:

- 1: Use soapy water to remove condensation on interior surfaces.
- 2: Use clean water to remove soap residue.
- 3: Use disinfection solution*. Airways should remain completely filled for at least 10 minutes.
- 4: Use clean water to remove disinfection solution.

After each step lift the head clear of the liquid to allow drainage.

- Allow to dry completely before reconnecting head to manikin.

* A sodium hypochlorite solution freshly prepared for single class use is an effective and recommended disinfectant. This solution should have at least 500 ppm free available chlorine, i.e. 20 ml of 2.5 - 5.0 % household bleach per 500 ml of water. Do not use alcohol to clean the upper airways or the Face Skin.

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