Instructions for Use

VWR Collection Sterile Cleanroom Mask with Earloop
113-6222

Intended use

The VWR Cleanroom face masks are designed to reduce process and products contamination in cleanroom and controlled environments.
Masks rigorously tested and manufactured in an ISO-certified environment; the tests are validated by an independent control laboratory
For use in industrial environment only

Read and Understand before Use

1. **This product is not a medical device in context of Regulation (EU) 2017/745/Directive 93/42/EEC nor personal protective equipment in context of Regulation (EU) 2016/425!**
2. **This product does not prevent Coronavirus transmission and is not intended for general public!**
3. Before use, the wearer should always check the mask is in good condition - not dirt, no damage to head strap attachment etc.
4. Facial hair and certain facial characteristics will prohibit effective use of this product.
5. Please make sure the mask fit properly before entering working environment.
6. Those masks are for single use. Do not reuse nor clean!
7. Dispose the products according to applicable local législations

General Specifications

Item	Test method	Standard
Material	-	21TB+25MB+24Film
Color	-	White
BFE (Bacterial Filtration Efficiency)	ASTM F2101	95.4%
Delta P	MIL-M-36954C	< 3.5
PFE(Particle Filtration Efficiency)	ASTM F2299	>99%
Particle Shedding	IEST-RP-CC003.3 (≥0.5 μ m/cu. ft.)	Level I 196 count/ft3
Earloop breakage test	-	14N for 20seconds
Sterilization (Gamma)	ISO11137	SAL-6

Packaging and Storage Specifications

Package/Quantity	10pcs/pack, 5packs/Inner Bag, 10Inner Bag/Case
Best storage Conditions	10° - 30°C
	40% - 75% humidity
	dry
Shell Life	3 years
Country of Origin	China

Instructions for Use

FITTING INSTRUCTION

(English)
1. Put loops around ears.
2. Press forefingers down on both sides of nose.
3. Adjust mask over nose & under chin.

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