VWR INCU-Line Shake
SHAKING INCUBATORS

MODELS: ILS4, ILS6

INSTRUCTION MANUAL

EUROPEAN CATALOGUE NUMBERS:

444-0732	INCU-Line ILS 4
444-0733	INCU-Line ILS 6

VERSION: 3
ISSUED: 09-18- 2014

CE

Legal Address of Manufacturer

United States
VWR International, LLC
100 Matsonford Rd
Radnor, PA 19087
800-932-5000
http://www.vwr.com

Europe
VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
http://be.vwr.com

Country of Origin
United States

TABLE OF CONTENTS

SECTION 1.0	INTRODUCTION
SECTION 2.0	SYMBOLS AND CONVENTIONS
SECTION 3.0	SPECIFICATIONS
SECTION 4.0	SCOPE OF DELIVERY
SECTION 5.0	INSTALLATION
SECTION 6.0	CONTROLS OVERVIEW
SECTION 7.0	OPERATION
SECTION 8.0	CLEANING AND ROUTINE CARE
SECTION 9.0	SERVICE AND TROUBLESHOOTING
SECTION 10.0	PARTS LIST
SECTION 11.0	WARRANTY
SECTION 12.0	EQUIPMENT DISPOSAL

Section

1

INTRODUCTION

The VWR Shaking incubators have been designed for general purpose incubations (up to 80ºC). These units are high performance shaking incubators for professional, industrial or educational use where the preparation or testing of materials is done at approximately atmospheric pressure and no flammable, volatile or combustible materials are being heated. These units are not intended for hazardous or household locations or use.

Temperature within the shaking incubator is controlled and regulated by a microprocessor. Temperature is set and displayed digitally. An independent, user settable safety thermostat protects samples and the shaking incubator from overheating in the unlikely event of a primary controller failure.

A fan, located in the chamber, assists in providing an exceptionally uniform environment. The fan also speeds temperature recovery after the door of the shaking incubator has been opened. For applications where the flow of air in the chamber is undesirable, the fan may be switched off.

The interiors of the ILS4 and ILS6 shaking incubators are constructed of stainless steel. The exteriors are made of cold rolled steel to resist corrosion and provide strength. The door of the shaking incubator consists of two parts. The inner door is glass to allow for viewing the chamber contents without disrupting the temperature environment. The outer door is solid and insulated. It is also heated to help in maintaining the chamber temperature. The heated door prevents condensation from forming on the inner glass viewing door.

Your satisfaction and safety require a complete understanding of this unit. Read this instruction manual thoroughly before attempting to operate the shaking incubator. All operators should be given adequate training before using this shaking incubator. Note that any alterations or modifications to the shaking incubator may void the warranty.

Section

2

SYMBOLS AND CONVENTIONS

Your VWR ILS4 and ILS6 incubator shakers use internationally accepted graphic symbols to help convey information to the user and to call the users attention to important safety precautions and guides for using this equipment.

 Indicates that user should consult manual further description or discussion.

 Indicates AC “Power On”

 Indicates “Heating”

°C Indicates “Degrees Centigrade”

 Indicates “Temperature”

 Indicates Potential “Shock Hazard”

 Indicates Protective “Earthground”

SAFETY PRECAUTIONS

Use of this product in any manner not specified by the manufacturer may impair the safety protection provided by the equipment and may result in physical damage and/or personal injury. Please read all operating instructions in this manual prior to use of this equipment.

1. Do not operate this unit in an explosive or flammable environment.
2. Do not incubate or shake flammable or explosive materials or highly reactive chemicals.
3. **Lifting / Handling**: These units are heavy and care should be taken to use appropriate lifting devices. Units should only be lifted from their bottom surfaces and not by doors, handles or knobs.
4. **Leveling:** The unit must sit level and solidly on the four leveling feet.
5. **Load testing:** Test all loads to be shaken under observation to insure load and unit stability.
6. **ILS6 Units are not stackable.** Do not stack one unit on top of another unit.

Section

3

SPECIFICATIONS

Temperature Range	Ambient +5° to 80°C
Temperature Uniformity	
ILS4	±0.2° @ 37°C
ILS6	±0.25° @ 37°C
Temperature Accuracy	±0.2°
Temperature Set / Display	Digital
Temperature Control	Microprocessor
Overtemperature Safety	Independent, settable
Door	Inner – glass
	Outer - solid, heated
Chamber Volume	
ILS4	0.048 cubic meters
ILS6	0.078 cubic meters
Exterior Construction	Cold rolled steel
Interior Construction	Stainless steel
Exterior Dimensions	
ILS4 (cm)	68.60W x 50.80D x 73.66H
ILS6 (cm)	78.74W x 66.04D x 86.36H
Interior Dimensions	
ILS4 (cm)	43.18W x 35.60D x 50.80H
ILS6 (cm)	53.34W x 50.80D x 63.50H
Standard Accessories	Adjustable leveling feet 2 Stainless Steel Shelves Adjustable Shelving Points
Electrical Ratings	230V, 50/60Hz, 4 amps
Fuse	5x20MM, T6.3A-250V

Section

4

SCOPE OF DELIVERY

4.1 **Inspection:** The carrier, when accepting shipment, also accepts responsibility for safe delivery and is liable for loss or damage. On delivery, inspect for visible exterior damage, note and describe on the freight bill any damage found, and enter your claim on the form supplied by the carrier.

4.2 Inspect for concealed loss or damage on the unit itself, both interior and exterior. If any, the carrier will arrange for official inspection to substantiate your claim.

4.3 **Return Shipment:** Save the shipping crate until you are sure all is well. If for any reason you must return the unit, contact VWR International, Inc. for authorization and supply the data plate information.

4.4 **Accessories:** Verify that all of the equipment indicated on the packing slip is included with the unit. Carefully check all packaging before discarding. This unit is equipped with 2 large shelves, 8 shelf clips and 4 leveling feet. If any of these items is missing, contact VWR International's Customer Service Department immediately.

Section

5

INSTALLATION

Under normal circumstances this unit is intended for use indoors, at room temperatures between 5° and 40°C, at no greater than 80% Relative Humidity (at 25°C) and with a supply voltage that does not vary by more than 10%. Customer service should be contacted for operating conditions outside of these limits.

5.1 **Location:** In selecting a location, consider all conditions that might affect performance, such as heat from radiators, ovens, autoclaves, etc. Avoid direct sun, fast-moving air currents, heating/cooling ducts and high-traffic areas. Allow a minimum of 5cm between the unit and walls or partitions that might obstruct free airflow.

5.2 **Power Source:** Check the data plate for voltage, cycle, phase and ampere requirements. If matched to your power source, plug the power cord into a grounded outlet. VOLTAGE SHOULD NOT VARY MORE THAN 10% FROM THE DATA PLATE RATING. These units are intended for 50/60Hz applications. A separate circuit is recommended to preclude loss of product due to overloading or circuit failure. ***NOTE: Electrical supply to the unit must conform to all national and local electrical codes.***

5.3 **Lifting / Handling:** These units are heavy and care should be taken to use appropriate lifting devices that are sufficiently rated for these loads. Units should only be lifted from their bottom surfaces. Doors, handles and knobs are not adequate for lifting or stabilization. The unit should be completely restrained from tipping during lifting or transport. All moving parts, such as shelves and trays should be removed and doors need to be positively locked in the closed position during transfer to prevent shifting and damage.

5.4 **Leveling:** The unit must sit level and solidly. Leveling feet are supplied and must be installed in the four holes in the bottom corners of the unit. With the feet installed and the unit standing upright, each foot can be raised by turning it in a counterclockwise direction. Adjust the foot at each corner until the unit stands level and solid without rocking. If the unit must be moved, turn the leveling feet all the way clockwise to prevent damage while moving.


5.5 **Cleaning:** The unit chamber should be cleaned and disinfected before use. Remove all of the interior parts, and clean thoroughly with a disinfectant that is

appropriate to your application. DO NOT use chlorine-based bleaches or abrasives as this will damage the chamber interior. Regular periodic cleaning is required.

Section

6

CONTROLS OVERVIEW

From left to right

6.1 SHAKER CONTROL SYSTEM:
Run green signal light: lit when shaker is running.

Time yellow signal light: lit when the shaker is set on time function.

RPM yellow signal light: lit when the shaker is set on rpm function.

6.2 Control Knob (Push to Shift / Turn to Set): by rotating the control knob right (+) or left (-) you can change the TIME or RPM. Push the control knob to switch between TIME and RPM values. If you rotate knob fast then the values go up or down on display very quickly.

6.3 Start / Stop button: Press to start or stop operation.

6.4 Light Switch: On turns internal light on.

6.5 Temperature Control: This control consists of the digital display and Up/Down buttons for inputting set point temperatures and calibration.

6.6 Heating Indicator: When light is on incubator is heating

6.7 Safety Activated: This pilot light comes ON when the Safety Thermostat is activated. Under normal operating conditions this pilot light should never be on.

6.8 Power Switch: The main power ON/OFF switch controls all power to the unit. It must be on before any systems are operational. The green switch will be lighted when ON.

Section

7

OPERATION

7.1 **Main Power:** Press the power switch to on. The unit will automatically detect supply frequency of 50 or 60Hz and displays F50 or F60. Then the display switches to show the set time and illuminates the time light.

7.2 **Set Main Temperature Control:** The controller has a 3-digit display for displaying chamber temperature or set point information. Up and Down arrow pads are used to change the set point and controller mode of operation. To enter the set point mode of operation, press either the Up or Down pad one time. The display will start to blink, going from bright to dim. While blinking, the display is showing the set point. To change the setpoint use the Up and Down arrow pads. If the arrow pads are not pressed for five (5) seconds, the display will stop blinking and will read the chamber temperature. After setting temperature, allow at least one hour for the chamber temperature to stabilize and 24 hours for maximum stabilization.

7.3 **Calibration:** The Main Temperature Controller was calibrated at the factory at 37°C; however, it is recommended that calibration is done once the unit is installed in its working environment and has been stable at set point for several hours. Place a certified reference thermometer in the chamber. Be certain the thermometer is not touching any shelving. Allow the temperature to stabilize again until the thermometer reads a constant value for one hour. Compare the digital display with the reference thermometer. If there is an unacceptable difference, put the display into calibration mode by pressing both the Up and Down arrow pads at the same time until the two outside decimal points begin to flash. While the decimal points are flashing the display can be calibrated by pressing the Up or Down arrow pads until the display reads the correct value. Allow the shaking incubator temperature to stabilize again, and recalibrate if necessary.

7.4 **Set Safety Thermostat:** The safety thermostat is located on back of unit and is manually set. It is completely independent of the Main Temperature Controller. As mentioned in step 7.3, the Safety Thermostat should be initially set to its maximum position, to allow the unit to stabilize. Once the shaking incubator is stable at the desired set point, turn the Safety Thermostat counterclockwise until the Safety Light turns on. Next, turn the Thermostat clockwise just until the Safety Light turns off. Then turn the Thermostat clockwise two (2) of its smallest divisions on its scale past the point where the Safety Light went out. This will set the Safety Thermostat at approximately 1^{o}C above main temperature set point.

7.5 Shaker Control System

- Time signal light is on. By rotating the Control Knob right (+) or left (-) the set time can be selected from 30 sec to 99 min 50 sec :

 99.5 ⇒ 99 min 50 sec
 9.59 ⇒ 9 min 59 sec
 0.30 ⇒ 30 sec

- If you want to set the timer to hold – turn the knob left or right until HLd is displayed. You can reach the HLd function by setting to under 0.30 or above 99.5

- Push down on the Control Knob to switch modes between time set mode and RPM set mode.

- RPM signal light is on. By rotating the Control Knob right (+) or left (-) the desired RPM can be selected :

 34 ⇒ 340 Rpm
 120 ⇒ 1200 Rpm

- Press START/STOP button to start operation

- RUN and TIME signal lights on. Shaker runs and counts down time from set value.

- NOTE: If there is a power outage during operation, shaking incubator will restart automatically, as soon as the electricity returns. Display will be flashing (warning you that the power was cut). You can turn the flashing off, by pressing the Control Knob.

- If you want change the RPM during shaking, push the Control Knob to put the shaker into RPM set mode (RPM signal light will be on). Rotate the Control Knob right (+) or left (-) to set the desired value. In the meantime RPM signal light will flash. When you stop rotating the Control Knob the signal light for RPM will stop flashing after 2 sec.

- When the run time has elapsed or when you press the START/STOP button, the display will show “End” and the RUN light will flash. When shaker completely stops the set speed and time will return to the last programmed values.

7.6 ADDITIONAL OPERATIONS

- **Door Open Switch:** This switch senses when the incubator shaker door is open and sends a signal to microprocessor which causes the shaker to stop running as a safety measure. The timer will pause and re-start when door is closed.

- If you want to view the set value for RPM during shaking rotate the Control Knob for ONE CLICK right (+) or left (-). The RPM signal light will then flash for 2 sec and the display will show the set RPM. After 2 sec the RPM and the light will stop flashing.

 NOTE: Shaker must be in RPM mode – RPM signal light must be on.

- If you want to view the set value for TIME during shaking rotate the Control Knob ONE CLICK right (+) or left (-). The TIME signal light for 2 sec and the display will show the set TIME. After 2 sec the display will return to the count TIME and the signal light will stop flashing.

 NOTE: Shaker must be in TIME mode – TIME signal light must be on.

Section

8

CLEANING AND ROUTINE CARE

No routine maintenance is required for the electrical or mechanical components of the unit. Both the interior and exterior of the shaking incubator should be wiped down with a soft cloth periodically to prevent any build-up of dust and grime.

Any spills in the shaking incubator should be cleaned immediately. Cleaning may be performed with a soft, damp cloth. If required, the chamber may be disinfected. Before using any type of disinfectant or cleaning agent in the shaking incubator, contact VWR Customer Service Department to ensure that the product will not damage the stainless steel interior of the shaking incubator. After cleaning the shaking incubator, rinse the cleaned surface with a damp cloth, using water only. Dry with a clean cloth.

Do not pour liquid into the chamber or immerse the unit in liquid. DO NOT use chlorine-based bleaches or abrasives as these will damage the stainless steel.

Section

9

SERVICE AND TROUBLE-SHOOTING

9.1 Troubleshooting Guide

Problem	Possible Solution
Oven will not power up or not heat	Check power cord, outlet and unit circuit breaker
Shaking platform will not operate	Check time set to Hld or a time setting Check that nothing is blocking shaker platform Press Stop / Start button
Temperature too high	Check set point and readjust if necessary Check calibration
Chamber temperature goes above set point and settles back to set point	Normal operation in initial heat up of if door opened for long period of time.
Temperature will not remain stable or display shows “LO”	Check that set point is at least 5°C above ambient (or 10°C above ambient when shaker is on) which is minimum operating point.
Indicated temperature is unstable	A slight variation of +/- 0.1°C is normal. Larger fluctuations may be ambient variations from drafts, door opening and closing, a fan obstruction or failure or electrical noise from RFI (motors, etc..)
Temperature is too low	If door has opened, unit may not have recovered yet. Confirm temperature set point. Safety thermostat set too low.
Unit will not heat above temperature that is below set point.	Confirm set point. Check temperature of chamber with a thermometer and recalibrate if needed.
Temperature display and reference thermometer do not match.	Be sure that unit has been allowed to stabilize for 1 hour. Thermometer should be suspended in the chamber and not touching any surface. Only certified reference thermometers should be used.
Cannot adjust set point or calibration or display persists call service	Turn unit off for 5 seconds to reset. If problem persists, call for service.
Unit calibrated at one temperature and not at another	This can be a normal condition if temperatures or load vary widely. For best accuracy, calibrate at set point.

Power key does not light	Check the power source Check fuses
Shaker stalls - Message <Er1>	Power unit down, then restart
Message <Er2>	Contact VWR service
Message <Er3>	Contact VWR service
Message <Er4>	Contact VWR service

MECHANICAL

Mech-1 Door not sealing	
	1/ Check physical condition of gasket. **2/ Tighten door latch till it pulls glass in.** **3/ Assure that gasket is in original location.** **4/ Adjust hinge blocks or twist the door.** **5/ Confirm that unit has not been damaged and body is square.**
Mech-2 Motor doesn't move	
	1/ If shaft spins freely: check connections to motor and check voltage to motor. **2/ If shaft rubs or is frozen, relieve binding and retest.**
Mech-3 Motor makes noise	
	1/ If noise is from the motor, tap the top of motor shaft with ball peen hammer. Remove back panel to access fan motor. **2/ If the sound gets worse, tap the other end of the shaft – avoiding touching the fan blade.** **3/ If there is no change, call VWR Customer Service.** **4/ If noise is from shaft or fan blade, realign shaft.**

OTHER

Other-1	**Controller on at all times - "locked-up"**	
		1/ Turn unit off and on to reset. **2/ If cannot change any condition on the front panel, call VWR Customer Service.**
Other-2	**Front panel displays are all off**	
		Check for wire damage.
Other-3	**Unit or wall fuse/circuit breaker is blown**	
		1/ Check wall power source. **2/ Compare current draw and compare to specs on data plate.** **3/ See what other loads are on the wall circuit.**
Other-4	**Unit will not turn on**	
		1/ Check wall power source. **2/ Check fuse/circuit breaker on unit or in wall.** **3/ See if unit is on, e.g., fan or heater, and just controller is off.** **4/ Check all wiring connections, especially around the on/off switch.**
Other-5	**Unit is smoking – Out of box**	
		This is not an uncommon occurrence when first operating new units. Put unit under vent and run at full power for one hour. – smoking is normal during first cycle to temperature.
Other–6	**Contamination in chamber**	
		1/ See cleaning procedure in see section 8.0. **2/ Develop and follow standard operating procedure for specific application; include definition of cleaning technique and maintenance schedule.**

Technical service

Web Resources

Visit the VWR's website at **www.vwr.com** for:
• Complete technical service contact information
• Access to VWR's Online Catalogue, and information about accessories and related products
• Additional product information and special offers
Contact us For information or technical assistance contact your local VWR representative or visit. **www.vwr.com**.

Section
10

PARTS LIST

User replaceable accessories and spare parts

Description	Quantity	Cat. No.
SHELF FULL for ILS6	1	444-0831
SHELF FULL for ILS4	1	444-0832
SHELF HALF for ILS4	1	444-0833
ILS4 INCUBATOR STACKING ADAPTER	1	444-0913
Fuse	1	444-0828
Power Cord, EU	1	444-0827
Power Cord, UK	1	444-0829
Power Cord, Swiss	1	444-0830

Section

11

WARRANTY

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Section

12

EQUIPMENT DISPOSAL

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.
Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.
For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland
VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy
VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-467791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway
VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax. 058 32 38 205
E-mail: info@pl.vwr.com

Portugal
VWR International -
Material de Laboratório, Lda
Edificio Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain
VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden
VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland
VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey
Pro Lab Laboratuar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu Işmerkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK
VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 86 86
E-mail: uksales@uk.vwr.com

Australia
VWR International, Pty Ltd.
Unit 1/31 Archimedes Place
Murarrie, Queensland 4172
Tel.: 1300 727 696
Fax: 1300 135 123

China
VWR International China Co., Ltd
Rm.219, 2100 Dongming Road
Pudong New District
Shanghai 200123
Tel.:+86-21-5898 6888
Fax:+86-21-5855 8801
E-mail: info_china@vwr.com

India
VWR Lab Products Private Limited
135/12, Brigade Towers, 2nd Floor
Front wing, Brigade Road,
Bengaluru, India – 560 025
Tel.: +91-80-41117125/26 (Bengaluru)
Tel.: +91-2522-647911/922 (Mumbai)
Fax: +91-80-41117120
E-mail: vwr_india@vwr.com

New Zealand
Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Singapore
VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3760
E-mail: sales@sg.vwr.com

GO TO VWR.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS OF YOUR LOCAL VWR DISTRIBUTOR

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