Technical Product Specification

Product

Name: PIPETTE HOLDER VI-3 BOX OF 100 SIPHON WHITE
Item no: 140413
Version: 2022-06

Material: PETG clear / TPE white



Dimensions: ø10 x 60 mm

Packing

Box:	100 pcs./ box

Labelling

Internal:	Vitrex
External:	

Regulation

Intended use:	For holding during filling and emptying of End-to-end pipette. For single use
ISO:	
CE:	EU Regulation 2017/746 on in vitro medical devices.
Storage:	Dry, at room temperature

IFU

1. Loading: Insert the pipette into the pipette holder. The end of the pipette must be visible in the glass tube.
2. Filling: Hold the pipette horizontally enough to fill by capillary action. Be sure not to cover the hole of the pipette.
3. Emptying: Cover the hole of the pipette holder and compress the bulb to empty the pipette. Rinse out with reagent by filling and emptying as above.

Vitrex Medical a/s
Vasekaer 6-8 - DK 2730 Herlev
Phone +45 4494 7011
vitrex@vitrex.dk
www.vitrexmedical.com