USER MANUAL
MANUAL DE USUARIO
MANUEL D'UTILISATION
BEDIENUNGSANLEITUNG

Microbiological Safety Cabinets
Cabinas de Seguridad Microbiológica
Cabines de sécurité microbiologique
Mikrobiologische Sicherheitswerkbänke

BIO II ADVANCE PLUS

Technology. Science. Life.

USER'S MANUAL

Microbiological Safety Cabinet

BIO II ADVANCE PLUS

Technology. Science. Life.

Translation of the original manual

-

Failure to read, understand and follow the instructions in this manual may result in damage to the unit, injury to operating personnel and poor equipment performance.

Our company reserves the right to update, improve, amend or otherwise modify the design or specifications of any product at any time without prior notice.

SUMMARY

0. SYMBOLS AND ABBREVIATIONS 4

1. GENERAL DATA 5

2. DESCRIPTION OF THE SAFETY CABINETS 6
- 2.1. Working principle 6
- 2.2. Applications and principe diagram 6
- 2.3. Components description 6
- 2.4. Standard equipment 7
- 2.5. Certification 7
- 2.6. Technical features 8

3. INSTALLING INSTRUCTIONS 9
- 3.1. Transport 9
- 3.2. Placing on site 9
- 3.3. Precautions for first connection 11
- 3.4. Electrical connection and services 12
- 3.5. Cabinet Exhaust Ducting 13
 - 3.5.1. Outside ducting via a canopy 13
- 3.6. Installation tests 13

4. OPERATING INSTRUCTIONS 14
- 4.1. Start-up 14
- 4.2. Front panel positions 14
- 4.3. General working recommendations 16
 - 4.3.1. Ergonomics 17
- 4.4. Definition of the working area 20
- 4.5. Control keyboard-display 21
- 4.6. Disconnection 27
- 4.7. Prolonged stoppage 27
- 4.8. Optional accessories 28

5. MAINTENANCE 29
- 5.1. Maintenance table 29
- 5.2. Absolute filters 30
- 5.3. Spare parts 31

6. CLEANING AND DISINFECTION 32
- 6.1. Surface cleaning and disinfection 32
- 6.2. Front glass cleaning and disinfection 32
- 6.3. Disinfection with formaldehyde 32
- 6.4. Disinfection with vaporized hydrogen peroxide 32

7. MAINTENANCE TABLE 33

8. ELECTRIC DIAGRAMS 33

9. TESTS 33
- 9.1. Downflow filters integrity test 33
- 9.2. Exhaust filters integrity test 33
- 9.3. Downflow air velocity test 33
 - 9.3.1. Set point alarm downflow air velocity test 33
- 9.4. Inflow air velocity test 33
 - 9.4.1. Set point alarm inflow air velocity test 33
- 9.5. Smoke test 33

10. CERTIFICATES 33
- 10.1. Declarations of conformity 33

11. TROUBLESHOOTING 34

0. SYMBOLS AND ABBREVIATIONS

List of stickers used in the biological cabinets.

"EARTH WIRE PROTECTED" STICKER.	
"FUNCTIONAL EARTH WIRE" STICKER.	
"BIOHAZARD" STICKER.	
"ATTENTION, ELECTRIC SHOCK RISK" STICKER.	
"SOCKET VOLTAGE" STICKER.	
"SOCKET MAXIMUM CURRENT" STICKER.	

"SAFETY WARNING PRIOR TO SERVICE OR REPAIR ACTIONS" STICKER

"FRONTAL MECHANISM OPENING" STICKER

Warning: before lifting or lowering the front opening mechanism, always ensure the window is fully closed.
Failure to do so may result in window cracking.

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Abbreviations used in this manual:

V	Volts
A	Amperes
W, kW, µW	Watts, kilowatts, microwatts
mm, cm, m	Millimetres, centimetres, metres
kg	Kilogram

s, h	Seconds, hours
Hz	Hertz
dB	Decibel
ºC	Celsius
GS	"Tested safety" according to German standards

1. GENERAL DATA

Operating Instructions Manual for Microbiological Safety Cabinets Class II **Bio II Advance Plus**, according to EN 12469 standard, which ensure (applying correct techniques in laboratory) high level protection for the operator, the environment, and the sample.

WARNING ⚠

This Manual must be read and understood before its installation, maintenance, and operation of the biological safety cabinet Class II. All the safety instructions given in the Manual must be followed.


Any use or application failing to comply with the recommendations expressed in the Manual may cause danger and loss of rights to manufacturing guarantee as well as freeing the manufacturer from all responsibility arising.

Any maintenance must be done by qualified staff according to the laboratory working correct conditions. In case these actions are not carefully carried out, as a result the safety cabinet can work in abnormal conditions.

All maintenance actions on this equipment must be performed by a qualified technician who is familiar with the proper maintenance procedure required for this equipment.

The European Union issued the Directive 2012/19/EU on Waste of Electrical and Electronic Equipment (WEEE).

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive. Mark symbol is shown at the right.
Syntegon Telstar sells products through distributors throughout Europe. Contact your local distributor for recycling/disposal.

The name plate & CE mark sticker is located at one side of the cabinet.

MANUFACTURER: **SYNTEGON TELSTAR TECHNOLOGIES, S.L.U.**
Av. Font i Sagué, 55
Parc Científic i Tecnològic Orbital 40
08227 TERRASSA (SPAIN)
Tel. (+34) 937 361 600
Fax (+34) 937 861 380
E-mail**:** telstar@telstar.com
http://www.telstar.com

2. DESCRIPTION OF THE SAFETY CABINETS

2.1. WORKING PRINCIPLE

CLASS II Biohazard Safety Cabinets are designed for microbiology work with the most advanced microprocessor control, ensuring a high level of protection for the operator, environment, and sample.

Exterior of cabinet made of laminated, lacquered, oven-dried steel. Functional design with reduced external dimensions and front access for changing filters.

Spacious internal chamber with sliding-swinging front made from laminated safety glass. Segmented work surface and liquid collection tray made of polished stainless steel. Armrest is detachable to facilitate cleaning of liquids tray.

Two-fan system (four in the Bio II Advance Plus 6) providing laminar and exhaust flow by means of a shared plenum. In the case of a breakdown of any fan, the other ones continue to provide airflow.

As can be seen in **Annex A2**, the air propelled by the fans (1) is moved into the chamber or plenum (3) and, on the one hand via the Absolute HEPA filter (4), is filtered and flows according to laminar regime in the working area (6). On the other hand, part of the air flows through an Absolute HEPA filter (5) before being expelled to the outside (8).

System of recycled (70%), sterile vertical laminar flow in depression and exhaust (30%) sterile air, approximately.

2.2. APPLICATIONS AND PRINCIPE DIAGRAM

The microbiological safety cabinet models are appropriate for handling biologically pathogenic samples or products (up to microorganism risk BSL3), ensuring sterile air in the work zone, as well as a high level of protection for the operator and the environment.

Principle diagram in **Annex A1.**

WARNING ⚠

Under no circumstances, is possible handling in corrosive, flammable, or explosive gas ambient. It is recommended to expel the extracted air to the exterior. Right working conditions are obtained only when the display is green and without any alarm message on the display.

2.3. COMPONENTS DESCRIPTION

(See Annex A2)

2.4. STANDARD EQUIPMENT

- Two inlet sockets 2A - 460W each located on either side of the work zone. Insulation level IP44.

- UV kit

Easy-to-work microprocessor-controlled UV germicidal lamp with a timer function. It can only be enabled in the *Disinfection* screen (**Chapter 4.5 point 2.4**) or by means of the **UV light Auto Power on** screen (**Chapter 4.5 point 3.5.1.2**). For safety reasons, this button is electrically linked with the light button, and it switches off if the front window is opened.

Average disinfection intensity over the work zone over 40 μW/cm^2 with 20% depreciation after a working life of 8000 hours.

WARNING ⚠

In the event of using this accessory, you should consider the normal precautions in the use of these types of lamps. Their use is only suitable for previous external sterilization of materials introduced in the cabinet. So, this lamp accessory should be used before the start of work in the cabinet, and it is essential to keep it turned off when the operator is in front of the working area or working on it. The glasses provide protection from direct radiation, but UV rays are reflected on polished metal parts and can seriously damage operator's eyes. That is why the glass front must be placed into the closed position to be able to switch on the U.V light, then automatically the fans will disconnect.

2.5. CERTIFICATION

Quality control and tests results certification carried out in line with standard EN 12469. See Certificate inside the cabinet.

2.6. TECHNICAL FEATURES

Features	Units.	Bio II Advance Plus 3	Bio II Advance Plus 4	Bio II Advance Plus 6
External dimensions (LxWxH)	mm	1049 x 764 x 1260	1354 x 764 x 1260	1964 x 764 x 1260
Internal dimensions (LxWxH)	mm	954 x 602 x 580	1259 x 602 x 580	1869 x 602 x 580
Front opening height	mm	182		
Weight	Kg	180	200	280
Downflow velocity	m/s	0.35	0.35	0.35
Inlet velocity	m/s	0.55	0.55	0.50
Voltage *	V	120 or 230		
Frequency	Hz	50/60		
Lighting	Lux	≥ 1000		
Noise	dB (A)	≤ 58 (Attached table)		
Vibration	mm RMS	< 0.005		
HEPA/ULPA H 14 Filters		Efficiency 99.995%, EN 1822 by MPPS method filtration system in laminar and exhaust flow (99.999% D.O.P. test @ 0.3µm particles)		
Fans		Two high-efficiency fans (four in cabinet Advance Plus 6). Flows measured by a thermoanemometric probe.		

* To ensure proper installation, check the attached electrical diagrams.

DISSOCIATED DECLARED VALUES OF NOISE EMISSION Standard(s) to which conformity is declared ISO 4871			
	$V_{downflow}$ = 0.35m/s V_{inlet} = 0.55 m/s (230 V, 50 Hz)		$V_{downflow}$ = 0.35m/s V_{inlet} = 0.50 m/s (230 V, 50 Hz)
	Bio II Advance Plus 3	**Bio II Advance Plus 4**	**Bio II Advance Plus 6**
L_{pA} (dB)	≤58	≤58	≤58
K_{pA} (dB)	2	2	2
Values measured according to acoustic test in standard EN 12469:2000 following standard ISO 11201			

3. INSTALLING INSTRUCTIONS

For a correct start-up procedure of the microbiological safety cabinet, observe the following step-by-step instructions.

3.1. TRANSPORT

The cabinet is supplied in an individual package and in the suitable way to prevent it from being damaged.

This manual and certificate of individual Quality controls of each cabinet are included in the scope of supply.

3.2. PLACING ON SITE

a) Special care should be taken when unpacking the front glass and the control display.

b) Use the same packaging base for the transport to its final location and it is highly recommended doing it with a pallet jack. The cabinet has two M10 eyebolts to help to position it on a raised surface. If there are no mechanical means to position it, this can be done manually.

WARNING ⚠

We advise against the use of lift trucks for placing the cabinet directly by its base. In such case, it should be handled with care and it is also a must that:
- the useful length of **the forks of the lift-truck are over 750 mm**.
- the **cabinet is fully supported all along its width** (front (A) and rear (B) structural profiles). Direct visual inspection during its handling is required.


Fig. 1 X

Fig. 2 ✓

Fig. 3 ✓

WARNING

For any handling during assembly or dismantling of the cabinet, appropriate safety measures should be taken, such as wearing safety gloves or shoes.

c) Cabinets should be placed away from entrances, corridors, ventilation supply points, etc., which might affect airflow in the work area.

d) Ensure the cabinet is set perfectly on the table or workbench to prevent vibrations.

e) To prevent the work area from overheating, we recommend that the room temperature is always within the comfort margins for working.

The optimal comfort margins are:

- Temperature: 18-28ºC
- Relative humidity: 50%±20%

If any of these conditions is not fulfilled, optimum cabinet performance cannot be guaranteed.

f) When necessary, air vents in the room should be modified or speed of air entry reduced so that air is not aimed directly at the cabinet, as this may produce turbulence.
Side air vents are not appropriate for rooms where cabinets are installed.

g) Due to the ceiling temperature affects air distribution, sites should be selected with enough space. We recommend 250 mm minimum free space between the exhaust and the ceiling.

h) The location of the cabinet will directly influence the noise level. This level will rise considerably in small spaces with very reflective surfaces because of the proximity of the cabinets to the walls, and where fan speed has risen to compensate for air currents in the room or the resistance of clogged filters.

i) Correct cabinet performance should be tested before putting it into operation if any modification has been made to the room or to the air intake system which may considerably change the airflow characteristics in the room.

3.3. PRECAUTIONS FOR FIRST CONNECTION

- Before connecting it to the mains supply, proceed to make a general cleaning to remove any dust particles that may be accumulated during transportation. Follow the instructions on cleaning and disinfection.


- Before using it for the first time, run the cabinet for around 6 hours to purge the filters. New filters give off a characteristic smell which tends to disappear in a short period of time.

WARNING ⚠

The cabinet should work at the voltage shown on the name plate. It should be connected to a suitable socket base (minimum 16A) with an earthed plug. The cable section recommended for the base is 1.5 mm^2.

WARNING ⚠

The plug connection must be in such a place to prevent it from being unplugged intentionally or accidentally during operation, e.g., next to a wall or closed to another equipment.

WARNING ⚠

Unplug the cabinet from the electric current supply, in case of electrical maintenance.

The cabinet can only be activated from the button in the frontal control keyboard.

The installation which the cabinet is connected to should meet the safety requirements in force in each country.

Exhaust air should be evacuated outside the lab if working with volatile substances not neutralized by absolute filters.

During normal work in these safety cabinets and because of a basic principle of recirculation and depression in the work zone, no perforated surface in the work zone should be covered with paper or sanitary fabric.

3.4. ELECTRICAL CONNECTION AND SERVICES

The electrical characteristics of the cabinet are as follows:

- Electrical class: Class I.
- Insulation level: IP-20.

The cabinet is supplied with a plug cord according to European standard (other plug cords are available as optional accessories). The on/off button is in the frontal control keyboard-display. (**Annex A3 pos. 1**).


The power socket in the wall should be in a perimeter around the cabinet of:

- 1.25 m from the middle of the cabinet to the sides
- 2.5 m from the base of the cabinet upwards

On the top right-hand side of the cabinet and identified by name are:


<u>230V und 120V versions (except for 6 ft size)</u>

- Voltage input (1)
- Checking valve fans pressure (2)
- Checking valve for D.O.P. control (3)

<u>Bio II Advance Plus 6 120V</u>

- No use (1)
- Checking valve fans pressure (2)
- Checking valve for D.O.P. control (3)
- Main fuse (4)
- Voltage input (5)

3.5. CABINET EXHAUST DUCTING

Ducting of the exhaust may be preferred to exhaust in the lab to, for example:

- Prevent from lab temperature rise.
- When volatiles and toxic vapors are involved during the process.

WARNING ⚠

Ducting the exhaust to the outside of the lab creates a depression in the lab.

3.5.1. OUTSIDE DUCTING VIA A CANOPY

In case of direct ducting outside it must be installed an additional fan to the duct, a canopy, and adjust the velocity parameters. GS has not been tested. It must be tested after installation.

WARNING ⚠

When the cabinet exhaust is connected to an assisted/air forced exhaust system, it is necessary to install a canopy between the cabinet and the exhaust duct. This will prevent from the loose of operator protection in case of any decrease of the exhaust airflow system (duct obstruction, blower control failure…).

See description & installation in **Annex A4.4.**

3.6. INSTALLATION TESTS

According to standard DIN EN 12469 and to ensure correct cabinet operation, after the installation and connection of the cabinet, it is compulsory to perform the following tests (before the equipment is used for the first time):

- Impulsion filters integrity test.
- Exhaust filters integrity test.
- Downflow air velocity.
- Downflow set points alarm test.
- Inflow air velocity test.
- Inflow set points alarm test.
- Smoke test.

Syntegon Telstar also recommends making the following tests:

- Light test.
- Noise level test.

Syntegon Telstar, its representatives, and distributors put a qualified technical team at your disposition that uses control equipment with maximum health protection for carrying out this installation test under full safety guarantees.

4. OPERATING INSTRUCTIONS

4.1. START-UP

Once the above points have been controlled, start up the microbiological safety cabinet by pressing the On/Off button (**Annex A3, pos. 1**). After entering the password correctly (by default the password is **0000**) or immediately if this function has been disabled, the fans will automatically start working.

Ensure that the aspiration area and air exhaust outlet are free (**Annex A2, pos. 7 & 8**).

The electronic system of regulation by microprocessor via a thermoanemometer probe, laminar down flow and exhaust airflow, needs a few seconds to stabilize during connection in order to begin to work normally.

The cabinet can only be used by qualified staff who meets the correct work techniques in the laboratory. If these techniques are not fulfilled, faulty working of the equipment may happen.


WARNING ⚠

The EN 12469 certification does not allow the deactivation of the password enquiry.


4.2. FRONT PANEL POSITIONS

By means of the "A" handles, the front glass can be moved from working position to closed position or vice versa. To move the glass just slide down or up the glass.

FROM WORKING POSITION


FROM CLOSED POSITION

In case it is needed to clean the glass by the inner part (Annex A5.2) or if it is necessary to place in the working area any big product or device, the glass-frontal group may be separated from the cabinet, providing full access to both the inner glass area and the working area. All needed precautions should be taken as the user protection decreases during that operation.

WARNING ⚠

To open the frontal (after the glass has been moved to the lower position), always pull the lateral handles "B" fixed to the metal supports never pull the handles of the glasses, that are specifically designed for the vertical movement of the glass. Failure to do so might result in glass cracking.

WARNING ⚠

The glass is required to remain in the closed position prior to moving the frontal of the cabinet. In case this is not fulfilled, the glass will not be properly adjusted with the frame and it might result in glass cracking.

WARNING ⚠

Do not put the hands in the laterals of the frontal support, to avoid any entrapment risk.

4.3. GENERAL WORKING RECOMMENDATIONS

a) We recommend starting up the cabinet around 10 minutes before you begin to work. This will ensure particles are swept from the work zone and from material introduced.

b) Wash your arms, hands, and nails well, using an anti-bacterial soap, before and after work. Staff should avoid touching their mouths and eyes.

c) We recommend wearing lab coats with adjusted cuffs and, in special works, protection gloves (or sleeve cuffs). Both the coat and cuffs should be made from a fabric that minimizes the emission of fibers and particles.

d) The work area inside the cabinet should not be used to store lab equipment. This can produce an accumulation of dust that is unnecessary and hazardous for working under sterile conditions. Only the fewest pieces of equipment should be placed to work correctly and with the minimum number of interruptions for the user (inserting and removing arms from the work zone).

e) Although the work zone is under a clean and sterile laminar down flow, the area around the cabinet might be contaminated. The sources of contamination should be known to avoid them when inserting arms, material, etc.

f) All the material needed for the work should be free of particles and carefully cleaned before it is placed on the work area.

g) Materials such as paper, wood, cardboard, pencils, rubbers, etc. should not be placed in the work area as they shed a lot of particles.

h) In the event of using pipettes, only mechanical aspiration type should be used. Blowing the pipette with the mouth should never be done because it is easy to inhale aerosols caused by suction.

i) When platinum handles are used, it is advisable to use electric incinerators, preferably single-use ones.

j) If the work requires the use of a gas flame from a Bunsen burner or similar, we recommend using a type where the flame is maintained by holding down a button. Please note that the use of a constant flame burner creates significant turbulence. Bear in mind that a flame too large could burn the absolute filters.

k) Shadows and turbulence provoked by objects, equipment and materials placed in the work area should be studied before beginning to work to evaluate their possible effects. Bear in mind that the laminarity of the airflow is not recovered until 2.5 times the diameter of the object that causes the obstruction.

l) When using bottles and tubes, they should preferably have a screw top instead of cotton stoppers which easily shed a lot of particles.

m) When the work is finished, all disposable products (handles, Petri dishes, etc.) as well as means of cultivation, samples, tubes, bottles, etc. should be removed from the cabinet in impermeable bags suitable for sterilization if necessary.

n) When working in the cabinet and during cleaning operations, it is important to protect the HEPA filters from any damage by avoiding knocking them, projecting liquids, or splashing, etc.

o) When placing or taking out the arm rest, take care not to strongly hit the lateral glasses, to prevent any damage of them.

p) Any equipment connected to the internal sockets inside the cabinet should meet all the electrical specifications.

q) People traffic generates air movements that may cause unexpected loss of the containment barrier and, in the case of Class II cabinets, a reduced level of product protection. Warning signs and other means should be used to avoid traffic during critical operations.

 Cabinet locations can be divided to reduce the effects of air movement produced by other cabinets and people traffic.

r) Suddenly doors opening and closing may disturb airflows in the cabinet, and it should be avoided when cabinets are being used.

WARNING 

The UV light system can only be used when the front glass is closed. Never connect equipment whose consumption exceeds the maximum one indicated on the internal sockets of the cabinets.

4.3.1. ERGONOMICS

Microbiological security cabinet Bio II Advance Plus and its accessories have been designed following recommendations and specifications of the standard EN ISO 14738:2010.

Users should consider their position when working with the cabinet, for their own safety.

Below are examples of four positions **NOT** recommended and one that **IS** recommended.


Fig. 1 X


Fig. 2 X


Fig. 3 X


Fig. 4 X

Adjusting the seat and adopting a correct position enables the user to work safely and to avoid any possible physical damage (sight, neck, back, etc.)


Fig. 5 ✓

In addition, to protect both the product and the user, it is necessary that the cabinet is placed in a location without air currents, which may cause turbulence in the laminar down flow.


Fig. 6 X

4.4. DEFINITION OF THE WORKING AREA

To prevent from creating any turbulence in the laminar down flow and to protect the user and the product, it is necessary to establish a working area inside the cabinet.

- Do not work in any area with perforations and do not cover them up.
- Always work 10 cm from the limits of the walls and perforated area.


Fig. 7

4.5. CONTROL KEYBOARD-DISPLAY

The following elements are found in the control keyboard:

A. *On/Off Button*  *(**Annex 3 pos. 1**):*

- Pressing when the cabinet is turned off:
 It goes to the *Password Request* screen. When entering the correct password (**Chapter 4.5 point 1**), it goes to the *Main Operation* screen.

- Pressing when the cabinet is turned on:
 It goes to the *Password Request* screen. When entering the correct password (**Chapter 4.5 point 1**), the cabinet switches off.

B. *Navigation buttons*  *and*  *(**Annex 3 pos. 2**):*

They move the cursor through all different function icons and select them. The icon pointed inverts its colors. If the selected item is a numeric field, the navigation buttons increase/decrease and validate the value.


C. *Backlight graphic display (**Annex 3 pos. 3**):*

It shows alarm messages, operating graphs, function icons and all the information required in cabinet-user communication. The background color changes to:

- Green: Working mode. All the flows are working fine and there are no alarms activated.
- Red: Alarm mode. Any alarm has been activated or the cabinet is in unsafe mode. In case of alarm, an acoustic alarm is activated together with the red light.
- Yellow: Disinfection mode.
- White blue: Information mode. It is used in Technical Service screens.

Below is an explanation of the different system screens:

1. ***Password Request*** screen *(by default **0000**)*

 Increase the value of the selected field.


 Decrease the value of the selected field.

 Confirm the value of the selected field and move the selection to the next field. After pressing for the final time:

- Right Password: it accesses the *Main Operation* screen.
- Wrong Password: Cursor returns to the first field.

2. ***Main Operation*** screen.

It shows the direct access icons to the most common functions for the user.

Activate/Deactivate fans.

Activate/Deactivate lights.

Press for 5 seconds the button, to access the light intensity screen.

Activate/Deactivate sockets.

Activate/Deactivate* gas solenoid valve (if assembled).


Activate/Deactivate ECO mode: fans in low-velocity, lights switched off and acoustic alarm can be muted.

Go into Disinfection Menu.

Go to Configuration Menu.

***REMARK:** *It is not possible to operate gas tap solenoid valve under unsafe conditions.*

2.1. **Fan Deactivation Confirmation** screen.


Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field.

2.2. **Light intensity** screen.


Increase the value.

Decrease the value.

Confirm the value.

2.3. **ECO Mode Confirmation** screen.


Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field.

2.4. **Disinfection Confirmation** screen.


Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field. After pin confirmation it goes to Disinfection Menu screen.

2.4.1. ***Disinfection Menu*** screen.


Go to *U.V. Menu* screen.

Go to *Vaporized hydrogen peroxide* screen.

Go to *Particles counter* screen (when applicable).

Back to *Main Operation* screen.

2.4.1.1. ***U.V. Menu*** screen.


Begin disinfection, if glass is closed, for the time determined by the timer.

Go to *U.V. Time Adjust* screen.

Back to *Disinfection* Menu screen.

2.4.1.1.1. ***U.V. Lamp*** screen.

Switch off the U.V. lamp and access to *Disinfection* Menu screen.

2.4.1.1.2. ***U.V. Time adjust*** screen.


Confirm the selected value and save it. *

Reset the numeric value.

Back to *U.V. Menu* screen.

***** *Value 000 means infinite time: start and stop will be made manually.*

2.4.1.2. **Vaporized hydrogen peroxide** screen.


ENT Start Vaporized hydrogen peroxide cycle.

Activate/Deactivate the sockets of the work zone.

Back to *Disinfection* Menu screen.

2.4.1.2.1. **Vaporized hydrogen peroxide cycle** screen.


Switch on manually the fans.

Switch off the fans.

Back to Disinfection menu screen (it requires confirmation).

3. ***Configuration Menu*** screen.

It shows the direct access icons to the configuration functions allowed to the user.

Go to Contrast screen.

Go to Maintenance screen.

Go to Language screen.


Go to PIN Configuration screen.

Go to Date & Hour screen.

Go to Sensors Information screen.

Go to Main Operation screen.

3.1. ***Contrast*** screen.


Increase the value.

Decrease the value.

Confirms the value and returns to *Configuration Menu* screen.

3.2. ***Maintenance Menu*** screen*.*

It shows information about the running hours of cabinet, fans, U.V. lamp and filters, together with the last Service visit. Besides it allows the following:


Record a new revision.

Get maintenance recommendations.

Show the historical of revisions and alarms.

Back to *Configuration Menu* screen.

3.3. ***Language Selection*** screen*.*


Go up.

Go down.

Confirms the selected language and goes back to *Configuration Menu* screen.

3.4. ***PIN Configuration*** screen*.*


Enable/disable the password request to Connect/Disconnect the cabinet.

Enable the user to change the user password.

Coded password. In case of forgotten password, by giving this code to Syntegon Telstar the previous password will be recovered.

Back to *Configuration Menu* screen.

3.5. **Date-hour** screen*.*


Saves the value.
(format dd/mm/yy hh:mm).

Cancel and/or edit the entered value.

Auto Power-on menu.

Back to *Configuration Menu* screen.


3.5.1. **Auto Power on** screen.

Cabinet Auto Power on.

UV light Auto Power on.

Back to *Date-hour* screen.

3.5.1.1. **Cabinet Auto Power on** screen.


Saves the value.
(format dd/mm/yy hh:mm).

Cancel and/or edit the entered value.

Back to *Auto Power on* screen.

3.5.1.2. **UV light Auto Power on** screen.


Saves the value.
(format dd/mm/yy hh:mm).

Cancel and/or edit the entered value.

Back to *Auto Power on* screen.

3.6. **Sensors Information** screen.

It shows information about the values of the sensors of the cabinet in real time.

Back to *Configuration Menu* screen.

D. Microprocessor:

The microprocessor keeps total control over the exhaust flow and impulsion flow via the thermoanemometer probe.
The microprocessor control system consists of:

- Backlight graphic display.
- Alarms shown on the display:

ALARM	SHOWN ON THE DISPLAY
Excessive laminar down flow	HIGH LAM. FLOW
Insufficient laminar down flow	LOW LAM. FLOW
Excessive inflow	HIGH INFLOW
Insufficient inflow	LOW INFLOW
Front open / Glass out of working position	WORK POSITION
Redundant channel alarm	REDUNDANT FAIL.
Power supply failure	POWER FAILURE
Memory fail	MEMORY FAILURE

- Messages shown on the display:

MESSAGE	SHOWN ON THE DISPLAY
Fans at reduced speed	ECO MODE

E. Auxiliary microprocessor or redundant circuit.

Should the main microprocessor fail, the auxiliary microprocessor will activate the optical alarm, and also the acoustic alarm. In such a case, slide the glass to the closed position and cover the cabinet exhaust. Disconnect the equipment and contact the Syntegon Telstar Technical Service.

4.6. DISCONNECTION

To stop the cabinet, it is necessary to push the On/Off button (**Annex A3 pos. 1**) to access the *Password Request* screen (**Chapter 4.5 point 1**). After entering the correct password, the cabinet will disconnect.

4.7. PROLONGED STOPPAGE

In the event of a prolonged stoppage, avoid environmental dust to go inside the cabinet by turning the cabinet off with the glass in the closed position. When it is switched on again, proceed like the first time it was switched on. See **Chapters 3.3** and **4.1**.

4.8. OPTIONAL ACCESSORIES

- Support table (two models for 770 mm and 900 mm working area height).
- Extension legs for support stands (*)
- Height adjustable support stand (electrically operated) (except for ducted cabinets) (*)
- IT version: articulated VESA bracket (*)
- Gas tap kit. (*)
- Vacuum tap kit.
- Compressed air tap kit.
- Nitrogen tap kit.
- Oxygen tap kit.
- Water tap kit.
- Carbon dioxide (CO_2) tap kit.

(*) *GS has not been tested for marked options. It must be tested after installation.*

WARNING 

Gas is a toxic and flammable substance which in the event of an accident can cause serious damage to the user and the cabinet.

The presence of a flame within the cabinet creates turbulence which affects its control and proper functioning; therefore, we recommend you do **NOT** use burners unless strictly essential and if the user is prepared to assume this responsibility.

Gas connection should be done in line with the safety standards of each country, region, and company.

For these reasons, the electrically controlled valve that enables gas entry can only be activated when the fans are working, and it is automatically disconnected when the alarm indicating a lack of impulsion flow is activated.

- Vaporized Hydrogen Peroxide decontamination kit.
- One-piece stainless-steel working surface. (*)
- Exhaust charcoal filter (except for ducted cabinets) (**Annex A4.2**) (*)
- Doble exhaust HEPA filter (except for ducted cabinets) (**Annex A4.3**) (*)
- Ducted cabinets items (canopy -**Annex A4.4**-, regulation kit, remote EC fan and flexible hoses) (*)

(*) *GS has not been tested for marked options. It must be tested after installation.*

5. MAINTENANCE

WARNING ⚠

All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

Before starting any maintenance work or dismantling the cabinet, the area should be wiped down with vaporized hydrogen peroxide or formaldehyde (the latter not recommended) (**Annexes A5.3 and A5.4**). Once this has been done, disconnect the cabinet from the power socket.

In general, a regular control of the proper working of all cabinet components should be done once a year.

5.1. MAINTENANCE TABLE

The frequencies determined for each operation, and in particular about changing filters, which appear on the Maintenance Table below are indicative values recommended for normal cabinet service. Very severe work conditions may reduce the time intervals shown on the Table.

TEST *(See Annex A9)*	***REGULARITY***
Downflow filters integrity Test	*Minimum once a year or after filter changes*
Exhaust filters integrity Test	
Downflow air velocity test	
Set point alarm downflow air velocity test	
Inflow air velocity Test	
Set point alarm inflow air velocity test	
Smoke test	
Light test	*Recommended during the startup of the cabinet*
Noise level Test	

5.2. ABSOLUTE FILTERS

The average life of HEPA filters in these cabinets is 4000/5000 hours.

After this time, the propelled airflow falls because of the increased loss of load of the filter as it is a mechanical filtering, although this does not mean that filtering efficacy is reduced. When it reaches the minimum flow measured by the thermoanemometric probe, the microprocessor activates both visual and acoustic alarms.

WARNING 

Do not replace absolute filters or directly handle fans in the microbiological safety cabinets without having first wiped the area down with vaporized hydrogen peroxide or formaldehyde (the latter not recommended). Opening the cabinet without taking the proper precautions can endanger your health. We recommend it is carried out by staff especially trained for this work.

When using formaldehyde, the laws in force in each country must be followed.

To ensure correct installation and adjustment of the filter, as well as ensuring there are no leaks, filters should be replaced by specialized technical staff equipped with the corresponding monitoring instruments: anemometer, leak detector, aerosol generator, particle counter, etc.

The used clogged absolute filters from safety cabinets should be destroyed. This process is always the responsibility of the user. Our technical service does not have any responsibility for possibly contaminated parts.

5.3. SPARE PARTS

Only use original spare parts from the manufacturer.

Bio II Advance Plus 3

No.	Description
1	DOWNFLOW / EXHAUST FAN
2	ABSOLUTE FILTER HEPA 915X457X69 LPD
3	ABSOLUTE FILTER HEPA 457X457X69 LPD
4	LED BAR Y-3014B-144-600-T 24V
5	UV LAMP
6	CPU CARD

Bio II Advance Plus 4

No.	Description
7	DOWNFLOW / EXHAUST FAN
8	ABSOLUTE FILTER HEPA 1220X457X69 LPD
9	ABSOLUTE FILTER HEPA 610X457X69 LPD
10	LED BAR Y-3014B-216-900-T 24V
11	UV LAMP
12	CPU CARD

Bio II Advance Plus 6

No.	Description
13	DOWNFLOW / EXHAUST FAN
14	ABSOLUTE FILTER HEPA 1830X457X69 LPD
15	ABSOLUTE FILTER HEPA 915X457X69 LPD
16	LED BAR Y-3014B-360-1500-T 24V
17	UV LAMP
18	CPU CARD

6. CLEANING AND DISINFECTION

6.1. SURFACE CLEANING AND DISINFECTION

See **Annex A5.1**

6.2. FRONT GLASS CLEANING AND DISINFECTION

See **Annex A5.2**

6.3. DISINFECTION WITH FORMALDEHYDE

Total decontamination of the cabinet could be done using disinfecting gas emission, which ensures penetration in sealed and inaccessible areas of the Cabinet. It also guarantees the sterilization of the working area, fans, plenum, HEPA filters and evacuation ducts.

Decontamination is usually done by fumigation with formaldehyde vapor. However, this may not be suitable in certain cases, such as when it is necessary to follow specific regulations. Alternatives may be used, e.g., hydrogen peroxide.

WARNING 

Disinfection with formaldehyde is not recommended due to it is known to be a human carcinogen and it should only be done in line with the safety standards of each country, region, and company.

The formaldehyde gas is irritant and toxic when inhaled. When solutions are used for surface disinfection, the vapors given out from the surface can injure the respiratory mucosa. This can be nullified using ammonia vapor after disinfection has been completed. It is also a skin irritant.

See **Annex A5.3** for the formalin disinfection procedure.

6.4. DISINFECTION WITH VAPORIZED HYDROGEN PEROXIDE

It is recommended to do decontamination of the cabinet by using vaporized hydrogen peroxide emission, which ensures penetration in sealed and inaccessible areas of the Cabinet. It also guarantees the sterilization of the working area, fans, plenum, HEPA filters and evacuation ducts.

Vaporized hydrogen peroxide systems provide rapid, safe, low-temperature decontamination methods for any enclosed area that may be contaminated with micro-organisms, including spore-forming bacteria.

See **Annex A5.4** for the disinfection with vaporized hydrogen peroxide procedure.

7. MAINTENANCE TABLE

See **Annex A6**

8. ELECTRIC DIAGRAMS

See **Annexes A7 (230V) & A8 (120V)**

9. TESTS

9.1. DOWNFLOW FILTERS INTEGRITY TEST

See **Annex A9.1**

9.2. EXHAUST FILTERS INTEGRITY TEST

See **Annex A9.2**

9.3. DOWNFLOW AIR VELOCITY TEST

See **Annex A9.3**

9.3.1. SET POINT ALARM DOWNFLOW AIR VELOCITY TEST

See **Annex A9.3.1**

9.4. INFLOW AIR VELOCITY TEST

See **Annex A9.4**

9.4.1. SET POINT ALARM INFLOW AIR VELOCITY TEST

See **Annex A9.4.1**

9.5. SMOKE TEST

See **Annex A9.5**

10. CERTIFICATES

10.1. DECLARATIONS OF CONFORMITY

See **Annex A10.1**

11. TROUBLESHOOTING

No.	Type of error	Causes of error	Correction of error
1	**Insufficient exhaust flow** Alarm "LOW INFLOW" (see Chapter 4.5 point D)	Exhaust partially covered.	Visually check the exhaust to remove any type of obstruction.
		Front perforation in the working area totally/partially covered.	Visually check the perforation and remove any type of obstruction.
		Distance between cabinet exhaust and lab ceiling too small.	Increase space between cabinet exhaust and the ceiling. In some cases, it is possible to modify or eliminate the false ceiling. Install a lower support table. This means the working area will be also lower. Duct the cabinet (**Chapter 3.5**).
		Clogged exhaust filter.	Check the date of the last filter change and the hours it has been under operation. (**Chapter 4.5 point 3.2**). Notify the technical service to change the filter.
		Probe not working properly.	Check the probe position is correct.
		Ducted cabinet and ducting too long.	Fit an additional fan at the end of the duct (**Chapter 3.5**).
2	**Insufficient laminar down flow** Alarm "LOW LAM. FLOW" (see Chapter 4.5 point D)	Front and rear perforation in the work zone totally/partially covered.	Visually check the perforation and remove any type of obstruction.
		Clogged laminar down flow filter.	Check the date of the last filter change and the hours it has been under operation. (**Chapter 4.5 point 3.2**). Notify the technical service to change the filter.
		Probe not working properly.	Check the probe position is correct (sensor parallel to the flow). Check and remove any obstruction to the probe.

No.	Type of error	Causes of error	Correction of error
3	**Excessive laminar down flow** Alarm "HIGH LAM. FLOW" (see Chapter 4.5 point D)	Probe not working properly.	Check the probe position is correct (sensor parallel to the flow).
		Air currents in front aperture.	Check and remove any airflow in front of the cabinet that distorts the reading of the sensor.
4	**Laminar down flow not working properly** Alarm "LOW LAM. FLOW" (see Chapter 4.5 point D)	Probe not working properly.	Check the probe position is correct (sensor parallel to the flow). Check and remove any obstruction to the probe. Check that the connection in the control card is correct (**Annex A7/A8**).
		Laminar down flow fan stopped.	Check the fuse has not blown (**Annex A7/A8**). Check that the connection in the control card is correct (**Annex A7/A8**). Notify the Technical Service to change the fan.
5	**Exhaust not working properly** Alarm "LOW INFLOW" (see Chapter 4.5 point D)	Exhaust completely covered up.	Visually check the exhaust to remove any type of obstruction.
		Probe not working properly	Check the probe position is correct (sensor parallel to the flow). Check that the connection in the control card is correct (**Annex A7/A8**).
		Exhaust fan stopped.	Check the fuse has not blown (**Annex A7/A8**). Check that the connection in the control card is correct (**Annex A7/A8**). Notify the Technical Service to change the fan.
6	**Front does not move**	Front is blocked.	Visually check the front movement to remove any type of obstruction.
7	**Glass out of position** Alarm "WORK POSITION" (see Chapter 4.5 pto. D)	Glass position	Assemble the glass in working position.
		Working position micro is loose	Open the frontal and check that the working position micros are not loose. Notify the Technical Service to receive instructions for adjusting the micros.
		Main control is damaged	Call a specialized technician to replace the control.

No.	Type of error	Causes of error	Correction of error
8	**Open frontal** Alarm "WORK POSITION" (see Chapter 4.5 pto. D)	Frontal contact is damaged	Slide down the glass to closed position and check that the frontal contacts are not loose and remove any type of obstruction. Notify the Technical Service to receive instructions for adjusting the contacts.
		Main control is damaged	Call a specialized technician to replace the control.
9	**Light does not turn on**	Light not working properly	Check that the connection in the control card is correct (**Annex A7/A8**).
		Led bar is not working	Check that the lamp is properly connected. Change the led bar.
10	**UV germicidal lamp does not turn on (optional accessory)**	Glass not in closed position	Check that the front slides properly.
		UV lamp does not work	Check that the lamp holder is connected properly. Check that the connection in the control card is correct (**Annex A7/A8**). Change the UV lamp.
		Reactance is damaged	Change the reactance.
11	**Intermittent alarm activated**	Front position security micro activated	Open the front and check the security micro is not loose. Call the Technical Service to receive instructions on adjusting micros.
		Main control is damaged	Call a specialized technician to replace the control.
12	**Display does not turn on**	Power cord cable disconnected	Check that the power cord cable is properly connected in cabinet and in socket on the wall.
		No power	Check voltage on the socket of the wall.
		Main control is damaged	Call a specialized technician to replace the control.

WARNING 

All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

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EN

1. Centrifugal Blower (x4)
2. Exhaust flexible tube
3. Laminar flow chamber or plenum
4. Laminar flow HEPA filter
5. Exhaust HEPA filter
6. Working area
7. Outside air inlet and protection barrier
8. Exhaust outlet

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

EN This version is supplied as a factory built (it includes a fan).

1. Connecting pipe
2. Support
3. Active charcoal filter box
4. Safety Cabinet

GS has not been tested for Active Carbon Filters.

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
EN

Double HEPA filter in the exhaust for extra security when the unit is expelling air into the laboratory (according to standard BS 5726).

GS has not been tested for the double HEPA filter.

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

EN

Installation
A. Install the canopy (1) by using the tubes (2) and the screws (3), access cover (4) facing the front of the unit.
B. Connect the duct to the canopy outlet (Ø200mm) with a hose clamp.

Use
The access cover (4) allows access to the unit exhaust outlet for exhaust airflow checking and for outlet closure for disinfection process.

GS has not been tested for the canopy. It must be tested after installation.

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EN

For worktable and sides, use sterile or non-use disables cloths, which do not shed particles or fibres.

For disinfection, these cloths will be used, slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or affect the stainless steel and the glass. Cleaning can also be carried out with water and soap prior to disinfection.
To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabinet; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

WARNING

It is important keep the glass in working position and the cabinet in working mode when performing the cleaning and disinfection of the working area of the cabinet. In case that it is needed to open the frontal in order to clean areas with difficult access, all needed cautions must be taken because the user protection decreases when the frontal is open.
It is recommended NOT to use disinfectant products with chlorinated derivatives or containing halogens or their compounds, because they can damage the stainless steel and painted surfaces.
Should some bleach be poured, the surface must be always thoroughly rinsed and completely dried with a spatula (glass cleaner type).

It is required to dry completely all surfaces always after cleaning or disinfecting in order to preserve the original qualities of the surfaces.


Avoid spray any product over the frontal mechanism area, in the lateral ends of the frontal glass.

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EN

To accede to corners or difficult access areas, another type of tools will be used as shown en la image.

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EN

The perforated stainless-steel worktable is divided into three parts to facilitate cleaning of the interior container and provide for sterilisation of these parts into a vessel.
Remember that every segment or shelve must be cleaned and/or disinfected both in its lower and upper sides. The parts will be cleaned inside the cabinet while this is running.

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EN

Also be sure of carrying out the cleaning/disinfection in the lower part of the working area (surface and corners).

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EN

Inside, safety cabinet cleaning is recommended on the following occasions:

a) Before starting any work in the safety cabinet.
b) After working in the safety cabinet.
c) Whenever there is a change of work program.
d) In event of liquid spilling on the worktable.
e) Before carrying out mechanical or biological control tests in the working area.

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EN

For the cleaning of front glass, you must use sterile textures or disposable ones without particle or fibres loosing.

For disinfection, use the previous textiles slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or alter the glass. Cleaning can also be carried out with water and soap prior to disinfection.

To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabin; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

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EN

For the glass exterior part cleaning, it is recommended to move the glass into closed position. In this position a major part of the glass is accessible. For the interior part of the glass cleaning, it is recommended to leave the glass in its working position in order to assure the user's protection.

WARNING 

It is important to clean and disinfect the cabinet's working area, with the glass in working position and the cabinet in performance. If it is necessary to separate the frontal-glass mechanism to clean the areas of difficult access, bring the glass to its closed position to set free the frontal mechanism. By separating the frontal-glass, the inner surface of the glass is totally accessible. There must be taken the necessary precautions because the user's protection is lower with the glass separated from cabinet.

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WARNING 

Do not put the hands in the laterals of the frontal support, in order to avoid any entrapment risk.

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
EN

WARNING ⚠

To separate the frontal-glass mechanism (after the glass has been moved into the lower position), pull always only the handles "B" of the metal support, *never using the handles of the glass, that are specifically designed for the vertical movement of the glass. Failure to do so might result in glass cracking.*

FR	Annexe A5.2. NETTOYAGE ET DÉSINFECTION VITRE FRONTALE	GE	Anhang A5.2. REINIGUNG UND DESINFEKTION DER FRONTSCHEIBE

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EN

Gas decontamination of a cabinet that has handled biologically pathogenic products is necessary in the following cases:

- *Before undertaking any safety cabinet maintenance work.*
- *Before removal of HEPA filters for replacement.*
- *Before carrying out mechanical or biological control tests in the safety cabinet.*
- *Before moving the safety cabinet from one work area to another.*
- *Before changing a work program.*
- *Whenever liquid is spilled on the worktable.*

The gas used in this procedure for total decontamination is formaldehyde, made by depolymerisation of paraformaldehyde (solid) or formol (liquid).

The quantity of formaldehyde to be used is calculated according to Safety Cabinet total volume (see table 1).

WARNING ⚠

Disinfection with formaldehyde is not recommended due to it is known to be a human carcinogen and it should only be done in line with the safety standards of each country, region and company.

	Bio II Advance Plus 3	**Bio II Advance Plus 4**	**Bio II Advance Plus 6**
Volumen total (m^3) ***Total volume (m^3)***	0,75	1,00	1,50
Formaldehido (ml) ***Formaldehyde (ml)***	45	60	90

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EN

Sufficient formaldehyde should be used to produce a persistent airborne concentration at least 50 mg/m^3.

Formaldehyde vapour can be generated by:

- *Evaporation of formalin (a volume fraction of ± 36% in water).*
 NOTE: 60 ml of formalin plus 60 ml water should be vaporized per cubic meter of cabinet volume.
- *Depolymerization of paraformaldehyde by heating, provided adequate humidity is available.*

WARNING ⚠

Formaldehyde steam is explosive when the percentage in volume of the mixture with dry air reaches 7.75%; its flash point is 430ºC. The national limits on exposure in the workplace should be taken into account.

The necessary material to perform decontamination by gas is as follows:

- *Formaldehyde.*
- *A formolizer/neutralizer unit capable to perform automatically both formolization and neutralization in sequence.*
- *A hermetic plastic cover (PVC or equivalent)*
- *PVC sheet or PVC laminated aluminum sheet*
- *Duct tape*

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EN

WARNING ⚠

Formaldehyde penetrates with difficulty and his efficiency depends on the temperature and humidity. It is more effective over 20ºC and with relative humidity of 65%. The use of excessive quantities can cause polymer deposits inside the cabinet and can contribute to the filter blockade.

It is recommended that workers use proper protection to make decontamination operation, in order to protect them from an accidental formaldehyde gas exposure. Mainly, they must cover the eyes with leak tight safety goggles and, when necessary, use respiratory protection with safety hood.

Formaldehyde steams are very irritant for the eyes, throat, and nose.

Nevertheless, performing decontamination work with enough precautions means that the operator will not face any contact with the gas.

Steps to follow:

a) Refer to local regulation concerning the use of formaldehyde.
b) Fumigation must be performed by a qualified person with a good knowledge of the procedure and the required precautions. At any case, a visible warning should indicate the cabinet is being fumigated.
c) Measure the needed quantity of formaldehyde to be used and put it in the formolizer (see table 1).

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EN

d) The safety cabinet should be sealed before fumigation to ensure that no formaldehyde can leak into the laboratory or the other rooms. For that, cover the safety cabinet with the plastic cover. This plastic cover is completely airtight and has been tested against leaks with the ammonia test. It consists of testing the plastic cover at a pressure of 80 Pa with an ammonia atmosphere inside. Externally, an ammonia detector checks no leakage. Check also that the PVC sheet or PVC laminated aluminum sheet is continuous (one single piece without any joint)

e) Place the PVC sheet or PVC laminated aluminum sheet below the table support of the cabinet (by raising first one side and then the other side)

f) Switch on the cabinet, disable the socket inlets and switch off the fans

g) Separate the frontal of the cabinet and raise the glass to the maximum position

h) Plug in the vaporized with formol into the working surface while keeping the plug switched off by the display control

i) Cover totally the cabinet using the hermetic PVC bag and letting go out only the cabinet power cord by the lower part of the bag

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EN

j) Seal the cabling together with the PVC bag and the PVC sheet, with duct tape.
Sticking process details:
 a) Duct tape with Omega shape covering the cabling 360º (fig. 1 & 2).
 b) It is necessary to turn several times the cabling to guarantee a good sealing all around the 360º (fig. 3)
 c) Raise the cabling by the free exterior part (fig. 4) and put another piece of duct tape (in red in fig. 5)
 d) Stick this piece of duct tape against the former piece of duct tape with the Omega shape (fig. 6) and then turn several times the cabling to guarantee a good sealing (fig. 7)

Fig 1


Fig 2


Fig 3


Fig 4


Fig 5


Fig 6


Fig 7

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EN

k) Complete seal the plastic cover, against the PVC sheet base, with duct tape, taking special care about avoiding any folding that may lead to some possible leak (fig. 10-11-12-13)

l) Cover with a complete second round of duct tape the joint between the plastic bag and the PVC sheet base (fig. 14 & 15)

Fig 10

Fig 11

Fig 12

Fig 13

Fig 14

Fig 15

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EN

m) Go to the decontamination/VHP menu and switch on the sockets by means of the display control

n) Press ENT to begin the formolization procedure. Switch manually the fans after 10 minutes for 15 seconds, once the formol is partially vaporized

o) Once the formol is fully vaporized (about 15 minutes), switch again the fans for 15 seconds

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EN

p) For a suitable fumigation, it should be ensured that the downstream side of the HEPA filters and the interior of the cabinet are exposed to formaldehyde enough time to ensure inactivation of microorganisms that have penetrated into the filter. This requires an exposure time of over 6 hours. Fumigation is most conveniently done overnight.

q) Even after the treatment, the HEPA filters can only be considered to be 'safe to handle when using the appropriate protective clothing' instead of no sterile. They should be treated in an autoclave or incinerated on removal. Following the fumigation, the seals of the cabinet should be removed, and the extractor fan switched on for a suitable period of time to eliminate the residual formaldehyde steam before using the cabinet again or carrying out any maintenance task.

r) For cabinets that recycle air into the room, the fumigation system used should ensure that no member of personnel is exposed to higher levels than those specified by the relevant authorities. For example, by using a suitable neutralisation system or temporary ducts to extract the exhaust fumes into the exterior. Suitable protective respiratory equipment designed for being used with formaldehyde should be available during fumigation in case of an emergency.

s) Once the process has finished, switch off the equipment in the same way, through the control panel.

FR	Annexe A5.3. DÉSINFECTION AU FORMOL	GE	Anhang A5.3. DESINFEKTION MIT FORMALDEHYD

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The gas used in this procedure is the vaporized hydrogen peroxide.

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- *1 vaporized hydrogen peroxide decontamination kit accessory.*
- *1 Plastic bag 1800x600 mm.*
- *1 Roll of packaging tape.*
- *1 Roll of adhesive tape.*
- *1 Allen key.*
- *Vaporization hydrogen peroxide system (like STERIS VHP® 100 or STERIS VHP® X-10).*
- *Vaporized hydrogen peroxide sensor.*

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

EN

Steps to follow:

a) Fumigation must be performed by a qualified person with a good knowledge of the procedure and the required precautions. At any case, a visible warning should indicate the cabinet is being disinfected.

b) *Install the vaporized hydrogen peroxide accessory kit in the side glass (input port), and the exhaust cover on top (output port). Seal the exhaust cover with tape.*

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EN

c) Separate the frontal of the cabinet.

d) To guarantee the tightness of the equipment, seal the cabinet front opening by means of one plastic bag and packaging tape, covering any gap that could leak hydrogen peroxide to the ambient.

e) Close the front window and place the glass in working position.

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

EN

f) Connect the vaporization hydrogen peroxide system hoses to the proper port. Check that input of vapor is fixed in the side glass, and output is fixed on the top.

g) Go to Vaporized hydrogen peroxide cycle *screen. Screen will be red.*

h) Cabinet is ready to start with decontamination cycle.

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EN

Decontamination process.

It must be performed by specialized personnel of the vaporization hydrogen peroxide system.

After cycle is finished, high concentration may remain inside the cabinet. Once the plastic sealing the front and the exhaust cover are removed, concentration will be diluted with room air and be reduced rapidly to accepted levels.

ATTENTION

*During process, the hydrogen peroxide concentration in the room must be monitored for safety reasons, so it will be necessary to place a sensor for that purpose. Ambient limit is 1 ppm*****, and at that concentration 8 h of continuous work is allowed. With ambient concentration of 5 ppm***** 2h of continuous work is allowed.*

**** 0.5 ppm in Germany***

a) Even after the treatment, the HEPA filters can only be considered to be 'safe to handle when using the appropriate protective clothing' instead of no sterile. They should be treated in an autoclave or incinerated on removal.
b) Suitable protective respiratory equipment designed for being used with hydrogen peroxide should be available in case of an emergency.
c) Once the process has finished, switch off the equipment through the control panel.

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PROJECT LEADER:
Victor Lazaro

AUTOMATION LEADER:
X. Cañellas

DRAWN BY:	DATE:
X.Cañellas	24/08/2015
APPROVED:	DATE:
X. Cañellas	24/08/2015
REVIEWED:	DATE:
X. Cañellas	24/08/2015
AS BUILT:	DATE:

REVISIONS:
ID: Z DATE: 30/01/2024 DRAWN BY: X.Cañellas
ID: ZA DATE: 08/08/2024 DRAWN BY: X.Cañellas
ID: W DATE: 15/01/2024 DRAWN BY: X.Cañellas
ID: X DATE: 25/01/2024 DRAWN BY: X.Cañellas
ID: Y DATE: 25/01/2024 DRAWN BY: X.Cañellas

Telstar®

Date 24/08/2015 Report of Drawings 1609430 Page 1
BIO-II-ADVANCE PLUS CABINET

AZBIL TELSTAR
TECHNOLOGIES, S.L.


Drawing	Description	Creation	Modification
01	SUMMARY	24/08/2015	01/12/2020
02	POWER	24/08/2015	01/12/2020
03	CONTROL BOARD	24/08/2015	01/12/2020
04	BLOWERS	24/08/2015	15/04/2020
05	LIGHTING AND UV KIT	24/08/2015	01/12/2020
06	SOCKETS & GAS VALVE	24/08/2015	08/08/2024
07	WINDOW DETECTORS	24/08/2015	28/05/2020
08	AIR FLOW SENSORS	24/08/2015	30/01/2024
09	DISPLAY AND BUTTONS	24/08/2015	09/11/2020
10	ACCESSORIES BOARD (OPTIONAL PARTS)	11/03/2020	01/12/2020
11	BILL OF MATERIALS	24/08/2015	08/08/2024

PROPERTY OF TELSTAR
REPRODUCTION FORBIDDEN

REPLACES:

REPLACED BY:

PROJECT NAME:
1609430

ELECTRICAL CHARACTERISTICS			
MODEL	ADVANCE 6	ADVANCE 4	ADVANCE 3
VOLTAGE	230VAC	230VAC	230VAC
POWER	2.0 KVA	1.65 KVA	1.62 KVA
CURRENT	8.7 A	7.16 A	7.04 A
FREQUENCY	50..60 HZ	50..60 HZ	50..60 HZ

Revision	Date	Draw by	Date:	24/08/2015
T	09/11/2020	X. Cañellas	Draw by:	X. Cañellas
M	28/05/2019	X.Cañellas	Aproved:	X. Cañellas
R	06/08/2020	X.Cañellas		

Project description: BIO-II-ADVANCE PLUS

Situation:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Zone:

Drawing description: POWER

Drawings: 2 / 11

Previous: 1


Drawing number: 1609430

Next: 3

Project reference: 977068


Drawing number: 2

632679
-PCB1
ATLAS
CONTROL BOARD


OPTIONAL

627589
-PCB2
ATLAS
ACCESSORIES BOARD

CONTROL BOARD			
CONNECTORS			FUSES
J1 - DISPLAY	J10 - SOCKETS 230VAC 4A	J19 - not used	F1 - 250VAC 8A/5x20R CONTROL BOARD
J2 - not used	J11 - P. SUPPLY EXTERNAL EXH. FANS	J20 - not used	F2 - 250VAC 10A/5x20R DOWNFLOW FANS
J3 - ELECTROVALVE	J12 - EXTERNAL EXH. FAN CONT. SIGNAL	J21 - WORK POSITION W.SENSOR	F3 - 250VAC 2.5A/5x20R EXHAUST FAN
J4 - LED LIGHTING	J13 - FANS CONTROL SIGNAL	J22 - WORK POSITION W.SENSOR REDN.	F4 - 250VAC 5A/5x20R SOCKETS
J5 - EXTERNAL EXH. AIR FLOW SENSOR	J14 - POWER SUPPLY FANS	J23 - not used	F5 - 250VAC 1.6A/5x20R LIGHT
J6 - AIR FLOW SENSOR	J15 - MAIN POWER SUPPLY 230/110VAC	J24 - not used	F6 - 250VAC 1A/5x20R UV LIGHT
J7 - ACCESSORY BOARD CONNECTION	J16 - not used	J25 - CLOSE POSITION W.SENSOR	
J8 - UV LIGHT	J17 - 24VDC POWER SUPPLY	J26 - CLOSE POSITION W.SENSOR REDN.	
J9 - not used	J18 - not used		

Revision	Date	Draw by	Date:	24/08/2015	Project description: BIO-II-ADVANCE PLUS	Situation:	Zone:	Drawings: 3 / 11	Previous: 2	Next: 4
Q	02/06/2020	X. Cañellas	Draw by:	X. Cañellas		Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: CONTROL BOARD		Drawing number: 1609430	
G	26/01/2017	M.Díaz	Aproved:	X. Cañellas				Project reference: 977068	Drawing number 3	
N	11/03/2020	X.Cañellas								


AIRFLOW FANS


6 FEET CABINET


3 & 4 FEET CABINET

Revision	Date	Draw by	Date:	24/08/2015	Project description: BIO-II-ADVANCE PLUS	Situation:	Zone:	Drawings: 4 / 11	Previous: 3	Next: 5
K	16/05/2019	X. Cañellas	Draw by:	X. Cañellas		Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: BLOWERS		Drawing number: 1609430	
O	15/04/2020	X. Cañellas	Aproved:	X. Cañellas						
H	30/01/2017	M.Díaz						Project reference: 977068	Drawing number 4	

Telstar®

J4
CIF03005
SH4
0V4
VS4
CABLE: 1720466
SHIELD
BLACK
RED
RED
BLACK
xx
-LED1
24VDC
**
CABINET
CODE LED1
WATTS
ADVANCE 3
622654
20 W
ADVANCE 4
621658
23 W
ADVANCE 6
621659
38 W
LED STRIP
J8
CIF03003
L8
N8
BROWN
BLUE
CABLES: 1637601
1335
-L1
230VAC
30W
L8
2018
-CB1
-LT1
1341
UV LIGHT
N8
Revision
Date
Draw by
Date:
Draw by:
Aproved:
Q
02/06/2020
X. Cañellas
J
30/05/2018
X. Cañellas
O
15/04/2020
X. Cañellas
2408/2015
X. Cañellas
X. Cañellas
Project description:
BIO-II-ADVANCE PLUS
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
LIGHTING AND UV KIT
Drawings:
5 / 11
Previous:
4
Next:
6
Project reference:
977068
Drawing number:
1609430
Drawing number
5
Telstar

Revision
Date
Draw by
Date:
Draw by:
Aproved:
B
11/05/2016
X. Cañellas
U
21/12/2022
X. Cañellas
ZA
08/08/2024
X. Cañellas
24/08/2015
X. Cañellas
X. Cañellas
Project description:
BIO-II-ADVANCE PLUS
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
SOCKETS AND GAS VALVE
Drawings: 6 / 11
Previous: 5
Next: 7
Project reference:
977068
Drawing number:
1609430
Drawing number
6
Telstar
J10
CIF03003
GR/YEL
BLUE
BROWN
PE
N10
L10
CABLES:
ADVANCE 3: 1637612
ADVANCE 4: 1637611
ADVANCE 6: 1637613
SOCKETS
637516
-SK1
230VAC
2A
637516
-SK2
230VAC
2A
OPTIONAL
J3
CIF02005
CABLES: 1637628
RED
BLACK
VS3
0V3
GAS VALVE
-SV1
640621
-SV1
24VDC

Revision	Date	Draw by	Date:	24/08/2015
B	11/05/2016	X. Cañellas	Draw by:	X. Cañellas
C	23/05/2016	X. Cañellas	Aproved:	X. Cañellas
P	28/05/2020	X. Cañellas		

Project description: BIO-II-ADVANCE PLUS

Situation:

Zone:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Drawing description: WINDOW DETECTORS

Drawings: 7 / 11

Previous: 6


Next: 8

Drawing number: 1609430

Project reference: 977068

Drawing number: 7

Telstar®

AIR
VELOCITY
PROBE

Revision	Date	Draw by	Date:	24/08/2015
D	31/05/2016	X. Cañellas	Draw by:	X. Cañellas
V	15/02/2023	X. Cañellas	Aproved:	X. Cañellas
Z	30/01/2024	X. Cañellas		

Project description: BIO-II-ADVANCE PLUS

Situation:

Zone:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Drawing description: AIR FLOW SENSORS

Drawings: 8 / 11

Previous: 7


Next: 9

Drawing number: 1609430

Project reference: 977068


Drawing number 8

Telstar®


					Situation:	Zone:	Drawings: 9 / 11	Previous: 8	Next: 10		
Revision	Date	Draw by	Date:	24/08/2015	Project description: BIO-II-ADVANCE PLUS	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: DISPLAY		Drawing number: 1609430		
O	15/04/2020	X. Cañellas	Draw by:	X. Cañellas							
R	02/06/2020	M.Díaz	Aproved:	X. Cañellas				Project reference: 977068	Drawing number 9		
S	09/11/2020	X.Cañellas									

OPTIONAL


RS485 COMMPORT

J28
15EDGK-3.81-4


VHP CONTROL

J31
15EDGK-3.81-5


BUTTERFLY VALVE

J33
15EDGK-3.81-3

POSITION SIGNAL


J32
2EDGK-7.5-3P

RELAY ON/OFF CONTROL
(MAX. 230VAC 10A)


BOOSTER FAN

J34
15EDGK-3.81-3

VELOCITY CONTROL


J35
15EDGK-3.81-3

RELAY ON/OFF CONTROL
(MAX. 230VAC 10A)

ALARM SIGNAL


J36
15EDGK-3.81-3

RELAY ON/OFF ALARM
(MAX. 230VAC 10A)


TEMP./HUMID./PRESS. SENSOR

J37
15EDGK-3.81-4


PARTICLE COUNTER SENSOR

J38
15EDGK-3.81-4

EXTERNAL AIRFLOW SENSOR

J39
15EDGK-3.81-4

ACCESSORIES BOARD	
CONNECTORS	
J39 - EXTERNAL AIRFLOW SENSOR	J33 - BUTTERFLY VALVE INPUT SIGNAL
J38 - PARTICLE COUNTER SENSOR	J32 - BUTTERFLY VALVE RELAY
J37 - TEMPERATURE/PRESSURE/HUMIDITY	J36 - ALARM RELAY SIGNAL
J28 - RS485 COMMPORT	J35 - BOOSTER FAN POWER RELAY
J40 - TO CONTROL BOARD CONNECTION	J34 - BOOSTER FAN 0-10V SIGNAL
J31 - VHP	

ACCESSORIES
(OPTIONALS PARTS)

Revision	Date	Draw by	Date:		Project description:
N	11/03/2020	X.Cañellas	Draw by:	X. Cañellas	BIO-II-ADVANCE PLUS
T	09/11/2020	X.Cañellas	Aproved:	X. Cañellas	

Situation:
Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description: ACCESSORIES BOARD OPTIONAL PARTS
Drawings: 10 / 11
Previous: 9
Next: 11
Drawing number: 1609430
Project reference: 977068
Drawing number: 10

Telstar

Date 20/12*2022 Bill of Materials Page 1
BIO-II-ADVANCE PLUS

TELSTAR
TECHNOLOGIES, S.L.

SYmbol	Drawing	Column	Description	N. Refer	Catalogue	Reference
-FL1	02	D	EMI FILTER YANBIXINKEJI	YB12EA-20A-Q	ESPE	634174
-CN1	02	D	CONNECTOR C14+2 FUSES	0717-2-CQ	ESPE	619568
-PS1	02	G	Power Supply MEANWELL	RS-75-24	ESPE	621116
-PCB1	03	C	Control board ATLAS	-	ESPE	632679
-PCB2	03	F	ATLAS ACCESSORIES EXPANSION CARD	---	---	627589
-FAN1	04	B	Fan with capacitor FANS-TECH	DH146E5-AGT-00	ESPE	607949
-FAN2	04	C	Fan with capacitor FANS-TECH	DH146E5-AGT-00	ESPE	607949
-FAN3	04	C	Fan with capacitor FANS-TECH	DH146E5-AGT-00	ESPE	607949
-FAN4	04	D	Fan with capacitor FANS-TECH	DH146E5-AGT-00	ESPE	607949
-LED1	05	E	BIO-II-ADVANCE 3 LED LIGHT 24VDC	Y-3014B-144-600-O	ESPE	622654
	05	E	BIO-II-ADVANCE 4 LED LIGHT 24VDC	Y-3014B-216-900-T	ESPE	622658
	05	E	BIO-II-ADVANCE 6 LED LIGHT 24VDC	Y-3014B-360-1500-T	ESPE	622659
-LT1	05	G	UV Lamp 30W Philips	-	ESPE	1341
-L1	05	F	Reactance 30W 220V	FL22-28-30-20-18 CFL (W)	ESPE	1335
-CB1	05	G	Starter 4-65W	-	ESPE	2018
-BE2	06	E	Socket plug YUADON	YD-001S	ESPE	637516
-BE3	06	F	Socket plug YUADON	YD-001S	ESPE	637516
-SV1	06	F	SOLENOID VALVE BURKERT	266293	ESPE	640621
-S2	07	C	Microswitch OMRON	V-155-1A5	ESPE	71692
-S1	07	F	Microswitch OMRON	V-155-1A5	ESPE	71692
-S4	07	F	Microswitch OMRON	D2F-L2-D	ESPE	71429
-S3	07	F	Microswitch OMRON	D2F-L2-D	ESPE	71429
-SP1	08	C	Airflow probe Degree Controls	FS22500	ESPE	614060
-DIS1	09	C	DISPLAY PARA BSC ATLAS	-	ESPE	632700

Revision	Date	Draw by	Date:	15/04/2020
Z	15/03/2023	X. Cañellas	Draw by:	X. Cañellas
ZA	08/08/2024	X.Cañellas	Aproved:	X. Cañellas
Y	25/01/2024	X.Cañellas		

Project description: BIO-II-ADVANCE PLUS
Situation:
Zone:
Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.
Drawing description: BILL OF MATERIAL
Drawings: 11 / 10
Previous: 10
Next:
Drawing number: 1609430
Project reference: 977068
Drawing number 11

EQUIPMENT:

BIO-II-ADVANCE PLUS 120VAC

PROJECT REFERENCE:

977068

CUSTOMER:

AZBIL TELSTAR TECHNOLOGIES, S.L.

DRAWING TYPE:

ELECTRIC

DRAWINGS:

11

MATERIAL O.F.:

PROJECT LEADER:

Victor Lazaro

AUTOMATION LEADER:

X. Cañellas

DRAWN BY:	DATE:
X.Cañellas	22/08/2017
APPROVED:	DATE:
X. Cañellas	22/08/2017
REVIEWED:	DATE:
X. Cañellas	22/08/2017
AS BUILT:	DATE:

REVISIONS:

ID: F	DATE: 05/12/2022	DRAWN BY:	X. Cañellas
ID: G	DATE: 19/01/2023	DRAWN BY:	X. Cañellas
ID: H	DATE: 16/052023	DRAWN BY:	X. Cañellas
ID: I	DATE: 25/01/2024	DRAWN BY:	X. Cañellas
ID: J	DATE: 29/04/2024	DRAWN BY:	X. Cañellas

Telstar®

Date 22/08/2017 Report of Drawings 1708351 Page 1
BIO-II-ADVANCE PLUS CABINET

AZBIL TELSTAR
TECHNOLOGIES, S.L.


Drawing	Description	Creation	Modification
01	SUMMARY	22/08/2017	--
02	POWER	22/08/2017	29/04/2024
03	CONTROL BOARD	22/08/2017	29/04/2024
04	BLOWERS	22/08/2017	08/05/2020
05	LIGHTING AND UV KIT	22/08/2017	02/06/2020
06	SOCKETS & GAS VALVE	22/08/2017	29/04/2024
07	WINDOW DETECTORS	22/08/2017	28/05/2020
08	AIR FLOW SENSORS	22/08/2017	-
09	DISPLAY AND BUTTONS	22/08/2017	06/08/2020
10	ACCESSORIES BOARD (OPTIONAL PARTS)	08/05/2020	29/04/2024
11	BILL OF MATERIALS	22/08/2017	29/04/2024

PROPERTY OF TELSTAR
REPRODUCTION FORBIDDEN

REPLACES:

REPLACED BY:

PROJECT NAME: 1708351

ELECTRICAL CHARACTERISTICS				
MODEL	ADVANCE 6	ADVANCE 4	ADVANCE 3	
VOLTAGE	120 VAC	120 VAC	120 VAC	
POWER	1.44 KVA	1.1 KVA	1.1 KVA	
CURRENT	12 A	8.5 A	8.5 A	
FREQUENCY	50..60 HZ	50..60 HZ	50..60 HZ	

Revision	Date	Draw by	Date:	22/08/2017
G	19/01/2023	X. Cañellas	Draw by:	X. Cañellas
H	16/05/2023	X. Cañellas	Aproved:	X. Cañellas
J	29/04/2024	X. Cañellas		

Project description: BIO-II-ADVANCE PLUS

Situation:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Zone:

Drawing description: POWER

Drawings: 2 / 11

Previous: 1

Drawing number: 1708351


Next: 3

Project reference: 977068

Drawing number: 2


Telstar

632679
-PCB1
ATLAS
CONTROL BOARD


OPTIONAL


646403
-PCB2
ATLAS
ACCESSORIES BOARD

CONTROL BOARD			
CONNECTORS			FUSES
J1 - DISPLAY	J10 - SOCKETS 230VAC 4A	J19 - not used	F1 - 250VAC 8A/5x20R CONTROL BOARD
J2 - not used	J11 - P. SUPPLY EXTERNAL EXH. FANS	J20 - not used	F2 - 250VAC 10A/5x20R DOWNFLOW FANS
J3 - ELECTROVALVE	J12 - EXTERNAL EXH. FAN CONT. SIGNAL	J21 - WORK POSITION W.SENSOR	F3 - 250VAC 2.5A/5x20R EXHAUST FAN
J4 - LED LIGHTING	J13 - FAN CONTROL SIGNAL	J22 - WORK POSITION W.SENSOR REDN.	F4 - 250VAC 5A/5x20R SOCKETS
J5 - EXTERNAL EXH. AIR FLOW SENSOR	J14 - POWER SUPPLY FANS	J23 - not used	F5 - 250VAC 1.6A/5x20R LIGHT
J6 - AIR FLOW SENSOR	J15 - MAIN POWER SUPPLY 230/110VAC	J24 - not used	F6 - 250VAC 1A/5x20R UV LIGHT
J7 - ACCESSORY BOARD CONNECTION	J16 - not used	J25 - CLOSE POSITION W.SENSOR	
J8 - UV LIGHT	J17 - 24VDC POWER SUPPLY	J26 - CLOSE POSITION W.SENSOR REDN.	
J9 - not used	J18 - not used		

Revision	Date	Draw by	Date:	22/08/2017	Project description:	Situation:	Zone:	Drawings: 3 / 11	Previous: 2	Next: 4
A	08/05/2020	X. Cañellas	Draw by:	X. Cañellas	BIO-II-ADVANCE PLUS	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: CONTROL BOARD		Drawing number: 1708351	
C	02/06/2020	X. Cañellas	Aproved:	X. Cañellas						
J	29/04/2024	X. Cañellas						Project reference: 977068	Drawing number 3	

AIRFLOW FANS
6 FEET CABINET
3 & 4 FEET CABINET
J14
CIF03003
J13
CTF0300T
CABLES: ADVANCE 3-4: 1636571
ADVANCE 6: 1636573
CABLES: ADVANCE 3-4: 1636672
ADVANCE 6: 1636673
-FAN1
-FAN2
-FAN3
-FAN4
632026
135W
0.9A
yellow
blue
L14
N14
PE
SH13
0V13
VS13
BRDVN
BLUE
Shield
Tach
0-10V
GND
Revision
A
Date
08/05/2020
Draw by
X. Cañellas
Date:
22/08/2017
Draw by:
X. Cañellas
Aproved:
X. Cañellas
Project description:
BIO-II-ADVANCE PLUS
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
BLOWERS
Drawings:
4 / 11
Previous:
3
Next:
5
Project reference:
977068
Drawing number
1708351
Drawing number
4
Telstar

Revision	Date	Draw by	Date:	22/08/2017
A	08/05/2020	X. Cañellas	Draw by:	X. Cañellas
C	02/06/2020	X. Cañellas	Aproved:	X. Cañellas

Project description: BIO-II-ADVANCE PLUS

Situation:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Zone:

Drawing description: LIGHTING AND UV KIT

Drawings: 5 / 11

Previous: 4


Next: 6


Project reference: 977068

Drawing number: 1708351

Drawing number 5

Telstar®

Revision
I
J
Date
25/01/2024
29/04/2024
Draw by
X. Cañellas
X. Cañellas
Date:
Draw by:
Aproved:
22/08/2017
X. Cañellas
X. Cañellas
Project description:
BIO-II-ADVANCE PLUS
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
SOCKETS AND GAS VALVE
Drawings: 6 / 11
Previous: 5
Next: 7
Project reference:
977068
Drawing number:
1708351
Drawing number
6
J10
CIF03003
GR/YEL
BLUE
BROWN
PE
N10
L10
CABLES:
ADVANCE 3: 1637612
ADVANCE 4: 1637611
ADVANCE 6: 1637613
SOCKETS
606918
-SK1
125VAC
2A
606918
-SK2
125VAC
2A
OPTIONAL
J3
CIF02005
CABLES: 1637628
RED
BLACK
VS3
0V3
GAS VALVE
640621
-SV1
24VDC
-SV1

Revision	Date	Draw by	Date:	22/08/2017
A	08/05/2020	X. Cañellas	Draw by:	X. Cañellas
B	28/05/2020	X. Cañellas	Aproved:	X. Cañellas

Project description: BIO-II-ADVANCE PLUS

Situation:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Zone:

Drawing description: WINDOW DETECTORS

Drawings: 7 / 11

Previous: 6


Drawing number: 1708351

Next: 8

Project reference: 977068


Drawing number: 7

Telstar

AIR
VELOCITY
PROBE

Revision	Date	Draw by	Date:	22/08/2017	Project description:	Situation:	Zone:	Drawings: 8 / 11	Previous: 7	Next: 9
E	01/12/2020	X.Cañellas	Draw by:	X. Cañellas	BIO-II-ADVANCE PLUS	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: AIR FLOW SENSORS		Drawing number: 1708351	
			Aproved:	X. Cañellas				Project reference: 977068	Drawing number 8	

Revision	Date	Draw by	Date:	22/08/2017
A	08/05/2020	X. Cañellas	Draw by:	X. Cañellas
D	06/08/2020	X. Cañellas	Aproved:	X. Cañellas

Project description: BIO-II-ADVANCE PLUS

Situation:

Zone:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Drawing description: DISPLAY

Drawings: 9 / 10

Previous: 8

Next: 11


Drawing number: 1708351

Project reference: 977068

Drawing number: 9


Telstar

OPTIONAL


RS485 COMMPORT

J28
15EDGK-3.81-4


VHP CONTROL

J31
15EDGK-3.81-5


BUTTERFLY VALVE

J33
15EDGK-3.81-3


J32
2EDGK-7.5-3P


BOOSTER FAN

J34
15EDGK-3.81-3


J35
15EDGK-3.81-3


ALARM SIGNAL

J36
15EDGK-3.81-3


TEMP./HUMID./PRESS. SENSOR

J37
15EDGK-3.81-4


PARTICLE COUNTER SENSOR

J38
15EDGK-3.81-4

EXTERNAL AIRFLOW SENSOR

J39
15EDGK-3.81-4

ACCESSORIES BOARD	
CONNECTORS	
J39 - EXTERNAL AIRFLOW SENSOR	J33 - BUTTERFLY VALVE INPUT SIGNAL
J38 - PARTICLE COUNTER SENSOR	J32 - BUTTERFLY VALVE RELAY
J37 - TEMPERATURE/PRESSURE/HUMIDITY	J36 - ALARM RELAY SIGNAL
J28 - RS485 COMMPORT	J35 - BOOSTER FAN POWER RELAY
J40 - TO CONTROL BOARD CONNECTION	J34 - BOOSTER FAN 0-10V SIGNAL
J31 - VHP	

ACCESSORIES
(OPTIONALS PARTS)

Revision	Date	Draw by	Date:	08/05/2020	Project description: BIO-II-ADVANCE PLUS	Situation:	Zone:	Drawings: 10 / 11	Previous: 9	Next: 11
A	08/05/2020	X. Cañellas	Draw by:	X. Cañellas		Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: ACCESSORIES BOARD OPTIONAL PARTS		Drawing number: 1708351	
E	01/12/2020	X. Cañellas	Aproved:	X. Cañellas				Project reference: 977068	Drawing number 10	
J	29/04/2024	X. Cañellas								

Date 22/08/2017 Bill of Materials Page 1

BIO-II-ADVANCE PLUS

TELSTAR
TECHNOLOGIES, S.L.

SYmbol	Drawing	Column	Description	N. Refer	Catalogue	Reference
-FL1	02	D	EMI FILTER YANBIXINKEJI	YB12EA-20A-Q	---	634174
-F1	02	D	FUSIBLE 15A 6.3x32 LANBAO ELECT	RTI-30	---	655770
-CN1	02	D	CONNECTOR IEC C20 250VAC 20A Helishi	R-305(00)	---	611678
-CN2	02	D	CONNECTOR C14+2 FUSES	0717-2-CQ	---	619568
-PS1	02	G	Power Supply MEANWELL	RS-75-24	---	621116
-TC1	03	C	Control board ATLAS	-	---	632679
-FAN1	04	B	Fan with capacitor FANS-TECH 120VAC	DH146E5-ADT-00	---	632026
-FAN2	04	C	Fan with capacitor FANS-TECH 120VAC	DH146E5-ADT-00	---	632026
-FAN3	04	C	Fan with capacitor FANS-TECH 120VAC	DH146E5-ADT-00	---	632026
-FAN4	04	D	Fan with capacitor FANS-TECH 120VAC	DH146E5-ADT-00	---	632026
-LED1	05	C	BIO-II-ADVANCE 3 LED LIGHT 24VDC	Y-3014B-144-600-O	---	622654
			BIO-II-ADVANCE 4 LED LIGHT 24VDC	Y-3014B-216-900-T	---	621658
			BIO-II-ADVANCE 6 LED LIGHT 24VDC	Y-3014B-360-1500-T	---	621659
-CB1	05	F	UV Lamp 30W Philips	-	---	1341
-BL1	05	F	BALLAST T8 120-277V	ACT2P32ISUNVE	---	637554
-BE2	06	E	SOCKET USA WATERPROOF 125 15A	-	---	606918
-BE3	06	F	SOCKET USA WATERPROOF 125 15A	-	---	606918
-SV1	06	F	SOLENOID VALVE BURKERT	266293	---	640621
-S2	07	C	Microswitch OMRON	V-155-1A5	---	71692
-S1	07	F	Microswitch OMRON	V-155-1A5	---	71692
-S4	07	F	Microswitch OMRON	D2F-L2-D	---	71429
-S3	07	F	Microswitch OMRON	D2F-L2-D	---	71429
-SP1	08	C	Airflow probe Degree Control	FS22500	---	614060
-DIS1	09	C	DISPLAY PARA BSC ATLAS	-	---	632700

Revision	Date	Draw by	Date:	22/08/2017
I	25/01/2024	X. Cañellas	Draw by:	X. Cañellas
J	29/04/2024	X. Cañellas	Aproved:	X. Cañellas
H	16/05/2023	X. Cañellas		

Project description: BIO-II-ADVANCE PLUS

Situation: | Zone:

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Drawing description: BILL OF MATERIAL

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EN *1. PURPOSE*

Check the integrity of the downflow filters.

2. SCOPE

Bio II Advance Plus

3. TOOLS

- *Aerosol generator (cool or hot).*
- ATI PAO-4, *Shell Ondina fluid or Emery 3004 oils.*
- *Photometer (ATI 2H).*

4. METHODOLOGY

*Introduce fume by the generator in the working area (Figure 1). Concentration in cabinet u*sually within 20-70µg/l.

- *Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Figure 2)*


Fig. 1

Fig. 2

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EN

- *With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters, including the perimeter of each filter pack, by passing the photometer probe in slightly overlapping strokes over the entire surface.*

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EN

- *Also scan the entire periphery of the filter and the junction between filter and filter mounting frame at a scanning rate that does not exceed 5 cm/s.*

Fig. 3

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EN ***5. ACCEPTANCE CRITERIA***

- *The photometer signal after the filter should be smaller or equal to a 0,01% of the specific concentration before the filter.*
- *The concentration value around the DOP port should be smaller or equal to the concentration of environment.*

WARNING

Positive pressure DOP connection*. Proceed to connect and disconnect the ports* ***only*** *when the cabinet is* ***switched off.***

Fig. 4. Resultado correcto: no fuga
Correct result: filter not leaking

Fig. 5. Resultado incorrecto: fuga
Wrong result: filter leaking

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EN Annex A9.2. EXHAUST FILTERS INTEGRITY TEST

EN 1. PURPOSE

*Prove the **not** existence of leaks in the HEPA exhaust filters*

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EN

- *With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters, by passing the photometer probe in slightly overlapping strokes over the entire surface.*

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EN

- *Also scan the entire periphery of the filter and the joint between filter and filter mounting frame at a scanning rate that does not exceed 5cm/s.*

Fig. 3


Fig. 4

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EN ***CANOPY***

- *In case of canopy are installed and can't remove it, scan the filter by the exhaust removing the registry cover canopy.*

 GS has not been tested for the canopy. It must be tested after installation.

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$V_{promedio}$ (m/s) = 1/n x (V_1 + V_2 + V_3 + V_n)

EN *1. PURPOSE*

- *Check that the speed control in different points of the working area is suitable to achieve product protection. To determine the speed average as well as the uniformity range in the laminar flow zone.*

2. SCOPE

- *Bio II Advance Plus*

3. TOOLS

- *Omnidirectional or unidirectional anemometer.*

4. METHODOLOGY

- *Inside the considered working area according to each cabinet model divide the working area in squares (see attached table).*

Then switch on the cabinet and wait until the velocity and the temperature.

- *Place the anemometer at 270mm high of the shelves, over each square (take all readings with the anemometer in the same direction and without covering any holes, whenever possible).*


- *Take note of the readings and the position of the anemometer in each point.*

- *For this calculation, we proceed to the following formula:*

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

EN ***5. ACCEPTANCE CRITERIA***

- *The average of the airflow velocity at the work surface should be ≥ 0.34 m/s and ≤ 0.36 m/s*
- *Every single value should not differ more than ± 20% from the average of the airflow velocity.*

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EN *Annex A9.3.1. SET POINT ALARM DOWNFLOW AIR VELOCITY TEST*

EN *1. PURPOSE*

- *Check correct work of the alarms when flow speed is minor or major than the admissible.*

2. SCOPE

- *Bio II Advance Plus*

3. TOOLS

- *Omnidirectional or unidirectional anemometer*

4. METHODOLOGY

- *Accede to Technical Service menu through Password Petition screen (User's Manual chapter. 4.5 Point 1), using for that the random security password that you are asked for.*
- *Move to Inflow Test Set Points screen and wait some seconds until the stabilization of the cabinet.*
- *Increase slowly downflow velocity through function button*  *until "Excessive Laminar Flow" alarm activated. Then take note of the display velocity.*
- Taking the same points as in the previous test *(Annex A9.3)*, measure and take note of air velocity of all points and calculate the average.
- *Decrease slowly downflow velocity through the function button*  *until "Insufficient Laminar Flow" activated. Then take note of the display velocity.*
- Taking the same points as in the previous test *(Annex A9.3)*, measure and take note of air velocity of all points and calculate the average.

EN *5. ACCEPTANCE CRITERIA*

All below criteria should be satisfied:

- (High down flow display velocity) -0.02 ≤ (Measured average down flow velocity at high alarm) ≤ (High down flow display velocity) +0.02
- (Down flow velocity at low alarm) -0.02 ≤ (Measured average down flow velocity at low alarm) ≤ (Displayed down flow velocity at low alarm) +0.02
- No individual measurement (max. or min. point) should differ from its average more than 20% when alarm activated.

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NOTA: Superficie aire de entrada:

Bio II Advance Plus 3: 0,149 m^2
Bio II Advance Plus 4: 0,203 m^2
Bio II Advance Plus 6: 0,313 m^2

EN *1. PURPOSE*

- *The purpose of this test is to calculate the inflow velocity through the front aperture of the equipment to obtain the Class II to the operator protection.*

2. SCOPE

- *Bio II Advance Plus*

3. TOOLS

- *The tool to use is an anemometer.*

4. METHODOLOGY

The following process shall be performed while the equipment is working:

- *Take 5 velocity readings at 2.5 cm of the protective grills over the exhaust filter.*
- *Calculate the average velocity: $V_{average}$ (m/s) = 1/5 x (V_1 + V_2 + V_3 + V_4 + V_5)*
- *Calculate the flow exhaust: $Q_{exhaust}$ (m^3/h) = ($V_{average}$ x Exhaust surface) x 3600*

NOTE: Exhaust surface: 0.078 m^2.

- *Finally, calculate the inflow velocity by the following formula: V_{inflow} (m/s) = $Q_{exhaust}$ / (Front Aperture surface x 3600)*

NOTE: Inflow surfaces:

Bio II Advance Plus 3: 0.149 m^2
Bio II Advance Plus 4: 0.203 m^2
Bio II Advance Plus 6: 0.313 m^2

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 (0,50 m/s ≤ $V_{entrada}$ ≤ 0,60 m/s).

EN *5. ACCEPTANCE CRITERIA*

- *Advance Plus 3 & 4: Test results are considered to be satisfactory when the inflow velocity is suitable to achieve operator protection (0.55 m/s ≤ V_{inflow} ≤ 0.66 m/s).*
- *Advance Plus 6: Test results are considered to be satisfactory when the inflow velocity is suitable to achieve operator protection (0.50 m/s ≤ V_{inflow} ≤ 0.60 m/s).*

Punto / *Point*	1	2	3	4	5
	V (m/s)				
Rejilla Extracción / *Exhaust Grill*					
Velocidad media (m/s) / *Average velocity (m/s)*					
Q extracción (m^3/h) / *Exhaust Q (m^3/h)*					
Velocidad entrada (m/s) / *Inflow Velocity (m/s)*					

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EN

1. PURPOSE

- Check the correct operation of inflow air velocity alarm when flow is minor than the admissible.

2. SCOPE

- *Bio II Advance Plus*

3. METHODOLOGY

- Accede to Technical Service menu through Password *Petition* screen (User's Manual Chapter 4.5 Point 1), using for it random security password that you ask for.
- Move to *Test Set Points* screen and wait some seconds until the stabilization of the cabinet
- Decrease slowly entry velocity through the function button  until "Insufficient Extraction flow" alarm is activated. Take note:
 - Display entry velocity.

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EN ***1. OBJECT***

- *Verify the cabinet for personnel and product protection.*

2. SCOPE

- *Bio II Advance Plus*

3. TOOLS

- *Smoke generator*

4. METHODOLOGY

- *Product protection: smoke is produced 5-10 cm from the protective glass (Fig. 1 and 2). The behaviour of the smoke is checked visually.*
- *Smoke is produced inside the cabinet, at about 5-10 cm over the lower edge of the glass (Fig. 3). The behaviour of the smoke is checked visually.*
- *User protection: smoke is produced inside the cabinet, as close as possible to the lower edge of the glass (Fig. 4). The behaviour of the smoke is checked visually.*


Fig. 1


Fig. 2

Fig. 3

Fig. 4

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EN ***5. ACCEPTATION CRITERIA***

Product protection:

- *The generated smoke in the outside of the front aperture will be dragged to the inflow area. No fume sign must reach the working surface.*
- *Turbulences in the working area shall not be created. The smoke generated should be dragged by the laminar flow to the surface of the working area.*

User protection:

- *No smoke shall be refluxed out of the cabinet once drawn in, guaranteeing the protection of the operator and environment.*

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UK DECLARATION OF CONFORMITY	
Manufacturer:	SYNTEGON TELSTAR TECHNOLOGIES, S.L.U.
Manufacturer Address:	Av. Font i Sagué, 55 Parc Científic i Tecnològic Orbital 40 08227 Terrassa (Barcelona), Spain
N.I.F.:	B-63797559
Equipment Description:	Biological Safety Cabinet **BIO II ADVANCE PLUS**
UK DIRECTIVES	**DECLARATION**
Supply of Machinery (Safety) Regulations 2008 ***(S.I. 2008/1597)*** *Electromagnetic Compatibility Regulations 2016* ***(S.I. 2016/1091)*** *Electrical Equipment (Safety) Regulations 2016* ***(S.I. 2016/1101)*** *The Restriction of the use of Certain Hazardous substances in Electrical and Electronic Equipment Regulations 2012* ***(S.I. 2012/3032)*** Specific standards: *Safety requirements for electrical equipment for measurement, control, and laboratory use* ***(BS EN 61010:2010+A1)*** *Performance criteria for microbiological safety cabinets* ***(BS EN 12469:2000)***	The present declaration of conformity is issued under the sole responsibility of Manufacturer. This document hereby certifies that the above equipment complies with the UK Directives listed, and their amendments. The equipment must be installed, tested, and operated in strict accordance with Manufacturer instructions and used by fully trained operators. Only approved Syntegon Telstar Technologies S.L.U. spares must be used for maintenance and repair and fitted in accordance with Manufacturer instructions. UK CA
Body Retaining Technical File:	Retained by Syntegon Telstar Technologies, S.L.U.
Signed by:	
Name:	Jordi Valls
Title:	Global Quality Manager
On Behalf of:	SYNTEGON TELSTAR TECHNOLOGIES, S.L.U.
Date:	January 2025

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