UV Lamp Installation Guide

RephiLe Catalogue Number: **2650110C1**
Description: UV Lamp, 254nm, 1/pk
Corresponding Millipore Part Numbers: **ZFRES00UV**

This product is a direct replacement of Millipore corresponding UV lamp. Please refer to Millipore's User Manual for the detailed instruction of installation.

Introduction

Where used: ASM (1997 Range)
In the Shipping Box: 1 UV Lamp with cable; 1 pair of gloves.

Disposing of the Used Lamp

Attention!

The used lamp contains metallic mercury and should therefore be disposed of by the owner in agreement with local regulations.

It is very important to keep the UV lamp perfectly vertical at all times that it is being moved into or out of the housing. This will avoid the possibility of breaking the UV lamp.

Warning!

The light from an UV lamp is very damaging to your eyes and to your skin. Do not turn on the Water System or attempt to turn on the UV lamp in any way until the UV lamp is inside its Quartz sleeve housing inside the Water System.

Supplier declines all responsibilities for such incidents linked to non compliance with this warning.

Procedure: A brief description of the installation process is described below:

Step	Action
1	Switch off the control box power.
2	Remove the protective cap on the reservoir.
3	Wear the supplied gloves so fingerprints do not get on the UV lamp.
4	Disconnect the lamp connector. Carefully remove the lamp from its quartz sleeve. Fit the new lamp into the quartz sleeve by holding it at the top. Connect the cable to the ASM system.
5	Verify that the protective cap is correctly fitted to avoid any UV radiation to the exterior.

Certificate of Quality

This product is manufactured in a RephiLe facility. This plant's Quality Management System is certified to the ISO9001:2008 Quality Systems Standard by an accredited registering body.

RephiLe certifies that the product described below meets RephiLe product quality release criteria.

UV Lamp
Catalogue number: **2650110C1**
Description: UV Lamp, 254nm, 1/pk
Lot Number: S5PAC5601
Mfg Date: 24-Mar-2015



Quality Assurance Manager

RephiLe is a trademark of RephiLe Bioscience, Ltd.
© 2010. RephiLe Bioscience, Ltd. All rights reserved.
All other trademarks are trademarks of their respective manufacturer's.
Rev. 03/10