USER MANUAL
MANUAL DE USUARIO
MANUEL D'UTILISATION
BEDIENUNGSANLEITUNG

Microbiological Safety Cabinets
Cabinas de Seguridad Microbiológica
Cabines de sécurité microbiologique
Mikrobiologische Sicherheitswerkbänke

BiOptima Exa

Telstar®
Life Science solutions

Green line

www.telstar-lifesciences.com

USER'S MANUAL

Microbiological Safety Cabinet

BiOptima Exa

Failure to read, understand and follow the instructions in this manual may result in damage to the unit, injury to operating personnel and poor equipment performance.

www.telstar-lifesciences.com

SUMMARY

0. SYMBOLS AND ABBREVIATIONS 4

1. GENERAL DATA 5

2. DESCRIPTION OF THE SAFETY CABINETS 6
2.1. Working principle 6
2.2. Applications 6
2.3. Components description 7
2.4. Standard equipment 7
2.5. Certification 7
2.6. Technical features 8

3. INSTALLING INSTRUCTIONS 9
3.1. Transport 9
3.2. Placing on site 9
3.3. Precautions for first connection 10
3.4. Electrical connection and services 11
3.5. Cabinet exhaust ducting 12
3.5.1. Additional fan connection 13
3.6. Installation tests 15

4. OPERATING INSTRUCTIONS 16
4.1. Start-up 16
4.2. Front panel positions 17
4.3. General working recommendations 19
4.3.1. Ergonomics 20
4.4. Definition of the working area 23
4.5. Control keyboard-display 24
4.6. Disconnection 30
4.7. Prolonged stoppage 30
4.8. Optional accessories 31

5. MAINTENANCE 32
5.1. Maintenance table 32
5.2. Absolute filters 33
5.3. Spare parts 34

6. CLEANING AND DISINFECTION 35
6.1. Surface cleaning and disinfection 35
6.2. Front glass cleaning and disinfection 35
6.3. Disinfection with vaporized hydrogen peroxide 35

7. MAINTENANCE TABLE 36

8. ELECTRICAL DIAGRAMS 36

9. TESTS 36
9.1. Downflow filters integrity test 36
9.3. Inflow filters integrity test 36
9.4. Downflow air velocity test 36
9.4.1. Set point alarm downflow air velocity test 36
9.5. Inflow air velocity test 36
9.5.1. Set point alarm inflow air velocity test 36
9.6. Smoke test 36

10. CERTIFICATES 37
10.1. EU & UKCA declarations of conformity 37

11. TROUBLESHOOTING 38

Translation of the original manual

0. SYMBOLS AND ABBREVIATIONS

List of stickers used in the biological cabinets.

"EARTH WIRE PROTECTED" STICKER.	
"FUNCTIONAL EARTH WIRE" STICKER.	
"BIOHAZARD" STICKER.	
"CYTOTOXIC" STICKER.	
"ATTENTION, ELECTRIC SHOCK RISK" STICKER.	
"SOCKET VOLTAGE" STICKER.	
"SOCKET MAXIMUM CURRENT" STICKER.	

"SAFETY WARNING PRIOR TO SERVICE OR REPAIR ACTIONS" STICKER.	
"DOP PORTS" STICKER.	

Abbreviations used in this manual:

V	Volts
A	Amperes
W, kW, µW	Watts, kilowatts, microwatts
mm, cm, m	Millimetres, centimetres, metres
kg	Kilogram

s, h	Seconds, hours
Hz	Hertz
dB	Decibel
ºC	Celsius
MSC	Microbiological Safety Cabinets

1. GENERAL DATA

Operating Instructions Manual for Microbiological Safety Cabinets Class II **BiOptima Exa** Class II B2 (100% exhaust air) according to NSF49 standards, which ensure (applying correct techniques in laboratory) high level protection for the operator, the environment and the sample. It is designed for tasks that involves volatile, toxic non-flammable and non-explosive chemicals or radionuclides in biological-chemical applications.

WARNING

This Manual must be read and understood before its installation, maintenance and operation of the biological safety cabinet Class II. All the safety instructions given in the Manual must be followed.


Any use or application failing to comply with the recommendations expressed in the Manual may cause danger and loss of rights to manufacturing guarantee as well as freeing the manufacturer from all responsibility arising.

Any maintenance must be done by qualified staff according to the laboratory working correct conditions. In case these actions are not carefully carried out, as a result the safety cabinet can work in abnormal conditions.

All maintenance actions on this equipment must be performed by a qualified technician who is familiar with the proper maintenance procedure required for this equipment.

The European Union issued the Directive 2012/19/EU on Waste of Electrical and Electronic Equipment (WEEE).

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive. Mark symbol is shown at the right.
Azbil Telstar sells products through distributors throughout Europe. Contact your local distributor for recycling/disposal.

The name plate & CE mark sticker is located at one side of the cabinet.

MANUFACTURER: **AZBIL TELSTAR TECHNOLOGIES, S.L.U**
Av. Font i Sagué, 55
Parc Científic i Tecnològic Orbital 40
08227 – Terrassa (Spain)
Tel (+34) 93 736 16 00
Fax (+34) 93 786 13 80
E-mail: Telstar@telstar.com
http://www.telstar.com

2. DESCRIPTION OF THE SAFETY CABINETS

2.1. WORKING PRINCIPLE

BiOptima Exa Biohazard Safety Cabinets are designed for microbiology work with the most advanced microprocessor control, ensuring a high level of protection for the operator, environment and sample.

Exterior of cabinet made of laminated, lacquered, oven-dried steel. Functional design with reduced external dimensions and front access for changing filters.

Spacious internal chamber with sliding-swinging front made from laminated safety glass. Segmented work surface and liquid collection tray made of polished stainless steel. Armrest is detachable.

Principle diagrams in **Annex A2.**

As can be seen in **Annex A2**, the air goes through the inlet in the upper part of the cabinet (4), then the air propelled by the fans (1) is moved into the chamber or plenum (2) and, via the Absolute HEPA filter (3), is filtered and flows according to laminar regime into the working area (5). The frontal air inlet (6) ensures the operator protection. Then, both frontal inflow and laminar flow pass through a first filter step with several HEPA filters (9). Access to these filters is completely safe because they can be replaced meanwhile user protection is working avoiding any risk of contamination.

Afterwards, the air is expelled to the outside without recirculation (7).

The unit is equipped with inlet (4) and aspiration prefilters (8) to protect the interior of the unit against pollution by coarse materials and to collect spilled fluids. Furthermore, the pre-filter will enhance the life span of the HEPA filters significantly.

2.2. APPLICATIONS

The MSC models are appropriate for handling biologically pathogenic samples or products (up to microorganism risk BSL3), ensuring sterile air in the work zone, as well as a high level of protection for the operator and the environment. **BiOptima Exa** may provide additional protection for the user when working with vapours and gases where aseptic conditions are required because total exhaust construction provides greater dilution of any volatile chemicals within the BSC work area, preventing the formation of explosive atmospheres.

WARNING 

Under no circumstances, is possible handling in corrosive, flammable or explosive gas ambient, where traditional fume cupboard must be used.

BiOptima Exa must always be ducted and air must be expelled to the exterior by means of an exhaust fan located in at the end of the duct.

Right working conditions are obtained only when the display is green and without any alarm message on the display.

2.3. COMPONENTS DESCRIPTION

(See Annex A2)

2.4. STANDARD EQUIPMENT

- Two inlet sockets 2A - 460W each located on either side of the work zone. Insulation level IP44.
- UV kit

Easy-to-work UV germicidal lamp with a timer function. It can only be enabled in the *Disinfection* screen (**Chapter 4.5 point 2.4**) or by means of the **UV light Auto Power on** screen (**Chapter 4.5 point 3.5.1.2**). For safety reasons, this button is electrically linked with the light button, and it switches off if the front window is opened.
Average disinfection intensity over the work zone of 200 µW/cm^2 with 20% depreciation after a working life of 8000 hours.

WARNING ⚠

In the event of using this accessory, you should consider the normal precautions in the use of these types of lamps. Their use is only suitable for previous external sterilisation of materials introduced in the cabinet. So, this lamp accessory should be used before the start of work in the cabinet and it is essential to keep it turned off when the operator is in front of the working area or working on it. The glasses provide protection from direct radiation, but UV rays are reflected on polished metal parts and can seriously damage operator's eyes. That is why the glass front must be placed into the closed position to be able to switch on the U.V light, then the fans will automatically disconnect.

2.5. CERTIFICATION

Quality control and tests results certification carried out in line with the applicable standards EN 12469 & NSF49. See Certificate inside the cabinet.

2.6. TECHNICAL FEATURES

Features	Units.	BiOptima Exa 4	BiOptima Exa 6
External dimensions (LxWxH) (Width including armrest)	mm	1338 x 832 x 2212	1948 x 832 x 2212
Internal dimensions (LxWxH)	mm	1190 x 548 x 666	1800 x 548 x 666
Front opening height	mm	200	
Weight	Kg	285	380
Downflow velocity	m/s	0.35	0.35
Inlet velocity / Exhaust flow	m/s m^3/h	0.6 1303	0.6 1971
Voltage *	V	230	
Frequency	Hz	50/60	
Lighting	Lux	≥ 1000	
Vibration	mm RMS	< 0.005	
HEPA/ULPA H 14 Filters		Efficiency 99.995%, EN 1822 by MPPS method filtration system (99.999% D.O.P. test @ 0.3µm particles)	
Fans		High-efficiency fans. Flows measured by thermoanemometric probes.	

* To ensure proper installation, check the attached electrical diagrams.

DISSOCIATED DECLARED VALUES OF NOISE EMISSION Standard(s) to which conformity is declared ISO 4871		
	$V_{downflow}$ = 0.35m/s V_{inlet} = 0.6 m/s (230 V, 50 Hz)	
	BiOptima Exa 4	**BiOptima Exa 6**
L_{pA} (dB)	≤58	≤58
K_{pA} (dB)	2	2
Values measured according to acoustic test in standard EN 12469:2000 using standard ISO 11201		

3. INSTALLING INSTRUCTIONS

For a correct start-up procedure of the MSC, observe the following step-by-step instructions.

3.1. TRANSPORT

The cabinet is supplied in an individual package and in the suitable way to prevent it from being damaged.

This manual and the individual Quality Control certificate of each cabinet are included in the scope of supply.

3.2. PLACING ON SITE

a) Special care should be taken when unpacking the front glass and the control display.

b) Use the same packaging base for the transport to its final location and it is highly recommended doing it with a pallet jack.

WARNING 

Handle the pallet jack carefully to avoid turn over risks due to the high gravity center.
Never open the frontal before the cabinet is fully fastened to its support stand.

WARNING 

We advise against the use of lift trucks for placing the cabinet directly by its base.

WARNING 

For any handling during assembly or dismantling of the cabinet, appropriate safety measures should be taken, such as wearing safety gloves or shoes.

c) Cabinets should be placed away from entrances, corridors, ventilation supply points, etc., which might affect airflow in the work area.

d) Ensure the cabinet is set perfectly on the table or workbench in order to prevent vibrations.

e) To prevent the work area from overheating, we recommend that the room temperature is always within the comfort margins for working.

The optimal comfort margins are:

- Temperature: 18-28ºC
- Relative humidity: 50% ± 20%

If any of these conditions is not fulfilled, optimum cabinet performance cannot be guaranteed.

f) When necessary, air vents in the room should be modified or speed of air entry reduced so that air is not aimed directly at the cabinet, as this may produce turbulence.
Side air vents are not appropriate for rooms where cabinets are installed.

g) The location of the cabinet will directly influence the noise level. This level will rise considerably in small spaces with very reflective surfaces because of the proximity of the

cabinets to the walls, and where fan speed has risen to compensate for air currents in the room or the resistance of clogged filters.

h) Correct cabinet performance should be tested before putting it into operation if any modification has been made to the room or to the air intake system which may considerably change the airflow characteristics in the room.

3.3. PRECAUTIONS FOR FIRST CONNECTION

- Before connecting it to the mains supply, proceed to make a general cleaning to remove any dust particles that may be accumulated during transportation. Follow the instructions on cleaning and disinfection.


- Before using it for the first time, run the cabinet for around 6 hours to purge the filters. New filters give off a characteristic smell which tends to disappear in a short period of time.

WARNING ⚠

The cabinet should work at the voltage shown on the name plate. It should be connected to a suitable socket base (minimum 16A) with an earthed plug. The cable section recommended for the base is 1.5 mm^2.

WARNING ⚠

The plug connection must be located in such a place to prevent it from being unplugged intentionally or accidentally during operation, e.g., next to a wall or closed to another equipment.

3.4. ELECTRICAL CONNECTION AND SERVICES

The electrical characteristics of the cabinet are as follows:

- Electrical class: Class I.
- Insulation level: IP-20.


The cabinet is supplied with a plug cord according to European standard (other plug cords available as optional). The on/off button is in the frontal control keyboard-display. (**Annex A3 pos. 1**).

The power socket in the wall should be located in a perimeter around the cabinet of:

- 1.25 m from the middle of the cabinet to the sides
- 2.5 m from the base of the cabinet upwards

On the top right-hand side of the cabinet are:

- Main fuse (1)
- Power supply (2)

WARNING

Unplug the cabinet from the electric current supply, in case of electrical maintenance.

The cabinet can only be activated from the button in the frontal control keyboard.

In the event of a power failure the cabinet will stop, showing alarm message and providing enough energy to close the front by continuously pressing the additional keyboard key. Once the cabinet is closed, it can be disconnected.

The installation which the cabinet is connected to should meet the safety requirements in force in each country.

The cabinet must be ALWAYS connected to a hard-ducted exhaust system with external fan to evacuate the air outside the laboratory.

During normal work in these safety cabinets no perforated surface in the work zone should be covered with paper or sanitary fabric.

3.5. CABINET EXHAUST DUCTING

The **BiOptima Exa** has no dedicated exhaust fan built in. It needs to be ducted to avoid the toxic vapours, prevent from the formation of explosive atmospheres when handling volatile compounds, assure an aseptic and particle-free working environment, and prevent from lab temperature rise.

WARNING

Ducting the exhaust to the outside of the lab creates a depression in the lab.

WARNING

An additional fan is required to duct the air outside of the laboratory, it must be connected to the CPU Board in order to adjust velocity and send an alarm if there is a breakdown.


The following conditions must be fulfilled in order not to disturb the normal operation of the unit and its performance:

- Exclusively intended hard duct of at least 300 mm diameter.
- Outlet protected from back airflow, closure, entry of any element in the duct which may result in air flow restriction (animals…).
- External fan, recommended extraction capacity greater than:

Model	Extraction flow (m^3/h)
BiOptima Exa 4	1300
BiOptima Exa 6	2000


Pressure drop ΔP = 370 Pa (ΔP of the cabinet with new filters = 185 Pa)

- Hermetic damper or manual valve to isolate the cabinet during the decontamination process.

3.5.1. ADDITIONAL FAN CONNECTION

Inside the upper box of the cabinet, a terminal block is provided for the external fan connection, that should follow this diagram:

PIN	COLOR	FUNCTION
1	Red	0..10V SIGNAL
2	Black	GND
3	White	NC
4	Brown	NO
5	Green	COM

Velocity signal: 0..10V and GND signal to control the fan velocity


FAN ON/OFF: COM and NC pins are a free potential contact used to give a signal to start the fan power. *Maximum values of free potential contacts are 250V 10A*.

Customer facilities installation examples:


<u>Manually regulated fans:</u>

- When the regulation of the additional fan speed must be done manually, the cabinet will only provide the signal to activate or deactivate the fan supply. The speed regulation will be carried out by means of a dimmer or similar regulator.
- For easy fan velocity adjustment, the dimmer should be installed in the same room as the cabinet.

- Manually regulated AC fans:


- EC fan manually controlled:


AC fans regulated from the cabinet:

o When the regulation of the velocity must be carried out from the cabinet, the 0-10V signal can be used to regulate the fan velocity by means of a frequency driver or a solid-state relay, and a free potential signal to turn on or turn off the device (relay or contactor) that provides the power supply to the fan.

EC fan controlled from the cabinet:

o When the regulation of the velocity must be carried out from the cabinet, the 0-10V signal can be used to regulate the fan velocity and a free potential signal to turn on or turn off the device (relay or contactor) that provides the power supply to the fan.

3.6. INSTALLATION TESTS

According to standard DIN EN 12469 and in order to ensure correct cabinet operation, after the installation and connection of the cabinet, it is compulsory to perform the following tests (before the equipment is used for the first time)

- Impulsion filters integrity Test.
- Inflow filters integrity test.
- Downflow air velocity.
- Downflow set points alarm test.
- Inflow air velocity Test.
- Inflow set points alarm test.
- Smoke test.

TELSTAR also recommends making the following tests:

- Light Test.
- Noise Level Test

Telstar, its representatives and distributors put a qualified technical team at your disposition that uses control equipment with maximum health protection for carrying out this installation test under full safety guarantees.

4. OPERATING INSTRUCTIONS

4.1. START-UP

Once the above points have been controlled, start up the MSC by pressing the On/Off button (**Annex A3, pos. 1**). After entering the password correctly (by default the password is **0000**) or immediately if this function has been disabled, the fans will automatically start working.

Ensure that the inlet area, aspiration area, and air exhaust outlet are free (**Annex A2, pos. 4, 6 & 7**).

The electronic system of regulation by microprocessor via two thermoanemometer probe, laminar down flow and exhaust airflow, needs a few seconds to stabilise during connection in order to begin to work normally.


The cabinet can only be used by qualified staff who meets the correct work techniques in the laboratory. If these techniques are not fulfilled, faulty working of the equipment may happen.


WARNING 

The EN 12469 certification does not allow the deactivation of the password enquiry.

4.2. FRONT PANEL POSITIONS

By means of the central buttons, the front glass can be moved.


 : It opens the front glass. Keep it pressed for some seconds to activate movement. Then release the button. When activated below working position, it will automatically stop just above it and go down till reaching the precise working position.

 : It moves down the front glass. Due to security reasons, it is necessary to keep it pressed continuously to allow down direction of the glass. Depending on the initial glass position, it will automatically stop at the precise working position or in the closed position.

Note: Lights will blink during glass movement.

Closed position | Working position | Open position

WARNING

The movement of the glass exerts a strong force. Check that no one has their hands, any part of body or obstacles below the glass when it is moving down, as it can cause injuries. Warn other people working near the cabinet.

In case it is needed to clean the glass by the inner part (Annex A5.2) or if it is necessary to place in the working area any big product or device, the front glass may be separated from the cabinet (after fully lifting it), providing full access to both the inner glass area and the working area. All needed precautions should be taken as the user protection decreases during that operation.

REMARK

In case of power failure, user is able to close the glass thanks to the battery that is installed for that purpose.

WARNING

The glass is required to be in full lifted position before closing the front glass towards the cabinet. In case this is not fulfilled, the glass may be damaged.

4.3. GENERAL WORKING RECOMMENDATIONS

a) We recommend starting up the cabinet around 10 minutes before you begin to work. This will ensure particles are swept from the work zone and from material introduced.

b) Wash your arms, hands and nails well, using an anti-bacterial soap, before and after work. Staff should avoid touching their mouths and eyes.

c) We recommend wearing lab coats with adjusted cuffs and, in special works, protection gloves (or sleeve cuffs). Both the coat and cuffs should be made from a fabric that minimises the emission of fibres and particles.

d) The work area inside the cabinet should not be used to store lab equipment. This can produce an accumulation of dust that is unnecessary and hazardous for working under sterile conditions. Only the fewest pieces of equipment should be placed in order to work correctly and with the minimum number of interruptions for the user (inserting and removing arms from the work zone).

e) Although the work zone is under a clean and sterile laminar down flow, the area around the cabinet might be contaminated. The sources of contamination should be known in order to avoid them when inserting arms, material, etc.

f) All the material needed for the work should be free of particles and carefully cleaned before it is placed on the work area.

g) Materials such as paper, wood, cardboard, pencils, rubbers, etc. should not be placed in the work area as they shed a lot of particles.

h) In the event of using pipettes, only mechanical aspiration type should be used. Blowing the pipette with the mouth should never be done because it is easy to inhale aerosols caused by suction.

i) When platinum loops are used, it is advisable to use electric incinerators, and preferably single-use (disposable ones) loops.

j) If the work requires the use of a gas flame from a Bunsen burner or similar, we recommend using a type where the flame is maintained by holding down a button. Please note that the use of a constant flame burner creates significant turbulence. Bear in mind that a flame too large could burn the absolute filters.

k) Shadows and turbulence provoked by objects, equipment and materials placed in the work area should be studied before beginning to work in order to evaluate their possible effects. Bear in mind that the laminarity of the airflow is not recovered until a distance of 2.5 times the diameter of the object that causes the obstruction.

l) When using bottles and tubes, they should preferably have a screw top instead of cotton stoppers which easily shed a lot of particles.

m) When the work is finished, all disposable products (handles, Petri dishes, etc.) as well as means of cultivation, samples, tubes, bottles, etc. should be removed from the cabinet in impermeable bags suitable for sterilisation if necessary.

n) When working in the cabinet and during cleaning operations, it is important to protect the HEPA filters from any damage by avoiding knocking them, projecting liquids or splashing, etc.

o) Any equipment connected to the internal sockets inside the cabinet should meet all the electrical specifications.

p) People traffic generates air movements that may cause unexpected loss of the containment barrier and, in the case of Class II cabinets, a reduced level of product

protection. Warning signs and other means should be used to avoid traffic during critical operations.

Cabinet locations can be divided to reduce the effects of air movement produced by other cabinets and people traffic.

q) Suddenly doors opening and closing may disturb airflows in the cabinet, and it should be avoided when cabinets are being used.

WARNING ⚠

The UV light system can only be used when the front glass is closed. Never connect equipment whose consumption exceeds the maximum one indicated on the internal sockets of the cabinets.

4.3.1. ERGONOMICS


Microbiological security cabinet BiOptima and its accessories have been designed following recommendations and specifications of the standard EN ISO 14738:2010.

Users should consider their position when working with the cabinet, for their own safety.


Adjusting properly the seat and adopting a correct position enables the user to work safely and to avoid any possible physical damage (sight, neck, back, etc.)

Below are examples of four positions **NOT** recommended and one that **IS** recommended.


In addition, in order to protect both the product and the user, it is necessary that the cabinet is placed in a location without air currents, which may cause turbulence in the laminar down flow.

4.4. DEFINITION OF THE WORKING AREA

To prevent from creating any turbulence in the laminar down flow and to protect the user and the product, it is necessary to establish a working area inside the cabinet.

- Do not work in any area with perforations and do not cover them up.
- Always work within the limits according to the following image.

4.5. CONTROL KEYBOARD-DISPLAY

The following elements are found in the control keyboard**:**

A. *On/Off Button*  *(**Annex 3 pos. 1**):*

- Pressing when the cabinet is turned off:
 It goes to the *Password Request* screen. When entering the correct password (**Chapter 4.5 point 1**), it goes to the *Main Operation* screen.

- Pressing when the cabinet is turned on:
 It goes to the *Password Request* screen. When entering the correct password (**Chapter 4.5 point 1**), the cabinet switches off.

B. *Navigation buttons*  *,*  *and*  *(**Annex 3 pos. 2**):*

They move the cursor through all different function icons and select them. The icon pointed inverts its colours. If the selected item is a numeric field, the navigation buttons increase/decrease and validate the value.


C. *Backlight graphic display (**Annex 3 pos. 3**):*

It shows alarm messages, operating graphs, function icons and all the information required in cabinet-user communication. The background colour changes to:

- Green: Working mode. All the flows are working fine and there are no alarms activated.
- Red: Alarm mode. Any alarm has been activated or the cabinet is in unsafe mode. In case of alarm, an acoustic alarm is activated together with the red light.
- Yellow: Disinfection mode.
- White-blue: Information mode. It is used in Technical Service screens.

Below is an explanation of the different system screens:

1. ***Password Request*** screen *(by default **0000**).*

 Increase the value of the selected field.


 Decrease the value of the selected field.

 Confirm the value of the selected field and move the selection to the next field. After pressing for the final time:

- Right Password: it accesses the *Main Operation* screen.
- Wrong Password: Cursor returns to the first field.

2. ***Main Operation*** screen.

It shows the direct access icons to the most common functions for the user.

Activate/Deactivate fans.

Activate/Deactivate lights.

* Press for 5 seconds the button, to access the light intensity screen.

Activate/Deactivate sockets.


Activate/Deactivate gas solenoid valve (if assembled).*

Go into Disinfection Menu.

Go to Configuration Menu.

***REMARK:** *It is not possible to operate gas tap solenoid valve under unsafe conditions.*

2.1. **Fan Deactivation Confirmation** screen.


Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field.

2.2. **Light intensity** screen.


Increase the value.

Decrease the value.

Confirm the value.

2.3. **Disinfection Confirmation** screen.


Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field. After pin confirmation it goes to Disinfection Menu screen.

2.3.1. ***Disinfection Menu*** screen.


Go to *U.V. Menu* screen.

Go to *Vaporized hydrogen peroxide* screen.

Go to *Particles counter* screen (when applicable).

Back to *Main Operation* screen.

2.3.1.1. ***U.V. Menu*** screen.


Begin disinfection, if glass is closed, for the time determined by the timer.

Go to *U.V. Time Adjust* screen.

Back to *Disinfection* Menu screen.

2.3.1.1.1. ***U.V. Lamp*** screen.

Switch off the U.V. lamp and access to *Disinfection* Menu screen.

2.3.1.1.2. ***U.V. Time adjust*** screen.


Confirm the selected value and save it.*

Reset the numeric value.

Back to *U.V. Menu* screen.

* *Value 000 means infinite time: start and stop will be made manually.*

2.3.1.2. **Vaporized hydrogen peroxide** screen.


Start Vaporized hydrogen peroxide cycle.

Activate/Deactivate the sockets of the work zone.

Back to *Disinfection* Menu screen.

2.3.1.2.1. **Vaporized hydrogen peroxide cycle** screen.


Switch on manually the fans.

Switch off the fans.

Back to Disinfection menu screen (it requires confirmation).

3. ***Configuration Menu*** screen.

It shows the direct access icons to the configuration functions allowed to the user.

Go to Contrast screen.

Go to Maintenance screen.

Go to Language screen.


Go to PIN Configuration screen.

Go to Date & Hour screen.

Go to Sensors Information screen.

Go to Main Operation screen.

3.1. ***Contrast*** screen.


Increase the value.

Decrease the value.

Confirms the value and returns to *Configuration Menu* screen.

3.2. ***Maintenance Menu*** screen*.*

It shows information about the running hours of cabinet, fans, U.V. lamp and filters, together with the last Service visit. Besides, it allows the following:


Record a new revision.

Get maintenance recommendations.

Show the historical of revisions and alarms.

Back to *Configuration Menu* screen.

3.3. ***Language Selection*** screen*.*


Go up.

Go down.

Confirms the selected language and goes back to *Configuration Menu* screen.

3.4. ***PIN Configuration*** screen*.*


Enable/disable the password request to Connect/Disconnect the cabinet.

Enable the user to change the user password.

Coded password. In case of forgotten password, by giving this code to Telstar the previous password will be recovered.

Back to *Configuration Menu* screen.

3.5. **Date-hour** screen*.*


Saves the value.
(format dd/mm/yy hh:mm).

Cancel and/or edit the entered value.

Access to *Auto Power-on* menu.

Back to *Configuration Menu* screen.


3.5.1. **Auto Power on** screen.


Cabinet Auto Power on.

UV light Auto Power on.

Back to *Date-hour* screen.

3.5.1.1. **Cabinet Auto Power on** screen.


Saves the value.
(format dd/mm/yy hh:mm).

Cancel and/or edit the entered value.

Back to *Auto Power on* screen.

3.5.1.2. **UV light Auto Power on** screen.


Saves the value.
(format dd/mm/yy hh:mm).

Cancel and/or edit the entered value.

Back to *Auto Power on* screen.

3.6. **Sensors Information** screen.

It shows information about the values of the sensors of the cabinet in real time.

Move window up.

Move window to the work-position.

Move window down.

Back to *Configuration Menu* screen.

D. Microprocessor:

The microprocessor keeps total control over the exhaust flow and impulsion flow via the thermoanemometer probes.
The microprocessor control system consists of:

- Backlight graphic display.
- Alarms shown on the display:

ALARM	SHOWN ON THE DISPLAY
Excessive laminar down flow	HIGH LAM. FLOW
Insufficient laminar down flow	LOW LAM. FLOW
Excessive inflow	HIGH INFLOW
Insufficient inflow	LOW INFLOW
Glass out of working position	WORK POSITION
Front opened	LIFTED FRONTAL
Redundant channel alarm	REDUNDANT FAIL.
Power supply failure	POWER FAILURE
Memory fail	MEMORY FAILURE
Exhaust duct obstructed	CLOGG.EXH.DUCT

- Messages shown on the display:

MESSAGE	SHOWN ON THE DISPLAY
Uv light switched on	GERMICIDE
Battery needs replacement	LOW BATTERY
Glass in upper position	WINDOW OPEN
Window is moving	WINDOW MOVING

E. Auxiliary microprocessor or redundant circuit.

Should the main microprocessor fail, the auxiliary microprocessor will activate the optical alarm, and also the acoustic alarm. In such a case, move the glass to the closed position and cover the cabinet exhaust. Disconnect the equipment and contact the Telstar Technical Service.

4.6. DISCONNECTION

To stop the cabinet, it is necessary to push the On/Off button (**Annex A3 pos. 1**) to access the *Password Request* screen (**Chapter 4.5 point 1**). After entering the correct password, the cabinet will disconnect.

4.7. PROLONGED STOPPAGE

In the event of a prolonged stoppage, avoid environmental dust to go inside the cabinet by turning the cabinet off with the glass in the closed position. It is recommended to disconnect the battery to avoid its full discharge. When it is switched on again, proceed like the first time it was switched on. See **Chapters 3.3** and **4.1**.

4.8. OPTIONAL ACCESSORIES

- Gas tap kit.
- Vacuum tap kit.
- Compressed air tap kit.
- Nitrogen tap kit.
- Oxygen tap kit.
- Water tap kit.

WARNING 

Gas is a toxic and flammable substance which in the event of an accident can cause serious damage to the user and the cabinet.

The presence of a flame within the cabinet creates turbulence which affects its control and proper functioning; therefore, we recommend you do **NOT** use burners unless strictly essential and if the user is prepared to assume this responsibility.

Gas connection should be done in line with the safety standards of each country, region and company.

For these reasons, the electrically controlled valve that enables gas entry can only be activated when the fans are working, and it is automatically disconnected when the alarm indicating a lack of impulsion flow is activated.

- Vaporized Hydrogen Peroxide decontamination kit.
- I.V. bar with hooks.
- Waste port (for solids, at the left side of the cabinet).
- Humidity or Temperature sensor.
- IT configuration (option for the integration of software solutions). With or without supports -for PC, keyboard and printer- at the right side of the cabinet
- Weighing stone.
- External roof fan.

5. MAINTENANCE

WARNING ⚠

All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

Before starting any maintenance, work or dismantling the cabinet, the area should be wiped down with vaporized hydrogen peroxide (**Annex A5.4**). Once this has been done, disconnect the cabinet from the power socket.

In general, a regular control of the proper working of all cabinet components should be done once a year.

NOTICE ⚠

It is recommended to change the batteries every 2-3 years to guarantee the proper function. The disposal of exhausted batteries must be done in accordance with local requirements. The batteries constitute a chemical residue.

5.1. MAINTENANCE TABLE

The frequencies determined for each operation, and in particular with regard to changing filters, which appear on the Maintenance Table below are indicative values recommended for normal cabinet service. Very severe work conditions may reduce the time intervals shown on the Table.

TEST *(See annex A9)*	***REGULARITY***
Downflow filters integrity Test	*Minimum once a year or after filter changes*
Aspiration filters integrity Test	
Downflow air velocity test	
Set point alarm downflow air velocity test	
Inflow air velocity Test	
Set point alarm inflow air velocity test	
Smoke test	
Light test	*Recommended during the start-up of the cabinet*
Noise level Test	

5.2. ABSOLUTE FILTERS

Under normal circumstances, the life span of a HEPA filter comprises several years. However, it will become clogged up gradually, until there is no longer a proper flow-through and the filter has to be replaced.

After this time, the propelled airflow falls because of the increased pressure drop of the filter as it is a mechanical filtering, although this does not mean that filtering efficacy is reduced. When it reaches the minimum flow measured by the thermoanemometric probes, the microprocessor activates both visual and acoustic alarms.

WARNING 

Do not replace absolute filters or directly handle fans in MSCs without having first wiped the area down with vaporized hydrogen peroxide. Opening the cabinet without taking the proper precautions can endanger your health. We recommend it is carried out by staff especially trained for this work.

To ensure correct installation and adjustment of the filter, as well as ensuring there are no leaks, filters should be replaced by specialised technical staff equipped with the corresponding monitoring instruments: anemometer, leak detector, aerosol generator, particle counter, etc.

The used clogged absolute filters from safety cabinets should be destroyed. This process is always the responsibility of the user. Our technical service does not have any responsibility for possibly contaminated parts.

WARNING 

The replacement of the aspiration filters must be done with safety conditions keeping the cabinet working and following the instructions described in the Service Manual.

The inlet pre-filter is assembled on the top of the cabinet. It can be cleaned regularly using a vacuum cleaner or easily replaced.
The aspiration pre-filter is assembled under the worktop. It can be replaced regularly.

According to EN12469:2000 Annex K the pre-filters need to be replaced on every maintenance occasion. Replacement can be done easily following the instructions described in Annexes A4.1 & A4.2.

5.3. SPARE PARTS

Only use original spare parts from the manufacturer.

Please do not hesitate to contact the Customers Service pointing out the bench model, operating frequency (Hz), voltage (V) and serial number of the unit.

BiOptima Exa all models

No.	Description
1	UV LAMP PHILIPS TUV 36T5 4P SE UNP
2	CPU CARD
3	BATTERY (2 UNITS)
4	SPEED PROBE UP TO 1.0 M/S (1 UNIT) (Downflow)
5	SPEED PROBE UP TO 10.0 M/S (1 UNIT) (Extraction)

BiOptima Exa 4

No.	Description
6	DOWNFLOW (2 UNITS)
7	PRE-FILTER BIOPTIMA 4 G3 (EN 779)
8	PRE-FILTER BIOPTIMA EXA 4 G3 (EN 779)
9	ABSOLUTE FILTER H14 1220X457X93
10	ABSOLUTE FILTER H14 610X457X93
11	LED LIGHTING BIOPTIMA 4
12	V-SHAPE FILTER (9 UNITS)

BiOptima Exa 6

No.	Description
13	DOWNFLOW (3 UNITS)
14	PRE-FILTER BIOPTIMA 6 G3 (EN 779)
15	PRE-FILTER BIOPTIMA EXA 6 G3 (EN 779)
16	ABSOLUTE FILTER H14 1830X457X93
17	ABSOLUTE FILTER H14 1220X457X93
18	LED LIGHTING BIOPTIMA 6
19	V-SHAPE FILTER (14 UNITS)

6. CLEANING AND DISINFECTION

6.1. SURFACE CLEANING AND DISINFECTION

See **Annex A5.1**

6.2. FRONT GLASS CLEANING AND DISINFECTION

See **Annex A5.2**

6.3. DISINFECTION WITH VAPORIZED HYDROGEN PEROXIDE

It is recommended to do decontamination of the cabinet by using vaporized hydrogen peroxide emission, which ensures penetration in sealed and inaccessible areas of the Cabinet. It also guarantees the sterilisation of the working area, fans, plenum, HEPA filters and evacuation ducts.

Vaporized hydrogen peroxide systems provide rapid, safe, low-temperature decontamination methods for any enclosed area that may be contaminated with micro-organisms, including spore-forming bacteria.

See **Annex A5.4** for the disinfection with vaporized hydrogen peroxide procedure.

7. MAINTENANCE TABLE

See **Annex A6**

8. ELECTRICAL DIAGRAMS

See **Annexes A7 (230V) and Annexes A8 (120V)**

9. TESTS

9.1. DOWNFLOW FILTERS INTEGRITY TEST

See **Annex A9.1**

9.3. INFLOW FILTERS INTEGRITY TEST

See **Annex A9.3**

9.4. DOWNFLOW AIR VELOCITY TEST

See **Annex A9.4**

9.4.1. SET POINT ALARM DOWNFLOW AIR VELOCITY TEST

See **Annex A9.4.1**

9.5. INFLOW AIR VELOCITY TEST

See **Annex A9.5**

9.5.1. SET POINT ALARM INFLOW AIR VELOCITY TEST

See **Annex A9.5.1**

9.6. SMOKE TEST

See **Annex A9.6**

10. CERTIFICATES

10.1. EU & UKCA DECLARATIONS OF CONFORMITY

See **Annex A10.1**

11. TROUBLESHOOTING

No.	Type of error	Causes of error	Correction of error
1	**Insufficient exhaust flow** Alarm "LOW INFLOW" (see Chapter 4.5 point D)	Front perforation in the working area totally/partially covered.	Visually check the perforation and remove any type of obstruction.
		Clogged aspiration filter.	Change the Prefilters (**Annex A4.1**). Check the date of the last filter change and the hours it has been under operation. (**Chapter 4.5 point 3.2**). Notify the technical service to change the filter.
		Cabinet ducting too long.	Fit an additional fan at the end of the duct (**Chapter 3.5**).
2	**Insufficient laminar down flow** Alarm "LOW LAM. FLOW" (see Chapter 4.5 point D)	Front and rear perforation in the work zone totally/partially covered.	Visually check the perforation and remove any type of obstruction.
		Clogged laminar down flow filter or aspiration filter	Change the Prefilter (**Annex A4.1**). Check the date of the last filter change and the hours it has been under operation. (**Chapter 4.5 point 3.2**). Notify the technical service to change the filter.
3	**Excessive laminar down flow** Alarm "HIGH LAM. FLOW" (see Chapter 4.5 point D)	Air currents in front aperture.	Check and remove any airflow in front of the cabinet that distorts the reading of the sensor.

No.	Type of error	Causes of error	Correction of error
4	**Laminar down flow not working properly** Alarm "LOW LAM. FLOW" (see Chapter 4.5 point D)	Probe not working properly.	Check that the connection in the control card is correct (**Annex A7**).
		Laminar down flow fan stopped.	Check the fuse has not blown (**Annex A7**). Check that the connections in the control card are correct (**Annex A7**). Notify the Technical Service to change the fan.
5	**Exhaust not working properly** Alarm "LOW INFLOW" (see Chapter 4.5 point D)	Probe not working properly.	Check that the connection in the control card is correct (**Annex A7**).
		External exhaust fan stopped.	Check the fuse has not blown (**Annex A7**). Check that the connections in the control card are correct (**Annex A7**). Notify the Technical Service to change the external fan.
6	**Window does not move**	Front is blocked.	Visually check the front movement to remove any type of obstruction.
			Check the fuse has not blown (**Annex A7**).
		Motor of the window does not work.	Notify the Technical Service to change the motor of the window.
		Inductive sensors unadjusted.	Notify the Technical Service to adjust the Inductive sensors position.
		Main control is damaged.	Call a specialised technician to replace the control.
7	**Glass out of position** Alarm "WORK POSITION" (see Chapter 4.5 pto. D)	Glass position.	Move the glass into the working position.
		Working position inductive sensors unadjusted.	Open the frontal and check that the working position inductive sensor is not loose. Notify the Technical Service to receive instructions to adjust the inductive sensor.
		Main control is damaged.	Call a specialised technician to replace the control.

No.	Type of error	Causes of error	Correction of error
8	**Open frontal** Alarm "LIFTED FRONTAL" (see Chapter 4.5 pto. D)	Lifted frontal inductive sensor unadjusted.	Open the frontal and check that the open frontal position inductive sensor is not loose. Notify the Technical Service to receive instructions to adjust the inductive sensor.
		Main control is damaged.	Call a specialised technician to replace the control.
9	**Light does not turn on**	Light not working properly.	Check that the connection in the control card is correct (**Annex A7**).
			Call a specialised technician to replace the Led bar.
10	**UV germicidal lamp does not turn on**	Closed frontal inductive sensor unadjusted.	Open the frontal and check that the closed frontal position inductive sensor is not loose. Notify the Technical Service to receive instructions to adjust the inductive sensor.
		Glass not in closed position.	Check the front slides properly.
		UV lamp does not work.	Check the lamp holders are properly connected. Check that the connection in the control card is correct (**Annex A7**). Change the UV lamp.
		Ballast is damaged.	Change the ballast.
11	**Intermittent alarm activated**	Inductive sensors are unadjusted.	Open the frontal and check that the inductive sensors are not loose. Notify the Technical Service to receive instructions to adjust the inductive sensors.
		Main control is damaged.	Call a specialised technician to replace the control.
12	**Display does not turn on**	Power cord cable disconnected.	Check that power cord cable is properly connected in socket on the wall.
		No power.	Check voltage on the socket of the wall.
		Main control is damaged.	Call a specialised technician to replace the control.

WARNING ⚠

All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

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EN

This procedure must be done with the cabinet working (external extraction fan working). Then the negative pressure in the working area will stay. This is an extra safety aspect. Since this procedure can be done in a short time, the extra pollution of the HEPA filters will be low.


This procedure is described for a right-hand orientated person. The described left hand may be changed for a right hand and right for left hand if this is done consistently (this procedure needs consistence due to contamination danger).


Supplies needed:


- Plastic bag (minimum measurement 75x50 cm with closure).
- Two pair of gloves.
- Tape.
- Possible disposing requirements.
- New pre-filter.


Procedure:


a) Switch on the unit.
b) Make sure that the window is in the working position and that there are no alarms.
c) Put on two pair of gloves.
d) Lift up the segmented worktop to the left of the chamber, in order to reach the prefilter. (Fig. 1)
e) Fold the edge of the plastic bag ±10 cm and place the opening of the plastic bag inside the front aperture with your left hand. (Fig. 2)
f) Position the plastic bag on the armrest.
g) Pull the pre-filter from its frame with your right hand.
h) Roll the filter up (still with your right hand). (Fig. 3)
i) Put the filter into the plastic bag. (Fig. 4)
j) Fold back the edge of the plastic bag, taking off your right glove and introducing it in the bag, not touching the inner side of the plastic bag. (Fig. 5)
k) Press the air out of the bag (with its opening still in the working space).
l) Seal the plastic bag airtight.
m) Remove the plastic bag and place it in an appropriate container.
n) Put the segmented worktop at the right and, using the right hand, introduce and place the new pre-filter in the frame. (Fig. 6)
o) Put the segmented worktop back in its position.
p) Take off the gloves and dispose them away in the appropriate container.

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
For the installation of the new prefilter, all steps should be carried out in reverse order, checking that it is perfectly located.


Fig. 1


Fig. 2

EN

WARNING ⚠

First, ensure that the fans are fully stopped.

For the replacement of the upper inlet prefilter, these steps should be followed:

a) Unplug the cord plug of the cabinet. The prefilter is located at the upper part of the cabinet (Fig. 1).
b) Remove the 4 screws that hold the upper grid of the prefilter by means of a 3mm allen key. (Fig. 2)
c) Remove the upper grid of the prefilter. Then you can access to the prefilter and replace it.

For the installation of the new prefilter, all steps should be carried out in reverse order, checking that it is perfectly located.


Fig. 1


Fig. 2

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EN

	**CYTOTOXIC DANGER** Before starting work in cytotoxic production the hazard risks of cytotoxic handling need to be evaluated and documented (industrial safety act, hazardous substances regulations). Based on these findings the employees must be instructed. For several cytotoxic, alternatives cleaning procedures could be applied.
	**BIOLOGICAL DANGER** Never clean a used unit without safety precautions

The unit must be cleaned regularly. Preferably, each time after the unit has been used.
- Put on gloves and sleeves.
- Use glasses and mouth protection.

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EN

For worktable and sides, use sterile or non-use disable cloths, which do not shed particles or fibres.

For disinfection, these cloths will be used, slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or affect the stainless steel and the glass. Cleaning can also be carried out with water and soap prior to disinfection.
To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabinet; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

WARNING

It is important keep the glass in working position and the cabinet in working mode when performing the cleaning and disinfection of the working area of the cabinet. In case that it is needed to open the frontal to clean areas with difficult access, all needed cautions must be taken because the user protection decreases when the frontal is open.
It is recommended NOT to use disinfectant products with chlorinated derivatives or containing halogens or their compounds, because they can damage the stainless steel and painted surfaces.
Should some bleach be poured, the surface must be always thoroughly rinsed and completely dried with a spatula (glass cleaner type).


It is required to dry completely all surfaces always after cleaning or disinfecting to preserve the original qualities of the surfaces.

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EN

To accede to corners or difficult access areas, another type of tools will be used as shown in the image.

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EN

The perforated stainless-steel worktable is divided into several parts to facilitate cleaning of the interior container and provide for sterilisation of these parts into a vessel.
Remember that every segment or shelve must be cleaned and/or disinfected both in its lower and upper sides. The parts will be cleaned inside the cabinet while this is running.

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
EN

Also be sure of carrying out the cleaning/disinfection in the lower part of the working area (surface and corners).

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EN

Inside, safety cabinet cleaning is recommended on the following occasions:

a) Before starting any work in the safety cabinet.
b) After working in the safety cabinet.
c) Whenever there is a change of work program.
d) In event of liquid spilling on the worktable.
e) Before carrying out mechanical or biological control tests in the working area.

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EN

For the cleaning of front glass, you must use sterile textures or disposable ones without particle or fibres loosing.

For disinfection, use the previous textiles slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or alter the glass. Cleaning can also be carried out with water and soap prior to disinfection.

To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabin; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

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EN

For the glass exterior part cleaning, it is recommended to move the glass into closed position. In this position a major part of the glass is accessible.
For cleaning the most accessible interior parts of the glass, it is recommended to leave the glass in its working position to assure the user's protection.

WARNING 

It is important to clean and disinfect the cabinet's working area, with the glass in working position and the cabinet in performance. If it is necessary to separate the frontal-glass mechanism to clean the areas of difficult access, bring the glass to its fully lifted position to set free the frontal mechanism. By separating the frontal-glass, the inner surface of the glass is totally accessible. There must be taken the necessary precautions because the user's protection is lower with the glass separated from cabinet.

WARNING 

The movement of the glass exerts a strong force. Check that no one has their hands, any part of the body or obstacles below the glass, as it can cause injuries. Warn other people working near the cabinet.

FR	Annexe A5.2. NETTOYAGE ET DÉSINFECTION VITRE FRONTALE	GE	Anhang A5.2. REINIGUNG UND DESINFEKTION DER FRONTSCHEIBE

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EN

Gas decontamination of a cabinet that has handled biologically pathogenic products is necessary in the following cases:

- *Before undertaking any safety cabinet maintenance work.*
- *Before removal of HEPA filters for replacement.*
- *Before carrying out mechanical or biological control tests in the safety cabinet.*
- *Before moving the safety cabinet from one work area to another.*
- *Before changing a work program.*
- *Whenever liquid is spilled on the worktable.*

The gas used in this procedure is the vaporized hydrogen peroxide.

The necessary material to perform decontamination by gas is as follows:

- *1 vaporized hydrogen peroxide decontamination kit accessory.*
- *1 Plastic bag 1800x600 mm.*
- *1 Roll of packaging tape.*
- *1 Roll of adhesive tape.*
- *1 Allen key.*
- *Vaporization hydrogen peroxide system.*
- *Vaporized hydrogen peroxide sensor.*

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

EN


Steps to follow:

a) Fumigation must be performed by a qualified person with a good knowledge of the procedure and the required precautions. At any case, a visible warning should indicate the cabinet is being disinfected.

b) *The vaporized hydrogen peroxide accessory kit needs to be installed on the side glass (inlet port), and on the bottom manhole cover (outlet port). Seal the upper pre-filter with a bag and packing tape. If you do not have a butterfly valve to seal the extraction duct, it is necessary to disassemble the hard duct and seal the cabinet extraction opening with a bag and duct tape.*

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EN

c) Separate the frontal of the cabinet.

d) To guarantee the tightness of the equipment, seal the cabinet front opening by means of one plastic bag and packaging tape, covering any gap that could leak hydrogen peroxide to the ambient.

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EN

e) Close the front window and place the glass in working position.

f) Connect the vaporization hydrogen peroxide system hoses to the proper port. Check that input of vapor is fixed in the side glass, and output is fixed on top.

g) Cabinet is ready to start with decontamination cycle.

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EN

Decontamination process.

It must be performed by specialized personnel of the vaporization hydrogen peroxide system.

After cycle is finished, high concentration may remain inside the cabinet. Once the plastic sealing the front and the exhaust cover are removed, concentration will be diluted with room air and be reduced rapidly to accepted levels.

ATTENTION

During process, the hydrogen peroxide concentration in the room must be monitored for safety reasons, so it will be necessary to place a sensor for that purpose. Ambient limit is 1 ppm, and at that concentration 8 h of continuous work is allowed. With ambient concentration of 5 ppm* 2h of continuous work is allowed.*

**** 0.5 ppm in Germany***

a) Even after the treatment, the HEPA filters can only be considered to be 'safe to handle when using the appropriate protective clothing' instead of no sterile. They should be treated in an autoclave or incinerated on removal.
b) Suitable protective respiratory equipment designed for being used with hydrogen peroxide should be available in case of an emergency.
c) Once the process has finished, switch off the equipment through the control panel.

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REVIEWED: X. Cañellas DATE: 12/11/2018

AS BUILT: DATE:

REVISIONS:
ID: F DATE: 12/11/2021 DRAWN BY: X.Cañellas
ID: G DATE: 01/03/2022 DRAWN BY: X.Cañellas
ID: DATE: DRAWN BY:
ID: DATE: DRAWN BY:
ID: DATE: DRAWN BY:

Telstar®

Date 23/04/2015 Report of Drawings 1756262 Página 1
BIOPTIMA EXA
AZBIL TELSTAR TECHNOLOGIES, S.L.


Drawing	Description	Creation	Modification
01	SUMMARY	12/11/2021	---
02	POWER	12/11/2018	---
03	CONTROL BOARD	12/11/2018	---
04	BLOWERS	12/11/2018	---
05	LIGHTING AND UV KIT	12/11/2018	01/03/2022
06	WINDOW MOTOR, SOCKETS & GAS VALVE	12/11/2021	01/03/2022
07	WINDOW DETECTORS	12/11/2018	---
08	AIR FLOW SENSORS	12/11/2018	---
09	DISPLAY AND BUTTONS	12/11/2018	---
10	ACCESSORIES BOARD (OPTIONAL PARTS)	12/11/2018	---
11	BILL OF MATERIALS	12/11/2021	---


PROPERTY OF TELSTAR
REPRODUCTION FORBIDDEN

REPLACES:

REPLACED BY:

PROJECT NAME:
1756262

CABLES: 1014825
CABLES: 1014825
CABLES: 1014783
-F1 16 A
L
PE
N
1L
1N15
1L15
1N15
-T1
644186
R12.5
12.5A
634174
-FL1
YB12E4-20A-Q
Note:
Adjust the output voltage to 28.5Vdc
0V17
27V17
2L15
2N15
620999
-PS1
230/110 VAC
27 VDC
1.3/2.5 A
200 W
CABLES: 1014783
YV/GR
BLUE
BROWN
PE
JN5
2EDGK-7.5-5P
J15
BATTERY
12 Vdc
1.3 Ah
-BT1
638420
-BT2
638420
BLACK
RED
CABLES: 1014826
CABLES: 1014783
0B16
24B16
27V17
0V17
J16
J17
2EDGK-5.08-2P
2EDGK-5.08-2P
ELECTRICAL CHARACTERISTICS
MODEL
EXA 6
EXA 4
VOLTAGE
230VAC
230VAC
POWER
2.2 KVA
1.8 KVA
CURRENT
9.6 A
8.2 A
FREQUENCY
50..60 HZ
50..60 HZ
Revision
Date
Draw by
Date:
Drawn by:
Aprovet:
12/11/2018
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
POWER
Drawings: 2 / 11
Previous: 1
Next 3
Project reference:
977083
Drawing number:
1756262
Drawing number
2
Telstar



CONTROL BOARD


CONNECTORS			FUSES
J1 - DISPLAY	J10 - SOCKETS 230VAC 4A	J19 - FRONT HOOD SENSOR	
J2 - WINDOW MOTOR	J11 - not used	J20 - not used	F1 - 250VAC 8A/5x20R DC POWER
J3 - ELECTROVALVE	J12 - not used	J21 - WORK POSITION W.SENSOR	F2 - 250VAC 10A/5x20R IMP FAN
J4 - LED LIGHTING	J13 - DOWNFLOW FAN CONTROL SIGNAL	J22 - not used, connected to J21	F3 - 250VAC 2.5A/5x20R EXT FAN
J5 - EXHAUST AIR FLOW SENSOR	J14 - POWER SUPPLY DOWNFLOW FANS	J23 - UP POSITION W.SENSOR	F4 - 250VAC 5A/5x20R SOCKETS
J6 - DOWN AIR FLOW SENSOR	J15 - MAIN POWER SUPPLY 230/10VAC	J24 - not used	F5 - 250VAC 1.6A/5x20R LIGH" (not used)
J7 - ACCESSORY BOARD CONNECTION	J16 - BATTERY	J25 - CLOSE POSITION W.SENSOR	F6 - 250VAC 1A/5x20R UV LIGHT
J8 - UV LIGHT	J17 - 28.5VDC POWER SUPPLY	J26 - CLOSE POSITION W.SENSOR REDN.	
J9 - not used	J18 - WINDOW PUSHBUTTONS		


Revision	Date	Draw by	Date:	12/11/2018	Project description:	Situation:	Zone:	Drawings: 3 / 11	Previous: 2
			Draw by:	X. Cañellas	BIOPTIMA EXA	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: CONTROL BOARD		Drawing number: 1756262
			Aproved:	X. Cañellas				Project reference: 977083	Drawing number 3


Next 4


DOWNFLOW FANS
6 FEET CABINET
4 FEET CABINET
J14
2DEGK-7.5-3P
CABLES: 1014777
L14
N14
PE
BROWN
BLUE
YV/GR
J13
15EDGK-3.81-3
CABLES: 1653380
SH13
0V13
VS13
yellow
blue
-FAN1
-FAN2
-FAN3
EC
Tach
+10V
0-10V
GND
L
N
611044
170W
1.4A
CABLES:
EXA 4: 1014763
EXA 6: 1014762
CABLES:
EXA 4: 1014828
EXA 6: 1014830
J34
15EDGK-3.81-3
CABLES: 1793929
BLACK
RED
J35
2DEGK-7.5-3P
GREEN
BROWN
WHITE
BLACK
RED
EXTERNAL FAN CONNECTIONS
(CUSTOMER FACILITIES)
COM
NO
NC
GND
0-10V
FREE POTENTIAL
CONTACT
(ON/OFF)
VELOCITY
SIGNAL
Revision
Date
Draw by
Date:
Drawn by:
Aprovet:
12/11/2018
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
BLOWERS
Drawings: 4 / 11
Previous: 3
Next 5
Project reference:
977083
Drawing number:
1756262
Drawing number
4
Telstar


J4
2EDGK-5.08-3
SH4
0V4
VS4
CABLES: 1653390
CABLES: 1014781
SHIELD
BLACK
RED
CABINET
EXA 4
EXA 6
LED STRIP
741339
741347
RED
BLACK
-LED1
24VDC
19W/m
LED STRIP
J8
2EDGK-7.5-3P
L8
N8
PE
CABLES: 1752382
BROWN
BLUE
PE
-L1
BLACK
WHITE
BLUE
BROWN
24
23
22
21
N8
L8
620816
-BL1
QTP-OPTIMAL 1x18-40
CABLES: 1653391
CABLES: 1014781
CABINET
EXA 4
EXA 6
UV LIGHTING
32073
32073
UV LIGHT
Revision
G
Date
01/03/2022
Draw by
X. Cañellas
Date:
23/04/2015
Draw by:
X. Cañellas
Aprovet:
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
LIGHTING AND UV KIT
Drawings:
5 / 11
Previous:
4
Next:
6
Project reference:
977083
Drawing number:
1756262
Drawing number
5
Telstar

J2
2EDGK-5.08-2P
J10
2EDGK-7.5-3P
J3
2EDGK-5.08-2P
WINDOW MOTOR
651474
-M1
24VDC
VS2
0V2
BLUE
BROWN
SOCKETS
510032
-SK1
230VAC
2A
510032
-SK2
230VAC
2A
PE
N
L
N10
L10
GR/YEL
BLUE
BROWN
CABLES: 1014758
CABLES: 1014772
EXA 4: 1014772
EXA 6: 1014774
GAS VALVE
640621
-SV1
24VDC
8W
VS3
0V3
RED
BLACK
CABLES: 1014775
CABLES: 1014779
Revision
F
G
Date
12/11/2021
01/03/2022
Draw by
X.Cañellas
X.Cañellas
Date:
Draw by:
Aprovet:
12/11/2018
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
WINDOW MOTOR
SOCKETS
GAS VALVE
Drawings: 6 / 11
Previous: 5
Next: 7
Project reference:
977083
Drawing number:
1756262
Drawing number
6

UP POSITION SENSOR
WORK POSITION SENSOR
CLOSE POSITION SENSOR
FRONT HOOD SENSOR
CLOSE POSITION SENSOR REDUNDANT
J23
J21
J22
J25
J19
J26
15EDGK-3.81-3
0V23
VS23
24V23
0V21
VS21
24V21
0V25
VS25
24V25
0V19
VS19
24V19
0V26
VS26
24V26
BLUE -
BLACK S
BROWN +
-IS1 64287 TL-V5M1B1 DMRDN
-IS2 64287 TL-V5M1B1 DMRDN
-IS3 64287 TL-V5M1B1 DMRDN
-IS4 64287 TL-V5M1B1 DMRDN
-IS5 64287 TL-V5M1B1 DMRDN
Revision
Date
Draw by
Date:
Drawn by:
Aprovel:
12/11/2018
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
WINDOW DETECTORS
Drawings: 7 / 11
Previous: 6
Next 8
Project reference:
977083
Drawing number:
1756262
Drawing number
7
Telstar

J5
15EDGK-3.81-4
1 2 3 4
24V5
0V5
VS5
1 RED
2 BLACK
3 WHITE
CABLES: 1014831
CABLES: 1014769
621613
-SP1
FS22503
0.2-10m/s
EXHAUST
VELOCITY
PROBE
J6
15EDGK-3.81-4
1 2 3 4
24V6
0V6
VS6
1 RED
2 BLACK
3 WHITE
CABLES: 1014831
CABLES: 1014769
614060
-SP2
FS22500
0.2-1m/s
DOWNFLOW
VELOCITY
PROBE
Revision
Date
Draw by
Date:
Draw by:
Aproved:
12/11/2018
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
AIR FLOW SENSORS
Drawings: 8 / 11
Previous: 7
Next: 9
Project reference:
977083
Drawing number:
1756262
Drawing number
8
Telstar

J1
15EDGK-3.81-5P
CABLES: 1014782
Shield
YELLOW
GREEN
BROWN
WHITE
J2
632700
-DIST
J18
15EDGK-3.81-4
CABLES: 1014832
GREEN
BLUE
RED
BLUE
RED
GREEN
J2
TELSTAR
J1
-VB1
1384557
Revision
Date
Draw by
Date:
Drawn by:
Aprovet:
12/11/2018
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
DISPLAY AND BUTTONS
Drawings: 9 / 11
Previous: 8
Next 10
Project reference:
977083
Drawing number:
1656262
Drawing number
9

627589
-PCB2
ATLAS
ACCESSORIES BOARD

RS485 COMMPORT

J28
15EDGK-3.81-4

4 SHIELD
3 GND
2 B -
1 A +

VHP CONTROL

J31
15EDGK-3.81-5

5 NC
4 NO
3 COM
RELAY ON/OFF MAX. 230VAC 10A
2 GND
1 ALARM (IN)
ALARM INPUT SIGNAL

BUTTERFLY VALVE

J33
15EDGK-3.81-3

3 SHIELD
2 GND
1 POSITION SIGNAL
POSITION SIGNAL

J32
2EDGK-7.5-3P

3 NC
2 NO
1 COM
RELAY ON/OFF CONTROL (MAX. 230VAC 10A)

BOOSTER FAN

J34
15EDGK-3.81-3

3 SHIELD
2 GND
1 0-10V OUT
VELOCITY CONTROL

J35
15EDGK-3.81-3

3 NC
2 NO
1 COM
RELAY ON/OFF CONTROL (MAX. 230VAC 10A)

ALARM SIGNAL

J36
15EDGK-3.81-3

3 NC
2 NO
1 COM
RELAY ON/OFF ALARM (MAX. 230VAC 10A)

TEMP./HUMID./PRESS. SENSOR

J37
15EDGK-3.81-4

4 SHIELD
3 4..20mA IN
2 GND
1 24VDC

PARTICLE COUNTER SENSOR

J38
15EDGK-3.81-4

4 SHIELD
3 4..20mA IN
2 GND
1 24VDC

EXTERNAL AIRFLOW SENSOR

J39
15EDGK-3.81-4

4 SHIELD
3 0-5V IN
2 GND
1 24VDC

ACCESSORIES BOARD	
CONNECTORS	
J39 - EXTERNAL AIRFLOW SENSOR	J33 - BUTTERFLY VALVE INPUT SIGNAL
J38 - PARTICLE COUNTER SENSOR	J32 - BUTTERFLY VALVE RELAY
J37 - TEMPERATURE/PRESSURE/HUMIDITY	J36 - ALARM RELAY SIGNAL
J28 - RS485 COMMPORT	J35 - BOOSTER FAN POWER RELAY
J40 - TO CONTROL BOARD CONNECTION	J34 - BOOSTER FAN 0-10V SIGNAL
J31 - VHP	

ACCESSORIES
(OPTIONALS PARTS)

Revision	Date	Draw by	Date:	12/11/2018	Project description: BIOPTIMA EXA
			Draw by:	X. Cañellas	
			Aproved:	X. Cañellas	

Situation:
Zone:
Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.
Drawing description: ACCESSORIES BOARD OPTIONAL PARTS

Drawings: 10 / 11
Previous: 9
Next: 11
Drawing number: 1756262
Project reference: 977083
Drawing number: 10

Telstar

Date 24/04/2014 Bill of Materials Page 1

BIOPTIMA EXA

AZBIL TELSTAR
TECHNOLOGIES, S.L.

Symbol	Drawing	Column	Description	Refer. Numb.	Catalogue	Reference
-FL1	02	B	POWER LINE FILTER EMC 20A YANBIXINKEJI	YB12E4-20A-Q	---	634174
-F1	02	B	SLOW BLOW FUSE 16A SCHURTER	0001.2516.T	---	616733
-I1	02	B	SINGLE PHASE LINE INDUCTOR POLYLUX	R12.5	---	644186
-PS1	02	E	POWER SUPPLY 27VDC 200W MEANWELL	RSP-200-27	---	620899
-BT1	02	F	BATERIA 12V 1,3AH BAOTE	BT-12M1.3AT	---	636420
-BT2	02	G	BATERIA 12V 1,3AH BAOTE	BT-12M1.3AT	---	636420
-PCB1	03	C	ATLAS BSC CONTROL CARD	---	---	632679
-PCB2	03	F	ATLAS ACCESSORIES EXPANSION CARD	---	---	627589
-FAN1	04	C	EBM EC FAN	K3G225-RE07-03	---	611044
-FAN2	04	E	EBM EC FAN	K3G225-RE07-03	---	611044
-FAN3	04	F	EBM EC FAN	K3G225-RE07-03	---	611044
-LED1	05	C	LED STRIP (BIOPTIMA EXA 4)	---	---	1741339
			LED STRIP (BIOPTIMA EXA 6)	---	---	1741347
-BL1	05	G	BALASTRO T5 OSRAM	QTP-OPTIMAL 1X18-40	---	620816
-L1	05	G	UV TUBE PHILIPS TUV 36W	TUV 36T5 4P SE UNP	---	32073
-M1	06	F	DC TUBULAR MOTOR 24VDC	FT45DS-20/8	---	651474
-SK1	06	C	ELECTRIC SOCKET	B-012H	---	510032
-SK2	06	D	ELECTRIC SOCKET	B-012H	---	510032
-SV1	06	F	ELECTROVALVE 24VDC 8W	266293 TYPE 6013	---	640621
-IS1	07	C	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS2	07	C	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS3	07	C	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS4	07	G	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS5	07	G	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-SP1	08	C	AIR VELOCITY SENSOR 0-10m/s 0-5V DEGREE	FS22503	---	621613
-SP2	08	E	AIR VELOCITY SENSOR 0-1m/s 0-5V DEGREE	FS22500	---	614060
-DIS1	09	D	DISPLAY PARA BSC ATLAS	---	---	632700
-WB1	09	F	WINDOW BUTTONS ATLAS	---	---	1384557

Revision	Date	Draw by	Date:	12/11/2018	Project description:	Situation:	Zone:	Drawings: 11 / 11	Previous: 10	Next:
F	12/11/2021	X.Cañellas	Draw by:	X. Cañellas	BIOPTIMA EXA	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: BILL OF MATERIALS		Drawing number: 1756262	
			Aproved:	X. Cañellas				Project reference: 977083	Drawing number 11	

SP	**Anexo A8. ESQUEMA ELÉCTRICO 120V 50/60 Hz**
EN	**Annex A8. ELECTRICAL DIAGRAM 120V 50/60 Hz**
FR	**Annexe A8. SCHÉMA ÉLECTRIQUE 120V 50/60Hz**
GE	**Anhang A8. SCHALTPLAN 120V 50/60 Hz**

EQUIPMENT:
BIOPTIMA EXA 120V

PROJECT REFERENCE:
977083

CUSTOMER:
AZBIL TELSTAR TECHNOLOGIES, S.L.

DRAWING TYPE:
ELECTRIC

DRAWINGS:
11

MATERIAL O.F.:

PROJECT LEADER:
Victor Lazaro

AUTOMATION LEADER:
10/06/2022

DRAWN BY:
X. Cañellas
DATE: 10/06/2022

APPROVED:
X. Cañellas
DATE: 10/06/2022

REVIEWED:
X. Cañellas
DATE: 10/06/2022

AS BUILT:
DATE:

REVISIONS:
ID: DATE: DRAWN BY:
ID: DATE: DRAWN BY:
ID: DATE: DRAWN BY:
ID: DATE: DRAWN BY:
ID: DATE: DRAWN BY:

Telstar

Date 10/06/2022 Report of Drawings 1917194 Página 1
BIOPTIMA EXA
AZBIL TELSTAR TECHNOLOGIES, S.L.


Drawing	Description	Creation	Modification
01	SUMMARY	10/06/2022	---
02	POWER	10/06/2022	---
03	CONTROL BOARD	10/06/2022	---
04	BLOWERS	10/06/2022	---
05	LIGHTING AND UV KIT	10/06/2022	---
06	WINDOW MOTOR, SOCKETS & GAS VALVE	10/06/2022	---
07	WINDOW DETECTORS	10/06/2022	---
08	AIR FLOW SENSORS	10/06/2022	---
09	DISPLAY AND BUTTONS	10/06/2022	---
10	ACCESSORIES BOARD (OPTIONAL PARTS)	10/06/2022	---
11	BILL OF MATERIALS	10/06/2022	---

PROPERTY OF TELSTAR
REPRODUCTION FORBIDDEN

REPLACES:


REPLACED BY:

PROJECT NAME:
1917194

ELECTRICAL CHARACTERISTICS		
MODEL	EXA 6	EXA 4
VOLTAGE	120VAC	120VAC
POWER	1.6 KVA	1.4 KVA
CURRENT	13 A	11.5 A
FREQUENCY	50..60 HZ	50..60 HZ


Revision			
Date			
Draw by			
Date:	10/06/2022		
Draw by:	X. Cañellas		
Aproved:	X. Cañellas		
Project description:	BIOPTIMA EXA		
Situation:			
Customer:	AZBIL TELSTAR TECHNOLOGIES, S.L.		
Zone:			
Drawing description:	POWER		
Drawings:	2 / 11		
Previous:	1		
Project reference:	977083	Drawing number:	1917194
		Drawing number:	2
Next:	3		

CONTROL BOARD

CONNECTORS		FUSES	
J1 - DISPLAY	J10 - SOCKETS 230VAC 4A	J19 - FRONT HOOD SENSOR	F1 - 250VAC 8A/5x20R DC POWER
J2 - WINDOW MOTOR	J11 - not used	J20 - not used	F2 - 250VAC 10A/5x20R IMP FAN
J3 - ELECTROVALVE	J12 - not used	J21 - WORK POSITION W.SENSOR	F3 - 250VAC 2.5A/5x20R EXT FAN
J4 - LED LIGHTING	J13 - DOWNFLOW FAN CONTROL SIGNAL	J22 - not used, connected to J21	F4 - 250VAC 5A/5x20R SOCKETS
J5 - EXHAUST AIR FLOW SENSOR	J14 - POWER SUPPLY DOWNFLOW FANS	J23 - UP POSITION W.SENSOR	F5 - 250VAC 1.6A/5x20R LIGHT (not used)
J6 - DOWN AIR FLOW SENSOR	J15 - MAIN POWER SUPPLY 230/110VAC	J24 - not used	F6 - 250VAC 1A/5x20R UV LIGHT
J7 - ACCESSORY BOARD CONNECTION	J16 - BATTERY	J25 - CLOSE POSITION W.SENSOR	
J8 - UV LIGHT	J17 - 28.5VDC POWER SUPPLY	J26 - CLOSE POSITION W.SENSOR REDN	
J9 - not used	J18 - WINDOW PUSHBUTTONS		

Revision	Date	Draw by	Date:	Project description:	Situation:	Zona:	Drawings: 3 / 11	Previous: 2	Next: 4
			10/06/2022	BIOPTIMA EXA	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: CONTROL BOARD		Drawing number: 1917194	
		Draw by:	X. Cañellas				Project reference: 977083		Drawing number: 3
		Aproved:	X. Cañellas						

Revision	Date	Draw by	Date:	10/06/2022
			Draw by:	X. Cañellas
			Aproved:	X. Cañellas

Project description: BIOPTIMA EXA

Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.

Situation:

Zone:

Drawing description: BLOWERS

Drawings: 4 / 11


Previous: 3


Next: 5

Project reference: 977083


Drawing number: 1917194

Drawing number: 4


J4
2EDGK-5.08-3
SH4
0V4
VS4
CABLES: 1653390
SHIELD
BLACK
RED
CABLES: 1014781
CABINET
LED STRIP
EXA 4
1741339
EXA 6
1741347
RED
BLACK
-LED1
24VDC
19W/m
LED STRIP
J8
2EDGK-7.5-3P
PE
N8
L8
BLACK
WHITE
32073
-L1
BLACK
WHITE
BLUE
BROWN
RED
RED
BLUE
BLUE
YELLOW
YELLOW
UV LIGHT
WHITE
BLACK
637566
-BL1
ACT2S39PSUNV
Revision
Date
Draw by
Date:
Draw by:
Aproved:
10/06/2022
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Situation:
Zona:
Drawing description:
LIGHTING AND UV KIT
Drawings: 5 / 11
Previous: 4
Next: 6
Project reference:
977083
Drawing number:
1917194
Drawing number
5

WINDOW MOTOR
SOCKETS
GAS VALVE
J2
2EDGK-5.08-2P
J10
2EDGK-7.5-3P
J3
2EDGK-5.08-2P
VS2
0V2
BLUE
BROWN
651474
-M1
24VDC
PE
N10
L10
GR/YEL
BLUE
BROWN
CABLES 1014768
CABLES 1014772
EXA 4: 1014772
EXA 6: 1014774
510032
-SK1
230VAC
2A
510032
-SK2
230VAC
2A
VS3
0V3
RED
BLACK
CABLES 1014775
CABLES 1014779
640621
-SV1
24VDC
8W
Revision
Date
Draw by
Date:
Draw by:
Aproved:
10/06/2022
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zona:
Drawing description:
WINDOW MOTOR
SOCKETS
GAS VALVE
Drawings: 6 / 11
Previous: 5
Next: 7
Project reference:
977083
Drawing number
1917194
Drawing number
6
Telstar


UP POSITION SENSOR


FRONT HOOD SENSOR


WORK POSITION SENSOR


CLOSE POSITION SENSOR



CLOSE POSITION SENSOR REDUNDANT

Revision	Date	Draw by			Project description:	Situation:	Zone:	Drawings: 7 / 11	Previous: 6	Next: 8
			Date:	10/06/2022	BIOPTIMA EXA	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: WINDOW DETECTORS		Drawing number: 1917194	
			Draw by:	X. Cañellas						
			Aproved:	X. Cañellas				Project reference: 977083	Drawing number 7	

Revision
Date
Draw by
Date:
Draw by:
Aproved:
10/06/2022
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Situation:
Zona:
Drawing description:
AIR FLOW SENSORS
Drawings: 8 / 11
Previous: 7
Next: 9
Project reference:
977083
Drawing number
1917194
Drawing number
8
J5
15EDGK-3.81-4
1 RED
2 BLACK
3 WHITE
24V5
0V5
VS5
CABLES: 1014769
CABLES: 1014831
EXHAUST VELOCITY PROBE
621613 -SP1 FS22503 0.2-10m/s
J6
15EDGK-3.81-4
1 RED
2 BLACK
3 WHITE
24V6
0V6
VS6
CABLES: 1014769
CABLES: 1014831
DOWNFLOW VELOCITY PROBE
614060 -SP2 FS22500 0.2-1m/s

J1
15EDGK-3.81-5P
CABLES: 1014782
0V1
12V1
B-1
A+1
SHIELD
YELLOW
GREEN
BROWN
WHITE
J2
632700
-DIS1
J18
15EDGK-3.81-4
CABLES: 1014832
GREEN
BLUE
RED
BLUE
RED
GREEN
J2
J1
TELSTAR
-VB1
1384557
Revision
Date
Draw by
Date:
Draw by:
Aproved:
10/06/2022
X. Cañellas
X. Cañellas
Project description:
BIOPTIMA EXA
Situation:
Customer:
AZBIL TELSTAR TECHNOLOGIES, S.L.
Zona:
Drawing description:
DISPLAY AND BUTTONS
Drawings: 9/11
Previous: 8
Next: 10
Project reference:
977083
Drawing number
1917194
Drawing number
9

627589
-PCB2
ATLAS
ACCESSORIES BOARD

RS485 COMMPORT

J28
15EDGK-3.81-4

4 SHIELD
3 GND
2 B -
1 A +

VHP CONTROL

J31
15EDGK-3.81-5

5 NC
4 NO
3 COM
2 GND
1 ALARM (IN)

RELAY ON/OFF MAX. 230VAC 10A
ALARM INPUT SIGNAL

BUTTERFLY VALVE

J33
15EDGK-3.81-3

3 SHIELD
2 GND
1 POSITION SIGNAL

POSITION SIGNAL

J32
2EDGK-7.5-3P

3 NC
2 NO
1 COM

RELAY ON/OFF CONTROL
(MAX. 230VAC 10A)

BOOSTER FAN

J34
15EDGK-3.81-3

3 SHIELD
2 GND
1 0-10V OUT

VELOCITY CONTROL

J35
15EDGK-3.81-3

3 NC
2 NO
1 COM

RELAY ON/OFF CONTROL
(MAX. 230VAC 10A)

ALARM SIGNAL

J36
15EDGK-3.81-3

3 NC
2 NO
1 COM

RELAY ON/OFF ALARM
(MAX. 230VAC 10A)

TEMP./HUMID./PRESS. SENSOR

J37
15EDGK-3.81-4

4 SHIELD
3 4..20mA IN
2 GND
1 24VDC

PARTICLE COUNTER SENSOR

J38
15EDGK-3.81-4

4 SHIELD
3 4..20mA IN
2 GND
1 24VDC

EXTERNAL AIRFLOW SENSOR

J39
15EDGK-3.81-4

4 SHIELD
3 0-5V IN
2 GND
1 24VDC

ACCESSORIES BOARD	
CONNECTORS	
J39 - EXTERNAL AIRFLOW SENSOR	J33 - BUTTERFLY VALVE INPUT SIGNAL
J38 - PARTICLE COUNTER SENSOR	J32 - BUTTERFLY VALVE RELAY
J37 - TEMPERATURE/PRESSURE/HUMIDITY	J36 - ALARM RELAY SIGNAL
J28 - RS485 COMMPORT	J35 - BOOSTER FAN POWER RELAY
J40 - TO CONTROL BOARD CONNECTION	J34 - BOOSTER FAN 0-10V SIGNAL
J31 - VHP	

ACCESSORIES
(OPTIONALS PARTS)

Revision	Date	Draw by	Date:	10/06/2022	Project description:	Situation:	Zone:	Drawings: 10 / 11	Previous: 9	Next: 11
			Draw by:	X. Cañellas	BIOPTIMA EXA	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: ACCESSORIES BOARD OPTIONAL PARTS		Drawing number: 1917194	
			Aproved:	X. Cañellas				Project reference: 977083	Drawing number 10	

Date 10/06/2022 Bill of Materials Page 1

BIOPTIMA EXA

AZBIL TELSTAR
TECHNOLOGIES, S.L.

Symbol	Drawing	Column	Description	Refer. Numb.	Catalogue	Reference
-FL1	02	B	POWER LINE FILTER EMC 20A YANBIXINKEJI	YB12E4-20A-Q	---	634174
-F1	02	B	SLOW BLOW FUSE 16A SCHURTER	0001.2516.T	---	616733
-I1	02	B	SINGLE PHASE LINE INDUCTOR POLYLUX	R12.5	---	644186
-PS1	02	E	POWER SUPPLY 27VDC 200W MEANWELL	RSP-200-27	---	620899
-BT1	02	F	BATERIA 12V 1,3AH BAOTE	BT-12M1.3AT	---	636420
-BT2	02	G	BATERIA 12V 1,3AH BAOTE	BT-12M1.3AT	---	636420
-PCB1	03	C	ATLAS BSC CONTROL CARD	---	---	632679
-PCB2	03	F	ATLAS ACCESSORIES EXPANSION CARD	---	---	627589
-FAN1	04	C	EBM EC FAN	K3G225RE1904	---	1629383
-FAN2	04	E	EBM EC FAN	K3G225RE1904	---	1629383
-FAN3	04	F	EBM EC FAN	K3G225RE1904	---	1629383
-LED1	05	C	LED STRIP (BIOPTIMA EXA 4)	---	---	1741339
			LED STRIP (BIOPTIMA EXA 6)	---	---	1741347
-BL1	05	G	BALAST T5	ACT2S39PSUNV	---	637566
-L1	05	G	UV TUBE PHILIPS TUV 36W	TUV 36T5 4P SE UNP	---	32073
-M1	06	F	DC TUBULAR MOTOR 24VDC	FT45DS-20/8	---	651474
-SK1	06	C	ELECTRIC SOCKET	B-012H	---	510032
-SK2	06	D	ELECTRIC SOCKET	B-012H	---	510032
-SV1	06	F	ELECTROVALVE 24VDC 8W	266293 TYPE 6013	---	640621
-IS1	07	C	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS2	07	C	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS3	07	C	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS4	07	G	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-IS5	07	G	INDUCTIVE SENSOR OMRON	TL-W5MB1	---	64287
-SP1	08	C	AIR VELOCITY SENSOR 0-10m/s 0-5V DEGREE	FS22503	---	621613
-SP2	08	E	AIR VELOCITY SENSOR 0-1m/s 0-5V DEGREE	FS22500	---	614060
-DIS1	09	D	DISPLAY PARA BSC ATLAS	---	---	632700
-WB1	09	F	WINDOW BUTTONS ATLAS	---	---	1384557

Revision	Date	Draw by	Date:	10/06/2022	Project description:	Situation:	Zone:	Drawings: 11 / 11	Previous: 10	Next:
			Draw by:	X. Cañellas	BIOPTIMA EXA	Customer: AZBIL TELSTAR TECHNOLOGIES, S.L.	Drawing description: BILL OF MATERIALS		Drawing number: 1917194	
			Aproved:	X. Cañellas				Project reference: 977083	Drawing number 11	

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EN

1. PURPOSE

Check the integrity of the downflow filters.

2. SCOPE

- *BiOptima Exa*

3. TOOLS

- *Aerosol generator (cool or hot).*
- ATI PAO-4, *Shell Ondina fluid or Emery 3004 oils.*
- *Photometer (ATI 2H).*

4. METHODOLOGY

*Introduce fume by the generator into the plenum on the top of the cabinet centered on the upper prefilter (Figure 1). Concentration in cabinet u*sually within 20-70µg/l.

- *Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Figure 2)*


Fig. 1

Fig. 2

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EN

- *With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters (Fig. 3a), or CG screen (Fig. 3b), including its perimeter, by passing the photometer probe in slightly overlapping strokes over the entire surface.*

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EN

- *Also scan the entire periphery of the filter and the junction between filter and filter mounting frame at a scanning rate that does not exceed 5 cm/s.*


Fig. 3a

Fig. 3b

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EN ***5. ACCEPTANCE CRITERIA***

- *The photometer signal after the filter should be smaller or equal to a 0,01% of the specific concentration before the filter.*
- *The concentration value around the DOP port should be smaller or equal to the concentration of environment.*

Fig. 4. Resultado correcto: no fuga
Correct result: filter not leaking

Fig. 5. Resultado incorrecto: fuga
Wrong result: filter leaking

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EN

1. PURPOSE

- *Prove the **not** existence of leaks in the HEPA inflow filters.*

2. SCOPE

- *BiOptima Exa*

3. TOOLS

- *Aerosol generator (cool or hot).*
- ATI PAO-4, *Shell Ondina fluid or Emery 3004 oils.*
- *Photometer (ATI 2H).*

4. METHODOLOGY

- *The DOP photometer sample port should be placed in the exhaust duct. Ideally, at 10 times the duct diameter (3 meters) from the first stage filtering to insure an adequate smoke mixing. Photometer tube must be inserted inside the duct in opposite direction of the flow.*

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EN

Introduce fume by the generator into the cabinet using soft tube. Concentration in cabinet usually within 20-70µg/l.

- *Introduce the smoke to cover each plenum of V filters, all kinds of corners around the V-filters, seams above the V-Filters fixing plate, seams around the side glasses and sockets, etc. (Figure 3).*
- *Check the photometer display reading during the whole process (Figure 4).*

Fig. 3

Fig. 4

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EN ***5. ACCEPTANCE CRITERIA***

- *The photometer signal after the filter should be smaller or equal to a 0.01% of the specific concentration before the filter.*
- *The concentration value around the DOP port should be smaller or equal to the concentration of environment.*

Fig. 5. Resultado correcto: no fuga
Correct result: filter not leaking

Fig. 6. Resultado incorrecto: fuga
Wrong result: filter leaking

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$V_{promedio}$ (m/s)= 1/n x (V_1 + V_2 + V_3 + V_n)

EN *1. PURPOSE*

- *Check that the speed control in different points of the working area is suitable to achieve product protection. To determine the speed average as well as the uniformity range in the laminar flow zone.*

2. SCOPE

- *BiOptima Exa*

3. TOOLS

- *Omnidirectional or unidirectional anemometer.*

4. METHODOLOGY

- *Inside the considered working area according to each cabinet model divide the working area in squares (see attached table).*

Then switch on the cabinet and wait until the velocity and the temperature becomes stable.


- *Place the anemometer at 270mm high of the shelves, over each square (take all readings with the anemometer in the same direction and without covering any holes, whenever possible).*


- *Take note of the readings and the position of the anemometer in each point.*

- *For this calculation, we proceed to the following formula:*

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EN *5. ACCEPTANCE CRITERIA*

- *The average of the airflow velocity at the work surface should be ≥ 0.34 m/s and ≤ 0.36 m/s*
- *Every single value should not differ more than ± 20% from the average of the airflow velocity.*

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EN *Annex A9.4.1. SET POINT ALARM DOWNFLOW AIR VELOCITY TEST*

EN *1. PURPOSE*

- *Check correct work of the alarms when flow speed is lower or higher than the admissible.*

2. SCOPE

- *BiOptima Exa*

3. TOOLS

- *Omnidirectional or unidirectional anemometer*

4. METHODOLOGY

- *Accede to Technical Service menu through Password Petition screen (User's Manual chapter. 4.5 Point 1), using for that the random security password that you are asked for.*
- *Move to Inflow Test Set Points screen and wait some seconds until the stabilization of the cabinet.*
- *Increase slowly downflow velocity through function button*  *until "Excessive Laminar Flow" alarm activated. Then take note of the display velocity.*
- Taking the same points as in the previous test *(Annex A9.4)*, measure and take note of air velocity of all points and calculate the average.
- *Decrease slowly downflow velocity through the function button*  *until "Insufficient Laminar Flow" activated. Then take note of the display velocity.*
- Taking the same points as in the previous test *(Annex A9.4)*, measure and take note of air velocity of all points and calculate the average.

EN *5. ACCEPTANCE CRITERIA*

All below criteria should be satisfied:

- (Displayed down flow velocity at high alarm) -0.02 ≤ (Measured average down flow velocity at high alarm) ≤ (Displayed down flow velocity at high alarm) +0.02
- (Down flow velocity at low alarm) -0.02 ≤ (Measured average down flow velocity at low alarm) ≤ (Displayed down flow velocity at low alarm) +0.02
- No individual measurement (max. or min. point) should differ from its average more than 20%.

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EN ***1. PURPOSE***

- *The purpose of this test is to calculate the inflow velocity through the front aperture of the equipment to obtain the Class II to the operator protection.*

2. SCOPE


- *BiOptima Exa*

3. TOOLS

- *Anemometer with a range till 7 m/s.*
- *Pressure gauge.*


4. METHODOLOGY

- *Switch on the cabinet and wait until the velocity and the temperature are stable.*
- *Use the anemometer to measure 9 points divided in 3 horizontal rows.*

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
- *Probe must be held perpendicular to the direction of the airflow. This can be checked by moving the probe horizontally and vertically. The highest reading shows that the probe is perpendicular. You should also be helped by smoke.*
- *First row (A) at 10 mm below the top of the working opening.*

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

- *Second row at the centre of the working opening.*

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

- *Once you get the 9 values, calculate the average per row and overall average.*
- *Take note of the display voltage of the exhaust fan and the downflow fans.*
- *Measure the exhaust plenum pressure.*
- *Measure the downflow plenum pressure (left port + right port / 2).*

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 BiOptima Exa 6: 86±20 Pa

EN *5. ACCEPTANCE CRITERIA*

- *The average inflow velocity must fulfill:*

 BiOptima Exa 4: 0.6 -0.02/+0.04 m/s

 BiOptima Exa 6: 0.6 -0.02/+0.04 m/s

- *The voltage of the downflow fan must fulfill:*

 BiOptima Exa 4: 5.55±0.5 V

 BiOptima Exa 6: 5.55±0.5 V

- *The downflow plenum pressure must fulfill:*

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EN ***1. PURPOSE***

- *Check the correct operation of inflow air velocity alarm when flow is lower or higher than the admissible.*

2. SCOPE


- *BiOptima Exa*

3. TOOLS

- *Anemometer with a range till 7 m/s.*


4. METHODOLOGY

- *This test has to be performed after calibration of the cabinet.*
- *Accede to Technical Service menu through Password Petition screen (User's Manual chapter. 4.5 Point 1), using for that the random security password that you are asked for.*
- *Move to Inflow Set Points Test screen and wait some seconds until the stabilization of the cabinet.*
- *Decrease slowly downflow velocity through the function button*  *until "Insufficient inflow" activated. Take note:*
 - *Anemometer velocity.*
 - *Display velocity.*
- *Use the anemometer to measure 9 points divided in 3 horizontal rows.*

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



- *Once you get the 9 values, calculate the average per row and overall average.*
- *The same inflow velocity average must be calculated for the alarm "Excessive inflow". Repeat the process increasing slowly the inflow velocity through the function button*  *until the alarm is activated.*

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EN 5. ACCEPTANCE CRITERIA

"Insufficient inflow" alarm:

- *The inflow average value should be the same as shown in the display ± 0.05 (measuring error).*

 Display v.-0.05 < Real inflow < Display v.+0.05

- *Display entry velocity ≥ 0.55 m/s as soon as the alarm is activated.*

"Excessive inflow" alarm:

- *The inflow average value should be the same as shown in the display ± 0.05 (measuring error).*

 Display v.-0.05 < Real inflow < Display v.+0.05

- *Display entry velocity ≤ 0.72 m/s as soon as the alarm is activated.*

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EN Annex A9.6. SMOKE TEST

EN 1. OBJECT

- *Verify the cabinet for personnel and product protection.*

2. SCOPE

- *BiOptima Exa*

3. TOOLS

- *Smoke generator*

4. METHODOLOGY

- *Product protection: smoke is produced 5-10 cm from the protective glass (Fig. 1 and 2). The behaviour of the smoke is checked visually.*
- *Smoke is produced inside the cabinet, at about 5-10 cm over the lower edge of the glass (Fig. 3). The behaviour of the smoke is checked visually.*
- *User protection: smoke is produced inside the cabinet, as close as possible to the lower edge of the glass (Fig. 4). The behaviour of the smoke is checked visually.*


Fig. 1


Fig. 2

Fig. 3

Fig. 4

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EN ***5. ACCEPTATION CRITERIA***

Product protection:

- *The generated smoke in the outside of the front aperture will be dragged to the inflow area. No fume sign has to reach the working surface.*
- *Turbulences in the working area shall not be created. The smoke generated should be dragged by the laminar flow to the surface of the working area.*

User protection:

- *No smoke shall be refluxed out of the cabinet once drawn in, guaranteeing the protection of the operator and environment.*

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
UK DIRECTIVES	DECLARATION
Supply of Machinery (Safety) Regulations 2008 ***(S.I. 2008/1597)*** *Electromagnetic Compatibility Regulations 2016* ***(S.I. 2016/1091)*** *Electrical Equipment (Safety) Regulations 2016* ***(S.I. 2016/1101)*** *The Restriction of the use of Certain Hazardous substances in Electrical and Electronic Equipment Regulations 2012* ***(S.I. 2012/3032)*** Specific standards: *Safety requirements for electrical equipment for measurement, control, and laboratory use* ***(BS EN 61010:2010+A1)*** *Performance criteria for microbiological safety cabinets* ***(BS EN 12469:2000)***	The present declaration of conformity is issued under the sole responsibility of Manufacturer. This document hereby certifies that the above equipment complies with the UK Directives listed, and their amendments. The equipment must be installed, tested, and operated in strict accordance with Manufacturer instructions and used by fully trained operators. Only approved Azbil Telstar Technologies S.L.U. spares must be used for maintenance and repair and fitted in accordance with Manufacturer instructions. UK CA
Body Retaining Technical File:	Retained by AZBIL TELSTAR TECHNOLOGIES, S.L.U.

Signed by:	
Name:	Jordi Valls
Title:	Global Quality Manager
On Behalf of:	AZBIL TELSTAR TECHNOLOGIES, S.L.U.
Date:	December 2022

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