Instruction manual

VWR® POLARIZING MICROSCOPE, VISISCOPE® Series 100POL

ML113 POL
BL113 POL

EU cat. no 630-2977
630-2976

Legal address of Manufacturer

Europe

VWR International bv,
Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001
Leuven
+3216385011
be.vwr.com

UK Importer:
VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
uk.vwr.com

Country of Origin

ITALY

Table of Contents

Warning 4
Safety Information 4
Package Contents 5
ML113 POL 5
BL113 POL 5
Unpacking 6
Installation 7
Microscope assembling 7
Intended use 8
Symbols and conventions 8
Instrument description 9
General Operating Instructions 10
Use of the microscope 10
Light intensity adjustment 10
Coarse focus tension adjustment 10
Stage 10
Adjust the interpupillary distance 10
Diopter adjustment 11
Aperture diaphragm 11
Use with rechargeable batteries 11
Use in brightfield 12
Use in polarized light 12
Troubleshooting 13
Repair and maintenance 14
Technical service 14
Warranty 14
Compliance with local laws and regulations 14
Equipment disposal 15
Appendix 16
Local VWR offices in Europe and Asia Pacific 20

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.
We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.
We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.
Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Package Contents

ML113 POL

① Frame
② Monocular observation head
③ Eyepiece
④ Tension adjustment tool
⑤ Objectives (4X / 10X / 40X)
⑥ Dust cover
⑦ Green filter
⑧ Power supply

BL113 POL

① Frame
② Binocular observation head
③ Eyepiece
④ Tension adjustment tool
⑤ Objectives (4X / 10X / 40X)
⑥ Dust cover
⑦ Green filter
⑧ Power supply

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

Installation

Microscope assembling

1. Remove the dust cap from the microscope frame and from the bottom of the observation head.
2. Insert the optical head above the stand and tighten the screw. (Fig. 1)

- **Hold the head with one hand during the locking in order to avoid that the head falls.**

Fig. 1

3. Insert both eyepieces into the tubes of the optical head. (Fig. 2)

- **For model ML113 POL: only one eyepiece is used as this model is monocular.**

Fig. 2

4. Insert the power supply jack in the socket placed at the rear side of the microscope. (Fig. 3)

Fig. 3

Intended use


For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

	CAUTION This symbol indicate the possibility of electrical shock
	CAUTION This symbol indicates a potential risk and alerts you to proceed with caution

Instrument description

General Operating Instructions

Use of the microscope

Light intensity adjustment

Operate on the light intensity adjustment dial ① to turn ON / OFF the microscope and to increase / decrease the illumination voltage. (Fig. 4)

Fig. 4

Coarse focus tension adjustment

- **Adjust the tension using the provided tool**.

The coarse knob tension is pre-setted in the factory.
To modify the tension according to personal's needs, rotate the ring ② using the provided tool. (Fig. 5)
Clockwise rotation increases the tension.

- If the tension is too loose, the stage could drop by itself or the focus can be easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

Fig. 5

Stage

Stage is equipped with two clips ③ for the sample placement. Standard specimen is glass slide, thickness 1,2 mm with coverslide 0,17 mm. (Fig. 6)
Specimen can be moved with hands.

1. Lift the clip, place the specimen and gently release the clip until the specimen is secured over the stage.(Fig. 6)
2. When used in polarized light stage can be rotated on 360°, be grabbing the knurled ring ④ placed on the outer side of the stage.


Fig. 6

Adjust the interpupillary distance

Only BL113 POL

Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig. 7)

- The graduation on the interpupillary distance indicator ⑤, pointed by the spot "." on the eyepiece holder, shows the distance between the operator's eyes.

The range of the interpupillary distance is 48- 75 mm.

Fig. 7

Diopter adjustment

Only BL113 POL

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptric adjustment ring ① to compensate. (Fig. 8)

- **The adjustment range is ±5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptric correction.**

Fig. 8

Aperture diaphragm

The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
Move the diaphragm lever ② (Fig. 9) toward left or right to decrease or increase the N.A. value.

Fig. 9

With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 10.


Fig. 10

Use with rechargeable batteries

When the microscope is plugged with the power supply, the LED indicator for the battery recharge ③ is lit. (Fig. 11)

- LED red: battery under charge
- LED green: battery fully charged.

When the microscope is unplugged, the LED is off.

During the normal use with batteries, LED is always OFF.


Fig. 11

Use in brightfield

Move the analyzer slider ① to the right to remove analyzer filter from the light path. (Fig. 12)
In this way it is possible to work in brightfield.


Fig. 12

Use in polarized light

1. Move the analyzer slider ① to the left to insert analyzer filter into the light path. (Fig. 12)
2. Insert 10X objective into the light path.
3. Remove the specimen from the stage.
4. Rotate the polarizer filter ② (Fig. 13) until a complete dark filed of view can be achieved. Now the "crossed-Nicol" position is obtained and it is possible work in polarized light.
5. Put a specimen on the stage.
6. Insert the desired objective.
7. Focus the specimen.
8. Begin the observation.



Fig. 13

- Polarizer can be removed from the light path simply swinging the polarizer unit toward right. (Fig. 14)

Fig. 14

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark	Power supply is unplugged	Connect
	Brightness is too low	Set brightness to a proper level
	Batteries are uncharged	Fully charge the batteries
Dirt or dust is visible in the field of view	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
Visibility is poor • Image is not good • Contrast is poor • Details are indistinct • Image glares	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or to open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	Polarizing filters are not in "crossed-Nicol" position	Set them until "crossed-Nicol" position can be achieved
One side of the image is out of focus	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Mechanical Section:		
The coarse focus knob is hard to turn	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric section:		
The LED doesn't turn on	No power supply	Check the power cord connection
	Batteries are uncharged	Fully charge the batteries
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
IV. Observation tube:		
Field of view of one eye does not match that of the other	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope

Repair and maintenance

Microscopy environment
This microscope is recommended for use in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

The microscope should always be kept vertical when moving it so that no moving parts, such as the eyepieces, fall out.
Never mishandle or impose unnecessary force on the microscope.
Never attempt to service the microscope yourself.
After use, turn down the illumination intensity control and turn the light off. Cover the microscope with included the dust cover, and keep it in a dry and clean place.

Electrical safety precautions

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
Users should observe all safety regulations of the region. The equipment has acquired the CE marking label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics
If the optical parts require cleaning first use compressed air.
If that is not sufficient use a soft lint-free cloth with water and a mild detergent.
And as a final option use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.

- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**

Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics. Do not disassemble objectives or eyepieces attempting to clean them.
For the best results, use the VWR cleaning kit (630-2081).
If you need to send the microscope to manufacturer for maintenance, please use the original packaging.

Technical service

Web Resources

Visit the VWR website at www.vwr.com for:

- **Complete technical service contact information.**
- **Access to the VWR Online Catalogue, and information about accessories and related products.**
- **Additional product information and special offers.**

Contact us **For information or technical assistance contact your local VWR representative or visit. www.vwr.com.**

Warranty


VWR warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you.

Appendix

Memo

Memo

Memo

Local VWR offices in Europe and Asia Pacific

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Ireland
VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
sales.ie@vwr.com

Italy
VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
info.it@vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
info.nl@vwr.com

Norway
VWR International AS
Brynsalleen 4,
0667 Oslo
Tel.: +47 22 90 00 00
info.no@vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
info.pl@vwr.com

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India
Avantor Performance Materials India
Private Limited
19th Floor, Building No. 5, Tower C,
DLF Cyber City, Phase- III,
Gurugram - 122002, Haryana, India
Tel.: +91-124-4656700
csindia@avantorsciences.com

Korea
Avantor Performance Materials Korea Ltd
2F ACE Gwanggyo Tower I, Daehak 4ro 17
Yeongtong-gu Suwon, Korea 16226
Tel.: +82 31 645 7256
saleskorea@avantorsciences.com

Middle East & Africa
VWR International FZ-LLC
Office 203, DSP Lab Complex,
Dubai Science Park,
Dubai, United Arab Emirates
Tel: +971 4 5573271
info.mea@vwr.com

Singapore
VWR Singapore Pte Ltd
The Metropolis
Tower 1, #05-03
9 North Buona Vista Drive
Singapore 138588
Tel: +65 6505 0760
sales.sg@vwr.com

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