RFM-T Series Refractometer

Bellingham + Stanley, a Xylem brand, Xylem Analytics Germany GmbH, AM Achalaich 11, 82362 Weilheim, Germany
Email: sales.bs.uk@xylem.com

QUICK START GUIDE
19-406'05 REV A, OCT-2023

Thank you for purchasing this instrument. In order to ensure that this product provides many years of service please follow the guidance in this document. To register your instrument with the manufacturer and print a warranty certificate visit:
www.bellinghamandstanley.com

Unpacking

Carefully remove all of the packing material. It is recommended that the box and other packing materials are retained so that, should the need arise, the instrument can be safely returned to the supplier for service. Check that all parts listed below are present and that no transit damage has occurred. If any are damaged or missing, contact the supplier immediately.

Contents list

Quantity	Item	Code
1	RFM-T Series Refractometer	See below
1	Mains lead set	55-109
1	Power supply	55-105
1	Quick Start Guides	19-406
2	Card type ID tag	37-530
2	Spare air vent filter	22-482
1	Touchscreen protector	19-204

RFM-T Refractometer code numbers:

Model	Complete refractometer	Refractometer module only
RFM330-T	19-30	19-330
RFM340-T	19-40	19-340
RFM960-T	19-60	19-960
RFM970-T	19-70	19-970
RFM960-TF	19-66	19-966
RFM970-TF	19-76	19-976
RFM990-T AUS32	19-73	19-973

Positioning the system

- Dry and indoors.
- Capable of holding the instrument's weight.
- Away from draughty or hot equipment like fans or heaters.
- Out of direct sunlight or strong ambient light.
- Away from potential sources of interference, such as RFI generating equipment.
- Within reach of a power point.
- Not using a power circuit that also has large motors or noise generating equipment connected to it.

Do not block the air vents.

Mains connection

The power supply adapter is supplied with a moulded mains cord and plug to suit one of several socket types. For UK leads, replace the fuse only with the type indicated on the plug.

- Voltage 110 to 230 V ~ ±10%
- Maximum current 2 A
- Frequency 50 to 60 Hz

RISK OF ELECTRIC SHOCK:

- For indoor use only.
- Must be kept dry.
- Disconnect the equipment from the mains supply before unplugging the mains lead from the instrument.
- Do not cover, designed to operate with free air convection.
- No cleaning required.
- Do not jet wash

Switching on and off

Plug in the power supply and switch on the mains supply, after startup the instrument will show a loading screen indicating the software version. Once loaded, the instrument will be in standby mode. To switch on, swipe a finger or stylus across the display. After a few moments the instrument will show the Measurement Display (or the Setup Wizard if switching on for the first time).


To switch to standby mode press and hold the Menu Button for several seconds and select Yes on the confirmation screen. It is then safe to turn off the mains supply and disconnect the power lead.

Manoeuvring through menus

The instrument features a capacitive, high-resolution touchscreen display. On-screen buttons can be easily used by simply touching the display. Use the supplied stylus and screen protector to keep the touchscreen clean and avoid damage from improper use. The touchscreen should only be used with a finger or stylus designed for use with a capacitive touchscreen.

Alternatively, a USB keyboard, available from Bellingham + Stanley, can be used to select items and navigate menus.

Many buttons feature a quick select option, indicated by the two arrows on the right of the button. Simply touch either arrow to cycle through the possible settings.

When in a menu, the Quit and Home buttons become available. Press the Quit Button to return to the previous menu. Any changes made will be saved. The Home Button will exit to the Measurement Display, again any changes made will be saved.

Setup Wizard


When switching on for the first time, after the instrument has completed its startup procedures, the Setup Wizard will appear. First select the desired language. A list of different options will then appear:

Copy user guides and PC software to USB: Makes copies of the full User Guides and PC software compatible with the instrument to a removable USB storage device. Insert a USB removable storage device before selecting. Full User Guides provide detailed instructions on instrument operation.
Set clock: Set the instrument to local time and set the date format.
Set reading mode: Make changes to the measurement settings.
Set record settings: Specify how reading results will be printed or saved.
Set user settings: Alter the security settings to a level suitable for the area of operation.

Select each in turn and alter the settings as desired, then press Quit to return to the Setup Wizard Menu. Each menu will be marked by a tick after the settings are altered. Once all required settings are changed press Quit on the Setup Wizard Menu to show the Measurement Display.

Auto Zero

When the instrument is turned on a Zero calibration will be requested. This is to ensure that readings are as accurate as possible. Whilst the temperature stabilises the Reading Result section of the screen will display a warning message.

Follow the onscreen instructions – ensure that the prism is thoroughly clean and that the presser is then lowered.

The instrument stores the details of the last sample used for Zero calibration – typically distilled water. If a different sample is to be used press the Change Sample button and enter the data as requested, the sample value must be entered in the scale used for calibration.

Once the sample details are correct apply the Zero sample to the instrument prism. The complete prism surface must be covered, 1ml is generally sufficient to achieve this. Lower the presser and press ok. Continuous readings will then be taken until the sample stabilises The Zero point will then be set.

This symbol is an internationally agreed indicator that the product bearing it should not be disposed of as general waste or garbage which might end up in landfill sites, but should instead be sent for special processing and/or recycling in those countries where appropriate legislation and facilities are in place.

Full User Guides can be downloaded from the instrument's Help Menu or Setup Wizard

Caution or warning.

Electrical hazard.

RFM-T Series Refractometer


QUICK START GUIDE

19-406'05 REV A, OCT-2023

Bellingham + Stanley, a Xylem brand, Xylem Analytics Germany GmbH, AM Achalaich 11, 82362 Weilheim, Germany
Email: sales.bs.uk@xylem.com

Measurement Display

From the Measurement Display readings can be taken, viewed and saved. Quick access to all the frequently used features of the instrument is also available.

1. Reading result: After taking a reading, the final result will be displayed here. Warning and alert messages will also appear here. For example, after the temperature set-point is changed, a Zero calibration will be requested.
2. Reading Configuration Panel: Displays the measurement settings for the currently displayed reading or, if a result is not currently displayed, the settings to be used for the next reading. It also allows for quick changes to the Basic reading settings.
3. Mode button: Gives access to edit all Basic and Advanced reading settings.
4. Read button: Initiates a reading. The complete prism surface must be covered, 1ml is generally sufficient to achieve this
5. Record button: Pushing the button will record the result to the Reading Log or printer depending on the print settings.
6. Zero button: Initiates a Zero calibration routine. A Zero can be performed with any sample. Follow on screen instructions.


Holding down the above buttons for several seconds will open an associated Menu.

7. Menu button: Opens the instrument's Main Menu from where all instrument settings can be changed. Hold for several seconds to switch the instrument to standby.
8. Data button: Opens the Data Menu, this includes viewers for saved results and the System Log as well as the ability to download the results and any saved PDFs.

Taking a reading

- Raise the presser and thoroughly clean the prism surface.
- Place a small amount of the sample to be measured on the prism. The complete prism surface must be covered, 1ml is generally sufficient.
- Lower the presser.
- Readings will start or press the Read button.

The instrument will display a progress bar. The reading will then be displayed using the scales and temperature compensations specified in the Mode Menu or current Method as shown on the Reading Configuration Panel.

Menu structure

- **1. Calibration**
 - 1. Zero
 - 2. Span (Top Calibration)
 - 3. Undo
 - 4. Report
- **2. Data**
 - 1. Saved Results
 - 2. System Log
 - 3. Copy Files to USB
 - 4. Move Files to USB
 - 5. Export Logs to USB
- **3. Setup**
 - 1. Setup Wizard
 - 2. Methods
 - N. Add new…
 - S. Settings…
 - 1. Use Methods
 - 2. Hide preset
 - 3. Selection Type
 - 3. Users
 - N. Add new…
 - 1. administrator
 - 2. calibrator
 - 3. operator
 - S. Settings…
 - 1. Protection
 - 1. basic 5. cfr
 - 2. low
 - 3. medium
 - 4. high
 - 2. User Types
 - N. Add new…
 - 1. administrator
 - 2. calibrator
 - 3. operator
 - 4. Reading
 - 1. Recurrence
 - 1. single
 - 2. continuous
 - 3. auto single
 - 2. Appl. Correction
 - 1. none
 - 2. offset
 - 3. citric acid
 - 3. Scales
 - N. Add new…
 - 1. From Library
 - 2. From Data
 - 4. Temp. Comps.
 - N. Add new…
 - 1. From Data
 - 5.Recording
 - 1. Mode
 - 1. none 5. statistics
 - 2. single 6. multi
 - 3. save
 - 4. single + save
 - 2. Printer
 - 1. none 5. network
 - 2. serial
 - 3. usb
 - 4. pdf
 - 3. LIMS Output
 - S. Select Fields
 - 1. none
 - 2. csv
 - 3. xml
 - 1. serial
 - 2. network
 - 3. usb file
 - 4. ftp file
 - 4. Record After Read
 - 1. no
 - 2. yes
 - 5. Traceability
 - 1. Location
 - 2. Batch Codes
 - 1. off
 - 2. number
 - 3. text
 - 4. dateinc
 - 3. Auto Increment
 - 1. off
 - 2. on
 - 4. Log Users
 - 1. off
 - 2. number
 - 3. text
 - 4. login
 - 5. Remember Text
 - 1. none
 - 2. batch
 - 3. user
 - 4. batch & user
 - 6. Ports
 - 1. Serial Port
 - 1. Baud Rate
 - 2. Word Length
 - 2. Network
 - 1. DHCP
 - 2. IP Address
 - 3. Subnet Mask
 - 4. Default Gateway
 - 6. System
 - 1. Language
 - 2. Remote Display
 - 1. off
 - 2. on
 - 3. Display Brightness
 - 4. Screensaver
 - 5. Time / Date
 - 1. Date Format
 - 1. dd/mm/yy
 - 2. mm/dd/yy
 - 2. Summer Time
 - 1. no
 - 2. yes
 - 3. Time Server
 - 4. Set Clock…
 - 4. UTC Offset
 - 6. Temp Control
 - 1. Mode
 - 1. off
 - 2. on
 - 2. Disable Presser
 - 1. no
 - 2. yes
 - 7. Clone From Backup
- **4. Maintenance**
 - 1. Report
 - 2. Filter Change Reminder
 - 1. New Filters Fitted
 - 2. Frequency
 - 3. Last Changed
 - 4. Status
 - 3. Backup Settings
 - 1. Backup Now
 - 2. Schedule
 - 3. Destination
 - 4. Last backup
 - 4. Clean Screen
 - 5. Update Software
- **5. Help**
 - 1. Information
 - 2. Contact
 - 3. User Guide
 - 4. PC software

This symbol is an internationally agreed indicator that the product bearing it should not be disposed of as general waste or garbage which might end up in landfill sites, but should instead be sent for special processing and/or recycling in those countries where appropriate legislation and facilities are in place.

Full User Guides can be downloaded from the instrument's Help Menu or Setup Wizard

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