USER'S MANUAL
MANUAL DE USUARIO
MANUEL D'UTILISATION

Vertical Laminar Flow Benches
Cabinas de Flujo Laminar Vertical
Hottes à flux laminaire vertical

AEOLUS V

Technology. Science. Life.

USER'S MANUAL

Vertical Laminar Flow Benches

AEOLUS V

Technology. Science. Life.

Translation of the original manual

-

Failure to read, understand and follow the instructions in this manual may result in damage to the unit, injury to operating personnel and poor equipment performance.

Our company reserves the right to update, improve, amend or otherwise modify the design or specifications of any product at any time without prior notice.

INDEX

0. SYMBOLS AND ABBREVIATIONS ____ 4

1. GENERAL DATA ____ 5

2. DESCRIPTION OF THE CABINET ____ 6
- 2.1. Working principle ____ 6
- 2.2. Applications ____ 7
- 2.3. Standard equipment ____ 7
- 2.4. Certification ____ 7
- 2.5. Technical features ____ 8

3. INSTALLING INSTRUCTIONS ____ 9
- 3.1. Transport ____ 9
- 3.2. Placing on site ____ 9
- 3.3. Precautions for first connection ____ 10

4. ELECTRICAL CONNECTION AND SERVICES ____ 12
- 4.1. Installation tests ____ 13
- 4.2. Stopping for long periods ____ 13

5. INSTRUCTIONS FOR USE ____ 14
- 5.1. General recommendations ____ 14
- 5.2. Ergonomics ____ 15
- 5.3. Definition of the working area ____ 16
- 5.4. Starting up ____ 17
- 5.5. Switching off ____ 17
- 5.6. Front glass positions ____ 17
- 5.7. Control panel ____ 19
 - 5.7.1. Components description ____ 19
 - 5.7.2. Working mode menu functions ____ 19
 - 5.7.3. Signal lights & alarms ____ 20
 - 5.7.4. User configuration menu functions ____ 21
- 5.8. Optional accessories ____ 23

6. CLEANING AND DISINFECTING ____ 24
- 6.1. Front Glass cleaning and disinfection ____ 25

7. MAINTENANCE INSTRUCTIONS ____ 26
- 7.1. Regular validation and control ____ 26
- 7.2. Maintenance Table ____ 26
- 7.3. Absolute filters ____ 27
- 7.4. Prefilters ____ 27

8. LIST OF SPARE PARTS ____ 28

9. TROUBLESHOOTING ____ 29

10. MAINTENANCE CHART ____ 31

11. ANNEXES ____ 32
- 11.1. Electrical diagram 120V / 230V 50/60 Hz ____ 32
- 11.2. Impulsion filter integrity test ____ 38
- 11.3. Air velocity test ____ 40
- 11.4. Smoke pattern test ____ 43
- 11.5. Declarations of conformity ____ 44

Translation of the original manual

0. SYMBOLS AND ABBREVIATIONS

List of stickers used in the Aeolus Cabinets.

	"EARTH WIRE PROTECTED" STICKER.
	"FUNCTIONAL EARTH WIRE" STICKER.
	"ATTENTION, ELECTRIC SHOCK RISK" STICKER.
	"SOCKET VOLTAGE" STICKER.
	"SOCKET MAXIMUM CURRENT" STICKER.

Abbreviations used in this manual:

V	Volts
A	Amperes
W, kW, µW	Watts, kilowatts, microwatts
mm, cm, m	Millimeters, centimeters, meters
kg	Kilogram

s, h	Seconds, hours
Hz	Hertz
dB	Decibel
ºC	Celsius
LFB	Laminar Flow Bench

1. GENERAL DATA

Operating Instructions Manual for vertical laminar flow sterile benches, which ensures high level protection for the sample:


Models: **Aeolus V 3** **Aeolus V 4** **Aeolus V 5** **Aeolus V 6**

WARNING

This Manual must be read and understood before its installation, maintenance, and operation of the LFB. All the safety instructions given in the Manual must be followed.
Any use or application failing to comply with the recommendations expressed in the Manual may cause danger and loss of rights to manufacturing guarantee as well as freeing the manufacturer from all responsibility arising.
Any maintenance must be done by qualified staff according to the laboratory working correct conditions. In case these actions are not carefully carried out, as a result the LFB can work in abnormal conditions.
All maintenance actions on this equipment must be performed by a qualified technician who is familiar with the proper maintenance procedure required for this equipment.

The European Union issued the Directive 2012/19/EU on Waste of Electrical and Electronic Equipment (WEEE).

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive. Mark symbol is shown at the right.
Syntegon Telstar sells products through distributors throughout Europe. Contact your local distributor for recycling/disposal.

The name plate & CE mark sticker is located at one side of the cabinet.

MANUFACTURER: **SYNTEGON TELSTAR TECHNOLOGIES, S.L.U.**
Av. Font i Sagué, 55
Parc Científic i Tecnològic Orbital 40
08227 TERRASSA (SPAIN)
Tel (+34) 93 736 16 00
Fax (+34) 93 786 13 80
E-mail**:** telstar@telstar.com
http://www.telstar.com

2. DESCRIPTION OF THE CABINET


2.1. WORKING PRINCIPLE

The vertical LFB are equipped with high efficiency centrifugal EC fans (flow measured by a thermoanemometric probe) and designed with the most advanced microprocessor control, ensuring a high level of protection for the sample.

Exterior of cabinet is made of laminated, stove-enamel coated steel. Functional design with reduced external dimensions and front access for changing filters.

Spacious internal chamber with sliding front made from tempered and laminated safety glass. Single work surface made of polished stainless steel and liquid collection tray of polyester coated steel.

The pre-filtered air (1) is blown by the fans (2) into chamber/plenum (3) where it is filtered through the HEPA absolute filter (4) and then arrives at the working area (5) providing laminar flow and clean air in the working area.

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2.2. APPLICATIONS

This model of vertical laminar flow sterile bench is designed for handling samples or products when the quality of sterile air must be guaranteed in the working area, provided that no protection to the user or the environment is required.

WARNING ⚠

Under no circumstances is allowed handling of products or samples in corrosive, flammable, or explosive gas ambient.

Never handle products or samples biologically pathogenic or likely to harm in any way the user or the environment.

Never handle cytostatic (chemotherapy) or radiological samples.

Right working conditions are obtained only when there is no alarm message on the display.

2.3. STANDARD EQUIPMENT

- Two inlet sockets 2A - 460W each located on both corners of the work zone. Insulation level IP44.

2.4. CERTIFICATION

Quality control and tests results certification carried out in line with current standards. Certificate inside the cabinet.

2.5. TECHNICAL FEATURES

Features	Unit	Aeolus V 3	Aeolus V 4	Aeolus V 5	Aeolus V6
External dimensions (without support stand)	WxDxH (mm)	1048x798x1220	1353x798x1220	1658x798x1220	1963x798x1220
Internal dimensions	WxDxH (mm)	925x693x615	1230x693x615	1535x693x615	1840x693x615
Weight	Kg	140	160	185	210
Air velocity	m/s	0.36	0.36	0.36	0.36
Lighting	Lux	≥1200	≥1500	≥1500	≥1700
Noise level**	dB (A)	≤ 55	≤ 55	≤ 55	≤ 55
Vibration	Mm RMS	≤0.005	≤0.005	≤0.005	≤0.005
Voltage*	V	120 or 230	120 or 230	120 or 230	120 or 230
Frequency	Hz	50/60	50/60	50/60	50/60
HEPA/ULPA H14 Filter	Efficiency	Efficiency 99.995%, EN 1822 by MPPS method (99,999% D.O.P. test @ 0.3µm particles)			
Fans		EC High efficiency Centrifugal.			

* To ensure proper installation, check the attached electrical diagrams.
** Installation conditions may result in different noise values. Declared values measured in anechoic chamber at 1 m distance (K_{pA}= 2 dB).

3. INSTALLING INSTRUCTIONS

For a correct vertical LFB start-up procedure, observe these instructions step by step.


3.1. TRANSPORT

The cabinet is supplied in an individual package and in the suitable way to prevent it from being damaged.

This User Manual and a certificate of individual Quality controls of each cabinet are included in the scope of supply.

3.2. PLACING ON SITE

a) Use the same packaging base for the transport to its final location and it is highly recommended doing it with a pallet jack. The cabinet has two screwed holes intended to fasten two M10 eyebolts to help position it on a raised surface. If there are no mechanical means to position it, this can be done manually.

WARNING ⚠

We advise against the use of lift trucks for placing the cabinet directly by its base. In such case, it should be handled with care, and it is also a must that:

- the useful length of **the forks of the lift-truck are over 750 mm**.
- direct visual inspection during its handling is done.

Fig. 1 X **Fig. 2** ✓

Fig. 3 ✓

WARNING ⚠

For any handling during assembly or dismantling of the cabinet, appropriate safety measures should be taken, such as wearing safety gloves or shoes.

b) Cabinets should be placed away from entrances, corridors, ventilation supply points, etc., which might affect airflow in the work area.

c) Ensure the cabinet is set perfectly on the table or workbench to prevent vibrations.

d) We recommend that the room temperature is always within the comfort margins for working. That will prevent the work area from overheating and assure optimum cabinet performance.

The optimal comfort margins that match with speed probe optimal range are:

- Temperature: 20-30ºC.
- Relative humidity: less than 90%.

If any of these conditions is not fulfilled, optimum cabinet performance cannot be guaranteed.

e) When necessary, air vents in the room should be modified or speed of air entry reduced so that air is not aimed directly at the cabinet, as this may produce turbulence.
Side air vents are not appropriate for rooms where cabinets are installed.

f) Due to the ceiling temperature affects air distribution, sites should be selected with enough space. We recommend 250 mm free space between the top of the cabinet and the ceiling.

g) The location of the cabinet will directly influence the noise level. This level will rise considerably in small spaces with very reflective surfaces because of the proximity of the cabinets to the walls, and where fan speed has risen to compensate for air currents in the room or the resistance of clogged filters.

h) Correct cabinet performance should be tested before putting it into operation if any modification has been made to the room or to the air intake system which may considerably change the airflow characteristics in the room.


3.3. PRECAUTIONS FOR FIRST CONNECTION

Before connecting it to the mains supply, proceed to make a general cleaning to remove any dust particles that might be accumulated during transportation. Follow the instructions about cleaning and disinfection.

Before using it for the first time, run the cabinet for around 6 hours to purge the filters. New filters give off a characteristic smell which tends to disappear in a short period of time.

WARNING ⚠

The cabinet should work at the voltage shown on the name plate. It should be connected to a suitable socket base (minimum 16A) with an earthed plug. The cable section recommended for the base is 1.5 mm^2.

4. ELECTRICAL CONNECTION AND SERVICES


The cabinet is supplied with a plug cord according to European standard (other plug cords available as optional accessories). The general on/off switch is on the top right-hand of the cabinet.

The power socket in the wall should be in a perimeter around the cabinet of:

- 1.25 m from the middle of the cabinet to the sides
- 2.5 m from the base of the cabinet upwards

On the top side of the cabinet are:

- Voltage input (1)
- Switch (2)
- Main fuse (3)
- Checking valve for D.O.P. control (4)

WARNING ⚠

Unplug the cabinet from the electric current supply, in case of electrical maintenance. The cabinet can only be activated by means of the switch located on the top right-hand of the cabinet. The installation to which the cabinet is connected should meet the safety requirements in force in each country.

4.1. INSTALLATION TESTS

In order to ensure correct cabinet operation, after the installation and connection of the cabinet, TELSTAR recommends performing the following tests (before the equipment is used for the first time). See Chapter 11 for detailed instructions about.

- Filters integrity
- Impulsion air velocity
- Smoke pattern

Telstar, its representatives, and distributors put a qualified technical team at your disposition that uses control equipment with maximum health protection for carrying out this installation test under full safety guarantees.

4.2. STOPPING FOR LONG PERIODS

Turn off the general switch on the top of the cabinet.

In the event of a long stop, avoid as much as possible the accumulation of dust by turning the cabinet off with the glass in the closed position. When using the bench again, proceed as if it were started for the first time. See chapter 3.3.

5. INSTRUCTIONS FOR USE

5.1. GENERAL RECOMMENDATIONS

a) We recommend starting up the cabinet around 10 minutes before you begin to work. This will ensure particles are swept from the work zone and from material introduced.
b) Wash your arms, hands, and nails well, using an anti-bacterial soap, before and after work. Staff should avoid touching their mouths and eyes.
c) We recommend wearing lab coats with adjusted cuffs and, in special works, protection gloves (or sleeve cuffs). Both the coat and cuffs should be made from a fabric that minimises the emission of fibres and particles.
d) The work area inside the cabinet should not be used to store lab equipment. This can produce an accumulation of dust that is unnecessary and hazardous for working under sterile conditions. Only the fewest pieces of equipment should be placed to work correctly and with the minimum number of interruptions for the user (inserting and removing arms from the work zone).
e) Although the work zone is under a clean and sterile laminar down flow, the area around the cabinet might be contaminated. You should know the sources of contamination should be known to avoid them when inserting arms, material, etc.
f) All the material needed for the work should be free of particles and carefully cleaned before it is placed on the work area.
g) Materials such as paper, wood, cardboard, pencils, rubbers, etc. should not be placed in the work area as they shed a lot of particles.
h) In the event of using pipettes, only mechanical aspiration type should be used. Blowing the pipette with the mouth should never be done because it is easy to inhale aerosols caused by suction.
i) When platinum handles are used, it is advisable to use electric incinerators, preferably single-use ones.
j) If the work requires the use of a gas flame from a Bunsen burner or similar, we recommend using a type where the flame is maintained by holding down a button. Please note that the use of a constant flame burner creates significant turbulence. Bear in mind that a flame too large could burn the absolute filters and damage other parts.
k) Shadows and turbulence provoked by objects, equipment and materials placed in the work area should be studied before beginning to work to evaluate their possible effects. Bear in mind that the laminarity of the airflow is not recovered until a distance of 2.5 times the diameter of the object that causes the obstruction.
l) When using bottles and tubes, they should preferably have a screw top instead of cotton stoppers which easily shed a lot of particles.
m) When the work is finished, all disposable products (handles, Petri dishes, etc.) as well as means of cultivation, samples, tubes, bottles, etc. should be removed from the cabinet in impermeable bags suitable for sterilisation if necessary.
n) When working in the cabinet and during cleaning operations, it is important to protect the HEPA filters from any damage by avoiding knocking them, projecting liquids, or splashing, etc.
o) Any equipment connected to the internal sockets inside the cabinet should be fitted with an earth wire and meet all the electrical specifications.
p) People traffic generates air movements that may cause unexpected loss of the laminarity airflow and a reduced level of product protection. Warning signs and other means should be used to avoid traffic during critical operations.
 Cabinet locations can be divided to reduce the effects of air movement produced by other cabinets and people traffic.
q) Suddenly doors opening and closing may disturb airflows in the cabinet, and it should be avoided when cabinets are being used.

WARNING ⚠

The UV light system can only be used when the front glass is closed. Avoid working while the UV light is on. Never connect equipment whose consumption exceeds the maximum one indicated on the internal sockets of the cabinet.

5.2. ERGONOMICS

Aeolus Series and its accessories have been designed following recommendations and specifications of the standard EN ISO 14738:2010.

Users should consider their position when working with the cabinet, for their own safety.

Below are examples of four positions **NOT** recommended and one that **IS** recommended. (Fig. E is correct).


Fig. A – X **Fig. B – X** **Fig. C – X** **Fig. D – X** **Fig. E - √**

Adjusting the seated position to adopt a correct posture will allow the user to work within safety conditions and moreover to avoid possible physical injury (view, cervical, back, etc.)

As well as besides protecting the product as of the user, it is necessary to locate the safety cabinet in a place where no air currents can be produced that could cause turbulence in the laminar flow.



Fig. F – X

5.3. DEFINITION OF THE WORKING AREA

To prevent from creating any turbulence in the laminar down flow and to protect the product, it is necessary to establish a working area inside the cabinet.

- Do not work cover up the area with perforations.
- Always work 50 mm from the limits of the lateral glasses and the rear perforated area.

5.4. STARTING UP

After dealing with the above points, turn on the key switch (located at top-right of the bench) to start up the laminar flow bench. Ensure that the air intake at the top is clear.

The electronic system of flow monitoring by microprocessor via a thermoanemometer probe, needs a few seconds to stabilise during connection to begin to work normally.

5.5. SWITCHING OFF


To stop the bench, turn the key switch (located at top-right hand of the bench).

5.6. FRONT GLASS POSITIONS


By means of the handles, the front glass can be moved from working position to closed position or vice versa. To move the glass just slide down or up the glass.

The frontal glass is held in working position (250mm approx.) by the action of magnets.


FROM WORKING POSITION

FROM CLOSED POSITION


In case it is needed to clean the glass by the upper part or the inner part (Chapter 6.1), frontal can be open, and glass may be risen at the appropriate height, providing access to the upper part of both inner and outer glass areas. Lower part of the inner glass area should be cleaned on working position of glass.

WARNING ⚠

Always assure that open frontal is securely held by its lateral supports before accessing to the upper part of the glass for cleaning or for any other required task.

If it is necessary to place in the working area any big product or device, the glass may be risen till the top, allowing a maximum span of about 480 mm.

5.7. CONTROL PANEL

5.7.1. COMPONENTS DESCRIPTION

LCD display with 3x5.5 mm size characters that shows all the information required in cabinet-user communication.

5.7.2. WORKING MODE MENU FUNCTIONS

button (1): operates the fan (on/off). By default, the fan speed is set as standard.

- By pressing the  down button (5) fans will run in eco mode (slower).
- By pressing  up button (6) they will run at normal speed again.


button (2): operates the lighting (on/off), linked electrically with the UV germicidal lamp button (3).

button (3): operates the UV germicidal lamp (on/off) (optional), linked electrically with the lighting switch (2). In case there has not been yet programmed a specific time (time= 00:00, see Programming the U.V. timer value in chapter 5.7.4), it will work until switching it off.

button (4): operates the AUX output for possible optional control relay.

5.7.3. SIGNAL LIGHTS & ALARMS

- Signal light (7) for air-impulsion system functioning, green.
- Signal light (8) for lighting switched on, blue.
- Signal light (9) for UV germicidal lamp switched on, yellow.
- Signal light (10) for auxiliary output activated, green.
- UV germicidal lamp on (12), red.
- Alarm light (11):
 - Insufficient laminar down flow alarm → red - steady.
 - Eco mode → red – intermittent.
 - Fans stopped → red – intermittent.

5.7.4. USER CONFIGURATION MENU FUNCTIONS

To enter, it should be pressed down (5) and up (6) buttons at the same time for some seconds.

- Button (1): shows the working hours of the fan or the filter (by pressing the button once or twice respectively)


- Button (3): shows the working hours of the U.V. light

Programming the U.V. timer value:

- Button (5): decreases the U.V. timer value

- Button (6): increases the U.V. timer value


Note: Timer value = 00:00 means timer is off

Changing the language:

By pressing 

button (4) again, language can be changed.

To see the available languages, press the arrow buttons (5) or (6).

By pressing 

button (4) again, it confirms the language and U.V. timer values and it goes back to the working mode menu.

5.8. OPTIONAL ACCESSORIES

The vertical laminar flow benches can include several additional accessories according to the work requirements.

These optional accessories are the following:

- Germicidal lamp (UV) kit

Easy-to-work microprocessor-controlled UV germicidal lamp with a timer function.

The "UV" button on the control panel switches it on. For safety reasons, this button is electrically linked with the light button, and it switches off if the front window is opened.

Average disinfection intensity over the work zone over 40 μW/cm^2 with 20% depreciation after a working life of 8000 hours.

WARNING ⚠

When this accessory is used, you should consider the normal precautions in the use of these types of lamps. It should only be used for the previous external sterilization of the elements introduced into the bench. Therefore, this lamp is an accessory to be used previously to the operation with the bench, and it is essential to have it turned off when the user is in front of the table or working at it.
The glass provides protection from direct radiation, but UV rays are reflected on polished metal parts and can seriously damage an operator's eyes. That is why the glass front must be into the closed position to be able to switch on the U.V. light, after the fans are disconnected.
The use of this accessory does not replace the need of proper manual disinfection and cleaning, which is explained in chapter 6.

- Vacuum stop-cock/tap
- Water stop-cock/tap
- Gas stop-cock/tap

WARNING ⚠

Gas is a toxic and flammable substance which in the event of an accident can cause serious damage to the user and the cabinet.

Gas connection should be done in line with the safety standards of each country, region, and company.

- Other gasses stop-cock/tap (compress air, nitrogen, oxygen)

WARNING ⚠

When compressed air or nitrogen are used, use an absolute filter in the line to avoid the introduction of particles inside the working area.

- Support tables (two heights available)
- Extension legs for support stands.
- Height adjustable support stand (electrically operated)
- IT version: articulated VESA bracket
- Hook fluid therapy rod

6. CLEANING AND DISINFECTING

The normal procedure consists of a previous cleaning prior to disinfection.

Cleaning can be carried out with water and soap.
For Stainless Steel parts, use sterile or non-use disables cloths, which do not shed particles or fibers.

For disinfection, these cloths will be used, slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or affect the stainless steel or painted parts neither the glasses.

To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabinet; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

Also be sure of carrying out the cleaning/disinfection in the lower part of the working area (surface and corners).

WARNING 

It is recommended NOT to use disinfectant products containing chlorine or halogens because they can damage the stainless steel and painted surfaces.


WARNING 

It is required to dry completely all surfaces always after cleaning or disinfecting in order to preserve the original qualities of the surfaces.

WARNING 

It is advised not to spray the filters when performing the cleaning/disinfection tasks. Not following that could result in damaging them.

To accede to corners or other areas with difficult access areas, another type of tools will be used as shown in the figures.

Inside, bench cleaning is recommended on the following occasions:

a) Before starting any work in the bench.
b) After working in the bench.
c) Whenever there is a change of work program.
d) In event of liquid spilling on the worktable.
e) Before carrying out mechanical or biological control tests in the working area.

6.1. FRONT GLASS CLEANING AND DISINFECTION

For the cleaning of front glass, you must use sterile textures or disposable ones without particle or fibers loosing.

For disinfection, use the previous textiles slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or affect the glass. Cleaning can also be carried out with water and soap prior to disinfection.

To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabin; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

For the glass exterior part cleaning, it is recommended to keep the glass in closed position. In this position a major part of the glass is accessible. For the upper part of the glass and interior part of the glass cleaning, see chapter 5.6.

WARNING ⚠

It is required to dry completely all surfaces always after cleaning or disinfecting in order to preserve the original qualities of the surfaces.

7. MAINTENANCE INSTRUCTIONS

In general, a regular control of the proper working of the bench components should be done once a year.

WARNING 

Turn off the general switch on the top of the cabinet and disconnect the plug cord of the bench from the power socket of the wall before attempting any Maintenance or disassembly work. All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

7.1. REGULAR VALIDATION AND CONTROL

Syntegon Telstar recommends regular validation and control (at least once a year). The tests to be performed must be the same as those done at the last manufacturing stage, issuing a result of certification.

7.2. MAINTENANCE TABLE

The frequencies determined for each operation, and in particular with regard to changing filters, which appear on the Maintenance table below are indicative values recommended for normal use of the bench. Very severe work conditions may reduce the time intervals shown on the table.

Maintenance table

OPERATION	**PERIODICITY**
Impulsion filter leakage test*	*Minimum once a year or after filter changes*
Impulsion air velocity test*	
Smoke visualization of product protection*	

(*) See Chapter 11 for detailed instructions about.

Replacing the absolute filter	Every 4000/5000 working hours**
Cleaning the prefilters	Every 500 working hours**
Replacing the prefilters	Every 1000 working hours**

(**) Depending on quality of ambient air.

7.3. ABSOLUTE FILTERS

The working hours of the HEPA filter are logged and indicated by the machine's elapsed time meter.
After HEPA filter has reached its average life, the propelled airflow falls because of the increased pressure drop of the filter as it is a mechanical filtering, however this does not mean that filtering efficacy is reduced. When it reaches the minimum flow measured by the thermoanemometric probe, the microprocessor activates a visual alarm on the control panel.

WARNING 

To ensure correct installation and adjustment of the filter, as well as ensuring there are no leaks, filters should be replaced by specialized technical staff equipped with the required monitoring instruments: anemometer, leak detector, aerosol generator, particle counter, etc.

7.4. PREFILTERS

The vertical laminar flow sterile benches are provided with a G3 prefilter in the air intake.

WARNING 

Never start up the bench without having previously installed the prefilter correctly, since this would cause an early HEPA filter clogging.

8. LIST OF SPARE PARTS

Only the use of original parts and accessories will guarantee correct operation of the vertical laminar flow bench.

Please do not hesitate to contact the Customers Service pointing out the bench model, operating frequency (Hz), voltage (V) and serial number of the unit.

Aeolus V 3

No.	Description
1	ABS.FILTER 915X515X69 LPD H14 + FILM
2	PREFILTER G3 AEOLUS V3
3	LED BAR Y-3014B-216-900-T 24V
4	UV LAMP 15W T8

Aeolus V 4

No.	Description
5	ABS.FILTER 1220X515X69 LPD H14 + FILM
6	PREFILTER G3 AEOLUS V4
7	LED BAR Y-3014B-288-1200-T 24V
8	UV LAMP 30W T8

Aeolus V 5

No.	Description
9	ABS.FILTER 1525X515X69 LPD H14 + FILM
10	PREFILTER G3 AEOLUS V5
11	LED BAR Y-3014B-360-1500-T 24V
12	UV LAMP 30W T8

Aeolus V 6

No.	Description
13	ABS.FILTER 1830X515X69 LPD H14 + FILM
14	PREFILTER G3 AEOLUS V6
15	LED BAR Y-3014B-432-1800-T 24V
16	UV LAMP 30W T8

9. TROUBLESHOOTING

- Alarm light (11):
 - Insufficient laminar down flow alarm → red - steady.
 - Eco mode → red – intermittent.
 - Fans stopped → red – intermittent.

No.	Type of error	Causes of error	Correction of error
1	**Insufficient flow** Red-steady ALARM	Prefilter on the top, partially covered	Visually check the prefilter on the top of the bench and remove any type of obstruction.
		Distance between top of cabinet and lab ceiling too small	Increase space between top of cabinet and the ceiling. In some cases, it is possible to modify or eliminate the false ceiling. Install a lower support table. This means the working area will be also lower.
		Clogged impulsion filter	Check the working hours that the filter has been under operation. (**Chapter 5.7.4**) Notify the technical service to change the filter.
		Main control is damaged	Call a specialized technician to replace the control.
		Fans are not working properly	Call a specialized technician to replace the fans.
		Probe not working properly	Call a specialized technician to replace the probe.
2	**Fans stopped or at eco-mode** Red-intermittent ALARM	Fans are switched off	Check airflow. If no airflow, you can switch on the fans.
		Bench working at eco-mode	Select normal speed or keep at eco-mode.
		Main control is damaged	Call a specialized technician to replace the control.
3	**Light does not turn on**	Led bar light not working	Call a specialized technician to replace the led bar light.
		Some connection is wrong	Notify the technical service to revise connections.
		Main control is damaged	Call a specialized technician to replace the control.

No.	Type of error	Causes of error	Correction of error
4	**UV germicidal lamp does not turn on (optional accessory)**	UV lamp does not work	Check that the lamp & lamp holders are properly connected. Change the UV lamp. Notify the technical service to revise cable connections.
		Glass not in closed position	Assure that the glass is fully closed (lower position).
		Ballast is damaged	Notify the Technical Service to change the ballast.
5	**Display does not turn on**	Power cord cable disconnected	Check that power cord cable is properly connected in the socket on the wall.
		Cabinet not switched on yet	Turn on the key switch (located at top-right of the bench).
		No power	Check voltage on the socket of the wall.
		Display is damaged	Notify the Technical Service to replace the display.
		Main control is damaged	Notify the Technical Service to replace the control.

WARNING 

All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

10. MAINTENANCE CHART

MAINTENANCE CHART			
DATE	**WORK DESCRIPTION**	**NAME**	**SIGNATURE**

11. ANNEXES

11.1. ELECTRICAL DIAGRAM 120V / 230V 50/60 HZ

EQUIPMENT:
Aeolus V Cabinet

PROJECT REFERENCE:
77020

CUSTOMER:
TELSTAR TECHNOLOGIES, S.L.

DRAWING TYPE:
ELECTRICO

DRAWINGS:
06

MATERIAL O.F.:

PROJECT LEADER:
Victor Lazaro

AUTOMATION LEADER:
X. Cañellas

DRAWN BY:
X. Cañellas

DATE:
03/07/2013

APPROVED:
X. Cañellas

DATE:
03/07/2013

REVIEWED:
X. Cañellas

DATE:
03/07/2013

AS BUILT:

DATE:

REVISIONS:
ID: P DATE: 02/06/2023 DRAWN BY: X.Cañellas
ID: Q DATE: 14/12/2023 DRAWN BY: X.Cañellas
ID: R DATE: 12/08/2024 DRAWN BY: X.Cañellas
ID: N DATE: 05/12/2022 DRAWN BY: X.Cañellas
ID: O DATE: 16/05/2023 DRAWN BY: X.Cañellas

Telstar®

Date 03/07/2013 Report of Drawings 1014191 Página 1

AEOLUS V

TELSTAR TECHNOLOGIES, S.L.


Drawing	Description	Creation	Modification
01	SUMMARY	05/02/2013	04/09/2014
02	POWER	05/02/2013	14/12/2023
03	CONTROL BOARD AND LIGHTING	05/02/2013	19/07/2021
04	MOTORS AND SOCKETS	05/02/2013	12/08/2024
05	OPTIONAL PARTS	05/02/2013	14/06/2017
06	BILL OF MATERIALS	05/02/2013	14/12/2023

PROPERTY OF TELSTAR
REPRODUCTION FORBIDDEN

REPLACES:

REPLACED BY:

PROJECT NAME:
1014191

ELECTRICAL CARACTERISTICS

MODEL	AEOLUS V 3	AEOLUS V 4	AEOLUS V 5	AEOLUS V 6
VOLTAGE	120 / 230 VAC	120 / 230 VAC	120 / 230 VAC	120 / 230 VAC
CURRENT	9.7 / 7.3 A	9.7 / 7.3 A	9.7 / 7.3 A	11.9 / 8.4 A
P.F.	0.6	0.6	0.6	0.6
*POWER	1164VA / 1680 VA	1164VA / 1680 VA	1164VA / 1680 VA	1430 / 1932 VA
FREQUENCY	50..60 HZ	50..60 HZ	50..60 HZ	50..60 HZ

* Electrical socket's maximum power included

Revision	Date	Draw by	Date:
O	16/05/2023	X.Cañellas	Draw by:
P	02/06/2023	X.Cañellas	Aproved:
Q	14/12/2023	X.Cañellas	

	Date:	03/07/2013
	Draw by:	X. Cañellas
	Aproved:	X. Cañellas

Project description: AEOLUS V CABINET

Situation:

Customer: TELSTAR TECHNOLOGIES, S.L.

Zone:

Drawing description: POWER

Drawings: 02 / 06

Previous: 01


Next: 03


Project reference: 77020

Drawing number: 1014191

Drawing number: 02


Telstar®

02/H1
1L
1N
PE
02/H1
02/H1
04/A1
04/A1
04/A1
02/G4
02/G4
23
24
X3
04/A4
X2
05/A6
CABLES: 1014206
CABLES: 1014208
AEOLUS 3,4,5: 1014208
AEOLUS 6: 1014208
Black
White
Red
White
Yellow
Blue
SENSOR
FAN
X1
X2
X3
X4
24VDC
DISPLAY
KEYPAD
CABLES: 1014232
-CTL1
1014196
AEOLUS PCBA
-DPL1
128B151
5A
5A
5A
K1
K2
K3
9X
1X6
2X6
3X6
4X6
5X6
6X6
24
25
1L
4L
05/A1
CABLES: 1014210
Red
Black
-LD1
**
24VDC
**
CABINET
CODE LD1
WATTS
AEOLUS V3
621658
23 W
AEOLUS V4
621427
30 W
AEOLUS V5
621659
38 W
AEOLUS V6
621660
45 W
Revision
Date
Draw by:
Date:
K
19/07/2021
X.Cañellas
B
04/09/2014
X.Cañellas
Draw by:
I
22/03/2018
X. Cañellas
Aproved:
03/07/2013
X. Cañellas
X. Cañellas
Project description:
AEOLUS V CABINET
Situation:
Customer:
TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
CONTROL BOARD AND LIGHTING
Drawings: 03 / 06
Previous: 02
Next: 04
Project reference:
77020
Drawing number:
1014191
Drawing number
03
Telstar

Revision	Date	Draw by:	Date:
N	05/12/2022	X.Cañellas	Draw by:
R	12/08/2024	X.Cañellas	Aproved:
L	20/10/2021	X.Cañellas	

Project description:	AEOLUS V CABINET
Situation:	
Customer:	TELSTAR TECHNOLOGIES, S.L.
Zone:	
Drawing description:	MOTORS AND SOCKETS
Drawings:	04 / 06
Previous:	03
Next:	05
Project reference:	77020
Drawing number:	1014191
Drawing number	04

PE
1N
4L
04/H1
04/H1
03/F4
03/B5
X2
OPTIONAL
-SP1
621613
AIRFLOW PROBE
OPTIONAL
1N
4L
637554
-BL1
120/230VAC
0.3/0.13A
White
Black
Red
Blue
Blue
N.C.
1341/1533
-L1
TB
30/15W
NO
COM
NC
-S1
500632
HF-D02
1A-3
NO
COM
NC
-S2
500632
HF-D02
1A-3
CABINET SIZE
UV KIT CODE
UV LAMP
4, 5, 6
1695432
1341
3
1679973
1533
Revision
Date
Draw by
Date:
B
04/09/2014
X.Cañellas
03/07/2013
C
10/02/2016
X.Cañellas
Draw by:
X. Cañellas
I
22/03/2018
X.Cañellas
Aproved:
X. Cañellas
Project description:
AEOLUS V CABINET
Situation:
Customer:
TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
OPTIONAL PARTS
Drawings: 05 / 06
Previous: 04
Next: 06
Project reference:
77020
Drawing number:
10114191
Drawing number
05
Telstar

Date 03/07/2013 Bill of Materials Page 1

AEOLUS V CABINET

TELSTAR
TECHNOLOGIES, S.L.

Symbol	Drawing	Column	Description	Refer. Numb.	Catalogue	Reference
-F1	02	B	FUSE 10A 6.3x32 LANBAO ELECT (230V Version)	RTI-30T10A250V	-	655980
			FUSE 12.5A 6.3x32 LANBAO ELECT (120V Version)	RTI-30T12.5A250V	-	655979
-FL1	02	B	EMC FILTER 20A YANBIXINKEJI	YB12E4-20A-Q	-	634174
-SW1	02	B	ROCKER SWITCH (250V 16A, 125V 16A)	-	-	657415
-PS1	02	E	Power Supply MEANWELL	RS-75-24	-	621116
-CTL1	03	C	Control Board AEOLUS	-	-	1014196
-DPL1	03	C	DiSPLAY AEOLUS	-	-	1288161
-LD1	03	G	LED BAR 900MM 800LUX 24V (Aeolus V 3)	Y-3014B-216-900-T	-	621658
			LED BAR 1200MM 800LUX 24V (Aeolus V 4)	Y-3014B-432-1800-T	-	621427
			LED BAR 1500MM 800LUX 24V (Aeolus V 5)	Y-3014B-288-1200-T	-	621659
			LED BAR 1800MM 800LUX 24V (Aeolus V 6)	Y-3014B-432-1800-T	-	621660
-F2	04	C	Fuse 5A 5x20mm	-	-	606806
-F3	04	G	Fuse 5A 5x20mm	-	-	606806
-FAN1	04	B	FAN BELONG 230VAC Version	B2ED146LB2M-1500RPM	-	634135
			FAN BELONG 115VAC Version	PB1ED146LB2M-1500RPM	-	637650
-FAN2	04	C	FAN BELONG 230VAC Version	B2ED146LB2M-1500RPM		634135
			FAN BELONG 115VAC Version	PB1ED146LB2M-1500RPM	-	637650
-FAN3	04	D	FAN BELONG 230VAC Version	B2ED146LB2M-1500RPM		634135
			FAN BELONG 115VAC Version	PB1ED146LB2M-1500RPM	-	637650
-BE2	04	F	SOCKET EUROPE YUADON	YD-001S	-	637516
-BE3	04	G	SOCKET EUROPE YUADON	YD-001S	-	637516
-SP1	05	C	AirFlow Probe Degree Controls 0-10M/S	FS22503	-	621613
-	05	F	KIT UV AEOLUS V 4, 5, 6	-	-	1695432
			KIT UV AEOLUS V 3	-		1679973

Revision	Date	Draw by	Date:	03/07/2013
Q	14/12/2023	X.Cañellas	Draw by:	X. Cañellas
R	12/08/2024	X.Cañellas	Aproved:	X. Cañellas
P	02/06/2023	X.Cañellas		

Project description: AEOLUS V CABINET

Situation:

Zone:

Customer: TELSTAR TECHNOLOGIES, S.L.

Drawing description: BILL OF MATERIALS

Drawings: 06 / 06

Previous: 05

Next:

Drawing number: 1014191

Project reference: 77020

Drawing number: 06

Telstar®

11.2. IMPULSION FILTER INTEGRITY TEST

PURPOSE

Check the integrity of the filters.

SCOPE

Aeolus V.

TOOLS

- Aerosol generator (cool or hot).
- ATI PAO-4, *Shell Ondina fluid or Emery 3004 oils*.
- *Photometer (ATI 2H).*

METHODOLOGY

Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Fig. 1).

Fig. 1

Inject aerosol into the flow before the filter, through the upper prefilter (Fig. 2).


Fig. 2

With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters (including the perimeter of each filter pack) by passing the photometer probe in slightly overlapping strokes over the entire surface.

Also scan the entire periphery of the filter and the junction between filter and filter mounting frame at a scanning rate that does not exceed 5 cm/s (Fig. 3).

Fig. 3

ACCEPTANCE CRITERIA

The photometer signal after the filter should be smaller or equal to a 0,01% of the specific concentration before the filter.

Figure 4. Correct result: filter not leaking

Figure 5. Wrong result: filter leaking

11.3. AIR VELOCITY TEST

PURPOSE

To check that the speed control in different points of the working area is suitable to achieve product protection. To determine the speed average as well as the uniformity range in the laminar flow zone.

SCOPE

Aeolus V.

TOOLS


The tool to use is an anemometer.


METHODOLOGY

Switch on the cabinet and wait until the velocity and the temperature becomes steady.


Depending on the model of the cabinet, take as many readings as indicated in the following images.

Place the omnidirectional anemometer at 450 mm from the working surface over each point according to images below.



Aeolus V 3

Aeolus V 4

Aeolus V 5

Aeolus V 6

For this calculation, we proceed to the following formula:

$$V_{average}\ (m/s) = 1/n \times (V_1 + V_2 + V_3 + V_n)$$

ACCEPTANCE CRITERIA

a) The average of the airflow velocity should be between the values shown below.
0.35m/s ≤ $V_{average}$ ≤ 0.37 m/s

Note: Measurement uncertainty of the anemometer used should be considered and could result in wrong measurement or comparison between measurements.

b) No individual measurement should differ from the average more than ± 20 %.

11.4. SMOKE PATTERN TEST

OBJECT

To verify the cabinet for product protection.

SCOPE

Aeolus V.


TOOLS

Smoke generator.

METHODOLOGY

Proceed as follows:

a) Generate smoke (uniformly) throughout the cabinet's working area.
b) Verify the behavior of the smoke visually.

ACCEPTANCE CRITERIA

Turbulences in the working area shall not be created. The smoke generated should be dragged by the laminar flow to the surface of the working area.

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UK DECLARATION OF CONFORMITY	
Manufacturer:	SYNTEGON TELSTAR TECHNOLOGIES, S.L.U.
Manufacturer Address:	Av. Font i Sagué, 55 Parc Científic i Tecnològic Orbital 40 08227 Terrassa (Barcelona) Spain
N.I.F.:	B-63797559
Equipment Description:	Laminar Flow bench **AEOLUS H, AEOLUS V**
UK DIRECTIVES	**DECLARATION**
Supply of Machinery (Safety) Regulations 2008 ***(S.I. 2008/1597)*** *Electromagnetic Compatibility Regulations 2016* ***(S.I. 2016/1091)*** *Electrical Equipment (Safety) Regulations 2016* ***(S.I. 2016/1101)*** *The Restriction of the use of Certain Hazardous substances in Electrical and Electronic Equipment Regulations 2012* ***(S.I. 2012/3032)*** Specific standards: *Safety requirements for electrical equipment for measurement, control, and laboratory use* ***BS EN 61010:2010+A1***	The present declaration of conformity is issued under the sole responsibility of Manufacturer. This document hereby certifies that the above equipment complies with the UK Directives listed, and their amendments. The equipment must be installed, tested, and operated in strict accordance with Manufacturer instructions and used by fully trained operators. Only approved Syntegon Telstar Technologies S.L.U. spares must be used for maintenance and repair and fitted in accordance with Manufacturer instructions. UK CA
Body Retaining Technical File:	Retained by Syntegon Telstar Technologies, S.L.U.
Signed by:	
Name:	Jordi Valls
Title:	Global Quality Manager
On Behalf of:	SYNTEGON TELSTAR TECHNOLOGIES, S.L.U.
Date:	January 2025

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

REPLACED BY:

PROJECT NAME:
1014191


CABLES: 1014501
CABLES: 1014501
-F1
230V 10A
120V 12.5A
L
N
PE
-SW1
657415
250VAC 16A
125VAC 16A
634174
-FL1
YB12E4-20A-Q
115/250V 20A
1L
1N
CABLES: 1014198
CABLES: 1014199
CABLES: 1014206
621116
-PS1
230 Vac
75 W
24 Vdc
23
24
03/A1
PE
ELECTRICAL CARACTERISTICS
MODEL
VOLTAGE
CURRENT
P.F.
*POWER
FREQUENCY
AEOLUS V 3
AEOLUS V 4
AEOLUS V 5
AEOLUS V 6
120 / 230 VAC
9.7 / 7.3 A
0.6
1164VA / 1680 VA
50..60 HZ
120 / 230 VAC
9.7 / 7.3 A
0.6
1164VA / 1680 VA
50..60 HZ
120 / 230 VAC
9.7 / 7.3 A
0.6
1164VA / 1680 VA
50..60 HZ
120 / 230 VAC
11.9 / 8.4 A
0.6
1430 / 1932 VA
50..60 HZ
* Electrical sockets maximum power included
Revision
O
P
Q
Date
16/05/2023
02/06/2023
14/12/2023
Draw by
X.Cañellas
X.Cañellas
X.Cañellas
Date:
Draw by:
Aproved:
03/07/2013
X. Cañellas
X. Cañellas
Project description:
AEOLUS V CABINET
Situation:
Customer:
TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description:
POWER
Drawings: 02 / 06
Previous: 01
Next: 03
Project reference:
77020
Drawing number:
1014191
Drawing number
02
Telstar


**			
CABINET	CODE LD1	WATTS	
AEOLUS V 3	621658	23 W	
AEOLUS V 4	621427	30 W	
AEOLUS V 5	621659	38 W	
AEOLUS V 6	621660	45 W	

Revision	Date	Draw by	Date:	
K	19/07/2021	X.Cañellas	03/07/2013	Project description: AEOLUS V CABINET
B	04/09/2014	X.Cañellas	Draw by: X. Cañellas	
I	22/03/2018	X. Cañellas	Aproved: X. Cañellas	

Situation:
Customer: TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description: CONTROL BOARD AND LIGHTING
Drawings: 03 / 06
Previous: 02
Next: 04
Project reference: 77020
Drawing number: 1014191
Drawing number: 03

Revision	Date	Draw by	Date:	Project description:	Situation:	Zone:	Drawings: 04 / 06	Previous: 03	Next: 05
N	05/12/2022	X.Cañellas	03/07/2013	AEOLUS V CABINET	Customer: TELSTAR TECHNOLOGIES, S.L.	Drawing description: MOTORS AND SOCKETS		Drawing number: 1014191	
R	12/08/2024	X.Cañellas	Draw by: X. Cañellas				Project reference: 77020		Drawing number: 04
L	20/10/2021	X.Cañellas	Aproved: X. Cañellas						

04/H1
04/H1
03/F4
PE
1N
4L
03/B5
X2
OPTIONAL
-SP1
621613
AIRFLOW PROBE
OPTIONAL
1341/1533
-L1
-T8
30/15W
N.C.
Blue
Blue
Red
637554
-BL1
120/230VAC
0.3/0.13A
Black
White
1N
4L
NO
COM
NC
-S2
600632
HF-D02
1A-3
-S1
600632
HF-D02
1A-3
CABINET SIZE
UV KIT CODE
UV LAMP
4, 5, 6
1695432
1341
3
1679973
1533
Revision
B
C
I
Date
04/09/2014
10/02/2016
22/03/2018
Draw by
X.Cañellas
X.Cañellas
X.Cañellas
Date:
Draw by:
Aproved:
03/07/2013
X. Cañellas
X. Cañellas
Project description:
AEOLUS V CABINET
Situation:
Customer:
TELSTAR
TECHNOLOGIES, S.L.
Zone:
Drawing description:
OPTIONAL PARTS
Drawings: 05/06
Previous: 04
Next: 06
Project reference:
77020
Drawing number:
1014191
Drawing number:
05
Telstar

Date 03/07/2013 Bill of Materials Page 1
AEOLUS V CABINET

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Project description: AEOLUS V CABINET
Situation:
Zone:
Customer: TELSTAR TECHNOLOGIES, S.L.
Drawing description: BILL OF MATERIALS
Drawings: 06 / 06
Previous: 05
Next:
Drawing number: 1014191
Project reference: 77020
Drawing number: 06

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DATE:

REVISIONS:

ID: P	DATE: 02/06/2023	DRAWN BY: X.Cañellas
ID: Q	DATE: 14/12/2023	DRAWN BY: X.Cañellas
ID: R	DATE: 12/08/2024	DRAWN BY: X.Cañellas
ID: N	DATE: 05/12/2022	DRAWN BY: X.Cañellas
ID: O	DATE: 16/05/2023	DRAWN BY: X.Cañellas

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



Revision	Date	Draw by	Date:	Project description: AEOLUS V CABINET	Situation:	Zone:	Drawings: 04 / 06	Previous: 03	Next: 05
N	05/12/2022	X.Cañellas	03/07/2013		Customer: TELSTAR TECHNOLOGIES, S.L.	Drawing description: MOTORS AND SOCKETS		Drawing number: 1014191	
R	12/08/2024	X.Cañellas	Draw by: X. Cañellas				Project reference: 77020		Drawing number: 04
L	20/10/2021	X.Cañellas	Aproved: X. Cañellas						

04/H1
04/H1
03/F4
PE
1N
4L
03/B5
X2
OPTIONAL
-SP1
621613
AIRFLOW PROBE
OPTIONAL
1341/1533
-L1
-T8
30/15W
N.C.
Blue
Blue
Red
637554
-BL1
120/230VAC
0.3/0.13A
Black
White
1N
4L
NO
COM
NC
-S2
600632
HF-D02
1A-3
-S1
600632
HF-D02
1A-3
CABINET SIZE
UV KIT CODE
UV LAMP
4, 5, 6
1695432
1341
3
1679973
1533
Revision
B
C
I
Date
04/09/2014
10/02/2016
22/03/2018
Draw by
X.Cañellas
X.Cañellas
X.Cañellas
Date:
Draw by:
Aproved:
03/07/2013
X. Cañellas
X. Cañellas
Project description:
AEOLUS V CABINET
Situation:
Customer:
TELSTAR
TECHNOLOGIES, S.L.
Zone:
Drawing description:
OPTIONAL PARTS
Drawings: 05/06
Previous: 04
Next: 06
Project reference:
77020
Drawing number:
1014191
Drawing number:
05
Telstar

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Project description: AEOLUS V CABINET
Situation:
Customer: TELSTAR TECHNOLOGIES, S.L.
Zone:
Drawing description: BILL OF MATERIALS
Drawings: 06 / 06
Previous: 05
Next:
Drawing number: 1014191
Project reference: 77020
Drawing number: 06

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