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Printed document is not subject to the change service

ENGLISH VERSION

1. Generally Admitted

This product specification is the benchmark for quality assurance in the DWK Life Sciences and for the customer's acceptance inspection. The product properties and product characteristics are defined in the specification and the article-related drawing.

The articles described in this specification are manufactured using a quality system in accordance with DIN ISO 9001. They fulfil the quality standard stipulated in the specification and drawing.

In addition, DURAN® fulfils the requirements of DIN ISO 3585 (borosilicate glass 3.3).

DURAN® Laboratory Glass articles are manufactured by machine production, manufacturing methods as follows:

- press-blow process
- press process
- blow process

2. Scope

This product specification is the basis for quality assurance in the DWK Life Sciences and for customers. It will apply, from the validity date of this specification, to all manufactured and supplied Laboratory Glass articles made from DURAN® borosilicate glass 3.3.

3. Validity

This product specification is valid from the release date and replaces all previous specifications

Printed document is not subject to the change service

4. Product description

DURAN® borosilicate glass 3.3 is highly resistant to water, neutral and acid solutions, strong alkaline solutions and their compositions, chlorine, bromine, iodine and organic substances.

Only hydrofluoric acid, boiling phosphoric acid and strong alkaline solutions strip the glass surface (glass corrosion) when high temperatures (>100°C) also occur.

General technical values	
Thermal conductivity λw at 100°C	1.2 W • m^{-1} • K^{-1}
Mean specific heat capacity Cp (20°C; 100°C)	0.8 J • g^{-1} • K^{-1}
Modulus of elasticity E (Young's modulus)	64 • 10^3 N • mm^{-2}
Poisson's ratio μ	0.20
Stress-optical coefficient K	4.0 • 10^{-6} mm^2• N^{-1}
Refractive index (λ = 587.6 nm) n_d	1.473
Dielectric properties (1 MHz 25°C) dielectric constant ε	4.6
Dielectric loss factor tan δ	37 • 10^{-4}
Temperature for the specific electrical resistance of 10^{-6} Ω • cm (DIN 52 326) t_k 100	250°
Coefficient of mean linear thermal expansion α (20°C; 300°C) according to DIN ISO 7991	3.3 • 10^{-6} • K^{-1}
Transformation temperature T_g	525°C
Density ρ at 25°C	2.23 g • cm^{-3}
Temperature of the glass at viscosity η in dPa • s: ➢ $1,0^{13}$ (annealing point) ➢ $10^{7.6}$ (softening point) ➢ 10^4 (working point)	 560°C 825°C 1260°C

Printed document is not subject to the change service

4.1 Chemical composition

DURAN® has the following approximate composition (in % by weight):

SiO_2	B_2O_3	Na_2O & K_2O	Al_2O_3
81	13	4	2

4.2 Chemical resistance

4.2.1 General

Thanks to the almost inert behaviour of DURAN® borosilicate glass, interaction, for example ion exchange, between the medium and the glass can be effectively excluded.

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Printed document is not subject to the change service

4.2.2 Hydrolytic Resistance


DURAN® corresponds to “Class 1” of 5 different hydrolytic resistance classes in accordance with ISO 719 (98 °C). The amount of Na_2O/g glass grain leached out after 1 hour in water at 98 °C is measured.

For DURAN® the quantity of Na_2O leached out is less than 31 µg/g of glass grain.

DURAN® also corresponds to “Class 1” of 3 different hydrolytic resistance classes in accordance with ISO 720 (121 °C). The quantity of Na_2O leached out after 1 hour in water at 121 °C is less than 62 µg/g of glass grain.

Due to its good hydrolytic resistance, DURAN® meets the requirements of the USP, EP and JP for neutral glass that corresponds to glass type 1.

It can therefore be used in an almost unrestricted way in pharmaceutical applications and in contact with foodstuff.


Example: Hydrolytic attack on DURAN® as a function of time (100°C).

Printed document is not subject to the change service

4.2.3 Acid Resistance

DURAN® corresponds to “Class 1” of 4 different acid classes in accordance with DIN 12116. As the surface removal after boiling for 6 hours in normal HCl is less than 0.7 mg/100 cm^2, DURAN® is classed as acid-resistant borosilicate glass. The quantity of alkaline metal oxides leached out in accordance with ISO 1776 is less than 100 µg Na_2O/100 cm^2.

Example


Acid attack on DURAN® as a function of acid concentration.

Printed document is not subject to the change service

4.2.4 Alkali Resistance

DURAN® corresponds to “Class 2” of 3 different alkali classes in accordance with ISO 695.

The surface erosion after 3 hours boiling in a mixture of equal volume fractions of sodium hydroxide solution (concentration 1 mol/l) and sodium carbonate solution (concentration 0.5 mol/l) is only 134 mg/100 cm^2.

Example

Alkali attack on DURAN® as a function of temperature

Printed document is not subject to the change service

4.3 Physical properties

4.3.1 Thermal properties

Temperature resistance when heated

The highest permissible continuous usage temperature of DURAN® glassware is 500°C. At temperatures of 525°C and above, the glass will begin to melt and, from temperatures of 860°C, it will change into a fluid state.

Temperature resistance when frozen

DURAN® can be cooled down to the maximum possible negative temperature and is also suitable for use in liquid nitrogen (approx. -196°C). During use/freezing special attention should be given to the expansion of the content. In general, use down to -70°C is recommended for DURAN® products. In addition to the geometry of the products the properties of the additional components used (e.g. screw caps) also need to be considered. When cooling and thawing, it is important to make sure that the difference in temperature does not exceed 100 K.

In practice, therefore, cooling and heating in stages is recommended. When freezing substances in DURAN® bottles, the receptacle must only be filled to a maximum level of 3/4. Furthermore, it should be frozen in a tilted position of 45° (surface enlargement) and subject to the properties of the screw caps or other components.

Thermal shock resistance

As it has a very low coefficient of linear expansion ($\alpha=3.3 \times 10^{-6}\ K^{-1}$), a feature of DURAN® is its high thermal shock resistance up to $\Delta T = 100$ K. This means that for a temperature change of 1 K, the glass changes by only 3.3×10^{-6} relative length units. The thermal shock resistance is influenced by wall thickness and product geometry.

4.3.2 Optical properties

In the spectral range of approximately 310 to 2,200 nm, the absorption of DURAN® is negligible. It is transparent and colourless. Larger coating thicknesses (axial inspection for tubes) appear slightly yellowish / greenish.

4.3.3 Pressure resistance

Laboratory Glass is not suitable for application with pressure or under vacuum, except of the DURAN® pressure plus+ bottles.

4.3.4 Cleaning

Automatic cleaning in a cleaning machine is generally considered more gentle than manual cleaning in a soaking bath. During machine cleaning, the glass surface only comes into contact with the alkalinity of the cleaning agent for a short, defined time interval, therefore minimising glass corrosion. The machine cleaning system accessories provide special attachments and locks so that glassware is securely fastened and protected against breakage.

Laboratory glass articles are suitable for manual cleaning as well as automated cleaning (minimum resistance: 250 cycles in a household dishwasher – at a maximum of 65°C). The ability to easily see and read the printed scale is the assessment criterion.

Printed document is not subject to the change service

4.3.5 Sterilisation

Laboratory glass articles are suitable for autoclaving as well as dry heat and plasma sterilization (H_2O_2).
In respect of sterilisation the following information should be observed:

- To avoid overpressure all containers should kept open during autoclaving.
- When sterilising media it is recommended to use the membrane cap that allows pressure equalisation through a PTFE membrane and hence the cap can be tightly closed. Therefore the risk of glass breakage is nearly excluded.
- DURAN® is no cer-stabilised glass therefore the gamma-ray sterilisation is only partially suitable. DURAN® may appear brown coloured after gamma-ray sterilisation.

5. Quality

Laboratory glass articles are manufactured in accordance with the specification documents. Specification documents are the specification and the drawing. The geometrical design is in accordance with the valid drawing. All tests are carried out according to the following conditions:

- Without magnifying glass
- Colour temperature > 3600 K (neutral white/daylight-white)
- 800 lux lighting
- Background black/white matt
- The viewing distance is one arm length (approx. 40-60 cm)
- The frequencies always refer to an article.
- All tests are carried out in random sampling procedures according to AQL DIN/ISO 2859-1.
- The basis for this is the simple sample of the general inspection level of stage II and the valid sampling plan.

Qualitative characteristics are divided into measurable and non-quantifiable characteristics. For all measurable characteristics, corresponding limits are defined in the form of a maximum or minimum permissible quality tolerance. The characteristics that cannot be clearly quantified are defined by means of good/bad limit samples.

Characteristics which do not impair the usability, but are visually recognizable and which are not contained in the specification, are not recognized as defects.

The basis for assessing the size of a characteristic is the equivalent diameter. Longish (oval, elliptical) features deviating from the circular shape are determined by the equivalent diameter:


Definition of equivalent diameter:

$$L = \frac{a+b}{2}$$

L = equivalent diameter of the characteristic

a = length of the glass characteristic


b = width of the glass characteristic

Printed document is not subject to the change service

Equivalent diameter =

All other characteristics are determined according to their actual size or according to the enclosing circle (illustrative figure):

5.1 Quantitative characteristics

The evaluation of quantitative characteristics is carried out according to the currently valid palletizing standard.

5.2 Quantifiable quality characteristics

5.2.1 Residual stress

After processing, all DURAN products are stress-relieved by means of a cooling furnace in order to avoid residual thermal stresses. The residual stress is therefore always measured after the cooling process. The limit value for the residual stress is: $\sigma \leq 150$ nm/cm optical difference (technically relaxed)
(residual stress constant K for DURAN = 4.0 x 10-6 mm^2/N)

5.2.2 Dimensions

All dimensions tolerated in the drawing or dimensions marked with min./max. are relevant inspection dimensions .

Printed document is not subject to the change service

5.2.3 Glass Characteristics: Pressing/Blowing Process and Blowing Process

Acceptable characteristics / characteristic values related all to one item.

Characteristic	Extent	Max. acceptable per article (ml)				
		≤ 100	≤ 500	≤ 1000	≤ 2000	≤ 5000
Bubbles, open, inside	All sizes	Not admissible				
Bubbles, open, outside	≤ 0,5 mm ∅ > 0,5 mm ∅	admissible not admissible				
Bubbles, closed	≤ 0,5 mm ∅	admissible				
	> 0,5 - 2,0 mm ∅	2	3	3	4	6
	> 2,0 – 4,0 mm ∅	0	0	1	2	2
	> 4,0 mm ∅	0	0	0	0	0
Knots	≤ 0,5 mm ∅	admissible				
	> 0,5 - 2,0 mm ∅	3	3	4	5	7
	> 2,0 ∅	0	0	0	0	0
Stone, not surrounded by glass	All sizes	Not admissible				
Stone, surrounded by glass	≤ 1,0 mm ∅	2	3	3	4	4
	>1,0 mm ∅	0	0	0	0	0
Dirt removable	All sizes	admissible				
Dirt not removable	>0,3 inside	Not admissible				
Dirt not removable, outside	>0,5 - ≤4mm ∅	2	3	3	4	4
	>4mm ∅ outside	Not admissible				
Metal / cinder	≤ 0,5 mm ∅	admissible				
	> 0,5 - <1.0 mm ∅	2	3	3	4	4
	≤1,0 mm ∅	0	0	0	0	0
	All sizes	Not admissible				
Mouth and thread	Axial and concentric runout; slanted, offset	Not admissible				
Zirconia marks	All sizes	Not admissible				
Crack	All sizes	Not admissible				
Breakage	All sizes	Not admissible				
Maximum number of characteristics per article	Bubbles (closed) and Knots up to a max. 0,5 mm ∅ are not regarded as a characteristic	3	3	5	5	7
Minimum distance between admissible characteristics (mm)		5	10	10	10	15

If there are several defects per article, the characteristic of the sharpest AQL group is evaluated.

Printed document is not subject to the change service

5.2.4 Glass Characteristics: Pressing Process

Acceptable characteristics / characteristic values related all to one item.

Charcteristic	Extent	Max. acceptable per article ∅ (mm)						
		≤ 50	≤ 75	≤ 100	≤ 150	≤ 200	≤ 250	≤ 300
Bubble, open	All sizes	Not admissible						
Bubble, closed	≤ 0,5 mm ∅	admissible						
	> 0,5 - 2,0 mm ∅	2	2	3	4	6	8	12
	> 2,0 - 4,0mm ∅	0	0	1	1	1	2	3
	> 4,0 mm ∅	0	0	0	0	0	0	0
Knots	≤ 0,5 mm ∅	admissible						
	> 0,5 - 2,0 mm ∅	1	2	3	4	5	6	8
	> 2,0 - 4,0 mm ∅	0	1	2	3	4	5	7
	> 4,0 mm ∅	0	0	0	0	0	0	0
Stone, not surrounded by glass	All sizes	Not admissible						
Stone, surrounded by glass	≤ 1,0 mm ∅	1	2	3	4	5	6	10
	≤ 2,0 mm ∅	0	0	0	0	0	1	1
	>2,0 mm ∅	0	0	0	0	0	0	0
Ziconium streak	All sizes	Not admissible						
Crack	All sizes	Not admissible						
Breakage	All sizes	Not admissible						
Dirt removable	All sizes	admissible						
Dirt not removable	≤ 0.5 mm ∅	admissible						
Dirt not removable, outside	>0,5 - < 4 mm ∅	1	2	3	4	5	6	10
	≤4 mm ∅	Not admissible						
Metal / cinder	≤ 0,5 mm ∅	admissible						
	> 0,5 - <1.0 mm ∅	1	2	3	4	5	6	10
	≤1,0 mm ∅	Not admissible						
Maximum number of characteristics per article	Bubbles (closed) and Knots up to a max. 0,5 mm ∅ are not regarded as a characteristic	2	3	4	5	6	8	10
Minimum distance between admissible characteristics (mm)		5	5	5	5	10	10	20

If there are several defects per article, the characteristic of the sharpest AQL group is evaluated.

Printed document is not subject to the change service

5.2.5 Printing features

The printing is done with the screen printing process. The printing ink is pressed with a squeegee through a fine-meshed fabric onto the glass to be printed. In places where no ink is to be printed, the mesh openings of the fabric are made colour-impermeable with the aid of a stencil. If the printing stencil is destroyed by small particles, additional, unintentional printing may occur. The printed image shows additional colour stains that are not part of the production marking. Depending on the size and number of colour spots, the colour spots are permissible or not. To evaluate the size of the colour stains, the stencil "colour stains" is used. The printed image should be clearly visible despite the colour stains.

The basis for assessing the size of the characteristic is the equivalent diameter. Longish (oval, elliptical) features deviating from the circular shape are determined by the equivalent diameter. To assess the size of the colour stain/fault, use the template "Colour stains/faults". The printed image should be clearly visible despite colour stains or imperfections.

Characteristic	**Extent**	**Max. acceptable per article (ml)**				
		≤ 100	**≤ 500**	**≤ 1000**	**≤ 2000**	**≤ 5000**
Missing print	≤ 0,25mm ∅	Admissible				
	0,25mm ∅< bis ≤ 0,5mm ∅	2	2	2	3	3
	0,5mm ∅< bis ≤ 0,75 mm ∅	1	1	1	1	1
	>0,75 ∅	0	0	0	0	0
	Maximum number of characteristics	2	2	3	4	4
Smeared print image	No longer legible or print shadow larger than 0.5 mm	Not admissible				
Frayed print edges	Fringes larger than 0.3 mm	Not admissible				
Adhesion Print image	Removal of particles during scrub test	Not admissible				

5.3 Non-quantifiable / Optical Quality Characteristics

The non-quantifiable characteristics are defined using good/bad limit samples. However, limit samples cannot be described or photographed in detail. For a better understanding, we have described a few typical features here:

5.3.1 Melting off characteristics

The melting edge is finally heated in the process in order to equalize the glass threads that have formed at the last point of contact. However, it can happen that the heating is not sufficient to

Printed document is not subject to the change service

completely equalize the glass thread. This can result in so-called melting nodes or sharp spots on the melting edge. Sharp spots that can lead to injuries are generally classified as rejects.


Figure 1: Example of a (tactile) melt-off filament (permissible)

Figure 2: Example of a melted top (not permitted)

Melt-off nodes, which lead to unevenness up to a size of less than 0.5 mm, are permissible; due to the measurement inaccuracy, unevenness up to 1 mm is not accepted as a reason for complaint.

In addition to this, fine bubbles can form in the manufacturing process due to an excess of oxygen in the melting process. The bubbles do not lead to measurable material weakening and can therefore only be classified as an optical feature. Articles with a bubble chain longer than or equal to 1/4 of the article circumference and bubbles larger than 0.4mm in diameter are classified as reject.

5.3.2 Adhesive glass splinters

During hot forming, glass breakage can occur, which adheres to the hot glass and thus poses a risk to the user. Due to the risk of injury, this characteristic is regarded as rejects. Due to the heat in the process, splinters are deformed in most cases and occur thus only as an elevation or bump on the glass surface. In that case, this feature is classified as a surface or down-melting feature.

Printed document is not subject to the change service

5.3.3 Metal/Cinder

When the article is blown out, metal particles and cinder can be applied to the surface through the blow heads or torches (pressing process only). These particles can be surrounded by glass or appear on the surface.

5.3.4 Scratches

Scratches can occur on the surface as a result of an impact or stop.

5.3.5 Thread & Mouth Characteristics (Press-Blowing only)

In the press-blowing process, the mechanical stress in the area of the mouth is particularly high. This can lead to surface characteristics, mouth cracks or deviation of the contour. Due to the process, this characteristic is unavoidable and permissible within the framework of the defined limit samples and gauge. Features such as axial runout and concentricity are described in the technical drawing and are permitted within the given tolerances. Leaning or staggered mouths are not admissible

5.3.6 Surface characteristic tires / roughness

The rotation blowing process can create unevenness in the mould which, in combination with other process parameters, can lead to radial depressions. This phenomenon is called a tyre and can be seen as a single tyre or as overall tyre maturity. Depending on the light source, the tires are difficult to see and cannot be measured directly. Shading of the tires is used as a criterion for better evaluation.




Figure 4: Single tyre (not admissible)

Overall tyre maturity (admissible)

Overall tyre maturity (not admissible)

The extent depends on the mould and the evaluation is in accordance with the limit samples.

In the pressing process, rough areas can occur due to the ingress of glass dust onto the surface. Admissibility is regulated by limit samples.

Printed document is not subject to the change service

5.3.7 Contour characteristics

The contour of the glass container is measured using a contour gauge. The dimensions of the contour gauge correspond to the dimensions determined for the shaping. If there is a gap between the contour gauge and the glass article, the article will be discarded if the gap dimension is greater than or equal to 2 mm. Depending on the application and the later assembly, a closer tolerance is possible. However, any requirements must be marked accordingly in the drawing.

Printed document is not subject to the change service

5.4 Acceptable quality level

DWK Life Sciences differentiates according to critical, major and minor characteristics. These characteristics are valuated according to different AQL classes. The batch size is based on the number of articles per pallet and varies between 1.201 and 3.200. The sample size is typically 125. The characteristics listed below are tested according to these specifications and the measurement methods also listed:

5.4.1 AQL blowing process

Critical	
Residual stress	AQL 0.65
Melted top	
Adhering glass fragments	

Major	
Cracks	AQL 1.0
Breakage	
Melt-off filament	
Zirconium streak	
Metal / cinder	

Minor	
Tolerated dimensions or dimensions marked with min./max.	AQL 2.5
Spout form characteristics	
melting off characteristics (melting off bubbles)	
Print characteristics	
Glass characteristics (bubbles, knots, stones)	
Contour characteristics	
Surface characteristics (tires)	
Scratches	
Dirt	
Complete contents within packaging	

Printed document is not subject to the change service

5.4.2 AQL pressing-blowing process

Critical	
Residual stress	AQL 0.65
Adhering glass fragments (glass in the mould)	

Major	
Cracks	AQL 1.0
Breakage	
Zirconium streak	
Sharp edge mouth	
Metal / cinder	

Minor	
Tolerated dimensions or dimensions marked with min./max.	AQL 2.5
Print characteristics	
Contour characteristics	
Glass characteristics (bubbles, knots, stones)	
Surface characteristics	
Container content	
Scratches	
Dirt	
Thread and Mouth characteristics	
Complete contents within packaging	
Folds/waves	

Printed document is not subject to the change service

5.4.3 AQL pressing process

Critical	
Residual stress	AQL 0.65

Major	
Cracks	AQL 1.0
Breakage	
Zirconium streak	
Sharp edge	
Metal / cinder	

Minor	
Tolerated dimensions or dimensions marked with min./max.	AQL 2.5
Glass characteristics (bubles, knots, stones)	
Contour characteristics	
Surface characteristics	
Scratches	
Dirt	
Edge not full	
Complete contents within packaging	

Printed document is not subject to the change service

6. Packaging

Items will be packed and checked in accordance with the valid DWK Life Sciences palletising standard. This ensures that given proper handling neither damage (within the agreed packaging means) nor loss (within the transport route: production/packaging/storage at DWK Life Sciences – delivery to customer) will occur.

7. Sales Packaging Labelling

Each packaging unit is provided with a label and a DWK Life Sciences pallet document containing the following information:

- Article drawing
- Palett number (DWK – pallet document)
- Article number (DWK –pallet document)
- Order number / Delivery number
- Manufacturing date (DWK –palett document)
- Quantity
- Mould Number
- EAN-Barcode
- Manufacturer

Individual boxes or trays which belong to a particular pallet load are labelled in accordance with the palletising standard relevant to the article.

8. Storage

Recommended storage conditions: protected from moisture, strong temperature fluctuations and direct sunlight.

9. Waste Disposal

Unused DURAN® glass can be classified as AVV No. 200102. DURAN® glass can be disposed of with household waste and can be deposited.

10. Complaints

The contact in the event of complaints is the Quality Management. Complaints must be sent to the following email address: complaint@dwk-lifesciences.com

Complaints can only be made about quality limit values or characteristics not maintained in this DWK Life Sciences specification. Agreement will be necessary between the customer and the DWK Life Sciences in each case prior to processing, as costs may be incurred through reworking, sorting, destruction or return of faulty units if applicable. The following information will be necessary, as well as the reason for the complaint:

- Order number / delivery number
- Complaint quantity in relation to the corresponding delivery quantity
- Receipt date at the customer
- Manufacturing date
- Delivery note number
- Article description
- Article number

Printed document is not subject to the change service

- Pallet number

A sample (free of contamination) should also be provided.

11. Attachments / Associated Documents

- Valid product drawing with corresponding item number
- Valid packaging norm with corresponding item number