Instruction manual

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China

Table of Contents

Introduction 4
Safety information 4
Intended Use 4
Package Contents 5
Installation 5
Environmental Conditions 5
Specifications
Mini Shaker 20 6
Mini Shaker 7
Microplate Shaker 8
Device Setup 9
Control Panel 10
Operating Instructions
Getting Ready 11
Idle Mode 11
Shake Mode 11
Pause Mode 12
Setting Menu 12
RS232 14
RS232 Command Set 15
Troubleshooting 16
Maintenance 16
Standard Tray Configurations 17
Standard Tray Accessories 17
Installation of Sticky Mat 17
Installation of Universal Harness 18
Installation of Adjustable Platform 18
Microplate Shaker Tray Setup 18
Microplate Shaker Tray Configurations 18
Microplate Shaker Tray Accessories 19
Installation of Deep Plate Kit 19
Universal Platform and Accessory Configurations - Mini Shakers 19
Installation of Universal Platform 19
Flask Clamp, Media Bottle, and Tube Rack Configurations for Universal Platforms 20
Installation of Flask Clamps and Media Bottles 20
Installation of Tube Racks 20
Accessories 21
Universal Harness 17
Adjustable Platform 18
Test Tube Racks 21
Flask Clamps 21
Technical Service 22
Warranty 22
Compliance 22
Disposal 23

Introduction

This manual contains installation, operation, and maintenance instructions for the VWR Open Air Shaker Series. Please read the manual completely before using.

Safety Information

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions, and false results.

WARNING	For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
ATTENTION	For important information about the product. May lead to equipment damage if not avoided.
NOTE	For useful information about the product.

Warning Symbols

General hazard

Electrical shock hazard

Safety Precautions

CAUTION: Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

Always operate unit on a level surface for best performance and maximum safety

DO NOT lift unit by the tray.

Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly.

DO NOT immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

WARNING! The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

WARNING! DO NOT use the unit in a hazardous atmosphere or with hazardous materials for which the unit was not designed.

CAUTION: To avoid electrical shock, completely cut of power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.

Earth Ground – Protective Conductor Terminal

Alternating Current

Intended Use

This instrument is intended for use in laboratories. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from VWR, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions,

the intended protection provided by the instrument may be impaired.

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Description	Quantity
Shaker	1
Power Cords: EU cord, UK cord, and a Swiss Plug Adapter	3
Power Supply	1
Non-Skid Rubber Mat *(Mini Shaker only)*	1
Microplate Clips to secure Microplates *(Microplate Shaker Only)*	4

Installation

Upon receiving the VWR Open Air/Microplate Shaker check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage, the carrier must be notified immediately.

After unpacking, place the Open Air/Microplate Shaker on a level bench or table, away from explosive vapors. Secure to an immovable work surface by pressing down on the four (4) corners of the unit, creating a strong suction to the work surface (**DO NOT** place on a bench mat or press on the tray). Do not position the equipment such that it is difficult to disconnect the power cord during use. Always place the unit on a clean and sturdy work surface.

Plug the power supply provided with the Open Air/Microplate Shaker into the power jack on the back of the unit. Only use the power supply provided with the unit. The Open Air/Microplate Shaker is supplied with a 3 conductor, grounded power cord. Connect to the provided power supply to the power cord and plug the power cord into a standard grounded outlet.

If the cord supplied does not meet your needs, please use an approved power cord that has ratings equal to or exceeding those of the originally provided cord and that complies with the local/national regulations of the country in which the equipment is to be used.

Environmental Conditions

Indoor Use Only:

Altitude: 0 to 2000m

Operating Temperature: -10°C to 40°C

Storage Temperature: -20°C to 65°C

Operating Humidity: 20 to 80% relative humidity, non-condensing

Storage Humidity: 20 to 80% relative humidity, non-condensing

Electrical Supply: 24VDC, 1.25A. (For use with certified or approved power supply, which must have a SELV and limited energy output.).

Voltage fluctuations: Mains supply voltage fluctuations up to ±10% of the nominal voltage.

Overvoltage category (Installation category): II

Pollution degree: 2

Avoid cold starts. Unit is not designed to start after being in a cold room environment. Bring unit into cold room from a room temperature environment, operate and remove unit from cold room as soon as operation is complete.

Product Specifications


Mini Shaker 20

Environmentally Preferable	Waste Reducing	
Motion	Orbital	
Orbit	20mm	
Display	75 x 20mm White Character LCD	
Weight Capacity	4kg	
Flask Clamp/ Tube Rack Capacity	**Standard Tray** 7 x 10 mL Flasks 6 x 125 mL Flasks 6 x 250 mL Flasks 4 x 500 mL Flasks 2 x Tube Racks	**Universal Tray** 39 x 10 mL Flasks 12 x 125 mL Flasks 8 x 250 mL Flasks 6 x 500 mL Flasks 3 x 1 L Flasks*** 2 x Tube Racks
Speed	50 – 300 RPM	
Speed Increments	1 RPM	
Speed Accuracy**	Above 100 rpm ± 2% of set speed, below 100 rpm ± 2 rpm	
Overall Dimensions (L x W x H)	43,2 x 26,2 x 11,9 cm	
Tray/ Platform Dimensions	30 x 22,2 cm	
Tray Construction	Aluminum tray with 11 mounting points for accessories. Non-slip rubber mat included	
Tray/ Platform Options	Universal Platform & Adjustable Platform	
Timer	Up to 99 hours 59 minutes	
Recall Feature	Retrieve last 5 recently used settings	
Data Logging & Remote Control	Built in RS232	
Audible Alarm	Yes	
Overload Protection	Yes	
Ingress Protection	IP21	
Operating Temperature Range	-10 °C to 40 °C, 80% RH, non-condensing	
Altitude	2000m	
Motor Type	Brushless DC Motor	
CO_2 Environment	Yes	
Safety Certification	CE Marked and TÜV SÜD listed	
Electrical	24VDC/ 1,25A	
Net Weight	8,5 kg	

**Conditions permitting. Variations in measurement process, vessel(s), and sample(s) will impact actual performance.

***Universal Platform 444-1425 required to mount 3 x 1 L flasks capacity. Note: 3 x1L flask without clamps are compatible with rubber mat.


Mini Shaker

Environmentally Preferable	Waste Reducing	
Motion	Orbital	
Orbit	3mm	
Display	75 x 20 mm White Character LCD	
Weight Capacity	4 kg	
Flask Clamp/ Tube Rack Capacity	**Standard Tray** 7 x 10 mL flasks 6 x 125 ml flasks 6 x 250 ml flasks 4 x 500 ml flasks 2 x tube racks	**Universal Platform** 39 x 10 mL flasks 12 x 125 ml flasks 8 x 250 ml flasks 6 x 500 ml flasks 3 x 1 L flasks*** 2 x Test tube racks
Speed	100 -1200 RPM	
Speed Increments	1 RPM	
Speed Accuracy**	Above 100 rpm ± 2% of set speed	
Overall Dimensions (L x W x H)	43,2 x 26,2 x 11,9 cm	
Tray/ Platform Dimensions	30 x 22,2 cm	
Tray Construction	Aluminum tray with 11 mounting points for accessories. Non-slip rubber mat included	
Tray/ Platform Options	Universal Platform & Adjustable Platform	
Timer	Up to 99 hours and 59 minutes	
Recall Feature	Retrieve last 5 recently used settings	
Data Logging & Remote Control	Built in RS232	
Audible Alarm	Yes	
Overload Protection	Yes	
Ingress Protection	IP21	
Operating Temperature Range	-10 °C to 40 °C 80% RH, non-condensing	
Altitude	2000 m	
Motor Type	Brushless DC Motor	
CO_2 Environment	Yes	
Safety Certification	CE Marked and TÜV SÜD listed	
Electrical	24VDC / 1,25A	
Net Weight	8,5 kg	

**Conditions permitting. Variations in measurement process, vessel(s), and sample(s) will impact actual performance.
***Universal Platform 444-1425 required to mount 3 x 1 L flasks capacity. Note: 3 x1L flask without clamps are compatible with rubber mat.


Microplate Shaker


Environmentally Preferable	Waste Reducing
Motion	Orbital
Orbit	3mm
Display	75 x 20 mm White Character LCD
Weight Capacity	1kg
Volume Capacity	4 x Microplates 2 x Microtube Racks 4 x Deep well Plates
Speed	100 -1200 RPM*
Speed Increments	1 RPM
Speed Accuracy**	Above 100 rpm ± 2% of set speed
Overall Dimensions	43,2 x 26,2 x 11,9 cm
Tray/ Platform Dimensions	28 x 19,7 cm
Tray Construction	Aluminum tray with hardware to secure Microplates and Deep well plates
Timer	Up to 99 hours 59 minutes
Recall Feature	Retrieve last 5 recently used settings
Data Logging & Remote Control	Built in RS232
Audible Alarm	Yes
Overload Protection	Yes
Ingress Protection	IP21
Operating Temperature Range	-10 °C to 40 °C, 80% RH, non-condensing
Altitude	2000 m
Motor Type	Brushless DC Motor
CO_2 Environment	Yes
Safety Certification	CE Marked and TÜV SÜD listed
Electrical	24VDC / 1,25A
Net Weight	8,5 kg

*Microplates and deep well plates are available with a variety of skirt styles. Depending on the skirt style the maximum speed may need to be reduced.
**Conditions permitting. Variations in measurement process, vessel(s), and sample(s) will impact actual performance.

Device Setup

Mini / Microplate Shaker

A. Time Display

B. Speed Display

C. Timer Status Indicator

D. Run/Speed Status Indicator

E. Timer Adjustment Knob

F. Speed Adjustment Knob

G. Power Switch

H. RS232 Port

I. DC Power Connector

Digital Control Panel


Figure 1: Mini and Microplate Digital Control Panel

A. Time Display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 99 hours 59 minutes. The display will indicate minutes and seconds until the timer reaches 59 minutes and 59 seconds (59:59), then the display will automatically display hours and minutes up to 99:59.

B. Speed Display: Displays the speed of the shaker. Turn the right knob (E) for setpoint control. Press knob (E) starts/stops shaking function. The run/speed indicator light (C) will be illuminated when the unit is shaking.

C. Time Status Indicator Light

Timed Mode: Indicator light will be illuminated while the time display is counting down. If the unit is paused while shaking with the timer set, then the indicator light will blink.

Continuous Mode: Timer is not set; the indicator light will be off during shaking while elapsed time counts up.

D. Run Status Indicator Light

E. Timer/ Menu Adjustment Knob

F. Speed Adjustment Knob

Operating Instructions

1.Getting Ready:

1. Connect the provided power supply to the unit.
2. Plug the female end of the provided power cord into the power supply and plug the male end of the power cord into a matching standard grounded outlet.
3. Toggle the power switch to the ON (•) position.
4. Allow unit to complete the start-up sequence:
 a. LCD Test – all display icons illuminate
 b. Model Type, software version
 c. Function Units Timer [Hours: Minutes], Shake Speed [RPM]
5. The main operating screen will return, and the previous timer and shaker settings will be displayed.

2.Idle Mode

Timer

1. By default, the timer will be set to 00:00 and count upwards when the shaking function is turned on.
 a. The timer begins in MM:SS mode
 b. Once the timer reaches 59:59 it will automatically switch to HH:MM mode at 01:00.
2. When the timer is set to any time between 00:01 and 99:59, the timer will count downwards from the set time once the shaking function is turned on and the ***Timer Status indicator Light (A)*** will be illuminated.
 a. The countdown timer can only be set in HH:MM mode when counting down.
 b. Once the timer reaches 01:00 in HH:MM mode, it will automatically switch to MM:SS mode at 59:59.
 c. Once the timer reaches 00:00 in the countdown mode shaking functions will shut off, the unit will beep 4 times and "End" will appear on the display (**A, Figure 1**).

Adjust Timer

1. Rotate the ***Time Adjustment Knob*** (**E, Figure 1**).
 a. The minutes will blink on the display.
 b. Rotate the knob to adjust the minutes value
 c. Quick press the ***Time Adjustment Knob*** (**E**) to select the minutes value
 d. The hours will now blink on the display.
 e. Rotate the ***Time Adjustment Knob*** (**E**) to adjust the hours value
 f. Quick press the ***Time Adjustment Knob*** (**E**) to confirm hours value and the timer setting

Reset Timer to Zero

1. Ensure shaker is not running, then rotate or quick press ***Time Adjustment Knob*** (**E**) so either minutes or hours are blinking
2. Long press ***Time Adjustment Knob*** (**E**) until 00:00 is on timer display

Adjust Set Speed

1. Rotate the ***Speed Adjustment Knob*** (**F, Figure 1**) to the desired setting.

3. Shake Mode

• Start Shake Function
 1. Long Press the ***Speed Adjustment Knob*** (**F**) until the unit beeps and the Run Status Indicator (D, Figure 1) illuminates.

• Adjust Speed while shaking
 1. Rotate the ***Speed Adjustment Knob*** (**F**)
 2. The ***Run Status Indicator Light*** (**D**) will blink with the "set" icon on the display indicating the unit is ramping to the new speed setting.
 3. Once the unit reaches the new speed setting the ***Run Status Indicator Light*** (**D**) will be illuminated and stop blinking. The "set" icon with disappear.

• Adjust Remaining Time *(count down only when timer is set)*
 1. Rotate the ***Time Adjustment Knob*** (**E**) to adjust the minutes value. The minutes value will blink while making adjustment.
 2. Quick press the ***Time Adjustment Knob*** (**E**) to confirm the minutes value, and now the hours value will blink.
 3. Rotate the ***Time Adjustment Knob*** (**E**) to adjust the hours value
 4. Quick press the ***Time Adjustment Knob*** (**E**) to confirm hours value and the new timer setting

NOTE: When making adjustments and the minute or hour is blinking; if the time setting remains idle without confirmation for 6 seconds it will reset to the current setting.

- Stop Shake

1. Long press the ***Speed Adjustment Knob*** (**F**) until the unit beeps and the ***Run Status Indicator Light*** (**D**) turns off

4. Pause Mode

The intent of Pause mode is to allow the user to stop the shaking function while maintaining the runtime/remaining time value. The display will show "PAUS", the ***Run Status Indicator*** **(D)** remains off while the ***Timer Status Indicator*** (**C**) blinks to indicate a paused timer. Pause mode must be enabled within the Settings Menu to allow access to this function.

- Pause Shake (when Pause feature is enabled in the Menu)

1. Quick press the ***Speed Adjustment Knob*** (**F**)
2. Unit will stop shaking and "PAUS" will appear on the display and stop the timer

- Adjust Remaining Time (count down only)

1. Rotate the ***Time Adjustment Knob*** (**E**) to adjust the minutes value. The minutes value will blink while making adjustment.
2. Quick press the ***Time Adjustment Knob*** (**E**) to confirm the minutes value, and now the hours value will blink.
3. Rotate the ***Time Adjustment Knob*** (**E**) to adjust the hours value
4. Quick press the ***Time Adjustment Knob*** (**E**) to confirm hours value and the new timer setting

- Resume Shake Function

1. Quick press the ***Speed Adjustment Knob*** (**F**)
2. The unit will start shaking, the ***Run Status Indicator Light*** (**D**) and the ***Timer Status Indicator Light*** **(C)** will illuminate.

- Stop Shake

1. Long press the ***Speed Adjustment Knob*** (**F**) until the unit beeps and the ***Run Status Indicator*** (**D**) turns off

Operating Tips

- Centering your sample and even weight distribution on the tray or platform helps with balance and stability.
- Ensure the unit is secure on a level work surface by pressing down on the four (4) corners of the unit, creating a strong suction to the work surface.
- As a safety feature, a built-in program will shut power off to the motor if the tray is prevented from rotating, or the unit is overloaded beyond its recommended weight capacity and display an E4 error.
- If the load is within specifications and centered on the tray and an E4 error persists, disabling the Overload Detection feature in the Menu may resolve this fault.
- Built-in memory maintains the last used speed and time settings during a power interruption.
- **Avoid cold starts.** Unit is not designed to start after being in a cold room environment. Bring unit into cold room from a room temperature environment.

5.Setting Menu

Menu Access

- While in Idle Mode, press and hold the ***Time Adjustment Knob*** (**E**) until "MENU" is displayed

Navigation

- Rotate either the ***Time Adjustment Knob*** (**E**) or ***Speed Adjustment Knob*** (**F**)
- Quick press either the ***Time Adjustment Knob*** (**E**) or ***Speed Adjustment Knob*** (**F**) to access submenu/select option
- To exit the Menu at any time, toggle the power switch.


Figure 2: Mini/ Microplate Shaker Settings Menu Structure; default settings underlined.

Recall

- While in Idle Mode, press and hold the ***Time Adjustment Knob*** (**E**) until "MENU" is displayed
- After "MENU" is displayed "RECALL" will appear on the display.
- Recall previous run conditions
 - Quick press either the ***Time Adjustment Knob*** (**E**) or ***Speed Adjustment Knob*** (**F**).
 - Rotate either knob to navigate through the previous 5 unique run conditions
 - To select a recall run condition, press, and hold the ***Speed Adjustment Knob*** (**F**) until the unit beeps and the **Run Status Indicator** (**D**, Figure 1) illuminates.
 - Unit will now run at the recall timer and speed setting that was selected.

Menu – System Settings

Pause Enable/Disable

- Navigate to the System Settings submenu ("SYS")
- Navigate to the Pause setting ("PAUS")
- Quick press to select value
- Rotate knob to select the desired setting
- Quick press to confirm value (En = Enabled or dIS = Disabled)

Beeper Enable/Disable

- Disabling the Beeper will silence and prevent beeps in the following scenarios: o Starting and stopping shaking
 - When the timer reaches 00:00 in timer countdown mode
- Navigate to the System Settings submenu ("SYS")
- Navigate to the Pause setting ("BEEP")
- Quick press to select value
- Rotate knob to select the desired setting (En = Enabled or dIS = Disabled)
- Quick press to confirm value

Power Recovery Enable / Disable

- Power Recovery is an optional feature that allows the unit to automatically restart the shaking functions when power is returned to the unit after a disconnect.
- Navigate to the System Settings submenu ("SYS")
- Navigate to the Power Recovery setting ("PWr")
- Quick press to select value
- Rotate knob to select the desired setting (En = Enabled or dIS = Disabled)
- Quick press to confirm value

Overload Detection Enable / Disable

- Optional load sensing feature which provides protection against improper positioning and excess load conditions.
- The unit will automatically slow to a safe operating speed and display "E4" if a hazardous load is detected. Error may be cleared by interacting with the unit; shake and timer functions will continue operating as standard.
- Navigate to the System Settings submenu ("SYS")
- Navigate to the Power Recovery setting ("LOAD")
- Quick press to select value
- Rotate knob to select the desired setting (En = Enabled or dIS = Disabled)
- Quick press to confirm value

Display Settings

- Adjust LCD brightness and configure the display timeout feature, which dims the LCD to set brightness after a selected period of idle time.
- Navigate to the System Settings submenu ("SYS")
- Navigate to the Power Recovery setting ("dISP")
- Adjust Backlight Brightness
 - Adjust the brightness percentage ("bL P")
 - Range: 0-100 %
 - Quick press to select the back light brightness percentage value
 - Rotate the knob to select the desired setting
 - Quick press to confirm value
- Backlight Timeout Enable/Disable
 - Quick press to select value ("bLtO")

- Rotate knob to select the desired setting (En = Enabled or dIS = Disabled)
- Quick press to confirm value

- When Backlight Timeout Brightness is Enabled, adjust Backlight Timeout Period
 - Quick press to select value (“tO t”)
 - Range: 1-4 hours
 - Rotate knob to select the desired setting (select number of hours)
 - Quick press to confirm value
- Adjust backlight brightness while timeout is active:
 - Adjust the brightness percentage (“tO P”)
 - Range: 0-100 %
 - Quick press to select the back light brightness percentage value
 - Rotate the knob to select the desired setting
 - Quick press to confirm value

Factory Reset

- Resetting the unit to Factory Default Settings will do the following:
 - Clear all recall parameters
 - Clears all System settings to the default value
- Navigate to the System Settings submenu (“SYS”)
- Navigate to the Reset setting (“rSt”)
- Quick press to select value (no or yes)
- Rotate knob to select the desired setting
- Quick press to confirm

6. RS232

RS-232 serial port provides two-way communications for data logging and unit control by means of a PC and a suitable application program.

- Interface connections chosen in accordance with EIA standard RS-232
- Device Connector: RS-232 DB9 Female
- Recommended Cables:
 - DB9 M/M Straight Through Serial Cable
 - USB-A to RS-232 DB9 Straight Through Serial Adapter Cable

- Transmission Procedure
 - Asynchronous character transmission in start-stop mode

- Transmission Type
 - Full Duplex

- Character Format
 - Start Bit(s): 1
 - Character Bits: 8
 - Parity Bit(s): None
 - Stop Bit(s): 1

- Transmission Speed (Baud Rate): 9600
- Data Flow Control: Xon/Xoff
- Syntax
 - Instructions and parameters separated by space (0x20)
 - CR LF termination (0x0D, 0x0A)
 - Max length: 80 characters
- ‘<command> A’ returned if the command is recognized; ‘L’ returned otherwise.

RS232 Command Set

ID <XXXX>	Set ID 1 ≤ XXXX ≤ 9999
ID	Read ID value
MODE	Read mode: 0: idle 1: shake 2: paused 99: error
START_SHAKE	Start shake function and timer
PAUSE_SHAKE	Pause shake function and timer
STOP_SHAKE	Stop the shake function; reset timer
TARGET_SPEED <xxxx>	Set the target speed value
TARGET_SPEED	Read target speed value
MEASURED_SPEED	Read measured speed value
TIMER	Read current timer value
TIMER 0	Set Timer Mode = Count Up
TIMER <HH>:<MM>:<SS>	Set Timer Mode = Count Down; Timer = HH:MM:SS
TIMER RESET	Reset timer to set value; timer continues counting. No change to shake function.
LOCK	Lock user interface
UNLOCK	Unlock user interface
PARAM <XXXX>	Enable output of comma separated list of parameters, formatted as follows: <TIMER [HH:MM:SS]>, <ID>, <MODE>, <TARGET SPEED>, <MEASURED SPEED>, <ERROR CODE>, Output can be periodic, where XXXX represents the desired time interval: XXXX = 0 Single output 1 ≤ XXXX ≤ 9999 Periodic output every XXXX seconds

Troubleshooting

Review the information in the table below to troubleshoot operating problems.
During operation, any rattling or ticking sounds may indicate a loose screw on the tray, a tray attachment, or an accessory. All accessories should be sufficiently tightened in place before starting the unit.
The following table lists error codes associated with common problems and possible causes and remedies. If the problem persists, contact VWR or your authorized dealer.

Problem	Cause	Solution
E3	Speed Error No/ Low Tray Speed, excess motor speed	Broken drive belt, malfunctioning speed sensor, seized bearing, mechanical obstruction, tray overload. Reset the power switch to clear this error. Check for mechanical obstructions. If the E3 error persists, remove power from the unit and contact a VWR representative for repair. End user should NOT attempt repair
E4	Load Error Uneven / excessive tray loading, high vibration, motor overload	Uneven / excessive loading, loose suction cup foot*, mechanical obstruction Toggle the rocker switch to clear this error. Reduction of load / even load distribution may resolve this error. **Alternatively, disabling the Overload Detection feature may resolve this fault, please see the Instruction Manual for** Entering the Menu and disabling this feature. If the E4 error persists, remove power from the unit and contact a VWR representative for repair. End user should NOT attempt repair.

***In the event a foot (suction cup) has come loose from the bench top, the unit may register an errant E4 error message due to the instability of the unit. Reset the power switch to clear this error. Firmly press down on the four (4) corners of the unit, creating a strong suction to the work surface (do not place on bench mat).**

Maintenance

The Open Air/Microplate Shaker is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. It needs no user maintenance beyond keeping the unit clean. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

The Open Air/Microplate Shaker should be given the care normally required for any electrical appliance. At least every three (3) months you should:

- Unplug the unit.
- Remove any accumulated dirt from the work surface and tray.
- Check all accessory items to make sure they are properly tightened.


Attention. Do not use solvents, harsh chemicals, ammonia, or abrasive cleaning agents The housing may be cleaned with a cloth dampened with a mild detergent if necessary.
Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly after the unit has cooled down. Before using any cleaning or decontamination method except as noted in this section, users should check with the manufacturer that the proposed method will not damage the equipment. The user is responsible for carrying out appropriate decontamination if hazardous material is spilled onto or into the equipment.


Mini Shaker Set-Up Standard Tray Configurations

The Open-Air Shaker is supplied with a tray designed to hold a variety of accessories (available separately).


1. Flat containers can be shaken by placing them on the non-skid rubber mat provided with the unit.
2. The tray also has mounting holes ready for use with flask clamps or test tube racks. See below examples of tray configurations

(6) Erlenmeyer Flask Clamps
125 ml – 250 ml

(4) 500 mL Erlenmeyer Flask Clamps

(2) Tube Racks

Standard Tray Accessories for Mini Shakers (available separately)

Universal Platform
444-1425

Sticky Mat
444-1417

Universal Harness OA
(4 short bands, 2 long bands)
444-1421

Adjustable Platform
444-0282

Installation of Sticky Mat

Sticky mats are recommended for flat bottom vessels with a wide base at low speeds. Sticky mat performance may diminish over time due to spills, dirt, and debris.

1. Turn unit off with power switch.
2. Remove non-skid rubber mat.
3. Place sticky mat directly on tray, then place sample on sticky mat. The sticky mat can be repositioned simply by peeling from the tray.
4. Vessels can be removed by tilting to the side

Note: Recommended maximum speed not exceed 250rpm.

Installation of Universal Harness

The Universal Harness attaches to tray to secure low-profile plates.

1. Turn unit off with power switch.
2. Remove non-skid rubber mat.
3. To use the short bands, install the center post thumbscrew on the tray. Hand tighten into place **DO NOT** over tighten,
4. Place sample(s) on tray.
5. Place the short band under the bent down corner of the tray and stretch the short band over the sample(s) and then place on the center post thumbscrew.
6. To use the long bands, place sample on tray, place the long band under the bent down corner of the tray and stretch diagonally over the sample and secure to the opposite bent down corner.

Not for use with the Microplate Shaker.

Installation of Adjustable Platform

The Adjustable Platform has a non-skid rubber surface with adjustable clamping bars to accommodate various vessel types. Stainless steel construction.

1. Turn unit off with power switch.
2. Remove non-skid rubber mat.
3. Line up the holes in the Universal Platform to the holes in the tray. Using the screws provided, hand tighten into place using a Phillips screwdriver. **DO NOT** over tighten.
4. Add glassware. Adjust rollers by loosening each of the thumbscrews. Slide the padded rollers to the desired position and retighten the thumbscrews.

Not for use with the Microplate Shaker.

Microplate Shaker Tray Setup


The Microplate Shaker is designed to hold two (2) or four (4) microplates, or two (2) Micro-Tube Racks.

1. Turn unit off with power switch
2. Place two (2) microplates or deep well blocks diagonally on the tray, or place four (4) microplates or deep well blocks on the tray. The plates/blocks do not have to be filled.
3. Slide plate/block into place. You are ready to use.


NOTE: The tray on the Microplate Shaker is not designed to hold flask clamps

Microplate Shaker Tray


2.


3.


Configurations

(2) Microplates


(4) Microplates

(2) Micro-Tube Racks

Note: The tray on the Microplate Shaker is not designed to hold Flask Clamps.

Microplate Shaker Tray Accessories (available separately)

Deep Plate Kit, set of 4
444-1416

Installation of Deep Plate Kit, set of 4

1. Turn unit off with power switch
2. With a screwdriver remove V-Clips from top plate
3. Line up the Deep Plate Adapter with the openings on the top plate and screw in with the thumbscrew.
4. Place deep well plates on tray. If more tension is required tighten the thumbscrew to secure the deep well pate.

Universal Platform and Accessory Configurations for Mini Shakers

The optional Universal Platform 30 x 22 cm 444-1425 accessory is required for configurations listed below. All flask clamps, medial bottle clamps and tube racks are sold separately.


Installation of Universal Platform

1. Turn unit off with power switch.
2. Remove non-skid rubber mat.
3. Loosen the four (4) thumbscrews on the shake tray.
4. Place the universal platform on the shaker tray by aligning the slots on the two sides of the platform over the four (4) screw knobs.
5. Allow the platform to rest on the screws.
6. When aligned and seated on the screws hand tighten the four (4) thumbscrews. **DO NOT** overtighten.


Flask Clamp, Media Bottle and Tube Rack Configurations for Optional Universal Platform


Turn unit off with power switch before installing accessories on the Universal Platform. See below examples of Universal Platform configurations.


(39) 10mL Erlenmeyer Flask Clamps


(24) 25mL Erlenmeyer Flask Clamps


(24) 50mL Erlenmeyer Flask Clamps


(12) 125mL Erlenmeyer Flask Clamps


(12) 125mL Erlenmeyer PVC Flask Clamps


(8) 250mL Erlenmeyer Flask Clamps


(8) 250mL Erlenmeyer PVC Flask Clamps


(6) 500mL Erlenmeyer Flask Clamps


(6) 500mL Media Bottle Clamps

(4) 500mL Erlenmeyer PVC Flask Clamps

(3) 1L Erlenmeyer Flask Clamps

(2) Tube Racks

Installation of Flask Clamps and Media Bottles (Stainless Steel)

1. Turn unit off with power switch
2. Line up the opening(s) in the clamp to the openings in the platform. Hand tighten screw(s) with a Phillips head screwdriver.
3. Insert the flask or media bottle into the clamp and wrap the supplied spring around the neck of the clamp to hold secure

Installation of Flask Clamps (PVC)

1. Turn unit off with power switch
2. Line up the opening(s) in the clamp to the openings in the platform. Hand tighten the screw(s) on top of the washer with a flat head screwdriver.
3. Insert the flask into the clamp.

Installation of Tube Racks on Universal Platform

1. Turn unit off with power switch.
2. Place one (1) adaptor bracket on the Universal Platform with the screw hole openings side facing the side of the platform. The flat wall side of the adapter bracket will be on the inside.
3. Place the tube rack with the clip side facing the flat wall side of the adaptor bracket.
4. Align the two (2) adaptor bracket screw hole openings and the two-tube rack screw hole openings with four (4) compatible openings in the platform.
5. Ensure there is tension on the tube rack so it will stay in position. If the tension is too tight or too loose, bend the clip accordingly.
6. Using the four (4) 3/8” screws, first secure the adaptor bracket.
7. Ensure there is tension on the test tube rack clip and then install the two (2) screws to secure the tube rack to the platform.

Accessories

Cat. No.	Description	Standard Tray (3mm/20mm)	Universal Platform (444-1425)	Microplate Shaker
444-7040	10mL Erlenmeyer Flask Clamp	7	39	-
444-7041	25mL Erlenmeyer Flask Clamp	7	24	-
444-7042	50mL Erlenmeyer Flask Clamp	7	24	-
444-7043	125mL Erlenmeyer Flask Clamp	6	12	-
444-7044	250mL Erlenmeyer Flask Clamp	6	8	-
444-7045	500mL Erlenmeyer Flask Clamp	4	6	-
444-7030	1 Liter Erlenmeyer Flask Clamp	-	3	-
444-7035	500mL Media Bottle Clamp	4	6	-
444-2960	125mL Erlenmeyer Flask Clamp, PVC	6	12	-
444-2961	250mL Erlenmeyer Flask Clamp, PVC	6	8	-
444-2962	500mL Erlenmeyer Flask Clamp, PVC	-	4	-
444-0276	1.5 to 2mL Micro-Tube Rack	2	2	2
444-0278	10 to 13mm Test Tube Rack	2	2	-
444-0279	14 to 16mm Test Tube Rack	2	2	-
444-0280	18 to 20mm Test Tube Rack	2	2	-
444-0281	22 to 25mm Test Tube Rack	2	2	-
444-0283	15mL Centrifuge Tube Rack	2	2	-
444-0284	50mL Centrifuge Tube Rack	2	2	-
444-1425	Universal Platform	1	-	-
444-0282	VWR Adjustable Platform	1	-	-
444-1417	Sticky Mat 22 x 18 cm	1	-	-
444-1416	Deep Plate Kit, set of 4	-	-	1
444-1421	Universal Harness	1	-	-

Test Tube Racks

Installation of Test Tube Racks - Mini Shaker

1. Unit should be in standby mode.
2. Remove non-skid rubber mat.
3. Align the clip on the rack to the side of the tray.
4. Line up the holes in the test tube rack to the holes in the tray. Using the screws provided, hand tighten into place using a Phillips screwdriver. **DO NOT** over tighten.

Installation of 1.5 to 2mL Micro-Tube Rack - Microplate Shaker

1.Unit should be in standby mode.
2. Remove V-Clips from tray.
3. Align the clip on the rack to the side of the tray.
4. Line up the holes in the micro-tube rack to the holes in the tray. Using the screws provided, hand tighten into place using a Phillips screwdriver. **DO NOT** over tighten.

Flask Clamps

Installation of Flask Clamps - Mini Shaker

1. Unit should be in standby mode.
2. Remove non-skid rubber mat.
3. Line up the hole(s) in the flask clamp to the hole(s) in the tray. Using the screw(s) provided, hand tighten into place using a Phillips screwdriver. **DO NOT** over tighten.
4. Insert flask/media bottle into clamp and wrap the spring around the neck of the clamp to hold secure.

Not for use with the Microplate Shaker.

Technical service

Web Resources

Visit the VWR website at **www.vwr.com** for:

- Complete technical service contact information
- Access to the VWR Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty

VWR warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation unless any refusal is due to a defect of the product.

Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
CE	VWR Corporation declares that the Open Air Mini and Microplate Shakers comply with directives 2011/65/EU, (EU) 2015/863, 2014/30/EU, 2014/35/EU, and standards EN 50581, EN 61010-2-010, EN 61010-2-051, EN 61326-1. The full text of the EU declaration of conformity is available at the following internet address: https://vwr.com
	This product complies with directive 2012/19/EU. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to https://vwr.com
	EN 61326-1
TÜV SÜD C US	CAN/CSA C22 261010-1, CAN/CSA C22 261010-2-010, CAN/CSA C22 261010-2-051 UL 61010-1, UL 61010-2-010, UL 61010-2-051

Global Notice

Warning: This is a Class B product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

Canada Notice


This Class B digital apparatus complies with Canadian ICES-003.

FCC Notice

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

—Reorient or relocate the receiving antenna.
—Increase the separation between the equipment and receiver.
—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
—Consult the dealer or an experienced radio/TV technician for help.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.
Instead, it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.
For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.
By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.
Thank you

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INDIA
Avantor Performance Materials India Pvt. Ltd.,
19 th Floor, Building No. 5, Tower C,
DLF Cyber City Phase - III,
Gurgaon - 122002, Haryana
Tel.: +91-124-4656700
csindia@avantorsciences.com

IRELAND
VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
sales.ie@vwr.com

ITALY
VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
info.it@vwr.com

KOREA
Avantor Performance Materials Korea Ltd
2F ACE Gwanggyo Tower I, Daehak 4ro 17
Yeongtong-gu Suwon, Korea 16226
Tel.: +82 31 645 7250
customerservicekr@avantorsciences.com

THE NETHERLANDS
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
info.nl@vwr.com

MEXICO
VWR International, S.de R.L. de C.V.
Km. 14.5 Carretera
Tlalnepantla-Cuautitlán
Col. Lecheria
Tultitlán Edo. de México
CP 54940
Tel.: +52 (55) 5005 0100
vwrmx@vwr.com

MIDDLE EAST & AFRICA
VWR International FZ-LLC
Office 203, DSP Lab Complex,
Dubai Science Park,
Dubai, United Arab Emirates
Tel: +971 4 5573271
Info.mea@vwr.com

NORWAY
VWR International AS
Brynsalleen 4,
0667 Oslo
Tel.: +47 22 90 00 00
info.no@vwr.com

POLAND
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
info.pl@vwr.com

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GO TO AVANTORSCIENCES.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS FROM YOUR LOCAL TEAM

VWR® Mini and Microplate Shakers
Instruction Manual
31018502

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VWR International GmbH
VWR International, LLC