VWR®

Instruction Manual
Analog Dry Block Heaters
Digital Dry Block Heaters

EN - English 1
FR - Français 17
ES - Español 34
IT - Italiano 51
DE - Deutsch 68
PT - Português 85
NL - Nederlands 102
NO - Norsk 107
DA - Dansk 112
SV - Svenska 117
FI - Suomi 122
HU - Magyar 127
PL - Polski 132

European Catalogue Numbers:

Analog Dry Block Heaters

Plug	1 Block	2 Block	3 Block	4 Block	6 Block
Euro	460-3249	460-3252	460-3255	460-3258	460-3261
UK	460-3250	460-3253	460-3256	460-3259	460-3262
Swiss	460-3251	460-3254	460-3257	460-3260	460-3263

Digital Dry Block Heaters

Plug	1 Block	2 Block	3 Block	4 Block	6 Block
Euro	460-3264	460-3267	460-3270	460-3273	460-3276
UK	460-3265	460-3268	460-3271	460-3274	460-3277
Swiss	460-3266	460-3269	460-3272	460-3275	460-3278

Digital Dry Block Heaters with Heated Lid

Plug	2 Block
Euro	460-0183
UK	460-0184
Swiss	460-0185

Revision 4
09/04/13

CE

Table of Contents

Package Contents 1
Warranty 1
Installation 2
Maintenance & Servicing 2
Intended Use 2
Equipment Disposal 2
Environmental Conditions 2
Safety Instructions 3
Warning 3
Standards & Regulations 3
Specifications - Analog Dry Block Heater 4
Operating Instructions - Analog Dry Block Heater. 5
Digital Control Panel 6
Specifications - Digital Dry Block Heater 7
Operating Instructions - Digital Dry Block Heater. 8-12
Technical Service 13
Troubleshooting 13
Accessories 14-16
EC Declaration of Conformity 137-141

Package Contents

Dry Block Heater
234cm detachable power cord
Instruction manual

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, VWR will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Legal Manufacturer
VWR International bvba • Researchpark Haasrode 2020 • Geldenaaksebaan 464 • B-3001 Leuven • + 32 16 385011 • http://be.vwr.com

Installation

Upon receiving the VWR Dry Block Heater, check to ensure that no damage has occurred in shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the Dry Block Heater on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of 15cm from vertical surfaces. Always place the unit on a sturdy work surface.

The Dry Block Heater is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The 230V unit plugs into a 230 volt, 50/60 Hz source.

It is necessary to fill the Dry Block Heater well(s) with modular heating blocks because empty block locations will affect performance. Place filled tubes in modular block(s), then place the modular block(s) into the Dry Block Heater well(s).

Maintenance & Servicing

The Dry Block Heater is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. It needs no user maintenance beyond keeping the surfaces clean. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly after the unit has cooled down. **DO NOT** use a cleaning agent or solvent on the front panel which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

Intended Use

These Dry Block Heaters are intended for general laboratory use.

Equipment Disposal

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Environmental Conditions

Operating Conditions: Indoor use only.

Temperature:	18 to 33°C
Humidity:	20% to 80% relative humidity, non-condensing
Altitude:	0 to 2000 M above sea level

Non-Operating Storage:

Temperature:	-20 to 65°C
Humidity:	20% to 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

Safety Instructions

Please read the entire instruction manual before operating the Dry Block Heater.

WARNING! DO NOT use the Dry Block Heater in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly, after the unit has cooled down. **DO NOT** immerse the unit for cleaning. **DO NOT** operate the unit if it shows signs of electrical or mechanical damage.

The Dry Block Heaters are designed to be operated in dry conditions. **DO NOT** put water, oil or other fluids in the wells of the units. The chamber that the modular heating blocks and baths sit in is not designed to be filled with liquid or other fluids. **DO NOT** place anything other then the appropriate heating block(s) in this cavity.

CAUTION! Units can reach temperatures of up to 150˚C. Use caution at all times.

WARNING! Dry Block Heaters are not explosion proof. Use caution when unit is on or when heating volatile materials.

Earth Ground - Protective Conductor Terminal

Alternating Current

Standards & Regulations

VWR international hereby declares under its sole responsibility that the products comply with the requirements of the following directives and associated standards:

Associated EU directives:

EMC directive	2004/108/EC
LVD directive	2006/95/EC

Safety standards:

EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use. Part I: General Requirements.
EN 61010-2-010	Part II: Particular requirements for laboratory equipment for the heating of materials.

EMC standards:

IEC 61326-1	IEC 61000-4-4
IEC 61000-3-2	IEC 61000-4-5
IEC 61000-3-3	IEC 61000-4-6
IEC 61000-4-2	IEC 61000-4-11
IEC 61000-4-3	


Analog Dry Block Heater Specifications

	1 Block	2 Block	3 Block	4 Block	6 Block
Cat. No. with Euro plug	460-3249	460-3252	460-3255	460-3258	460-3261
Cat. No. with UK plug	460-3250	460-3253	460-3256	460-3259	460-3262
Cat. No. with Swiss plug	460-3251	460-3254	460-3257	460-3260	460-3263
Dimensions (L x W x H):	31,5 x 20,3 x 8,9cm	39,1 x 20,3 x 8,9cm	46,7 x 20,3 x 8,9cm	42,9 x 20,3 x 8,9cm	53,1 x 20,3 x 8,9cm
Electrical 230V 50/60 Hz:	0,5 amps, 110 watts	0,9 amps, 210 watts	0,9 amps, 210 watts	1,3 amps, 310 watts	1,7 amps, 410 watts
Fuses:	5mm x 20mm, 5 amp quick acting				
Temperature low range:	ambient +5°C to 100°C				
Temperature high range:	75°C to 150°C				
Stability @ 37°C:	+/-1,5°C	+/-2°C	+/-2°C	+/-2,5°C	+/-2,5°C
Uniformity w/in the block @ 37°C:	+/-0.4°C				
Uniformity across similar blocks @ 37°C:	N/A	+/-0,5°C	+/-0,6°C	+/-0,7°C	+/-0,9°C
Stability @ 60°C:	+/-3°C	+/-4°C	+/-4°C	+/-5°C	+/-5°C
Uniformity w/in the block @ 60°C:	+/-0,6°C				
Uniformity across similar blocks @ 60°C:	N/A	+/-0,8°C	+/-1°C	+/-1,2°C	+/-1,4°C
Heat-up time to 100°C:	45 min.	50 min.	60 min.	70 min.	75 min.
Controls:	rocker switch heat indicator light low temperature knob, variable 1 to 10 dial markings high temperature knob, variable 1 to 10 dial markings				
Ship weight:	2,6kg	2,9kg	3,6kg	3,9kg	4,5kg


Analog Dry Block Heater Operating Instructions

These multi-purpose units are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting/boiling points and a wide variety of other laboratory procedures.

1. Switch the three (3) position rocker power switch from the center off position to the low range or high range position as desired. The controls are divided into two separate heating ranges, both thermostats having graduations to assist in setting the desired temperatures. The left hand 'low temperature adjust' controls from slightly above ambient to approximately 100°C. The right hand 'high temperature adjust' controls from approximately 75°C to 150°C.
2. The rocker switch has a center off position and is used to select the desired operating range. When operating at the point where the two thermostats overlap in temperature range, the proper thermostat must be chosen for the task being performed. Move the power switch to the desired operating range and turn the matching range temperature control knob clockwise to increase the temperature within the range selected. The heat indicator light will illuminate during the operation of the heater.
3. The temperature may be verified by placing a calibrated thermometer in the test solution or by insertion in the modular block thermometer well provided. This hole fits regular glass bulb thermometers or small diameter digital probes. Due to air currents and radiation losses, the temperature in the test solution will be lower than the temperature in the block itself. For the most accurate readings a thermometer should be placed in a sample test tube with solution matching the samples being tested. If the temperature is too high or too low, adjust clockwise to increase temperature, counter-clockwise to decrease temperature. Slight adjustments will usually suffice to correct the temperature setting. When the heat indicator light flashes on and off intermittently, check the temperature again. Allow sufficient time for the temperature to stabilize before re-adjusting. This procedure should be followed until the desired temperature is reached.

VWR Analog 1 Block Heater with modular block and tubes

Digital Dry Block Heater Control Panel

The front panel of the Dry Block Heater contains all the controls and displays needed to operate the unit.

A. Standby button/standby indicator light: The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to start the temperature and time functions. The standby indicator light will shut off. Press the standby button again and the unit will once again be in standby mode.

B. Caution hot top indicator light: Illuminates when the plate temperature is above 40°C.

C. Plate indicator light: Illuminates when the optional external RTD probe is not being used. The temperature displayed is the plate temperature.

D. Probe indicator light: Illuminates when the optional external RTD probe is plugged in. The temperature displayed is the probe temperature, NOT the plate temperature.

E. Actual indicator light: Illuminates when the temperature displayed is the actual temperature of the plate/RTD probe.

F. Set-point indicator light: Illuminates when the set-point temperature is displayed.

G. Temperature display: Displays the actual/set-point temperatures in conjunction with the actual/set-point indicator lights. **H.** Up/down arrows for set-point control. On/off button starts/stops the heating function.

I. Time display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9999. **J.** Up/down arrows for set-point control. On/off button starts/stops the timer function.

Digital Dry Block Heater Specifications

	1 Block	2 Block	3 Block	4 Block	6 Block	2 Block with Lid
Cat. No. with Euro plug:	460-3264	460-3267	460-3270	460-3273	460-3276	460-0183
Cat. No. with UK plug:	460-3265	460-3268	460-3271	460-3274	460-3277	460-0184
Cat. No. with Swiss plug:	460-3266	460-3269	460-3272	460-3275	460-3278	460-0185
Dimensions (L x W x H):	31,5 x 20,3 x 8,9cm	39,1 x 20,3 x 8,9cm	46,7 x 20,3 x 8,9cm	42,9 x 20,3 x 8,9cm	53,1 x 20,3 x 8,9cm	39,1 x 20,3 x 17,8cm
Electrical 230V 50/60 Hz:	0,5 amps, 110 watts	0,9 amps, 210 watts	0,9 amps, 210 watts	1,3 amps, 310 watts	1,7 amps, 410 watts	1,65 amps, 400 watts
Fuses:	5mm x 20mm, 5 amp quick acting					
Temperature Range:	ambient +5°C to 120°C					ambient +5°C to 100°C
Stability @ 37°C:	+/-0,2°C					
Uniformity w/in the block @ 37°C:	+/-0,2°C					
Uniformity across similar blocks @ 37°C:	N/A	+/-0,1°C	+/-0,15°C	+/-0,2°C	+/-0,3°C	+/-0,1°C
Stability @ 60°C:	+/-0,4°C					
Uniformity w/in the block @ 60°C:	+/-0,4°C					
Uniformity across similar blocks @ 60°C:	N/A	+/-0,5°C	+/-0,7°C	+/-0,8°C	+/-1°C	+/-0,5°C
Heat-up time:	45 min.	50 min.	60 min.	70 min.	75 min.	50 min.
Controls:	See page 6					
Ship Weight:	2,6kg	2,9kg	3,6kg	3,9kg	4,5kg	3,2kg

Digital Dry Block Heater Operating Instructions

Designed for applications that require repeatable results and superior temperature stability. These multi-purpose units are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting/boiling points, and a wide variety of other laboratory procedures.

For best performance, the Dry Block Heater should be used in a stable environment. The unit's environment should have no air currents, drafts or temperature changes and it cannot be placed in direct sunlight. The unit requires a stable electrical supply that is free of voltage fluctuations. An unstable environment will adversely affect the performance of the unit. For example, even minor air drafts or temperature changes will adversely affect the unit's ability to maintain a stable temperature.

1. **Getting ready:**
 a. When using the external RTD probe, plug the RTD probe into the three (3) pin DIN connector at the back of the unit and place the thermometer portion in the thermometer well of the modular block. When using multiple blocks, place the RTD probe in the front right modular block.
 b. Press the standby button to change the unit from standby mode. The temperature display, time display and probe indicator light will illuminate. When not using the RTD probe, the plate indicator light will illuminate. The temperature display will alternate between the actual and set-point temperatures.
2. **Setting temperature:**
 a. Press the up/down arrows to the right of the temperature display until you reach the desired temperature. Holding down either the up or down arrow will cause the set temperature to change rapidly, a single pressing of either key will move the set temperature by 0.1°C. When you release the button, the display will blink off and then on, indicating the new set temperature has been accepted. Once the set-point has been programmed and the keys are not being pressed, the temperature display and temperature indicator light will illuminate. The actual and set-point indicator lights will alternate between set and actual temperatures. There are three (3) audible beeps to indicate the set-point temperature has been reached.
 b. Set-point temperature adjustments can be made without interrupting heating using the up/down arrows to the right of the temperature display. After the change has been made and you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
 c. To stop heating, press the on/off button to the right of the temperature display.
 d. Allow time for the temperature to stabilize. The actual temperature displayed is the temperature at the bottom of the modular block or of the RTD probe. Once the displayed actual temperature agrees with the set temperature, several minutes should be allowed for the temperature to stabilize throughout the block evenly.

Overshoot protection: If the unit exceeds the set temperature by 10°C, the unit will automatically stop heating.

3. **Setting timed mode:** Programmed time.
 a. Press the up/down arrows to the right of the time display until you reach the desired time.
 b. Start this function by pressing the on/off button to the right of the time display, the unit will run for the selected time. When using the timer in conjunction with the heating function, when the time display reaches zero (0:00), four (4) audible beeps will indicate the time down function is complete. Both the time and heating functions will shut off automatically and the time display will default back to the set time. To repeat for the same time, simply depress the on/off button again.
 c. To interrupt an automatic timing cycle before it is completed, press the on/off button to the right of the time display. The time display will flash until you resume the time function by pressing the on/off button again. This interrupt will not stop the heating function, the heating function will stop only when the timer reaches zero (0:00).

4. **Setting time to zero (0:00) and continuous mode:** Accumulated time.
 a. Press and hold the on/off button to the right of the time display. After three (3) seconds, the display will indicate the previous set time.
 b. Simultaneously press both the up and the down arrows, the display will indicate zero (0:00). The unit time is now set to zero (0:00) minutes. Alternately, you can use the up/down arrows to get to zero (0:00).
 c. Press the on/off button to the right of the time display, the display will indicate actual running time. The up/down arrows will become inactive. To stop timer, press the on/off button again. **IMPORTANT:** This will not affect the heating function. Press the on/off button to the right of the heat display to interrupt the heating function.
 d. To reset, press and hold the on/off button to the right of the time display. After three (3) seconds the display will indicate the previous set time, which was zero (0:00).
5. **Turning unit off:**
 a. To turn the unit off, press the standby button, the temperature and time displays will be blank, the standby indicator light will turn on.

SINGLE POINT CALIBRATION

This procedure is used to fine tune and calibrate the block temperature at a specific temperature setting. It will only be active without an external temperature probe connected. This process may be repeated for up to three (3) separate set-points. If a fourth calibration set-point is entered, the first set-point entered will be overwritten.

1. Turn unit on.
2. Set desired temperature.
3. Stabilize twenty (20) minutes or more, measuring the block temperature with a calibrated precision instrument or thermometer.
4. Press and hold standby button, then press the temperature up button once. The unit will beep two (2) times, confirming calibration mode. The display will now be flashing.
5. Press the temperature up/down arrows until the display matches the temperature probe/thermometer.
6. Press standby button to exit calibration mode and return to normal heating.

This process may be repeated at the same set-point multiple times for fine tuning if desired.

The unit will now use the biased offset for that specific temperature setting and increase or decrease temperature accordingly to bring the block temperature to set temperature. The decimal point of the display will flash to indicate a biased offset is being used. All other temperature settings will use the standard internal calibration. This offset will be stored in memory and retained until reset.

To restore unit to factory setting:

Press and hold standby button while pressing the temperature down button once. The reset will be confirmed with two (2) beeps. Press the standby button to exit calibration mode and return to normal heating.

BEEPER PREFERENCE

To silence beeper operation (except for error codes), with the unit in standby mode, press and hold the time on/off button and press the standby button. To restore normal beeper operation, remove AC power to unit for 10 seconds and then restore. Alternately, you may have to turn the unit on and press and hold the standby button and press and hold the time on/off button simultaneously.

MODULAR HEATING BLOCKS AND SAMPLES

Only VWR Modular Heating Blocks are to be used in the Dry Block Heaters.

1. Select the proper blocks that will be a precise fit for the tubes, plates, or vials that you will be using for your application.
2. Only plastic or glass tubes, plates, or vials are to be used in the Modular Heating Blocks. Metal vessels will adversely affect the temperature performance of the unit. Metal vessels will dissipate too much heat into the air, thus adversely affecting the temperature readings of the unit.
3. Blocks need to be in all positions on the unit so the heating plate is not exposed to the environment.
4. To ensure proper heating, the tubes, plates, or vials that are used must be the correct size for the Modular Heating Block. The tubes, plates or vials must fit securely in the hole with no air gaps and maintain as much contact as possible with the wall of the block. This will ensure good thermal contact between the Modular Heating Block and the tube, plate or vial you are using for your application.
5. For proper heating, the fluid level within your tube, plate or vial should not exceed the height of the Modular Heating Block. If your application requires the fluid level to be above the height of the block, then a temperature cover is recommended.
6. When using a temperature measuring device in a sample, the end of the probe should be placed at the bottom of the sample and the height of the liquid should not exceed the height of the Modular Heating Block. Ensure that your temperature measuring device is designed for immersion in liquids.
7. Contact your VWR representative for information on ordering Modular Heating Blocks and other Dry Block Heater accessories.

OPTIONAL EXTERNAL RTD PROBE KIT (FOR DIGITAL UNITS)

1. If your application requires a high level of accuracy, the optional external RTD probe should be used with the Dry Block Heater.
2. Follow the 'Operating Instructions' for the proper installation of the optional external RTD probe. With the RTD probe plugged into the back of the unit, place the thermometer portion in the thermometer well of the Modular Heating Block. With the external RTD probe in place, the RTD probe is now driving the displayed temperature setting for the operation of the unit, not the heating plate of the unit. Once the optional external RTD probe is properly installed, the probe indicator LED above the temperature display will be illuminated.

TEMPERATURE TESTING

1. To ensure good thermal conductivity, select the proper Modular Heating Block for your application. Select the proper tubes or vials which are the proper size for the Modular Heating Block. Close contact, with no air gaps, must be maintained between the walls of the Modular Heating Block well and the sides of the tubes or vials.
2. Fill the tubes or vials so that the fluid level does not rise above the top surface of the Modular Heating Block.
3. Select a temperature measuring device that is designed for immersion in liquids. Place your calibrated temperature measuring device in one of the samples so it reaches the bottom of the tube or vial. Once the temperature measuring device is placed in the fluid sample, ensure the fluid level is still below the top surface of the Modular Heating Block.
4. Set the desired temperature on the unit, allow the unit to reach this temperature, and let the unit stabilize for an additional twenty (20) minutes or more before taking any temperature readings.
5. The temperature can also be tested by utilizing the thermometer hole in the Modular Heating Block. A calibrated temperature measuring device can be inserted into the thermometer hole where there is a snug fit and close contact between the walls of the block and the temperature device. The temperature device must reach the bottom of this hole with no air gaps. Follow the procedure above for allowing the unit to reach temperature and stabilize before taking any temperature readings.

6. If the measured temperature on your temperature measuring device does not match the actual temperature on the display of the unit (for Digital units only), then the single point calibration procedure can be used. By doing this, the unit will now be more accurate at that set-point for your specific application.

STABILITY TESTING

1. The manufacturer has performed temperature stability tests on Dry Block Heaters. The manufacturer used calibrated temperature measuring devices for the stability tests. Units were set-up with the proper number of Modular Heating Blocks so the heater plates were not exposed to the environment. A calibrated temperature measuring device was inserted into the thermometer well of one block on each unit. The testing temperature was set, the units were allowed to heat up and stabilize for a minimum of twenty (20) minutes and then temperature readings were recorded at regular intervals for four (4) hours. These tests confirmed the temperature stability of the units.
2. The recommended procedure for testing the unit's stability is as follows:
 a. Set-up the unit in a stable environment.
 b. Set-up the unit with the proper number of Modular Heating Blocks for the unit. Place a calibrated temperature measuring device into the thermometer well of the Modular Heating Block. The temperature device should have a snug fit with no air gaps in the thermometer well. Set the temperature of the unit. Allow the unit to reach temperature and stabilize for twenty (20) minutes or more and then take temperature readings for stability.
 c. Another method to test for stability is to set up the unit with the proper number of Modular Heating Blocks and then place the proper fitting tubes or vials into the blocks. Fill the tubes or vials with liquid, where the liquid level is lower than the top surface of the Modular Heating Block. Use a calibrated temperature measuring device that is designed for immersion in liquids, and place the temperature probe in the bottom of one of the tubes or vials in one block. Set the temperature of the unit. Allow the unit to reach temperature and stabilize for twenty (20) minutes or more and then take temperature readings for stability. Please note the characteristics of the liquid used for the test and the fit of the tube or vial in the Modular Heating Block can affect the results of a stability test for the Dry Block Heaters. Tubes and vials must be the proper size for the Modular Heating Block so there is a secure fit with no air gaps to ensure good thermal contact.

UNIFORMITY TESTING

1. The manufacturer has performed uniformity tests on Dry Block Heaters to ensure even heating across the entire heater plate. For the one block heaters, the manufacturer used a specially designed Modular Heating Block with five (5) temperature holes. Five independent calibrated temperature measuring devices were used during the test. For multiple block units, Modular Heating Blocks were placed in all positions to cover the entire heating plate and then an independent calibrated temperature measuring device was used in the thermometer well of each block. For all tests, the tip of the temperature measuring device was an exact fit for the temperature well in the Modular Heating Block with no air gaps. The temperature was set and the units were allowed to heat up and stabilize for twenty (20) minutes or more. Readings were then taken at regular intervals to monitor the temperature uniformity for four (4) hours. These tests were repeated after rotating the temperature measuring devices to different thermometer well positions to confirm the temperature uniformity of the units.
2. The recommended procedure for testing the unit's temperature uniformity is as follows:
 a. Set-up the unit in a stable environment.
 b. Set-up the unit with the proper number of Modular Heating Blocks to cover the entire heating plate. Then set up the blocks with the proper fitting tubes or vials. Place liquid samples in tubes or vials where the liquid level is lower than the height of the top surface of the Modular Heating Block.

DIGITAL DRY BLOCK HEATER OPERATING INSTRUCTIONS (CONT'D)

c. Select multiple calibrated temperature measuring devices that are designed for immersion in liquids. Use these temperature measuring devices at the same time across various positions in the Modular Heating Block(s). Place the temperature probes in the bottom of several of the filled tubes or vials. Ensure the liquid levels in the tubes or vials with the temperature probes does not rise above the height of the top surface of the Modular Heating Block.

d. Set the temperature of the unit. Allow the unit to reach temperature and stabilize for twenty (20) minutes or more and then take temperature readings from all the temperature devices to test for temperature uniformity.

e. Please note the characteristics of the liquid used for the test and the fit of the tube or vial in the Modular Heating Block can affect the results of the uniformity test for the Dry Block Heaters. Tubes and vials must be the proper size for the Modular Heating Block so there is a secure fit with no air gaps to ensure good thermal contact.

TECHNICAL SERVICE

For more information or technical assistance contact your local VWR representative or visit www.vwr.com.

Web resources: Visit VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

TROUBLESHOOTING

To clear error press the standby button

Problem	Cause	Solution
Unit fails to power on	Missing or blown fuse	Add or replace fuse as necessary. If problem persists, please contact your VWR representative for repair.
E1	Faulty temperature sensor	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E2	Thermocouple failure or Heating element failure	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E3	Unit cannot reach set-point or Probe not in thermometer well	If using probe, verify probe is in thermo well and follow Single Point Calibration instructions on page 9 If problem persists, please contact your VWR representative for repair.

ACCESSORIES

Microcentrifuge Tube Blocks
Single block.

Description	Cat. No.
0,5 mL Tube	460-3211
1,5mL Tube	460-3235
1,5mL Eppendorf™ Tube	460-3212
2mL Eppendorf™ Tube	460-3245
2mL Corning™ Tube	460-3246

Titer Plate Block
Double block. Fits 2/4/6 block Dry Block Heaters. Ideal for 96-well or 384-well titer plates. Recessed well for better stability, flat surface good for flat and round bottom plates.

Description	Cat. No.
Titer Plate	460-3238

Conical Bottom Centrifuge Tube Blocks
Single block.

Description	Cat. No.
15mL Tube	460-3221
50mL Tube	460-3223

Standard Test Tube Blocks
Single block.

Description	Cat. No.
6mm Tube	460-3213
10mm Tube	460-3214
12/13mm Tube (20 wells)	460-3216
12/13mm Tube (16 wells)	460-3215
15/16mm Tube	460-3217
17/18mm Tube	460-3243
20mm Tube	460-3218
25mm Tube	460-3219
35mm Tube	460-3222

Combination Blocks
Single block. These blocks have been designed for variable sized samples.

Description	Cat. No.
Test Tube Combination: 6mm, 12/13mm, 25mm	460-3220
Centrifuge Tube Combination: 1,5mL, 15mL, 50mL	460-3247
Micro-Tube Combination 0,5mL, 1,5mL, 2mL	460-3248

ACCESSORIES (CONT'D)

Vial Blocks

Single block. Designed for sample/serum and scintillation vials.

Description	Cat. No.
12mm Vial	460-3279
15mm Vial	460-3280
16mm Vial	460-3287
17mm Vial	460-3281
19mm Vial	460-3282
21mm Vial	460-3283
23mm Vial	460-3284
25mm Vial	460-3285
28mm Vial	460-3286

PCR-Plate, Tube, Strip Blocks

Single or double block. Tapered tube wells for 0,2mL tubes. Spaced for easy access and removal.

Description	Cat. No.
96-well PCR plate	460-3226
10 x 8 Tube Strips	460-3224
Individual Tubes	460-3225

Cuvette Block

Single block. Two parallel slots fit 6 cuvettes in each slot, side-by-side.

Description	Cat. No.
(12) 12,5mm Cuvettes	460-3237

Solid Blocks

Single and double block. For use as a low temperature hotplate, slide drying, or for custom drilling to make a custom block.

Description	Cat. No.
Single	460-3236
Double	460-3239

Low Temperature Cover

Plexiglass cover reduces air flow for additional temperature stability in low-temperature applications. Two sides are open 6,4mm.

Description	Cat. No.
2 Block Cover	460-3240
4 Block Cover	460-3241
6 Block Cover	460-3242

ACCESSORIES (CONT'D)

Stainless Steel Sand Baths
Ideal for irregular shaped vessels. Stainless steel construction for superior corrosion resistance. Designed to hold sand, stainless steel shot, or non-volatile liquids.

For Unit	Cat. No.
1 Block Dry Block Heater	460-3227
2 Block Dry Block Heater	460-3228
3 Block Dry Block Heater	460-3229
4 Block Dry Block Heater	460-3230
6 Block Dry Block Heater	460-3231
Accessories	**Cat. No.**
Sand, 0,45kg	460-3232
Stainless Steel Shot, 0,45kg	460-3233

External Temperature Probe Kit
Enables the unit to read actual block or sample temperature and display that temperature on the control panel. The optional External Temperature Probe Kit includes a stainless steel RTD probe, 45,7cm stainless steel support rod, thermometer / temperature probe extension clamp and hook connector. The RTD probe is designed to fit perfectly into the thermometer well of each modular block.

Description	Cat. No.
Optional External Temperature Probe Kit	460-3288

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Title of Council Directive on harmonization of the laws of the member states
Council Directive 73/23/EEC of 19 February 1973 on the harmonisation of the laws of Member States relating to Electrical Equipment designed for use within certain voltage limits
Council Directive 89/336/EEC of 3 May 1989 on the approximation of the laws of the Member States relating to Electromagnetic Compatibility

and are in conformity with the following standard(s) and/or other normative document(s):
et sont conformes à la (aux) norme(s) ou autre(s) document(s) normatif(s) :
(den) folgenden Normen oder anderer normativer Dokumente entsprechen:

Reference number	Title of the applied standards:
IEC 61010-1	**Safety requirements for electrical equipmentfor measurement, control and laboratory use – Part 1: General requirements – First Edition**
IEC 61010-2-010	**Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-10: Particular requirements for laboratory equipment for heating of materials**
EN 45501; 55022A EN 61000-3-2;-3-3 EN 61000-4-2; -4-3; -4-4 EN 61000-4-5; -4-6; -4-11	**Electromagnetic compatibility (EMC)**

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Yves Van Damme
European Regulatory Affairs Manager

Frankie Vangeel
Managing Director

VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium
Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

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EN61010-1:2001	Relating to safety requirements for electrical equipment for measurement, control and laboratory use – part 1: general requirements – second edition, 2001, including amendment 1 (1992) and amendment 2 (1995)
EN-61010-2-010:2005	Relating to safety requirements for electrical equipment for measurement, control and laboratory use – part 2-010: particular requirements for laboratory equipment for heating of materials

VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium
Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

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Ireland / Northern Ireland
VWR International Ltd / VWR International
(Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy
VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com
info@internationalpbi.it

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway
VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax. 058 32 38 205
E-mail: labart@pl.vwr.com

Portugal
VWR International - Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain
VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden
VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland
VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

UK
VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

China
VWR (Shanghai) Co., Ltd
2nd Floor, Building 4,
Lane 998, Halei Rd,
Zhangjiang Hi-tech Park
Shanghai, 201203
China
Tel.: +86-21 589 868 88
Fax: +86-21 585 588 01
E-mail: info_china@vwr.com

India
VWR Lab Products Private Limited
135/12, Brigade Towers, 2nd Floor
Front wing, Brigade Road,
Bengaluru , India – 560 025
Tel.: +91-80-41117125/26 (Bengaluru)
Tel.: +91-2522-647911/922 (Mumbai)
Fax: +91-80-41117120
E-mail: vwr_india@vwr.com

Singapore
VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Legal Manufacturer:
VWR International bvba • Researchpark Haasrode 2020 • Geldenaaksebaan 464 • B-3001 Leuven • + 32 16 385011 • http://be.vwr.com
549010-00 (Rev 4)