VWR

Microscope VisiScope® series 384

INSTRUCTION MANUAL

Model	European Catalogue Number
TL384 POL	630-1936

Version: 1
Issued: 11, 03, 2015

CE

Legal Address of Manufacturer

Europe
VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
http://be.vwr.com

Country of origin: **ITALY**

Table of Contents

Warning
Safety information
Package contents
Unpacking
Intended use
Symbols and conventions
Product specifications
Overview
Using the microscope
Troubleshooting
Repair and maintenance
User replaceable accessories and spare parts
Technical service
Warranty
Compliance with local laws and regulations
Disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.
We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.
We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.
Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Package Contents

DESCRIPTION	QUANTITY
Microscope stand with head, focus and rotating stage	1
WF10x/20 eyepiece	2
Bertrand lens attachment	1
Objective E-PLAN IOS POL 4x	1
Objective E-PLAN IOS POL 10x	1
Objective E-PLAN IOS POL 40x	1
Objective E-PLAN IOS POL 60x	1
Condenser with rotating polarizer	1
Retardation plates (lambda, lambda/4, quartz wedge)	3
Allen wrench for objective centering	2
Power supply 6Vdc	1
Dust cover	1

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take care to ensure that the optical items (objectives and eyepieces) do not fall out and get damaged. Using both hands (one around the neck and one around the base), lift the microscope from the container and put it on a stable desk.

1. First of all remove all the components from the original package and store the box in a safe and dry place.
2. Place the main body of the microscope on a flat and sturdy surface.


Don’t connect the power supply cable at this stage.

3. Place the attachment with analyzer /Bertrand Lens onto the main body and turn the head locking screw (1).

4. Insert the optical head and lock the screw.

5. After removing the protection caps, fully insert both eyepieces.

6. Connect the external power supply cable to the socket on the back of the microscope.

 If there is no need to install additional accessories, the instrument is now ready for use.

Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Symbols and conventions


The following chart is an illustrated glossary of the symbols that are used in this manual.

⚠	**CAUTION** This symbol indicates a potential risk and alerts you to proceed with caution


Product specifications

Head	Trinocular, 30° inclined, 360° rotating.
Eyepieces	WF10X/20mm with cross hair.
Bertrand lens	Swing-out type; centrable.
Polarizing attachment	0°-90° rotating analyzer filter. Tint plates included: 1° order red (λ), λ/4, quartz wedge.
Nosepiece	4-positions with centering mechanism for all objectives.
Objectives	E-PLAN IOS POL (Strain-free) 4x/0.10, 10x/0.25, 40x/0.65, 60x/0.80.
Focusing system	Coaxial coarse and fine.
Stage	160mm dia.; 360° rotating with stop knob and 0.1° Vernier.
Condenser	N.A. 1.25, with iris diaphragm, focusable and centrable. With rotating polarizer filter.
Illumination	P-LED^3, with manual brightness control.
Power supply	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A

Overview

Overview

PHOTO/VIDEO PORT
EYEPIECES (ONE WITH CROSSLINE)
DIOPTRIC ADJUSTMENT RING
BERTRAND LENS (13)
ANALYZER ROTATION SCALE (8)
ANALYZER INSERTION LEVER (12)
RETARDATION PLATE (11)
HOLES FOR NOSEPIECE CENTERING (7)
OBJECTIVE
ROTATING STAGE
STAGE ROTATION LOCK SCREW
CONDENSER CENTERING SCREWS (2)
POLARIZER

Using the microscope

Observation head
Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.


Place the specimen on the stage
Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

Illumination system settings
The microscope is fitted with a white LED illuminator. Turn the brightness adjustment knob to a brightness suitable for observation.

Adjust interpupillary distance
Hold the right and left parts of the observation head with both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.
The white dot (°) placed on the right eyepiece shows the set interpupillary distance. Just remember this value to help on later settings.

Focus tension adjustment
The tension of the coarse focusing knob is preset by factory.
To change the tension according to your preference, just rotate the knob (3) clockwise in order to increase it.
Excessive tension could damage the mechanism of focus. A too loosed tension causes the descent of the stage by gravity or a sudden loss of focus. In this case, rotate the knob (3) to increase the tension.

Diopter adjustment
Turn the dioptric adjustment ring on the left eyepiece halfway. Turn the coarse focus knob in order to focus the slide with an objective with low magnification. Adjust the fine focus knob until you obtain a clear and defined picture observing with the right eye, and then repeat the operation with the left dioptric compensation ring and the left eye. When the image appears in focus, choose the necessary objective with the revolving nosepiece.

Condenser
The condenser, placed below the stage, is equipped with the following features:

Raise or lower the condenser through the knob (5) to obtain a clear and uniform illumination of the sample.
To center the condenser: completely close the iris diaphragm (4). Using the condenser centering screws (2), move the diaphragm in the center of the field of view. Then gradually expand the diaphragm until it is tangent to the edges of the field of view. If necessary, you can perform an additional adjustment.
The condenser is centered when the edges of iris diaphragm are tangent to the field of view.

In the bottom part of the condenser there is a polarizer filter, which can be rotated through 360°, with a reference notch for a precise angle measurement.

Numerical aperture setting
The value of the numerical aperture (N.A.) of the diaphragm is an indication of the contrast of the illumination system. Matching the value of illumination system's N.A. with that of the objective ensures the best results in terms of contrast and image quality.
To set the numerical aperture of the illuminator, adjust the opening of the iris diaphragm (4). In this way you control contrast and image resolution.
For samples with low contrast set the iris to about 75% of the value of the objective's numerical aperture.


System centering
In order to correctly center the whole optical system, first remove the Bertrand lens out of the path by rotating its disc. Remove also the lower polarizer filter.

The rotating stage is factory pre-centered and so it doesn't need any additional operation.

It's necessary to center the objectives following this procedure:

1. Place a specimen on the stage.
2. Insert 10x objective.
3. Insert the two little wrenches into the holes for nosepiece centering (7).
4. Focus your sample while continuously rotating the stage clockwise and then counter-clockwise by a little angle (e.g. 30° or 45°). During these oscillations you should spot a point on the sample that doesn't rotate on a circumference but only revolving about itself (fixed point). With the two objective's centering wrenches, bring this point in the center of the crossline.
5. Insert the 20x objective.
6. Insert the two little wrenches into the holes for nosepiece centering (7), above the 20x objective.
7. Rotate the two wrenches until the same point of the sample is brought in the center of the crossline.
8. Repeat the operation for the other objectives.

Polarized light observation
In order to correctly use the polarized light, it's necessary to perform a setting for "light extinction":
1. Remove from the optical path the sample, retardation plates, etc...
2. Insert the 10x objective.
3. Insert in the optical path the analyzer filter and rotate the rotation scale (8) in correspondence to position "0".
4. Put the angular scale (9) of the polarizer filter to "0" position.
5. Loosen the lock screw of the polarizer and, looking into the eyepieces, rotate the scale of the polarizer until you obtain the complete extinction of the light (totally dark).
6. Then lock the polarizer screw.

Bertrand lens
Bertrand lens is employed in conoscopic observation mode (lens inserted) or orthoscopic mode (lens not inserted).

Orthoscopic observation
The path of the polarized light is parallel to the optical axis, thus allowing the observation of the optical characteristics of the specimen. You can use objectives from 4x to 100x.
1. Remove the Bertrand lens out from optical path.
2. Insert the analyzer filter using the proper lever (12) and look at the sample.
3. Insert a retardation plate (11) pushing it completely inside the slot under the head.

To remove the plate from the optical path, pull it outward till the first click.

Conoscopic observation
1. Use objectives from 20X to 100X.
2. Insert the analyzer using the proper lever (12) and put it in the extinction position.
3. Insert the Bertrand lens.
4. Insert the objective of choice (from 20X to 100X).
5. Open the aperture diaphragm.
6. It's possible to center the conoscopic image adjusting the screws (10) placed under of the head.


IMAGE CONOSCOPIC

In this manner it's possible to look at the conoscopic images that form on the back focal plane of the objective. If the peripheral part of the image is too dark, you can just move the condenser vertically to find the optimal position for the illlumination.

Retardation plates
Retardation plates (1 lambda, corresponding to 530nm or Red 1° order, ¼ lambda corresponding to 137nm, and quartz wedge) are used to make the polarization analysis more precise and refined.
Insert or remove the plate in the optical path of the microscope by pushing it into the slot (11) under the head.
Retardation plates have no visible effect when used in normal brightfield light.

Repair and maintenance

Microscopy environment
This microscope is recommended for use in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

After using the microscope

The microscope should always be kept vertical when moving it so that no moving parts, such as the eyepieces, fall out.
Never mishandle or impose unnecessary force on the microscope. Never attempt to service the microscope yourself. Turn down the illumination intensity control and turn the light off. Cover the microscope with the included dust cover, and keep it in a dry and clean place.

Electrical safety precautions

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics
If the optical parts require cleaning first use compressed air.
If that is not sufficient use a soft lint-free cloth with water and a mild detergent.
And as a final option use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics. Do not disassemble objectives or eyepieces in an attempt to clean them.

For the best results, use the VWR cleaning kit (see catalogue number below).

If you need to send the microscope to manufacturer for maintenance, please use the packaging if possible.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
LIGHT DOESN'T TURN ON	Power supply not connected	Check that the 6Vdc power supply jack is well inserted on the rear of the microscope
		Rotate the brightness adjustment control and check if there is an increase in the light output
IMAGE CANNOT BE SEEN OR IS DARK.	The iris diaphragm aperture is not completely opened	Completely open the iris diaphragm aperture
	Brightness level is low	Rotate the brightness adjustment potentiometer
	Objective is not aligned with the optical axis	Rotate the nosepiece until an objective is well inserted in the optical path (it "clicks")
IMAGE IS UNCLEAR, BLURRED OR HAS INSUFFICIENT CONTRAST.	Objectives or filters are dirty	Wipe them clean
	The iris diaphragm aperture is not opened correctly	Open the iris diaphragm aperture completely
	Condenser at wrong height	Rotate the condenser knob until you see a uniform illumination

User replaceable accessories and spare parts

DESCRIPTION	QUANTITY	CAT. NO.
Eyepiece WF20x	2	630-2000
Eyepiece WF15x	2	630-2001
Eyepiece WF20x	2	630-2002
Micrometer eyepiece WF10x/20mm	1	630-2003
Slide micrometric 26x76mm 1mm	1	630-1650
Objective E PLAN IOS 4X for serie 200&300	1	630-1639
Objective E PLAN IOS 10X for serie 200&300	1	630-1640
Objective E PLAN IOS 20X for serie 200&300	1	630-1641
Objective E PLAN IOS 40X for serie 200&300	1	630-1642
Objective E PLAN IOS 100X for serie 200&300	1	630-1643
Objective 60x/0,80 IOS E-PLAN for serie 300	1	630-2049
Darkfield condenser for dry objectives for serie 300	1	630-2060
Phototube adapter for APS-C sensor	1	630-1645
Camera adapter CCD 0,45X	1	630-1644
Battery solar pack	1	630-1637

Technical service

Web Resources

Visit the VWR's website at **www.vwr.com** for:
• Complete technical service contact information
• Access to VWR's Online Catalogue, and information about accessories and related products
• Additional product information and special offers
Contact us For information or technical assistance contact your local VWR representative or visit. **www.vwr.com.**

Warranty


VWR International warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.
Instead it is your responsibility to correctly dispose of your equipment the end of its life cycle by handling it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.
For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland
VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy
VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

New Zealand
Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Norway
VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax. 058 32 38 205
E-mail: info@pl.vwr.com

Portugal
VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore
VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain
VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden
VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland
VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey
Pro Lab Laboratuar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu Işmerkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK
VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

GO TO VWR.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS OF YOUR LOCAL VWR DISTRIBUTOR

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