Instruction manual

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Operating Instructions

Incubating Mini Shaker Set-Up 8

Incubating Microplate Shaker Set-Up 8

Incubating Microplate Shaker Tray Configurations 8

Getting Ready 10

Idle Mode 10

Shake Mode 10

Heating Mode 11

Turning Unit Off 11

Setting Menu 11

System Settings 12

RS232 13

Troubleshooting 14

Accessories 15

Technical Service 15

Warranty 15

Compliance 16

Disposal 16

Introduction

This manual contains installation, operation, and maintenance instructions for the VWR Incubating Mini and Microplate Series. Please read the manual completely before using.

Safety Information

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions, and false results.

WARNING	For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
ATTENTION	For important information about the product. May lead to equipment damage if not avoided.
NOTE	For useful information about the product.

Warning Symbols

General hazard

Electrical shock hazard

CAUTION! HOT

Safety Precautions

Please read the entire instruction manual before operating the Incubating Mini/Microplate Shaker.

WARNING! DO NOT use the Incubating Mini/Microplate Shaker in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

CAUTION! The caution hot indicator light warns that the temperature of the chamber air is above 40°C. The light will illuminate and remain lit when the temperature of the tray reaches approximately 40°C. When the heat is turned off, the caution hot indicator light will stay lit until the temperature of the chamber air is less than 40°C

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the tray or lid.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly. **DO NOT** immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground – Protective Conductor Terminal. Protective earthing of the equipment is achieved via connection of the provided power cord to a compatible grounded power outlet.

Alternating Current

Intended Use

This instrument is intended for use in laboratories. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from VWR, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

Package Contents

Incubating Mini Shaker or Incubating Microplate Shaker

Description	Quantity
Shaker	1
Power Cords 230V: EU cord, UK cord, and a Swiss Plug Adapter	3
Microplate Clips *(Microplate Shakers ONLY)*	4

Installation

Upon receiving the VWR Incubating Mini/Microplate Shaker, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the Incubating Mini/Microplate Shaker on a level bench or table, away from explosive vapors. Secure to an immovable work surface by pressing down on the four (4) corners of the unit, creating a strong suction to the work surface (DO NOT place on a bench mat). Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit. Always place the unit on a sturdy work surface.

The Incubating Mini/Microplate Shaker is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The 120V unit plugs into a 120-volt, 50/60 Hz source. The 230V unit plugs into a 230-volt, 50/60 Hz source.

If the cord supplied does not meet your needs, please use an approved power cord that has ratings equal to or exceeding those of the originally provided cord and that complies with the local/national regulations of the country in which the equipment is to be used.

Maintenance & Servicing

The Incubating Mini/Microplate Shaker is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. However at least every three (3) months you should:

- Unplug the unit.
- Remove any accumulated dirt from the base and tray.
- Check all accessible items to make sure they are properly tightened.

The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. **DO NOT** use a cleaning agent or solvent on the front panel or lid which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

Environmental Conditions

Operating Conditions: Indoor use only.

Temperature: 5 to 40°C

Humidity: maximum 80% relative humidity, non-condensing

Altitude: 2000 M above sea level

Non-Operating Storage:

Temperature: -20 to 65°C

Humidity: maximum 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.


Product Specifications

Incubating Mini Shaker

Overall dimensions (L x W x H):	43.2 x 26.2 x 29 cm
Interior Dimensions (L x W x H):	27.9 x 19.7 x 15.3 cm
Tray Dimensions (L x W):	28 x 19.7 cm
Height with Open Lid:	41.9 cm
Tray Construction:	Aluminum Tray with 77 mounting points for accessories
Motion:	Orbital
Orbit:	3 mm
Display:	115 x 30 mm, White Character LCD
Electrical:	230 volts, 2.0 amps, 460 watts
Heating Power:	350 W at 230V
Fuses:	5mm x 20mm, 5-amp quick acting, 250V
Temperature Range:	Ambient +5°C - 70°C
Temperature Set Resolution:	0.1°C
Temperature Uniformity:	±0.5 °C at 37 °C
Temperature Calibration (SPC) Set Points:	3 temperature points
Hot Indicator:	"HOT" texts and Yellow Triangle
Over Temperature Protection:	Yes, shut off the heater but keep the motor running as it is.
Operating Temperature Range:	5 °C to 40 °C, 80% RH, non-condensing
Speed Range:	100-1200 rpm
Speed Steps:	1 rpm increments
Speed Accuracy:**	± 2% of set speed
Motor Overload Protection:	Yes, motor shuts off and unit remains heating
Weight Capacity:	4 kg
Volume Capacity:	35 x 10 mL flasks 20 x 25 mL flasks 12 x 50 mL flasks 8 x 125 ml flasks 5 x 250 ml flasks 2 x Microtube Racks 1 x Test tube Rack
Motor Type:	Brushless DC Motor
Driving System:	Triple Eccentric
Timer:	Up to 99 hours 59 minutes (increased in 1 second increments)
Data Logging & Remote Control:	built in RS232
Programming:	On-board programming for 8 separate, 5-step programs
Recall Feature:	Retrieve last 5 recently used settings
Set Point Retention:	Yes – displays last setting
Controls:	See on page 9
Ingress Protection:	IP21
Altitude:	2000m
Shipping Weight:	13.9 kg
Safety Certification:	CE Marked and TÜV SÜD listed

**Conditions permitting. Variations in measurement process, vessel(s), and sample(s) will impact actual performance.


Incubating Microplate Shaker

Overall Dimensions (L x W x H):	43.2 x 26.2 x 20.8 cm
Interior Dimensions (L x W x H):	28 x 19.7 x 7 cm
Tray Dimensions:	28 x 19.7 cm
Height with Open Lid:	41.9 cm
Tray Construction:	Aluminum tray with hardware to secure Microplates and Deep Well Plates
Motion:	Orbital
Orbital:	3 mm
Display:	115 x 30 mm White Character LCD
Electrical:	230 volts, 2.0 amps, 460 watts
Heating Power:	350W at 230V
Fuses:	5mm x 20mm, 5-amp quick acting, 250V
Temperature Range:	Ambient +5°C - 70°C
Temperature Resolution:	0.1°C
Temperature Uniformity:	±0.5 °C at 37 °C
Temperature Calibration (SPC) Set Points:	3 temperature points
Hot Indicator:	"HOT" texts and Yellow Triangle
Over Temperature Protection:	Yes, shut off the heater but keep the motor running as it is.
Operating Temperature Range:	5 °C to 40 °C, 80% RH, non-condensing
Speed Range:	100-1200 rpm
Speed Steps:	1 rpm increments
Speed Accuracy:**	± 2% of set speed
Motor Overload Protection:	Yes, motor shuts off and unit remains heating
Weight Capacity:	1kg
Volume Capacity:	4 x Microplates 2 x Microtube Racks 4 x Deep Well Plates
Motor Type:	Brushless DC Motor
Driving System:	Triple Eccentric
Timer:	Up to 99 hours 59 minutes (increased in 1 second increments)
Data Logging & Remote Control:	built in RS232
Programming:	On-board programming for 8 separate, 5-step programs
Recall Feature:	Retrieve last 5 recently used settings
Set Point Retention:	Yes, displays last setting
Controls:	See on page 9
Ingress Protection:	IP21
Altitude:	2000 m
Shipping Weight:	13.9 kg
Safety Certification:	CE Marked and TÜV SÜD listed

**Conditions permitting. Variations in measurement process, vessel(s), and sample(s) will impact actual performance.

Incubating Mini Shaker Set-Up

The Incubating Mini Shaker is supplied with a tray designed to hold a variety of accessories.

1. Flat containers can be shaken by placing them on the tray.


See page 14, for optional accessories and unit capacities.

Incubating Microplate Shaker Set-Up


The Incubating Microplate Shaker is designed to hold two (2) or four (4) microplates, or two (2) Micro-Tube Racks.

1. Place two (2) microplates or deep well blocks diagonally on the tray, or place four (4) microplates or deep well blocks on the tray. The plates/blocks do not have to be filled.
2. Place the corner of the plate/block under the spring located at each corner of the tray.
3. Slide plate/block into place. You are ready to use.


2.


3.


Incubating Microplate Shaker Tray Configurations

(2) microplates


(4) microplates

(2) Micro-Tube Racks

Note: The Tray on the Incubating Microplate Shaker is not designed to hold flask clamps.

Control Panel

The front panel of the Incubating Mini/Microplate Shaker contains all the controls and displays needed to operate the unit.

A. **Actual Temperature display:** Displays the actual temperature.

B. **Set-Point Temperature display:** Displays Set-Point temperature.

C. **Heat indicator light:** Illuminates when the heat function is heating.

D. **Timer:** Counts up when unit is running or counts down from set-time.

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F. **Set-Point Speed display:** Displays the set-point speed selected.

G. **Temperature Adjustment Knob:** for set-point control. Click knob for on/off button starts/stops the heating function.

H. **Speed Adjustment Knob:** for set-point control. Click knob for on/off button starts/stops the time function.

I. **Speed indicator light:** will be illuminated when the unit is shaking

J. **Caution hot indicator light:** Illuminates when the air temperature of the chamber is above 40°C (104°F).

Incubating Mini and Microplate Shaker Operating Instructions

The Incubating Mini/Microplate Shakers have been designed for the temperature, speed, and time functions to work independently of one another. The temperature and speed can be reset without resetting the timer and the timer can be stopped and started without interrupting the heating and shaking functions.

1. Getting Ready:

1. Connect the provided power supply to the unit.
2. Plug the female end of the provided power cord into the power supply and plug the male end of the power cord into a matching standard grounded outlet.
3. Toggle the power switch to the ON (•) position.
4. Allow unit to complete the start-up sequence:
 a. LCD Test – all display icons illuminate
 b. Model Type, software version
 c. Function Units Timer [Hours: Minutes], Shake Speed [RPM]
5. The main operating screen will return, and the previous timer and shaker settings will be displayed.

2. Idle Mode

Timer

1. By default, the timer will be set to 00:00 and count upwards when the shaking function is turned on.
 a. The timer begins in MM:SS mode
 b. Once the timer reaches 59:59 it will automatically switch to HH:MM mode at 01:00.
2. When the timer is set to any time between 00:01 and 99:59, the timer will count downwards from the set time once the shaking function is turned on and the ***Timer Status indicator Light (C)*** will be illuminated.
 a. The countdown timer can only be set in HH:MM mode.
 b. Once the timer reaches 01:00 in HH:MM mode, it will automatically switch to MM:SS mode at 59:59.
 c. Once the timer reaches 00:00 in the countdown mode shaking functions will shut off, the unit will beep 1 time and "End" will appear on the display (**C, Figure 1**).

Adjust Timer

1. Click the **Speed *Adjustment Knob*** (**B, Figure 1**).
 a. The minutes will blink on the display.
 b. Rotate the knob to adjust the minutes value
 c. Quick press the ***Speed Adjustment Knob*** (**B**) to select the minutes value
 d. The hours will now blink on the display.
 e. Rotate the ***Speed Adjustment Knob*** (**B**) to adjust the hours value
 f. Quick press the ***Speed Adjustment Knob*** (**B**) to confirm hours value and the timer setting

Reset Timer to Zero

1. Ensure unit is not running, click then rotate the **Speed *Adjustment Knob*** (**B**) so either minutes or hours are blinking
2. Press and hold ***Speed Adjustment Knob*** (**B**) until 00:00 is on timer display

Adjust Set Speed

1. Rotate the ***Speed Adjustment Knob*** (**B, Figure 1**) to the desired setting then hold

3. Shake Mode

- Start Shake Function
 1. Quick press the ***Speed Adjustment Knob*** (**B**) until the unit beeps and the Run Status Indicator (D, Figure 1) illuminates.
- Adjust Speed while shaking
 1. Rotate the ***Speed Adjustment Knob*** (**B**)
- Adjust Remaining Time (count down only when timer is set)
 1. Ensure the unit is not running then quick press and rotate the ***Speed Adjustment Knob*** (**B**) to adjust the minutes value. The minutes value will blink while making adjustment.
 2. Quick press the ***Speed Adjustment Knob*** (**B**) to confirm the minutes value, and now the hours value will blink.
 3. Rotate the ***Speed Adjustment Knob*** (**B**) to adjust the hours value
 4. Quick press the ***Speed Adjustment Knob*** (**B**) to confirm hours value and the new timer setting

 NOTE: When making adjustments and the minute or hour is blinking; if the time setting remains idle without confirmation for 6 seconds it will reset to the current setting.
- Stop Shake

1. Press and hold the ***Speed Adjustment Knob*** (**J**) until the unit beeps and the ***Run Status Indicator Light*** (**H**) turns off

4. Heating Mode:

- Setting Temperature
 1. Rotate the ***Temperature Adjustment Knob*** (**G**) to desired temperature then press and hold for unit to begin heating
 2. ***Set-Point Temperature*** (**B**) and ***Actual Temperature*** (**A**) will display, ***Heating indicator light*** (**C**) will illuminate indicating the unit is heating.
- Changing Set Temperature
 1. Using the ***Temperature Adjustment Knob*** (**G**) adjust the set-point to the desired temperature.
 2. Press and hold the ***Temperature Adjustment Knob*** (**G**) until ***Heating indicator light*** (**C**) blinks and ***Set-Point Temperature*** (**B**) displays new set temperature.

CAUTION HOT indicator:

The caution hot indicator light warns that the temperature of the air in the chamber is above 40°C (104°F). The light will illuminate and remain lit when the temperature of the air in the chamber reaches approximately 40°C (104°F). When the heat is turned off, the caution hot indicator light will stay lit until the temperature of the air in the chamber is less than 40°C (104°F).

5. Turning unit off:

a. To turn the unit off, toggle the powerswitch. The temperature, speed and time displays will be blank. The Incubating Mini/Microplate Shaker should be kept in standby mode when not in use. To completely cut off power to the unit, disconnect the power cord from the unit or unplug from the wall outlet.

Operating Tips

As a safety feature, a built-in program will shut power off to the motor if the tray is prevented from rotating, or the unit is overloaded beyond its recommended weight capacity.

Built-in memory maintains the last used temperature, speed, and time settings during a power interruption.

6. Setting Menu

Menu Access

- While in Idle Mode, press and hold the ***Temperature Adjustment Knob*** (**G**) until "MENU" is displayed

Navigation

- To exit the Menu at any time, toggle the power switch.


Figure 2: Mini/ Microplate Shaker Settings Menu Structure; default settings underlined.

Recall

- While in Idle Mode, press and hold the ***Temperature Adjustment Knob*** (**G**) until "MENU" is displayed
- After "MENU" is displayed "RCL" will appear on the display.
- Recall previous run conditions
 - Rotate the **Temperature Adjustment Knob (G)** to navigate through the previous 5 unique run conditions

- To select a recall run condition, press, and hold the ***Speed Adjustment Knob*** (**H**) until the unit beeps and the **Run Status Indicator** (**H**, Figure 1) illuminates.
- Unit will now run at the recall timer and speed setting that was selected.

User Programs

- The User Programs feature allows the user to create, store, delete, and run a series of shaking programs for repeat applications in the lab.
- Each program can have a maximum of 5 steps and a time must be entered for each step. To keep the shaking off during a step enter "0" for the speed.
- The unit will store up to 8 individual 5-step programs.
- To access the Program Feature:
 - While in Idle Mode, press and hold the ***Temperature Adjustment Knob*** (**G**) until "MENU" is displayed, after "MENU" is displayed "RECALL" will appear on the display.
 - Rotate the ***Temperature Adjustment Knob*** (**G**) until "PROG" appears on the display
 - Quick press either the ***Temperature Adjustment Knob*** (**G**) to access the Program menu and Quick Press to select an option
 - Run program
 - Navigate to the User Program submenu ("PROG")
 - Navigate to the Run option
 - Navigate existing User Programs; quick press the speed knob.
 - Press and hold the ***Temperature Adjustment Knob*** (**G**) until the unit beeps and the **Run Status Indicator** (**C**, Figure 1) illuminates
 - **NOTE**: While a program is running, the end user cannot adjust the speed or time settings.
 - To stop a program, press and hold the ***Speed Adjustment Knob*** (**H**) until the unit beeps and the ***Run Status Indicator*** (**C**) turns off or turn off the power switch
 - Create new User Program
 - Navigate to the User Program submenu ("PROG")
 - Select the Edit option
 - Select new program ("Edit PrGn")
 - Select new step ("PrG<X> StPn")
 - Adjust Hours parameter
 - Quick press to select Minutes parameter
 - Adjust Minutes parameter
 - Quick press to select RPM parameter
 - Adjust RPM parameter
 - Quick press to confirm step parameters
 - Rotate either knob to add a new step ("PrG<X> StPn") to the program
 - Repeat above steps to set Hours, Minutes, Speed
 - To edit an existing step while creating a new program, rotate either knob and select an existing step to further adjust following the above steps
 - Select Back to complete program creation and save the program
 - Edit saved program
 - Navigate to the User Program submenu ("PROG")
 - Select the Edit option
 - Select program to be edited ("Edit PrG<X>")
 - Select new step ("PrG<X> StPn") to add a step to a program or select the existing step to be edited ("PrG<X>" "StP<Y>")
 - Adjust step parameters
 - Select Back to complete program edit and save the program
 - Clear program step
 - Navigate to the User Program submenu ("PROG")
 - Select the Clear ("CLr") option
 - Select the appropriate program ("Edit PrG<X>")
 - Select the step to be cleared ("PrG<X>" "StP<Y>")
 - Select the "yes" option to verify clear request
 - Clear program
 - Navigate to the User Program submenu ("PROG")
 - Select the Clear ("CLr") option
 - Select the appropriate program ("Edit PrG<X>")
 - Select the ALL option ("PrG<X>" "ALL")
 - Select the "yes" option to verify clear request

Single Point Calibration Procedure

This procedure is used to fine tune and calibrate the Incubating Mini-Microplate Shaker at a specific temperature setting. This process may be repeated for up to three (3) separate set-points. If a fourth calibration set-point is entered, the first set-point entered will be overwritten.

1. Turn unit on.
2. Press and hold the ***Temperature Adjustment Knob*** (**G**) to open the MENU. Navigate through to "CAL" then select "SPC". You will see the SPC Icon illuminate.
3. Select "New", then set desired temperature. Long press the ***Temperature Adjustment Knob*** (**G**) to start temperature regulation, the SPC Icon will blink. Allow the unit to stabilize.
4. Once the unit is stabilized, ***Set-Point Temperature Display*** (**B**) will flash. Using the ***Temperature Adjustment Knob*** (**G**) to enter the external temperature to match the temperature probe/thermometer, then quick press.
5. The unit will begin to calibrate to the set-point temperature.
6. Once the unit has reached the set-point temperature the display will blink "SAVE". If you are done with the calibration quick press save.
7. To *adjust* the SPC point turn the ***Temperature Adjustment Knob*** (**G**) clockwise to select "ADJ" and quick press. Enter the externally measured temperature and allow unit to stabilize.
8. To *clear* the SPC point turn the ***Temperature Adjustment Knob*** (**G**) clockwise to select "CLR" and quick press. Select Yes/ No. This process may be repeated at the same set-point, multiple times for fine tuning if desired.

The unit will now use the biased offset for that specific temperature setting and increase or decrease temperature accordingly to bring the temperature to set temperature. All other temperature settings will use the standard internal calibration. This offset will be stored in memory and retained until reset.

System Settings

Beeper Enable/Disable

- Disabling the Beeper will silence and prevent beeps in the following scenarios:
 - Starting and stopping shaking
 - When the timer reaches 00:00 in timer countdown mode
- Navigate to the System Settings submenu ("SYS")
- Navigate to the beep setting ("BEEP")
- Quick press to select value (En = Enabled or dIS = Disabled)
- Rotate ***Temperature Adjustment Knob*** (**G**) knob to select the desired setting
- Quick press to confirm value

Power Recovery Enable / Disable

- Power Recovery is an optional feature that allows the unit to automatically restart the shaking functions when power is returned to the unit after a disconnect.
- Navigate to the System Settings submenu ("SYS")
- Navigate to the Power Recovery setting ("PWr")
- Quick press to select value (En = Enabled or dIS = Disabled)
- Rotate ***Temperature Adjustment Knob*** (**G**) knob to select the desired setting
- Quick press to confirm valued

Overload Detection Enable / Disable

- Optional load sensing feature which provides protection against improper positioning and excess load conditions.
- The unit will automatically slow to a safe operating speed and display "E4" if a hazardous load is detected. Error may be cleared by interacting with the unit; shake and timer functions will continue operating as standard.
- Navigate to the System Settings submenu ("SYS")
- Navigate to the Power Recovery setting ("LOAD")
- Quick press to select value (En = Enabled or dIS = Disabled)
- Rotate knob to select the desired setting
- Quick press to confirm valued

Display Settings

- Adjust LCD brightness and configure the display timeout feature, which dims the LCD to set brightness after a selected period of idle time.
- Navigate to the System Settings submenu ("SYS")

- Navigate to the Power Recovery setting (“dISP”)
- Adjust Backlight Brightness
 - Adjust the brightness percentage (“bL P”)
 - Range: 0-100 %
 - Quick press to select the back light brightness percentage value
 - Rotate the knob to select the desired setting
 - Quick press to confirm value
- Backlight Timeout Enable/Disable
 - Quick press to select value (“bLtO”)
 - Rotate knob to select the desired setting (En = Enabled or dIS = Disabled)
 - Quick press to confirm value

Factory Reset

- Resetting the unit to Factory Default Settings will do the following:
 - Clear all programs
 - Clear all recall parameters
 - Clears all System settings to the default value
- Navigate to the System Settings submenu (“SYS”)
- Navigate to the Reset setting (“rSt”)
- Quick press to select value (no or yes)
- Rotate ***Temperature Adjustment Knob*** (**G**) knob to select the desired setting
- Press and hold to confirm

7. RS232

RS-232 serial port provides two-way communications for data logging and unit control by means of a PC and a suitable application program.

- Interface connections chosen in accordance with EIA standard RS-232
- Device Connector: RS-232 DB9 Female
- Recommended Cables:
 - DB9 M/M Straight Through Serial Cable
 - USB-A to RS-232 DB9 Straight Through Serial Adapter Cable
- Transmission Procedure
 - Asynchronous character transmission in start-stop mode
- Transmission Type
 - Full Duplex
- Character Format
 - Start Bit(s): 1
 - Character Bits: 8
 - Parity Bit(s): None
 - Stop Bit(s): 1
- Transmission Speed (Baud Rate): 9600
- Data Flow Control: Xon/Xoff
- Syntax
 - Instructions and parameters separated by space (0x20)
 - CR LF termination (0x0D, 0x0A)
 - Max length: 80 characters
- ‘<command> A’ returned if the command is recognized; ‘L’ returned otherwise.

RS232 Command Set

ID <xxxx>	Set ID 1 ≤ xxxx ≤ 9999
ID	Read ID value
MODE	Read mode: 0: idle 1: shake 2: heating 3: shaking + heating 99: error
START_SHAKE	Start shake function and timer
STOP_SHAKE	Stop the shake function; reset timer
TARGET_SPEED <xxxx>	Set the target speed value
TARGET_SPEED	Reads target speed value
MEASURED_SPEED	Reads measured speed value
START_HEAT	Start heat function and timer
STOP_HEAT	Stop heat function; reset timer
TARGET_TEMP<XX>	Set the target temperature value
TARGET_TEMP	Reads target temperature value
MEASURED_TEMP	Reads measured temperature value
TIMER	Reads current timer value
TIMER 0	Set Timer Mode = Count Up
TIMER <hh>:<mm>:<ss>	Set Timer Mode = Count Down; Functional shut down @ timer = HH:MM:SS
TIMER RESET	Reset timer to set value; timer continues counting from reset value. No change to shake functions.
LOCK	Lock user interface
UNLOCK	Unlock user interface
PARAM <xxxx>	Enable output of comma separated list of parameters, formatted as follows: XXXX = 0: Single parameter dump 0 < XXXX < 9999: Parameter dump every XXXX seconds. <TIMER [HH:MM:SS]>, <ID>, <MODE>, <TARGET SPEED>, <MEASURED SPEED>, <TARGET_TEMP>, <MEASURED_TEMP>, <ERROR CODE>

Troubleshooting

Review the information in the table below to troubleshoot operating problems. During operation, any rattling or ticking sounds may indicate a loose screw on the tray, a tray attachment, or an accessory. All accessories should be sufficiently tightened in place before starting the unit.


Problem	Description	Cause	Solution
N/A	Unit fails to power on	Missing or blown fuse	Add or replace fuse as necessary. If the problem persists, please contact your VWR representative for repair.
E1	RTD Open	RTD fault or over temp	The E1 error should **NOT** be addressed by the end user. Switch the unit off and contact your VWR representative for repairs.
E2	RTD Short	RTD fault or under temp	The E2 error should **NOT** be addressed by the end user. Switch the unit off and contact your VWR representative for repairs.
E3	Speed Error	ceased bearing. drive belt broken mechanical obstruction loose foot (suction cup) *	Flip the standby switch to clear this error and remove the mechanical obstruction. If the E3 error persists the reason may be a ceased bearing or broken drive belt and should **NOT** be addressed by the end user. Switch the unit off and contact your VWR representative for repairs.
E4	Unit overloaded	maximum load exceeded loose foot (suction cup) *	Flip the standby button to clear this error. Be sure the load is within the maximum load capacity before restarting the unit. If the E4 error persists, switch the unit off and contact your VWR representative for repairs

* In the event a foot (suction cup) has come loose from the bench top, the unit will register an errant E04 or E03 error message due to the instability of the unit. Flip the standby switch to clear this error. Firmly press down on the four (4) corners of the unit, creating a strong suction to the work surface (**DO NOT** place on bench mat). Press the standby button to resume operation.)

Accessories

Test Tube Racks

Description	Tube Capacity	Platform Capacity	Cat. No.
10 to 13mm Test Tube Rack	63	1	444-0278
14 to 16mm Test Tube Rack	48	1	444-0279
15mL Centrifuge Tube Rack	35	1	444-0283
50mL Centrifuge Tube Rack	12	1	444-0284
1.5 to 2mL Micro-Tube Rack	70	2	444-0276

Installation of Test Tube Racks - Incubating Mini Shaker

1. Unit should be in standby mode.
2. Align the clip on the rack to the right side of the tray.
3. Press the test tube rack in place, you will hear the test tube rack clip into position. No mounting screws are required.

Installation of 1.5 to 2mL Micro-Tube Rack - Incubating Microplate Shaker

1. Unit should be in standby mode.
2. Align the Micro-Tube Rack side-to-side on the tray.
3. Press the Micro-Tube Rack in place, you will hear the Micro-Tube Rack clip into position. No mounting screws are required.

Flask Clamps

Description	Material	Platform Capacity	Cat. No.
10mL Erlenmeyer Flask Clamp	Stainless Steel	35	444-7040
25mL Erlenmeyer Flask Clamp	Stainless Steel	20	444-7041
50mL Erlenmeyer Flask Clamp	Stainless Steel	12	444-7042
125mL Erlenmeyer Flask Clamp	Stainless Steel	8	444-7043
250mL Erlenmeyer Flask Clamp	Stainless Steel	5	444-7044
125mL Erlenmeyer Flask Clamp	PVC	8	444-2960
250mL Erlenmeyer Flask Clamp	PVC	4	444-2961

Installation of Flask Clamps - Incubating Mini Shaker

1. Unit should be in standby mode.

2. Line up the hole(s) in the flask clamp to the hole(s) in the tray. Using the screw(s) provided, hand tighten into place using a flathead screwdriver. **DO NOT** over tighten.
3. Insert flask into clamp and wrap the spring around the neck of the clamp to hold secure.

Not for use with the Incubating Microplate Shaker.

Technical service

Web Resources

Visit the VWR website at **www.vwr.com** for:

• Complete technical service contact information

• Access to the VWR Online Catalogue, and information about accessories and related products

• Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty

VWR warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation unless any refusal is due to a defect of the product.

Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
CE	VWR Corporation declares that the Incubating Mini and Microplate Shakers comply with directives 2011/65/EU, (EU) 2015/863, 2014/30/EU, 2014/35/EU, and standards EN 50581, EN 61010-2-010, EN 61010-2-051, EN 61326-1. The full text of the EU declaration of conformity is available at the following internet address: https://vwr.com
	This product complies with directive 2012/19/EU. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to https://vwr.com
	EN 61326-1
TÜV SÜD C US	CAN/CSA C22 261010-1, CAN/CSA C22 261010-2-010, CAN/CSA C22 261010-2-051 UL 61010-1, UL 61010-2-010, UL 61010-2-051

Global Notice

Warning: This is a Class B product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

Canada Notice


This Class B digital apparatus complies with Canadian ICES-003.

FCC Notice

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

—Reorient or relocate the receiving antenna.
—Increase the separation between the equipment and receiver.
—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
—Consult the dealer or an experienced radio/TV technician for help.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead, it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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INDIA
Avantor Performance Materials India Pvt. Ltd.,
19 th Floor, Building No. 5, Tower C,
DLF Cyber City Phase - III,
Gurgaon - 122002, Haryana
Tel.: +91-124-4656700
csindia@avantorsciences.com

IRELAND
VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
sales.ie@vwr.com

ITALY
VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
info.it@vwr.com

KOREA
Avantor Performance Materials Korea Ltd
2F ACE Gwanggyo Tower I, Daehak 4ro 17
Yeongtong-gu Suwon, Korea 16226
Tel.: +82 31 645 7250
customerservicekr@avantorsciences.com

THE NETHERLANDS
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
info.nl@vwr.com

MEXICO
VWR International, S.de R.L. de C.V.
Km. 14.5 Carretera
Tlalnepantla-Cuautitlán
Col. Lecheria
Tultitlán Edo. de México
CP 54940
Tel.: +52 (55) 5005 0100
vwrmx@vwr.com

MIDDLE EAST & AFRICA
VWR International FZ-LLC
Office 203, DSP Lab Complex,
Dubai Science Park,
Dubai, United Arab Emirates
Tel: +971 4 5573271
Info.mea@vwr.com

NORWAY
VWR International AS
Brynsalleen 4,
0667 Oslo
Tel.: +47 22 90 00 00
info.no@vwr.com

POLAND
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
info.pl@vwr.com

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GO TO AVANTORSCIENCES.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS FROM YOUR LOCAL TEAM

VWR

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VWR International AB
VWR International GmbH
VWR International, LLC