Instruction Manual

VWR® Digital & Analogue Heavy Duty Vortex Mixers & Microplate Vortex Mixers

EU cat. no

Vortex Mixer Standard Heavy Duty

Euro plug and Swiss adapter: 444-7133

UK plug: 444-7134

Advanced Heavy Duty Vortex Mixer

Euro plug and Swiss adapter: 444-7130

UK plug: 444-7131

Analogue Microplate Vortex Mixers

Euro plug and Swiss adapter: 444-0986

UK plug: 444-0993

Digital Microplate Vortex Mixer

Euro plug and Swiss adapter: 444-0985

UK plug: 444-0992

EN	-	English		1
FR	-	Français		14
ES	-	Español		27
DE	-	Deutsch		40
IT	-	Italiano		54
PT	-	Português		67
BG	-	Български		80
CZ	-	Čeština		84
DK	-	Dansk		88
EE	-	Eesti		92
FI	-	Suomi		96
GR	-	Ελληνικά		100
HR	-	Hrvatska		104
HU	-	Magyar		108
IS	-	Íslenska		112
LT	-	Lietuvių		116
LV	-	Latviešu		120
NL	-	Nederlands		124
NO	-	Norsk		128
PL	-	Polski		132
RO	-	Română		136
SI	-	Slovenščina		140
SK	-	Slovenčina		144
SV	-	Svenska		148

Revision C Issue 05/2026

Legal address of Manufacturer

Europe

VWR International bv,
Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001
Leuven
+3216385011
be.vwr.com

UK Importer:
VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
uk.vwr.com

Country of Origin

United States

Table of Contents

Introduction	**4**
Safety Information	**4**
Package Contents	**4**
Installation	**4**
Maintenance & Servicing	**5**
Enviornment Conditons	**5**
Specifications	**6**
Vortex Mixer Control Panel	**6**
Analogue Control Panel	**6**
Digital Control Panel	**6**
Touch Mode Operation	**7**
Operating a Vortex Mixer	**7**
Operating an Analogue Vortex Mixer	**7**
Operating a Digital Vortex Mixer	**7**
Accessory Installation	**8**
Troubleshooting	**8**
Heavy Duty Vortex Mixer Accessories	**9**
Microplate Vortex Mixer Accessories	**10**
Technical Service	**11**
Warranty	**11**
Compliance with Local Laws and Regulations	**11**
Equipment Disposal	**12**
Distributor Contact Information	**13**

Introduction

This manual contains installation operation and maintenance instructions. Please read the manual completely before using.

Safety Information

Please read the entire instruction manual before operating the Vortex Mixer.

WARNING! DO NOT use the Vortex Mixer in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift the Vortex Mixer by the head. All heads, including cup head, are removable. They can pop off if lifted by them.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly. **DO NOT** immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground - Protective Conductor Terminal

Alternating Current

Package Contents

Heavy Duty Vortex Mixer Analogue or Digital

Description	Quantity
183cm detachable power cord (EU or UK) and Swiss Adapter*	1
Cup Head	1
Universal Holder	1
Universal Holder Cover	1
Foam Insert (1, 5-2,0 mL microtubes, holds 38)	1

Microplate Vortex Mixer Analogue or Digital

Description	Quantity
183cm detachable power cord (EU or UK) and Swiss Adapter*	1
Cup Head	1
Microplate Holder, Single	1

** Dependent on model*

Installation

Upon receiving the VWR Vortex Mixer, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, remove the protective coverings from the feet and place the Vortex Mixer on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is clean and free of dust. Always place the unit on a sturdy work surface.

The Vortex Mixer is supplied with a power cord that should be plugged into a properly grounded outlet. If the supplied power cord does not meet your needs, please use an approved power cord that suits local codes and electric supply. Replacement of the plug must be made by a qualified electrician.

Maintenance & Servicing

The Vortex Mixer is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. **DO NOT** use a cleaning agent or solvent on the front panel which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

Environmental Conditions

Operating Conditions: Indoor use only.

Temperature: 18 to 33°C

Humidity: 20% to 85% relative humidity, non-condensing

Altitude: 0 to 2000 M above sea level

Non-Operating Storage:

Temperature: -20 to 65°C

Humidity: 20% to 85% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.


Heavy Duty Vortex Mixer Specifications

Overall dimensions (L x W x H): 24,1 x 16,8 x 16,5cm
Electrical (50/60 Hz): 0,13 amps | 30 watts
Fuses: 5mm x 20mm, 5 amp 250 volt quick acting
Speed range:
On Mode: 300-2500rpm
Touch Mode: 300-3500rpm
Orbit: 4,9mm
Duty rating: Continuous duty
Controls:
Analogue: Touch/off/on - rocker switch
Speed knob: variable 1 to 10 dial marks
Digital: See page 6
Capacity: 1,1kg
Ship weight: 6,8kg

Microplate Vortex Mixer Specifications

Overall dimensions (L x W x H): 26,6 x 13,7 x 11,4cm
Electrical (50/60 Hz): 0,13 amps | 30 watts
Fuses: 5mm x 20mm, 5 amp 250 volt quick acting
Speed range:
On Mode: 300-2500rpm
Touch Mode: 300-3500rpm
Orbit: 3,5mm
Duty rating: Continuous duty
Controls:
Analogue: Touch/off/on - rocker switch
Speed knob: variable 1 to 10 dial marks
Digital: See page 6
Capacity: 2 Microplates
Ship weight: 5,4kg

Analogue Vortex Mixer Control Panel


The front panel of the Analogue Vortex Mixer contains all the switches and control knob needed to operate the unit.

A. 3-way rocker switch: touch/off/on rocker switch starts/stops the vortexing function.

B. Speed control: The speed control knob controls the vortexing speed. User should turn knob clockwise until the vortexer reaches the desired speed.

Digital Vortex Mixer Control Panel

The front panel of the Digital Vortex Mixer contains all the switches, controls and displays needed to operate the unit.

A. 3-way rocker switch: touch/standby/on rocker switch starts/stops the vortexing function.

B. Time display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9,999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9,999.

C. Up/down arrows for set-point control.

D. Speed display: Displays the speed of the Vortex Mixer.

E. Up/down arrows for set-point control.

Vortex Mixer Introduction

Your Vortex Mixer is ready for most one-handed applications. The Vortex Mixer operates by mixing samples just prior to testing.

Touch Mode Operation

NOTE: Due to the robust mechanical design of this unit, more downward force is required to operate this unit in Touch Mode to mix single tubes as compared to other models. Another option is to set your unit to the ON mode for mixing a single tube.

To ensure the stability of your vortex mixer, it is very important to start with a smooth, clean work surface and the mixers (4) feet. If the mixer starts to “walk”, remove the power to the unit and clean the work surface and the (4) feet with isopropyl alcohol. Press the unit down firmly on the work surface to ensure the feet are properly bonded.

Analogue Vortex Mixer Operating Instructions

1. Make sure the 3-way rocker switch is in the center, off position. Plug the power cord into a properly grounded outlet.
2. For continuous operation, push the rocker switch to the right, on position. For intermittent/touch operation, push the rocker switch to the left, touch position.
3. In either case, turn the speed knob to the desired setting. In the on position, you will see the head in motion immediately. In the touch position, the vessel you are using must be pushed down on the head attachment to achieve motion. To get the desired mix, vary the speed, using the speed knob, and/or vary the angle of contact and pressure against the head.
4. When finished with either continuous or intermittent/touch operation, return rocker switch to the center, off position.
5. The continuous, on mode is for the accessories.
6. The intermittent, touch mode is intended for short mixing times of one (1) minute or less at full speed.
7. To completely cut off power to the unit, disconnect the power cord from the wall outlet.

Digital Vortex Mixer Operating Instructions

1. Make sure the 3-way rocker switch is in the center, standby position. Plug the power cord into a properly grounded outlet. The unit is now on and in a standby mode.
2. For continuous operation, push the rocker switch to the right, on position. For intermittent/touch operation, push the rocker switch to the left, touch position.
3. In either case, set the speed. In the on position, you will see the head in motion immediately. In the touch position, the vessel you are using must be pushed down on the head attachment to achieve motion. To get the desired mix, vary the speed, using the speed up/down arrows, and/or vary the angle of contact and pressure against the head.

Setting speed:

To set the speed press the up/down arrows below the speed display until you reach the desired speed.

- Speed while running in the on position is 300 to 2500rpm
- Speed while running in the touch position is 300 to 3500rpm

Setting time:

Time operation works in both touch and on positions.

- To run in timed mode (programmed time), start by pressing the rocker switch to the left, touch position. Press the up/down arrows below the time display until you reach the desired time. While running unit in timed mode, the time display will show time remaining, counting down to zero (0:00).
- To run in continuous mode (accumulated time), reset the timer to zero (0:00) before running. Time display will show the accumulated time and run until you press the rocker switch to the center, standby position.

4. When finished with either continuous or intermittent/touch mode, return rocker switch to the center, standby position.
5. The continuous, on mode is for accessories.
6. The intermittent, touch mode is intended for short mixing times of one (1) minute or less at full speed.
7. To completely cut off power to the unit, disconnect the power cord from the wall outlet.

Beeper Reference:

To silence

the beeper operation (except for error codes), with the unit in standby mode, remove accessories or cup head. Press and hold the head mount. Press the time down arrow. Release the time down arrow. Release the head mount. To restore to factory setting, repeat procedure.

Accessory Installation Instructions (All Units)

Installing the cup head, microplate trays, universal holder, or other head attachments:

1. Move the 3-way rocker switch to the center, off/standby position.
2. Place accessory on head mount of the unit, aligning the square on the bottom of the accessory with the square on the top of the unit. Lower accessory until it is sitting on the head mount.
3. Turn the accessory counterclockwise until it snaps into place (approximately 1/8 of a turn). Add samples to attachments after the accessory is secured onto the Vortex Mixer. Be sure samples are secure before turning unit on.
4. To remove an installed accessory, move the 3-way rocker switch to the center, off/standby position. Turn the accessory **clockwise** (approximately 1/8 of a turn) and then lift the accessory straight up. Remove samples before removing the accessories.
 - To load a microplate into the microplate tray, position the microplate against the foam bumpers first, then push down the other side to ensure a secure fit. To release a microplate, push it against the microplate tray's foam bumper and lift to remove.
 - Maximum speed may vary when vortexing microplates. Speed may be dependent on the height of the skirt on the microplate or material of the microplate.
 - Rubber vibration dampeners on the bottom of the microplate tray and universal holder may need to be reinserted into the proper position when vortexing at high speeds. Squeeze firmly as you maneuver it back into position.

Accessory Installation Instructions (Heavy Duty Vortex Mixer)

Installing the vessel holders:

1. Attach the Universal Holder to the vortex mixer.
2. Place the Grip Mat in the bottom of the Universal Holder.
3. Place the flask on top of the Grip Mat, in the center of the holder.
4. Pull the center hole of the Vessel Holder over the sample container.
5. With the arrow on the Vessel Holder close to the front of the unit, the straps should line up with the posts of the Universal Holder.
6. Push the hole on the strap over the posts on the holder. Secure the strap so that it fits into the slot on the holder. Repeat this for the remaining two straps.
7. Be sure the sample is secure before turning unit on.

Installing the Universal Holder Cover:


1. Place the Foam Insert into the Universal Holder.
2. Lay the Universal Holder Cover on top of the foam and align the edges.
3. Push the hole on the cover down over the post on the holder. Squeeze the strap so that it fits into the slot on the holder, leaving the cover sitting flat on the foam. Repeat this process for the remaining two posts.
4. To remove the accessory cover, secure the straps and pull straight up.

Troubleshooting (All Units)

Review the information in the table below to troubleshoot operating problems.

Error	Cause of Error	How to Fix
Analogue Models Motor does not run in the touch or on mode **Digital Models** E3 Error	**Analogue and Digital Models** There is no motion on the motor or the motor is not working properly	**Analogue and Digital Models** 1. Turn the rocker switch to the center off/standby position 2. Check to verify the load has not been exceeded 3. Check and clear any possible obstructions of the head mount 4. If the unit still does not work, please contact your VWR representative

Accessories for the Heavy-Duty Vortex Mixer

Cup Head

Replacement for cup head supplied with Heavy Duty Vortex Mixer. Designed for mixing one tube at a time.

Cat. No. 444-7136

Universal Holder

Replacement for Universal Holder supplied with Heavy Duty Vortex Mixer to mix irregularly shaped objects when Universal Holder Cover is in place. Also used as a holder for the foam tube inserts. Measures 15,5x15,0x7,6cm.

Cat. No. 444-7139

Universal Holder Cover

Replacement for Universal Holder Cover supplied with Heavy Duty Vortex Mixer. The cover allows for mixing irregularly shaped objects.

Cat. No. 444-7140

Microplate Holder (single)

Designed to hold one microplate.

Cat. No. 444-7137

Microplate Holder (double)

Designed to hold two standard microplates.

Cat. No. 444-7138

Microplate Holder (quad)

Designed to hold four standard microplates.

Cat. No. 444-7151

Stackable Microplate Holder Four

Designed to maximize the capacity of the Heavy Duty Vortex Mixer to eight microplates by stacking the tray on top of the four plate holder.

Cat. No. 444-7152

Accessories for the Heavy-Duty Vortex Mixer

Foam Insert for 0,5mL microtubes

Foam insert holds (52) 0,5mL microtubes.

Cat. No. 444-7141

Foam Insert for 1,5 to 2,0mL microtubes

Foam Insert for 1,5 to 2,0mL microtubes
Foam insert holds (38) 1,5 to 2,0mL microtubes.

Cat. No. 444-7142

Foam Insert for 12-13mm test tubes

Foam insert holds (34) 12-13mm diameter test tubes.

Cat. No. 444-7143

Foam Insert for 15-18mm test tubes

Foam insert holds (20) 15-18mm diameter test tubes. Ideal for 15mL centrifuge tubes.

Cat. No. 444-7144

Foam Insert for 19-21mm test tubes

Foam insert holds (18) 19-21mm diameter test tubes.

Cat. No. 444-7145

Foam Insert for 22-25mm test tubes

Foam insert holds (13) 22-25mm diameter test tubes.

Cat. No. 444-7146

Foam Insert for 26-29mm test tubes

Foam insert holds (4) 26-29mm diameter test tubes. Ideal for 50mL centrifuge tubes.

Cat. No. 444-7147

Foam Insert Blank

Foam insert without tube openings.

Cat. No. 444-7153

Small Vessel Holder

Rubber holder secures 125 and 250mL Erlenmeyer flasks. Vessel holder also includes a grip mat.

Cat. No. 444-7148


Large Vessel Holder

Rubber holder secures 500 and 1000mL Erlenmeyer flasks. Vessel holder also includes a grip mat.

Cat. No. 444-7149

Accessories for the Microplate Vortex Mixer

Cup Head

Replacement for cup head supplied with Microplate Vortex Mixer. Designed for mixing one tube at a time.

Cat. No. 444-7136

Microplate Holder (single)

Replacement for Microplate Holder supplied with Microplate Vortex Mixer. Designed to hold one microplate.

Cat. No. 444-7137

Microplate Holder (double)

Designed to hold two standard microplates.

Cat. No. 444-7138

Technical service

Web Resources

Visit the VWR website at **www.vwr.com** for:

• Complete technical service contact information

• Access to the VWR Online Catalogue, and information about accessories and related products

• Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty


VWR International warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with Local Laws and Regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization unless any refusal is due to a defect of the product.

Equipment disposal

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycleend by handing it over to an authorized facility for separate collection and recycling. It is also you responsibility to decontaminate the equipment in case of biological, chemical, and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.
For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.
Thank you

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India
VWR Lab Products Private Limited
No.139. BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India – 560058
Tel.: +91 80 28078400
Email: vwr_india@vwr.com

Ireland
VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 01 88 22 222
Email: sales.ie@vwr.com

Italy
VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02-3320311
Email: info.it@vwr.com

Korea
Avantor Performance Materials Korea Ltd
2F ACE Gwanggyo Tower I
Daehak 4ro 17
Yeongtong-gu Suwon
Korea 16226
Tel.: +82 31 645 7250
Email: saleskorea@avantorsciences.com

Middle East & Africa
VWR International FZ-LLC
Office 203, DSP Lab Complex,
Dubai Science Park
Dubai, United Arab Emirates
Tel.: +971 4 5573271
Email: info.mea@vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 020 4808 400
Email: info.nl@vwr.com

Norway
VWR International AS
Brynsalleen 4
0667 Oslo
Tel.: +47 22 90 00 00
Email: info.no@vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 058 32 38 200
Email: info.pl@vwr.com

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