NUCLEOSIL® Carbohydrate

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NUCLEOSIL® Carbohydrate

Note: All HPLC columns from MACHEREY-NAGEL are supplied with a certificate, which contains specifications and test results of the column.
The NUCLEOSIL® Carbohydrate column is a quality product based on the robust silica NUCLEOSIL®. This column is specifically developed for HPLC analysis. Due to the chemical nature of the phase, the lifetime of the column highly depends on the measurement and the treatment after measurement. Consequently, prior to column installation, you should familiarize yourself with the contents of this instruction leaflet. If carefully and properly used excellent chromatographic results and long column lifetime can be achieved. This product can be used for the separation of mono-, di- and trisaccharides and their quantitative determination. All HPLC columns must exclusively be used in accordance with universally accepted laboratory regulations and HPLC working methods. Before running the column the entire analytical system (column and equipment) has to be carefully checked by the operator. Chromatographic conditions (mobile phase, flow, temperature etc.) must be adapted to the analytical task. MACHEREY-NAGEL does not give any warranty and is not liable for the success of a separation or application. If you have any questions after reading this leaflet, please call our service / technical support.

Table of contents

- Safety indication
- Description of the column
- Application note
- Installation
- Guard columns
- Sample
- Eluent
- Flow rate and pressure
- Temperature
- Detection
- Equilibration
- Column storage
- Troubleshooting
- Column regeneration
- Abstract

Safety indication
Follow the general safety instructions for handling of HPLC solvents used as mobile phases (e.g., acetonitrile, methanol) and take precautions against any kind of injuries or damage to health (e.g., skin and eye protection in case of broken capillaries). Disposal of used HPLC columns must follow international, national and local environmental protection regulations. The use of HPLC columns is only permitted to staff members, who are qualified in their field. Keep HPLC columns away from children. MACHEREY-NAGEL disclaims and excludes all warranties of any kind or nature whatsoever and MN shall not be liable for any damages (whether direct, indirect, foreseeable, incidental, compensatory, consequential or special), whether based upon warranty, contract, tort or strict liability, if damages and / or losses occur caused by improper use, maintenance, neglect or improper treatment (especially opening of the column and exposure of the column bed).


Description of the column
As stationary phase the NUCLEOSIL® Carbohydrate column contains an amino phase, based on spherical silica and modified by a special procedure. This column is almost exclusively operated under aqueous organic elution conditions. Separation is achieved by polar interactions with the stationary phase. Probably hydrogen bonding effects are involved. Especially for reversed phase applications like the determination of sugars the column is supplied with the eluent acetonitrile – water (79:21, v/v).

Application note
The figure below shows the rapid separation of fructose, glucose, saccharose, maltose and lactose. Such separations are often required for the quality control of honey, chocolate or ice cream.

Separation of sugars

Column: EC 250/4 NUCLEOSIL® Carbohydrate
REF 720905.40
Eluent: acetonitrile – water (79:21, v/v)
Flow rate: 2 mL/min
Temp.: 25 °C
Detection: RI
Injection: 10 µL

Peaks:
1. Fructose
2. Glucose
3. Saccharose
4. Maltose
5. Lactose

MN Appl. No. 102480

Installation
The column should be installed in the flow direction indicated on the column label. It is connected with 1/16" capillaries and fittings, typical for HPLC instruments.

Guard columns
For protection and an extension of column lifetime the column should always be used with a guard column. The filter elements and the adsorbent in the guard column retain contaminants from the sample or the eluent. Connection of the guard column with the separation column is made by a suitable guard column holder (see ***www.mn-net.com*** or the MN chromatography catalog). Cartridge replacement is required when increased column pressure and / or loss of performance is observed.

Sample
Sample solutions should be passed through a syringe filter (e.g., CHROMAFIL® Xtra PET, 0.45 µm, 25 mm, REF 729220) before entering the column. If injected sample solutions are still turbid even after filtration, the lifetime of the column may be significantly reduced. The sample volume should be as small as possible to achieve an optimal resolution.

Eluent
NUCLEOSIL® Carbohydrate columns are supplied with the eluent acetonitrile – water (79:21, v/v). Important sugars can be isocratically determined with this eluent composition (see application note). For further applications the composition of the organic aqueous eluent can be varied. A use in gradient mode is also possible. Thus, oligosaccarides with longer chains can be separated using a gradient with an increasing aqueous part. In certain cases addition of buffer (e.g., 10 mmol/L Tris-phosphate, Na citrate, Na acetate) improves the separation by adjustment to an optimum pH value. Also, small concentrations of organic amines (0.01 – 0.1 %) can enhance the separation. A pH value below 2 and above 9 should always be avoided. Strong acidic or basic conditions can result in dissolution of the column bed or the organic modification. The amount of buffer salts should be kept as low as possible. Note the solubility limit of the buffer in the eluent. An increase of the organic portion can result in precipitation of buffer salts and plugging of the column. Before start of operation with an eluent containing a buffer the column should be first preconditioned with a minimum of 10 column volumes acetonitrile – water (25:75, v/v). Furthermore, depending on the pH value of the eluent, the NH_2 group can be protonated. Therefore, it acts as a weak anion exchanger, which alters the retention behavior. Especially bulky buffer salts like phosphate can develop ionic interactions, thereby effectively shielding the polar NH_2 group. These phenomena lead to a retention time shift. After finishing measurements with buffer containing eluents, or if problems with a retention shift happen after this application, the column should be regenerated (see column regeneration). A changing to normal phase mode is not recommended. If necessary, it should only be made with an intermediate flushing step with THF. Eluents should be always filtered through a 0.2–0.45 µm membrane filter and degassed.

Flow rate and pressure
Flow rate (recommended for analytical columns with 2–4.6 mm ID: 0.2–2.0 mL/min) influences the required time, the resolution and the column lifetime. It is limited by the back pressure, which should not exceed the maximum of 400 bar. In mixtures of methanol and water viscosity reaches a maximum at about 40 % methanol. For this reason a reduced flow rate is recommended, when changing the eluent composition. We recommend controlling back pressure regularly. If a high pressure results from the use of the column at nominal flow rates, this usually indicates that some deposition of contaminants has occured on the packing material, which must be removed (see troubleshooting).

Temperature
Column temperatures up to 60 °C are possible; for a long lifetime 30–40 °C is recommended. However, they should be at least 30 °C below the boiling temperature of the eluent, in order to ensure proper detection. Variation of the temperature influences retention times and the peak shape. Optimum temperatures for successful separations should be determined empirically.

Detection
Spectrophotometers, refractometers and electrochemical detectors can be used with the column.

Equilibration
Prior to measurement of samples the column must be rinsed with the eluent at the same flow rate and temperature as the method to be applied. Column equilibration is finished, when the baseline of the detector no longer shows a drift (generally after 10 column volumes).

Column storage
The original eluent acetonitrile – water (79:21, v/v) is recommended for storage. For long-term storage mobile phases containing inorganic salts are not recommended (see regeneration). Methanol is also not recommended for a longer storage, because of a possible impurity with metal ions (e.g., iron(III)). For column storage be sure the end fittings are tightly sealed using column end plugs, because storage without these seals can result in drying of the packing material. Under these circumstances rinse the column with approx. 10 column volumes of the eluent of storage at a flow rate of max. 0.2 mL/min.

Troubleshooting
The following outline describes the symptoms of performance loss and their cause. All columns are subject to the strict regulation and control of our quality assurance system. Columns based on silica are robust and hold their separation efficiency for long periods by correct maintenance and treatment. According to experience, column failures are mostly a result of injection of contaminants to the sorbent bed. The usage of a guard column, as well as an appropriate sample pretreatment will help to minimize these risks.
Use the outline below to help determine the cause of a possible performance loss:

Symptom / Error / Cause	Prevention / Remedy
Baseline drift	
· insufficient period for equilibration with the eluent	longer or better equilibration
· contaminated eluent	use freshly prepared solvents and reagents
· temperature	column temperature control
Broad peaks	
· mixing and / or diffusion before / behind the column	keep length and ID of capillaries at a minimum
· too large sample volume	smaller injection volume
Peak interference; too fast elution	
too fast elution and / or insufficient separation by:	
· improper column temperature or flow rate	optimize concerned parameter
· elution power of eluent is too high	optimize eluent system
Increasing back pressure; degradation of the separation performance	
contamination of sorbent by:	
· particulate accumulation on frit or sorbent bed from sample, eluent or system	prepare fresh eluent; prefilter samples and eluent, use in-line filter / rinse LC system, clean the sorbent beforehand check solubility of buffer salts / remove them by rinsing (see column regeneration)
· precipitation of buffer salts	
Insufficient separation; degradation of the separation with regular column pressure	
contamination with:	
· fats, oils, lipids from sample (coating of sorbent surface) and other organic substances from improperly prepared eluent or matrices	remove organic substances by sample preparation / clean the sorbent (see column regeneration)
protonation of NH_2 group by:	
· usage of acidic buffers as eluent	deprotonation by rinsing with weakly basic solutions (see column regeneration)
Double peaks (dead volume)	
· faulty fittings (capillaries, ferrules, nuts)	use “PEEK Fingertight Fittings”, REF 718770 or REF 718778 / replace fittings
· dissolution of silica by too high pH value of eluent	consider pH range of column / replace column

Column regeneration
In some cases the performance of the column can be restored by removing contaminants from the sorbent bed or by regeneration of the phase. It is important, however, to locate the source of contamination before using the column for the analysis of samples again.
1. Prepare fresh eluent: Sometimes the performance loss is traced to eluent contamination. Therefore, prepare fresh eluent and flush all liquid lines before using the column again. The eluent should be filtered through a 0.2–0.45 µm membrane and degassed prior to use.
2. Cleaning of sorbent: To remove contamination rinse the column with a minimum of 10 column volumes (see table below) at the original flow rate and temperature as follows:
 - acetonitrile – water or methanol – water (10:90, v/v) for removal of the buffer
 - 100 % methanol to remove polar organic compounds
 - 100 % acetonitrile to remove medium polar organic compounds (possibly T= 40 °C)
 - 100 % tetrahydrofuran to remove nonpolar organic compounds
 - if necessary, 100 % tetrahydrofuran with inverse flow direction at 1/5 of original flow rate
 - convert column to storage condition with acetonitrile – water (79:21, v/v) at original flow rate

 An adequate indicator for a clean column is a constant baseline. At constant temperature you should observe less than 2–3 mAU drift during a running time of 5 minutes with an isocratic run.
3. Regeneration: After finishing measurement and for a possible retention shift after the usage of buffer, rinse as follows (30–50 min/step, flow rate of max. 1 mL/min):
 - acetonitrile – 20 mmol/L ammonium acetate or formate, pH 8–8.5 (95:5, v/v)
 - acetonitrile – water (10:90, v/v)
 - gradually increase the part of acetonitrile in steps of 20 % to the storage conditions
4. Column replacement: The above procedures will restore performance only in certain cases. Some organic contaminants are particularly refractory and may not respond to treatment. Also dead volume, due to column compression can generally not be repaired. Under these circumstances, column replacement is necessary. It is highly advisable to locate the cause of the problem before installing a new column.

1 column volume (250 mm length x 4 mm ID column) ≙ 3,15 mL

Abstract
To extend column lifetime, please keep in mind the following:
1. As eluents organic aqueous eluent systems (e.g., acetonitrile – water or buffer) are recommended. Eluents should be filtered through a 0.2–0.45 µm membrane and degassed.
2. Filter samples through a 0.2–0.45 µm CHROMAFIL® Xtra PET syringe filter before injection.
3. Use a guard column for samples with difficult matrices.
4. The recommended flow rate for analytical columns (ID 2–4.6 mm) is 0.2–2.0 mL/min.
5. Adjust flow rate to keep column pressure below 400 bar.
6. Store the column in acetonitrile – water (79:21, v/v).
7. Use analytical grade reagents and HPLC grade solvents for all work. Discard any solutions that show evidence of bacterial growth.

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MACHEREY-NAGEL GmbH & Co. KG
MACHEREY-NAGEL AG
MACHEREY-NAGEL Inc.