VWR

Microscope VisiScope® series 300

INSTRUCTION MANUAL

Model	European Catalogue Number
TL384 M	630-1939

Version: 1
Issued: 27, 05, 2014

CE

Legal Address of Manufacturer

Europe
VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
http://be.vwr.com

Country of origin: **ITALY**

Table of Contents

Warning
Safety Information
Package Contents
Intended use
Symbols and conventions
Unpacking
Product specifications
Overview
Using the microscope
Troubleshooting
Repair and maintenance
User replaceable accessories and spare parts
Technical service
Warranty
Compliance with local laws and regulations
Disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.
We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.
We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.
Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Package Contents

DESCRIPTION	QUANTITY
Microscope stand with nosepiece, stage, condenser	1
Optical head (trinocular)	1
Objective PLAN IOS MET 4x	1
Objective PLAN IOS MET 10x	1
Objective PLAN IOS MET 20x	1
Objective PLAN IOS MET 50X	1
Eyepiece WF10X/20mm	2
Metallographic attachment	1
Polarizer & Analyzer filters	2
Glass insert for sample support	1
Dust cover	1
Power supply output 6Vdc	1

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take care to ensure that the optical items (objectives and eyepieces) do not fall out and get damaged. Using both hands (one around the neck and one around the base), lift the microscope from the container and put it on a stable desk.

Place the observation head onto the top of the neck and tighten the lock-screw. Insert the eyepieces into the eye tubes.

Connect the provided power supply to the power supply input jack on the rear of the microscope.

Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Symbols and conventions


The following chart is an illustrated glossary of the symbols that are used in this manual.


⚠	**CAUTION** This symbol indicates a potential risk and alerts you to proceed with caution

Product specifications

Head:	Trinocular, 30° inclined, 360° rotating. Interpupillary adjustment 48-75 mm.
Eyepieces:	WF10X/20mm.
Nosepiece:	5-position reversed revolving nosepiece. Ball bearing rotation.
Objectives:	PLAN IOS MET 4x/0.10, 10x/0.25, 20x/0.40, 50x/0.75.
Focusing system:	Coaxial coarse and fine.
Stage:	Double layer mechanical sliding stage, dimensions 216x150mm, moving range 78x54mm. Belt-drive in X direction. Glass insert for sample support.
Condenser:	Abbe condenser, sliding-in, N.A. 1.25 with centering system.
Illumination:	Transmitted light: P-LED^3, with manual brightness control.
Power supply:	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A.
Metallurgical Attachment:	High-power white P-LED^3; field and aperture diaphragms.
Filters:	Blue, polarizer & analyzer.

Overview

PHOTO PORT
DIOPTRIC ADJUSTMENT RING
EYEPIECE
LOCK SCREW
NOSEPIECE
OBJECTIVE
STAGE
TRANSLATION KNOBS
CONDENSER
IRIS DIAPHRAGM
CONDENSER-CENTRING SCREWS
FINE FOCUSING KNOB
LED ILLUMINATOR
FOCUS TENSION

LED COLLIMATING
LENS
LED HOUSE
ANALYZER
FILTER HOLDER
POLARIZER
FILTER HOLDER
FIELD DIAPHRAGM
APERTURE
DIAPHRAGM

CONDENSER
ADJUSTMENT KNOB
COARSE
FOCUSING KNOB
FOCUS STOP KNOB
BRIGHTNESS
ADJUSTMENT KNOB

Using the microscope

Adjust the observation head
Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

Place the specimen on the stage
Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage. A glass plate is provided in order to support metallurgical samples.

LED lamp settings – Transmitted illumination
The microscope is fitted with high brightness white LED illuminator. Press the selector switch on the back of the main body to I position. Turn the brightness adjustment knob to a brightness suitable for observation.

LED lamp settings – Epillumination
The microscope is fitted with high brightness white LED illuminator. Press the selector switch on the back of the main body to II position. Turn the brightness adjustment knob to a brightness suitable for observation.

Adjust interpupillary distance
Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

Focus and diopter adjustment
Turn the dioptric adjustment ring on the left eye tube to the “0” position. Loosen the focus-stop knob, rotate the coarse focusing knob to bring the slide into focus using an objective with a low magnification, and then lock the focus-stop knob. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left dioptric adjustment ring until a sharp and clear image can be observed also with the other eye. Turn the tension-adjust-knob to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.

Condenser adjustment
Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser centring screws can be used to centre the condenser.

Set the numerical aperture
Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution. There are two aperture diaphragms on the microscope: for transmitted light and for epillumination light.

Adjust the field diaphragm for epillumination
Move the lever on the metallographic attachment to fully close the field diaphragm. Then slowly open it and stop when the whole field of view is illuminated. This condition prevents the light that falls outside the field of view from being scattered and producing noise.

Polarisation and colour filters
The blue, yellow and frosted glass filters can be inserted in the slot in the metallographic attachment, just in front of the Led house. In the same slot you can insert the polarizer filter. The analyzer filter can be inserted in the slot below the optical head, after removing the protection cap.
Turn the graduated wheel on the analyzer filter to perform polarisation analysis of your sample.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
LIGHT DOESN'T TURN ON	Power supply not connected	Check that the 6Vdc power supply jack is well inserted on the rear of the microscope
		Rotate the brightness adjustment control and check if there is an increase in the light output
IMAGE CANNOT BE SEEN OR IS DARK.	The iris diaphragm aperture is not completely opened	Completely open the iris diaphragm aperture
	Brightness level is low	Rotate the brightness adjustment potentiometer
	Objective is not aligned with the optical axis	Rotate the nosepiece until an objective is well inserted in the optical path (it "clicks")
IMAGE IS UNCLEAR, BLURRED OR HAS INSUFFICIENT CONTRAST.	Objectives or filters are dirty	Wipe them clean
	The iris diaphragm aperture is not opened correctly	Open the iris diaphragm aperture completely
	Condenser at wrong height	Rotate the condenser knob until you see a uniform illumination

Repair and maintenance

Microscopy environment
This microscope is recommended for use in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

After using the microscope

The microscope should always be kept vertical when moving it so that no moving parts, such as the eyepieces, fall out.
Never mishandle or impose unnecessary force on the microscope. Never attempt to service the microscope yourself. Turn down the illumination intensity control and turn the light off. Cover the microscope with the included dust cover, and keep it in a dry and clean place.

Electrical safety precautions

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics
If the optical parts require cleaning first use compressed air.
If that is not sufficient use a soft lint-free cloth with water and a mild detergent.
And as a final option use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics. Do not disassemble objectives or eyepieces in an attempt to clean them.

For the best results, use the VWR cleaning kit (see catalogue number below).

If you need to send the microscope to manufacturer for maintenance, please use the packaging if possible.

User replaceable accessories and spare parts

DESCRIPTION	QUANTITY	CAT. NO.
Eyepiece WF10x/20mm for BL/TL 254 and 300	2	630-2000
Eyepiece WF15x for BL/TL 254 and 300	2	630-2001
Eyepiece WF20x for BL/TL 254 and 300	2	630-2002
Micrometer eyepiece WF10x/20mm for BL/TL 254 and 300	1	630-2003
Slide micrometric 26x76mm 1mm for 300	1	630-1650
Objective E PLAN IOS 4x for serie 200&300	1	630-1639
Objective E PLAN IOS 10x for serie 200&300	1	630-1640
Objective E PLAN IOS 20X for serie 200&300	1	630-1641
Objective E PLAN IOS 40X for serie 200&300	1	630-1642
Objective E-PLAN IOS 60x for serie 300	1	630-2049
Objective E PLAN IOS 100x for serie 200&300	1	630-1643
Objective PLAN IOS MET 50x/0,75	1	630-2159
Objective PLAN IOS MET 4x/0.10	1	630-2160
Objective PLAN IOS MET 10x/0.25	1	630-2161
Objective PLAN IOS MET 20x/0.40	1	630-2162
Phototube adapter for APS-C sensor	1	630-1645
Heating stage	1	630-1980
Battery solar pack	1	630-1637

Technical service

Web Resources

Visit the VWR's website at **www.vwr.com** for:
• Complete technical service contact information
• Access to VWR's Online Catalogue, and information about accessories and related products
• Additional product information and special offers
Contact us For information or technical assistance contact your local VWR representative or visit. **www.vwr.com.**

Warranty


VWR International warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.
Instead it is your responsibility to correctly dispose of your equipment the end of its life cycle by handling it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.
For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland
VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy
VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

New Zealand
Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Norway
VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax. 058 32 38 205
E-mail: info@pl.vwr.com

Portugal
VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore
VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain
VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden
VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland
VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey
Pro Lab Laboratuar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu Işmerkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK
VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

GO TO VWR.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS OF YOUR LOCAL VWR DISTRIBUTOR

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