REF 945023 / 945024 05.23

de

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Method:
Determination of the acute toxicity of **liquid-dried** luminescent bacteria according to DIN EN ISO 11348-2. The measuring unit is the natural light output (luminescence) of the microorganism *Vibrio fischeri* NRRL B-11177. Obtained will be the level of inhibition of the light intensity in the sample in comparison to an uninhibited control. The level of inhibition of the light intensity in the sample is a measure of its toxicity.

Area of application:
Waste water, aqueous extracts and leachates, fresh water (surface water and ground water), seawater, brackish water, eluates of sediments (fresh water brackish water and seawater), interstitial pore water, and single substances dissolved in water.

Measurement range:
0 – 100 % inhibition

Indication of the results:
- **% Inhibition** of the light intensity in the sample in comparison to an uninhibited control
- **G_L-Value** (according to DIN EN ISO 11348-2: **LID-Value = Lowest Ineffective Dilution**): Reciprocal value of the first dilution of a sample that causes an inhibition of the light intensity of less than 20 %.
- **EC_{xx}**: Concentration of a sample that causes an inhibition of the light intensity of exactly xx % (e.g. EC_{50}-Value: Concentration of the sample that causes an inhibition of the light intensity of 50 %).
- **TU-Values** (**T**oxicity **U**nits according to the definition of the American Environmental Protection Agency U.S. EPA): 100 divided by EC_{50}-Value.

Contents:	**REF 945023**
sufficient for	200 determinations
	10 vials with liquid-dried luminescent bacteria according to DIN EN ISO 11348-2
	10 vials with "BioFix® *Lumi* Reactivation Solution" according to DIN EN ISO 11348-2, portioned ready-to-use
	1 flask with ready-to-use sodium chloride standard solution
	REF 945024
sufficient for	400 determinations
	20 vials with liquid-dried luminescent bacteria according to DIN EN ISO 11348-2
	20 vials with "BioFix® *Lumi* Reactivation Solution" according to DIN EN ISO 11348-2, portioned ready-to-use
	1 flask with ready-to-use sodium chloride standard solution

Safety precautions:
This test kit does not contain any hazardous substances for which labelling is obligatory.
The luminescent microorganism *Vibrio fischeri* NRRL B-11177 has never appeared as a pathogen up to now. In accordance with the code of practice B006 1/92 ZH 1/346 of the employer's liability insurance association of the chemical industry *Vibrio fischeri* is assigned to the risk group 1 i.e. there is no risk for humans and vertebrates.

Storage:
The test reagents can be kept at -15 °C to -21 °C until the date indicated on the packaging. The liquid-dried luminescent bacteria should only be removed from the freezer compartment immediately prior to their reactivation. Reactivated luminescent bacteria should be used within 4 hours! Intermediate storage of reactivated luminescent bacteria must only be in undiluted state (see step 1) in a fridge at +2 °C to +8 °C.
Upon storage of the reactivated luminescent bacteria a reduction of the natural light emission and a shift of the sensitivity spectrum can be possible. It is not recommended to refreeze the reactivated luminescent bacteria and is not subject to the guarantee obligations of MACHEREY-NAGEL.

Reactivation and dosage of the liquid-dried BioFix® *Lumi* luminescent bacteria:
The **reactivation** of liquid-dried BioFix® *Lumi* luminescent bacteria is always done using the same method regardless of the test. The **dosage** of the resuspended luminescent bacteria depends on the optional test method of the luminescent bacteria test according to DIN EN ISO 11348-2 without measuring GL1 (step 2.1) or with measuring GL1 (step 2.2).

Important note:
Further details to the test procedures can be derived from the instructions for use of the corresponding tests, test systems and handbooks of the luminometers used.

Step 1: Reactivation

1. Remove a deep frozen vial of **"BioFix® *Lumi* Reactivation Solution"**, thaw out, and maintain a temperature of **+15 °C** in a suitable thermoblock (or water bath).
2. Remove a deep frozen vial with **BioFix® *Lumi* luminescent bacteria** from the freezer immediately before the reactivation and add **1 mL "BioFix® *Lumi* Reactivation Solution"**.
3. **Defrost and dissolve** the luminescent bacteria in a water bath (room temperature) for **2 min**.
4. Then place the vial with the resuspended luminescent bacteria in a suitable thermoblock (or water bath) and maintain a temperature of **+15 °C** for **15 min**.

Step 2.1: Dosage for measurement without GL1

1. **Mix** the resuspended luminescent bacteria **with the remaining "BioFix® *Lumi* Reactivation Solution"**.
2. Transfer **0.5 mL** of the **reactivated luminescent bacteria suspension** to each provided cuvette.
3. Maintain at **+15 °C** for **15 min**.
4. Start the test by **measuring** the **initial light intensity I_0**.

Important! After the measurement of the initial luminescence and prior to the addition of the control or sample solution it is important to keep the cuvettes at +15 °C.

Details for further performance of the test procedure has to be found in the instructions and manuals of the luminometer used.

or

Step 2.2: Dosage for measurement with GL1

1. **Add 4 mL of "BioFix® *Lumi* Reactivation Solution"** to the resuspended luminescent bacteria and mix well.
2. Place **0.2 mL** of this suspension in the **test cuvettes for GL1 calculation** (control and sample mixtures).
3. Mix the remaining resuspended luminescent bacteria with the remaining "BioFix® *Lumi* Reactivation Solution".
4. In the remaining test cuvettes provided, transfer **0.5 mL** of the **reactivated luminescent bacteria suspension**.
5. Maintain at **+15 °C** for **15 min**.
6. Start the test by **measuring** the **initial light intensity I_0**.

Important! After the measurement of the initial luminescence and prior to the addition of the control or sample solution it is important to keep the cuvettes at +15 °C.

Details for further performance of the test procedure has to be found in the instructions and manuals of the luminometer used.

Analytical quality assurance:
The compliance of certain validity criteria is required by DIN EN ISO 11348-2. The compliance is guaranteed by MACHEREY-NAGEL by means of the enclosed test certificate.
According the above mentioned norm the sensitivity of the luminescent bacteria has to be controlled by the end user. The required information such as reference substances of concentrations can be derived from the enclosed test certificate.
A test independent of the standard according to DIN EN ISO 11348-2 can be carried out using the accompanying ready-to-use sodium chloride standard solution. Mixing 0.5 mL of this standard solution with 0.5 mL resuspended luminescent bacteria (in accordance with step 2.1) results in a 40 % to 60 % inhibition of the light intensity after 30 min incubation at +15 °C

pH-value:
According to the DIN EN ISO 11348-2 no correction of the pH is required if the pH value of the sample is between 6.0 and 8.5. It may be noted that pH-values of less than 6.0 and greater than 8.5 can lead to pH-related inhibitions of luminescence.

Interferences:
Non or hardly soluble or volatile substances or compounds which are difficult to dissolve, or substances which react with the diluent or the test suspension or which change their state during the test can falsify the test result or make it difficult to reproduce the test results.
Strong colored or turbid samples can cause losses of light due to light adsorption or scattering. In accordance with the procedure described in DIN EN ISO 11348-2, annex A, these effects can be minimized by an appropriate sample treatment (centrifugation or sedimentation) or by using color correction cuvettes (REF 940006).
Oxygen is required for the bioluminescence (> 3 mg/L) and can be limiting in case of samples with a strong oxygen consumption. This can finally result in an inhibition of light emission.
If the water sample contains easily biodegradable nutrients a loss of the light intensity which is not caused by contaminants can occur.
As the test organism *Vibrio fischeri* is a marine organism the use of seawater samples can lead to an increase of luminescence masking potential inhibitory effects of the sample (see DIN EN ISO 11348-2 annex D).
Salt concentrations of the sample higher than 30 g/L NaCl or concentrations of other compounds with similar osmolarity can in conjunction with the salts in the test solutions lead to hyperosmotic effects. To prevent these effects the overall salt concentration should not exceed the osmolarity corresponding to 35 g/L NaCl.
Clorine containing samples effect the test procedure so that a dechlorination is mandatory, e.g., by addition of a 1 % sodium thiosulfate solution.

Disposal:
Luminescent bacteria and test preparations can be disposed of without difficulty via the sink. Restrictions can arise if the sample contains hazardous or toxic substances which are subject to waste disposal regulations. The user him / herself is responsible for the orderly disposal of such test preparations in accordance with the respectively valid guidelines and regulations.

Literature:
DIN EN ISO 11348-2:2009-05
Determination of the inhibitory effect of water samples upon the light emission of *Vibrio fischeri* (luminescent bacteria test); part 2; procedure with liquid-dried bacteria

AQS Data Sheet for Quality Assurance in Water, Waste Water and Silt Tests by the Federal Work Community Water of Germany (LAWA) P-9/8

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0 · info@mn-net.com · ***www.mn-net.com***

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MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germania