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DURAN WHEATON KIMBLE

II DURAN® Glass

Scope
Material data sheet / specification valid for all articles made of DURAN® Glass that are manufactured by the DWK Life Sciences GmbH. DURAN® Glass is conforming to the requirements of ISO 3585.

Trade name: DURAN®

Material name: Borosilicate glass 3.3 (special glass)

Material use: Manufacture of laboratory, industrial and special industrial glassware, electrical home appliance, housewares and of glass for application in pharmaceutical manufacturing processes.

Composition
DURAN® Glass has the following approximate composition (weight %)

SiO_2	B_2O_3	Na_2O & K_2O	Al_2O_3
81 %	13 %	4 %	2 %

Chemical and physical properties
DURAN® Glass is practically inert. However an ion exchange between the medium and the glass can not be excluded.

Coefficient of expansion	Acid resistance	Acid resistance	Hydrolytic resistance	Alkalis resistance
3.3 x 10-6 K-1	Class S 1 (DIN 12116)	≤ 100 µg Na_2O per dm² (ISO 1776)	Class HGB 1 (ISO 719) / Class HGA 1 (ISO 720)	Class A 2 (ISO 695)

GOST 13905-2005 Glass containers
(Method of testing the water resistance of inner surface)

Test	Description	Units	Test Result
GOST 13905-2005	Hydrochloric acid consumption	cm3	< 0.01
GOST 13905-2005	Water resistance	cm3	< 0.01

ASTM-Compliance
DURAN® Glass is compliant with “Type I, Class A low-expansion borosilicate glass” as defined by ASTM-E438-1992 2.

Pharmacopoeial compliance
DURAN® glass is classified as a “Type I neutral glass” as defined by the current editions of the following pharmacopoeias: USP <660> Containers – Glass Ph. Eur. (3.2.1) Glass Containers for Pharmaceutical Use JP <7.01> Test for Glass Containers for Injections.

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Spezifikation /specification / spécificités / specifikacija DURAN®

T_g DURAN®: 525 °C
If heated above 525 °C, the material softens.

Recommended storage conditions

Protect from humidity, and large temperature fluctuations.
Storage under normal conditions in accordance with standard climatic conditions +23°C and 50% relative humidity in accordance with ISO 554:1976.

Technical properties	
Thermal conductivity λ_W 90 °C	1,2 W • m^{-1} • K^{-1}
Mean specific heat capacity C_p (20 °C; 100 °C)	0,8 J • g^{-1} • K^{-1}
Young`s modulus	64 GPa
Poisson's ratio μ	0,20
Stress-optical coefficient K (DIN 52314:1977-11)	4,0 • 10^{-6} mm^2 • N^{-1}
Refractive index (λ = 587,6 nm) n_d	1,473
Dielectric constant (1 MHz 25 °C)	4,6
Loss factor tan δ	37 • 10^{-4}
Temperature for specific electrical resistance 10^{-6} Ω • cm (DIN 52326) t^k 100	250 °C
Thermal expansion α 20 °C; 300 °C (ISO 7991:1997)	3,3 • 10^{-6} K^{-1}
Transformation temperature T_g	525 °C
Density ρ bei 25 °C	2,23g • cm^{-3}
Temperatures at the viscosities η in dPa • s: 10^{13} (upper cooling point) $10^{7,6}$ (softening point) 10^4 (processing temperature)	 560 °C 825 °C 1260 °C
CAS no.	65997-17-3

Possible hazard

DURAN® is no dangerous material and DURAN® has no health-threatening effects. DURAN® has no dangerous ingredients corresponding to legal chemical requirements and to hazardous materials regulations in actual issuance. DURAN represents no danger to human health or the environment.

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