Instruction manual

VWR® Professional Hotplate-Stirrer

EU Cat No.

444-1510	Ceramic 18 x 18 cm Hotplate Stirrer
444-1511	Ceramic 18 x 18 cm Hotplate Stirrer Cert
444-1512	Ceramic 25 x 25 cm Hotplate Stirrer
444-1513	Ceramic 25 x 25 cm Hotplate Stirrer Cert

Revision A 11/2025

EN

Legal address of Manufacturer

Europe

VWR International bv,
Researchpark Haasrode 2020
Geldenaaksebaan
464 B-3001 Leuven
+3216385011
be.vwr.com

UK Importer:
VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
uk.vwr.com

United States
VWR International LLC
100 Matsonford Rd
Radnor, PA 19087
+1 800-932-5000
www.vwr.com

Country of Origin

USA

EN

Table of Contents

1. Introduction
 - 1.1 Safety Information 4
 - 1.2 Intended Use 5
 - 1.3 Package Contents 5
 - 1.4 Installation 5
2. Overview
 - 2.1 Specifications 6
 - 2.2 Control Panel 8
 - 2.3 Interface 8
 - 2.4 Display 9
3. Operation
 - 3.1 Heating 10
 - 3.1.1 Heating Operation Tips 10
 - 3.1.2 Heating Operations 11
 - 3.2 Stirring 12
 - 3.2.1 Stirring Operation Tips 13
 - 3.3 Timer Control 13
 - 3.4 Single Point Calibration 14
4. Menu
 - 4.1 Enter or Exit the Menu 15
 - 4.2 Menu Map 15
 - 4.3 Heat Plate Temperature Limitation – TLIM 15
 - 4.4 Single Point Calibration – SPC 15
 - 4.5 System – SYS 16
5. RS232 16

Troubleshooting 18
Technical Service 19
Warranty 19
Equipment Disposal 20

EN

1. Introduction

This manual contains installation, operation and maintenance instructions for VWR Professional Hotplate-Stirrer. Please read the manual completely before using.

1.1. Safety Information

Definition of Signal Warnings and Symbols

WARNING	**For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.**
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
Attention	For important information about the product. May lead to equipment damage if not avoided.
Note	For useful information about the product.

Warning Symbols

Caution

Explosion Hazard

Electrical Shock

Caution, hot surface

Protective Conductor Terminal

Alternating Current

Safety Precautions

WARNING! The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

- Always operate unit on a level surface for best performance and maximum safety.
- DO NOT lift unit by the top plate.
- Disconnect unit from the power supply prior to maintenance and servicing.
- Spills should be removed promptly, after the unit has cooled down.
- Alkalis spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.
- DO NOT immerse the unit for cleaning.
- DO NOT operate the unit at high temperatures without a vessel/sample on the top plate.
- DO NOT operate the unit if it shows signs of electrical or mechanical damage.
- Protective earthing of the equipment is achieved via connection of the provided power cord to a compatible grounded power outlet.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.

WARNING! unit is not explosion proof. Use caution when unit is on or when heating volatile materials.

WARNING! DO NOT use the unit in explosive atmospheres or with materials that could cause a hazardous environment from processing. Keep in mind the material blink point relative to the target temperature that has been set.

CAUTION! The top plate can reach up-to 550°C, DO NOT touch the heated surface. Use caution at all times. Keep the unit away from explosive vapors and clear of papers, drapery and other flammable materials. Keep the power cord away from the heater plate.

Beware of the following risks when heating:

- flammable materials
- low boiling point combustible substances
- glass breakage as a result of mechanical shaking power
- incorrect container size
- too much medium
- unsafe condition of container

Earth Ground - Protective Conductor Terminal

Alternating Current

1.2. Intended Use

This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from VWR, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations: however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

1.3. Package Contents

Check that there's no damage occurred during shipment. Notify the carrier for the damages found at the time of unpacking.
Packing list:

- Hotplate-Stirrer
- (20.3cm) Stainless Steel PT1000 Temperature Probe
- Power Cord
- Stir Bar

1.4. Installation

Upon receiving the VWR Hotplate-Stirrer, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage, the carrier must be notified immediately.

After unpacking, place the Hotplate-Stirrer on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of six (6) inches from vertical surfaces. Do not position the equipment such that it is difficult to disconnect the power cord during use. Always place the unit on a sturdy work surface.

The Hotplate-Stirrer is supplied with a 3 conductor, grounded power cord that should be plugged into a matching standard grounded outlet. If the cord supplied does not meet your needs, please use an approved power cord that has ratings equal or exceeding those of the originally provided cord and that complies with the local/national regulations of the country in which the equipment is to be used. Replacement of the plug must be made by a qualified electrician.

2. Product Overview

2.1. Specifications


Model	Professional 18 x 18 cm Ceramic Hotplate Stirrer
Function	Heating and stirring
Top Plate Size	18 x 18 cm
Top Plate Material	Ceramic
Control	Digital
Display	115 x 30 mm White Character LCD
Timer	1 minute to 99 hours 59 minutes
Temperature Range	Ambient +5°C – 550°C
Temperature Resolution	0.5°C
Temperature Readability	0.5°C
Temperature Stability of Top Plate*	±1% (> 100 °C), ± 1 °C (≤ 100 °C)
Temperature Stability with Temperature Probe**	± 0.5 % (> 100 °C), ± 0.5 °C (≤ 100 °C)
SPC Point	3 points for plate temperature 3 points for probe temperature
Probe Connection	PT1000 (Class A)
Wattage	1050W
Speed Range	50-1800 RPM
Speed Resolution	10 RPM
Speed Stability**	+/-2%
Stirring Capacity	20L H2O
Weight Capacity	23kg
Data Logging and Remote Control	Built-in RS232
IP Rating	IP32
Housing Material	Painted Aluminum
Dimensions (L x W x H)	32.9 x 22.1 x 11.6 cm
Net Weight	4.6 kg
Electrical (50/60 Hz)	220-240V, 7A

Note:
***The parameter in the specification table is applicable to 5 cm diameter center of the top plate.**
**** Conditions permitting. Variations in temperature and speed measurement processes, vessels, ambient and sample will impact actual performance. To improve temperature accuracy of the system, use the Single Point Calibration feature.**

Model	Professional 25 x 25 cm Ceramic Hotplate Stirrer
Function	Heating and stirring
Top Plate Size	25 x 25 cm
Top Plate Material	Ceramic
Control	Digital
Display	115 x 30 mm White Character LCD
Timer	1 minute to 99 hours 59 minutes
Temperature Range	Ambient +5°C – 500°C
Temperature Resolution	0.5°C
Temperature Readability	0.5°C
Temperature Stability of Top Plate*	±1% (> 100 °C), ± 1 °C (≤ 100 °C)
Temperature Stability with Temperature Probe**	± 0.5 % (> 100 °C), ± 0.5 °C (≤ 100 °C)
SPC Point	3 points for plate temperature 3 points for probe temperature
Probe Connection	PT100, PT1000 (Class A)
Heating Power	2066W
Speed Range	50 – 1800 RPM
Speed Resolution	10 RPM
Speed Stability**	+/-2%
Stirring Capacity	22L H2O
Weight Capacity	23kg
Data Logging and Remote Control	Built in RS232
IP Rating	IP32
Housing Material	Painted Aluminum
Dimensions (L x W x H)	40.2 x 28.7 x 11.8 cm
Net Weight	6.7 kg
Electrical (50/60 Hz)	220-240V, 10A


Note:
***The parameter in the specification table is applicable to (5 cm) diameter center of the top plate.**
**** Conditions permitting. Variations in temperature and speed measurement processes, vessels, ambient and sample will impact actual performance. To improve temperature accuracy of the system, use the Single Point Calibration feature.**

2.2. Control Panel


#	Controls	Functions
A	**Left Knob**	Temperature Knob
B	**Right Knob**	Speed Knob
C	**Display Screen**	Show operation status, parameters, and menu settings

2.3. Interface

#	Description
D	**External Temperature Probe Port**
E	**RS232 Port**
F	**Power Entry Module**
G	**Rod Holders**

2.4. Display

Heating Display		General Display		Stirring Display	
A	Real-Time Temperature Display	**J**	Timer Display	**K**	Stirring Speed Display
B	Heating Indicator			**L**	Stirring Indicator
C	Set Temperature Display			**M**	Set Speed Icon
D	°C Icon			**N**	Speed Unit – RPM Icon
E	Set Temperature Icon				
F	Single Point Calibration Icon				
G	External Probe Connection Icon				
H	Set Temperature Limit Icon				
I	Hot Top Indicator				

Attention:
The **Hot Top Indicator Icon** will illuminate when the top plate temperature is >= 40°C.
This icon will continue to light up when the instrument is switched Off.

3. Operations

This section covers topics and instructions of the instrument operation.
The topics would include:

Topics
Heating
Stirring
Timer Control
Single Point Calibration (SPC)

3.1. Heating

This section covers topics, including heating operation tips and instructions for heating operation. Proper heating operation is crucial in laboratory settings, where precise temperature control is essential for accurate experimentation.

Topics
Heating Operation Tips
Heating Operations

3.1.1. Heating Operation Tips

Measuring the Sample Temperature

The temperature display on the unit represents the estimated top plate temperature, not the sample temperature. The vessel contents being heated may be at a lower temperature depending on the size and thermal conductivity of the vessel. It may be beneficial to monitor the temperature of the vessel contents and adjust the setpoint temperature accordingly.
When external probe is in use, the temperature display on the unit represents the sample temperature measured by the probe.
For precise sample temperature control, we recommend using the External Temperature Probe.
For details of working with External Temperature Probe to measure sample temperature, please refer to *Measuring the Sample Temperature with External Temperature Probe*

Heating Overshoot

The unit may overshoot the temperature before stabilizing at the setpoint. The two methods to minimize overshoot are:

- Metal containers minimize overshoot.

 Caution! When heating metal containers on a ceramic top plate, it is recommended to use the lowest temperature setting possible to limit thermal stress to the ceramic top plate.

- If a glass vessel is used, start with a temperature setpoint 5 to 10°C below the desired temperature. When the temperature stabilizes at this lower setting, increase the heater to the final temperature.

Typical Time to Boil Water

The chart below is an example of an approximate time to boil for the specified amount of water in a specific vessel. These values are only approximate and can vary from unit to unit. Values are based on 23°C water in an ambient environment of 23°C.

Unit Size	Heater Temp. Limit	Volume of Water	Typical Time to Boil
18 x 18 cm	500°C	500mL in 1L beaker	≈ 15 min
25 x 25 cm	500°C	500mL in 1L beaker	≈ 17 min

3.1.2. Heating Operations

This section covers the following topics:

Topics
Start Heating
Adjust Temperature while Heater is On
Turn off the Heater
Measure the Sampe Temperature with External Temperature Probe
Set a Heating Plate Temperature Limitation
Disable the Heating Function

3.1.2.1. Start Heating

1. Rotate the **Left Knob** to adjust the target temperature.
2. Long press the **Left Knob** for 2 seconds to start heating.

When the heater is turned on:

- The **Heating Indicator** will be illuminated.
- You will hear a beep sound (if the beeper is turned on).
- The **Heating Indicator** will blink when the temperature is ramping.
- The **Heating Indicator** will stop blinking when the temperature has reached the target.

3.1.2.2. Adjusting Temperature while Heater is On

1. When the heater is on, rotate the **Left Knob** to adjust the temperature set-point.
2. Confirm the edit of temperature set-point. After rotating the **Left Knob**, the **Set Temperature Display** will start to blink. Short press the **Left Knob** to confirm the edit.

3.1.2.3. Turn off the Heater

Long press the **Left Knob** and release the knob when you hear a beep sound, or release in 1.5 second. When heating is stopped, the green light will fade.

3.1.2.4. Measuring the Sample Temperature with External Temperature Probe

1. Connect the external probe to the External RTD Probe Port on the rear side of the unit.
2. Adjust the temperature set-point, then long press the **Left Knob** to start.

Once the External Temperature Probe is connected:

- The External **Probe Icon** on the screen will illuminate.
- The display will show the temperature of the sample measured by the external probe instead of the heating plate temperature.

Note: The hot top plate icon  will still illuminate once the heater temperature reaches 40°C.

Attention: DO NOT plug or unplug the external temperature probe while the heater is on.

3.1.2.5. Set a Heating Plate Temperature Limitation

1. Press and hold the **Left Knob** to enter the Menu.
2. Rotate the **Left Knob** to navigate to **TLIM**.
3. Short press the **Left Knob** to edit the temperature limitation. The digits will blink while editing.
4. Short press the **Left Knob** again to confirm the setting.

When the temperature limitation is adjusted:

- The temperature limitation icon  will illuminate.
- The temperature limitation icon will blink when setting the temperature up to the limitation.

3.1.2.6. Disable the Heating Function

Heating function can be disabled for heating prohibited applications. Set Temperature Limitation (TLIM) to 0 to disable the heating function.

1. Press and hold the **Left Knob** to enter the Menu.
2. Rotate the knob to navigate to TLIM.
3. Press the knob to edit the temperature limitation to “0”. The digits will blink while editing.
4. Press the **Left Knob** again to confirm the setting.

When the heating function is disabled:

- The temperature limitation icon will illuminate.
- The temperature setting digits set would display **0**.
- The heating function is always disabled even when the actual temperature is below Zero.
- When user attempts to start heating, both the temperature setting digits and the TLIM symbol will blink. And the timer won’t work.

3.2. Stirring

This section covers topics, including stirring operation tips and instructions for stirring operation. Proper stirring operation is crucial in laboratory settings, where precise speed control is essential for accurate experimentation.

Topic
Stirring Operation Tips
Stirring Operation

3.2.1. Stirring Operation Tips

The stirrer increases speed at a steady rate unit the set-point is reached. The stirrer may not be able to reach its setpoint for the following reasons:

- The stir bar may be too large.
- The liquid may be too viscous.
- The magnetic strength of stir bar has reduced over time.

As the stirring speed will vary according to liquid viscosity, stir bar length, stir bar magnetic strength, and distance from top plate. Adjust one or all of these to achieve the desired stirring speed. For example: the closer the reaction vessel is the top plate, the strong magnetic connection between the unit and the stir bar.

3.2.2. Stirring Operation

This section covers the following topics:

Topics
Start Stirring
Adjust Speed while the Stirrer is On
Turn off Stirrer

3.2.2.1. Start Stirring

1. Rotate the **Right Knob** to adjust the target speed
2. Long press the **Right Knob** for 1.5 seconds to start stirring.

When the stirrer is turned on:

- The **Stirring Indicator** will be illuminated.
- The unit will beep if the beeper is turned on.
- The **Stirring Indicator** will blink when the speed is ramping.
- The **Stirring Speed Display** will show the current speed and the target speed in the repeating cycle.
- The **Stirring Indicator** will stop blinking when the speed has reached the target.

3.2.2.2. Adjust Speed while the Stirrer is On

1. When the stirrer is on, rotate the **Right Knob** to adjust the speed set-point.
2. Confirm the edit of speed set-point. After rotating the **Right Knob**, the **Stirring Speed Display** will start to blink. Short press the **Right Knob** to confirm the adjustment.

3.2.2.3. Turn off Stirrer

Long press the **Right Knob** and release the knob when you hear a beep sound, or release in 1.5 seconds.
When stirring is stopped, the green light will fade.

3.3. Timer Control

By default, the Timer will be set to 00:00 and count upwards when the heating or stirring functions are turned on.
Users can also set a Count-Down Timer to specify a time for heating or stirring the samples.
Once the set time is reached, the device will automatically switch off, ensuring that the samples are not overheated or over-stirred.
This section covers the following topics:

Topics
Set a Timer
Reset the Timer
Timer Starting Criterion

3.3.1. Set a Timer

To set a timer:
1. Short press the knob to enter timer configuration.
 The screen will display “HH:MM” and then minutes setting will start to blink.
2. Rotate the knob to set the “minutes” of timer. Then short press the knob to confirm.
3. Rotate the knob to set the “hours” of timer. Then short press the knob to confirm.
 The timer is now adjusted. The screen will display the value of the timer.

3.3.2. Reset a Timer

To reset the timer
1. Short press the knob to enter timer configuration.
2. Press and hold the knob for 2-3 seconds, until the timer is reset to 00:00.

3.3.3. Timer Starting Criterion

Start the Timer as soon as Heating or Stirring is On

Set the Timer Start Setting (TMDE) to STD to start Timer when heating or stirring is switched on.
1. Press and hold the **Left Knob** to enter Menu.
2. Rotate the **Left Knob** and navigate to **SYS** > **TMDE**
3. Select **STD**
 The Timer will start as soon as heating or stirring is switched on.

Start the Timer as soon as Heating Reached the Target Temperature

Set the Timer Start Setting (TMDE) to Temperature Dependent Mode (TEMP) to start Timer when reaching the target heating temperature.
1. Press and hold the **Left Knob** to enter Menu.
2. Rotate the **Left Knob** and navigate to **SYS** > **TMDE**
3. Select **TEMP**

The timer will start as soon as heating or stirring is switched on.

3.4. Single Point Calibration (SPC)

3.4.1 Start a Single Point Calibration (SPC)

1. To enter the **Single Point Calibration**, the unit must first be at the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top-level menu. Then wait for the screen to show **SPC**.
2. Briefly press the left knob to enter the **SPC** settings menu.
3. Briefly press the left knob again to change the calibration temperature. The temperature will begin to blink to indicate that it can be modified.
4. Rotate the left knob to scroll to the desired temperature.
5. Press and hold the left knob until the unit beeps and the illuminates to begin calibration at the set temperature.
 a) The unit will begin to heat to the set temperature.
 b) The **Single Point Calibration Icon** on the display will blink to indicate that the calibration is running.
 c) The left and right knobs will be disabled until the calibration is completed.
 d) If the external probe is connected, the stirrer will turn on simultaneously.
 e) To cancel the calibration while it is running, flip the standby switch off to turn off the unit.
6. Once the unit has reached the calibration temperature, the **Single Point Calibration Icon** and the calibration temperature will blink.
7. With a secondary temperature measurement device, measure the temperature of the top plate, or the heated sample at the location of the external probe if using probe control.
8. Rotate the left knob to scroll to the measured temperature from the secondary temperature measurement device, and then briefly press the left knob to select the new temperature.
9. The unit will begin to regulate temperature with compensated error. When this is complete, **SAVE** will appear and blink on the screen.

10. Measure the temperature at the same location as step 7.
11. Rotate the left knob to select from:
 a) **SAVE** – to retain calibration (stirring will stop if using probe control).
 b) **ADJ** – to prompt fine tuning of calibration (return to step 8).
 c) **RJCt** – to cancel the single point calibration and return to the initial calibration menu.
12. If you select **SAVE**, the screen will return to the top level of the settings menu.
13. Rotate the left knob to scroll to **ESC** (EXIT)
14. Briefly press the left knob to return to the main operating screen.
15. If you select **RJCt** in step 11, the screen will return to the top level of the settings menu. Rotate the left knob to **ESC** (EXIT) and briefly press the left knob to return to the main operating screen.

3.4.2 Adjust a Calibrated Temperature

1. To adjust a calibrated temperature, the unit must first be in the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top level menu. Then wait for the screen to show **SPC**.
2. Briefly press the left knob to enter the **SPC** settings menu.
3. Rotate the left knob to select which calibration temperature to adjust if two calibration records are stored. The unit can store 1 point (Plate Controlled Heating Temperature) and 1 point (Probe Controlled Heating Temperature) calibration record. For calibration performed with an external probe, the screen will display **External Probe Icon**.
4. Press and hold the left knob until the unit beeps and the **Heater Indicator** illuminates to begin calibration. Then repeat steps in the previous **Start a Single Point Calibration** Section from step 5 to the end to finish the process.

3.4.3 Clear a Calibrated Temperature

1. To clear a calibrated temperature, the unit must first be at the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top-level menu. Then wait for the screen to show **SPC**.
2. Briefly press the left knob to enter the **SPC** settings menu.
3. Rotate the left knob to scroll to **CLr**.
4. Briefly press the left knob to enter the menu. If there is no stored SPC point, select **BAcK** to return to the previous screen.
5. Rotate left knob to scroll to the calibration point you want to clear. The calibration points are stored in ascending order by temperature.
6. Briefly press the left knob to confirm your selection.
7. Rotate the left knob to select from:

a) **YES** – to clear the selected calibration temperature.

b) **NO** – to return to **SPC** settings menu.

8. Briefly press the left knob to confirm your selection and return to the Calibration settings menu.

3.4.4 Use a Calibrated Temperature

1. To heat to a calibrated temperature, the unit should return to the main operating screen.
2. Rotate the left knob to scroll the heat setting to the stored calibration temperature. The **SPC** icon will appear once the temperature is correct.
3. Press and hold the left knob until the unit beeps and the illuminates. The unit will heat to the stored calibration temperature.

4. Menu

4.1. Enter or Exit the Menu

Enter the Menu

Note: Heating and Stirring need to be turned off to enter the Menu.

1. Press and hold the knob until **MENU** is displayed.
 MENU would be displayed for 2 seconds, and then the first sub-menu will be displayed for configuration.
2. Rotate the knob to select sub-menu. Press the knob to enter the sub-menu.
3. Configure the settings.
4. Press the knob to confirm the adjustment, and back to the Menu.

Exit the Menu

1. If in the sub-menu, rotate the knob to **BACK**, then press the knob to go back to the Menu.
2. From the Menu, rotate the knob to **EXIT**, then press the knob.

4.2. Menu Map

Menu	Sub-Menu	Refer to
TLIM	Max °C – 0°C adjustable	*Heat Plate Temperature Limitation – TLIM*
SPC	• **NEW** • **CLR**	*Single Point Calibration – SPC*
SYS	• Beeper **BEEP** • Timer Start Setting **TMDE** • Power Recovery **PWRR** • Reset **RSET** • System Version **V 1.01** (Example display of system version)	*System – SYS*

4.3. Heat Plate Temperature Limitation – TLIM

TLIM allows users to limit the maximum temperature of the heating plate for safety concerns regarding sample blink-point. Heating function can be disabled by setting the **TLIM** limit value at 0°C. For details, please refer to *Set a Heating Plate Temperature Limitation* Temperature Limit Icon will be displayed when a temperature limit is activated.	

4.4. Single Point Calibration – SPC

Menu Options	Definition
NEW	Set a new SPC point. Storage of SPC points: • Top Plate: 3 • External Probe: 3
CLR	Clear SPC point

Note: For calibration procedures, please refer to *Single Point Calibration (SPC).*

4.5. System – SYS

SYS Sub-Menu	Definition	Options
BEEP	Set to enable or disable beeper.	• Enable **EN** (default) • Disable **DIS**
TMDE	Timer Start Setting controls the initiation of the timer in count or countdown situations	• **STD** (default): timer start when heating or stirring is on • Temperature dependent **TEMP:** the timer will start once the heater temperature or probe temperature reaches the set temperature in count or countdown situations.
PWRR	If power has been disconnected, Power Recovery allows the unit to automatically start heater and stirrer functions again when power is returned.	• Enable **ON** • Disable **OFF** (default)
RSET	Restore the instrument to factory settings.	• **YES** (long press the knob to confirm selection) • **NO**
V 1.01[1]	Display of the current system version.	-

Note 1: Example display of system version

5. RS232

RS-232 serial port provides two-way communications for data logging and unit control by means of a PC and a suitable application program.

Setup

- Interface connections chosen in accordance with EIA standard RS-232
- Device Connector: RS-232 DB9 Female
- Recommended Cables:
 - DB9 M/M Straight Through Serial Cable
 - USB-A to RS-232 DB9 Straight Though Serial Adapter Cable
- Transmission Procedure
 - Asynchronous character transmission in start-stop mode
- Transmission Type
 - Full Duplex
- Character Format
 - Start Bit(s): 1
 - Character Bits: 8
 - Parity Bit(s): None
 - Stop Bit(s): 1
- Transmission Speed (Baud Rate): 9600
- Data Flow Control: Xon/Xoff
- Syntax
 - Instructions and parameters separated by space (0x20)
 - CR LF termination (0x0D, 0x0A)
 - Max length: 80 characters
- '<command> A' returned if the command is recognized; 'L' returned otherwise.

User Commands

RS232 User Commands	
ID <XXXX>	Set ID [1-9999] Returns ID value [1-9999] if <XXXX> is blank
MODEL	Returns equipment model ID
SERIAL	Return production serial number
VERSION	Returns software revision
MODE	Returns value: 0: Idle 1: Heating, Plate Controlled 2: Heating, Probe Controlled 3: Stirring 4: Heating (Plate) & Stirring 5: Heating (Probe) & Stirring 99: Error
START_HEAT	Start heat function and timer
STOP_HEAT	Stop the heat function, reset timer if stir is not active
START_STIR	Start stir function and timer
STOP_STIR	Stop the stir function, reset timer if heat is not active
TARGET_TEMPERATURE <XXXX>	Set the target temperature value Returns the target temperature value if <XXXX> is blank
TARGET_SPEED <XXXX>	Set the target speed value Returns the target speed value if <XXXX> is blank
MEASURED_TEMPERATURE	Returns measured temperature value Returns both plate and probe values if in probe mode
MEASURED_SPEED	Returns measured speed value
TIMER	Returns current timer value
TIMER <HH>:<MM>:<SS>	Timer Mode = Count Down Run end @ timer = HH: MM: SS
TIMER_RESET	**Mode = Idle:** Reset timer to 00:00:00 (count up) **Mode = running (1-5):** Reset timer to set value; timer continues counting from reset value
LOCK	Lock user interface
UNLOCK	Unlock user interface
PARAM XXXX	XXXX = 0: Single parameter dump 0 < XXXX < 9999: Parameter dump every XXXX seconds <Timer [hh:mm:ss]>, <ID>, <MODE>, <TARGET TEMP>, <MEASURED TEMP>, <TARGET SPEED>, <MEASURED SPEED>, <ERROR CODE>.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

Problem	Cause	Solution
N/A	Unit fails to power on Missing or blown fuse	Not fixable by user, please contact A VWR representative.
E1	Plate RTD Disconnected	Not fixable by user, please contact A VWR representative.
E2	Plate RTD Error	Not fixable by user, please contact A VWR representative.
E3	Stir Error	Flip the standby switch off and back on. If the problem persists, please contact a VWR representative.
E4	Probe RTD Disconnected	Flip the standby switch off and back on. If the problem persists, please contact a VWR representative.
E5	Probe RTD Error	•Flip the standby switch off, and remove the probe from the unit. Then, flip the standby switch on to return the normal operating mode. •If this does not solve the issue, replace the external probe or the main PCBA and then try again. •If the problem persists, please contact a VWR representative.
E7	User probe error. An External Probe is connected to the unit while heating is on.	Flip the standby switch off and back on. If the problem persists, please contact A VWR representative.
E8	Plate over temperature	Not fixable by user, please contact A VWR representative.
E9	Heating not working	Not fixable by user, please contact A VWR representative.
E10	Triac Fault	Not fixable by user, please contact A VWR representative.
Ac Err	Mains frequency out of the range 40~55Hz (rating 50Hz) or 55~70Hz (rating 60Hz)	Regulate mains frequency within the range.

Technical Data

Equipment Ratings:
Indoor use only

Altitude:	0 to 2000m
Operating temperature:	5°C to 40°C
Storage temperature:	-20°C to 65°C
Operating humidity:	20 to 80% relative humidity, non-condensing
Storage humidity:	20 to 80% relative humidity, non-condensing
Voltage fluctuations:	Mains supply voltage fluctuations up to ±10% of the nominal voltage.
Overvoltage category (Installation category):	II
Pollution degree:	2

COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
CE	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), (EU)2015/863 (RoHs), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at www.VWR.com/
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.VWR.com/
	EN 61326-1
TÜV SÜD C US	CAN/CSA-C22.2 No. 61010-1-12, CAN/CSA-C22.2 No. 61010-2-010, CAN/CSA-C22.2 No. 61010-2-051
	UL 61010-1, UL 61010-2-010, UL 61010-2-051

Global Notice
Warning: This is a Class B product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

ISED Canada Compliance Statement:
CAN ICES-003(B) / NMB-003(B) FCC Notice

ISO 9001 Registration
The management system governing the production of this product is ISO 9001 certified.

Technical service

Web Resources

Visit the VWR website at **www.vwr.com** for:

- Complete technical service contact information
- Access to the VWR Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty


VWR warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.
Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.
For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.
By doing so, you will help to conserve natural and environmental resources, and you will ensure that your equipment is recycled in a manner that protects human health.
Thank you

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VWR International AB
VWR International GmbH
VWR International, LLC