Cole-Parmer®

Cole-Parmer® Paraffin Wax Dispenser - WD-200

- ‘Wax on tap’ for histology and pathology labs
- Over-temperature cut-out at 120°C
- Non-drip heated tap to prevent blockages
- Anti-microbial coating (Ag+) to inhibit growth of bacteria
- Slimline design to maximise bench space



coleparmer.com

antylia SCIENTIFIC



Cole-Parmer® Paraffin Wax Dispenser - WD-200

The Cole-Parmer® paraffin wax dispenser is designed for on-demand delivery of molten wax. It heats wax to a liquid state and provides convenient 'wax on tap' for Histology and Pathology laboratories. The paraffin wax dispenser has a large 7.5 L capacity tank, for melting up to 6 kg of pelletised wax at a time.

The innovative design of the instrument ensures that it takes up minimal bench space in the laboratory. The wax dispenser has digital temperature control which provides accurate temperature control of the molten wax, and is coupled with an ultra-fast heating system for rapid melting of pelletised wax.

An excellent feature is that the wax dispenser tap is heated via a dedicated heating system which prevents blockages and ensures an even flow of wax on-demand. Delivery of molten wax is via a non-drip lever tap with three settings (off, on-demand pull dispense and continuous flow). The inner tank has a filter screen (0.5 mm mesh) fitted to prevent coarse particles from blocking the delivery tap.

The wax dispenser is controlled via an intuitive digital control panel. The set temperature can be controlled using the up and down arrows on the control panel. The heater indicator will illuminate to show heater activity. The wax dispenser is designed to melt bulk quantities of pelletised wax in as quick a time as possible. When set at 65˚C the instrument will take ~90 minutes to fully melt 5 kg of pelletised wax, with significant quantities of molten wax available for use after as little as 30 minutes.

The equipment is protected by a miniature circuit breaker and a temperature cut out. If an over-temperature condition occurs, the wax dispenser will automatically turn off at 120°C.

Technical Specification

Specification	WD-200
Capacity	7.5 L
Material	Stainless Steel
Temperature	Ambient to 70°C
Dimensions (d x w x h)	500 x 181 x 455mm
Weight	7.0 kg

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Wax Dispenser 230V	85230-07	WD-200	WD-200D	MH8524
Wax Dispenser 120V	85230-08	WD-200	WD-200D-120	MH8524X1

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

AZ-EMEA-LT-00245 1121

Cole-Parmer®
essentials
an antylia scientific company

Find out more at antylia.com