# Short user's manual

## Star-Blender 400 Silent

EU cat. No 710-2143

Manufactured in France

### Warning

Before use, please read these instructions thoroughly. Any damage incurred due to non-respect of these instructions will render the warranty void.
The unit must be placed on its feet, on a stable and horizontal surface. NEVER operate the **Star-Blender 400 Silent** in a different position.
The unit must be used in a well-ventilated room in normal conditions of temperature, pressure and humidity. If it has been stored at a temperature below 10°C, allow the unit to return to room temperature before operating.
The housing should only be removed by a qualified and authorized technician. If the housing must be removed, be sure to disconnect the power cord first.
Each unit is delivered in a specific box. Please keep the box and the inside packing for future use.
For best operation, it is important to maintain and clean regularly the **Star-Blender 400 Silent.** Please refer to the "MAINTENANCE" section for the cleaning procedure and the list of non-recommended products.

### Safety information

To avoid electrical shock:
The unit must always be powered with a grounded lead in compliance with electrical standards.

### Getting started

1. Connect the **Star-Blender 400 Silent** to the power mains (100-240 VAC). Put the switch in position "I".
2. Set the blending time using the time keys ( + **or** − ) or the ∞ for open running.
3. Set the blending speed 1, 2, 3, 4 (4, 6, 8, 10 strokes/s) with the speed keys ( + **or** − )
4. Pull the lever down to open the door. Place the bag with the sample in the unit on the bagrest.
5. Close the door by raising the lever until it clicks shut.
6. The **Star-Blender 400 Silent** starts automatically and the lighting under the front side flashes orange slowly.
7. At the end of the cycle, the **Star-Blender 400 Silent** stops automatically. The lighting turns green, flashes quickly then more slowly. At the end, the unit beeps.
8. Open the door and remove the bag. The **Star-Blender 400 Silent** is ready for the next cycle.
9. After 30 mn of inactivity, **Star-Blender 400 Silent** goes automatically to stand-by mode. To reactivate it, simply press the ⏻ key.

### Maintenance

**Encountering problems**

Do not use sodium hypochlorite (bleach) or quaternary ammonium to clean the unit.
It is important to keep the unit clean, in particular the blending chamber and the door.

**Daily maintenance of the housing and blending chamber**

The housing of the blender can be cleaned with a wet cloth or a product to clean stainless steel. The inside of the unit, the paddles and the drip tray can be cleaned and disinfected with 70° ethanol. Make sure that the mains switch is in position "0".
To clean the blending chamber:
1. Unplug the power cord.
2. Put the unit on the edge of the table.
3. Open the door flat.
4. Remove the drip tray and clean it.
5. Click the paddles and remove them.
6. Clean the blending chamber.
7. Reposition the paddles by clicking them and put back the drip tray in the blending chamber.

8. Lift up the door and plug the power cord.
9. Position the switch at the back of the unit on "I".
10. Close the door. Test the good functioning of the unit.

**Cleaning of door seals**
Regularly lubricate the door seals with a lubricant containing silicone.

**For more information, please consult the full manual"**