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Contents	Page
Safety information	3
German	4
English	7
wiring diagramm	10
spare parts	11

Technische Daten / Technical Data		**PZ 44**
Abmessung Platte / Plate size	mm	440 x 290
Temp.-Vorgabe / -Anzeige Range of temp.	°C	20..450
Schaltdifferenz / Const. Temp.	K	± 0,5
Leistung / Power	W	3300
Gewicht	kg	26 kg

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Safety information

- Read operating instruction before use.

- Always make sure the voltage specified on the type plate corresponds to your local mains voltage.

- use the device only with the protective grounding conductor connected.

- Never use the device in explosion-prone environments.

- Do not expose the device to high air humidity.

- Never insert objects through slots and vent holes.

- Have repairs carried out only by qualified personnel.

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The appliance is constructed with a cold casing and can be placed on heat sensitive surfaces. Height adjustment for stable, horizontal placement can be achieved by rotating the four appliance feet.

1. OPERATION: AT FRONT

The appliance is switched on with the rocker switch (illuminated when in operation). The digital panel indicates the temperature after a few seconds.
Keep the SET key pressed. Use the buttons UP and DOWN to select the desired temperature. You can increase the temperature by pressing the UP key and accordingly decrease it by pressing the DOWN key. You only need to press the SET key to control the set temperature value at any time.

2. SELECTION OF HEATING CAPACITY

The appliance has a wide control range of 20-450°C and is therefore fitted with a power controlling unit by means of which the wattage can again be set from 10-100%. This control is situated on the right side of the appliance.
Required capacity is dependent on required warmth, temperature and heat loss. A higher capacity than is absolutely necessary for required warmth should not be set.
The blinking light-emitting diode indicates the electronic controlling function. Illumination denotes power supply to the heater according to requirements. If the capacity is set at 100%, the diode illuminates continuously during the heating phase. The red light diode at the lower right edge of the last figure indicates that heat has fallen below the required value and the appliance is in a heating phase.

3. OPERATION: RIGHT SIDE

The following wattages can be set with the rotating knob: 825 and 3300 watts heat the entire surface of the plate. At 0 the heater is disconnected from the power supply, for instance to observe cooling down.

SPECIAL SETTING 1650 WATTS

If this wattage is set, the right half of the plate only is heated. The temperature falls away from the right to the left half of the plate by creepage. The temperature is measured by a built-in thermostat on the right side of the plate which permits function at all three settings. Depending on the digital panel temperature setting, a greater or lesser temperature gradient occurs from the right to the left side. Since with this special setting several hundred varying temperature values are involved, in general the temperature gradient cannot be given. Rule of thumb

estimates are based on a temperature gradient of circa 40% from the right plate half to the lef t edge.

4. Contacter - Outlet

This appliance can either be used with a normal electronic contact thermometer (e.g. EBRO GFX 393, Ikatron ETS 3, etc.) or a mercury-contact thermometer.

If a contact thermometer is not being used then the Dummy Connector (which is included) must be plugged in.

The Dummy Connector acts as a safety contact between 3 and 5) This impedes further heating should the plug accidentally be removed while being used with a contact thermometer.

When using the hotplate with a contact thermometer, the hotplate temperature should be accordingly higher adjusted. Should the required temperature shown on the contact-thermometer not be reached then this is normally due to the fact that the hotplate temperature has been set too low. Please ensure that the temperature-limiter is also accordingly set at a higher temperature.

a. Operating with an electronic contact thermometer

An electronic C-Thermometer can be directly connected to the 5-pole socket at the back of the appliance. The necessary power for the contact thermometer will be supplied over this socket.
The required temperature should now be adjusted on the contact thermometer (see the contact thermometer operating instructions).

b. Operating with a mercury contact thermometer

When operating with a mercury contact thermometer, our cable (Order No.. SK 85) must be connected to the thermometer.
Please make sure that the heating output approximately matches the heat-requirement, this can be adjusted with the power regulator (turning knob 10-100%), High surplus heating output always results in post-heating. The more precise the heating output adjustment is, the more exact is the temperature tolerance.

c Operating with a temperature cut-out, Timeswitch, etc..

Instead of the bridge between 3 and 5, other contact-giving devices may be connected that will switch off the heating when the contact has been opened.
When using a temperature cut-out please make sure that the cut-out has a thermal contact to the medium being heated

Please note that **no** external voltage should be used on the contact socket.

Information: the contact-thermometer socket is fitted with 12 V (galvanically separated).

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Harry Gestigkeit GmbH