Allgemeine Produktdaten

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General Product Information

ratiolab® cuvets *with clearly improved photometric properties, optimized form and the small wall thickness of the cuvets provides increased heat transfer resulting in more constant sample temperatures, very low variation of extinction values, exellent optical transmission rage, cavity-sorted production and measurements even in the UV range, path length 10 mm, outer dimensions* 12,5x12,5x45mm

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
UV cuvette

- ☐ Developed specially for photometric identification of proteines, ssDNA, dsDNA, RNA und oligo-nukleotides in UV range
- ☐ Ideal for measurings at **260 nm, 280 nm and in visible wave length range**
- ☐ Coat thickness 10 mm

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\+ = **Very good chemical resistance**
Continuous exposure (test up to 30 days) = no damage.

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o = **Good to limited chemical resistance**
Continuous exposure (tested within the period from the 7th to the 30th day) = minor damages that are partially reversible (softening, swelling, loss of mechanical strength, discoloration).

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\- = **Low chemical resistance**
Not suitable for permanent exposure to the listed materials. Damage can occur immediately (reduction of mechanical strength, deformation, discoloration, cracking, dissolution).

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- PS Polystyrol / polystyrene
- PMMA Polymethylmethacrylat / polymethylmethacrylate

Chemisches Material / chemical material	PS-Halbmikro/ semi-micro	PS-Makro/ macro	PMMA-Halbmikro semi- micro	PMMA-Makro/ macro	UV-Küvette Halbmikro UV-cuvettes semi-micro	Lösemittel-beständige Küvette solvent resistant cuvettes	UV-Küvette Makro UV-cuvettes macro
Aceton / Acetone	-	-	-	-	+	x	+
Aldehyde	-	-	o	o	o	o	o
Alkohole, aliphatisch / alcohol, aliphatic	+	+	-	-	+	+	+
Ammoniak / ammonia	+	+	+	+	+	+	+
Benzaldehyd / benzaldehyde	-	-	-	-	+	+	+
Butanon / butanone	-	-	-	-	+	+	+
Chloroform / chloroform	-	-	-	-	-	-	-
Dioxan / dioxane	-	-	-	-	+	+	+
DMF	-	-	-	-	+	+	+
Essigsäure / acetic acid	-	-	-	-	+	+	+

Chemisches Material / chemical material	PS-Halbmikro/ semi-micro	PS-Makro/ macro	PMMA-Halbmikro semi- micro	PMMA-Makro/ macro	UV-Küvette Halbmikro UV-cuvettes semi-micro	Lösemittel-beständige Küvette solvent resistant cuvettes	UV-Küvette Makro UV-cuvettes macro
(100%)							
Ester / ester	-	-	o	o	o	o	o
Ethylacetat / ethyl acetate	-	-	-	-	+	+	+
Ether (offenkettige) / ether (open-chain)	-	-	-	-	-	-	-
Fluss-Säure / hydro-fluoric acid (10%)	+	+	+	+	+	+	+
Hexan / hexane	-	-	+	+	-	+	-
Isopropanol / isopropanol	+	+	+	+	+	+	+
Kohlenwasserstoffe-aliphatisch (Methan, Propan, Ethan) / hydrocarbons - aliphatic (methane, propane, ethane	-	-	+	+	+	+	+
Kohlenwasserstoffe-aromatisch (Benzol) / hydrocarbons – aromatic (benzol)	-	-	-	-	-	-	-
Kohlenwasserstoffe-halogene (Fluor, Chrom, Brom) / hydrocarbons – halogene (fluor, chrome, bromine)	-	-	-	-	-	-	-
Natronlauge / lye of soda	+	+	+	+	+	+	+
Oxidierende Säuren, Oxidationsmittel / oxidising acids	-	-	-	-	+	+	+
Salpetersäure / azotic acid (65%)	-	-	-	-	+	+	+
Salzsäure / hydrochloric acid (36%)	+	+	-	-	+	+	+

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The above mentioned data was proved by criteria stabilities, thermic stabilities and physical characteristics and summarised according to specialist literature and the information given by the basic material manufacturers. A careful evaluation of application characteristics by the applicator ranks first.

The careful posed recommendations of specialist literature respectively basic material manufacturers shall inform and consult – though the aptitude test by the applicator under respective application conditions can not be replaced.

Terms, application fiels, dimensions according to DIN 58 963 "Cuvettes for photometric measurings – single right-angle cuvettes made of plastics, dimensions, requirements"