VWR®
We Enable Science

Instruction Manual
Mini Shakers
Microplate Shakers

EN - English 1
FR - Français 14
ES - Español 28
IT - Italiano 42
DE - Deutsch 56
PT - Português 70
NL - Nederlands 84
NO - Norsk 88
DA - Dansk 92
SV - Svenska 96
FI - Suomi 100
HU - Magyar 104
PL - Polski 108

European Catalogue Numbers:

Mini Shaker:
Euro plug: 444-0268
UK plug: 444-7093
Swiss plug: 444-0269

Mini Shaker 15:
Euro plug: 444-7106
UK plug: 444-7107
Swiss plug: 444-7108

Microplate Shaker:
Euro plug: 444-0270
UK plug: 444-7094
Swiss plug: 444-0271

Revision 5
15/04/14

CE

Table of Contents

Package Contents	1
Warranty	1
Installation	2
Maintenance & Servicing	2
Intended Use	2
Equipment Disposal	2
Environmental Conditions	2
Safety Instructions	3
Warning	3
Standards & Regulations	3
Specifications - Mini Shaker	4-5
Specifications - Microplate Shaker	6-7
Control Panel	8
Operating Instructions	9-10
Technical Service	11
Troubleshooting	11
Accessories	12-13
EC Declaration of Conformity	112-116

Package Contents

Mini Shaker or Microplate Shaker
Non-skid rubber mat (Mini Shaker only)
234cm detachable power cord
Instruction manual

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Installation

Upon receiving the VWR Mini/Microplate Shaker, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the Mini/Microplate Shaker on a level bench or table, away from explosive vapors. Secure to an immovable work surface by pressing down on the four (4) corners of the unit, creating a strong suction to the work surface (**DO NOT** place on a bench mat). Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit. Always place the unit on a sturdy work surface.

The Mini/Microplate Shaker is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The 230V unit plugs into a 230 volt, 50/60 Hz source.

Maintenance & Servicing

The Mini/Microplate Shaker is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. However at least every three (3) months you should:

- Unplug the unit.
- Remove any accumulated dirt from the base and tray.
- Check all accessible items to make sure they are properly tightened.

The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. **DO NOT** use a cleaning agent or solvent on the front panel which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

Intended Use

The Mini/Microplate Shaker is intended for use in research laboratories and for general laboratory use.

Equipment Disposal

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Environmental Conditions

Operating Conditions: Indoor use only.

* For use in CO2 environments, incubators or cold rooms.
* Temperature: -10 to 60°C
Humidity: maximum 80% relative humidity, non-condensing
Altitude: 0 to 2000 M above sea level

Non-Operating Storage:

* Temperature: -20 to 65°C
Humidity: maximum 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

* **Avoid cold starts:** Unit is not designed to start after being in a cold room environment. Bring unit into cold room from a room temperature environment, operate and remove unit from cold room as soon as operation is complete.

Safety Instructions

Please read the entire instruction manual before operating the Mini/Microplate Shaker.

WARNING! DO NOT use the Mini/Microplate Shaker in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the tray.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly. **DO NOT** immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground - Protective Conductor Terminal

Alternating Current

Standards & Regulations

VWR international hereby declares under its sole responsibility that the products comply with the requirements of the following directives and associated standards:

Associated EU directives:

EMC directive	2004/108/EC
LVD directive	2006/95/EC
ROHS directive	2002/95/EC
WEEE directive	2002/96/EC

Safety standards:

IEC 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use. Part I: General Requirements.

CAN/CSA C22.2 No. 61010-2-051-04
CAN/CSA C22.2 No. 0-M91
CAN/CSA C22.2 No. 61010-1-04
UL 61010-1

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
Mini Shaker Set-Up


The Mini Shaker is supplied with a tray designed to hold a variety of accessories.

1. Flat containers can be shaken by placing them on the non-skid mat provided with the unit.
2. The tray also has mounting holes ready for use with flask/media bottle clamps or test tube racks. See below for tray configurations.


Mini Shaker Tray Configurations


(35) 10mL Erlenmeyer Flask Clamps


(20) 25mL Erlenmeyer Flask Clamps


(15) 50mL Erlenmeyer Flask Clamps


(12) 125mL Erlenmeyer Flask Clamps


(6) 250mL Erlenmeyer Flask Clamps

(4) 500mL Erlenmeyer Flask Clamps

(3) 500mL Media Bottle Clamps

(2) Test Tube Racks

MICROPLATE SHAKER SPECIFICATIONS

Microplate Shaker with plasticware and 1.5 to 2mL Micro-Tube Rack


Overall dimensions (L x W x H):	43,2 x 27,9 x 10,2cm
Tray dimensions (L x W):	27,9 x 19,7cm
Electrical (50/60 Hz):	230 volts, 5 amps, 25 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Speed range:	100 to 900rpm
Speed accuracy:	±2%
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	3mm
Capacity:	4 microplates or 2 Micro-Tube Racks
Controls:	see page 8
Tray material:	aluminum
Ship weight:	11,4kg
European Catalogue Numbers:	Euro plug: 444-0270 UK plug: 444-7094 Swiss plug: 444-0271

Microplate Shaker Set-Up


The Microplate Shaker is designed to hold two (2) or four (4) microplates, or two (2) Micro-Tube Racks.

1. Place two (2) microplates or deep well blocks diagonally on the tray, or place four (4) microplates or deep well blocks on the tray. The plates/blocks do not have to be filled.
2. Place the corner of the plate/block under the spring located at each corner of the tray.
3. Slide plate/block into place. You are ready to use.


2.


3.


4.


Microplate Shaker Tray Configurations

(2) microplates


(4) microplates

(2) Micro-Tube Racks

NOTE: The tray on the Microplate Shaker is not designed to hold flask clamps.

Control Panel

The front panel of the Mini/Microplate Shaker contains all the controls and displays needed to operate the unit.

A. Standby button/standby indicator light: The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to activate the speed and time functions. The standby indicator light will shut off and the speed display and time displays will illuminate. Press the standby button again and the unit will once again be in standby mode.

B. Speed display: Displays the speed of the shaker. **C.** Up/down arrows for set-point control. On/off button starts/stops shaking function. **D.** The speed indicator light will be illuminated when the unit is shaking.

E. Time display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9999. **F.** Up/down arrows for set-point control. On/off button starts/stops the time function.

OPERATING INSTRUCTIONS

The Mini/Microplate Shakers have been designed for the speed and time functions to work independently of one another. The speed can be reset without resetting the timer and the timer can be stopped and started without interrupting the shaking function.

1. **Getting ready:**
 a. Plug the power cord into a properly grounded outlet. The standby indicator light will illuminate, verifying power to the unit.
 b. Press the standby button to move the unit from standby mode. The standby indicator light will turn off and the speed and time displays will illuminate, displaying the previously used settings.
2. **Setting speed:**
 a. Press the up/down arrows below the speed display until you reach the desired speed. When you release the button, the display will blink off and then on indicating the new set speed has been accepted.
 b. Press the on/off button to start the shaking function. The indicator light below the speed display will illuminate to indicate the shaking function is in use and remain lit until shaking has ceased. The microprocessor controlled ramping feature slowly increases speed until the set-point is reached which helps to avoid splashing, and provides excellent low end control.
 c. Speed adjustments can be made without interrupting shaking by using the up/down arrows below the speed display. After the change has been made and you release the button, the display will blink off and then on indicating the new set speed has been accepted.
 d. To stop the shaking function, press the on/off button below the speed display. The speed indicator light will turn off.
3. **Setting time to zero (0:00) and continuous mode:** Accumulated time.
 a. Press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time.
 b. Simultaneously press both the up and the down arrows, the display will indicate zero (0:00). The unit time is now set to zero (0:00) minutes. Alternately, you can use the up/down arrows to get to zero (0:00).
 c. Press the on/off button below the time display. The display will indicate accumulated time. The up/down arrows will become inactive. To stop timer, press the on/off button again. **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function.
 d. To reset, press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time, which was zero (0:00).
4. **Setting timed mode:** Programmed time.
 a. Press the up/down arrows below the time display until you reach the desired time.
 b. Start this function by pressing the on/off button below the time display. The unit will run for the selected time, the up/down arrows will become inactive while the timer is running. The unit will stop shaking when time display reaches zero (0:00). Four (4) audible beeps will indicate the count down function is complete. The time display will default back to the set time. To repeat for the same time, simply press the on/off button again.
 c. To interrupt an automatic timing cycle before it is completed, press the on/off button below the time display. The display will flash off and on to indicate the time function is on “hold”. **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function. Restart the timer by pressing the on/off button below the time display. Unit will continue counting down to zero (0:00). When the display reaches zero (0:00), you will hear the four (4) audible beeps that indicate the count down function is complete and the shaking function will cease.

OPERATING INSTRUCTIONS (CONT'D)

5. **Turning unit off:**

 a. To turn the unit off, press the standby button. The speed and time displays will be blank, the standby indicator light will illuminate. The Mini/Microplate Shaker should be kept in standby mode when not in use. To completely cut off power to the unit, disconnect the power cord form the unit or unplug from the wall outlet.

OPERATING TIPS

As a safety feature, a built-in program will shut power off to the motor if the tray is prevented from rotating, or the unit is overloaded beyond its recommended weight capacity.

Built-in memory maintains the last used speed and time settings during a power interruption.

Technical Service

For more information or technical service contact your local VWR representative or visit www.vwr.com.

Web resources: Visit VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional products information and special offers

Troubleshooting

To clear error press the standby button

Problem	Cause	Solution
Unit fails to power on	Missing or blown fuse	Add or replace fuse as necessary. If problem persists, please contact your VWR Representative for repair.
Unit is excessively noisy	Rattling or ticking sounds may indicate a loose screw on the tray	Ensure that tray is secured tightly. If problem persists, please contact your VWR Representative for repair.
E3	Drive system failure Loose suction cup foot Ceased bearing Broken drive belt	Remove the mechanical obstruction. In the event a foot (suction cup) has come loose from the bench top, firmly press down on the four (4) corners of the unit creating a strong suction to the work surface (**DO NOT** place on bench mat). If problem persists, please contact your VWR Representative for repair.
E4	Unit overload Or Loose suction cup foot	Ensure the load is evenly distributed and does not exceed the maximum load capacity for the unit. In the event a foot (suction cup) has come loose from the bench top, firmly press down on the four (4) corners of the unit creating a strong suction to the work surface (**DO NOT** place on bench mat). If problem persists, please contact your VWR Representative for repair.

ACCESSORIES

TEST TUBE RACKS


DESCRIPTION	TUBE CAPACITY	PLATFORM CAPACITY	CAT. NO.
10 to 13mm Test Tube Rack	63	2	444-0278
14 to 16mm Test Tube Rack	48	2	444-0279
18 to 20mm Test Tube Rack	35	2	444-0280
22 to 25mm Test Tube Rack	24	2	444-0281
15mL Centrifuge Tube Rack	42	2	444-0283
50mL Centrifuge Tube Rack	12	2	444-0284
1,5 to 2mL Micro-Tube Rack	70	2	444-0276

Installation of Test Tube Racks - Mini Shaker

1. Unit should be in standby mode.
2. Remove non-skid rubber mat.
3. Align the clip on the rack to the right side of the tray.
4. Line up the holes in the test tube rack to the holes in the tray. Using the screws provided, hand tighten into place using a flathead screwdriver. **DO NOT** over tighten.

Installation of 1,5 to 2mL Micro-Tube Rack - Microplate Shaker

1. Unit should be in standby mode.
2. Align the Micro-Tube Rack side-to-side on the tray.
3. Press the Micro-Tube Rack in place, you will hear the Micro-Tube Rack clip into position. No mounting screws are required.

FLASK CLAMPS

DESCRIPTION	PLATFORM CAPACITY	CAT. NO.
10mL Erlenmeyer Flask Clamp	35	444-7040
25mL Erlenmeyer Flask Clamp	20	444-7041
50mL Erlenmeyer Flask Clamp	15	444-7042
125mL Erlenmeyer Flask Clamp	12	444-7043
250mL Erlenmeyer Flask Clamp	6	444-7044
500mL Erlenmeyer Flask Clamp	4	444-7045
500mL Media Bottle Clamp	3	444-7035

Installation of Flask Clamps - Mini Shaker

1. Unit should be in standby mode.
2. Remove non-skid rubber mat.
3. Line up the hole(s) in the flask clamp to the hole(s) in the tray. Using the screw(s) provided, hand tighten into place using a flathead screwdriver. **DO NOT** over tighten.
4. Insert flask/media bottle into clamp and wrap the spring around the neck of the clamp to hold secure.

Not for use with the Microplate Shaker.

ACCESSORIES (CONT'D)

UNIVERSAL HARNESS

DESCRIPTION	MICROPLATE CAPACITY	PLATFORM CAPACITY	CAT. NO.
Universal Harness	2 or 4	1	444-0277

Installation of Universal Harness - Mini Shaker

1. Unit should be in standby mode.
2. Remove non-skid rubber mat.
3. Insert one band through the side of one "U" in the U-shaped cut-outs on the center post. Pull band through to secure to the U-shaped cut-out. Repeat for the three remaining bands.
4. After the bands are installed in the center post, line up the holes in center post with the two holes in middle of the tray. Using the screws provided, hand tighten into place using a flat-head screwdriver. DO NOT over tighten.
5. Place samples(s) on tray.
6. Stretch the harness bands over the sample(s) and then under each of the bent-down corners of the tray.

Not for use with the Microplate Shaker.

UNIVERSAL PLATFORM

The Universal Platform has a non-skid rubber surface with adjustable clamping bars to accommodate various vessel types. Stainless steel construction.

DESCRIPTION	SIZE (L x W)	PLATFORM CAPACITY	CAT. NO.
Universal Platform	29,7 x 22,1cm	1	444-0282

Installation of Universal Platform - Mini Shaker

1. Unit should be in standby mode.
2. Remove non-skid rubber mat.
3. Line up the holes in the Universal Platform to the holes in the tray. Using the screws provided, hand tighten into place using a flathead screwdriver. **DO NOT** over tighten.
4. Add glassware. Adjust rollers by loosening each of the thumbscrews. Slide the padded rollers to the desired position and retighten the thumbscrews.

Not for use with the Microplate Shaker.

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Title of Council Directive on harmonization of the laws of the member states
Directive 2006/95/EC of the European Parliament and of the Council of 12 December 2006 on the harmonisation of the laws of Member States relating to Electrical Equipment designed for use within certain voltage limits
Directive 2004/108/EC of the European Parliament and of the Council of 15 December 2004 on the approximation of the laws of the Member States relating to electromagnetic compatibility

and are in conformity with the following standard(s) and/or other normative document(s):
et sont conformes à la (aux) norme(s) ou autre(s) document(s) normatif(s) :
(den) folgenden Normen oder anderer normativer Dokumente entsprechen:

Reference number	Title of the applied standards:
IEC 61010-1	**Safety requirements for electrical equipmentfor measurement, control and laboratory use – Part 1: General requirements – First Edition**
IEC 61010-2-010	**Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-010: Particular requirements for laboratory equipment forheating of materials**
IEC 61010-2-051	**Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-051: Particular requirements for laboratory equipment for mixing and stirring**
EN 45501; 55022A **EN 61000-3-2; -3-3** **EN 61000-4-2; -4-3; -4-4** **EN 61000-4-5; -4-6; -4-11**	**Electromagnetic compatibility (EMC)**

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Yves Van Damme
European Regulatory Affairs Manager

Frankie Vangeel
Managing Director

VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium
Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

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Ireland / Northern Ireland
VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy
VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway
VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax. 058 32 38 205
E-mail: info@pl.vwr.com

Portugal
VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain
VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden
VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland
VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey
Pro Lab Laboratuar Teknolojileri Ltd.□ti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu I□merkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK
VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

Australia
VWR International, Pty Ltd.
Unit 1/31 Archimedes Place
Murarrie, Queensland 4172
Tel.: 1300 727 696
Fax: 1300 135 123

China
VWR (Shanghai) Co., Ltd
2nd Floor, Building 4,
Lane 998, Halei Rd,
Zhangjiang Hi-tech Park
Shanghai, 201203
Tel.:+86-21-5898 6888
Fax:+86-21-5855 8801
E-mail: info_china@vwr.com

India
VWR Lab Products Private Limited
135/12, Brigade Towers, 2nd Floor
Front wing, Brigade Road,
Bengaluru, India – 560 025
Tel.: +91-80-41117125/26 (Bengaluru)
Tel.: +91-2522-647911/922 (Mumbai)
Fax: +91-80-41117120
E-mail: vwr_india@vwr.com

New Zealand
Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Singapore
VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

GO TO VWR.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS OF YOUR LOCAL VWR DISTRIBUTOR