VWR®
We Enable Science

VWR

Microscope VisiScope® series 200

INSTRUCTION MANUAL

Model	European Catalogue Number
ML214PL	630-3295
ML224PL	630-3294
BL214PL	630-3292
BL224PL	630-3291
TL224PL	630-3293
BL224PL T1	630-3290

Version: 1.0
Issued: 23, 12, 2020

CE

Legal Address of Manufacturer

Europe
VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
http://be.vwr.com

Country of origin: **ITALY**

Table of Contents

1. Warning 3
2. Safety Information 3
3. Package Contents 3
4. Unpacking 3
5. Intended use 3
6. Symbols and conventions 4
7. Product specifications 4
8. Instrument description 5
8.1 ML214PL / ML224PL / BL214PL / BL224PL / TL22PL 5
8.2 BL224PL T1 6
8.3 Assembling the microscope 7
8.3.1 ML214PL/ML224PL/BL214PL/BL224PL/TL22PL 7
8.3.2 BL224PL T1 8
8.4 Polarizing set (optional) 10
9. Use of the microscope 11
9.1 Switching on the microscope 11
9.2 Light intensity adjustment 11
9.3 Coarse focus tension adjustment 11
9.4 Stage 11
9.5 Adjust the interpupillary distance 12
9.6 Diopter adjustment 12
9.7 Use of oil immersion objective 12
9.8 Aperture diaphragm 13
9.9 Use of the polarizer (optional) 13
10. Microphotography 14
10.1 Cameras with projection lens 14
10.2 Reflex camera 14
11. Use of software and digital head 15
12. Micrometric Slide VWR005 15
13. Troubleshooting 16
14. Repair and maintenance 17
15. Technical service 17
16. Warranty 17
17. Compliance with local laws and regulations 17
Disposal 18

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.
We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.
We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.
Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package Contents

DESCRIPTION	QUANTITY
Microscope stand with 4 position nosepiece, stage, condenser	1
Optical head (ML: monocular; BL: binocular; TL: trinocular)	1
Objective 4x	1
Objective 10x	1
Objective 40x	1
Objective 100x (or Objective 60x for ML214PL and BL214PL)	1
Eyepiece WF10x/18mm	1 (ML SERIES) 2 (BL and TL SERIES)
Objective vials	4
Green filter, 30mm diameter	1
Immersion oil (not in ML214PL and BL214PL)	1
Dust cover	1
Power supply output 6V	1

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take care to ensure that the optical items (objectives and eyepieces) do not fall out and get damaged. Using both hands (one around the neck and one around the base), lift the microscope from the container and put it on a stable desk.

5. Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION
This symbol indicates a potential risk and alerts you to proceed with caution


ELECTRICAL SHOCK
This symbol indicates a risk of electrical shock.

7. Product specifications


Model	Head	Eyepiece(s)	Nosepiece	Objectives	Stage	Focusing	Illuminator	Condenser	Power supply
ML224PL	Monocular, 360° rotating, 30° inclined	Widefield 10x/18mm	Quadruple, reversed	N-PLAN 4x/0.10, 10x/0.25, 40x/0.65, 100x/1.25	Double layer, 125x115 mm with X-Y 70x30mm moving range	Coaxial coarse and fine focusing	P-LED^2 system with adjustable intensity	1.25 N.A. Abbe type with adjustable height and iris diaphragm	100-240 Vac 50-60 Hz external power adapter 6Vdc output
ML214PL	Monocular, 360° rotating, 30° inclined	Widefield 10x/18mm	Quadruple, reversed	N-PLAN 4x/0.10, 10x/0.25, 40x/0.65, 60x/0.85	Double layer, 125x115 mm with X-Y 70x30mm moving range	Coaxial coarse and fine focusing	P-LED^2 system with adjustable intensity	1.25 N.A. Abbe type with adjustable height and iris diaphragm	100-240 Vac 50-60 Hz external power adapter 6Vdc output
BL224PL	Binocular 360° rotating, 30° inclined with adjustable interpupillary distance and dioptric compensation	Widefield 10x/18mm	Quadruple, reversed	N-PLAN 4x/0.10, 10x/0.25, 40x/0.65, 100x/1.25	Double layer, 125x115 mm with X-Y 70x30mm moving range	Coaxial coarse and fine focusing	P-LED^2 system with adjustable intensity	1.25 N.A. Abbe type with adjustable height and iris diaphragm	100-240 Vac 50-60 Hz external power adapter 6Vdc output
BL214PL	Binocular 360° rotating, 30° inclined with adjustable interpupillary distance and dioptric compensation	Widefield 10x/18mm	Quadruple, reversed	N-PLAN 4x/0.10, 10x/0.25, 40x/0.65, 60x/0.85	Double layer, 125x115 mm with X-Y 70x30mm moving range	Coaxial coarse and fine focusing	P-LED^2 system with adjustable intensity	1.25 N.A. Abbe type with adjustable height and iris diaphragm	100-240 Vac 50-60 Hz external power adapter 6Vdc output
TL224PL	Trinocular 360° rotating, 30° inclined with adjustable interpupillary distance and dioptric compensation	Wide Field 10x/18mm	Quadruple, reversed	N-PLAN 4x/0.10, 10x/0.25, 40x/0.65, 100x/1.25	Double layer, 125x115 mm with X-Y 70x30mm moving range	Coaxial coarse and fine focusing	P-LED^2 system with adjustable intensity	1.25 N.A. Abbe type with adjustable height and iris diaphragm	100-240 Vac 50-60 Hz external power adapter 6Vdc output
BL224PL T1	Binocular Digital observation head, inclined 30° and rotatable 360°.	Wide Field 10x/18mm	Quadruple, reversed,	N-PLAN 4x/0.10, 10x/0.25, 40x/0.65, 100x/1.25	Double layer, 125x115 mm with X-Y 76x30mm moving range	Coaxial coarse and fine focusing	P-LED^2 system with adjustable intensity	1.25 N.A. Abbe type with adjustable height and iris diaphragm	100-240 Vac 50-60 Hz external power adapter 6Vdc output

8. Instrument description

8.1 ML214PL / ML224PL / BL214PL / BL224PL / TL22PL

8.2 BL224PL T1

8.3 Assembling the microscope

8.3.1 ML214PL/ML224PL/BL214PL/BL224PL/TL22PL

1. Remove the dust cap from the microscope frame and from the bottom of the observation head.
2. Insert the optical head above the stand and tighten the screw. (Fig. 1)

- **Hold the head with one hand during the locking in order to avoid that the head falls.**

Fig. 1

3. Insert both eyepieces into the tubes of the optical head. (Fig. 2)

Fig. 2

4. Insert the power supply jack in the socket placed at the rear side of the microscope. (Fig. 3)

Fig. 3

Only for TL224PL

5. Unscrew the protection cap mounted on the photo port and screw the photo tube. (Fig. 4)

Fig. 4

8.3.2 BL224PL T1

1. Remove the dust cap from the microscope frame and from the bottom of the observation head.
2. Insert the optical head above the stand and tighten the screw. (Fig. 5)

- **Hold the head with one hand during the locking in order to avoid that the head falls.**

Fig. 5

3. Insert both eyepieces into the tubes of the optical head. (Fig. 6)
4. Insert the power supply jack in the socket placed at the rear side of the microscope. (Fig. 3)

Fig. 6

5. Fix the rotating part of the junction using the black wing-nut □. (Fig. 7)


Fig. 7

6. Then hook the Tablet PC onto the 4 screws of the junction and pull toward down to firmly lock the Tablet PC in the holder. (Fig. 8)

- To unlock the Tablet PC proceed with the opposite operation: push toward up and remove it from the holder.

Fig. 8

7. Plug one side of the cable named *CAMERA CONNECTION (USB + OTG)* ② to the digital head and the other side to the Tablet PC using the connector ③. (Fig. 9-10).
8. Plug the cable named *POWER SUPPLY CONNECTION* to the Tablet PC for battery recharge using the connector ④. (Fig. 10)

- **The connector ④ has a label on the back side of the Tablet PC ⑤ (Fig. 11). This connector is used for battery recharge and for camera connection.**
- **The connector ③ is used ONLY for the camera connection.**
- **DO NOT use the connector ③ for battery recharge. This could cause malfunctions / damages.**

- The Tablet has been set with the Rotation function disabled: this prevents any flipping of the Live View in order to get a continuous and as large as possible view of your slide also when the Tablet is removed from the holder.
- To enable this function again: you can activate the Rotation by swiping the screen on his bottom right side and selecting Settings + Screen. Anyway, it's not suggested to activate the function when the camera is in Live View mode as it may give troubles when the camera runs at high resolutions.

Fig. 9


Fig. 10

Fig. 11

8.4 Polarizing set (optional)

1. Place the polarizer on the light exit ① at the base of the microscope. (Fig. 12)

Fig. 12

2. Loosen the head fixing knob ② and remove the head from the microscope frame. (Fig. 13)


Fig. 13

3. Insert the analyzer into the hole inside the frame ③. (Fig. 14)
4. Put back the head into its original position and lock the fixing knob.


Fig. 14

9. Use of the microscope

9.1 Switching on the microscope

Operate on the main switch ① placed in the rear side of the microscope, moving the selector on “I” (Fig. 15)


Fig. 15

9.2 Light intensity adjustment

Operate on the light intensity dial to increase or decrease the illumination intensity. (Fig. 16)

Fig. 16

9.3 Coarse focus tension adjustment

- **Adjust the tension using the provided tool**.

The coarse knob tension is preset in the factory.

1. To modify the tension according to personal’s needs, rotate the ring using the provided tool. (Fig. 17).

- Clockwise rotation increases the tension.
- If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

Fig. 17

9.4 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17 mm. (Fig. 18)

1. Open the spring arm of the slide holder ② and place the slide from the front on the stage.
2. Gently release the spring arm of the slide holder.

- **A sudden release of the spring arm could cause the falling of the slide.**


Fig. 18

9.5 Adjust the interpupillary distance

Except ML series
Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig. 19)

- The graduation on the interpupillary distance indicator ①, pointed by the spot "." on the eyepiece holder, shows the distance between the operator's eyes.

The range of the interpupillary distance is 48- 75 mm.


Fig. 19

9.6 Diopter adjustment

Except ML series
1. Look into the right eyepiece with your right eye only, and focus on the specimen.
2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ② to compensate. (Fig. 20)

- **The adjustment range is ±5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's diopter correction.**

Fig. 20

9.7 Use of oil immersion objective

ML224PL / BL224PL / TL22PL / BL224PL T1
1. Focus the specimen with a low power objective.
2. Lower the stage.
3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 21)

- **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality**.
- To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be circular and bright).
- To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.

4. Insert immersion objective.
5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).

- **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult (even not impossible) due to the presence of an additional thickness on the objective.**

Fig. 21

9.8 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Move the diaphragm ring ① (Fig. 22) on the value corresponding to the objective in use. In this case the optimal setting of the condenser is achieved.

- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 23.

Fig. 22


Fig. 23

9.9 Use of the polarizer (optional)

1. Remove the specimen from the stage.
2. Looking inside the eyepieces, rotate the polarizer until the darkest position is achieved.
3. Once the dark is achieved ("extinction" or "Crossed Nicol" position) it is possible to begin the observation.

10. Microphotography

10.1 Cameras with projection lens

1. Remove dust caps from camera and projection lens.
2. Screw the projection lens to camera thread. (Fig. 24)


Fig. 24

3. Insert the projection lens into the photo tube. (Fig. 25)

Fig. 25

10.2 Reflex camera

1. Screw the “T2” ring (not provided) at the end of the projection lens (M-173), then install the whole onto the reflex camera. (Fig. 26)

Fig. 26

2. Insert the projection lens into the photo tube. (Fig. 27)

Fig. 27


11. Use of software and digital head

The camera inside the digital head is driven by PROVIEW software.

For the instructions about the use of the software, please refer to the specific instruction manual.


12. Micrometric Slide VWR005

Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereomicroscopes)

1 DIV=0.01mm

For biological microscopes calibration

1 DIV=0.1mm

For stereo microscopes calibration

13. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
Visibility is poor. • Image is not good. • Contrast is poor. • Details are indistinct. • Image glares	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or to open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage.
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Mechanical Section:		
The coarse focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric section:		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
IV. Observation tube:		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
V. Microphotography:		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder.	Cover eyepieces and viewfinder with a dark cloth

14. Repair and maintenance

Microscopy environment
This microscope is recommended for use in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

The microscope should always be kept vertical when moving it so that no moving parts, such as the eyepieces, fall out. Never mishandle or impose unnecessary force on the microscope.
Never attempt to service the microscope yourself.
After use, turn down the illumination intensity control and turn the light off. Cover the microscope with included the dust cover, and keep it in a dry and clean place.

Electrical safety precautions

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics
If the optical parts require cleaning first use compressed air.
If that is not sufficient use a soft lint-free cloth with water and a mild detergent.
And as a final option use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room**.

Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
Do not disassemble objectives or eyepieces attempting to clean them.
For the best results, use the VWR cleaning kit.
If you need to send the microscope to manufacturer for maintenance, please use the original packaging.

15. Technical service

Web Resource

Visit the VWR's website at **www.vwr.com** for:
- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. **www.vwr.com.**

16. Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

17. Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.
Instead it is your responsibility to correctly dispose of your equipment the end of its life cycle by handling it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.
For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland
VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy
VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

New Zealand
Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Norway
VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax. 058 32 38 205
E-mail: info@pl.vwr.com

Portugal
VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore
VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain
VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden
VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland
VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey
Pro Lab Laboratuar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu Işmerkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK
VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

GO TO VWR.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS OF YOUR LOCAL VWR DISTRIBUTOR

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