Instruction Manual
Incubating Orbital Shakers

EN - English 1
FR - Français 15
ES - Español 31
IT - Italiano 47
DE - Deutsch 63
PT - Português 79
NL - Nederlands 95
NO - Norsk 99
DA - Dansk 103
SV - Svenska 107
FI - Suomi 111
HU - Magyar 115
PL - Polski 119

European Catalogue Numbers:

Large Capacity Shaker
EU/UK/CH Plug: 444-7084

Model 5000I Incubating Orbital Shaker
EU/UK/CH Plug: 444-7109

Model 5000IR Incubating/Refrigerating Orbital Shaker
EU/UK/CH Plug: 444-7112

Revision 13
23/02/2019

Table of Contents

Package Contents	1
Warranty	1
Installation	2
Maintenance & Servicing	2
Intended Use	2
Equipment Disposal	2
Environmental Conditions	2
Safety Instructions	3
Warning	3
Standards & Regulations	3
Control Panel	4
Specifications	5-7
Operating Instructions	8-11
Technical Service	11
Troubleshooting	12
Accessories	13-14
EC Declaration of Conformity	123-126

Package Contents

Incubating Orbital Shaker
Non-skid rubber mat
234cm detachable power cord
Instruction manual

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations
The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Installation

Upon receiving the VWR Incubating Orbital Shaker, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, it requires two (2) people to lift the Incubating Orbital Shaker, from the bottom, to place on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit. Always place the unit on a sturdy work surface.

The Incubating Orbital Shaker is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The 230V unit plugs into a 230 volt, 50/60 Hz source.

Maintenance & Servicing

The Incubating Orbital Shaker is built for long, trouble-free, dependable service.No lubrication or other technical user maintenance is required. However at least every three (3) months you should:

- Unplug the unit.
- Remove any accumulated dirt from the base and tray.
- Check all accessible items to make sure they are properly tightened.

The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. **DO NOT** use a cleaning agent or solvent on the front panel or lid which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

Intended Use

The Incubating Orbital Shaker is designed for a variety of shaking applications.

Equipment Disposal

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Environmental Conditions

Operating Conditions: Indoor use only.

Temperature:	15 to 32°C
Humidity:	maximum 80% relative humidity, non-condensing
Altitude:	0 to 2000 M above sea level

Non-Operating Storage:

Temperature:	-20 to 65°C
Humidity:	maximum 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

Safety Instructions

Please read the entire instruction manual before operating the Incubating Orbital Shaker.

WARNING! DO NOT use the Incubating Orbital Shaker in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Exercise extreme caution when moving the unit, there is a risk of pinching or crushing fingers.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the tray, front bezel or lid.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly. **DO NOT** immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

CAUTION! The caution hot indicator light warns that the temperature of the air in the chamber is above 40°C. The light will illuminate and remain lit when the temperature of the air in the chamber reaches approximately 40°C. When the heat is turned off, the caution hot indicator light will stay lit until the temperature of the air in the chamber is less than 40°C.

Earth Ground - Protective Conductor Terminal

Alternating Current

Standards & Regulations


VWR international hereby declares under its sole responsibility that the products comply with the requirements of the following directives and associated standards:

Associated EU directives:

EMC directive	2014/30/EU
LVD directive	2014/35/EU
ROHS directive	2011/65/EU
WEEE directive	2012/19/EU

Safety standards:
EN 61010
EN 61010-2-051
EN 61010-2-010
EN 50581

EMC standards:
EN 61326-1

Control Panel

The front panel of the Incubating Orbital Shaker contains all the controls and displays needed to operate the unit.

A. Standby button/standby indicator light: The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to activate the temperature, speed and time functions. The standby indicator light will shut off and the temperature, speed and time displays will illuminate. Press the stand-by button again and the unit will once again be in standby mode.

B. Temperature display: Displays the actual/set-point temperatures in conjunction with the actual/set-point indicator lights. **C.** Up/down arrows for set-point control. On/off button starts/stops the heating function. **D.** The heat indicator light will be illuminated when the unit is heating.

E. Actual indicator light: Illuminates when the temperature displayed is the actual temperature of the air in the chamber.

F. Set-point indicator light: Illuminates when the set-point temperature is displayed.

G. Caution hot indicator light: Illuminates when the air temperature of the chamber is above 40°C (104°F).

H. Speed display: Displays the speed of the shaker. **I.** Up/down arrows for set-point control. On/off button starts/stops shaking function. **J.** The speed indicator light will be illuminated when the unit is shaking.

K. Time display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9999. **L.** Up/down arrows for set-point control. On/off button starts/stops the time function.

Specifications - Large Capacity Incubating Shaker

Overall dimensions (L x W x H):	64,8 x 35,6 x 40,6cm
Interior dimensions (L x W x H):	34 x 30,5 x 24,1cm
Tray dimensions (L x W):	33 x 27,9cm
Electrical (50/60 Hz):	230 volts, 5 amps, 450 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Temperature range:	ambient +5°C to 65°C
Temperature uniformity:	±0,5°C at 37°C
Speed range:	15 to 500rpm
Speed accuracy:	
Above 100rpm	±1% of set speed
Below 100rpm	±1rpm
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	19mm
Maximum weight capacity:	15,9kg
Controls:	see page 4
Tray material:	aluminum
Ship weight:	37,7kg
European Catalogue Number:	EU/UK/CH plug: 444-7084

Specifications - Model 5000I Incubating Shaker

Overall dimensions (L x W x H):	81,5 x 67,6 x 59,7cm
Interior dimensions (L x W x H):	52,3 x 63 x 43,2cm
Tray dimensions (L x W):	45,7 x 45,7cm
Electrical (50/60 Hz):	230 volts, 8 amps, 800 watts
Fuses:	5mm x 20mm, 8 amp quick acting
Temperature range:	ambient +5°C to 65°C
Temperature uniformity:	±0,5°C at 37°C
Speed range:	15 to 500rpm
Speed accuracy:	
Above 100rpm	±1% of set speed
Below 100rpm	±1rpm
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	25mm
Maximum weight capacity:	22,7kg
Controls:	see page 4
Tray material:	aluminum
Ship weight:	97,5kg
European Catalogue Number:	EU/UK/CH plug: 444-7109

Specifications - Model 5000IR Incubating/Refrigerating Shaker

Overall dimensions (L x W x H):	104,4 X 67,6 X 59,7cm
Interior dimensions (L x W x H):	52,3 x 63 x 43,2cm
Tray dimensions (L x W):	45,7 x 45,7cm
Electrical (50/60 Hz):	230 volts, 10 amps, 800 watts
Fuses:	5mm x 20mm, 10 amp quick acting
Temperature range:	ambient -15°C to 65°C
Temperature uniformity:	±0,5°C at 37°C
Speed range:	15 to 500rpm
Speed accuracy:	
Above 100rpm	±1% of set speed
Below 100rpm	±1rpm
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	25mm
Maximum weight capacity:	22,7kg
Controls:	see page 4
Tray material:	aluminum
Ship weight:	117,9kg
European Catalogue Number:	EU/UK/CH plug: 444-7112

5000IR Upper Ambient Control Limit

Please follow the procedure below if setting the 5000IR unit's temperature at/near the room's ambient temperature.

For optimum temperature control in the chamber of the 5000IR unit, the upper ambient control limit is to be set a minimum of 3°C above the measured room temperature.

Example: Measured (ambient) temperature is 22°C, the upper ambient control limit should be set at a minimum of 25°C.

With the uinit at the heat/cool setting as described below, and in cool mode ("cool" on time display), a temperature setting is displayed in the temperature window. The default temperature setting of the 5000IR is 28°C. This is the factory set upper ambient control limit. This setting can be adjusted between 20°C and 32°C by the end user.

- Any unit set point above the upper limit setting will heat only
- Any unit set point below 20°C will cool only
- Any set point between 20°C and the selected upper limit setting, up to 32°C, the chiller and heater will both control chanber temperature

Heat/Cool Setting Procedure to Adjust Upper Ambient Control Limit

1. Put unit in standby mode.
2. Press and hold the temperature down, speed down and time down time arrows simultaneously, now press the standby button. Release the standby button and then release all 3 of the down arrows simultaneously.
3. Time display will not read "cool" and the temperature display will read 28°C. If the time display reads "heat" press the time up/down arrow until the word "cool" appears in the display.
4. Now you can set the upper ambient control limit. Using the temperature up/down arrows adjust this temperature to be 3°C above the measured ambient temperature.
5. Once the desired temperature is selected, press the standby button to return to normal operation.
6. This procedure does not require recalibrating the unit.

Additional 5000IR Notes:

- 5000IR units have a toggle switch on the left side. This switch must be in the "on" position (depress "I") for the unit to function properly.
- For 5000IR units, make sure the chiller's hose is positioned to drain freely (into a sink when possible).
- Do not allow the end of the chiller hose to become submerged.
- If chiller hose is positioned in a drain container, the container should have a volume of at least 8 liters if the unit will be unattended for 48 hours.
- There are two strips of adhesive on the bottom of the rubber mat. To prevent the mat from sliding on the tray during operation, remove the plastic film to expose the adhesive and firmly press the mat down onto the tray to secure.

Operating Instructions

The Incubating Orbital Shaker has been designed for the temperature, speed and time functions to work independently of one another. The temperature and speed can be reset without resetting the timer and the timer can be stopped and started without interrupting the heating and shaking functions.

1. **Getting ready:**
 a. Plug the power cord into a properly grounded outlet. The standby indicator light will illuminate, verifying power to the unit.
 b. Press the standby button to move the unit from standby mode. The standby indicator light will turn off and the temperature, speed and time displays will illuminate, displaying the previously used settings.

 Note: 5000IR units have a toggle switch on the left side. This switch must be in the "on" position (depress "I") for the unit to function properly.

 Note: See the procedure on page 8 for 5000IR units that will have the chamber temperature set at/near the ambient temperature.

2. **Setting temperature:**
 a. Press the up/down arrows below the temperature display until you reach the desired temperature. When you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
 b. Press the on/off button to start the heating function. The indicator light below the temperature display will illuminate to indicate the heating function is in use and remain lit until heating has ceased.
 c. Temperature adjustments can be made without interrupting heating by using the up/down arrows below the temperature display. After the change has been made and you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
 d. To stop the heating function, press the on/off button below the temperature display. The heat indicator light will turn off.

CAUTION HOT indicator:
The caution hot indicator light warns that the temperature of the air in the chamber is above 40°C. The light will illuminate and remain lit when the temperature of the air in the chamber reaches approximately 40°C. When the heat is turned off, the caution hot indicator light will stay lit until the temperature of the air in the chamber is less than 40°C.

3. **Setting speed:**
 a. Press the up/down arrows below the speed display until you reach the desired speed. When you release the button, the display will blink off and then on indicating the new set speed has been accepted.
 b. Press the on/off button to start the shaking function. The indicator light below the speed display will illuminate and blink until the set-point is reached. Once the set-point is reached the light will stop blinking and remain lit until shaking has ceased. The micro-processor controlled ramping feature slowly increases speed until the set-point is reached which helps to avoid splashing, and provides excellent low end control.
 c. Speed adjustments can be made without interrupting shaking by using the up/down arrows below the speed display. After the change has been made and you release the button, the display will blink off and then on indicating the new set speed has been accepted.
 d. To stop the shaking function, press the on/off button below the speed display. The speed indicator light will turn off.

OPERATING INSTRUCTIONS (CONT'D)

4. **Setting time to zero (0:00) and continuous mode:** Accumulated time.
 a. Press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time.
 b. Simultaneously press both the up and the down arrows, the display will indicate zero (0:00). The unit time is now set to zero (0:00) minutes. Alternately, you can use the up/down arrows to get to zero (0:00).
 c. Press the on/off button below the time display. The display will indicate accumulated time. The up/down arrows will become inactive. To stop timer, press the on/off button again. **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function.
 d. To reset, press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time, which was zero (0:00).
5. **Setting timed mode:** Programmed time.
 a. Press the up/down arrows below the time display until you reach the desired time.
 b. Start this function by pressing the on/off button below the time display. The unit will run for the selected time, the up/down arrows will become inactive while the timer is running. The unit will stop shaking when time display reaches zero (0:00). Four (4) audible beeps will indicate the count down function is complete. The time display will default back to the set time. To repeat for the same time, simply press the on/off button again.
 c. To interrupt an automatic timing cycle before it is completed, press the on/off button below the time display. The display will flash off and on to indicate the time function is on "hold". **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function. Restart the timer by pressing the on/off button below the time display. Unit will continue counting down to zero (0:00). When the display reaches zero (0:00), you will hear the four (4) audible beeps that indicate the count down function is complete and the shaking function will cease.
6. **Turning unit off:**
 a. To turn the unit off, press the standby button. The temperature, speed and time displays will be blank, the standby indicator light will illuminate. The Incubating Orbital Shaker should be kept in standby mode when not in use. To completely cut off power to the unit, disconnect the power cord from the unit or unplug from the wall outlet.

OPERATING TIPS

Opening the lid on the Incubating Orbital Shaker will cause the unit to pause shaking and/or heating. Close the lid and the unit will resume shaking and/or heating at current settings.

Centering your sample and even weight distribution on the tray helps with balance and stability.

As a safety feature, a built-in program will shut power off to the motor if the tray is prevented from rotating, or the unit is overloaded beyond its recommended weight capacity.

The shaker will automatically restart after a power interruption. Built-in memory maintains the last used temperature, speed and time settings during a power interruption.

LOAD SENSING FUNCTION

The Incubating Orbital Shaker is equipped with a load sensing function that can be activated by the user. This function provides protection against improper positioning of load and maximum load being exceeded. When activated, the unit will automatically sense improper load conditions and slow to a safe running speed, then display that speed followed by an E04 error message on the speed display. The unit will also beep three (3) times every 60 seconds until the error is reset by pressing the speed on/off button. To activate the load sensing function use the following steps:

1. Place the unit in standby mode.
2. Press and hold the speed on/off button and press the standby button. The unit will beep two (2) times, confirming the function is enabled.

3. To restore normal operation, remove AC power to the unit for ten (10) seconds and then restore.

If the E04 error occurs be sure the load is within the maximum specification and properly balanced (centered on tray) and/or reduce sample size/speed before restarting the unit. If the E04 error occurs due to acceptable sample vibration or another vibration source, the vibration sensing function can be disabled as described above.

BEEPER PREFERENCE (muting audible alarm)

To silence beeper operation (except for error codes), with the unit in standby mode, press and hold the time on/off button and press the standby button. Release the standby button first, and then release the on/off button. To restore normal beeper operation, remove AC power to unit for ten (10) seconds and then restore.

SINGLE POINT CALIBRATION PROCEDURE

This procedure is used to fine tune and calibrate the Incubating Orbital Shaker at a specific temperature setting. This process may be repeated for up to three (3) separate set-points. If a fourth calibration set-point is entered, the first set-point entered will be overwritten.

1. Turn unit on.
2. Set desired temperature.
3. Stabilize one (1) hour or more, measuring the temperature with a temperature probe/thermometer.
4. Press and hold standby button, then press the temperature up arrow once. The unit will beep two (2) times, confirming calibration mode. The display will now be flashing.
5. Press the temperature up/down arrows until the display matches the temperature probe/thermometer.
6. Press standby button to exit calibration mode and return to normal heating.

This process may be repeated at the same set-point, multiple times for fine tuning if desired.

The unit will now use the biased offset for that specific temperature setting and increase or decrease temperature accordingly to bring the temperature to set temperature. The decimal point of the display will flash to indicate a biased offset is being used. All other temperature settings will use the standard internal calibration. This offset will be stored in memory and retained until reset.

To restore unit to factory setting:

Press and hold the standby button while pressing the temperature down arrow once. The reset will be confirmed with two (2) audible beeps. Press the standby button to exit calibration mode and return to normal heating.

Speed Calibration Procedure:

This procedure is used to self calibrate the incubating shakers. The tray should be free of any samples, vessels, and accessories prior to calibrating.

1. Turn unit on. Speed and time displays will be illuminated.
2. Press and hold the standby button and momentarily press the speed on/off button. The speed display should read "CAL".
3. The unit will be run for approximately one (1) minute and automatically calibrate.

RS-232 SERIAL PORT

The RS-232 serial port provides two-way communications for data logging and unit control. **NOTE:** To comply with EMI/EMC requirements, the length of the external cable shall not exceed 3 meters.

TECHNICAL SERVICE

For more information or technical assistance contact your local VWR representative or visit www.vwr.com.

Web resources: Visit VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

TROUBLESHOOTING


To clear error press the standby button

Problem	Cause	Solution
Unit fails to power on	Missing or blown fuse	Check and replace fuse if necessary If problem persists, please contact your VWR representative for repair.
Unit is excessively noisy	Rattling or ticking sounds may indicate a loose screw on the tray	Ensure that tray is secured tightly. If problem persists, please contact your VWR representative for repair.
Unit not shaking at proper speed	-	Perform speed calibration on page 11 If problem persists, please contact your VWR representative for repair.
Unit not heating to proper temperature	-	Perform single point calibration on page 11 If problem persists, please contact your VWR representative for repair.
E1	RTD open or ±1ºC temperature deviation from set-point (after unit has stabilized)	This error connot be fixed by the end user. If problem persists, please contact your VWR representative for repair.
E2	RTD shorted or temperature below 0ºC	This error cannot be fixed by the end user. If problem persists, please contact your VWR representative for repair.
E3	Mechanical obstruction Drive system failure Ceased bearing Drive belt broken	Remove mechanical obstruction. If problem persists, the reason may be the drive system. Please contact your VWR representative for repair.
E4	Improper positioning of load or Maximum load exceeded	Ensure the load is evenly distributed and does not exceed the maximum load capacity for the unit. See “Load Sensing Function” on page 10-11. Verify load has not exceeded the maximum allowable weight. Reduce load if necessary. If error code resumes, contact your VWR representative for repair.
E6	Over temperature error	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E7	Default sensor Error (5000IR units only)	Error will shut unit off if set temperature is 20°C or below. Unit is still operational with a set temperatue above 20°C. Error will be indicated every 5 minutes until repaired. Please contact your VWR representative for repair.

ACCESSORIES

PLATFORMS

Platform can be used for mounting Erlenmeyer flask/media bottle clamps and test tube racks. Slides over-top included tray. Tighten with adjustment screws. Stainless steel construction.

DESCRIPTION	DIMENSIONS (L x W)	CAT. NO.
Platform	33 x 28cm	444-7087
Platform	47,5 x 47,5cm	444-0967

TEST TUBE RACKS - *Requires Platform accessory.*

DESCRIPTION	TUBE CAPACITY	33 x 28CM PLATFORM CAPACITY	47,5 x 47,5 PLATFORM CAPACITY	CAT. NO.
Half Size Test Tube Racks				
10 to 13mm Test Tube Rack	63	2	4	444-0278
14 to 16mm Test Tube Rack	48	2	4	444-0279
18 to 20mm Test Tube Rack	35	2	4	444-0280
22 to 25mm Test Tube Rack	24	2	4	444-0281
15mL Centrifuge Tube Rack	35	2	4	444-0283
50mL Centrifuge Tube Rack	12	2	4	444-0284
1,5 to 2mL Micro-Tube Rack	70	2	4	444-0276
Full Size Stationary Test Tube Racks				
10 to 14mm Test Tube Rack	48	N/A	3	444-7037
16 to 20mm Test Tube Rack	33	N/A	3	444-7038
21 to 25mm Test Tube Rack	21	N/A	3	444-7039
Full Size Pivoting Test Tube Racks				
13mm Test Tube Rack	90	N/A	2	444-2955
16mm Test Tube Rack	60	N/A	2	444-2956
20mm Test Tube Rack	40	N/A	2	444-2957
25mm Test Tube Rack	24	N/A	2	444-2958
30mm Test Tube Rack	21	N/A	2	444-2959

Installation of Half Size Test Tube Racks:

1. Unit should be in standby mode.
2. Remove non-skid rubber mat. Slide Platform over-top included tray.
3. Using the hardware provided, attach brackets to the Platform. Position brackets as shown. The proper spacing between the brackets is 19,7cm.
4. Align the test tube rack side-to-side with the brackets on the Platform.
5. Press the test tube rack in place, you will hear the test tube rack clip into position. No mounting screws are required.

Installation of Full Size Pivoting Test Tube Racks *(Model 5000I/5000IR Only)***:**

1. Unit should be in standby mode.
2. Line up the 4 openings on the bottom of the outer rack with the four compatible openings on the Universal Platform. Secure with the four screws provided.
3. If a tilt angle is desired, loosen wing knobs, tilt inner holder and tighten knobs.
4. To maximize the platforms capacity for holding racks in the stationary position, remove the outer holder and secure the inner holder to the Universal Platform using the four openings on the bottom of the inner holder.


Installation of Full Size Stationary Test Tube Racks *(Model 5000I/5000IR Only)***:**

1. Line up the two openings on the sides of the test tube rack with the two compatible openings on the Universal Platform.
2. Place the washer on top of the opening on the test tube rack, place the clamp knob through the washer and rack and hand tighten the clamp knob to the platform.

Accessories

PLATFORMS

Platform can be used for mounting Erlenmeyer flask/media bottle clamps and test tube racks. Slides over-top included tray. Tighten with adjustment screws. Stainless steel construction.

Description	Dimensions (L x W)	Cat. No.
Platform	33 x 28cm	444-7087
Platform	47,5 x 47,5cm	444-0967

FLASK CLAMPS - *Requires Platform accessory*

Description	33 X 28cm Platform Capacity	47,5 x 47,5 Platform Capacity	Cat. No.
Stainless Steel Flask Clamps			
10mL Erlenmeyer Flask Clamp	60	113	444-7040
25mL Erlenmeyer Flask Clamp	25	64	444-7041
50mL Erlenmeyer Flask Clamp	13	32	444-7042
125mL Erlenmeyer Flask Clamp	10	20	444-7043
250mL Erlenmeyer Flask Clamp	9	20	444-7044
500mL Erlenmeyer Flask Clamp	7	13	444-7045
1L Erlenmeyer Flask Clamp	4	8	444-7030
2L Erlenmeyer Flask Clamp	N/A	5	444-7031
2.8L Erlenmeyer Flask Clamp	N/A	2	444-7032
4L Erlenmeyer Flask Clamp	N/A	4	444-7033
6L Erlenmeyer Flask Clamp	N/A	2	444-7034

Description	33 X 28cm Platform Capacity	47,5 x 47,5 Platform Capacity	Cat. No.
Stainless Steel Media Bottle Clamps			
500mL Media Bottle Clamp	5	16	444-7035
1L Media Bottle Clamp	2	10	444-7036
PVC Flask Clamps			
125mL Erlenmeyer Flask Clamp	10	20	444-2960
250mL Erlenmeyer Flask Clamp	8	18	444-2961
500mL Erlenmeyer Flask Clamp	5	12	444-2962
1L Erlenmeyer Flask Clamp	2	8	444-2963
2L Erlenmeyer Flask Clamp	N/A	4	444-2964

Installation of Flask Clamps:

1. Unit should be in standby mode.
2. Remove non-skid rubber mat. Slide Platform over-top included tray.
3. Line up the hole(s) in the flask clamp to the hole(s) in the Platform. Using the screw(s) provided, hand tighten into place using a flathead screwdriver. **DO NOT** over tighten.
4. Insert flask/media bottle into clamp and wrap the spring around the neck of the clamp to hold secure.

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The undersigned states that this declaration is issued under the sole responsibility of the manufacterer, VWR international bvba, and that the object(s) of the declaration described in annex 1 are in conformity with the relevant Union harmonization legislation.
In the event of unauthorized modification of any of the products listed in the annex 1, this declaration becomes invalid.

Object Name: VWR Incubating Orbital Shakers

Relevant EC Directives:	Low Voltage Directive	2014/35/EU
	EMC Directive	2014/30/EU
	RoHS	2011/65/EU
Harmonized Standards:	EN 61326-1:2013	
	EN 61010-1:2001	
	EN 61010-2-010:2003	
	EN 61010-2-051:2003	
	EN 50581:2012	

Document Number: 14871.2017.02

Place & date of issue: Leuven, 3/03/2017

Signed for and on behalf of VWR International bvba:

Yves Van Damme
European Regulatory Affairs Manager

Kris Nijs
Managing Director

VWR International bvba
Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven
http://be.vwr.com

Annex 1

Cat. No.	Description
444-7084	VWR Incubating Orbital Shaker, Digital, Large Capacity, EU-Plug
444-7085	VWR Incubating Orbital Shaker, Digital, Large Capacity, UK-Plug
444-7086	VWR Incubating Orbital Shaker, Digital, Large Capacity, CH-Plug

VWR
We Enable Science

VWR International bvba
Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven
http://be.vwr.com

EU Declaration of Conformity

Manufacturer: VWR International bvba, Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven

The undersigned states that this declaration is issued under the sole responsibility of the manufacterer, VWR International bvba, and that the object(s) of the declaration described in annex 1 are in conformity with the relevant Union harmonization legislation.
In the event of unauthorized modification of any of the products listed in the annex 1, this declaration becomes invalid.

Object Name: VWR Incubating Orbital Shakers 5000

Relevant EC Directives:	Low Voltage Directive	2014/35/EU
	EMC Directive	2014/30/EU
	RoHS	2011/65/EU
Harmonized Standards:	EN 61326-1:2013	
	EN 61010-1:2001	
	EN 61010-2-010:2003	
	EN 61010-2-051:2003	
	EN 50581:2012	

Document Number: 14871.2017.05

Place & date of issue: Leuven, 3/03/2017

Signed for and on behalf of VWR International bvba:

Yves Van Damme
European Regulatory Affairs Manager

Kris Nijs
Managing Director

VWR International bvba
Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven
http://be.vwr.com

Annex 1

Cat. No.	Description
444-7109	VWR Incubating Orbital Shaker Large Capacity, 5000I
444-7112	VWR Incubating Refrigerating Orbital Shaker Large Capacity, 5000IR

Austria
VWR International GmbH
Graumanngasse 7
1150 Vienna
Tel.: +43 01 97 002 0
Email: info.at@vwr.com

Belgium
VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 016 385 011
Email: vwr.be@vwr.com

China
VWR International China Co., Ltd
Shanghai Branch
Room 256, No. 3058 Pusan Road
Pudong New District
Shanghai 200123
Tel.:+86-21-5898 6888
Fax:+86-21-5855 8801
Email: info_china@vwr.com

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India
VWR Lab Products Private Limited
No.139. BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India – 560058
Tel.: +91-80-28078400
Email: vwr_india@vwr.com

Ireland / Northern Ireland
VWR International Ltd / VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 01 88 22 222
Email sales.ie@vwr.com

Italy
VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02-3320311
Email: info.it@vwr.com

The Netherlands
VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 020 4808 400
Email: info.nl@vwr.com

Norway
VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: +47 22 90 00 00
Email: info.no@vwr.com

Poland
VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 058 32 38 200
Email: info.pl@vwr.com

Portugal
VWR International –
Material de Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, nº 6
2610-088 Alfragide
Tel.: +351 21 3600 770
Email: info.pt@vwr.com

Singapore
VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Email: sales.sg@vwr.com

Spain
VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
Email: info.es@vwr.com

Sweden
VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: +46 08 621 34 00
Email: info.se@vwr.com

Switzerland
VWR International AG
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 044 745 13 13
Email: info.ch@vwr.com

UK
VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
Email: uksales@vwr.com

Legal Manufacturer:
VWR International bvba • Researchpark Haasrode 2020 • Geldenaaksebaan 464 • B-3001 Leuven • + 32 16 385011 • http://be.vwr.com
715085-00 (Rev 13)