MD

avantor™
delivered by VWR™

EN Q Path® NOVABRUSH® DETACHABLE BRUSHES
Instruction for use

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INTENDED PURPOSE

The Q Path NovaBrush family of devices are single-use cervical brushes intended for the collection of cells from the endocervix (Q Path® NovaBrush® 1) or from the ecto- and endocervix (Q Path® NovaBrush® 2, 3 and 4).
These cells are subsequently fixed and assessed for morphological changes associated with cervical cancer using conventional or liquid-based cytology or subjected to testing for human papillomavirus.
For professional use.

INDICATIONS AND CONTRAINDICATIONS

Indications

The Q Path® NovaBrush® 1, 2, 3 and 4 are indicated for the collection of endo- and ectocervical cells for subsequent conventional cytology, liquid-based cytology, and human papillomavirus testing.

Contraindications

Q Path® NovaBrush® family is contraindicated for use on pregnant women after 10 weeks of pregnancy.

Potential side effects

No potential side effects.

Warnings and precautions

- Brushes are intended for single use only. The re-use can be source of viral, bacterial, and fungal contaminations. Furthermore, the re-use of a contaminated brush carries the risk of a false diagnosis.
- Factors that may decrease the quality of the smear:
 - Menstruations, important losses of blood.
 - Infection and\or vaginal inflammation.
 - Severe genital atrophy (menopause).
 - Preliminary physical action or chemical irritation as:
 - Vaginal touch before the smear or the sexual intercourse within hours which precedent
 - Disinfecting liquid, lubricant.
 - Vaginal ovules (within the limits of 48 h).
 - Colposcopy with acetic acid (within the limits of 24 h).
 - Smear or cervical surgery in the 3 months.
 - Radiotherapy
- Check that brushes are clean and without debris before use.
- Collect the specimen according to the protocol.

Patient Population

The Q Path® NovaBrush® Family is intended to be used for adult female patients, including pregnant women before 10 weeks of pregnancy.

User Profile

The Q Path® NovaBrush® Family is intended for professional use by clinicians or other adequately qualified healthcare professionals.

Instruction for use

MD For professional use only

English

Q Path® NOVABRUSH® DETACHABLE BRUSHES

Use environment

The Q Path® NovaBrush® Family is intended to be used in professional healthcare environments, for example hospitals, doctor's offices and diagnostic laboratories.

COMPOSITION

Polypropylene and polyethylene

STORAGE AND SHELF LIFE

Brushes Q Path® NovaBrush® must be stored in a dry location, away from any dust.
Do not use if the package is damaged.
Do not use after the expiry date indicated on the package.


ADDITIONAL SPECIAL EQUIPMENT

- Speculum.
- Vial of Q Path® Easyfix® cellular conservative.
- Slide and cytology fixative Q Path® Labofix®.
- Gloves.

PROTOCOL

Liquid-based cytology

1. Insert the NovaBrush® into the cervix for a specimen sampling on the exo-endocervical junction.
2. Execute the sampling by a gentle pressure and a rotation of 5 times in a clockwise direction and then pull out the brush with precaution avoiding the contact with the speculum.
3. Insert the NovaBrush® in the prefilled vial with cytology fixative, shake the NovaBrush®.
4. To release the brush, use a gentle pressure on the handle downwards, holding firmly the vial: the brush will release easily.
5. Remove the handle, close the vial, identify the patient.

Conventional cytology

1. Insert the NovaBrush® into the cervix for a specimen sampling on the exo-endocervical junction.
2. Execute the sampling by a gentle pressure and a rotation of 5 times in a clockwise direction and then pull out the brush with precaution avoiding the contact with the speculum.
3. Spread out the sample on the slide, taking care to cover the largest surface area possible.
4. Fix the sample using a fixative e.g. Q Path® Labofix®.
5. Allow to dry and clearly identify the slide.
6. Discard the brush in a disposal container provided for this purpose.
7. Ensure that the sample reaches the cytopathology laboratory as quickly as possible.

METHODOLOGY

Brushes Q Path® NovaBrush® allow to do a representative sampling of cells of the cervix for a cytological examination.

Instruction for use

MD For professional use only

English

Q Path® NOVABRUSH® DETACHABLE BRUSHES

Exo- and endocervicals cells are simultaneously taken: the long strands collect the endocervical cells and the courts, the exocervical cells.
An optimal sampling with Q Path® NovaBrush® 2 and the NovaBrush® 3 requires 5 clockwise rotations of the brush between the thumb and the index finger.

DISPOSAL

Used brush can be a source of infection or microbial hazard. Dispose of the product in accordance with current applicable legislation.

NOTE FOR THE USER

Any serious incident that has occurred in relation to this product must be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or the patient is established. Please contact your local VWR partner through the vwr.com website.

We pay the utmost attention on the quality of our products. In case of any malfunction and dissatisfaction please use the complaints form on the contact page of your local vwr.com website and indicate the reference and batch number of the product concerned.

This instruction for use was established in 2021-07 and cancels all previous instructions. The information is based on our knowledge and our experience this day. The attention of the users is drawn on the eventual risks when the product is used for a different use than it is dedicated. The user is not exempt to know and implement all the texts regulating his activity. He will take under his responsibility the precautions bound in the use that he does of the product.

SYMBOLS USED (EN ISO 15223-1:2021)

Symbol	Symbol Title and Description
MD	Medical device: Indicates the item is a medical device.
REF	Catalogue number: Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Do not re-use: Indicates a medical device that is intended for one single use only.
LOT	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Caution: To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.

Symbol	Symbol Title and Description
	Manufacturer: Indicates the medical device manufacturer.
	Consult instructions for use or consult electronic instructions for use: Indicates the need for the user to consult the instructions for use.
	Keep dry: Indicates a medical device that needs to be protected from moisture.
	Use-by date: Indicates the date after which the medical device is not to be used.
	Do not use if package is damaged and consult instructions for use: Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.

VWR International S.A.S.
Immeuble Estréo, 1-3 rue d´Aurion
93110 Rosny-sous-Bois, France
+33 (0) 145 14 85 00

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