

Certificate of Conformity

Material Number: 85111.003
Material description: Melting Point Standard Set;
Vanillin (81 – 83°C)
Sulphanilamide (164 - 166°C)
Caffeine (235 – 238°C)

Lot Number: 240030
Use by *: 2026-Nov-28
Date of Quality Control: 2024-Jan-30



This Melting Point Standard Set consists of the following Components:

Material No.	Material Product Description	Lot No
85103.001	Melting Point Standard Vanillin 81 to 83°C	240018
85107.001	Melting Point Standard Sulphanilamide 164 to 166°C	240022
85108.001	Melting Point Standard Caffeine 235 to 238°C	240018

The Certificates of Analysis for the individual Melting Point Standard Set components are available under the above referenced Material No. and Lot No.

This certificate must not be reproduced except in full.

Signature

We certify that this batch conforms to the specifications listed above. This document has been produced electronically and is valid without a signature

Signed on behalf of VWR: 31st January 2024

A.Hoszowska

QC Technician

Other Information

* The use by date is the date to which VWR guarantees the above specification will be met for the product in an unopened container stored under our recommended conditions.