

PYREX® Reusable Media Bottles

PYREX reusable laboratory bottles are designed as storage vessels for reagents, sterile tissue culture media and components, biological fluids, aqueous and nonaqueous solutions. Choose Low Actinic bottles for light sensitive materials or PYREXPLUS® bottles for the added protection of a PVC coating.

Materials

PYREX bottles are made from 7740 glass type, a chemically resistant, low expansion, borosilicate glass. The bottles are supplied with clear Polypropylene pouring rings, and orange caps are fully autoclavable. A wide range of optional closures are available for the bottles including:

- ▶ Plug seal, GL45 color-coded caps to aid in sorting and identification
- ▶ High temperature, PBT caps with PTFE coated liners (Cat. No. **16157-464**), ETFE copolymer rings (Cat. No. **16157-497**), and Silicone septa (Cat. No. **11310-694**) are suitable for dry heat sterilization (180°C).
- ▶ Vented, GL45 caps with 0.22 µm PTFE hydrophobic membranes (Cat. No. **16157-412**)
- ▶ Open Top, PBT caps in GL25, GL32, and GL45 sizes
- ▶ Silicone or PTFE-faced Silicone Septa for the open-top caps

PBT = Polybutylene-teraphthalate; PTFE = Polytetrafluoro ethylene; ETFE = Ethylenetetrafluoro-ethylene



Ordering Information

1395 PYREX Laboratory Storage Bottles

VWR Cat. No.	Corning Cat. No.	Capacity (mL)	Approx. Diameter x Height (mm)	Grad. Range (mL)	Grad. Interval (mL)	Qty/ Cs
16157-071	1395-25	25	36.5 x 70	10-25	5	10
16157-068	1395-50	50	46 x 88	20-50	10	10
16157-103	1395-100	100	56 x 100	40-80	20	10
89427-656	1395-150	150	50 x 125	62-110	25	10
16157-136	1395-250	250	70 x 138	50-200	50	10
16157-169	1395-500	500	86 x 176	100-400	100	10
89427-658	1395-750	750	95 x 203	100-600	100	10
16157-191	1395-1L	1000	101 x 225	100-900	100	10
16157-227	1395-2L	2000	136 x 262	400-1800	200	10
89427-655	1395-3X	3500	160 x 295	500-3000	500	1
16157-250	1395-5L	5000	186 x 335	500-4500	500	1
16157-282	1395-10L	10000	234 x 410	2000-9000	1000	1

1396 PYREX Square Laboratory Storage Bottles

VWR Cat. No.	Corning Cat. No.	Capacity (mL)	Thread Size	Approx. Diameter x Height (mm)	Grad. Range (mL)	Grad. Interval (mL)	Qty/ Cs
11310-680	1396-100	100	GL32	50 x 105	40 - 80	20	10
11310-682	1396-250	250	GL45	64 x 138	50 - 200	50	10
11310-684	1396-500	500	GL45	78 x 176	100 - 400	100	10
11310-686	1396-1L	1000	GL45	94 x 217	100 - 900	100	10





1397 PYREX® GLS 80 Wide Mouth Bottles

VWR Cat. No.	Corning Cat. No.	Capacity (mL)	Thread Size	Approx. Diameter x Height (mm)	Grad. Range (mL)	Grad. Interval (mL)	Qty/Cs
89090-062	1397-500	500	GLS80	101 x 148	100 - 500	100	10
89090-054	1397-1L	1000	GLS80	101 x 218	100 - 1000	100	10
89090-060	1397-2L	2000	GLS80	136 x 248	400 - 2000	200	10

1399 PYREX Slim Line Storage Bottles

VWR Cat. No.	Corning Cat. No.	Capacity (mL)	Thread Size	Approx. Diameter x Height (mm)	Grad. Range (mL)	Grad. Interval (mL)	Qty/Cs
10127-908	1399-125	125	GL45	55 x 124	25 - 125	25	4
10127-910	1399-250	250	GL45	66 x 158	50 - 250	50	4
10127-912	1399-500	500	GL45	78 x 193	100 - 500	100	4
10127-914	1399-1L	1000	GL45	93 x 253	200 - 1000	100	4



51395 PYREX Low Actinic Bottles

VWR Cat. No.	Corning Cat. No.	Capacity (mL)	Thread Size	Approx. Diameter x Height (mm)	Grad. Range (mL)	Grad. Interval (mL)	Qty/Cs
33500-510	51395-25	25	GL25	36.5 x 70	10 - 25	5	4
33500-512	51395-50	50	GL32	46 x 88	20 - 50	10	4
16157-067	51395-100	100	GL45	56 x 100	40 - 80	20	4
16157-065	51395-250	250	GL45	70 x 138	50 - 200	50	4
16157-063	51395-500	500	GL45	86 x 176	100 - 400	100	4
16157-061	51395-1L	1000	GL45	101 x 225	100 - 900	100	4
89090-098	51395-2L	2000	GL45	136 x 262	400 - 1800	200	4
89090-104	51395-5L	5000	GL45	186 x 335	500 - 4500	500	1
89090-920	51395-10L	10000	GL45	234 x 410	2000 - 9000	1000	1



61626 PYREXPLUS® PVC-coated Storage Bottles

VWR Cat. No.	Corning Cat. No.	Capacity (mL)	Thread Size	Approx. Diameter x Height (mm)	Grad. Range (mL)	Grad. Interval (mL)	Qty/Cs
16157-057	61626-250	250	GL45	70 x 138	50 - 200	50	4
16157-055	61626-500	500	GL45	86 x 176	100 - 400	100	4
16157-053	61626-1L	1000	GL45	101 x 225	100 - 900	100	4
16157-051	61626-2L	2000	GL45	136 x 262	400 - 1800	200	4



Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.



vwr.com

Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy and additional disclaimers. © 2026 Avantor Inc. All rights reserved.

The information contained within is accurate as of the date of publication and subject to change without notice.
For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.
© 2014, 2026 Corning Incorporated. All rights reserved. 3/26 CLS-GL-041 REV1 VWR

0326 Lit. No. 030524WREV