

Certificate of Conformity

VWR Catalogue Number	See List
Description	Petri dishes, glass, with lid, VWR®, Non-Sterile
Country of Origin	Manufactured in Czech Republic
Date of Issue (yyyy-mm-dd)	2022-06-20

Catalogue Number	Diameter [mm]	Height [mm]
391-0650	40	12
391-0651	60	15
391-0652	70	15
391-0653	80	15
391-0654	90	15
391-0655	100	15
391-0656	100	20
391-0657	120	20
391-0658	150	25
391-0659	180	30
391-0660	200	30
391-0661	200	50
391-0665	40	12
391-0666	60	15
391-0667	80	15
391-0668	90	15
391-0669	100	15
391-0674	150	25
391-0671	100	20
391-0673	120	20

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Soda-lime glass.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.

Storage Conditions Room Temperature, dry place

Shelf Life No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.