



Quality Certificate

VWR* DISPOSABLE SPATULAS

This document certifies that this product meets the following criteria:

North American Catalog No:	80081-188
Batch Number:	WE20140401
Description:	VWR SPATULA PP DISP BLUE PK300

This document certifies that the above product is manufactured in a facility that also manufactures plastic products for food contact applications.

The manufacturing process includes the following procedures: high heat extrusion, a proprietary method for cutting both ends of the spatula to create the spoon/scrapper end and the narrow scoop end, and, lastly, treatment in a "blue light" room prior to packaging. All employees handling the product wear lab coats, gloves, masks and hairnets. However, the above referenced product is NOT sold as sterile.

The specification for the plastic used in this product is: 100% Virgin Food Grade Polypropylene.

No raw materials that contain, or are derived from, animals are specified in the manufacture of this product.

The manufacturing process meets FDA standards for plastics coming in contact with food.

Quality control for the product is performed upon receipt of each batch. This product is tested for adherence to size, weight and other manufacturing specifications.

The product referenced in this certificate, conforms to the manufacturer specifications.

Signed:

A handwritten signature in blue ink, appearing to read "Ken Crossley", is written over a light blue horizontal line.

Ken Crossley
Manager
Quality Assurance

Date: April 15, 2014