

Certificate of Conformity

VWR Catalogue Number	613-1062
Description	Pipette Tip, Reload System, Graduated, 100 – 1250 μ L Volume, Non-Sterile
Lot Number	21351632
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2021-12-02

Quality System Compliance

Products are manufactured under the **ISO 9001 & ISO 13485** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

PCR Inhibitor: Products tested to be considered free of PCR inhibitors.

Product Specifications

Material	High Grade Polypropylene.
Sterilization	Non-Sterile.
ATP Assay	This product has been tested and is ATP free. The assay detection limit is $<10^{-13}$ mg/ μ L.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination. The assay detection limit is $<10^{-7}$ Units/ μ L for DNase and $<10^{-9}$ Units/ μ L for RNase.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). Tallow derived additives may be used in the manufacturing of the raw material. The material used in the product comply with the Regulations 1069/2009/EC, and 142/2011/EC and the "Note for guidance EMEA 410/01, rev.3".
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	The acceptance level for product is <0.06 EU/mL.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in a room temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature.