



## Exterior


Front View



Side View



Top View


115VAC/60Hz/1PH  
NEMA-5-15R

Cat# 76537-060/VWR-HC-RFC20A

| Compartment  | Cu.Ft. | Defrost | Door    | Int Door | Shelves       | W" | D"     | H" | Refrigerant | Amps | Weight   |
|--------------|--------|---------|---------|----------|---------------|----|--------|----|-------------|------|----------|
| Freezer      | 5.15   | Auto    | 1 Solid | 0        | 1 Fixed       | 30 | 34 1/8 | 69 | R600a       | 6    | 171 lbs. |
| Refrigerator | 15.35  | Cycle   | 1 Solid | 0        | 2 Adj/1 Fixed |    |        |    |             |      |          |