

## Certificate of Conformity Finntip™ Filter Sterile

Product Code: 94060313  
Product Name: Finntip Filter 200ul, sterile  
Lot Number: 21201A0  
Exp: 2026-07

We hereby confirm that these high quality Finntip™ Pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of the following contaminants:

|           |                               |
|-----------|-------------------------------|
| Human DNA | < $1.9 \times 10^{-11}$ g/tip |
| DNase     | < $9.4 \times 10^{-4}$ U/tip  |
| RNase     | < $6.5 \times 10^{-8}$ U/tip  |
| Endotoxin | < 0.05 EU/ml                  |

The sterility (SAL  $10^{-6}$ ) is ensured by irradiating this product according to ISO 11137 guidelines.

In Joensuu, 16th September 2021



Heidi Salokangas, Quality Manager  
Laboratory Plastics Essentials *Kaisu Rannamäki*  
Thermo Fisher Scientific Oy